User Guide: Liquids NMR

Varian NMR Spectrometer Systems With VNMR 6.1C Software Pub. No. 01-999161-00, Rev. C1002



User Guide: Liquids NMR Varian NMR Spectrometer Systems With VNMR 6.1C Software Pub. No. 01-999161-00, Rev. C1002

Revision history:

A0800 – Initial release for VNMR 6.1C software. B0601 – Updated Chapter 13 for LC-NMR2000 and STAR 5.5 B0801 – Added setLP1 macro to chapter 3. B0202 – Corrected information in Table 6 in chapter 1. C1002 – Chapters 13 and 14 moved to LC-NMR Operation and Installation, updated Pbox section.

Applicability of manual:

UNITY INOVA, MERCURY VxWorks Powered NMR spectrometer systems (shortened to MERCURY-VX throughout this manual), MERCURY, UNITYplus, GEMINI 2000, UNITY, and VXR-S NMR superconducting spectrometer systems with VNMR 6.1C software installed.

Technical contributors: Greg Brissey, Steve Cheatham, Bayard Fetler, Phil Hornung, Dan Iverson, Boban John, Frits Vosman, Evan Williams Technical writers: Michael Carlisle, Everett Schreiber Technical editor: Dan Steele

Copyright ©2002 by Varian, Inc. 3120 Hansen Way, Palo Alto, California 94304 http://www.varianinc.com All rights reserved. Printed in the United States.

The information in this document has been carefully checked and is believed to be entirely reliable. However, no responsibility is assumed for inaccuracies. Statements in this document are not intended to create any warranty, expressed or implied. Specifications and performance characteristics of the software described in this manual may be changed at any time without notice. Varian reserves the right to make changes in any products herein to improve reliability, function, or design. Varian does not assume any liability arising out of the application or use of any product or circuit described herein; neither does it convey any license under its patent rights nor the rights of others. Inclusion in this document does not imply that any particular feature is standard on the instrument.

UNITY INOVA, MERCURY, Gemini, GEMINI 2000, UNITYplus, UNITY, VXR, XL, VNMR, VnmrS, VnmrX, VnmrI, VnmrV, VnmrSGI, MAGICAL II, AutoLock, AutoShim, AutoPhase, limNET, ASM, and SMS are registered trademarks or trademarks of Varian, Inc. Sun, Solaris, CDE, Suninstall, Ultra, SPARC, SPARCstation, SunCD, and NFS are registered trademarks or trademarks of Sun Microsystems, Inc. and SPARC International. Oxford is a registered trademark of Oxford Instruments LTD. Ethernet is a registered trademark of Xerox Corporation. VxWORKS and VxWORKS POWERED are registered trademarks of WindRiver Inc. Dell is a registered trademark of Dell Computer Corporation. Windows is either a registered trademark of Microsoft Corporation in the United States and or other countries.Other product names in this document are registered trademarks of their respective holders.

Overview of Contents

SAFETY PRECAUTIONS	23
Introduction	27
Chapter 1. Advanced 1D NMR	29
Chapter 2. 1D Experiments	71
Chapter 3. Multidimensional NMR	87
Chapter 4. Multidimensional and Advanced Experiments	135
Chapter 5. Indirect Detection Experiments	169
Chapter 6. Data Analysis	199
Chapter 7. Pulse Analysis	223
Chapter 8. Variable Temperature Operation	261
Chapter 9. Carousel, SMS, and NMS Automation	269
Chapter 10. VAST Accessory Operation	311
Chapter 11. PFG Modules Operation	371
Chapter 12. PFG Modules Experiments	393
Index	401

Table of Contents

SAFETY P	RECAUTIONS	23
Introductio	on	27
Chapter 1.	Advanced 1D NMR	29
-	rking with Experiments	
	Iti-FID (Arrayed) Spectra	
	Arrayed Parameters	30
	Multiple Arrays	31
	Setting Array Order and Precedence	32
	Interactively Arraying Parameters	32
	Resetting an Array	32
	Array Limitations	33
	Acquiring Data	33
	Processing	33
	Display and Plotting	33
	Saving and Retrieving	34
	Pulse Width Calibration Step-by-Step	35
$1.3 T_1$	and <i>T</i> ₂ Analysis	35
	Setting Up The Experiment	36
	Processing the Data	36
	Analyzing the Data	36
	Exponential Analysis Menu	37
	T ₁ Data Workup: Step-by-Step	37
1.4 Kir	etics	37
	Setting Up the Experiment	37
	Processing the Data	38
	Kinetics Step-by-Step	38
1.5 Dif	fusion Experiments/DOSY	38
	Pulsed Gradient Experiments	38
	Pulsed Gradient Experiment Setup	40
	Gradient Calibration	42
	Data Reduction	42
	Data Display	45
	Variations on the pge Pulse Sequence	46
	DOSY Experiments	46
	Filter Diagonalization Method	66
	Using FDM	66
Chapter 2.	1D Experiments	71
-	T—Attached Proton Test	72
	Applicability	72
	Parameters	72
	Technique	72
	References	72
	Related Commands and Macros	73
2.2 BIN	NOM—Binomial Water Suppression	73

Applicability	73
Parameters	73
Reference	73
2.3 CPMGT2—Carr-Purcell Meiboom-Gill <i>T</i> ₂ Measurement	74
Applicability	74
Parameters	74
T ₂ Measurement	74
Acquisition and Processing	74
2.4 CYCLENOE—Cycled NOE Difference Experiment	75
Applicability	75
Parameters	75
Technique	76
2.5 D2PUL—Standard Two-Pulse Using Decoupler as Transmitter	76
Applicability	76
Parameters	77
Technique	77
2.6 DEPT—Distortionless Enhancement by Polarization Transfer	77
•	77
Applicability	
Parameters	78
Technique	78
Potential Problems	79
Reference	79
Related Commands and Macros	79
2.7 INEPT—Insensitive Nuclei Enhanced by Polarization Transfer	80
Applicability	80
Parameters	80
Technique	81
Reference	81
2.8 JUMPRET—Jump-and-Return Water Suppression	81
Applicability	81
Parameters	81
Reference	81
2.9 NOEDIF—NOE Difference Experiment	82
Applicability	82
Parameters	82
Phase Cycling	82
Procedure	83
Reference	85
	85
2.10 PRESAT—1D Water Suppression	85 85
Applicability	
Parameters	85
2.11 S2PUL—Standard Two-Pulse Sequence	85
Applicability	85
Parameters	85
2.12 S2PULR—Standard Two-Pulse in Reverse Configuration	86
Applicability	86
Parameters	86
Technique	86

Chapter 3. Multidimensional NMR	. 87
3.1 Interferograms	87
3.2 2D NMR Step-by-Step	88
To Process Stored Data	88
To Acquire a Simple COSY	89
3.3 Phase-Sensitive 2D NMR	89
3.4 Data Acquisition: Arrayed 2D	
Hypercomplex Method	90
TPPI Method	91
Real-Time 2D	92
Macros for 2D Experiments	92
3.5 Weighting	93
Parameters	93
Setting Values	93
Interactive Weighting	94
3.6 Phasing Before the 2D Transform	95
3.7 Baseline Correction	
First-Point Multiplier	97
Baseline Correction	98
FID Drift Correction	98
Spectral Drift Correction	99
3.8 Processing Phase-Sensitive 2D and 3D Data	
Processing Programs	101
Common Coefficients for wft2d Processing	101
Sign of f ₁ Frequencies	102
2D Solvent Subtraction Filtering	103
Left Shift, Frequency Shift, Phase Rotation	104
2D Processing of 3D Data	104
3.9 2D and 3D Linear Prediction	104
3.10 Phasing the 2D Spectrum	106
3.11 Display and Plotting	107
Display Modes	107
Display and Plot Limits	107
Maximum Intensity	109
Axis Label and Direction	109
Display Scaling	110
Grid Lines	110
Color Maps and Contour Plots	110
Whitewashed Spectra	111
Label Display	111
Projection of 2D Data	111
2D Referencing	112
Rotating Homonuclear 2D-J Spectra	112
Symmetrizing Data	112
Setting Negative Intensities to Zero	112
Automatic Analysis	113
3.12 Interactive 2D Color Map Display	113
Interactive 2D Display Menus	115
Controlling the Display with the Mouse	116

Changing the Display	116
Treating 2D Traces as 1D Spectra	118
3.13 Interactive 2D Peak Picking	118
Interactive 2D Peak Picking Menus	121
Automatic 2D Peak Picking	125
Interactive Peak Picking or Editing	125
Automatic Integration	126
Interactive Integration and Editing	126
Labeling and Commenting Peaks	126
Displaying Peaks in dconi	127
Peak File Manipulations	127
3.14 3D NMR	127
3D Acquisition	128
3D Processing	129
3D Display	130
3D Pulse Sequences	130
Experiment Setup	130
Data Processing	131
3.15 4D NMR Acquisition	132
Chapter 4. Multidimensional and Advanced Experiments	
4.1 Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY	
Applicability	137
Parameters	137
Technique	137
Potential Problems	139
4.2 COSY—Correlated Spectroscopy	
Applicability	139
Parameters	139
Technique	139
4.3 COSYPS—Phase-Sensitive COSY	139
Applicability	139
Parameters	139
Phase Cycling	140
Technique	141
References	141
4.4 DQCOSY—Double-Quantum Filtered COSY	141
Applicability	141
Parameters	141
Phase Cycling	142
Technique	142
Potential Problems	143
References	143
4.5 HET2DJ—Heteronuclear 2D-J	143
Applicability	143
Parameters	143
Technique	144
Potential Problems	144
	177

	145
Applicability	145
Parameters	145
1	146
	147
References	147
4.7 HETCORPS—Absolute-Value and Phase-Sensitive HETCOR	147
Applicability	148
Parameters	148
Recommendations	148
4.8 HOM2DJ—Homonuclear J-resolved 2D	148
Applicability	148
Parameters	149
Technique	149
Potential Problems	150
References	150
4.9 INADEQUATE—Double-Quantum Transfer Experiment	150
	150
	150
	151
-	151
	151
	151
	152
	152
	153
	153
	153
	154
	154
*	155
	156
	156
	156
	156
	157
	157
	158
	158
	158
	158
	158
	158
	158
	158 159
	159
	159
	159
4.16 TNNOESY—NOESY with Water Suppression	159

Applicability	160
Parameters	160
4.17 TNROESY—ROESY with Water Suppression	160
Applicability	160
Parameters	
Technique	161
4.18 TNTOCSY—TOCSY with Water Suppression	
Applicability	
Parameters	
References	162
4.19 TOCSY—Total Correlation Spectroscopy	
Applicability	
Parameters	
Phase Cycling	
Technique	
References	
Related Macros	
4.20 TROESY—Transverse ROESY	
Applicability	
Parameters	
Reference	
4.21 HCCHTOCSY Pulse Sequence	
Applicability	
Parameters	
Technique	
4.22 HMQCTOCSY Pulse Sequence	
Parameters	
Technique	
-	
4.23 HMQC-TOCSY 3D Pulse Sequence	
Applicability Parameters	
Technique	
4.24 HSQC-TOCSY 3D Pulse Sequence	
Applicability	
Parameters	168
Chapter 5. Indirect Detection Experiments	169
5.1 Requirements for Indirect Detection Experiments	169
Probes	
RF System	
Pulse Sequences	
HMQC Pulse Sequence	
HMQCR Pulse Sequence	
5.2 The Basic HMQC Experiment	
Spin-Echo Difference Experiment	
BIRD Nulling	
Transmitter Presaturation for High-Dynamic Range Signals	
5.3 Phase-Sensitive Aspects of the Sequence	

5.4 Cancellation Efficiency	
5.5 Pros and Cons of Decoupling	
5.6 Specifications Testing	
5.7 Using the HMQC and HMQCR Sequences	
5.8 Recabling Single-Broadband Systems	
5.9 Recabling Dual-Broadband Systems	
5.10 Filters for Indirect Detection	
5.11 Tuning the Probe in the Reverse Mode	
Tune the ¹ H Channel	
Tune the X Channel	
5.12 Controlling Transmitter Power in the Reverse Mode	
5.13 Indirect Detection Calibration	
5.14 Typical Experimental Protocol for HMQC Experiments	
5.15 Differences for ¹⁵ N Indirect Detection	
5.16 HSQC Experiment	
Applicability	
Parameters	
hapter 6. Data Analysis	19
6.1 Spin Simulation	19
Spin Simulation Step-by-Step	
Spin Simulation Menus	
Entering a Spin System	
Spin Simulation Parameters	
Performing a Spin Simulation	
Iterative Mode	
Spin Simulation Files	
6.2 Deconvolution	
Deconvolution Step-by-Step	
Performing Deconvolution	
Display and Plotting	
Deconvolution Menu	
6.3 Reference Deconvolution	
Reference Deconvolution of 1D Spectra	
Reference Deconvolution of 2D Spectra	
References	
6.4 Addition and Subtraction of Data	
Add/Subtract Menu	
Noninteractive Add/Subtract	
Interactive Add/Subtract	
6.5 Regression Analysis	
Regression Commands and Menus	
Regression Analysis Step-by-Step	
Contents of "analyze.out" File	
Contents of "regression.inp" File	
6.6 Chemical Shift Analysis	
hapter 7. Pulse Analysis	
7.1 Pulse Shape Analysis	

Directory and File Operations	
Attribute Selection	. 225
Scale and Reference	. 225
Cursors	. 225
Simulation Overview	. 226
Simulation Parameters	. 226
Performing a Simulation	. 227
Creating a Pulse	. 228
7.2 Pandora's Box	228
Getting Started	. 229
Calibrating the RF Field	. 229
Creating Waveforms from Macros	
Creating Waveforms from UNIX	
Pbox File System	
Pbox Equation Evaluator	
Reserved Shape Names	
Pbox VNMR Parameters	
Wave String Variables	
Creating Waveforms Using Menus	
Pbox Macro Reference	
Pbox PSG Statements	
Pulse Shaping "On-Fly"	
Pbox_psg.h include Pulse Sequence Statements	
shonfly.c Sequence	
Pbox UNIX Commands	
Chapter 8. Variable Temperature Operation	. 261
8.1 Startup	. 261
8.2 Operating Procedures	262
8.3 Temperature-Related Command	264
8.4 Operating Recommendations	265
8.5 VT Controller Safety Circuits	265
8.6 VT Interlock Parameters	
0.0 V T Interfock T draineers	266
	266
Chapter 9. Carousel, SMS, and NMS Automation	266 267
Chapter 9. Carousel, SMS, and NMS Automation 9.1 Carousel Autosampler	266 267 . . 269
•	266 267 . . 269 269
9.1 Carousel Autosampler Configuring VNMR for the Carousel	266 267 269 269 . 270
9.1 Carousel Autosampler	266 267 269 269 . 270 . 271
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel	266 267 269 269 . 270 . 271 . 273
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel Adjusting the Eject Air	266 267 269 269 . 270 . 271 . 273 . 274
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel Adjusting the Eject Air Loading and Unloading Samples	266 267 269 269 . 270 . 271 . 273 . 274 . 275
 9.1 Carousel Autosampler	266 267 269 269 270 271 273 274 275 276
 9.1 Carousel Autosampler	266 267 269 269 269 270 271 273 274 275 276 277
 9.1 Carousel Autosampler	266 267 269 269 270 271 273 274 275 276 277 278
 9.1 Carousel Autosampler	266 267 269 269 270 271 273 274 275 276 277 278 279
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel Adjusting the Eject Air Loading and Unloading Samples Running NMR on One Sample at a Time Running Automated NMR on Up to Nine Samples Inserting Samples Manually with the Carousel Attached Carousel Error Codes and Recovery 9.2 SMS Autosampler	266 267 269 269 270 271 273 274 275 276 277 278 280
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel Adjusting the Eject Air Loading and Unloading Samples Running NMR on One Sample at a Time Running Automated NMR on Up to Nine Samples Inserting Samples Manually with the Carousel Attached Carousel Error Codes and Recovery 9.2 SMS Autosampler Configuring VNMR for the SMS Autosampler	266 267 269 269 270 271 273 274 275 276 277 278 280 281
 9.1 Carousel Autosampler	266 267 269 269 270 271 273 274 275 276 277 278 278 280 281 282
9.1 Carousel Autosampler Configuring VNMR for the Carousel Checking Out the Carousel Mounting and Removing the Carousel Adjusting the Eject Air Loading and Unloading Samples Running NMR on One Sample at a Time Running Automated NMR on Up to Nine Samples Inserting Samples Manually with the Carousel Attached Carousel Error Codes and Recovery 9.2 SMS Autosampler Configuring VNMR for the SMS Autosampler	266 267 269 269 270 271 273 274 275 275 276 277 278 279 280 281 282 282 282

SMS Error Codes and Recovery	284
9.3 NMS Autosampler	285
Before Using NMS	286
Configuring VNMR for the NMS Autosampler	286
Running NMR on One Sample at a Time	286
Running Automated NMR	287
9.4 General Automation Tasks For All Sample Changers	289
Preparing and Initiating an Automation Run	289
Setting Up an Automation Run for Multiple Users	290
Monitoring an Automation Run	291
Using Sample Changers in Continuous Walkup Mode	293
Adding Samples to an Automation Run in Progress	294
9.5 Changing Sample Changers or Serial Ports	295
9.6 Using Gradient Autoshimming with Automation	295
9.7 Automation Run Description	295
Basic Automation Run	296
Automation Behind the Scenes	297
While an Automation Run is in Progress	298
When an Automation Run is Finished	299
Parameters for Automation	299
Variable Temperature Control During Automation	299
9.8 Customizing the Sample Entry Window	300
9.9 Automated Data Acquisition	301
Optimizing Acquisition Macros	302
Customizing Macro Operation	302
Example of Customizing a Macro	304
9.10 Automated Data Processing	308
9.11 File Structures in an Automation Run	308
Chapter 10. VAST Accessory Operation	311
10.1 Using the VAST Accessory	311
To Prepare VAST for Use	312
To Set Up NMR Experiments for VAST	314
To Change Samples with VNMR	
To Shut Down a VAST System	317
10.2 Solvent Suppression in VAST	317
Setting Up Solvent Suppression	317
Troubleshooting Solvent Suppression	322
Evaluating Solvent Mixture Equilibration	322
Solvent Suppression: Background Information	322
10.3 Processing, Displaying, and Plotting VAST Data Sets	323
Creating a Pseudo 2D Data Set	323
Processing, Displaying, and Plotting Glued VAST Data	326
Defining a Custom Display Order with plate_glue	327
Examples of Plots of a VAST Data Set	328
Summary of VAST Display and Plot Options	329
10.4 Using CombiPlate to Analyze Data	331
Preparing VNMR Data For Analysis Using CombiPlate	331
Data Analysis Using CombiPlate And VNMR	332

Analyzing Data Using CombiPlate Without VNMR.	334
Checking And Fixing The Color Map	336
10.5 Vast Process, Display, and Plot Macros	. 337
10.6 Preparing the Hardware and Configuring VNMR	341
Connecting the Transfer Tube	341
Connecting the Air Tubing	342
Connecting Signal and Power Cables	342
Configuring VNMR for VAST	344
10.7 Calibrating Volumes and Flow Rates	345
To Calibrate Probe Volume	345
To Calibrate Sample Volume	350
To Calibrate Flow Rate Parameters	351
To Calibrate XYZ Positions of the Arm	354
10.8 Acquiring Data on Standard Test Samples	354
10.9 Evaluating Carryover	355
10.10 VAST Interface Description	355
SAMPLE Def	356
Rack Def. Pane	360
Main Control	360
Calibrations	361
10.11 Customizing the enter Window for VAST	362
10.12 Files that Control VAST Operation	
10.13 Writing VAST Protocols	
·	
Chapter 11. PFG Modules Operation	. 371
-	
11.1 Configuring the Software	371
-	371
11.1 Configuring the Software11.2 PFG Amplifier Operation11.3 Shimming PFG Systems	371 372 374
 11.1 Configuring the Software	371 372 374 374
11.1 Configuring the Software11.2 PFG Amplifier Operation11.3 Shimming PFG Systems	371 372 374 374
 11.1 Configuring the Software	371 372 374 374 374
 11.1 Configuring the Software	371 372 374 374 374
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 374
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 374
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 375 375 376 377
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375 376 377 377 377
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375 376 377 377 377
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375 376 377 377 377
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 375 375 376 377 377 378
 11.1 Configuring the Software 11.2 PFG Amplifier Operation 11.3 Shimming PFG Systems Performa I and Performa II Performa XYZ 11.4 Setting Up Software for Imaging Pulse Sequences Calibrating the Gradients Creating a Gradient Table Setting the System Gradient Coil 11.5 Homospoil Gradient Type 11.6 Gradient Shimming Configuring Gradients and Hardware Control Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming System Menu General User Gradient Shimming How Gradient Shimming Works 	 371 372 374 374 374 374 374 374 374 374 375 375 375 376 377 377 378 379 379 379 379
 11.1 Configuring the Software	 371 372 374 374 374 374 374 374 374 374 375 375 375 376 377 377 378 379 379 379 379
 11.1 Configuring the Software 11.2 PFG Amplifier Operation 11.3 Shimming PFG Systems Performa I and Performa II Performa XYZ 11.4 Setting Up Software for Imaging Pulse Sequences Calibrating the Gradients Creating a Gradient Table Setting the System Gradient Coil 11.5 Homospoil Gradient Type 11.6 Gradient Shimming Configuring Gradients and Hardware Control Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming System Menu General User Gradient Shimming How Gradient Shimming Works References How Making a Shimmap Works 	 371 372 374 374 374 374 374 374 374 374 375 375 375 376 377 377 378 379 379 379 379
 11.1 Configuring the Software 11.2 PFG Amplifier Operation 11.3 Shimming PFG Systems Performa I and Performa II Performa XYZ 11.4 Setting Up Software for Imaging Pulse Sequences Calibrating the Gradients Creating a Gradient Table Setting the System Gradient Coil 11.5 Homospoil Gradient Type 11.6 Gradient Shimming Configuring Gradients and Hardware Control Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming Quitting the Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming Method Mapping the Shimming Quitting the Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming Method Mapping the Shimming Mow Gradient Shimming Works References 	 371 372 374 374 374 374 374 374 374 375 375 376 377 377 378 379 380
 11.1 Configuring the Software 11.2 PFG Amplifier Operation 11.3 Shimming PFG Systems Performa I and Performa II Performa XYZ 11.4 Setting Up Software for Imaging Pulse Sequences Calibrating the Gradients Creating a Gradient Table Setting the System Gradient Coil 11.5 Homospoil Gradient Type 11.6 Gradient Shimming Configuring Gradients and Hardware Control Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming Quitting the Gradient Shimming Quitting the Gradient Shimming How Gradient Shimming Works References How Making a Shimmap Works Shimmap Files and Parameters How Automated Shimming Works 	 371 372 374 374 374 374 374 374 374 374 375 375 376 377 377 378 378 379 380 381
 11.1 Configuring the Software 11.2 PFG Amplifier Operation 11.3 Shimming PFG Systems Performa I and Performa II Performa XYZ 11.4 Setting Up Software for Imaging Pulse Sequences Calibrating the Gradients Creating a Gradient Table Setting the System Gradient Coil 11.5 Homospoil Gradient Type 11.6 Gradient Shimming Configuring Gradients and Hardware Control Gradient Shimming Method Mapping the Shims Starting Gradient Shimming Quitting the Gradient Shimming Quitting the Gradient Shimming How Gradient Shimming Works References How Making a Shimmap Works Shimmap Files and Parameters 	 371 372 374 374 374 374 374 374 374 375 375 375 376 377 377 377 378 379 380 381 382 383 383

Full Deuterium Gradient Shimming Procedure for Lineshape	386
Setting Up Automation	386
Suggestions for Improving Results	387
Gradient Shimming Menus	388
Chapter 12. PFG Modules Experiments	393
12.1 GCOSY—PFG Absolute-Value COSY	393
Parameters	393
Processing	394
12.2 GHMQC—PFG HMQC	394
Parameters	394
Processing	394
12.3 GHMQCPS—PFG HMQC, Phase Sensitive	394
Processing	395
Recommendations	395
12.4 GHSQC—PFG HSQC, Absolute Value or Phase Sensitive	395
Parameters	395
Processing	396
12.5 GMQCOSY—PFG Absolute-Value MQF COSY	396
Parameters	396
Processing	396
12.6 GNOESY—PFG NOESY	397
Parameters	397
Processing	397
12.7 GTNNOESY—PFG TNNOESY	398
Parameters	398
Processing	398
12.8 GTNROESY—PFG Absolute-Value ROESY	398
Parameters	398
Processing	400
12.9 PFG Selective Excitation	400
Parameters	400
Reference	400
Index	401

Table of Contents

List of Figures

Figure 1. PGE Pulse Sequences	. 39
Figure 2. Data Processing Macros Flowchart	. 43
Figure 3. Sample analyze.inp File	
Figure 4. Gradient Stimulated Echo Element	. 47
Figure 5. Tcl-Tk Acquisition Panel of Doneshot Pulse Sequence	. 50
Figure 6. Dbppste Experiment	. 50
Figure 7. DgcsteSL Experiment	. 51
Figure 8. Oneshot DOSY Experiment	. 52
Figure 9. Tcl-Tk Process Panel for 2D_DOSY Pulse Sequences	. 55
Figure 10. Dgcstecosy (AV Mode) Experiment	. 57
Figure 11. Dgcstehmqc Experiment (AV Mode)	. 57
Figure 12. Tcl-Tk Process2 Panel for the 3D-DOSY Pulse Sequences	
Figure 13. fdm1.inparm File	. 69
Figure 14. APT Pulse Sequence	. 72
Figure 15. BINOM Pulse Sequence	. 73
Figure 16. CPMGT2 Pulse Sequence	. 74
Figure 17. D2PUL Pulse Sequence	. 76
Figure 18. DEPT Pulse Sequence	. 77
Figure 19. INEPT Pulse Sequence	. 80
Figure 20. NOEDIF Pulse Sequence	. 82
Figure 21. S2PULR Pulse Sequence	. 86
Figure 22. Data Flow in Phase-Sensitive 2D Transformation	. 99
Figure 23. Interactive 2D Contour Display (dconi Program)	114
Figure 24. Interactive 2D Peak Picking (ll2d Program)	121
Figure 25. Absolute-Value COSY Pulse Sequence	136
Figure 26. RELAY- COSY Pulse Sequences	
Figure 27. COSYPS Pulse Sequence	140
Figure 28. DQCOSY Pulse Sequence	141
Figure 29. HET2DJ Pulse Sequence	143
Figure 30. HETCOR Pulse Sequence	145
Figure 31. HOM2DJ Pulse Sequence	149
Figure 32. INADEQUATE Pulse Sequence	150
Figure 33. MQCOSY Pulse Sequence	152
Figure 34. NOESY Pulse Sequence	153
Figure 35. ROESY Pulse Sequence	156
Figure 36. TOCSY Pulse Sequence	162
Figure 37. HMQC-TOCSY 3D Pulse Sequence	
Figure 38. HMQC Pulse Sequence with null<>0 and mbond='n'	
Figure 39. HMQCR Pulse Sequence	
Figure 40. Heteronuclear Spin-Echo Difference Experiment	176

Figure 41. HMQC Pulse Sequence, Showing Movement of Attached Protons	176
Figure 42. Evolution Time Added Between X-Nucleus Pulses	177
Figure 43. HMQC with BIRD Pulse Nulling Effects	178
Figure 44. Basic HMQC Pulse Sequence	182
Figure 45. Normal ¹³ C Spectrum of ¹³ CH ₃ I	186
Figure 46. Normal ¹ H Spectrum of ¹³ CH ₃ I	186
Figure 47. HMQC Without and With X-Nucleus Pulses	187
Figure 48. Calibration of pwx, Coarse and Fine	188
Figure 49. HMQC 4-Pulse Cancellation without ² H Bandpass Filter	189
Figure 50. HMQC 4-Pulse Cancellation with ² H Bandpass Filter	189
Figure 51. HMQC Without and With ¹³ C Decoupling	190
Figure 52. Coupled HMQC Spectrum of ¹³ CH ₃ I	191
Figure 53. Verifying Cancellation with pwx=0,90	193
Figure 54. Optimizing the BIRD Nulling Time	194
Figure 55. Coupled HMQC Spectrum of 3-Heptanone	195
Figure 56. Expansion of Coupled 3-Heptanone HMQC Showing Multiplets	196
Figure 57. Decoupled HMQC Spectrum of 3-Heptanone	196
Figure 58. Display of Regression Fittings (expl Program)	218
Figure 59. Temperature Control Window	265
Figure 60. Carousel Autosampler Carousel and Driver	270
Figure 61. Carousel with Sensor Arm in Locked and in Run Positions	271
Figure 62. Lock Pin on Sensor Arm Disengaged and Engaged	271
Figure 63. Bottom View of Carousel, Showing Connector and Groove	272
Figure 64. Optical Sensor and Proper Sample Floating Height	274
Figure 65. Sample Entry Form Window for the Carousel Autosampler	277
Figure 66. Manually Loading and Unloading Samples Through Position 1	279
Figure 67. SMS Autosampler	281
Figure 68. Sample Entry Form Window (100 samples) for the SMS Autosampler	283
Figure 69. NMS Window	286
Figure 70. Sample Entry Form Window for the NMS Autosampler	288
Figure 71. Sample Status Window (status Program)	292
Figure 72. Locate Window (status Program)	292
Figure 73. Sample Entry Form Window for Walkup Operation	293
Figure 74. Gilson 215 Liquid Handler	312
Figure 75. Gilson Liquid Handler Window	312
Figure 76. VAST Sample Entry Form Window (enter program)	315
Figure 77. LC-NMR Pane for VAST	
Figure 78. VAST Sequence Pane	318
Figure 79. VAST Pane	319
Figure 80. VAST Acq & Obs Pane	320
Figure 81. VAST Sequence Pane and Frequency Suppression Options	320
Figure 82. Partial doneQ File	
Figure 83. Queue Name and Scout Directory Fields	325
Figure 84. plate_glue Window and e4x4 glue order.	327
Figure 85. Array of 1D Spectra from a VAST Data Set	
Figure 86. Plot of a Subset from a 96-Well Sample Plate	

Figure 87. Contour Plot of a VAST Data Set	329
Figure 88. Stacked Plot of a VAST Data Set	330
Figure 89. Spectrum of the Sample in Well F6	330
Figure 90. CombiPlate Started by combishow From vnmr	332
Figure 91. CombiPlate Display - Binary Display For Each Region	333
Figure 92. CombiPlate Save File Window	334
Figure 93. CombiPlate Column and Row Input Screen	335
Figure 94. CombiPlate Field Box	335
Figure 95. CombiPlate Window for VAST Data Analysis	336
Figure 96. Connection Between the Rheodyne Injector Valve and Transfer Tube	341
Figure 97. VAST Air Connections	343
Figure 98. VAST Signal Cable Connections	344
Figure 99. Microflow Probe and Transfer Tube Volumes	345
Figure 100. Injecting Solution into a Microflow Probe	346
Figure 101. Microflow Probe Flow Cell	348
Figure 102. Finding the Optimum Sample Volume	352
Figure 103. Sample Definition Pane	356
Figure 104. Rack Definition Window	360
Figure 105. Main Control Window	360
Figure 106. Calibrations Window	362
Figure 107. Gradient Shimming Pulse Sequence	379
Figure 108. Mapping the z1 Shim	380
Figure 109. Shimmap Plot	381
Figure 110. Curve Fit Plot	384
Figure 111. Display of Shim Adjustments for Each Iteration	384

List of Figures

List of Tables

Table 1. Commands and Macros for Working with Experiments	. 29
Table 2. Arrayed Spectra Commands and Parameters	31
Table 3. T_1 and T_2 Analysis Commands and Parameters	. 35
Table 4. Kinetics Analysis Commands and Parameters	. 38
Table 5. Diffusion Analysis Commands	39
Table 6. Gradient Control Parameter Set	40
Table 7. Tools for the DOSY Experiment	. 47
Table 8. Dbppste Parameters	. 51
Table 9. DgcsteSL Parameters	. 51
Table 10. Oneshot DOSY Parameters	. 52
Table 13. Dgcstecosy Parameters	57
Table 14. Dgcstehmqc Parameters	. 58
Table 16. Arrayed 2D & 3D Data Acquisition Commands and Parameters	. 91
Table 17. Macros for 2D Experiments	. 92
Table 18. Weighting Parameters for ni and ni2 Dimensions	. 93
Table 19. Commands and Macros for Setting 2D Weighting Values	. 94
Table 20. Commands and Parameters for Phasing Before the 2D Transform	. 96
Table 21. Baseline and Drift Correction Commands and Parameters	. 97
Table 22. Tools for Processing Phase-Sensitive 2D and 3D Data	100
Table 23. 2D and 3D Linear Prediction (LP) Commands and Parameters	105
Table 24. Commands and Parameters for Phasing the 2D Spectrum	106
Table 25. 2D Display and Plotting Commands and Parameters (Part 1 of 2)	108
Table 26. 2D Display and Plotting Commands and Parameters (Part 2 of 2)	109
Table 27. Interactive 2D Color Map Display Commands and Parameters	114
Table 28. Interactive 2D Peak Picking Commands and Parameters	118
Table 29. 3D NMR Commands and Parameters (Part 1 of 2)	128
Table 30. 3D NMR Commands and Parameters (Part 2 of 2)	129
Table 31. 4D NMR Acquisition Commands and Parameters	133
Table 32. Parameters for HMQC and HMQCR Pulse Sequences	182
Table 33. System Bandpass Filters for Indirect Detection Probes.	184
Table 34. Parameter Values for HMQCR on Natural Abundance Sample	195
Table 35. Spin Simulation Commands and Parameters	200
Table 36. Deconvolution Commands and Parameters	207
Table 37. Reference Deconvolution Commands	209
Table 38. Add/Subtract Experiment Commands and Parameters	212
Table 39. Regression Commands	216
Table 40. Curve Types.	
Table 41. Chemical Shift Analysis Commands	221
Table 42. Pulse Shape Analysis Commands and Parameters	223
Table 43. Pbox Commands and Parameters	259

261
296
300
304
313
349
353
372
376
377
378
385
393
394
396
397
397
398
399
400

SAFETY PRECAUTIONS

The following warning and caution notices illustrate the style used in Varian manuals for safety precaution notices and explain when each type is used:

- **WARNING:** Warnings are used when failure to observe instructions or precautions could result in injury or death to humans or animals, or significant property damage.
- **CAUTION:** Cautions are used when failure to observe instructions could result in serious damage to equipment or loss of data.

Warning Notices

Observe the following precautions during installation, operation, maintenance, and repair of the instrument. Failure to comply with these warnings, or with specific warnings elsewhere in Varian manuals, violates safety standards of design, manufacturing, and intended use of the instrument. Varian assumes no liability for customer failure to comply with these precautions.

WARNING: Persons with implanted or attached medical devices such as pacemakers and prosthetic parts must remain outside the 5-gauss perimeter from the centerline of the magnet.

The superconducting magnet system generates strong magnetic fields that can affect operation of some cardiac pacemakers or harm implanted or attached devices such as prosthetic parts and metal blood vessel clips and clamps.

Pacemaker wearers should consult the user manual provided by the pacemaker manufacturer or contact the pacemaker manufacturer to determine the effect on a specific pacemaker. Pacemaker wearers should also always notify their physician and discuss the health risks of being in proximity to magnetic fields. Wearers of metal prosthetics and implants should contact their physician to determine if a danger exists.

Refer to the manuals supplied with the magnet for the size of a typical 5-gauss stray field. This gauss level should be checked after the magnet is installed.

WARNING: Keep metal objects outside the 10-gauss perimeter from the centerline of the magnet.

The strong magnetic field surrounding the magnet attracts objects containing steel, iron, or other ferromagnetic materials, which includes most ordinary tools, electronic equipment, compressed gas cylinders, steel chairs, and steel carts. Unless restrained, such objects can suddenly fly towards the magnet, causing possible personal injury and extensive damage to the probe, dewar, and superconducting solenoid. The greater the mass of the object, the more the magnet attracts the object.

Only nonferromagnetic materials—plastics, aluminum, wood, nonmagnetic stainless steel, etc.—should be used in the area around the magnet. If an object is stuck to the magnet surface and cannot easily be removed by hand, contact Varian service for assistance.

Warning Notices (continued)

Refer to the manuals supplied with the magnet for the size of a typical 10-gauss stray field. This gauss level should be checked after the magnet is installed.

WARNING: Only qualified maintenance personnel shall remove equipment covers or make internal adjustments.

Dangerous high voltages that can kill or injure exist inside the instrument. Before working inside a cabinet, turn off the main system power switch located on the back of the console, then disconnect the ac power cord.

WARNING: Do not substitute parts or modify the instrument.

Any unauthorized modification could injure personnel or damage equipment and potentially terminate the warranty agreements and/or service contract. Written authorization approved by a Varian, Inc. product manager is required to implement any changes to the hardware of a Varian NMR spectrometer. Maintain safety features by referring system service to a Varian service office.

WARNING: Do not operate in the presence of flammable gases or fumes.

Operation with flammable gases or fumes present creates the risk of injury or death from toxic fumes, explosion, or fire.

WARNING: Leave area immediately in the event of a magnet quench.

If the magnet dewar should quench (sudden appearance of gasses from the top of the dewar), leave the area immediately. Sudden release of helium or nitrogen gases can rapidly displace oxygen in an enclosed space creating a possibility of asphyxiation. Do not return until the oxygen level returns to normal.

WARNING: Avoid liquid helium or nitrogen contact with any part of the body.

In contact with the body, liquid helium and nitrogen can cause an injury similar to a burn. Never place your head over the helium and nitrogen exit tubes on top of the magnet. If liquid helium or nitrogen contacts the body, seek immediate medical attention, especially if the skin is blistered or the eyes are affected.

WARNING: Do not look down the upper barrel.

Unless the probe is removed from the magnet, never look down the upper barrel. You could be injured by the sample tube as it ejects pneumatically from the probe.

WARNING: Do not exceed the boiling or freezing point of a sample during variable temperature experiments.

A sample tube subjected to a change in temperature can build up excessive pressure, which can break the sample tube glass and cause injury by flying glass and toxic materials. To avoid this hazard, establish the freezing and boiling point of a sample before doing a variable temperature experiment.

Warning Notices (continued)

WARNING: Support the magnet and prevent it from tipping over.

The magnet dewar has a high center of gravity and could tip over in an earthquake or after being struck by a large object, injuring personnel and causing sudden, dangerous release of nitrogen and helium gasses from the dewar. Therefore, the magnet must be supported by at least one of two methods: with ropes suspended from the ceiling or with the antivibration legs bolted to the floor. Refer to the *Installation Planning Manual* for details.

WARNING: Do not remove the relief valves on the vent tubes.

The relief valves prevent air from entering the nitrogen and helium vent tubes. Air that enters the magnet contains moisture that can freeze, causing blockage of the vent tubes and possibly extensive damage to the magnet. It could also cause a sudden dangerous release of nitrogen and helium gases from the dewar. Except when transferring nitrogen or helium, be certain that the relief valves are secured on the vent tubes.

WARNING: On magnets with removable quench tubes, keep the tubes in place except during helium servicing.

On Varian 200- and 300-MHz 54-mm magnets only, the dewar includes removable helium vent tubes. If the magnet dewar should quench (sudden appearance of gases from the top of the dewar) and the vent tubes are not in place, the helium gas would be partially vented sideways, possibly injuring the skin and eyes of personnel beside the magnet. During helium servicing, when the tubes must be removed, carefully follow the instructions and safety precautions given in the manual supplied with the magnet.

Caution Notices

Observe the following precautions during installation, operation, maintenance, and repair of the instrument. Failure to comply with these cautions, or with specific cautions elsewhere in Varian manuals, violates safety standards of design, manufacturing, and intended use of the instrument. Varian assumes no liability for customer failure to comply with these precautions.

CAUTION: Keep magnetic media, ATM and credit cards, and watches outside the 5-gauss perimeter from the centerline of the magnet.

The strong magnetic field surrounding a superconducting magnet can erase magnetic media such as floppy disks and tapes. The field can also damage the strip of magnetic media found on credit cards, automatic teller machine (ATM) cards, and similar plastic cards. Many wrist and pocket watches are also susceptible to damage from intense magnetism.

Refer to the manuals supplied with the magnet for the size of a typical 5-gauss stray field. This gauss level should be checked after the magnet is installed.

Caution Notices (continued)

CAUTION: Keep the PCs, (including the LC STAR workstation) beyond the 5gauss perimeter of the magnet.

Avoid equipment damage or data loss by keeping PCs (including the LC workstation PC) well away from the magnet. Generally, keep the PC beyond the 5-gauss perimeter of the magnet. Refer to the *Installation Planning Guide* for magnet field plots.

CAUTION: Check helium and nitrogen gas flowmeters daily.

Record the readings to establish the operating level. The readings will vary somewhat because of changes in barometric pressure from weather fronts. If the readings for either gas should change abruptly, contact qualified maintenance personnel. Failure to correct the cause of abnormal readings could result in extensive equipment damage.

CAUTION: Never operate solids high-power amplifiers with liquids probes.

On systems with solids high-power amplifiers, never operate the amplifiers with a liquids probe. The high power available from these amplifiers will destroy liquids probes. Use the appropriate high-power probe with the high-power amplifier.

CAUTION: Take electrostatic discharge (ESD) precautions to avoid damage to sensitive electronic components.

Wear a grounded antistatic wristband or equivalent before touching any parts inside the doors and covers of the spectrometer system. Also, take ESD precautions when working near the exposed cable connectors on the back of the console.

Radio-Frequency Emission Regulations

The covers on the instrument form a barrier to radio-frequency (rf) energy. Removing any of the covers or modifying the instrument may lead to increased susceptibility to rf interference within the instrument and may increase the rf energy transmitted by the instrument in violation of regulations covering rf emissions. It is the operator's responsibility to maintain the instrument in a condition that does not violate rf emission requirements.

Introduction

This manual is designed to help you perform liquids NMR experiments on your Varian NMR spectrometer system using VNMR software (VNMR is Varian's NMR application software package). The manual contains the following chapters:

- Chapter 1, "Advanced 1D NMR," covers working with 1D experiments at a more advanced level, such as handling multi-FID spectra, kinetics, and T₁ and T₂ analysis.
- Chapter 2, "1D Experiments," describes a number of common 1D pulse sequences for "everyday" use.
- Chapter 3, "Multidimensional NMR," describes data acquisition, processing, and display for 2D and 3D NMR. Acquisition of 4D NMR is also covered.
- Chapter 4, "Multidimensional and Advanced Experiments," describes a number of 2D, 3D, and advanced liquids experiments.
- Chapter 5, "Indirect Detection Experiments," covers indirect detection experiments, also known as *heteronuclear multiple-quantum coherence* (HMQC) experiments.
- Chapter 6, "Data Analysis," covers spin simulation, deconvolution, addition and subtraction of data, regression analysis, and pulse shape analysis.
- Chapter 7, "Pulse Analysis," describes the Pulsetool and Pandora's Box tools.
- Chapter 8, "Variable Temperature Operation," describes VT startup and operation.
- Chapter 9, "Carousel, SMS, and NMS Automation," covers automation run acquisition and operation, including user programming of the sample entry window.
- Chapter 10, "VAST Accessory Operation," covers setting up and using the VAST sample changer.
- Chapter 11, "PFG Modules Operation," describes operation using Varian Pulsed Field Gradient modules.
- Chapter 12, "PFG Modules Experiments," describes experiments using PFG modules.

Notational Conventions

The following notational conventions are used throughout all VNMR manuals:

- Typewriter-like characters identify VNMR and UNIX commands, parameters, directories, and file names in the text of the manual. For example: The shutdown command is in the /etc directory.
- Typewriter-like characters also show text displayed on the screen, including the text echoed on the screen as you enter commands. For example: Self test completed successfully.
- Text shown between angled brackets(<>) in a syntax entry is optional. For example, if the syntax is seqgen s2pul<.c>, entering the ".c" suffix is optional, and typing seqgen s2pul.c or seqgen s2pul is functionally the same.
- Lines of text containing command syntax, examples of statements, source code, and similar material might not fit the width of the page. In such cases, lines are broken at a convenient point (such as at a comma), a backslash (\) is inserted at the break, and the line is continued as the next line of text. This notation will be familiar to

C programmers. The backslash is not part of the line and, except for C source code, should not be typed when entering the line.

- Because pressing the Return key is required at the end of almost every command or line of text you type, the Return key is mentioned only in cases where it is *not* used. This convention avoids repetition of the instruction "press the Return key."
- Text with a change bar (like this paragraph) identifies material new to VNMR 6.1C that was not in the previous version of VNMR. Refer to the *Release Notes* for a description of new features of the software.

Other Manuals

This manual should be your basic source of information on intermediate and advanced liquids NMR. Other VNMR 6.1 manuals you might need to reference include:

- Getting Started
- VNMR Command and Parameter Reference
- Walkup NMR
- VNMR User Programming
- VNMR and Solaris Software Installation
- User Guide: Solid-State NMR
- User Guide: Imaging

All of these manuals are shipped with the VNMR software. Note that these manuals, other Varian hardware and installation manuals, and most Varian accessory manuals are also provided online so that you can view the pages on your workstation and print copies.

Types of Varian Spectrometer Systems

In parts of this manual, the type of system (UNITY INOVA, MERCURY-VX, MERCURY, GEMINI 2000, UNITY plus, UNITY, or VXR-S) must be considered in order to properly use the software.

- UNITY INOVA and MERCURY-VX are the current systems sold by Varian.
- UNITYplus, UNITY, and VXR-S are spectrometer lines that preceded the UNITYINOVA.
- *GEMINI 2000* is a separate line of spectrometers that preceded the *MERCURY* and *MERCURY-VX*.

Chapter 1. Advanced 1D NMR

Sections in this chapter:

- 1.1 "Working with Experiments," this page
- 1.2 "Multi-FID (Arrayed) Spectra" on page 30
- 1.3 "T₁ and T₂ Analysis" on page 35
- 1.4 "Kinetics" on page 37
- 1.5 "Diffusion Experiments/DOSY" on page 38

This chapter describes working with 1D NMR liquids experiments at a more advanced level than the manual *Getting Started*

1.1 Working with Experiments

Table 1 lists commands for working with experiments.

Table 1. Commands and Macros for Working with Experiments

```
Command
                                  Display experiment library
explib
jexp1, jexp2, ..., jexp9999
                                 Join existing experiment
md(<from_exp,>to_exp)
                                  Move display parameters between experiments
mf(<from_exp,>to_exp)
                                 Move FIDs between experiments
                                 Move parameters between experiments
mp(<from_exp,>to_exp)
Macros
cexp(<exp_dir,>exp_number)
                                  Create a VNMR experiment
delexp(exp_number)
                                  Delete an experiment
```

To obtain a view of the experiments on a system, type explib. The monitor displays the experiment library of the currently available experiment files (expl, exp2, ..., exp9999).

For each experiment, the following is displayed:

- The name of the experiment and its current size.
- The pulse sequence currently active in the experiment.
- The first 50 characters of the text file in the experiment.
- The names contained in subexperiments.

An experiment is created with the macro cexp(n), where n is a number from 1 to 9999 (e.g., cexp(4) creates the experiment exp4). An experiment is deleted by the macro delexp(n), where n is a number from 2 to 9999 (exp1 cannot be deleted). For example, entering delexp(4) deletes exp4.

Chapter 1. Advanced 1D NMR

When multiple experiments are created, an issue arises concerning how to individually work with each experiment. To handle this matter, only one experiment is allowed at a time to be currently active (i.e., in the foreground for manipulation), although background processing can be occurring in other experiments at the same time.

To make another experiment currently active (provided it has been created), you "join" it by entering the jexp command with the experiment number as an argument. For example, to join exp6, enter jexp6.

The macros jexp1 to jexp99999 have a number of convenient uses. For example, they clear the graphics window, display a parameter set appropriate to the newly joined experiment, and restart the main menu. To prevent any of these actions, use the jexp command with an experiment number as an argument. For example, perhaps you want to simultaneously display a second spectrum from another experiment, but do not want the currently displayed spectrum to be erased. In this case, enter jexp(6) rather than jexp6.

The mp, mf, and md commands move FIDs and parameters between experiments:

- mp (<n,>m) moves parameters from experiment n to experiment m, for example, mp (4,5). If n is omitted, parameters are moved from the currently active experiment to experiment m.
- mf (<n , >m) moves the last acquired FID and the associated parameters.
- md(<n,>m) moves only those "saved display" parameters associated with the commands s1 through s9.

The Workspace menu contains most of the commands for working with experiments:

Library Exp 2 Exp 3 Exp 4 Exp 5 Create New Delete

To enter this menu, click on the Main menu button, then select the Workspace button. You can also display this menu by entering the command menu('workspace'). The number of buttons in this menu depends on how many experiments have been created on your system. In the Workspace menu shown above, you see that five experiments have been created (remember that expl is always present by default and cannot be deleted).

1.2 Multi-FID (Arrayed) Spectra

Many experiments require obtaining a series of FIDs, related to each other through the variation of one or more parameters. For example, suppose it is necessary to run a series of spectra at four different temperatures: 30°C, 50°C, 70°C, and 90°C. Instead of acquiring four separate sets of data, it is possible to create an array in which the temp parameter is given four successively different values. These four subexperiments are now all treated as a single experiment. Entering go begins successive acquisition of all four experiments. One command can be used to transform all the spectra, one command to display all the spectra on the screen simultaneously, one command to plot all the spectra, and one command to save all the spectra.

Arrayed Parameters

Table 2 lists command and parameters used with arrayed parameters.

To create an array for a numeric parameter, enter the arrayed values separated by commas, (e.g., temp=30, 50, 70, 90 or pw=5, 10, 15, 20, 25). Alphanumeric parameters can

Commands		
array*	Easy entry of linearly spaced array values	
<pre>da<(param1<,param2,>)></pre>	Display acquisition parameter arrays	
<pre>ds<(index)></pre>	Display trace from arrayed 1D spectra	
dss*	Display stacked spectra	
dssa*	Display stacked spectra automatically	
dssan*	Display stacked spectra automatically without erasing	
dssh*	Display stacked spectra horizontally	
dsshn*	Display stacked spectra horizontally without erasing	
dssn*	Display stacked spectra without erasing screen	
ft*	Fourier transform 1D data	
pl*	Plot spectra	
<pre>rt<(file<,'nolog'>)></pre>	Retrieve FIDs	
<pre>svf<(file<, 'nolog'><, 'arch'>)></pre>	> Save FIDs in current experiment	
wft*	Weight and Fourier transform 1D data	
<pre>* array<(parameter<,num_steps,start,step_size)> dss<(start,finish<,step>)><,options>)> dssa<(start,finish<,step>)><,options>)> dssa<(start,finish<,step>)><,options>)> dssh<(start,finish<,step>)><,options>)> dssh<(start,finish<,step>)><,options>)> dssh<(start,finish<,step>)><,options>)> dssh<(start,finish<,step>)><,options>)> ft<(coptions,><'nf'><,start><,finish<,step>)>, ft('inverse', exp_num, expansion_factor) pl<(<start,finish<,step>><, 'int'><, 'all'><,options>)> wft<('options,><'nf'><,start><,finish<,step>)>, wft('inverse', exp_num, expansion_factor)</start,finish<,step></pre>		
Parameters		
array {string}	Parameter order and precedence	
ho {number}	Horizontal offset	
il {'y', 'n'}	Interleave arrayed and 2D experiments	
vo {number}	Vertical offset	

Table 2. Arrayed Spectra Commands and Parameters	Table 2.	Arrayed Spectra	Commands a	and Parameters
--	----------	-----------------	------------	----------------

also be arrayed. To perform two experiments in which the decoupler is off in one case and on in the other, for example, you can use dm = 'n', 'y'.

Not all parameters can be arrayed. Non-arrayable acquisition parameters include processing parameters, display parameters, and any parameter that changes the number of data points to be acquired, such as np, sw, dp, and at.

To display the values of the arrayed parameter, the da command is used. da displays all values of arrayed parameters if entered without an argument. If one or more parameters are listed as an argument, da displays only the specified parameters.

Multiple Arrays

Two or more parameters can be arrayed in an experiment. For example, an experiment to perform a series of decoupling experiments using an array of decoupler power levels and an array of decoupler frequencies might be set up with dpwr=17, 20, 23 and dof=295.1, 345.6, 507.2, 1245.5. In this example, *twelve* experiments are performed (i.e., three different values of decoupler power dpwr are used), and for each of those values, four different values of the decoupler offset dof are used.

Setting Array Order and Precedence

Whenever an array of one or more parameters is set up, the parameter array becomes important. This parameter tells the system the name of the parameter or parameters that are arrayed, and the order and precedence in which the arraying is to take place.

The string parameter array can have one of several forms:

- array='' means no parameter is arrayed (this value is two single quotation marks with no space between, not a double quotation mark).
- array='x' means parameter x is arrayed.
- array='y, x' means parameters x and y are arrayed, with x taking precedence. The order of the experiments is x_1y_1 , x_2y_1 , $...x_ny_1$, x_1y_2 , x_2y_2 , $...x_my_2$, $...x_my_n$, with a total of m×n experiments being performed.
- array='x,y' means parameters x and y are arrayed, with y taking precedence. The order of the experiments is x₁y₁, x₁y₂,...x₁y_n, x₂y₁, x₂y₂,...x₂y_n,...x_my_n, with a total of m×n experiments being performed.
- array='(x, y)' means parameters x and y are jointly ("diagonally") arrayed. The number of elements of the parameters x and y must be identical, and the order of experiments is $x_1y_1, x_2y_2, ... x_ny_n$, with n experiments being performed.

As you enter one or more arrayed parameters, array is automatically set for you. Only if you want to change the order or precedence is it necessary to enter array directly.

Interactively Arraying Parameters

Separate from the array *parameter* is the array *macro*. If you enter the array macro without an argument, an interactive mode is started in which you are asked for the following information, in this order:

- The name of the parameter to be arrayed.
- The number of values of the parameter.
- The starting value.
- The magnitude of the difference between elements in the array.

Using the information you provide, an arrayed parameter is set up. The restrictions are that only numeric parameters can be arrayed and all values of the array must satisfy the limits of the parameter.

Entering array with a parameter name as an argument, (e.g., array('pw')) still starts an interactive mode but the program only asks for the remaining three items of information.

If you enter the macro with all four pieces of information as arguments (in this order—parameter name, number of steps, starting value, and step size), array bypasses the interactive mode completely. For example, entering array('tof', 5, 1000, -50) sets the tof parameter to have 5 elements with the values in the order 1000, 950, 900, 850, 800.

Resetting an Array

Once an array is created, it is possible to change the value of a single element of the array by typing, for example, pw[2]=11.3, where the 2 enclosed in brackets indicates which element of the array to modify (array elements are counted starting at 1).

To reset an arrayed parameter to a single value, enter a single value for the parameter (e.g., pw=10). The array parameter is automatically modified to reflect this change.

Array Limitations

Regular multiple arrays can include up to 20 parameters, each of which can be a simple parameter or a diagonal array (a set of parameters), which can include up to 10 parameters. The total number of elements of all arrays is essentially unlimited $(2^{32}-1)$.

Acquiring Data

Once any parameter is an array, entering go (or related commands and macros) generates not just one, but an entire array of spectra. These spectra can then be examined either individually or as a group, as described below.

Autogain cannot be used in an arrayed experiment. You can either use gain='y', which sets the gain to the previously determined value, or set gain equal to a fixed value.

Arrayed acquisitions can be interleaved, in which a part of each experiment is done in turn rather than starting and finishing each experiment sequentially. The interleave function is controlled by the parameter il.

- If il='y', experiments are interleaved. bs transients are performed for each member of the array, followed by bs more transients for each member of the array, and so on until nt transients are collected for each member of the array. Thus, il is relevant only if bs (block size) is less than nt (number of transients).
- If il='n', all transients are acquired for the first experiment in the array, then all transients for the second experiment, etc.

Processing

The command ft or wft is used to transform all of the spectra. Both commands take the same arguments and options:

- 'acq' does not transform elements that have already been transformed.
- 'nodc' does not perform FID drift correction.
- 'nods' prevents an automatic spectral display (sane as ds command).
- 'zero' zeroes the imaginary channel of the FID before Fourier transform.

Phasing can be done on any spectrum. Only one set of phase correction parameters exists, so all spectra have the same phase at any one time (although the phase can of course be changed when examining different spectra).

Display and Plotting

The command ds(index) displays interactively the requested spectrum from the array. The index can have one, two, or three numbers, depending on the dimensionality of the spectral array. Spectra are always scaled according to the number of completed transients ct; if nt is arrayed (nt=1, 2, 4, 8), each spectrum is scaled by its *own* ct.

Other spectra display commands are dss, dssn, dssa, dssan, dssh, dsshn and dssl. These are not interactive like the ds command. They display stacked spectra in which each spectrum is offset with respect to the previous spectrum. The order of stacking can be left to right, right to left, top to bottom, or bottom to top, depending on whether the horizontal offset (ho) and vertical offset (vo) parameters are positive or negative. Some of these commands set ho and vo automatically.

The spectra display commands function as follows:

- dss displays stacked spectra using the current values of ho and vo to set the order of stacking.
- dssn displays stacked spectra the same as dss, but the graphics window is not erased before starting the display. This allows composite displays of many spectra to be created.
- dssa displays stacked spectra automatically (i.e., vo and ho are automatically adjusted to fill the screen in a lower left to upper right presentation).
- dssan displays stacked spectra automatically the same as dssa, but the graphics window is not erased before starting the display.
- dssh displays stacked spectra horizontally (i.e., vo is set to zero and ho is adjusted to fill the screen from left to right).
- dsshn displays spectra horizontally the same as dssh, but the graphics window is not erased before starting the display.
- dssl displays a label for each element in a set of stacked spectra. The label is an integer value starting with 1 and extending up to the number of spectra in the display.

The command pl plots stacked spectra with the same format as displayed by dss.

The argument syntax <(start,finish<,step>)><,options>)> is used by the dss command, variants of dss, and by the pl command. The arguments are the following:

- start is the index of the first spectra when displaying multiple spectra. It is also the index number of a particular trace to be viewed when displaying arrayed 1D spectra or 2D spectra.
- finish is the index of the last spectra when displaying multiple spectra. Because the parameter arraydim is automatically set to the total number of spectra, it can be used to set finish to include all spectra.
- step is the increment for the spectral index when displaying multiple spectra. The default step is 1.
- options can be any of the following:
 - 'all' is a keyword to display all of the spectra.

'int' is a keyword to only display the integral, independently of the value of the parameter intmod.

'top' or 'side' are keywords that cause the spectrum to be displayed either above or at the left edge, respectively, of a contour plot. This assumes that the parameters sc, wc, sc2, and wc2 are those used to position the contour plot. This option does not apply to dssa, dssan, dssh, or dsshn.

'dodc' is a keyword for all spectra to be drift corrected independently.

'red', 'green', 'blue', 'cyan', 'magenta', 'yellow', 'black', and 'white' are keywords that select a color. This option does not apply to dssa, dssan, dssh, dsshn, or pl.

'pen1', 'pen2', 'pen3', etc. specify a pen number on a plotter. This option does not apply to dss or any of its variants.

Saving and Retrieving

The commands rt and svf retrieve and save arrayed data, just like single 1D data sets. The entire data, consisting of a number of FIDs, is stored and retrieved together as a single file.

Pulse Width Calibration Step-by-Step

To illustrate using arrays, note how the following steps perform a pulse width calibration:

- 1. Set up parameters and obtain a normal spectrum of any sample. For best results, one or more intense signals should appear near the center of the spectrum.
- 2. Enter pw=5. You can use some other small value if you wish.
- 3. Enter **nt=1**.
- 4. Obtain a spectrum and phase it properly. Set **d1** to $5*T_1$.
- 5. Enter pw=5,10,15,20,25,30 ai go.

You can use some other set of suitable values for the pw array.

- 6. After the experiment finishes acquisition, enter wft dssh.
- 7. Find the experiment where the signal goes through its 180° or 360° null. Enter **da** to remind yourself of the values of the pw array.
- 8. To reset the array, enter **pw=10**.

1.3 T_1 and T_2 Analysis

One relatively common form of arrayed experiment is the inversion-recovery T_1 experiment. In this experiment, the nuclei are allowed to relax to equilibrium (d1), then inverted with a 180° pulse (p1), given a variable time to return to equilibrium (d2), and finally given a monitoring 90° pulse (pw) to measure their peak height as a function of d2. Under most circumstances, the behavior of the peak heights as a function of d2 will be exponential, and this exponential time is the T_1 .

Table 3 lists commands and parameters associated with the T_1 and T_2 analysis.

Commands		
autoscale	Resume autoscaling after limits set by scalelimits	
<pre>dels(index1<,index2,>)</pre>	Delete spectra from T_1 or T_2 analysis	
dot1*	Set up a T_1 experiment	
<pre>expl<(<options,>line1,line2,)></options,></pre>	Display exponential or polynomial curves	
fp*	Find peak heights	
pexpl*	Plot exponential or polynomial curves	
scalelimits*	Set limits for scales in regression	
t1	T_1 exponential analysis	
tls	T_1 exponential analysis with short output table	
t2	T_2 exponential analysis	
t2s	T_2 exponential analysis with short output table	
<pre>* dot1<(minimum_T1_estimate,maximu</pre>	<pre>um_T1_estimate,time)></pre>	
<pre>fp<(<'phase',><index1,index2,>)></index1,index2,></pre>		
<pre>pexpl<(<options,>line1,line2,)></options,></pre>		
<pre>scalelimits(x_start,x_end,y_start,y_end)</pre>		
Parameter		
<pre>npoint {1 to fn/4}</pre>	Number of points for fp peak search	

Table 3. T_1 and T_2 Analysis Commands and Parameters

Setting Up The Experiment

The standard two-pulse sequence is set up to perform the T_1 experiment. You can start if you wish by entering appropriate values for p1, pw, d1, and an array of values for d2.

Alternatively, you can use the dot1 macro. dot1 sets up all parameters to perform a T_1 experiment, including d1, pw, p1, nt, and an array of d2 values, based on information you enter. The three arguments that can be input are the minimum expected T_1 , the maximum expected T_1 , and the total time in hours the experiment should take. If no arguments are provided, dot1 prompts the user for the information.

Be sure that the parameter pw90 is set properly and contains the correctly calibrated 90degree pulse width, because dot1 uses this information.

Processing the Data

Once the data is acquired, process the data as follows:

- 1. Enter **wft ds(arraydim)** to display the last spectrum (or **ds(1)** for a T_2 experiment to display the first spectrum).
- 2. Phase this spectrum properly.
- 3. Select a threshold and adjust the threshold line position.
- 4. Enter dpf or dll to display a line list and locate lines for the system.
- Enter fp to measure the peak height of each peak in an array of spectra. If optional line indexes are supplied to fp as arguments (e.g., fp(1,3)), only the peak heights of the corresponding lines are measured.

The npoint parameter (if defined and set "on") determines the range of data points over which the fp command searches for a maximum for each peak.

Analyzing the Data

 T_1 and T_2 analysis is performed by the t1 and t2 macros, respectively. t1 and t2 measure relaxation times for all lines in the line listing and display an extended listing of observed and predicted peak intensities. t1s and t2s perform the same calculation as t1 and t2 but produce a shorter output, showing only a summary of the measured relaxation times.

The command expl displays exponential/polynomial curves resulting from T_1 , T_2 , or kinetic analysis. Optional input of line numbers as arguments allows displaying only selected lines. Similarly, the command pexpl plots the same curves.

The macro autoscale returns the command expl to autoscaling in which scale limits (set by scalelimits) are determined that will display all the data in the expl input file. The macro scalelimits causes the command expl to use typed-in scale limits. If no arguments are given, scalelimits asks for the desired limits. The limits are retained as long as an expl display is retained.

To delete spectra from the tl or t2 analysis (or from tls or t2s), enter dels(index1<, index2>...). This command deletes the spectra selected by the indexes from the output file fp.out of the fp command used by the tl or t2 analysis. Spectra can be restored by rerunning fp.

Exponential Analysis Menu

Most of the commands for working with T_1 and T_2 analysis are available by clicking on Main Menu button, followed by the Analyze button, and then the Exponential button. The following menu, called the Exponential Analysis menu, is displayed

T1 Proc T1 Analysis T2 Proc T2 Analysis Plot Print Return

This menu can also be displayed by entering the command menu('tlt2anal').

T₁ Data Workup: Step-by-Step

The following procedures accomplish the same result.

Using Menus

1. From the Permanent menu, click on Main Menu > File > Set Directory > Parent

The text output window shows a list of directories (entries with a slash as last character) and files (if any). The status window (at the top of the screen) shows the directory path you are in currently.

- 2. Click the **Parent** button as many times as necessary on until the message "Directory now /" appears in the status window.
- 3. Click the mouse on vnmr/ so it turns to inverse video, then click on the **Change** button. You are now in the vnmr directory.
- 4. Click the mouse on fidlib/ so it turns to inverse video, then click on the **Change** button. You are now in the fidlib directory.
- 5. Click the mouse on the tldata.fid so it turns to inverse video, then click on the buttons **Return** > **Load**.
- 6. Click on Analyze > Exponential > T1 Proc > T1 Analysis.

The T_1 analysis appears in the text output window.

Using Commands

- 1. Enter rt('/vnmr/fidlib/tldata.fid').
- 2. Enter wft dssh full ds(arraydim) aph.
- 3. Click on Next > Th. Use the left mouse button to set the threshold.
- 4. Enter dll fp t1 center expl.

1.4 Kinetics

The arraying capability of the VNMR software provides for the acquisition of data for the study of kinetics. Table 4 lists commands and a parameter for kinetics analysis.

Setting Up the Experiment

Usually, the best procedure is to array the preacquisition delay parameter pad. For example, if pad=0, 3600, 3600, 3600, 3600, the system acquires the first spectrum immediately (pad[1]=0), waits 3600 seconds (pad[2]=3600), acquires the second

Commands	
kind	Kinetics analysis, decreasing intensity
kinds	Kinetics analysis, decreasing intensity, short form
kini	Kinetics analysis, increasing intensity
kinis	Kinetics analysis, increasing intensity, short form
Parameter pad {number, in sec}	Preacquisition delay

Table 4. Kinetics Analysis Commands and Parameters

spectrum, waits another 3600 seconds, etc. Because 3600 seconds is 1 hour, this inserts a wait of one hour between acquisitions. After all the spectra have been obtained, they are processed much like T_1 or T_2 data.

Processing the Data

If the signal decreases exponentially with time, the output is matched to the equation I=A1*EXP(-T/TAU)+A3. The analysis is done by the macro kind, or by macro kinds if a short output table is desired.

If the signal increases exponentially with time, the output is matched to the equation I=-A1*EXP(-T/TAU)+A3-A1 with analysis done by the macro kini, or by the macro kinis for a short output table.

Kinetics Step-by-Step

The following steps are typical in processing a kinetics experiment:

- 1. Enter wft dssh full ds aph.
- 2. Click on Next > Th. Use the left mouse button to set the threshold.
- 3. Enter **dll fp**.
- 4. Enter kind, kini, kinds, or kinis, as appropriate (see Table 4).
- 5. If desired, adjust sc, wc, sc2, and wc2 by entering **center** or **full**.
- 6. Enter expl.


1.5 Diffusion Experiments/DOSY

Software for diffusion measurements using a pulsed Z-gradient probe is an optional feature of Varian NMR systems. This section describes the pulsed gradient experiment and its analysis. Table 5 summarizes the commands used in diffusion experiments. See the *PFG installation* manual concerning installing the Performa PFG module, and see the *VNMR and Solaris Software Installation* manual for instruction on loading the diffusion software.

Pulsed Gradient Experiments

The pulse sequences diagrammed in Figure 1 illustrate the two pulsed gradient experiments. These experiments are fully described in the literature by Stejskal and Tanner (Stejskal, E. O.; Tanner, J. E. *J. Chem. Phys.* **1965**, *42*, 288-292, and Tanner, J. E. *J. Chem. Phys.* **1970**, *52*, 2523–2526).

Commands		
analyze*	Generalized curve fitting	
expfit*	Make least-squares fit to polynomial or exp. curve (UNIX)	
expl	Display exponential or polynomial curves	
<pre>expladd(integral_region)</pre>	Add another diffusion analysis to current display	
<pre>pexpladd(integral_region)</pre>	Add another diffusion analysis to current plot	
pge	Convert parameter set to PGE pulse sequence	
pge_calib	Calibrate gradient strengths for PGE pulse sequence	
<pre>pge_data(array_index)</pre>	Extract data from single element of PGE pulse sequence	
pge_output	Output results from PGE pulse sequence	
pge_process	Automated processing of data from PGE pulse sequence	
pge_results*	Calculate diffusion constant for integral region	
<pre>pge_setup<('no')></pre>	Set up gradient control parameters for PGE pulse sequence	
<pre>* analyze('expfit",xarray,<,option,option,>)</pre>		
expfit options <analyze.inp>analyze.list (UNIX)</analyze.inp>		
pge_results(integral_region<,reference_region>)		



pge — Stimulated Echo Experiment (stim='y')

Figure 1. PGE Pulse Sequences

|

A pulse sequence named PGE and its associated help file are provided. A parameter set /vnmr/parlib/pge is used as the master parameter set. Table 6 lists the new gradient control parameters.

Name	Description	Units	Used In
stim	select normal echo or stimulated echo experiment	'n','y'	pulse sequence
d0-d5	acquisition delays	sec	pulse sequence
p1,p2,pw	RF pulse widths	msec	pulse sequence
g1,g2	gradient pulse widths	msec	pulse sequence
grad_cw	continuous-wave gradient amplitude	gauss/cm	pulse sequence
grad_p1	pulsed gradient amplitude	gauss/cm	pulse sequence
grad_p2	pulsed gradient amplitude	gauss/cm	pulse sequence
dac_cw	continuous-wave gradient amplitude	dac units	pulse sequence
dac_p1	pulsed gradient amplitude	dac units	pulse sequence
dac_p2	pulsed gradient amplitude	dac units	pulse sequence
g_max	maximum gradient value		macro pge_setu
g_min	minimum gradient value		macro pge_setu
g_steps	number of gradient values		macro pge_setu
g_array	type of gradient array		macro pge_setu
nt_array	type of nt array		macro pge_setu
nt_first	first nt value		macro pge_setu
nt_fract	fractional nt increment		macro pge_setur

Table 6. Gradient Control Parameter Set

The parameters p2 and d5 are only used for the stimulated echo experiment (stim='y'). The parameter p1 defines the first rf pulse, pw defines the last rf pulse. Delay d1 occurs before the first rf pulse.

With the current pulse sequence, grad_pl defines the amplitudes of both pulsed gradients: gl and g2. Unused parameters grad_p2 and dac_p2 are provided if you wish to modify the pulse sequence to use different gradients. These parameters are unused by the pge pulse sequence that is provided.

The PGE pulse sequence checks the average power output of the gradient coils for a safe value. Maximum value is 100 gauss/cm for a continuous wave gradient. This value is hard-coded into the pulse sequence by defining a numerical value in the pulse sequence for the variable dutycycle. The value is calculated from the following expression:

dutycycle =
$$\sqrt{\frac{(\text{grad}_{sw})^2 \times \text{time} + (\text{grad}_{p1})^2 \times (\text{g1} + \text{g2})}{\text{time}}}$$

where time=g1+g2+d0+d1+d2+d4+d5+p1+p2+pw+at.

The user should change this to fit the needs of a specific probe. Comments in the pulse sequence file /vnmr/psglib/pge.c contain more details about the sequence.

Pulsed Gradient Experiment Setup

There are two ways to retrieve a set of parameters for the pulsed gradient experiment:

- Enter rt('/vnmr/parlib/pge'). This command returns a complete parameter set into the active experiment.
- Enter pge(macro_name). This appends the required parameters to those already present in the active experiment.

For the data from a diffusion experiment to be analyzed properly, it is necessary to define integral regions in the Fourier transformed spectra. The simplest way to define these regions is to run the pulsed gradient experiment as a simple echo experiment $(grad_pl=0)$: first, to check that the values of τ and 2τ start the data acquisition at the top of the echo, and second, to define integral regions in the experiment where the pulsed field gradient experiment will be subsequently done. The phase parameters (rp and tp) should also be determined for this spectrum. If data acquisition was started at the top of the echo, the lp phase parameter should be zero. The easiest way to define integrals is by using the RESETS button of the ds command.

Once the parameters are in the experiment, those associated with gradient control can be set by entering pge_setup. This macro has a single optional argument to turn off interactive questioning; any argument can suffice (e.g., pge_setup('no')).

The macro pge_setup performs the following three tasks.

1. Sets the gradient (grad_p1) array.

If the interactive mode is used, pge_setup prompts for the values of g_max, g_min, g_steps, and g_array. These parameters are used to calculate the gradient amplitude array. Manual override is provided by typing in each value. The value of g_array can be set to 'linear' if the gradient values are equally spaced between g_min and g_max, or set to 'square' if the square of the gradient values is equally spaced between the square of g_min and the square of g_max.

2. Sets the number of transients (nt) array.

If the interactive mode is used, pge_setup prompts for the values of nt_array and nt_first. The values of nt_array are 'same' and 'fraction'. The parameter nt_first is used to set to the value of the first nt array element. If 'same' is selected, all elements of the nt array are set to nt_first. If 'fraction' is selected, nt_fract is set so that elements of the nt array are calculated according to the equation:

nt_i=nt_first x [1 - i x (1 - nt_fract) / g_steps]²

Manual override is provided by typing in each value.

3. Does necessary housekeeping.

The array parameter is set so that nt and dac_p1 form a "diagonal array." The time macro is executed to display the experiment duration. The average gradient level of the particular parameter combination is checked when acquisition is attempted. The wexp parameter is set to 'pge_process' to perform appropriate data processing at the end of the experiment.

Other acquisition parameters can be altered by typing in new values. Once a "good" set of parameters is entered, it may be saved for future recall with the svp (save parameters) command. The automatic processing of the diffusion data following data acquisition is initiated with the au command if wexp is set to 'pge_process'.

Gradient Calibration

Calibration constants, which relate DAC values (in DAC units) to resulting values (in gauss/cm), are stored in conpar with names grad_cw_coef and grad_p_coef. These coefficients are in units of (gauss/cm)/(dac unit). The DAC values are whole numbers while gradient values are real numbers that may contain fractional parts.

Whenever the dac_x or grad_x parameters are changed, where x is cw, p1, or p2, macros are available to adjust the dependent parameter, taking into account possible minimum and maximum values and housekeeping in case the parameters are arrayed. The mathematical relationships are defined as follows:

- When dac_x is changed: grad_x=dac_x*grad_x_coef
- When grad_x is changed: dac_x=grad_x/grad_x_coef grad_x=dac_x*grad_x_coef

The second step taken, when $grad_x$ is changed, is necessary because the calculation of dac_x is rounded to the nearest integer, which necessitates that $grad_x$ then be recalculated so that it corresponds to dac_x .

A macro pge_calib is provided to assist in the calibration of acquisition parameters that control the gradient power levels (i.e., the DAC values). After phasing and selecting the integral region of the standard sample, run a pulsed gradient experiment and process data with the pge_process macro. The calculated diffusion constant is displayed in the text window. Then run the pge_calib macro to recalculate the coefficient. This macro resets the set of data, followed by processing with the pge_process macro, and should now give the diffusion constant that was selected with the pge_calib macro. After running pge_calib, run pge_setup again to calculate DAC values with the new coefficient.

Data Reduction

The pge_process macro performs several tasks (see Figure 2), calling macros as appropriate.

- 1. Transforms the data with the specified weighting function.
- 2. For each array element (i.e., spectrum), writes into a text file the following information: the gradient value for that spectrum, the integral amplitudes, and the spectral parameters needed to reproduce the data.

The file is stored in the current experiment directory and its name reflects which element of the array this information pertains to (e.g., info_1 to info_n). The macro pge_data is called by pge_process to do this for each array element.

The pge_data macro has the element number passed as an argument. This feature allows an operator to manually adjust the spectral parameters for a single element and then invoke the pge_data macro to update the raw information in the appropriate text file.

A sample information file for the second element of an array is the following: Gradient amplitude for spectrum 2 is 63.9715 Spectral parameters: rp= 344.604 lp= 0 lvl= 0 tlt= 0 sp= -5000 wp= 10000 is= 691.5 ins= 1 fn= 32768 2 Integral Regions Value 7440186 -46.5117 4.58914 -325.582 -1604.65 36.8508




Figure 2. Data Processing Macros Flowchart

To manually correct the data from a single element of the data array, the following sequence can be used: select the desired array element with ds (e.g., ds(5)); manually adjust the phase or drift corrections; then, type pge_data(5) to write the new integral information into the info_5 text file.

3. The pge_process macro calls the pge_results macro.pge_results collects the necessary information from the info_n files and constructs a text file analyze.inp that is used to calculate the diffusion coefficient. pge_results requires a single argument that indicates which integral region to use when recalculating the diffusion coefficient with the information from info_n files.

The pge_process macro calls pge_results once for each of the defined integral regions. Figure 3 is a sample analyze.inp text file.

List of 8 x-y data	
pairs	
diffusion data for	\setminus
integral region 2	
Grad*Grad ln(Amp)	
0 4.23277	
4092.35 3.60688	
8192.57 3.01594	
12305.2 2.37857	
16369.4 1.72547	
20619.8 1.13322	

Figure 3.	Sample analyze.inp	
File		

The diffusion coefficient is calculated using Equation 1:

$$\begin{split} \mathbf{A}_{i} &= \mathbf{A}_{0} \exp -\gamma^{2} \mathbf{D} \bigg\{ \frac{2}{3} t^{3} \times \text{grad}_{cw}^{2} + \delta^{2} \bigg(\Delta - \frac{1}{3} \delta \bigg) \text{grad}_{pli}^{2} \\ &- \delta \bigg[(t1^{2} + t2^{2}) + \delta (t1 + t2) + \frac{2}{3} \delta^{2} - 2\tau^{2} \bigg] \text{grad}_{pli} \times \text{grad}_{cw} \bigg\} \end{split}$$

where A_i is the observed integral value of an NMR resonance for the *i*th element of the gradient array. A_0 is the integral value of an NMR resonance just after the first 90° pulse in the pulse sequence. Δ is the self-diffusion coefficient. In the case of the two-pulse echo sequence, the variables τ , δ , Δ , t1, and t2 (refer to Figure 1) are calculated from the pulse sequence variables as follows:

$$\tau = \frac{(p1)}{2} + d0 + g1 + d2 + \frac{pw}{2}$$

$$\delta = g1$$

$$\Delta = g1 + d2 + pw + d3$$

$$t1 = \frac{p1}{2} + d0$$

$$t2 = d4$$

In the case of the stimulated-echo sequence, the equation for Δ is as follows:

D = g1 + d2 + p2 + d5 + pw + d3

Equation 1 can be recast as the following, which becomes Equation 2:

 $ln(A_i) = ln(A_0) + D \times C0 + D \times C1 \times \texttt{grad_pl}i + D \times C2 \times \texttt{grad_pl}i \times \texttt{grad_pl}i$

where:

$$C0 = \frac{2}{3}\gamma^{2}\tau^{3} \operatorname{grad_cw}^{2}$$

$$C1 = \gamma^{2}\delta\left[(t1^{2} + t2^{2}) + \delta(t1 + t2) + \frac{2}{3}\delta^{2} - 2\tau^{2}\right] \times \operatorname{grad_cw}$$

$$C2 = -\gamma^{2}\delta\left(\Delta - \frac{1}{3}\delta\right)$$

The fitting program analyze accepts two arguments that instruct it to perform a polynomial fit. The selection of which polynomial to fit depends on whether grad_cw is zero. If it is zero, the second and third terms of Equation 2 vanish and a first-order polynomial $y = c0 + c1 \times x$ is used where:

$$y = In(A_i)$$

$$c0 = In(A_0) + D \times C0$$

$$c1 = D \times C2$$

$$x = \text{grad pl}i$$

Otherwise, a second-order polynomial $y = c0 + c1 \times x + c2 \times x \times x$ is used, where:

 $y = In(A_i)$ $c0 = In(A0) + D \times C0$ $c1 = D \times C0$ $c2 = D \times C2$ $x = grad_pli$

Another argument of analyze is the complete name of a text file (analyze.inp) that contains the x-y data pairs. The output of this calculation is written into a text file in the current experiment directory. The name of this text file reflects the integral region on which the analysis was performed. This name has the form regionx_results, where the x is the integral region number. Using the experimental delays, grad_cw and gamma, the pge_results macro calculates the diffusion constant and the time-zero integral amplitude from the fitting parameters c0, c1, and c2. The results of these calculations are appended to the text file that contains the least-squares analysis results.

The diffusion coefficients of both components of a two component mixture can be calculated assuming the following condition is met. It is possible to find one integral region where the NMR resonance is due to only one component of the mixture. The diffusion coefficient is calculated using that integral region with the processing already described. For integral regions, where the NMR intensity results from both components of a two component mixture, Equation 1 transforms to the following, which becomes Equation 3:

 $A_i = a0 \times exp[D \times (C0 + C1 \times grad_pli + C2 \times grad_pli \times grad_pli)] + a2 \times exp[a1 \times D \times (C0 + C1 \times grad_pli + C2 \times grad_pli \times grad_pli)]$

The diffusion coefficient *D* is available from the separate reference integral region. The constants *C*0, *C*1, and *C*2 are defined in Equation 2. The fitting parameters are *a*0, *a*1, and *a*2. In order to perform the non-linear least squares analysis of Equation 3, the $pge_results$ macro is supplied with two arguments (e.g., $pge_results(1,3)$). The first argument is the region on which to perform the analysis (just as for the single-component analysis case) and the second argument is the integral region used to get the value of *D*. The fitting parameter a0 corresponds to the time-zero integral amplitude of the reference component; *a2* corresponds to the time-zero amplitude of the other component; *a1* corresponds to the ratio of the two diffusion coefficients.

Data Display

The macro pge_output prints the experimental parameters and the results of the diffusion calculations. It also prints graphs of the line fitting results and the spectra.

As with any printing operation, the pge_output macro calls printon, does a cat of the regionx_results files, and then calls printoff. The plotting is done with the pexpl command. The analogous expl command displays graphs on the screen.

Two macros are supplied that add the results of the separate calculation to an existing graph. These are called expladd and pexpladd, for graphics display and plotting, respectively. Each requires a single argument that specifies the number of the region whose results are to be added to the existing plot or graph.

To plot or display the results of a two-component analysis, the commands pexpl and expl, respectively, are provided. For example, to plot the results of a two-component analysis, enter pexpl('square', 'log'). This command makes a plot of the square of the gradient versus the natural logarithm of the amplitude.

Variations on the pge Pulse Sequence

In addition to the basic pulse sequence pge for diffusion measurements, there are pulse sequence variations on pge:

pgeramp	Ramps gradients, unlike pge. Once the pge parameter set has been recalled, set $seqfil='pgeramp'$ to execute pgeramp. Use this pulse sequence when probe impedance is highly mismatched to the gradient amplifier output. pgeramp determines ramp length (defined in μ s). When executed, this pulse sequence determines the number of steps in ramping the gradient based on the value of tramp (default value is 200 μ s) and the gradient strength. In arrayed series of gradients, lower gradients have fewer steps and higher gradients have more steps.
g2pulramp	Analogous to g2pul except that the gradients are ramped and ramp time is determined by tramp.g2pulramp is executed by setting seqfil='g2pulramp'. It determines the number of steps in ramping the gradient based on the value of tramp (the default value is 200 µs) and the gradient strength.

DOSY Experiments

The DOSY (Diffusion Ordered SpectroscopY) application separates the NMR signals of mixture components based on different diffusion coefficients. Generally speaking, DOSY increases the dimensionality of an NMR experiment by one. In 2D DOSY the initial diffusion weighted NMR spectra are one-dimensional; adding diffusion weighting to a 2D NMR experiment such as COSY or HMQC gives 3D DOSY spectra.

The DOSY analyzes involves two steps. These steps are executed by the dosy macro.

- 1. Set up and acquire a DOSY spectrum.
- 2. Determine the diffusion coefficients for each line (or cross-peak) in the spectrum. Take line (or cross-peak) positions and diffusion coefficients and show the results in a DOSY plot.

Table 7 shows the available tools for DOSY.

DOSY Pulse Sequences

Previous DOSY pulse sequences used an unhelpful choice of parameter names. These names have been corrected, but compatibility with old data has largely been maintained. The dosy macro attempts to identify the relevant information from the parameters and the pulse sequence name; if it fails, dosy starts a dialog asking for three pieces of required information:

- The width of the gradient pulse(s) used for dephasing before diffusion.
- · The diffusion delay between dephasing and rephasing.
- · For bipolar sequences, the time between the positive and negative gradient pulses.

New sequences always start with "D" and are supplied with this version of the DOSY software. The sequences calculate the time portion of the exponent governing diffusional attenuation, storing the calculation as dosytimecubed, and the Larmor frequency of the diffusing spins, storing that calculation as dosyfrq.

The macros bppste and showdosy are obsolete in VNMR 6.1C.

Command	Function	
cleardosy	Delete any temporarily saved data in the current (sub) experiment.	
ddif	Synthesize and display DOSY plot.	
Dgcstecosy	Set up parameters for the Dgcstecosy.c pulse sequence.	
Dgcstehmqc	Set up parameters for the Dgcstehmqc.c pulse sequence.	
DgsteSL	Set up parameters for the DgsteSL.c pulse sequence.	
Doneshot	Set up parameters for the Doneshot.c pulse sequence.	
dosy	Process DOSY experiments.	
Dbppste	Set up parameters for the Dbppste.c pulse sequence.	
Dbppsteinept	Set up parameters for the Dbppsteinept.c pulse sequence.	
fbc	Apply baseline correction for each spectrum in the array.	
fiddle*	Perform reference deconvolution.	
makedosyparams	Create DOSY-related parameters (called by setup macros).	
makeslice	Synthesize 2D projection of a 3D DOSY spectrum in diffusion limits.	
redosy	Restore the previous 2D DOSY display from the subexperiment.	
setup_dosy	Start dialog to set up gradient levels for DOSY experiments.	
sdp	Show diffusion projection.	
setgcal	Set the gradient calibration constant.	
showoriginal	Restore the first 2D spectrum in a 3D DOSY experiment.	
undosy	Restore the original 1D NMR data from the subexperiment.	
* fiddle(option<,file><,option<,file>><,start><,finish><,increment>)		

Table 7. Tools for the DOSY Experiment

General Considerations

The DOSY experiments are probably the most demanding gradient sequences in NMR. In conventional coherence-pathway-selected-experiments, you can optimize the experimental conditions for a given gradient setting. However, in DOSY, very often the whole scale of available gradient power is used and the high-resolution NMR conditions still must be maintained. Convection, i.e., moving liquid columns along the sample axis (primarily due to temperature gradients), does not seriously hurt coherence-pathway-selected-experiments (apart from the obvious intensity losses), but, it can make the DOSY analysis completely useless.

DOSY pulse sequences use the gradient stimulated echo element (or one of its modifications), shown in Figure 4.




Figure 4. Gradient Stimulated Echo Element

In the DOSY experiments the strength of the diffusion encoding gradient is arrayed and the diffusion coefficients are calculated according to the Stejskal-Tanner formula:

$$S(G_{i}) = S(0) \exp(-D_i \gamma^2 \delta^2 (G_{i})^2 (\Delta - \delta/3))$$

where $S(G_{zi})$ and S(0) are the signal intensities obtained with the respective gradients strengths of G_{zi} and 0. D is the diffusion coefficient. γ is the gyromagnetic constant, δ is the gradient pulse duration and Δ is the diffusion delay.

The formula itself provides valuable hints on how to set DOSY-related parameters in different pulse sequences.

• $(\gamma^* \delta^* G_{zi})^2$ is the gradient pulse area.

Nuclei with higher γ are more sensitive to diffusion than the low- γ nuclei. If possible, observe ¹H, ¹⁹F, or at least do the diffusion encoding step on the high- γ nucleus. See "Dbppsteinept (DOSY Bipolar Pulse Pair Stimulated Echo INEPT) Experiment," page 52.

Shaping a gradient dramatically reduces its phase encoding efficiency. Although VNMR software can support gradient shaping on UNITY*plus* or ^{UNITY}*INOVA* spectrometers, no advantages are expected from using shaped gradients.

• δ is the gradient pulse duration.

During δ (and the subsequent gradient stabilization delay, gstab) the magnetization is transverse and subject to T₂ relaxation and homonuclear J-evolution. Do not use long δ values in the presence of large homonuclear couplings or short T₂ relaxation times ($\delta << T_2$ or 1/J).

- G_z is the gradient strength. Use as high values as possible provided high-resolution NMR conditions are still maintained (no phase, amplitude, and line shape distortions).
- Δ is the diffusion delay. Convection can always be an unwanted competitor to diffusion and T₁ relaxation attenuates the signal intensities. Do not use unnecessarily long diffusion delays (Δ <T₁).

Some of the previous recommendations might seem contradictory. Of course, in real cases you need to find an acceptable compromise between them.

The separation efficiency in the diffusion domain is determined by the accuracy of the measured diffusion coefficients. DOSY does not necessarily intend to get absolute diffusion coefficients (in mixtures, it is difficult to discuss "absolute" numbers anyway); the relative differences in the D values might be adequate for separation.

Note: Changing the solvent of a DOSY mixture might change the diffusion coefficients and the separation power of the method. The solvent might play a similar role in DOSY as the different columns in HPLC.

Diffusion coefficients errors can either be statistical or systematic. The most obvious source of statistical errors is inappropriate signal-to-noise ratio; therefore in DOSY experiments, relatively high S/N values must be reached even with the strongest phase encoding gradients. Systematic errors are primarily caused by instrumental imperfections (such as gradient nonlinearity over the active sample volume, phase distortions, changes in experimental lineshape as a function of gradient amplitude). Systematic errors can be minimized by careful pulse sequence design (see *Magn. Reson. Chem.* (1998), **36**: 706) and by adding a suitable internal reference to the sample (a component producing a strong, well isolated singlet peak in the spectrum) suitable for reference deconvolution (FIDDLE) when processing DOSY.

When setting up DOSY experiments, consider the following recommendations:

- Be sure that the probe parameter is set to the probe you intend to use and Probegcal has the right value (the setup macros extract the gradient strength, gcal, from the probe file and store it in the local parameter DAC_to_G.) Pulse power levels and pw90 values are also read from the probe calibration file.
- Set z0 precisely on resonance and carefully adjust the lock phase. Incorrect adjustment might cause progressive phase errors with increasing gradient power.
- Do not spin the sample.
- Use an adequate number of data points for proper spectral digitization.
- When running long experiments, use interleaved acquisition.
- To minimize temperature gradients (and convection), avoid using extreme (low and high) temperatures. For solutions with very low viscosity, you might prefer to completely switch off the VT controller.
- If you can find a substance suitable for reference deconvolution, add it to the mixture before running DOSY (in proton spectra, TMS might be an ideal candidate).

2D-DOSY Spectroscopy

The current DOSY package includes four 2D DOSY sequences:

- Dbppste
- DgcsteSL
- Doneshot
- Dbppsteinept

Setting Up 2D-DOSY Experiments

- 1. Start setting up any of the four experiments by recording a normal s2pul spectrum on the nucleus to be observed.
- 2. Calibrate (or check) pulse widths if necessary.
- 3. Before calling the setup macro, which always has the same name as the pulse sequence itself, reduce the spectral window to the region of interest.
- 4. Each sequence has a parameter called delflag. By setting delflag='y', the actual DOSY sequence is activated. Setting delflag='n' enables you to go back to the basic s2pul (Dbppste, DgcsteSL, Doneshot) or INEPT (Dbppsteinept) sequence without changing the experiment or the parameter set.
- 5. In all of the sequences, the phase encoding gradient duration is defined by the gtl parameter (the total defocusing time). Its strength is defined by the gzlvll parameter and the diffusion delay by the del parameter. The actual DOSY setup determines the proper relationship among these three parameters. The best setting primarily depends on the sample itself (e.g., solvent, viscosity, molecular size and shape, the isotope to be detected) and on the experimental conditions (e.g., temperature). Therefore, it is recommended that you use the DOSY sample to optimize the experimental parameters. For small or medium sized molecules, it might be useful to set gtl=0.002 and del=0.05 sec and to array the gradient strength:

gzlvl1=500, 5000, 15000, 20000, 25000, 30000 for Performa II gzlvl1=50, 500, 1000, 1500, 2000 for Performa I gradient systems

6. For the maximum gradient power used in the DOSY experiment, select the gzlvll value, which attenuates the signal intensities to 5% to 15% of the intensities obtained

Chapter 1. Advanced 1D NMR

with the weakest gradient pulse. If the intensity drop is not sufficient at the end of the array, you can increase del or gtl. If no signal is detected towards the end of the array, decrease del or gtl and repeat the procedure again.

- 7. Before the final setup, optimize the alfa delay to reach ideal baseline performance.
- 8. After having determined suitable values for gt1, de1, and the maximum gradient power, call the setup_dosy macro.

setup_dosy asks for the number of gradient levels, for the weakest and strongest gradient power to be used in the experiment and sets up a range of gzlvll values with their squares evenly spaced. The minimum gradient strength may be set to 0.3-0.5 G/cm. The number of different pulse areas to use depends on the range of diffusion coefficients to be covered and the balance between systematic and random errors but typically is in the range of 10 to 30. As in any quantitative experiment, there is a balance to be struck when choosing a repetition rate between signal-to-noise and accuracy. But in DOSY experiments, a delay of 1-2 T₁ suffices, provided that care is taken to establish a steady state before acquiring data. It is recommended to set ss < 0 to have steady-state pulses at every new array element and run the acquisition interleaved (il = 'y').

Each sequence is equipped with a Tcl-Tk acquisition panel, which provides direct access to parameters and setup related commands. Figure 5 shows the acquisition panel of the Doneshot sequence.

Acquisition Parameter	ers	DOSY flags	DOSY parameters	Acq & Obs
Observe nucleus	H1	👅 Phasecycle flag	Maximum selectable gradient level 32767	
Spect frequency	499.87	📕 Diffusion delay flag	Diffusion gradient level 1000	Decouplers
Spectral width	8000.0	_ _ Interleave	Diffusion gradient length (sec) 0.002	Sequence
Acquisition time	3.0		Gradient stabilization delay (sec) 0.0002	
Acq. complex points	24000	Sample	Diffusion delay (sec) 0.06	Flags & Cond
Recycle delay	1.5	Date Jul 16 1999	Unbalancing factor (~0.2) 0.2	Process
Transmitter offset	0.0	File exp	Spoiling gradient level 1000	
Obs. pulse power	61	temp. 30 C	Spoiling gradient length (sec) 0.0005	Process2
Observe pulse	5.55	Probe triax	Dosy freq.(MHz) (set by go) 499.869	Display
Block size	32	Solvent D20	dosytimecubed (set by go) 1.55039e-07	41 1 0
Steady state	-2		DAC_to_G (Gaus/cm) 0.0021872	display2
Transients	8	Channel : Nucl. Status	PFGON : mny	LCNMR/STARS
Completed transients	0	Observe : H1	Gradient Type: ttt	Text
Receiver gain	16	Decoupler : H1 - nnn		IGAC
				Spare
Clear subexperiment	Dis	splay Pulse Sequence Setup DOSY	Experiment Time Acquire	Setup EXP

Figure 5. Tcl-Tk Acquisition Panel of Doneshot Pulse Sequence

Dbppste (DOSY Bipolar Pulse Pair Stimulated Echo Experiment)




Figure 6. Dbppste Experiment

Parameters	Function	
calibflag	Correct systematic errors in DOSY experiments.	
DAC_to_G	Store gradient calibration value in DOSY sequences.	
del	Actual diffusion delay.	
delflag	y runs the Dbppste sequence.	
	n runs the normal s2pul sequence.	
fn2D	Fourier number to build up the 2D display in F2.	
gstab	Gradient stabilization delay (~200-300 us)fn2D.	
gt1	Total diffusion-encoding pulse width.	
gzlvll	Diffusion-encoding pulse strength.	

Table 8	. Dbppste	Parameters
---------	-----------	------------

DgcsteSL (DOSY Gradient Compensated Stimulated Echo with Spin Lock) Experiment




Figure 7. DgcsteSL Experiment

Table 9.	DgcsteSL	Parameters
----------	----------	------------

Parameters	Function	
calibflag	Correct systematic errors in DOSY experiments.	
DAC_to_G	Store gradient calibration value in DOSY sequences.	
del	Actual diffusion delay.	
delflag	y runs the DgcsteSL sequence.	
	n runs the normal s2pul sequence.	
fn2D	Fourier number to build up the 2D display in F2.	
gstab	Gradient stabilization delay (~200-300 us).	
gt1	Total diffusion-encoding pulse width.	
gz_alt	Flag to invert the gradient sign on alternating scans (default is n).	
gzlvl1	Diffusion-encoding pulse strength.	
prg_flg*	y selects purging trim pulse (default).	
	n omits purging pulse.	
prgpwr*	Power level for the purge pulse.	
prgtime*	Purging pulse length (in second).	
tweek	Tuning factor to limit eddy currents. It can be set between 0 and 1, usually set to 0.0.	
	pulse can effectively eliminate the dispersion signal components.	

Be careful not to create convection in the sample by the trim pulse.

The "Oneshot" Experiment

The total gradient power transmitted to the sample remains independent on the phase encoding gradient power. Although the sequence design makes phase cycling unnecessary and, unlike other DOSY sequences, the Doneshot sequence can be run with a single transient per array element, it is recommended to turn on the cyclops cycle: phasecycleflag='y'




Figure 8. Oneshot DOSY Experiment

The lock refocusing gradient is determined by kappa and gzlvll. The dummy heating gradients are automatically adjusted by the sequence. For the maximum gradient power available in the experiment, use:

gzlvl_max > gzlvl1*(1+kappa)

Table 10. Oneshot DOSY Parameters

Parameters	Function	
calibflag	Correct systematic errors in DOSY experiments.	
DAC_to_G	Store gradient calibration value in DOSY sequences.	
del	Actual diffusion delay.	
delflag	y runs the Doneshot sequence.	
	n runs the normal s2pul sequence.	
fn2D	Fourier number to build up the 2D display in F2.	
gstab	Gradient stabilization delay (~200-300 µs).	
gt1	Total diffusion-encoding pulse width.	
gt3	Spoiling gradient duration (sec).	
gzlvl1	Diffusion-encoding pulse strength.	
gzlvl3	Spoiling gradient strength.	
gzlvl_max	Maximum gradient strength accepted.	
	(32767 with Performa II or III, 2047 with Performa I)	
kappa	Unbalancing factor between bipolar pulses as a proportion of gradient strength (recommended:~0.2).	
phasecycleflag	Flag to turn on or off the phase cycle.	

Dbppsteinept (**D**OSY **B**ipolar **P**ulse **P**air **St**imulated **E**cho **INEPT**) Experiment

This sequence uses the higher "resolving power" of the wide 13 C chemical shift range, while the phase encoding and decoding step is more effectively done on the 1 H magnetization.

Parameters	Function
calibflag	Correct systematic errors in DOSY experiments.
DAC_to_G	Store gradient calibration value in DOSY sequences.
del	Actual diffusion delay.
delflag	y runs the Dbppsteinept sequence.
	n runs the normal INEPT sequence.
fn2D	Fourier number to build up the 2D display in F2.
gt0	Gradient duration of the spoiler gradient (for satflag='y').
gt1	Total diffusion-encoding pulse width.
gzlvl0	Spoiling gradient power.
gzlvl1	Diffusion-encoding pulse strength.
jlxh	One-bond H-X coupling.
mult	Carbon multiplicity.
	2 selects CHs (doublets).
	3 gives CH2s down and CHs CH3s up. 4 enhances all protonated carbons.
nn	90° hard ¹ H pulse.
qq pplvl	Decoupler power level for pp pulses.
satflag	y gives a grad-90(X) grad sequence during the diffusion delay to destroy the
Sacitay	X-magnetization not originated from INEPT transfer.

Table 11. Dbppsteinept Parameters

Processing 2D-DOSY Experiments

After DOSY data has been acquired, it must be processed to give a 2D DOSY spectrum. Processing data involves the following stages:

- 1. Reference deconvolution with the command fiddle (optional, but useful if there is a suitable reference line that slowly diffuses).
- 2. Baseline correction with the macro fbc (also optional, but strongly recommended).
- 3. Extraction of diffusion data from the spectra and synthesis of a 2D DOSY spectrum with the macro dosy.

fiddle

The fiddle program enables you to use reference deconvolution to correct the line shapes, frequencies, phases, etc. of the signals caused by instrumental imperfections. Full instructions for its use are given in section 6.2 "Deconvolution" on page 205 and in the *Command and Parameter Reference* manual. Reference deconvolution of DOSY spectra removes systematic errors caused by the disturbance of the magnetic field and field/ frequency lock caused by gradient pulses. It is best to use fiddle with the writefid option to store the corrected data, recall the corrected data, and set all the weighting parameters to n before Fourier transforming and proceeding to the next step.

fbc

The fbc macro applies bc type baseline correction to all spectra in an array. Use the partial integral mode to set integral regions to include all significant signals. Leave blank as large an area of baseline as you possible can; this minimizes systematic errors in diffusion coefficient fits caused by baseline errors.

dosy

This macro uses the commands dll and fp to determine the heights of all signals above the threshold defined by the parameter th. Then it fits the decay curve for each signal to a Gaussian using the program dosyfit, storing a summary of all diffusion coefficients and their estimated standard errors and various other results as follows:

In the directory \$HOME/vnmrsys/Dosy:

- diffusion_display.inp
- general_dosy_stats
- calibrated_gradients
- fit_errors
- diffusion_spectrum

In the current experiment:

• A second copy of diffusion_display.inp.

The dll program is limited to handling 512 lines, so very crowded spectra might need to be processed in sections by appropriately choosing sp and wp. dosy then runs the command ddif to synthesize the 2D DOSY spectrum.

The peak representation and the accuracy of the peak heights might increase with higher digital resolution, i.e., zero-filling the FIDs (fn>np) can occasionally improve the quality of the DOSY data. In extreme cases, even fn=512k is allowed by the software. Building up a 2D data set (and a 2D display) with this data size would not make sense; therefore a new parameter, fn2D (with a maximum limit of 64k), is introduced in the 2D-DOSY sequences, replacing fn when setting up the 2D display.

Note: The 2D DOSY display is set up in the same experiment where the data processing takes place.

The synthesized spectrum contains fn1/2 traces in the diffusion domain (f1), and fn2D real data points in the spectral domain (f2); fn1 is limited to the range 128-1024. Normally fn2D of 16k suffices. If fn2D*fn1 is too large, spectral synthesis and display will be slow and/or VNMR might run out of disk space.

Note: After displaying a 2D spectrum, the variable ni is set to fn1/2 (this setting is required by dconi). So if more data is to be acquired or the sequence is to be displayed (dps) you must set ni back to zero.

By default, dosy uses all the experimental spectra and covers the whole diffusion range seen in the experimental peaks. Either one or three arguments, shown in Table 12, can be supplied to dosy to change the defaults.

Command	Function
dosy('prune')	Start a dialog to allow one or more spectra to be omitted from the analysis.
dosy(d1,d2)	d1 and d2 are numbers causing the diffusion range of the synthesized spectrum to be limited to $d1*10^{-10}$ m ² /sec and $d2*10^{-10}$ m ² /sec.
dosy('prune',d1,d2)	Combine the previously described arguments.

Table 12. dosy Commands

The message Systematic Gz deviations indicates that the random errors in the data are sufficiently small. It might be worthwhile correcting for the small systematic errors in the field gradients, produced by the spectrometer hardware, by using the decay curves of selected signals to provide an internal calibration of the relative gradient strengths. To correct for systematic gradient errors, do the following procedure:

- 1. Set the display/threshold parameters to select a few strong, well-resolved signals, which are known to arise from single species (i.e., the signals are not composites of overlapping signals from species with different diffusion coefficients). Enter dosy to perform the analysis a first time, readjust the display and threshold to contain all the signals of interest.
- 2. Enter **undosy calibflag='y' dosy**. The second analysis uses the shapes of the decay curves in the first analysis to correct for systematic errors. Remember to set calibflag back to 'n' if you wish to stop using the internal gradient calibration.

WARNING: If the argument prune was used for the initial run of the dosy macro, you must ensure that the same increments are deleted in the second run. Use undosy calibflag='y' dosy('prune') and specify the same increment number(s).

The two-dimensional DOSY display (and plot) is constructed by taking the bandshape of a given signal from the first (lowest gradient area) spectrum and convoluting it in a second dimension with a Gaussian line centred at the calculated diffusion coefficient and with a width determined by the estimated error of the diffusion coefficient obtained from the fitting process.

To extract spectra of the mixture components separated along the diffusion axis, select the region of interest using the two cursors in the interactive 2D display (dconi) mode and click on Proj (projection) and Hproj(sum) (horizontal projection). The spectrum can be plotted with the Plot menu.

When processing 2D DOSY spectra, you might find the Tcl-Tk process panel, shown in Figure 9, useful..

Proces	ssing		Linear Prediction	L	Acq & Obs
Weighting	Interactive	•	🔶 back	Calculate full DOSY Spectrum	
line broad	2	F	🔶 forward		Decouplers
sinebell	0		coefs	Calculate full DOSY with dialog	Sequence
shift	0		basis pts		Flags & Cond.
gaussian	0.68		starting at		D
shift	0		predicted pts	Calculate partial DOSY Spectrum	Process
	Transform		starting at	Lower diffusion limit: 6 m2/sec	Process2
Acquired dat	ta points: 1	6384		Upper diffusion limit: 7 m2/sec	Display
Fourier No.	(1D, fn):	32768		Calculate partial DOSY with dialog	
Fourier No.	(2D, fn2D):	32768			display2
Fourier No.	(DIFF, fn1):	256	▲ ▼ □ Calibration f	Lag	LCNMR/STARS
				Recall Original NMR Data	Text
	Process		FIDDLE (TMS)		
Baseline	Correct ALL	Spectra	FIDDLE (NO-TMS) Recall Diffusion Display	Spare
					Setup EXP

Figure 9. Tcl-Tk Process Panel for 2D_DOSY Pulse Sequences

WARNING: Do not process the data with the dosy macro until the acquisition has been completed.

sdp

The sdp macro (show diffusion projection) displays the integral projection of a DOSY dataset onto the diffusion axis. sdp uses the file userdir+'/Dosy/diffusion_spectrum' as input for the sdp command. Only use sdp in an experiment in which data can be overwritten because it modifies parameters such as sw and at.

3D-DOSY Experiments

3D DOSY adds a diffusion domain to "conventional" 2D experiments such as COSY or HMQC. The package contains sequences for DOSY-COSY (Dgcstecosy) and DOSY-HMQC (Dgcstehmqc), but it is easy to add diffusion encoding to many other 2D experiments. The 3D DOSY sequences provide better resolving power than the 2D counterparts (the probability of overlapping cross-peaks in 2D is much lower than the probability of overlapping lines in 1D proton detected experiments), at the expense of data size and experiment time.

An arrayed set of 2D experiments is performed using different values of gradient strength (gzlvll). The data is doubly Fourier transformed, and the first 2D spectrum is used to manually define 2D integral regions. dosy analyzes then fits the integral volumes in successive increments to Gaussians and synthesizes 2D integral projections of the 3D dataset between defined diffusion limits. Full 3D display is not implemented; although with patience, you can achieve a similar effect by performing a series of projections.

Setting Up 3D-DOSY Experiments

- 1. Make sure that the "conventional" parameters of the COSY / HMQC experiment, such as pulse widths, transmitter offset, spectral window are correctly set.
- 2. As with 2D DOSY, try to find suitable lower and upper bounds for the gradient strength gzlvll. There is no need to run 2D experiments for this purpose; the first increment from a 2D run is normally adequate (ni=1).
- 3. In a COSY experiment with higher quantum filter (qlvl>1), the first increment does not contain signals. Set the incremented delay (d2) to 0.05-0.1 during the gradient optimization process. Set d2 back to zero when starting the DOSY-COSY experiment.
- 4. Use the setup_dosy macro to set up an array of 5 to 10 different gzlvll values.

The full 3D experiments then can be acquired. Note the total experiment time when choosing the number of gzlvll values, ni and nt. Both sequences are equipped with a Tcl-Tk acquisition panel.

Dgcstecosy (**D**OSY **G**radient **C**ompensated **S**timulated **E**cho **COSY**) Experiment (AV mode)





Figure 10. Dgcstecosy (AV Mode) Experiment

Using the ('t2dc') argument to wft2d can be useful.

Table 13. Dgcstecosy Parameters

Parameters	Function
calibflag	Correct systematic errors in DOSY experiments.
DAC_to_G	Store gradient calibration value in DOSY sequences.
del	Actual diffusion delay.
gt1	Total diffusion-encoding pulse width.
gzlvl1	Diffusion-encoding pulse strength.
gstab	Gradient stabilization delay (~200-300 us).
tweek	Tuning factor to limit eddy currents.
	(can be set between 0 and 1, usually set to 0.0).
gt2	Gradient duration for pathway selection.
gzlvl2	Gradient power for pathway selection.
qlvl	Quantum filter level (1= single quantum, 2=double quantum).

Dgcstehmqc (**D**OSY **G**radient **C**ompensated **S**timulated **E**cho **HMQC**) Experiment (AV mode)



Process N-type data with wft2d(1) to process the first 2D experiment.

Figure 11. Dgcstehmqc Experiment (AV Mode)

Parameters	Function	
del	Actual diffusion delay.	
gt1	Total diffusion-encoding pulse width.	
gzlvl1	Diffusion-encoding pulse strength.	
gt2	First coherence pathway selection gradient in HMQC.	
gzlvl2	Gradient power for gt2.	
gt3	Second coherence pathway selection gradient in HMQC.	
gzlvl3	Gradient power for gt3.	
gt4	Refocusing gradient in HMQC.	
gzlvl4	Gradient power for gt4.	
gstab	Gradient stabilization delay (~200-300 us).	
pwx	90° X pulse.	
pwxlvl	Power level for pwx.	
j	One-bond H-X coupling.	
mbond	Flag to select multiple-bond correlation (HMBC).	
taumb	Delay for magnetization transfer for mbond= 'y'.	
c180	Flag to use the 180° X pulse a composite pulse.	

Table 14.	Dgcstehmqc Parameters
-----------	-----------------------

If gzlvl4 has the same sign as gzlvl2 and gzlvl3 (N-type selection), use wft2d to process the data.

If gzlvl4 has opposite sign to gzlvl2 and gzlvl3 (P-type selection), use wft2d('ptype',1) to process the data. Using the ('t2dc') argument to wft2d can be useful.

Before using the sequence for the first time the coherence pathway selected gradients needs to be calibrated for a given probe and gradient amplifier.

The choice of decoupling method in the DOSY-HMQC experiment is crucial, as even relatively low values of dpwr can cause sufficient convection currents to invalidate DOSY results. On UNITY*plus* or UNITY*INOVA* systems equipped with a PPM module in the ¹³C channel adiabatic decoupling schemes (WURST, STUD) is recommended.

Processing 3D-DOSY Experiments

In order to analyze 3D results, it is necessary to manually define the individual signal regions in the 2D spectrum.

- 2D Fourier transform the first increment of the 3D data set (i.e., the increment with the lowest gzlvll value), using proper weighting functions in both dimensions: wft2d(1) for COSY wft2d('ptype',1) for HMQC
- 2. Correctly set vs2d and th, then define the signal regions in the first spectrum using the standard ll2d command and its options (e.g., 'reset', 'volume', 'clear', 'combine') The options are easily accessible via the dconi/Peak/ Edit menu set.
- 3. Include all the components of a given multiplet (cross-peak) in a single integral region, provided that there is no contamination by other signals. Grouping signals in this way maximizes the signal-to-noise ratio available for data fitting.

This step offers you the unique opportunity to exclude apparent spectral artifacts (t1-noise, decoupling sidebands, spurious peaks, etc. from the DOSY analysis.) Because the manual peak selection is probably the most boring and time-consuming step of the whole procedure, after it is completed, it is worth storing the file (using the command ll2dbackup) in the same directory where the FID is stored for later processing.

4. After the signal regions have been defined, enter the command dosy.

The macro extracts the volume of each region for every value of gzlvl1 (this step involves, among other things, as many 2D Fourier transforms as there are gzlvl1 increments). dosy then fits the volumes as functions of gzlvl1, returning with a display in which each signal region is labelled with its diffusion coefficient (10^{-10} m²/sec) and with its standard error in brackets. The coefficients are automatically displayed when the dosy macro is completed using the label facility of the ll2d command. Thus, 6.05(0.05) means a diffusion coefficient of $6.05*10^{-10}$ m²/sec (+/- $0.05*10^{-10}$ m²/sec). The 2D spectrum on which the display is based is that of the first 2D increment of the 3D experiment. A copy of the diffusion results is available from the file userdir+/Dosy/diffusion_display_3D.inp. This file contains three columns:

- The peak number (as obtained by ll2dmode='nynn')
- The diffusion coefficient
- · The standard error
- 5. The display of diffusion coefficients as numbers on the screen can result in very crowded display. You can change the type of information shown by using the ll2dmode parameter (for details see the *Command and Parameter Reference Manual* and Figure 12).
- 6. In order to make the analysis easier, use sdp to obtain the integral projection of the 3D data set onto the diffusion axis. You can use this diffusion spectrum to choose suitable diffusion regions for which to examine 2D projections of the 3D DOSY data.

WARNING: Warning: Be sure to use sdp in an experiment in which data can be overwritten!

 In the experiment containing the 3D data, enter the command: makeslice(d1,d2)

where d1 and d2 are the diffusion limits (in units of $10^{-10} \text{ m}^2/\text{s}$) between which the 2D projection of the 3D DOSY spectrum is required. The makeslice macro builds the slice and displays it after a few seconds. makeslice uses, among other things, the diffusion information in the file userdir+'/Dosy/diffusion_display_3D.inp'.

- 8. To return to the original spectrum, enter showoriginal. This command reverts to the original 2D spectrum for the first value of gzlvl1.
- 9. You can repeat the sequence makeslice showoriginal as needed with different diffusion values (or slice thicknesses), but you must use showoriginal in between the display of two slices.

Both sequences are equipped with a Tcl-Tk Process2 panel, shown in Figure 12, providing access to necessary functions and parameters to process 3D DOSY data.

Retransform the 1st 2D spectrum	Acq & Obs Decouplers
Calculate DOSY spectrum	Sequence
Labeling options (cross,number,box,Diff.coef.): mmyn	Flags & Cond.
Calculate DOSY spectrum with dialog Redisplay spectrum	Process
	Process2
Show 2D pojection within the specified limits:	Display
Lower diffusion limit: 4.1 m2/sec Upper diffusion limit: 4.7 m2/sec	display2
	LCNMR/STARS
	Text
Show ORIGINAL 2D (first spectrum)	Spare
	Setup EXP

Figure 12. Tcl-Tk Process2 Panel for the 3D-DOSY Pulse Sequences

Sample FIDs to Practice DOSY Processing

The package includes a few 2D and 3D FIDs (in /vnmr/fidlib/Dosy) to practice DOSY processing. Except for Doneshot, every pulse sequence has an example. Processing Doneshot data is not unique; it requires the same procedure as the Dbppste or the DgcsteSL sequences.

Dbppste.fid

The sample is a mixture of three dipeptides (Phe-Val, Phe-Glu, and Phe-Gly) and 3 (trimethylsilyl)-1-propane-sulfonic acid dissolved in D_2O .

- Load the data file into the experiment: cd('/vnmr/fidlib/Dosy') rt('Dbppste.fid')
- 2. Enter **ft** and adjust the phase of the first spectrum.
- 3. Set the cursor to the TSP singlet and enter **nl rl(0)**.
- 4. Change to a directory (cd) in which you have write permission.
- 5. Set the cursors 30 Hz either side of the TSP singlet, set **lb='n'** and **gf=0.75**, and enter the command

```
fiddle('satellites','TMS','writefid','temp')
```

to perform reference deconvolution on all the data, regularizing the lineshapes so that the peak heights in successive spectra accurately reflect the signal integrals.

6. Recall the FID created in the previous step and retransform it:

```
rt('temp') gf='n' lb='n' ft
```

- 7. The integral regions have already been set in the supplied parameters. Display the integral to see where the resets have been positioned.
- 8. Enter **fbc** to perform baseline correction.
- Set the threshold below the peaks of interest: (vs=500, th=3)
- 10. Enter **dosy**.
- 11. To zoom into the diffusion region of interest, enter undosy dosy(4,7).

The following examples describe how to process DOSY data in the command mode (Commands column) or by using the Tcl-Tk Process (2D-DOSY) or Process2 (3D-DOSY) panels (Buttons in the Process Window column). The middle column has instructions or comments about both types of operations.

DgcsteSL.fid

The sample is a mixture of adenosine mono-, di-, tri-phosphate (AMP, ADP, ATP) and K_2 HPO₄ in D₂O (pH=7). The data was acquired in a 3mm probe with direct ³¹P observe.

Commands	Comments, Instructions for Both	Buttons in the Process Window
	Recall the FID:	
	cd('/vnmr/fidlib/Dosy')	
	rt('DgcsteSL.fid')	
lb=2 wft	Fourier transform	Process
fbc	Do baseline correction	Baseline correct All spectra
dosy	Execute dosy.	Calculate full DOSY spectrum
	To have better diffusion resolution,	
	calculate a partial dosy spectrum:	
undosy		Recall original NMR data
dosy(6.1,7.1)		Calculate partial DOSY spectrum
	To display (and plot) the diffusion spectrum, join another experiment and execute sdp.	

Dbppsteinept.fid

The sample is a mixture of sucrose, methyl-alfa-D-glucopyranosid, 1,3,5,-O-methylidenemio-inosytol, and dioxane (as internal reference) in D_2O . The experiment was run using an AutoSwitchable gradient probe.

Command	Comments, Instructions for Both	Buttons in the Process Window
	Recall the FID: cd('/vnmr/fidlib/Dosy') rt('Dbppsteinept.fid')	
ft	Fourier transform	Process (unset lb and gf)
fbc	Do baseline correction	Baseline correct All spectra
lb=-0.4 gf=0.7	Set weighting functions. Expand the spectrum and put the two cursors around the most intense line dixonane \pm 15 Hz.	(Activate lb and gf)
fiddle	Execute fiddle.	FIDDLE (No TMS)
	Display full spectrum and set threshold.	
dosy	Execute dosy.	Calculate full DOSY spectrum
	To have better diffusion resolution, calculate a partial dosy spectrum:	
undosy		Recall original NMR data

Command	Comments, Instructions for Both	Buttons in the Process Window
dosy(2.0,5.0)		Calculate partial DOSY spectrum
	To display (and plot) the diffusion spectrum, join another experiment and execute sdp.	

Dgcstecosy.fid

The sample is a mixture of sucrose, methyl-alfa-D-glucopyranosid, and 1,3,5,-O-methylidene-mio-inosytol in D_2O . The experiment was run using an AutoSwitchable gradient probe.

Command	Comments, Instructions for Both	Buttons in the Process2 Window
	Recall the FID: cd('/vnmr/fidlib/Dosy') rt('Dgcstecosy.fid')	
wft2d(1)	Fourier transform.	Retransform the 1st 2D spectrum
	Signal regions for this file have already been saved. Recall 112d file: 112d 112d('read' 'Dgcstecosy.112d')	
112dmode='nnyn'	Check preset regions. Each cross peak of interest is boxed.	Labelling options: nnyn
dconi		Redisplay spectrum
dosy	Execute dosy.	Calculate DOSY spectrum
	When ready, COSY spectrum is displayed again with each cross peak labelled by its diffusion coefficient and its error.	
	Join another experiment and display diffusion projection: sdp	
	A set of signals appear:	
	4.1-4.8–1,3,5,-O-methylidene- mio-inosytol	
	3.6-3.9–methyl-alpha- D-glucopyranosid	
	2.8-3.1-sucrose	
	Other three lines between 3.2 and 3.6 D (10^{-10} m ² /sec) are overlapping diagonal peaks.	
	Rejoin DOSY experiment.	
112dmode='nnnn' dconi	Reset peak labels.	Labelling options: nnnn Redisplay spectrum
	To display inosytol spectrum:	

Command	Comments, Instructions for Both	Buttons in the Process2 Window
makeslice(4.1,4.8)		Low Lim.: 4.1, Up. Lim: 4.8
		Show 2D projection within
	To display glucopyranosid projection, recall original 2D.	
showoriginal		Show original 2D (first spectrum)
<pre>makeslice(3.6,3.9)</pre>	Display glucopyranosid.	Low. Lim.: 3.6, Up. Lim: 3.9
	To display the last projection, recall the original 2D.	
showoriginal		Show original 2D (first spectrum)
<pre>makeslice(2.8,3.1)</pre>	Display sucrose.	Low Lim.: 2.8, Up. Lim: 3.1
		Show 2D projection within

Note: By accident, this cosy spectrum was run with an unusual parameter setting (sw<>sw1). The setting was absolutely unintended and should not affect the DOSY processing. The operator (P.Sandor, Darmstadt) assures you that the sequence also gives proper results with adequate parametrization.

Dgcstehmqc.fid

The sample is a mixture of quinine, geraniol, and camphene (and TMS) in deuteromethanol. See *J. Magn. Reson.* (1998), **131:** 131-138.

Commands	Comments, Instructions for Both	Buttons in the Process2 Window
	Recall FID:	
	cd('/vnmr/fidlib/Dosy') rt('Dgcstehmqc.fid')	
wft2d('ptype',1)	Fourier transform.	Retransform the first 2D spectrum
	Signal regions for this file have already been saved. Recall ll2d file:	
	ll2d ll2d('read','Dgcstehmqc.ll 2d')	
ll2dmode='nnyn'	Check preset signal regions. Each cross peak of interest is boxed	Labelling options: nnyn
dconi		Redisplay spectrum
dosy	Execute dosy.	Calculate DOSY spectrum

Commands	Comments, Instructions for Both	Buttons in the Process2 Window
	When ready, HMQC spectrum is displayed again with each cross peak labelled by its diffusion coefficient and its error.	
	Join another experiment and display diffusion projection: sdp	
	A set of signals appears:	
	7.0-8.5 – quinine	
	10.0-11.6 – geraniol	
	14.0-15.4 – camphene	
	Other lines around $18 \text{ D} (10^{-10} \text{ m}^2/\text{sec})$ are methanol and TMS.	
	Rejoin DOSY experiment.	
ll2dmode='nnnn' dconi		Labelling options: nnnn Redisplay spectrum
	To display quinine spectrum:	
<pre>makeslice(7.0,8.5)</pre>		Low. Lim.: 7.0, Up. Lim: 8.5 Show 2D projection within
	To display another projection, you need to recall the original 2D.	
showoriginal	-	Show original 2D (first spectrum)
<pre>makeslice(10.0,11.6)</pre>	Display geraniol.	Low. Lim.: 10.0, Up. Lim: 11.6
		Show 2D projection within
	To display last projection, recall original 2D.	
showoriginal		Show original 2D (first spectrum)
makeslice(14.0,15.4)	Display camphene.	Low. Lim.: 14.0, Up. Lim: 15.4 Show 2D projection within

DOSY-Related Literature

Morris, K.F.; Johnson, C.S., Jr. "Resolution of Discrete and Continuos Molecular Size Distributions by Means of Diffusion-Ordered 2D NMR Spectroscopy," *J. Am. Chem. Soc.* (1993), **115**: 4291-4299.

Wider, G.; Dötsch, V.; Wütrich, K. "Self-Compensating Pulsed Magnetic-Field Gradients for Short Recovery Times," *J. Magn. Reson.* (1994), **108** (Series A): 255-258.

Barjat, H.; Morris, G.A.; Smart, S.; Swanson, A.G.; Williams, S.C.R. "High-Resolution Diffusion-Ordered 2D Spectroscopy (HR-DOSY) – A New Tool for the Analysis of Complex Mixtures," *J. Magn. Reson.* (1995), **108** (Series B): 170-172.

Wu, D.; Chen, A.; Johnson, C.S., Jr. "An Improved Diffusion-Ordered Spectroscopy Experiment Incorporating Bipolar-Gradient Pulses," *J. Magn. Reson.* (1995), 115 (Series A: 260-264.

Gozansky, E.K.; Gorenstein, D.G. "DOSY-NOESY: Diffusion-Ordered NOESY," J. Magn. Reson. (1996), **111**, (Series B): 94-96.

Wu, D.; Chen, A.; Johnson, C.S., Jr. "Three-Dimensional Diffusion-Ordered NMR Spectroscopy: The Homonuclear COSY-DOSY Experiment," *J. Magn. Reson.* (1996), **121 (Series A):** 88-91.

Wu, D.; Chen, A.; Johnson, C.S., Jr. "Heteronuclear-Detected Diffusion-Ordered NMR Spectroscopy through Coherence Transfer," *J. Magn. Reson.* (1996), **123 (Series A):** 215-218.

Jerschow, A.; Müller, N. "3D Diffusion-Ordered TOCSY for Slowly Diffusing Molecules," *J. Magn. Reson.* (1996), **123 (Series A):** 222-225.

Birlikaris, N.; Guittet, E. "A New Approach in the Use of Gradients for Size-Resolved 2D-NMR Experiments," *J. Am. Chem. Soc.* (1996), **118:** 13083-13084.

Jerschow, A.; Müller, N. "Suppression of Convection Artifacts in Stimulated Echo Diffusion Experiments. Double-Stimulated-Echo Experiments," *J. Magn. Reson.* (1997), **125:** 372-375.

Barjat, H.; Morris, G.A.; Swanson, A.G., "A Three-Dimensional DOSY-HMQC Experiment for the High-Resolution Analysis of Complex Mixtures," *J. Magn. Reson.* (1998) **131:** 131-138.

Pelta, M.D.; Barjat, H.; Morris, G.A.; Davis, A.L., Hammond, S.J. "Pulse Sequences for High Resolution Diffusion-Ordered Spectroscopy (HR-DOSY)," *Magn. Reson. Chem.* (1998), **36:** 706.

Tillett, M.L.; Lian, L.Y.; Norwood, T.J. "Practical Aspects of the Measurement of the Diffusion of Proteins in Aqueous Solution.," *J. Magn. Reson.* (1998), **133**: 379-384.

Gounarides, J.S.; Chen, A.; Shapiro, M.J. "Nuclear Magnetic Resonance Chromatography: Applications of Pulse Field Gradient Diffusion NMR to Mixture Analysis and Ligand-Receptor Interactions," *Journal of Chromatography B* (1999), **725**: 79-90.

DOSY Review Papers

Morris, G.A.; Barjat, H., "High Resolution Diffusion Ordered Spectroscopy," *Methods for Structure Elucidation by High Resolution NMR*, ed. K. Kövér, Gy. Batta, Cs. Szántay, Jr. (Amsterdam: 1997), pp. 209-226.

Morris, G.A.; Barjat, H.; Horne, T.J. "Reference Deconvolution Methods (FIDDLE)," *Progress in Nuclear Magnetic Resonance Spectroscopy* (1997), **31:** 197-257.

Johnson C.S. Jr., "Diffusion-Ordered Nuclear Magnetic Resonance Spectroscopy: Principles and Applications," *Progress in Nuclear Magnetic Resonance Spectroscopy* (1999), **34:** 203-256.

Filter Diagonalization Method

Filter Diagonalization Method (FDM) is a non-Fourier data processing method that extracts spectral parameters (peak positions, line widths, amplitudes, and phases) of Lorentzian lines directly from the time-domain signal by fitting FID data to a sum of damped complex sinusoids. The spectral parameters (saved in curexp/datdir/fdml.parm) are also called "line list" and are used to construct an "ersatz" spectrum of the NMR data.

FDM is slower than Fast Fourier Transform, but it offers better resolution in the case of truncated signals and the option of processing only a selected spectrum region. FDM has the potential to work well with corrupted data, and the potential to produce a line list with each line represents a true NMR peak.

FDM reads input parameters from a file created by the fdml macro, using default (optimal) values. You can change any of the parameters from the command line. Table 15 lists fdml parameters. If the spectrum is not referenced with rl, the reference rfl is also read from curpar. The section "Changing Local Variables," page 68 describes how you can override the default setting.

In most cases, you only need to decide the number of data points to be used and the spectrum window to be processed. By default, half of the FID data or 3000 data points, whichever is smaller, is used.

The window to be processed is determined by VNMR parameters sp and wp. If the data is already processed (using FT or FDM) and displayed, you can process the displayed region again with FDM by typing fdm1. If the data is not processed and displayed, sp and wp are read from the curpar file, but sp and wp might not be what you want.

Using FDM

The following steps describe how to do normal activities such as phasing, zooming in, zooming out, and processing a spectrum window with the fdml macro.

Parameter	Description	
cheat	No cheat if cheat=1, lines are narrower if cheat<1.	
cheatmore	No cheatmore if cheatmore=0.	
error	Error threshold for throwing away poles.	
fdm	1 for fdm, −1 for dft.	
Gamm*	Smoothing width (line broadening).	
Gcut	Maximum width for a pole.	
idat	-4 for ASCII complex FID file, -5 for VNMR FID file.	
kcoef	kcoef>0; use "complicated" dk (k)1 always preferred.	
Nb*	Number of basis function in a single window.	
Nbc*	Number of coarse basis vectors.	
Nsig*	Number of points to use, 3000 is ok.	
Nskip*	Number of points to skip.	
rho	rho=1 is optimal.	
SSW	A test parameter.	
t0	Delay of the first point.	
theta	Overall phase of FID (rp in radians).	
wmin	Minimum spectrum frequency in hertz.	
wmax	Maximum spectrum frequency in hertz.	
* Global; see "Global Parameters," page 67 for more information.		

 Table 15.
 fdm1 Parameters

- 1. Display the FID data and use the right mouse button to select the data points to be used by FDM.
- 2. Process the data with ft (it uses all FID points), then display and reference the spectrum.
- 3. Place the cursor on a region of interest, zoom in on it, then type fdml or select the **Process2** panel and click on **1D FDM** to process the data. If you select but do not zoom in on a region, the whole spectrum in display is processed.

A new menu appears with **Stop FDM** and **Display** buttons. The calculation might take a few seconds to a few minutes depending on the number of data points used and the size of spectrum window to be processed. To abort the process, click on **Stop FDM**. To check if the process is finished, click on **Display**. If the process is finished, display the spectrum.

Reprocessing a Spectrum

The **1D FDM** button is displayed on the **Process2** panel. Use this button to reprocess a spectrum.

Changing Parameters

Relevant fdml global parameters are displayed on the **Process2** panel with current values. You can change these parameters. The value of a global parameter is saved to curpar and it remains the same until you change it from the parameter panel or make a new assignment using the command line. You can also change the parameters from the fdml command line as described in the section "Changing Local Variables," page 68.

Global Parameters

The following FDM parameters are global.

- Nsig is the number of FID points to use. You initialize it with the right mouse button position (crf+deltaf) *sw. If Nsig=0, half of the FID data points or 3000, whichever is smaller, is used. Nsig can be changed from the parameter panel, the command line Nsig=nnnn, the right mouse button (when the FID is displayed), or the command line fdml('Nsig',nnnn). In general, the more peaks you have, the more data points it takes to fit the spectrum. To check the reliability of the FDM method, change Nsig a few times and reprocess the data to see if you get the same result.
- Nskip is the number of data points to skip at the beginning of a FID. By default, zero points are skipped. In some cases, you can improve baseline by skipping the first one or two points.
- Nb is the number of basis functions (poles) used to fit each of the windows in an FDM calculation. The default is 10. FDM breaks down the specified window into smaller windows. In general, bigger Nb gives better results, especially better baseline. Sensible values for Nb are between 10 and 50.
- Nbc is the number of additional poles (coarse basis functions) to be used. The default is zero, but setting Nbc to an integer larger than zero (typically 4-10) might improve the baseline.
- Gamm is the smoothing width (line broadening). The default is 0.2*sw/Nsig, which is about a tenth of the FT resolution. Typical values are 0.1 to 1.0.

Using bigger Nsig, Nb, Nbc, or a spectral window significantly slows down the calculation.

Changing Local Variables

FMD parameters that are not commonly used are set as fdm1 local variables. These parameters are listed with global parameters in Table 15. You can change local variables only from the fdm1 command line. Parameter values are lost after the completion of the macro. To use a value again, you must reenter it; otherwise, fdm1 sets the value to the default. To change more than one local variable, enter the variables from the same command line.

You can change any of the FDM parameters from the fdml command line and you can change both global and local variables. Values entered from the fdml command line override the default, the change from the **Process2** panel, and the value that you select with the cursor. Enter command line arguments by giving the parameter name in single quotation marks and a value separated by a comma, for example:

```
fdml('cheat',0.8)
fdml('Nsig',3000)
fdml('Nsig',3050)
fdml('Nb',20)
fdml('Nbc',10,'Nb',20)
fdml('Nsig',3000,'Nb',20,'Gamm',0.5)
fdml('wmin',-1600,'wmax',1600)
fdml('wmin',-1600,'wmax',0)
```

cheat is a factor multiplied to the line width. There is no cheat when cheat=1; lines are narrower when cheat<1.

wmin is the minimum spectrum frequency in Hz. The default is sp+rfl-sw/2. wmin is the upper field.

wmax is the maximum spectrum frequency in Hz. The default is wmin+wp. wmax is the lower field.

The center of the full spectrum is zero.

Seeing Parameter Values

Parameters are set to their default values. Normally, you do not need to change these parameters or you might change some of the global parameters. You cannot inquire values of local fdml parameters in the same way that you inquire global VNMR parameters with echo or ?. To see the values of all parameters used, look in the fdml.inparm file created by the fdml macro in the datdir directory of the current experiment. Figure 13 shows the format of the fdml.inparm file; the number of spaces and tabs is arbitrary.

References

J. Chen and V. A. Mandelshtam, J. Chem. Phys. (2000) 112: 4429-4437.

V. A. Mandelshtam, J. Magn. Reson. (2000) 144: 343-356.

A. A. De Angelis, H. Hu, V. A. Mandelshtam and A.J. Shaka, *J. Magn. Reson.* (2000) **144:** 357-366.

fid_filename	idat	
t0	theta	
fdm		
par		
fn_SplD spec	ctype axis	
wmin	wmax	
Nsig		
Nskip		
rho	Nb	
error		
Npower	Gamm Go	cut
cheat	cheatmore	
Nbc kcoef		
SSW SW		
fidmt specfmt		
i_fid		

Figure 13. fdm1.inparm File

Chapter 1. Advanced 1D NMR

Chapter 2. 1D Experiments

Sections in this chapter:

- 2.1 "APT—Attached Proton Test," page 72
- 2.2 "BINOM—Binomial Water Suppression," page 73
- 2.3 "CPMGT2—Carr-Purcell Meiboom-Gill T₂ Measurement," page 74
- 2.4 "CYCLENOE—Cycled NOE Difference Experiment," page 75
- 2.5 "D2PUL—Standard Two-Pulse Using Decoupler as Transmitter," page 76
- 2.6 "DEPT—Distortionless Enhancement by Polarization Transfer," page 77
- 2.7 "INEPT—Insensitive Nuclei Enhanced by Polarization Transfer," page 80
- 2.8 "JUMPRET—Jump-and-Return Water Suppression," page 81
- 2.9 "NOEDIF—NOE Difference Experiment," page 82
- 2.10 "PRESAT—1D Water Suppression," page 85
- 2.11 "S2PUL—Standard Two-Pulse Sequence," page 85
- 2.12 "S2PULR—Standard Two-Pulse in Reverse Configuration," page 86

This chapter describes a number of common 1D pulse sequences for "everyday" use. Each pulse sequence has a macro, with the same name as the pulse sequence, to help you set up the experiment. Many macros are written with the assumption that you have done a "normal" 1D experiment on the sample first. For example, after you have obtained a carbon spectrum, you can type apt to set up an APT experiment.

Each macro retrieves an associated file with parameters such as pw90, tpwr, dmf, etc. from a central location like /vnmr/probe or \$vnmruser/probe. Other specific parameters, such as mult for DEPT and d2 for APT, come from your user's parlib directory or, if the file is not found there, from the system directory /vnmr/parlib.

It is important that these parameters be correct. The first time the macro for an experiment is entered (by typing apt, for example), the system retrieves the default parameters and values.

To change any of these values (for example, the default d2 time in apt is 7 ms and you want the default to be 3 ms), make the appropriate change in the displayed parameters, and then save the modified parameters either in your user's parlib or in the system /vnmr/parlib. Notice that because files in the directory /vnmr/parlib are available to all users, only the system administrator vnmrl has permission to save the files in this directory.

To view complete listings of each pulse sequence, print out or look at the contents of the files in the directory /vnmr/psglib on your system disk. You can also enter dps to view a graphical representation.

2.1 APT—Attached Proton Test

The apt<(solvent)> macro converts a 13 C parameter set to an APT experiment, where solvent is the name of the solvent to be used. If solvent is not supplied, solvent either defaults to CDCl₃ or, if in automation mode, solvent is read from the sampleinfo file.

Use half as many transients as the normal carbon 13, assuring multiples of 4. This must be done before changing other parameters because as soon as any parameter is changed, the number of completed transients is zero (ct=0). Figure 14 is a diagram.




Figure 14. APT Pulse Sequence

Applicability

APT is available on all systems.

Parameters

d2 is the tau delay (in sec): 3 ms nulls aromatic CH; 4 ms nulls all else except J > 125 Hz and quaternaries; 6 to 8 ms gives quaternaries, CH₂'s up methyls; CH's down.

p1 is a 180° pulse for observe nucleus (composite pulses used).

pw is a normal observe pulse. It need not be 90°.

d3 is a short delay (typically 1 ms); if d3=0, the second 180° pulse is omitted and pw should be set to 90° or greater (θ =180– θ).

dm='yny' for an APT spectrum; dm='yyy' for a simple spin-echo spectrum.

rof1 is set to 0 to turn receiver off during the entire pulse sequence.

rof2 is the receiver dead time after last pulse in sequence.

Technique

To set up the experiment, enter apt. The apt macro sets up the experiment either by modifying an existing ¹³C parameter set or by using default ¹³C parameters. apt also sets the τ delay d2 to 7 ms, which gives CH's, CH₃'s down and CH₂'s up.

References

Rabenstein; Nakashima, T. Anal. Chem. 1979, 51, 14651A.
Lecocq, C.; Lallemand, J. J. Chem. Soc. Chem Comm. 1981, 150.
Patt, S.; Shoolery, J. J. Magn. Reson. 1982, 46, 535.

Related Commands and Macros

The aptaph macro automatically phases APT spectra.

2.2 BINOM—Binomial Water Suppression

The binom macro sets up parameters for the binomial water suppression sequence. Figure 15 is a diagram of the BINOM sequence with seq=1510.




Figure 15. BINOM Pulse Sequence

Applicability

BINOM is not supplied on MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

p1, p2, p3 are the pulses of the first half of the pulse sequence; all pulse sequences are symmetric (p3 is used only for 1-5-10-10-5-1). If p2=0 or p3=0, those pulses are derived from p1. If p1=0, p1 is derived from pw, which is assumed to be a 90° pulse.

offset is the frequency, in Hz, from the carrier at which maximum excitation occurs with the binomial sequences. If offset > 0, suppression is on-resonance. If offset < 0, suppression is off-resonance (a null at the offset specified in Hz from the transmitter offset position) and d2 is calculated by the program. If offset=0, a value of d2 entered gives off-resonance suppression at the corresponding offset from the transmitter position.

d2 is directly entered only if offset=0; otherwise, it is calculated.

seq is 11, 121, 1331, 146 (gives 1-4-6-4-1), or 1510 (gives 1-5-10-10-5-1).

rof2 is the receiver gating time, in µs, after last pulse in sequence.

rof1 is the receiver gating time, in μ s, before and after all other pulses. If rof1=0, the gate receiver is off during the entire sequence.

Reference

Hore, P. J. J. Magn. Reson. 1983, 55, 83-300.

Starcuk; Sklenar J. Magn. Reson. 1985, 61, 567-570.

2.3 CPMGT2—Carr-Purcell Meiboom-Gill T₂ Measurement

The cpmgt2 macro modifies a parameter set to perform a Carr-Purcell Meiboom-Gill T_2 measurement. Figure 16 is a diagram of the CPMGT2 sequence.




Figure 16. CPMGT2 Pulse Sequence

Applicability

CPMGT2 is available on all systems.

Parameters

pw is set for 90° pulse.

pl is set for 180° pulse.

d1 is set for 1 to 3 times the T_1 value.

d2 is the spin-echo cycle time; typical value is 1 to 10 ms.

bt ("big tau") is $4n \times d2$, the total time for T_2 relaxation. Typically, bt is arrayed for a T_2 experiment. If bt is not a multiple of $4 \times d2$, it is rounded automatically so that it is.

T₂ Measurement

1. Run a spectrum of the material. Then enter the macro cpmgt2.

This macro sets up the T_2 measurement parameters.

2. Create an array of bt. Values are entered in seconds. This array differs depending on the T_2 's to be measured. An example is bt=.020,.040,.080,.1,.120,.140

The computer determines the actual values used (selected as close as possible to the values you entered but still satisfying the equality $bt=4n \times d2$). In the sequence diagram shown in <Red>Figure 12, the part in parentheses is repeated *n* times. The length of bt is limited, usually to about 0.5 seconds, but the upper limit varies depending on a variety of factors.

Acquisition and Processing

- 1. Enter ga to start data collection. After data collection completes, enter the commands below to process the data.
- 2. Enter ds(1) to display the first spectrum,
- 3. Click the mouse on the th button to obtain a threshold line, then place the threshold below the tops of the peaks of interest.

- 4. Enter nll to find the peak frequency.
- 5. Enter fp to find the top of the peak.
- 6. Enter t 2 to calculate and print T_2 's of the peak.
- 7. If the exponential curve is not displayed, enter expl.
- 8. To plot the exponential curves, enter pexpl.
- 9. To plot the data, enter dssh pl('all') pap page.

2.4 CYCLENOE—Cycled NOE Difference Experiment

The cyclenoe macro sets up parameters for a cycled NOE difference experiment.

Applicability

CYCLENOE requires that the observe channel be equipped with direct synthesis rf and a linear amplifier. This experiment does not apply to *MERCURY-VX*, *MERCURY*, or *GEMINI 2000* systems; see the NOEDIF experiment instead.

Parameters

pw is a 90° excitation pulse (at power tpwr).

intsub='y' sets internal subtraction of data acquired by on-resonance and offresonance selective excitation on alternate scans. intsub='n' makes data acquired by on-resonance and off-resonance selective excitation to be stored separately; only 1D with satfrq is collected.

satfrq is the frequency of selective saturation (on-resonance).

control is the off-resonance selective saturation frequency (an inactive parameter if intsub='n'). If control is inactive, off-resonance spectra is not collected.

cycle='y' does on-resonance saturation using frequency cycling around the frequency satfrq given by spacing and pattern; cycle='n' does off-resonance saturation at control.

spacing is the spacing, in Hz, of the multiplet.

pattern is the pattern type (1 for singlet, 2 for doublet, etc.). Fit the desired pattern to some value, even if some frequencies do not fall on NMR lines

tau is the time spent on a single irradiation point during cycling.

satpwr is the power of selective irradiation (typical values are at lowest power).

sattime is the total length of irradiation at frequency satfrq.

mix is the mixing time.

sspul='y' does trim(x)-trim(y) before d1.

nt is a multiple of 16 if intsub= 'n', nt is a multiple of 32 if intsub= 'y'.

Technique

CYCLENOE does alternate scan subtraction of two FIDs in which the saturation frequency is moved on-resonance and off-resonance (for intsub='n'). Separate data tables are stored for intsub='n'. Power may be reduced from ordinary NOE experiments because the irradiation can be cycled (cycle='y') through the lines of a multiplet.

- 1. One way to set up satfrq is by entering dn='H1' and using sd the same as in homodecoupling. Then enter satfrq=dof and dn='C13'.
- 2. Adjust proper satpwr by setting nt=1 and arraying satpwr from 3 to -16 on UNITYINOVA and UNITY*plus* systems, and from 3 to 0 on UNITY and VXR-S systems. Enter proper values for pattern and spacing. tau is typically a few hundred milliseconds. sattime is usually several seconds.
- 3. Acquire the data and select the power necessary for 50 to 75% saturation. Then set nt to a large number (several hundred).
- 4. Set the control frequency as near as possible to satfrq to make the control and satfrq conditions as close as possible. The control frequency should be "in the noise," not on top of a multiplet. Several protons may be done simultaneously for one control frequency.
- 5. Enter appropriate frequencies for patterns and spacings.
- 6. Set array='(satfrq,pattern,spacing)'. This performs the proper number of experiments.
- 7. Run nonspin and temperature-regulated. A large number of transients result in better subtraction.

2.5 D2PUL—Standard Two-Pulse Using Decoupler as Transmitter

The d2pul macro sets up parameters for a two-pulse sequence with the decoupler configured as a transmitter. Figure 17 is a diagram of the D2PUL sequence.




Figure 17. D2PUL Pulse Sequence

Applicability

D2PUL is available on all UNITY-series systems. On *MERCURY-VX* and *MERCURY* systems, the high-band transmitter is used as observe and decoupler—use S2PUL. On *MERCURY-VX*

and *MERCURY* and *GEMINI 2000* systems, if tn is 'H1' and dn is not 'H1', the software now automatically uses the decoupler as the observe channel and the broadband channel as the decoupler channel. This channel swapping also makes the pulse sequence d2pul macro obsolete on the *GEMINI 2000*. If you want to run the equivalent of d2pul, set tn='C13' and dn='H1', and then run s2pul.

Parameters

dof (decoupler offset) should equal tof (transmitter offset) for proper signal detection.

homo must be set to 'n'.

tpwr, instead of dpwr, controls the decoupler power level during pulses on UNITY*INOVA*, UNITY*plus*, and UNITY systems equipped with a linear amplifier on the decoupler rf channel. dpwr controls the decoupler power level at all other times during a pulse sequence for this type of system configuration.

For systems with class C decoupler amplifiers, dhp is set to 255 by the pulse sequence statement declvlon during all pulses and is reset to its parameter value for all other times during the pulse sequence.

Technique

Acquire a normal spectrum, then enter d2pul to convert the parameter set.

2.6 DEPT—Distortionless Enhancement by Polarization Transfer

The dept macro sets up parameters for a DEPT pulse sequence, an improved version of the INEPT experiment. Figure 18 is a diagram of the DEPT sequence.




Figure 18. DEPT Pulse Sequence

Applicability

DEPT is available on all systems.

Parameters

pw is the 90° pulse (in μ s) on the X-nucleus (usually ¹³C) at power equal to tpwr.

pp is the proton 90° pulse width supplied from decoupler (refer to the decoupler 90° pulse width test using polarization transfer in the manual *Getting Started*). pp is set at the power level pplvl if linear amplifiers are present, or at full power if class C amplifiers are used. Optimum spectral editing requires a carefully calibrated value of pp. *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* broadband systems use linear amplifiers and the parameter pplvl must be set correctly. *GEMINI 2000* ¹H/¹³C spectrometers do not have linear amplifiers and pplvl has no meaning.

dmf sets the modulation frequency for broadband decoupling of protons at power equal to dhp (if class C amplifier) or dpwr (if linear amplifiers). On *GEMINI 2000*, UNITY, and VXR-S systems, the modulation mode is WALTZ. On UNITY INOVA, MERCURY-VX, MERCURY, and UNITY plus systems, a variety of modulation modes can be used. dmf sets the frequency for all modes.

j is an average X-H (usually C-H) coupling constant in Hz (tau is an internal parameter set to 1/2J).

mult is an arrayed parameter (0.5, 1, 1, 1.5) that leads to a value for the θ pulse of mult*pp: mult=0.5 (θ = 45°) gives approximately equal excitation of all protonated carbons; mult=1.0 (θ = 90°) excites CH's only (or mainly); mult=1.5 (θ = 135°) gives CH's, CH₃'s up, CH₂'s down.

d1 is a relaxation delay (1 to 3 times the value of t_1) for the protons connected to the NMR active nucleus, for example, ¹³C or ¹⁵N.

at is the acquisition time (t_2 period).

nt (number of transients) should be a multiple of 4. A multiple of 16 is suggested.

dm (decoupler modulation) is set to 'nny'.

satdly is an optional saturation delay. If satdly is greater than 0, PW(90)-HS-SATDLY ¹³C presaturation occurs immediately following the dl delay. If satdly=0, no PW(90)-HS-SATDLY ¹³C presaturation occurs. The recommended value is satdly=0.00. (satdly is not available on *MERCURY-VX*, *MERCURY*, or *GEMINI 2000.*)

hs='nn' for no homospoil is recommended. If homospoil is used, hst (homospoil time) is typically 10 ms (hs must be 'nn' on a *GEMINI 2000*).

Technique

To set up the experiment, enter dept. If a ${}^{13}C$ parameter set is present in the experiment, the dept macro modifies existing parameters; otherwise, it uses the default ${}^{13}C$ parameters. dept also sets mult=0.5,1,1,1.5.

Specific parameters for dept not found in the ¹³C parameter set are retrieved from the user's parlib directory, or if not found there, from the /vnmr/parlib system directory. In either case, it is important that the parameters in these directories be correct.

To initiate acquisition, enter au if you wish automatic processing, spectral editing, and plotting. For auto processing, the parameter wexp should be set to wexp='deptp'. The edited spectrum is constructed as follows:

1. ALL: 0.23*(s1 + s2) + 2*s1 2. CH: s2 + s3 3. CH2: s1 - s4 4. CH3: -0.77*(s1 + s3) + s1 + s4

where s1 represents spectrum 1, s2 represents spectrum 2, etc. Note that quaternary carbons are suppressed in the DEPT spectrum.

The dept experiment may be set up and run manually or by using the menu system. The following description uses the standard VNMR system menus and commands. Alternatively, the *GLIDE* system can be used to run the spectrum in a completely automated fashion.

- 1. A good sample to use the first time you run the experiment is 30% menthol in CDCl₃. Set up a carbon experiment by clicking on the **setup** menu button followed by the **C,CDCl3** button.
- 2. Set nt=4 and enter dept. This sets up the experiment.
- 3. Enter **ga** to acquire the data. Four spectra are acquired, which you can edit.
- 4. After data acquisition, display the first spectrum by entering **ds(1)**.
- 5. Click on the **th** button to obtain a threshold line. Select a threshold level below the top of all peaks.
- Enter adept dssa to analyze the spectra.
- 7. Enter **pldept** to plot the DEPT spectra.

Potential Problems

The most common failure of dept is poor subtraction in the edited spectra. Poor subtraction is usually caused by improper calibration of the decoupler 90° pulse pp. If the dept does not work, check the ¹³C 90° pulse width, and decoupler 90° pulse width calibrations. Other causes of poor cancellation are lock saturation, leading to unstable lock, improper vibration isolation of the system, or temperature change during the experiment.

Reference

Doddrell, D.; Pegg, D.; Bendall, M. J. Magn. Reson. 1982, 48, 323.

Related Commands and Macros

The following macros assist in processing, analyzing, and plotting DEPT data:

- The adept command automatically analyzes a set of four DEPT spectra and edits the spectra so the spectra is arrayed in the following order: CH₃ carbons only, CH₂ carbons, CH carbons only, and all protonated carbons. adept produces a text file dept.out in the current experiment directory that contains the result of the analysis. Refer to the description of adept in the VNMR Command and Parameter Reference for the arguments available with adept and other information about the command.
- The autodept macro process DEPT spectra, plots the unedited spectra, edits the spectra, plots the edited spectra, and prints out editing information.
- The deptgl macro sets up parameters for a DEPTGL sequence (not distributed with *MERCURY-VX*, *MERCURY*, or *GEMINI 2000* systems) for spectral editing and polarization transfer experiments. Refer to the description of deptgl in the *VNMR Command and Parameter Reference* for the parameters used with deptgl and a reference to the literature.

- The deptproc macro weights, Fourier transforms, and phases each spectrum in a DEPT data set.
- The padept command performs the adept analysis and plots the resulting spectra with a scale and the assigned line listing. Refer to the description of padept in the *VNMR Command and Parameter Reference* for the keyword arguments used with padept and other information.
- The pldept macro plots out DEPT data, edited or not edited.

2.7 INEPT—Insensitive Nuclei Enhanced by Polarization Transfer

The inept macro sets up parameters for an INEPT pulse sequence. Figure 19 is a diagram of an INEPT sequence with mult=3.




Figure 19. INEPT Pulse Sequence

Applicability

INEPT is available on all systems.

Parameters

j is a coupling constant (0 allows entry of d2, d3).

pw is a 90° pulse on the observe nucleus.

pp is a 90° pulse on protons using the decoupler.

mult selects multiplicity:

mult=0 gives a normal experiment (d1-pw-at sequence);

- mult=2 selects ch's (doublets);
- mult=3 gives ch2's down, ch's and ch3's up; and
- mult=4 enhances all protonated carbons.

dm= 'nny' gives decoupled spectrum; dm= 'nnn' gives coupled spectrum.

```
dmm is set to 'ccw' or 'ccf'.
```

focus='y' gives refocusing for coupled spectrum (decoupled spectra are always refocused).

normal='y' gives normal multiplets in coupled spectra.

If the decoupler channel uses a linear amplifier, pplvl is the power level for the proton decoupler pulse and dpwr is the power level for broadband proton decoupling. If the decoupler channel uses a class C amplifier, maximum power is used for the proton decoupler pulse and dhp specifies the power level for broadband proton decoupling.

dl is a relaxation delay (1 to 3 times the value of t_1) for the protons connected to the NMR active nucleus, for example, ¹³C or ¹⁵N.

Technique

To set up the experiment, enter inept. If a parameter set is present in the experiment, the inept macro modifies existing parameters; otherwise, it uses the default parameters.

Specific parameters for inept not found in the spectral parameter set are retrieved from the user's parlib directory, or if not found there, from the /vnmr/parlib directory. In either case, it is important that the parameters in these directories be correct.

Reference

Morris, G.A.; Freeman, R. J. Am. Chem. Soc. 1979, 101, 760.

2.8 JUMPRET—Jump-and-Return Water Suppression

The jumpret macro sets up parameters for the JUMPRET pulse sequence.

Applicability

JUMPRET is available on all systems except MERCURY-VX, MERCURY, and GEMINI 2000.

Parameters

pw is the 90° pulse width at power tpwr.

p1 is a correction to second pulse width, typically 0 to 1 µs. Try increments of 0.025 µs.

d2 is set for desired excitation maximum (typically 100 to 200 μ s).

lp should be able to phase spectrum with small lp (use calfa).

```
lfs, zfs: try solvent suppression digital filtering:
for wft('zfs'), try ssfilter=300, ssntaps=11, and ssorder=7.
for wft('lfs'), try ssfilter=100, ssntaps=200, and ssorder='n'.
```

Reference

Plateau, P.; Gueron, M. J. Amer. Chem. Soc. 1982, 104, 7310.

2.9 NOEDIF—NOE Difference Experiment

The noedif macro converts a ¹H parameter set to perform the NOE (Nuclear Overhauser Enhancement) difference experiment. Figure 20 is a diagram of the sequence. The NOEDIF experiment performs subtraction directly in the computer memory, avoiding the necessity for an add–subtract step.




Figure 20. NOEDIF Pulse Sequence

Applicability

NOEDIF is available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* only. For UNITY-Series systems, see CYCLENOE.

Parameters

d2 is set to 0 seconds.

f0 to f5 are individual decoupling frequencies through which to cycle (set unused frequencies to 0). f6, f7, f8, and f9 are also defined and can also be used.

You must start entering irradiation frequencies at f0 through f5. If f0 is left empty (0.0), then the entire irradiation period is skipped (as if dl=0).

tau is the time per individual decoupling (number of cycles equals $d1/(tau \times n)$, where n is the number of different frequencies).

dofoff is the position of decoupler during d1 of control experiment.

dof is the position of decoupler during d2 and acquisition.

ctrl is set to 'y' for control experiment; set to 'n' for decoupling cycling.

Phase Cycling

If ctrl='y': v1 = +x -x+y -y +x -x +y -v ... oph = +x -x +y -y +x -x+y If ctrl='n': v1 = +x +x+y +y -x -x -v -v ... oph = -x + x-y +y +x -x+v

Procedure

The following procedure performs the NOE difference experiment on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* spectrometers. The order of the first three steps is critical.

- Run a proton spectrum in the usual way, store the spectrum in exp1, then enter mf(1,2) to move the FID to exp2. Enter wft.
- 2. Enter **noedif**.

The noedif macro converts the pulse sequence from S2PUL to NOEDIF and creates the parameters to set up the experiment. The screen displays instructions.

- 3. Enter **ds**. Find the peak or multiplet to be irradiated in the NOE difference experiment and expand the spectral region around it so that it occupies about 20 percent of the width of the screen. Display a cursor and select the frequencies for cycling during irradiation as follows:
 - If the line is a single line, multiple irradiation points may be chosen by setting the cursor at the point halfway up the line on the left side, then enter sd followed by f0=dof. Reset the cursor at the peak of the line, then enter sd followed by f2=dof. Next, set the cursor at the point halfway up the line on the right side, enter sd followed by f3=dof. The decoupler now cycles over these three frequencies, remaining at each frequency for a time equal to the parameter tau seconds. Enter tau=0.1 unless it is already set to this value.
 - If the proton to be irradiated is a doublet, f1 through f4 are normally used, leaving f5 at 0. Set up the frequencies f1 to f4 as described above, setting the cursor in sequence on the side of the first line, then on the second, then on the other side of the first line, and finally on the other side of the second.
 - Triplets and higher multiplets are usually cycled over the frequencies of each of the lines in the multiplet. Set the frequencies fl to f# (where # is the number of lines in the multiplet) as described above, putting the cursor on the peak of each line.
- 4. Enter **ss=16** to set the steady-state parameter **ss** to perform the sequence 16 times before acquiring data, allowing the magnetization of the nuclei to reach equilibrium state.
- 5. Enter **d1=4** to set the delay d1 to 4 seconds. In samples with very long T_1 values, a longer delay is necessary to achieve equilibrium. Conversely, some samples with shorter T_1 's may give good results with a shorter delay.

The total number of steps equals d1/tau. Do not exceed 60 steps. If the number of steps exceeds 60, the acquisition parameters may be loaded into the acquisition computer by the host computer when the go command is given, but the acquisition computer may not be able to hold so many instructions and will not start to acquire data.

If the values of d1 and tau create more than 60 steps, increase the tau delay. For example, a d1 of 10 seconds with a tau=.1 produces 100 steps, which is too many. Increasing tau to 0.2 produces 50 steps, an acceptable number.

An integral number of frequency cycles (f1 through f#) is performed during d1. The computer changes tau if necessary to accomplish this.

- 6. Set dm='ynn' to turn on the decoupler during the delay period d1.
- 7. Set **dofoff** to determines where the decoupler is set to irradiate during alternate pulses when it does not irradiate the multiplet.

Chapter 2. 1D Experiments

The default value of doff is -10000 Hz, but it can be set anywhere in the spectrum that is not close to a line, even within the spectral range of the observed lines. The exact position should not matter, but more than one setting may be tried to determine experimentally whether it does or not.

8. Set ctrl='n' to cycle the decoupler during the NOE difference experiment. Alternate FIDs are subtracted, beginning with the second pulse for which the decoupler is set to the doff position. This results in positive peaks for positive NOEs and a large negative signal for the irradiated peak or multiplet.

If ctrl='y', the decoupler is not cycled and alternate FIDs are not subtracted, so the resulting spectrum is the normal proton spectrum, obtained with the same pulse width and delay time as the NOE experiment. This gives the control spectrum to compare with the NOE difference spectrum. Both spectra can be obtained in one experiment by entering ctrl='y', 'n' to array the value of ctrl.

- Set dpwr for *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* broadband systems, or dlp for ¹H/¹³C *GEMINI 2000*, to give a value of decoupling field strength at the sample of about 2 to 3 Hz. This requires that the homonuclear decoupler is calibrated.
 - For *MERCURY-VX* and *MERCURY*, dpwr controls decoupler power in 1-dB steps, from 0 to 63.
 - For *GEMINI 2000* broadband, dpwr controls the homonuclear decoupler field strength, set in 0.5 dB steps from 0 to 50. If equipped with PIN-diode attenuators, additional fine control can be used through dlp on BB systems.
 - For *GEMINI 2000* ¹H/¹³C, dlp controls the proton homodecoupler power level (if present) from 0 to 2047 in arbitrary units, where 2047 is full power.

The actual field strength can be calibrated exactly, but extreme accuracy is usually not necessary. Typical values vary widely from system to system:

- On MERCURY-VX and MERCURY, use dpwr set to 10.
- On *GEMINI 2000* broadband systems, start at dpwr=5 and increase or decrease as necessary.
- On *GEMINI 2000* ¹H/¹³C systems, start at dlp=1650 and increase in units of 50 until the correct power is obtained.

A simple method to calibrate the decoupler for NOEDIF is to irradiate the peak of interest with a small amount of power and partially saturate the peak. Compare the spectra with and without irradiation (set dm='ynn', 'nnn'). When the peak intensity has decreased by 50 to 75%, the power is optimum for NOE difference.

- 10. Set **bs='n'** to disable block size storage.
- 11. Turn off the spinner and control sample temperature.
- 12. Run the NOE difference experiment for 512 pulses (if time permits). In general, the number of repetitions of the pulse sequence should be large enough to build up good signal-to-noise, since the NOE difference peaks are only a few percent of the normal spectral intensity, and should also be large enough to statistically average unavoidable perturbations of the field or frequency. A value of nt between 256 and 1024 should be suitable for most work unless very small samples are studied.

Because the NOE difference experiment measures small differences between large signals, stability is very important. Be sure the spectrometer is protected from perturbing effects such as building vibrations, magnetic noise in the immediate vicinity (like moving iron objects and scanning magnetic fields), or strong air currents.

13. One possible cause of poor cancellation in this experiment is the noise in the lock channel. Set lock power to just below (3 dB) saturation and turn lock gain down as far as possible. The system remains locked at lock levels of 20%.

The optimum combination of lock gain and lock power must be determined for each sample. A series of experiments with nt set to 16 or 32 should allow the best combination to be found without investing too much time. Then the value of nt should be increased and the data collected.

Reference

Kinns, M.; Sanders, J. K. M. J. Magn. Reson., 1984, 56, 518.

2.10 PRESAT—1D Water Suppression

The presat macro sets up a standard two-pulse sequence with optional composite observe pulse: {trim(x)trim(y)}..dl..satdly..pl..d2..pw..at.

Applicability

PRESAT is available on all systems with a linear amplifier on the observe channel. It is not available on *GEMINI 2000* systems.

Parameters

satmode set to 'y' gives observe transmitter saturation at satfrq with power satpwr (use like dm, e.g., satmode='yyn' or satmode='ynn')

sspul='y' does trim(x)trim(y) to destroy all magnetization.

composit='y' uses composite 90° for pw (discriminates relative to B_1).

2.11 S2PUL—Standard Two-Pulse Sequence

The s2pul macro converts the current experiment into an experiment suitable for the S2PUL pulse sequence.

Applicability

S2PUL is available on all systems.

Parameters

A full description of S2PUL, including a diagram of the sequence, is provided in the manual *Getting Started*.

2.12 S2PULR—Standard Two-Pulse in Reverse Configuration

The s2pulr macro sets up parameters for a S2PULR pulse sequence. Figure 21 is a diagram of the sequence. The local oscillator (L.O.) signal must be taken from the decoupler board. No decoupling is supported in this sequence.




Figure 21. S2PULR Pulse Sequence

Applicability

S2PULR applies to VXR and UNITY systems only.

Parameters

The observe channel uses the decoupler hardware and is controlled by the parameters dn (which must be set to 'H1'), dof, dpwr (or dhp), p1, and pw.

The macros movetof and movesw cannot be used with this sequence except in the following way: tof=dof movetof (or movesw) dof=tof.

Technique

- 1. Acquire a ¹H spectrum in the normal manner.
- 2. Enter **s2pulr** to convert the parameter set.

Chapter 3. Multidimensional NMR

Sections in this chapter:

- 3.1 "Interferograms," this page
- 3.2 "2D NMR Step-by-Step," page 88
- 3.3 "Phase-Sensitive 2D NMR," page 89
- 3.4 "Data Acquisition: Arrayed 2D," page 90
- 3.5 "Weighting," page 93
- 3.6 "Phasing Before the 2D Transform," page 95
- 3.7 "Baseline Correction," page 96
- 3.8 "Processing Phase-Sensitive 2D and 3D Data," page 99
- 3.9 "2D and 3D Linear Prediction," page 105
- 3.10 "Phasing the 2D Spectrum," page 106
- 3.11 "Display and Plotting," page 107
- 3.12 "Interactive 2D Color Map Display," page 113
- 3.13 "Interactive 2D Peak Picking," page 118
- 3.14 "3D NMR," page 127
- 3.15 "4D NMR Acquisition," page 132

In some respects, 2D NMR is similar to an arrayed 1D experiment. In both, as a function of time (one of the time variables in the pulse sequence), we obtain a series of FIDs that we then transform to become a series of spectra. For 2D experiments, however, the times for each experiment are *not* explicitly specified. Instead, two new parameters are used: sw1, which describes our "2D" spectral width (to be discussed shortly), and ni, the number of increments, which sets the number of different experiments we will do. The implicit time variable will then be incremented from experiment to experiment as determined by sw1.

3.1 Interferograms

Once the data are obtained and transformed along the acquisition dimension, we have a series of spectra. If this were a 1D arrayed experiment, like an inversion-recovery T_1 experiment, we would see that the peak heights behave exponentially as a function of time. In 2D experiments, however, the peaks heights will oscillate as a function of time, and that oscillation is the information of interest. To unravel this information, we first transpose the matrix to form a series of *interferograms*.

Each interferogram contains a series of points that represent the peak height at a particular frequency in the original spectrum as a function of time. Of course, most of these interferograms contain only noise, because many of the frequencies in the original

Chapter 3. Multidimensional NMR

spectrum also contained noise. However, some interferograms, namely those corresponding to the peaks in the original spectrum, contain useful information.

The time that is varied in a 2D experiment is known as the *evolution time* or t_1 , because it is the first of two key time periods in the 2D experiment. Evolution time is controlled in VNMR software by the parameter d2. This time is normally calculated by setting the number of increments to the value of the parameter ni and the increment value to 1/swl. The value of ni determines if a 2D experiment will be run. Initially, d2 can be set to any value but is usually set to zero.

The d2 array does not appear in the display da (i.e., d2 is "implicitly" arrayed). Only the first value of d2 appears as the parameter value in the display dg. A minimum of eight increments must be used for ni to do a 2D transform. Typical numbers range from 32 to 512. The addpar('2d') macro creates the 2D acquisition parameters ni, sw1, and phase (par2d functions the same).

The time during which the signal is detected is known as the *detection time* or t_2 , because it is the second of the two key time periods. After transform of the signals detected during the time t_2 , the "normal" spectrum appears along the f_2 axis. The second transform reveals information about the frequencies of oscillations during the t_1 time period along the f_1 axis.

Many parameters that refer to the new f_1 axis in a 2D experiment are identified by the number 1 (e.g. sw1, lb1, fn1), whereas the normal 1D parameters control f_2 .

The process of transformation, transposition to interferograms, and second transformation may seem complicated, however, it can all be reduced to literally a single command, or even a single menu choice, that starts an acquisition of a 2D experiment and performs all the necessary processing when the experiment is done. So the process can be fairly simple.

3.2 2D NMR Step-by-Step

The examples in this section provide an excellent way to get started in 2D NMR. The 2D concepts involved are described in greater detail later in this chapter.

To Process Stored Data

1. In the VNMR menus, click on **Main Menu** > **Workspace**.

The Workspace menu appears on the lower row of menus.

- 2. Join the experiment in which you want to perform the 2D processing by clicking on the button for the experiment: Exp2, Exp3, etc. If the experiment does not exist, click on Create New, then click on the button to join the experiment.
- 3. Enter rt('/vnmr/fidlib/fid2d').

The text window displays the parameters values for the 2D sample data.

4. Click on Process > Select Params > Sinebell > Large > Return.

This establishes reasonable initial choices for the processing parameters for your experiment. A number of other choices are available for processing.

5. Click on Full Transform.

This initiates 2D processing. The processing ends with the display of a color map of the 2D frequency data.

- 6. Adjust the vertical scale vs2d and threshold th to display the data in an appropriate fashion. Each time you change vs using the keyboard, be sure to click the Redraw button to redisplay the color map.
- 7. To symmetrize the data, enter **foldt**.
- 8. To produce a 2D contour plot, enter **pcon page**.

To Acquire a Simple COSY

1. Perform a 1D ¹H experiment on a sample of interest. After processing, display two cursors and place them around a region containing only the peaks, omitting most of the baseline region at both ends of the spectrum.

Certain peaks, such as a small residual CHCl₃ peak in a spectrum otherwise containing only aliphatic resonances, can be ignored if they are far enough away not to "fold back into" the spectrum.

2. Enter movesw ga.

This moves the spectrum window and then obtains a narrowed spectrum.

- 3. Use the mp(<from,>to) command to move the 1D parameters to another experiment. For example, enter mp(1,4) if you performed the 1D experiment in exp1 and you want to perform the 2D experiment in exp4. Then join the 2D experiment by entering, for example, jexp4.
- 4. Enter **cosy**.

This is a macro that will modify your parameter set to perform an absolute-value COSY experiment. The total estimated time for the experiment should be displayed on the screen. If you wish, change nt and type time to alter the conditions. For your first experiment, nt=16, ss=1, and dl=1 is recommended.

5. Enter **wexp='do2d'**.

The value of the wexp parameter specifies the end of experiment process.

6. Enter **au**.

This starts the acquisition. When the data acquisition is complete, the data is automatically processed and displayed.

- 7. Adjust the vertical scale vs and threshold th to display the data in an appropriate fashion. Most of the screen should be black when vs and th are properly adjusted, with colored spots representing the diagonal and cross peaks. Each time you change vs using the keyboard, be sure to enter dconi to redisplay the contour plot. The parameter th can be adjusted interactively.
- 8. To symmetrize the data, enter **foldt**.
- 9. To plot the data, enter **pcon page**.

3.3 Phase-Sensitive 2D NMR

For many years, 2D NMR experiments were performed and displayed in an absolute- value mode. Just as in 1D NMR, phase-sensitive processing and display offers better sensitivity and resolution. Phase-sensitive 2D NMR by itself is simply the ability to display and plot *phased data*, as opposed to absolute-value data. There are four kinds of experiments in which a user might want to examine phase-sensitive data:

- A 2D experiment in which the data are not expected to appear in absorption mode in both directions, but in which it is nonetheless desirable to observe the data in a phase-sensitive presentation.
- A 2D experiment in which the data, processed in a suitable way, *are* expected to appear in absorption mode in both directions. Heteronuclear 2D-J is such an experiment.
- An experiment in which two different experiments are performed for each value of t₁, typically using different phase cycles, producing a full complex data set for the second transformation. We shall refer to this method, popularized by States, Haberkorn, and Ruben (*J. Magn. Reson.* **1982**, *48*, 286), as the *hypercomplex method*
- An experiment in which the phase of the excitation pulse is incremented as a function of t₁ (TPPI or Time Proportional Phase Incrementation, see Marion and Wuthrich, *Biochem. Biophys. Res. Commun.* **1983**, *113*, 967) and which produces real data along the t₁ axis.

Each of the experiments referred to above can be performed and processed with the VNMR software. Complex transforms are usually performed along t_1 , which is the ideal situation for the hypercomplex method. TPPI data can be processed along t_1 with either a complex FT or a real FT, depending upon the method of data collection.

In general, the hypercomplex method is the method of choice using VNMR software. A natural first reaction, since this method requires *two* data tables instead of one, is to assume that it requires twice as much storage as the TPPI method. This is untrue, however, for the same reason that a real 1D transform covering a given spectral width requires exactly as much data as a complex 1D transform— the sampling rate must be twice as high in the real case to produce the same result. In the same way, the TPPI method requires only one data table, but requires sampling to occur twice as frequently along t_1 , thereby incurring twice the data size per data table to produce the same real resolution. So in this sense the two experiments are equivalent in data storage requirements and experimental time.

3.4 Data Acquisition: Arrayed 2D

2D experiments have one implicitly "arrayed" parameter, d2. Like 1D experiments, however, 2D experiments can also have other parameters explicitly arrayed. This feature can be used, of course, for purposes that have nothing to do with phase-sensitive 2D, such as running a series of 2D-NOE experiments using different mixing times. As we shall see below in discussing the processing of such data, this feature alone opens up a variety of experiments, including addition/subtraction of two or more 2D experiments. Table 16 lists commands and parameters associated with arrayed 2D and 3D data acquisition.

Hypercomplex Method

The hypercomplex method of phase-sensitive 2D NMR requires the use of two data tables. Appropriate pulse sequences must be created (see below for more details of this point) which, as a function of some parameter, generate a different sequence of pulses or pulse phases suitable for generating the two component experiments of the hypercomplex method.

Any parameter may be used for this purpose. As a convention, we use the parameter phase, which takes on values of 0, 1, or 2:

• A value of phase=0 can be used to produce a phase cycle suitable for a conventional (non-phase-sensitive) 2D experiment.

Commands	
addpar('2d' '3d')	Add 2D or 3D parameters to the current experiment
par2d	Create 2D acquisition parameters
wft2d*	Weight and Fourier transform 2D data
* wft2d<(<options,>coe</options,>	efficients)>
Parameters	
d2*	Incremented delay for 1st indirectly detected dimension
d3 {number, in seconds}	Incremented delay for 2nd indirectly detected dimension
ni {number}	Increments in 1st indirectly detected dimension
ni2 {number}	Increments in 2nd indirectly detected dimension
phase {number}	Phase selection
phase2 {number}	Phase selection for 3D acquisition
sw1 {100 to 100000, in Hz}	Spectral width in 1st indirectly detected dimension
sw2 {100 to 100000, in Hz}	Spectral width in 2nd indirectly detected dimension
* d2 {0 to 4095, in sec} GEMIN	<i>I 2000</i> , {0 to 8190, in sec} other systems

Table 16. Arrayed 2D & 3D Data Acquisition Commands and Parameters

• Running instead an array of experiments with phase=1, 2 produces the two experiments suitable for the hypercomplex method.

You may ascertain the possible values of phase by reading the source code in the psglib directory for any particular pulse sequence.

TPPI Method

The TPPI method of phase-sensitive 2D NMR requires one data table when phase=3. The data must be processed along t_1 with a complex Fourier transform by setting proc1 (which sets the type of data processing to be performed on the t_1 interferogram) to 'ft'. This manner of implementing TPPI leads to a doubling of the f_1 frequency axis.

When an arrayed 2D experiment is run in this manner, there is in reality a double array: d2 (the evolution time) and phase. The order of these arrays is such that the phase array is cycled the most rapidly, so that the order of the experiments is, for example:

d2=0	phase=1	
d2=0	phase=2	States-Haberkorn
d2=1/sw1	phase=1	
d2=1/sw1	phase=2	
10 0		
d2=0	phase=3	
d2=0 d2=1/sw1	phase=3 phase=3	TPPI (non-arrayed)
	-	TPPI (non-arrayed)

Not all pulse sequences have the TPPI method incorporated.

When an experiment is in progress, the acquisition status window displays a count of the current FID and the number of completed transients (ct) in that FID. As indeed happens with a 1D arrayed experiment, the current FID number is actually the *total* count of the completed FIDs to this point, including all arrays. Since the phase parameter is cycling the most rapidly, and since typically phase is an array of two values, the current FID number is typically *twice* the number of the current increment. For example, when the counter reads FID 54, this means that 27 FIDs of the first type of experiment have been

completed, 26 of the second type, and the system is working on the 27th experiment of the second type.

Real-Time 2D

VNMR software can perform *real-time 2D*, that is, 2D actions while the experiment is still in progress. Once eight or more increments have been completed, you can enter the wft2d command to perform the full 2D transform on the data that exists up to that point.

For some experiments, such as heteronuclear chemical shift correlation and homonuclear 2D-J experiments, you will be surprised at how few increments are needed to resolve the resonances of interest. Others may require more increments. In any case, if you find you have sufficient data to solve the problem, you can always abort the experiment, so that the remaining increments are not performed, and proceed to the next problem.

Macros for 2D Experiments

Table 17 lists pulse sequences for 2D experiments that have a macro on the system to help you set up the experiment.

Macros	
cosy	Set up COSY, correlated spectroscopy
cosyps	Set up COSYPS, phase-sensitive COSY
dqcosy	Set up DQCOSY, double-quantum filtered COSY
het2dj	Set up HET2DJ, heteronuclear 2D-J experiment
hetcor<(exp_num)>	Set up HETCOR, heteronuclear chemical shift correlation
hmqc<(isotope)>	Set up HMQC, reverse detection heteronuclear multiple quantum
hmqcr	Set up HMQCR, HMQC with "reverse configuration"
hom2dj	Set up HOM2DJ, homonuclear J-resolved 2D experiment
inadqt**	Set up INADEQUATE experiment
mqcosy<(level)>	Set up MQCOSY, multiple quantum filtered COSY
noesy	Set up NOESY, laboratory frame Overhauser experiment
relayh	Set up absolute valueCOSY or single or double RELAY-COSY
roesy<(ratio)>**	Set up ROESY, rotating frame Overhauser experiment
tncosyps*	Set up TNCOSYPS, COSYPS with water suppression
tndqcosy*	Set up TNDQCOSY, DQCOSY with water suppression
tnmqcosy*	Set up TNMQCOSY, MQCOSY with water suppression
tnnoesy*	Set up TNNOESY, NOESY with water suppression
tnroesy*	Set up TNROESY, ROESY with water suppression
tntocsy*	Set up TNTOCSY, TOCSY with water suppression
tocsy**	Set up TOCSY, total correlation spectroscopy
	URY-VX, MERCURY, and GEMINI 2000 systems
** Not available on GEMI	NI 2000 system

Each macro displays suggestions on how to select values for the parameters in that experiment. Complete listings of each pulse sequence are found on your system disk in the directory /vnmr/psglib. Note that the INADEQUATE, MQCOSY, ROESY, and TOCSY sequences are not available on *GEMINI 2000* systems. None of the tn sequences are available on *GEMINI 2000, MERCURY-VX*, or *MERCURY* systems.

For further information, including diagrams and literature references, see Chapter 4, "Multidimensional and Advanced Experiments,".

3.5 Weighting

This section describes weighting functions for processing in the ni and ni2 dimensions.

Parameters

Each 2D weighting parameter (the parameters ending with "1") listed in Table 18 is analogous to a similarly named parameter for 1D experiments listed the manual *Getting Started* (the same parameter name but without the 1).

Parameters	
awc1 {'n', number}	Additive weighting const. in 1st indirectly detected dimension
awc2 {'n', number}	Additive weighting const. in 2nd indirectly detected dimension
fn1 {'n', number}	Fourier number in 1st indirectly detected dimension
fn2 {'n', number}	Fourier number in 2nd indirectly detected dimension
gf1 {'n', number in sec}	Gaussian function in 1st indirectly detected dimension
gf2 {'n', number in sec}	Gaussian function in 2nd indirectly detected dimension
gfs1 {'n',number in sec}	Gaussian shift constant in 1st indirectly detected dimension
gfs2 {'n',number in sec}	Gaussian shift constant in 2nd indirectly detected dimension
<pre>lb1 {'n',number in Hz}</pre>	Line broadening in 1st indirectly detected dimension
1b2 {'n',number in Hz}	Line broadening in 2nd indirectly detected dimension
<pre>sb1 {'n',number in sec}</pre>	Sinebell constant in 1st indirectly detected dimension
sb2 {'n',number in sec}	Sinebell constant in 2nd indirectly detected dimension
<pre>sbs1 {'n',number in sec}</pre>	Sinebell shift constant in 1st indirectly detected dimension
<pre>sbs2 {'n',number in sec}</pre>	Sinebell shift constant in 2nd indirectly detected dimension
<pre>wtfile1 {",file}</pre>	User-defined weighting in 1st indirectly detected dimension
wtfile2 {",file}	User-defined weighting in 2nd indirectly detected dimension

Table 18. Weighting Parameters for ni and ni2 Dimensions

Parameters for 2D experiments are used for processing the t_1 domain (the interferograms) or ni (first indirectly detected) dimension. Parameters ending with a "2" are used for processing the ni2 (second indirectly detected) dimension.

In non-phase-sensitive (absolute-value and power) 2D experiments, "pseudo-echo," sinebell, or sinebell-squared weighting is typically used to attenuate long dispersion tails. This weighting is often responsible for a significant loss in sensitivity in such 2D experiments.

In phase-sensitive 2D experiments, the key in using weighting functions is to ensure that the weighted FID or interferogram decays to zero by the end to avoid "truncation wiggles." The Gaussian function (gf and gfl parameters) is ideally suited for this; typical values might be gf=0.6*at, gfl=0.6*ni/swl(=0.6*atl). Resolution enhancement (using negative lb or negative lb1) may be helpful in cases of spectral overlap, but can also be dangerous, since the "dips" that it can induce around the sides of peaks show up as peaks of opposite sign in the 2D plot, complicating analysis.

Setting Values

Setting values for 2D weighting parameters can be done through various macros and an interactive program (described in the next section). Table 19 lists these tools.

The following macros are available to generate values for 2D weighting parameters:

Commands	
pseudo<(C1,C2,C3,C4)>	Set parameters for pseudo-echo weighting
<pre>sine<(shift<,npoints<,domain>)></pre>	Find values for a sine window function
sinebell	Select default parameters for sinebell weighting
<pre>sinesq<(shift<,npoints<,domain>)></pre>	Find values for a sine-squared window function
<pre>wti<(element_number)></pre>	Interactive weighting
<pre>wtia<(element_number)></pre>	Interactive weighting for 2D absorptive data
Macros	
gaussian	Set up unshifted Gaussian window
pi3ssbsq	Set up pi/3 shifted sinebell-squared window
pi4ssbsq	Set up pi/4 shifted sinebell-squared window
sqcosine	Set up unshifted cosine-squared window
sqsinebell	Set up unshifted sinebell-squared window

Table 19.	Commands	and Macros	for	Setting	2D	Weighting	Values

- gaussian sets up an unshifted Gaussian window function in 1, 2, or 3 dimensions. This macro checks whether data is 1D, 2D, or 3D. The argument t1_inc is the number of t1 increments; the default is ni. The argument t2_inc is the number of t2 increments; the default is ni2.
- pi3ssbsq and pi4ssbsq are macros that respectively set up pi3 and pi4 shifted sinebell-squared window functions in 1, 2, or 3 dimensions and check whether data is 1D, 2D, or 3D. Both macros use the t1_inc and t2_inc arguments, described in the previous description of the gaussian macro.
- pseudo<(C1, C2, C3, C4)> generates an initial guess at good values for lb, gf, lb1, and gf1 for absolute-value 2D experiments. To generate modified guesses, four coefficients (described in the VNMR Command and Parameter Reference) are available.
- sinebell sets sb and sb1 to one-half the acquisition time. Other weighting is turned off. Use sinebell in absolute-value 2D experiments only.
- sine(<shift><, npoints><, domain>) calculates appropriate values for sb and sbs (if argument domain is 'f2') or for sb1 and sbs1 (if domain is 'f1') to achieve a sine window function. If shift is greater than 0, the starting value for the window function is given by sine(pi/shift); otherwise, the starting value is 0. npoints specifies the number of real points the window function spans. If domain is not specified, the value of the parameter trace is used as the default.
- sinesq(<shift><,npoints><,domain>) calculates appropriate values for sb and sbs (if domain='f2') or for sb1 and sbs1 (if domain='f1') to achieve a sine squared window function. The arguments are used the same as the sine command.
- sqcosine and sqsinebell are macros that respectively set up unshifted cosinesquared and sinebell-squared window functions in 1, 2, or 3 dimensions and check whether data is 1D, 2D, or 3D. Both macros use the t1_inc and t2_inc arguments, described in the previous description of the gaussian macro.

Interactive Weighting

The wti<(element_number)> command allows interactive setting of weighting parameters for both t₂ FIDs and t₁ interferograms (both the ni and ni 2 dimension). The optional argument element_number specifies which FID element or interferogram

trace is to be used in adjusting the weighting parameters. The default value is the currently active element or trace. wti responds appropriately to phfid and lsfid for t₂ FIDs, phfid1 and lsfid1 for t₁ interferograms defined by ni, and lsfid2, and phfid2 for t₁ interferograms defined by ni2.

The following parameters are used with interactive weighting:

- awc, awc1, and awc2 set the additive weighting constant; added in to the weighting function after the lb and sb (sbs) contributions but before the gf (gfs) contributions.
- gf, gf1, and gf2 set the Gaussian apodization constant, in seconds.
- gfs, gfs1, and gfs2 set the Gaussian function shift, in seconds. This shifts the origin of the Gaussian function; active only if gf (or gf1) is active.
- 1b, 1b1, and 1b2 set the line broadening factor, in Hz; a positive value gives sensitivity enhancement; a negative value gives resolution enhancement.
- sb, sb1, and sb2 set the sinebell time period, in seconds; a negative value gives a sine squared bell.
- sbs, sbs1, and sbs2 set the sinebell shift, in seconds; shifts the origin of the sinebell; active only if sb (or sb1) is active.

These parameters can be typed in or changed with the left mouse button in the proper field. The right mouse button turns off the spectrum for a faster response to changes in the weighting function.

The wtia<(element_number)> command allows the same weighting parameters to be set interactively for 2D absorptive data. The argument element_number is used the same as in the wti command.

3.6 Phasing Before the 2D Transform

Table 20 summarizes the commands and parameters discussed in this section.

For a phase-sensitive 2D display, only the ph (phase-sensitive along f_2) command and the ph1 (phase-sensitive along f_1) command are relevant to ensure that phasing either is performed during the 2D FT or can be performed after the 2D FT along each dimension. ph and ph1 should be executed according to the following rules:

- pmode=' ' is both ph and ph1 must be executed before performing 2D FT.
- pmode='partial' is only ph must be executed before performing 2D FT.
- pmode='full' is either command must be executed before 2D FT.

To obtain pure 2D absorptive lineshapes requires a properly phased spectrum in the first dimension, that is, along f₂. To facilitate this operation, the wft and ft programs (actually the same program) allow normal 1D transforms on 2D data. After an experiment is complete, or while it is in progress, type wft(1) or ft(1). The first FID will be transformed. After the transform, the spectrum will be displayed. Phase this spectrum in the usual way. If pmode='full', all phasing along f₂ can be performed after the 2D FT if corrections in f₂ phasing are required. To prevent automatic spectral display (ds), type wft('nods').

One subtle point remains. In, for example, a COSY experiment acquired using the hypercomplex method, the first spectrum of the first array element contains data resulting from the summing of the signal from 90_X-90_X and 90_X-90_{-X} pulse sequences. But both of these experiments produce in principle *no* signal in the *xy* plane! In this case, phasing is

Commands					
ds*	Display a spectrum				
ft*	Fourier transform 1D data				
ft2d*	Fourier transform 2D data				
ph	Set phased mode along directly detected dimension				
phl	Set phased mode along 1st indirectly detected dimension				
wft*	Weight and Fourier transform 1D data				
* ds<(index)>, ds<(op	tions)>				
ft<(<options,><'nf'</options,>	<<,start><,finish><,step>)>				
ft('inverse',exp_nu	mber,expansion_factor),				
ft2d(array_element)					
ft2d<(<options,><pl< td=""><td colspan="5"><pre>ft2d<(<options,><plane_number,><coefficients>)>,</coefficients></plane_number,></options,></pre></td></pl<></options,>	<pre>ft2d<(<options,><plane_number,><coefficients>)>,</coefficients></plane_number,></options,></pre>				
ft2d('ni' 'ni2',ele	<pre>ft2d('ni' 'ni2',element_number,increment),</pre>				
('ni' 'ni2',increment, <coefficients>)</coefficients>					
wft<(<options,><'nf'><,start><,finish><,step>)>,</options,>					
<pre>wft('inverse',exp_number,expansion_factor)</pre>					
Parameters					
lp1 {-3600 to 3600, in deg}	First-order phase in 1st indirectly detected dimension				
1p2 {-3600 to 3600, in deg} First-order phase in 2nd indirectly detected dimensio					
pmode {",'partial','full'}	Processing mode for 2D data				
rp1 {-360 to 360, in deg}	Zero-order phase in 1st indirectly detected dimension				
rp2 {-360 to 360, in deg} Zero-order phase in 2nd indirectly detected dimension					

Table 20. Commands and Parameters for Phasing Before the 2D Transform

performed on the first spectrum of the *second* experiment, which is displayed by ds(2). Different experiments may require different first domain phasing procedures. With the standard pulse sequences, both array elements can be phased for pure absorption simultaneously. For cases where one data set is "in phase" and the other is "out of phase" for a given set of phase parameters (lp and rp), the ft2d program can be instructed to extract the properly phased data from each experiment.

After the full transform, f_1 phasing is possible using the lpl and rpl parameters if pmode has been set to 'partial' or 'full'. f_1 phasing for the ni2 dimension is done using the parameters lp2 and rp2. Once satisfactory f_1 and f_2 phasing has been obtained, future retransforms may be done with pmode='' (two single quotes with no space between the quotes). This results in a faster initial display of the processed data.

3.7 Baseline Correction

Table 21 lists commands and parameters covered in this section.

The ft and ft2d commands (and related commands) multiply the first point of each FID by fpmult (the default value is 1.0, except that if the processing involves backward extension of the time-domain data with linear prediction, the default value is then 0.5) and the first point of each interferogram by fpmult1 (default value is 0.5) for the ni dimension or fpmult2 (default value is 0.5) for the ni 2 dimension. fpmult attempts to compensate for the first point distortion caused by analog filters (see Otting, Widmer, Wagner and Wüthrich, *J. Magn. Reson.* **1986**, *66*, 187).

The effect of using the fpmult is to perform a linear baseline correction on all f_2 data, reducing negative-going ridges along f_2 in phase-sensitive 2D data. This correction is not needed in experiments such as COSY where the FID *starts* at zero and grows or in absolute-

Commands			
bc*	2D baseline correction		
calfa	Recalculate alfa so that the left phase is zero		
cdc	Cancel drift correction		
cfpmult	Calculate first-point multiplier for 2D experiments		
crof2<(alfa)>	Recalculate rof2 so that lp=0		
dc	Calculate spectral drift correction		
dc2d('f1' 'f2')	Apply drift correction to 2D spectra		
ft*	Fourier transform 1D data		
ft2d*	Fourier transform 2D data		
wft*	Weight and Fourier transform 1D data		
<pre>* bc(trace_direction<,or</pre>	der><,trace_start><,trace_end>)		
ft<(<options,><'nf'><,</options,>	<pre>start><,finish><,step>)></pre>		
	er,expansion_factor),ft2d(array_element)		
	e_number,> <coefficients>)>,</coefficients>		
ft2d('ni' 'ni2',elemen			
<pre>ft2d('ni' 'ni2',increment,<coefficients>) wft<(<options,><'nf'><,start><,finish><,step>)>,</options,></coefficients></pre>			
wft('inverse',exp_numb			
Parameters	,,,		
alfa $\{0 \text{ to } 1e8, \text{ in } \mu s\}$	Set alfa delay before acquisition		
dcg {'dc', 'cdc'} Drift correction group			
fpmult {'n', number} First point multiplier for np FID data			
fpmult1 {'n', number}	First point multiplier for ni interferogram data		
fpmult2 {'n', number}	First point multiplier for ni2 interferogram data		
$1p \{-3600 \text{ to } 3600, \text{ in deg}\}$	First-order phase along directly detected dimension		
rof2 {0 to 8190, in ms}			

Table 21. Baseline and Drift Correction Commands and Parameters

value mode presentation if pseudo-echo or sinebell processing is used, because the processing function goes to zero at $t_2=0$, forcing all FIDs to start at zero amplitude.

Unless lp is approximately zero, fpmult will affect both the dc offset and the curvature of the spectrum during 2D data processing. Obtaining a trial spectrum and phasing it to pure absorption will provide from the spectrum the current values of the parameters alfa and lp. Using these parameters, the calfa macro can calculate a new value for alfa so that lp is rendered approximately 0.

The crof2 macro recalculates a new value for rof2 (receiver gating time following a pulse) based upon the current rof2 and lp (first-order phase) values, so that lp is rendered approximately 0. For crof2 to work properly, a trial spectrum must be obtained and phased to pure absorption. This spectrum provides the current rof2 and lp values for crof2. The value of the alfa delay is left constant, provided rof2 does not become less than 1 μ s.

First-Point Multiplier

The best value of fpmult is a function of the filter setting and should be determined empirically. It can be determined before, during, or after the 2D experiment by using wft(1).

1. With a properly phased first increment spectrum on the screen, enter dc.

- 2. Position the mouse-controlled arrow at the right edge of spectrum baseline (to keep track of the ideal baseline position).
- 3. Enter cdc and observe the new position of the baseline. It typically drops.
 - If the baseline goes negative, set fpmult to greater then 1.0 (try 1.5) and enter wft(1).
 - If the baseline rises but does not return to the position indicated by the mouse arrow, increase the value of fpmult and enter wft(1) again. If in doing so the baseline rises above the ideal level, reduce fpmult and try again.

Only a few tries are required before the proper value of fpmult is found.

The macro cfpmult is provided that does the procedure above automatically. This applies only to data in t_2 . No equivalent macro for t_1 data is provided. Normally, no correction for fpmult1 is necessary.

Setting fpmult='n' and fpmult1='n' disables these features of the ft and ft2d programs. This would be the usual value for sinebell or pseudo-echo processing.

Baseline Correction

An alternative to the use of fpmult is the use of the baseline correction command bc. The implementation of bc in 2D processing uses the spline or second to twentieth order polynomial fitting of predefined baseline regions. These regions are set up prior to the use of bc by setting integral resets (set intmod='partial') so that integrals appear only over regions of the spectrum with signals present. These may be set after a wft(1), as described in the manual *Getting Started*. The quality of the baseline correction may be assessed by bc(1). In setting baseline regions near the ends of the spectrum, the bc operation does essentially the equivalent of fpmult because this represents a simple dc correction.

If this mode is used, the wft2d command must have an argument 'bc' and an optional order value such as wft2da('bc',1). After the first (t₂) transform, the bc program executes, and the interferograms are then calculated from the baseline-corrected data. After the second transform, bc may be used again, this time along f₁ using the command bc('f1'). Make sure that the resets are appropriate. They will automatically be so if sw=sw1 in a homonuclear experiment. Of course, in a heteronuclear shift correlation experiment no proper reset points may be set, in general.

FID Drift Correction

A dc offset in time-domain data transforms into a "center glitch" in the frequency spectrum. For 1D data, the ft program automatically applies a dc correction to the FID. Such a correction is not applied to 2D FIDs or interferograms unless explicitly requested. The command ft2d('t2dc') causes a dc correction to be applied to each t₂ FID before the first FT and ft2d('t1dc') causes a dc correction to be applied to each t₁ interferogram prior to the second FT. In both cases, the last one-sixteenth of the time-domain data is used to calculate the dc correction.

The dcg parameter contains the results of the dc or cdc command. This parameter cannot be set in the usual way but it can be queried (dcg?) to determine whether drift correction is active.

Spectral Drift Correction

The command dc2d('f1'|'f2') is the 2D equivalent of dc and is run only after the 2D transform. Use dc2d('f1') for corrections along f_1 and dc2d('f2') for corrections along f_2 . The drift correction calculation is done separately for each trace in the 2D data set.

3.8 Processing Phase-Sensitive 2D and 3D Data

After 2D data has been acquired, the complete 2D transformation can be performed with a single command, with or without weighting as appropriate. Table 22 lists commands and parameters associated with this processing.

The general flow of information during the process of phase-sensitive 2D transformation is depicted in Figure 22 for pmode='partial'.




Figure 22. Data Flow in Phase-Sensitive 2D Transformation

A series of complex FIDs, obtained as a function of t_1 , are transformed to become a series of spectra. Each spectrum consists of a real and imaginary part. Each spectrum is then phase rotated, according to the phase correction determined from an individual spectrum. We now have a series of spectra, each consisting of an absorptive and a dispersive part, formed as linear combinations of the original real and imaginary parts. Complex interferograms are

Commands			
addpar('2d' '3d')	Add selected 2D or 3D parameters to current experiment		
dg2	Display group of 3rd and 4th channel/3D parameters		
foldt<('symm' 'triang')>	Fold COSY-like spectrum along diagonal axis		
ftld*	Fourier transform of 2D data		
ftlda<(options)>	Fourier transform "halfway" for pure absorption 2D data		
ft2d*	Fourier transform 2D data		
ft2da<(options)>	Fourier transform for pure absorption 2D data		
ft3d*	Perform a 3D FT on a 3D FID data set (VNMR, UNIX)		
par3d	Create 3D acquisition, processing, display parameters		
parlp	Create parameters for linear prediction		
wftld*	Weight and Fourier transform f ₂ for 2D data		
wftlda<(options)>	Weight and FT "halfway" for pure absorption 2D data		
wft2d*	Weight and Fourier transform 2D data		
wft2da<(options)>	Weight and Fourier transform for pure absorption 2D data		
<pre>* ftld(element_number), ft</pre>	<pre>ld<(<options,><coefficients>)></coefficients></options,></pre>		
ft2d<(<options,><plane_n< td=""><td>umber,><coefficients>)>,</coefficients></td></plane_n<></options,>	umber,> <coefficients>)>,</coefficients>		
ft2d('ni' 'ni2',element_	number,increment)		
ft2d('ni' 'ni2',incremen			
	_files><,'nocoef'><,plane_type>)> (VNMR)		
ft3d -e exp_number -f -r	-		
wftld(element_number), w	ftld<(<options,><coefficients>)></coefficients></options,>		
wft2d<(<options,>coeffic</options,>			
wft2d<(<options,>coeffic Parameters</options,>	ients)>		
wft2d<(<options,>coeffic Parameters daslp</options,>	ients)> Increment for t1 dependent first-order phase correction		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number}</options,></pre>	ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off',partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number, in Hz} lsfrq1 {number, in Hz}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq1 {number, in Hz} lsfrq2 {number, in Hz}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz} lsfrq1 {number, in Hz} pmode {",'partial','full'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum Processing mode for 2D data</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz} lsfrq1 {number, in Hz} pmode {",'partial','full'} proc {'ft', 'rft', 'lp'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum Processing mode for 2D data Type of processing on np FID</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off',partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz} lsfrq1 {number, in Hz} lsfrq2 {number, in Hz} pmode {",'partial','full'} proc1 {'ft', 'rft', 'lp'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum Processing mode for 2D data Type of processing on np FID Type of processing on ni interferogram</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn1 {'n',number} intmod {'off','partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz} lsfrq1 {number, in Hz} lsfrq2 {number, in Hz} pmode {",'partial','full'} proc1 {'ft', 'rft', 'lp'} proc2 {'ft', 'rft', 'lp'}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum Processing mode for 2D data Type of processing on np FID Type of processing on ni2 interferogram</pre>		
<pre>wft2d<(<options,>coeffic Parameters daslp fn {'n',number} fn1 {'n',number} fn2 {'n',number} intmod {'off',partial','full'} lsfid {'n',number} lsfid1 {'n',number} lsfid2 {'n',number} lsfid2 {'n',number} lsfrq {number, in Hz} lsfrq1 {number, in Hz} lsfrq2 {number, in Hz} pmode {",'partial','full'} proc1 {'ft', 'rft', 'lp'} proc2 {'ft', 'rft', 'lp'} sw {100 to 100000, in Hz}</options,></pre>	<pre>ients)> Increment for t1 dependent first-order phase correction Fourier number along directly detected dimension Fourier number in 1st indirectly detected dimension Fourier number in 2nd indirectly detected dimension Integral display mode Number of complex points to left shift np FID Number of complex points to left-shift ni interferogram Number of complex points to left-shift ni2 interferogram Frequency shift of the fn spectrum Frequency shift of the fn1 spectrum Frequency shift of the fn2 spectrum Processing mode for 2D data Type of processing on np FID Type of processing on ni2 interferogram Spectra width</pre>		

Table 22. Tools for Processing Phase-Sensitive 2D and 3D Data

then formed out of corresponding points along the frequency axis from each of the spectra, and transformed to produce the final 2D spectrum.

The eight lines with arrows in the center of Figure 22 represent coefficients used during the process of transformation. The real and imaginary part of the interferograms can be formed from any linear combination of the real and imaginary parts of one or more spectral sets after the first Fourier transformation. We shall refer to these coefficients below according to the following scheme: RR1 is the coefficient used to multiply the real part (first R) of spectra in set 1 (the 1) before it is added to the real part (second R) of the interferogram.

IR2 would thus represent the contribution from the imaginary part of spectra in set 2 to the real part of the interferogram, and so on.

For pmode= 'full', another set of complex interferograms are formed from these two sets of f_2 spectra. This set of interferograms is 90° out-of-phase in f_2 to the previous set and can be constructed without any *additional* coefficients.

Different experiments will require different coefficients. Some, such as heteronuclear 2D-J experiments, consist of only one FID and spectral set, and hence there will be a total of four coefficients. Others, including hypercomplex 2D experiments, will consist of two original data sets and hence a total of eight coefficients. Other experiments are possible with three or even more data sets, requiring in each case four times as many coefficients as the number of data sets (see the macro wft2dac).

If there are n data sets to be transformed, as in typical phase-sensitive experiments, 4n coefficients must be supplied. The first 2n coefficients are the contributions to the real part of the interferogram, alternating between real and imaginary parts of the successive data sets. The next 2n coefficients are the contributions to the imaginary part of the interferogram, in the same order.

Thus, using the definition that the first letter refers to the source data set, the second letter refers to the interferogram, and the number identifies the source data set, we have the cases shown in the table on the right.

Data Sets	Coefficient Order
1	RR1, IR1, RI1, II1.
2	RR1, IR1, RR2, IR2, RI1, II1, RI2, II2.
3	RR1, IR1, RR2, IR2, RR3, IR3, RI1, II1, RI2, II2, RI3, II3.

The coefficients are generally

1, 0, or -1, but other coefficients are acceptable. Any *real* coefficient can be used, and as many coefficients can be non-zero as is desired. Up to 32 coefficients can be supplied, which at four per data set allows the addition, subtraction, etc., of eight 2D data sets (that is, eight different phase cycles). See the macro wft2dac for more information.

Processing Programs

A number of processing programs are available:

- ftld(coefficients) performs only the first Fourier transformation along the f₂ dimension (without weighting) and matrix transposition, allowing the display of interferograms with the wti, dcon, and dconi commands.
- wftld(coefficients) functions the same as ftld except weighting is included.
- ft2d<(<option,>coefficients)> performs a complete transformation in 2D, without weighting, after 2D data has been acquired. If the first Fourier transformation has already been done using ft1d, wft1d, ft1da, or wft1da, then ft2d performs only the second (t₁) transform. 'ptype' or 'ntype' can be used as the first argument to select P-type or N-type peak selection. The coefficients argument are discussed below.
- wft2d<(<option,>coefficients)> performs the same as ft2d except weighting is included. To perform a normal 2D transform on the n-th element in an arrayed 2D experiment, type wft2d(n).
- ft2da<('bc',polynomial_order)>runs complete phase-sensitive Fourier transform after the 2D FID data has been acquired. 'bc' is a keyword to perform a

baseline correction on the f_2 spectra prior to the Fourier transform along f_1 . polynomial_order is the order of the polynomial used in the baseline correction.

- wft2da<('bc',polynomial_order)> functions the same as ft2da except weighting is included.
- ftlda functions the same as ft2da except a Fourier transform along f_1 is omitted.
- wftlda functions the same as ftlda except weighting is included.

For some 2D data sets, you can save much time by selectively transforming the t_1 interferograms. ft2d('f2sel') allows only preselected f_2 regions to be transformed along t_1 ; the t_1 interferograms in the non-selected f_2 regions are zeroed but *not* transformed. The same mechanism used to select baseline regions for baseline correction (bc) is used to select the f_2 regions that are to be transformed along t_1 . Set intmod='partial' and partition the integral of the spectrum into several regions. The even numbered f_2 regions, e.g., 2, 4, etc., will be transformed along t_1 ; the odd numbered ones will not be transformed along t_1 .

Unreliable peak heights can be caused by Fourier transformation of truncated time-domain data, instead of Fourier numbers fn and fn1 being too low, as might be intuitively expected. To obtain properly defined signals, take one of the following steps:

- · Collect data until the signal has decayed to zero in the time domain, or
- Transform the data with zero-filling (fn>=2*np, fn1>=4*ni).

Taking one of these steps is particularly important in 2D spectra with antiphase or dispersive signals, where underdigitization can lead to signal cancellation.

Common Coefficients for wft2d Processing

A magnitude-mode transform, in which the real part of the interferogram is formed from the real part of the spectra and the imaginary part of the interferogram is formed from the imaginary part, would require wft2d(1,0,0,1). Changing the sign of the imaginary part of the interferogram serves to change the effective direction of the f1 frequency axis, as is required for data in which N-type peaks are detected. This can be done with wft2d(1,0,0,-1).

In some experiments, including heteronuclear 2D-J, the basic data are purely amplitude modulated, with a starting amplitude of +1. After the first transformation and phasing operation are complete, the dispersion part of each spectrum serves only to produce a phase-twist in the final spectrum without contributing any information. Setting the imaginary part of the second transform to zero produces a pure absorption display in both domains. To do this, type wft2d(1,0,0,0) or ft2d(1,0,0,0).

In the hypercomplex method for pure absorption 2D data, we have two complete sets of spectra and must therefore provide eight coefficients to specify the composition of the interferograms. A typical execution of the method described by States, Haberkorn, and Ruben, assuming that the first spectrum of the first data block has been phased for absorption, requires wft2d(1,0,0,0,0,1,0) to produce pure absorption spectra. For this coefficient set, the standard macro wft2da (or ft2da) is supplied; thus, in the most common case, you do not need to enter coefficients but simply type wft2da.

Other manipulations of two data blocks are formatted similarly. A magnitude-mode 2D experiment that is the sum of the two different experiments can be constructed by wft2d(1,0,1,0,0,1,0,1). For a COSY experiment, this would produce the P-type experiment. Subtracting data block two from block one, which for a COSY experiment gives the N-type COSY, would be accomplished by wft2d(1,0,-1,0,0,1,0,-1). Thus

two different absolute-value 2D experiments (P-type and N-type), and a phase-sensitive 2D experiment, can all be produced from the *same* data set, without acquiring the data again.

Different combinations of data sets with appropriate phase cycling might allow selection of various quantum orders in a *single* experiment. Note that since the coefficients may be different from one, it is possible essentially to phase shift each experiment *separately* (phase shift the receiver) *after* the experiment is done. For TPPI data with phase=3, only one data set is collected, and the imaginary part of the second transform is set to zero: wft2d(1,0,0,0).

The parameters proc and proc1 can be used to select the type of processing to be performed along t₂ and t₁, respectively. proc (or proc1) accepts the value 'ft' (complex FT, the default if proc or proc1 is not defined), 'rft' (real FT), or 'lp' (linear prediction processing on complex data):

- All Varian data and simultaneously sampled Bruker FID data should be processed with a complex FT along t₂, i.e., proc='ft'.
- Sequentially sampled Bruker FID data should be processed with a real FT along t₂, i.e., proc='rft'.
- Varian hypercomplex (phase=1, 2) and standard TPPI (phase=3) 2D FID data should be processed along t₁ with proc1='ft'.
- If 2D Bruker data is converted with the convertbru program or if the new method for acquiring TPPI data on a Varian system is used, a real FT should be performed along t₁, i.e., proc1='rft'. The command wft2da, which is equal to wft2d('ptype', 1, 0, 0, 0, 0, 0, 1, 0), then properly processes 2D Bruker data.
- If the 'lp' processing option is selected, additional parameters must be set to fully define how the time-domain data is to be manipulated. Refer to the description of the addpar or parlp macros in the *VNMR Command and Parameter Reference* for more information.

Occasionally, you may process a non-phase-sensitive data set in the phased (ph) mode by accident. Especially if pseudo-echo weighting has been used, the resulting data will not appear very pleasant! This problem cannot be cured if pmode=' ' without reprocessing the data. If pmode='partial', enter the absolute-value mode command av1 and redisplay the data with dcon. If pmode='full', *any* mode of display along *both* frequency dimensions is fully accessible without having to reprocess the data.

Sign of f₁ Frequencies

Different experiments have the potential to produce sign reversals along f_1 . Certain programs in the system, however, require that the sign of frequencies be the same in both dimensions. These include the axis labeling programs and the symmetrization program foldt. A simple method exists to overcome this problem without altering the experiment to be performed. Leaving the sign of the real part of a transform unchanged while reversing the sign of the imaginary part has exactly the effect of reversing the f_1 frequencies after transformation. Thus, multiplying all of the coefficients for the imaginary spectral part after the first transform by -1 will reverse the f_1 frequencies as desired. As an even simpler alternative, the keyword 'ptype' may be supplied as an argument to the transform program such as wft2d('ptype'). In the absence of either the 'ntype' or 'ptype' keyword, wft2d defaults to the 'ntype' option.

Once a particular set of coefficients is decided upon, a macro can be used to eliminate the necessity of typing four or eight coefficients each time. For example, if you perform the heteronuclear 2D-J experiment described above, you could create an appropriate macro.

wft2da is in fact a macro equivalent to wft2d(1,0,0,0,0,0,1,0) in some cases. Other macros for other combinations of coefficients can be created (see wft2da for an illustration).

2D Solvent Subtraction Filtering

In a 2D transform, solvent subtraction is invoked on t_2 FIDs in the same manner as 1D usage. The parameters ssfilter and ssorder select the processing option as follows:

- The zfs (zero-frequency suppression) option is selected if both ssfilter and ssorder are set to a value other than "Not Used."
- The lfs (low-frequency suppression) option is selected if ssfilter is set to a value other than "Not Used" and ssorder is set to "Not Used."
- The zfs and lfs options are both turned off if ssfilter is set to "Not Used."

These options are used with the ftld, wftld, ft2d, and wft2d commands, or with the ftlda, wftlda, ft2da, and wft2da macros.

Left Shift, Frequency Shift, Phase Rotation

If the parameter lsfidl is set to a value other than 'n', the interferogram is left-shifted by lsfidl complex (or hypercomplex) points before weighting and Fourier transformation are performed. The value of lsfidl must lie between 0 and ni-1.

The parameter lsfrq sets a frequency shift of fn spectral data, in Hz, with a negative value resulting in peaks being shifted upfield (to the right) and a positive value in peaks being shifted downfield (to the left). lsfrq is the time-domain equivalent of lp within VNMR. lsfrq operates on complex np FID data, referred to as the t₂ dimension in a 2D experiment. Similarly, the parameters lsfrq1 and lsfrq2 set a frequency shift of the fn1 and fn2 spectrum, respectively.

If the parameter phfid1 is set to a value other than 'n', the interferogram is phaserotated by phfid1 degrees (zero-order phase rotation) before weighting and Fourier transformation are performed.

If the parameter daslp exists, each interferogram is phase-rotated by daslp times (interferogram number) degrees before weighting and Fourier transformation are performed. This phase-rotation has the effect of "shearing" f1 traces of a 2D data set.

2D Processing of 3D Data

Acquisition and full processing of 3D data using VNMR is available provided the parameters ni2 and sw2 have been created (d3 is the incremented delay in the ni2 dimension). Also available is 2D processing of "slices" of the 3D data matrix, which can be performed as described below.

ft2d('ni2') transforms non-arrayed 2D data that have been collected with ni2 and sw2 (instead of ni and sw1). The addpar('3d') macro creates the necessary processing parameters for the ft2d('ni2') operation (par3d functions the same as addpar('3rf')).

ft2d('ni', #) is used to selectively transform a particular np-ni 2D plane within a non-arrayed 3D data set; # is an integer that can range from 1 to ni2 in this example.

ft2d('ni2', #) is used to selectively transform a particular np-ni2 2D plane within a non-arrayed 3D data set; # is an integer that can range from 1 to ni in this example.

If an arrayed 3D data set is to be selectively processed, the format of the arguments to ft2d changes. For example, ft2d('ni', #1, #2) performs a 2D transform along np and ni of the #2-th ni2 increment and the #1-th element within the explicit array. This yields a 2D np-ni frequency plane. #1 ranges from 1 to ni2; #2 ranges from 1 to [arraydim/(ni*ni2)].

Arrayed 3D data sets can also be subjected to 2D processing to yield 2D absorptive spectra. If the States-Haberkorn method is used along both f_1 (ni2 dimension) and f_2 (ni dimension), there will generally be four spectra per (ni,ni2) 3D element. In this case, the command ft2d('ni2', #1, <16 coefficients>) would perform a 2D transform along np and ni2 of the #1-th ni increment using the ensuing 16 coefficients to construct the 2D t₁-interferogram from appropriate combinations of the four spectra per (ni,ni2) 3D element. Use the proc2 parameter to specify the type of data processing to be performed on the ni2 interferogram (3D): 'ft' for complex FT, 'rft' for real FT, or 'lp' for linear prediction processing on complex data. The macro dg2 displays 3D processing parameters.

3.9 2D and 3D Linear Prediction

Just as in 1D linear prediction, the technique of linear prediction can be used in 2D and 3D data processing. Table 23 lists 2D and 3D linear prediction (LP) parameters and macros to create and display the parameters. For more information on the specific parameters, refer to the manual *VNMR Command and Parameter Reference*.

Table 23.	2D	and 3D	Linear	Prediction	(LP)	Commands and Parameters
-----------	----	--------	--------	------------	------	-------------------------

Commands	
addpar('lp')	Add LP parameters to current experiment
dglp	Display group of LP parameters
parlp	Create parameters for LP
setLP1	Sets F1 linear prediction parameters.
Parameters	
<pre>lpalg1 {'lpfft', 'lparfft'}</pre>	LP algorithm for ni dimension
lpalg2 {'lpfft', 'lparfft'}	LP algorithm for ni2 dimension
<pre>lpext1 {number}</pre>	LP data extension for ni dimension
<pre>lpext2 {number}</pre>	LP data extension for ni2 dimension
lpfilt1 {number}	LP coefficients to calculate, ni dimension
lpfilt2 {number}	LP coefficients to calculate, ni2 dimension
lpnupts1 {number}	LP number of data points, ni dimension
lpnupts2 {number}	LP number of data points, ni2 dimension
lpopt1 {'b', 'f' }	LP algorithm data extensions, ni dimension
lpopt2 {'b', 'f' }	LP algorithm data extension, ni2 dimension
<pre>lpprint1 {number}</pre>	LP print output for ni dimension
<pre>lpprint2 {number}</pre>	LP print output for ni2 dimension
lptrace1 {number}	LP output spectrum, ni dimension
lptrace2 {number}	LP output spectrum, ni2 dimension
proc1 {'ft', 'rft', 'lp'}	Type of processing on ni interferogram
proc2 {'ft', 'rft', 'lp'}	Type of processing on ni2 interferogram
strtext1 {1 to $ni/2$ }	Starting point for LP data extension, ni dimension
strtext2 $\{1 \text{ to } ni2/2\}$	Starting point for LP data extension, ni2 dimension
<pre>strlp1 {number}</pre>	Starting point for LP calculation, ni dimension
<pre>strlp2 {number}</pre>	Starting point for LP calculation, ni2 dimension

Linear prediction parameters are created for the evolution axes by entering addpar('lp', 1) to create parameters for the t_1 axis, or by addpar('lp', 2) to create parameters for the t_2 axis (3D experiments only). Macros parlp(1) and parlp(2) function the same as addpar('lp', 1) and addpar('lp', 2), respectively.

The dglp macro displays the linear prediction parameters for both (or all three) domains. The parameter procl controls the transformation process along t_1 , and proc2 controls the transformation process along t_2 . Using the same method of transformation is not necessary along two (or three axes). You might, for example, employ a backwards linear prediction in t_2 of a 2D experiment and a forwards linear prediction along t_1 , or perhaps a simple Fourier transformation along t_2 and a backwards linear prediction along t_1 .

3.10 Phasing the 2D Spectrum

Table 24 summarizes the commands and parameters discussed in this section.

Commands ds<(index)>,ds<(options)>	Display a spectrum
Parameters lp1 {-3600 to +3600, in deg}	Left phase in 1st indirectly detected dimension
pmode {", 'partial', 'full'}	Processing mode for 2D data
rp1 {-360 to +360, in deg} trace {'f1', 'f2', 'f3'}	Zero-order phase in 1st indirectly detected dim. Mode for <i>n</i> -dimensional data display

 Table 24. Commands and Parameters for Phasing the 2D Spectrum

The phase constants lpl and rpl control the phase correction along f_1 in phase-sensitive data. In most 2D experiments, these should be near zero, but because of finite pulse widths and delays present in the pulse sequence, they may be far from zero. If the pulse sequence properly compensates for these pulse widths and delays, it is possible to have zero lpl and rpl. Most of the setup macro set lpl and rpl to zero so that the first display will indicate the need (if any) for phase correction in f_1 . The same techniques as used in 1D phasing are employed here, with a minor difference.

- 1. Enter **f full** to display the full data matrix in a full chart display.
- 2. To phase the 2D spectrum, use the horizontal cursor present in the interactive display to identify a peak toward the right-hand edge of the spectrum. Note the trace number indicated at the top of the display (you can "memorize" this by setting **r1** equal to its value.)
- 3. Select one or more other traces at f₁ values more toward the center and left parts of the spectrum. If there is a diagonal in the spectrum with large peaks, these will be the most sensitive with which to work. Use **r2**, **r3**, etc. to "memorize" these trace values. A minimum of two is needed, one at the far right and one at the far left.
- 4. Enter **ds(r1)**. Phase this spectrum as you would a 1D spectrum using the Phase button in the displayed menu. Click the mouse on the peak displayed near the right edge of the spectrum. Phase up this spectrum (thus setting rp1). Do not "click" in the left part of the spectrum at this time.
- 5. Enter **ds**(**r2**). The second trace appears. Click the mouse near the right edge of the spectrum (to fix rpl at the previously determined value) and do not rephase. Move the mouse to the peak at the left, click and phase it (thus setting lpl).

6. Enter ds(r1) to recheck rp1. Repeat the process again if necessary.

In homonuclear correlation spectra (such as NOESY, TOCSY, and ROESY), use the diagonal peaks for phasing. It there are strong cross-peaks, you can phase an f_1 trace exactly like a 1D spectrum. Phase HMQC spectra by progressively working from right to left, with several peaks selected along the way to make sure that lpl does not go through an extra revolution that would induce some baseline roll.

Corrections in f₂ phasing may be obvious in the 2D data when they are not in the first increment 1D spectrum. If pmode='full' before the 2D transform, f₂ phasing may be corrected without retransforming by setting trace='f2' and using the same approach as described for f₁ phasing. Transformation of the data again is necessary if pmode='' or pmode='partial'. No f₁ phasing is possible after transformation if $pmode=''; f_1$ rephasing after the transform is possible (but not f₂ rephasing) if pmode='partial'. Do baseline corrections such as dc2d or bc only after data are properly phased in f₁ and f₂.

3.11 Display and Plotting

This section discusses noninteractive 2D display and plotting. Table 25 lists the many commands and parameters available. Interactive 2D color map display (the dconi program) and interactive 2D peak picking (the ll2d program) are covered on page 113 and page 118, respectively.

Display Modes

The commands for 2D display modes are analogous to the command and modes in f_2 :

- ph1 selects the phase-sensitive mode in f_1 . This is equivalent to the ph mode in f_2 .
- pal selects the phase angle mode in f_1 . This is equivalent to the pa mode in f_2 .
- av1 selects the absolute-value mode in f_1 . This is equivalent to the av mode in f_2 .
- pwr1 selects the power mode in f_1 . This is equivalent to the pwr mode in f_2 .

The dmgl parameter stores the display mode ('phl', 'pal', 'avl', or 'pwrl') for f_1 . It can be set with the commands phl, pa, avl, and pwrl, respectively. If dmgl is not set to one of these values, the display mode for f_1 is then selected based upon dmg, the display mode parameter both for f_2 in 2D data and for 1D data.

Display and Plot Limits

The center, left, right, and full commands set the parameters sc, wc, sc2, and wc2 as appropriate to produce a display (and subsequent plot) in the relevant portion of the screen (and page) described by the command. The parameter sc2 is the start of the chart in the second direction and the parameter wc2 is the width of the chart. Together, sc2 and wc2 control the dimensions of the second axis (or the y axis) in a 2D contour plot.

The asize macro adjusts the parameters sc, wc, wc2, and sc2 so that the displayed resolution along both f_1 and f_2 is the same. The use of asize is not suggested for heteronuclear experiments where the chemical shift spread of one nucleus is much greater than that of the other.

Another command, fullt, sets display limits for a full screen with room for traces. The parameter spl is the start of the plot in f_1 , and wpl is the width of the plot. Both are analogous to the corresponding sp and wp 1D parameters.

Command					
acosy	Automatic analysis of COSY data				
acosyold	Automatic analysis of COSY data, old algorithm				
asize	Make plot resolution along f_1 and f_2 the same				
avl	Select AV mode along 1st indirectly detected dimension				
center	Set display limits for center of screen				
centerswl	Move cursor to center of spectrum in 1st indirect dimension				
crll	Clear reference line in 1st indirectly detected dimension				
dcon<(options)>	Display noninteractive color intensity map				
df2d<(array_index)>	Display FIDs of 2D experiment				
dpcon*	Display plotted contours				
dpconn*	Display plotted contours without erasing screen				
ds<(index)>, ds<(options)>	Display a spectrum				
ds2d<(options)>	Display 2D spectra in whitewash mode				
ds2dn<(options)>	Display 2D spectra in whitewash mode without erasing				
dss*	Display stacked spectra				
dsww*	Display spectra in whitewash mode				
foldcc	Fold INADEQUATE data about 2-quantum axis				
foldj	Fold J-resolved 2D spectrum about f1=0 axis				
foldt	Fold COSY-like spectrum along diagonal axis				
full	Set display limits for a full screen				
fullt	Set display limits for a full screen with room for traces				
grid<(<spacing,><colors>)></colors></spacing,>	Draw a grid on a 2D display				
image	Display noninteractive gray scale image				
imageprint	Plot noninteractive gray scale image				
left	Set display limits to left half of screen				
nm2d<(noise_mult)>	Normalize intensity of 2D spectrum				
noisemult<(noise_mult)>	Control noise multiplier for automatic 2D processing				
pal	Set phase angle mode in 1st indirectly detected dimension (C)				
pa pacosy	Set phase angle mode in indirectly detected dimension (C) Plot automatic COSY analysis				
pcon*	Plot contours on a plotter				
peak2d*	Return information about maximum in 2D data				
phl	Select phased mode on 1st indirectly detected dimension				
pl*	Plot spectra				
pl2d*	Plot 2D spectra in whitewash mode				
plgrid*	Plot a grid on a 2D plot				
plww*	Plot spectra in whitewash mode				
<pre>* dpcon(<options,><levels,sp< pre=""></levels,sp<></options,></pre>	*				
dpconn(<options,><levels, spacing="">)</levels,></options,>					
dss<(<start,finish<,step>></start,finish<,step>	<,options>)>				
dsww<(<start,finish<,step>><,'int'>)></start,finish<,step>					
<pre>pcon<(<'pos' 'neg'><,'noaxis'><,levels><,spacing>)></pre>					
<pre>peak2d:\$maximum_intensity<,\$trace,\$point> </pre>					
<pre>pl<(<start,finish<,step>><,'int'><,'all'><,options>)> pl2dt<(upobagoul_fillul_fillpbu)></start,finish<,step></pre>					
pl2d*<('nobase' 'fill' 'fillnb')> plgrid<(<spacing><,><pen>)></pen></spacing>					
<pre>plgrid<(spacing><,>>pen>)> plgrid<(start_f2,increment_f2,start_f1,increment_f1<,pen>)></pre>					
<pre>plywe<(start,finish,step><,'all'>)></pre>					
Continued on next page					

Table 25.	2D 3	Display	and Plotting	Commands and	Parameters	(Part 1	of 2)
-----------	------	---------	--------------	--------------	------------	---------	-------

0	
Command proj*	Project 2D data
	5
pwr1	Select power mode on 1st indirectly detected dimension
right	Set display limits to right half of screen
rl1<(frequency)>	Set reference line in 1st indirectly detected dimension
rotate<(number_degrees)>	Rotate 2D data
setswl*	Set spectral width in 1st indirectly detected dimension
<pre>* proj(exp_number<,'sum'></pre>	
setsw1(nucleus,downfiel	dppm,upfieldppm):offset
Parameter	
axis*	Axis label for displays and plots
cr1 {frequency, in Hz}	Cursor position in 1st indirectly detected dimension
dmg {'pa','ph','av','pwr'}	Display mode for data
dmg1 {'pa1','ph1','av1','pwr1'}	Display mode along 1st indirectly detected dimension
sc {0 to wcmax, in mm}	Start of chart
sf {0 to the value of at, in sec}	Start of FID
sf1 {0 to (2xni)/sw1, in sec}	Start of interferogram in 1st indirectly detected dimension
sc2 {0 to wc2max, in mm}	Start of chart in second direction
sp {number, in Hz}	Start of plot in directly detected dimension
sp1 {number, in Hz}	Start of plot in 1st indirectly detected dimension
$scalesw1 \{'n', number > 0.0\}$	Scale spectra width in 1st indirectly detected dimension
trace {'f1','f2','f3'}	Mode for n-dimensional data display
vs2d	Vertical scale for 2D displays
wc {5 to wcmax, in mm}	Width of chart
wc2 {number, in mm}	Width of chart in second direction
wf {0 to the value of at, in sec}	Width of FID
wf1 {0 to $(2xni)/sw1$, in sec}	Width of interferogram in 1st indirectly detected dimension
wp{number, in Hz}	Width of plot in directly detected dimension
wp1 {number, in Hz}	Width of plot in 1st indirectly detected dimension
* axis {'1','2','3','c','d','h','k','m','n','p'	
	, - ,

Table 26. 2D Display and Plotting Commands and Parameters (Part 2 of 2)

Maximum Intensity

The peak2d command searches the area defined by sp, wp, sp1, and wp1 in a 2D data set for a maximum intensity. It returns the maximum intensity value found, the trace number of the maximum, and the data point number of the maximum on that trace.

Axis Label and Direction

For 2D, the parameter axis has two letters, with the first letter describing the detected spectral axis (f₂), and the second letter describing the indirectly detected axis (f₁). The special letter d is added to reference any indirectly detected axis to the parts per million of the decoupler channel. The special letter e is added to reference any indirectly detected axis to the parts per million of the second decoupler channel. The letter n is used to suppress the axis display on one or both axes.

The parameter trace selects the horizontal axis. trace='f1' displays f_1 axis horizontally so that f_1 traces can be displayed. trace='f2' displays f_2 axis horizontally so that f_2 traces can be displayed.

Display Scaling

The nm2d<(noise_mult)> macro sets up the parameters vs2d and th automatically for a 2D contour plot and color map display. nm2d measures the highest signal in the spectrum and sets vs2d such that the highest signal is in the range of the highest color level. It then measures the root-mean-square noise in the second trace, both in f_1 and f_2 (the f_1 trace may contain an f_1 axial signal if FAD was applied). From both traces, nm2d measures two regions, avoiding the center and both ends of the trace. From the four resulting root-mean-square noise figures, the lowest value is taken and multiplied with a noise multiplier, either the value specified by the argument noise_mult or the default (8 for ¹H, ¹⁹F and ³¹P [high-dynamic-range nuclei] and homonuclear spectra in general, or 3 for other spectra) if no argument is present. If the multiplied noise figure is below th=1, vs2d is scaled up; otherwise, th is increased to the desired level. nm2d works both with absolute-value and phase-sensitive spectra. trace can be set to 'f1' or 'f2'.

The macro noisemult<(noise_mult)> predetermines the noise multiplier used by nm2d when starting automatic 2D experiments. The default is 8 for homonuclear 2D spectra or 4 for other spectra. The argument noise_mult overrides the default.

Grid Lines

A grid of horizontal and vertical lines over a 2D display can be drawn by the grid macro. By default, grid lines are drawn in blue at approximately 1 cm intervals, rounded so that the intervals fall at a multiple of 1, 2, or 5 of Hz or ppm. To change the defaults, enter grid with a different spacing (in cm) or a different color ('red', 'green', etc.); for example, grid(2, 'white') gives white grid lines at 2 cm intervals.

The grid command also can define a grid, using the following syntax:

grid<(startf2,incrf2,startf1,incrf1,color)>

The arguments define the frequency and increments between grid lines in the f_2 and f_1 directions and the color of the grid lines.

The plgrid macro uses the same arguments as grid, but plots the grid instead.

Color Maps and Contour Plots

The main command for 2D color map display is dconi, but it is interactive and described separately on page 113. The rest of the 2D display commands are noninteractive.

- dcon produces a "contour plot" (actually a color intensity map) on the screen. dconn performs the same as dcon except dconn does not erase the screen before starting the display. Both dcon and dconn accept various options as arguments.
- image macro displays a dcon noninteractive display of an image using grayscale and linear scaling of the intensity. This macro is useful for adjusting the display while using the dconi command. The imageprint macro can send the same image to the plotter.
- dpcon(<mode,><levels, spacing>) produces a true contour plot display, and pcon<(mode,><levels, spacing)> produces a contour plot on the plotter.
 levels represents the maximum number of contour levels that will be shown; the default is 4. spacing represents the spacing of successive contour levels; the default is 2. dpconn(<mode,><levels, spacing>) produces the same display as dpcon but does not erase the screen before starting. In phase-sensitive spectra, mode can be 'pos', to display positive peaks only, or 'neg', to display negative peaks only.

df2d<(array_index)> produces a color intensity map of the raw 2D FIDs as a function of t₁ and t₂, where array_index is the index of the array to display. The display can be modified by subsequent display commands.

Whitewashed Spectra

The dsww<(start,finish,step)> command displays one or more spectra with whitewashing (traces in front "block" the view of traces behind them). Use the argument 'all' to display all spectra.plww<(start,finish,step)> plots the same spectra.

The ds2d<(options)> command produces a stacked display of 2D spectra in the whitewash mode. Certain options are available: 'nobase' activates the th parameter to suppress all intensity below th, 'fill' fills in the peaks, and 'fillnb' combines base suppression and peak filling. The ds2dn<(options)> command creates the same stacked display as ds2d but does not erase the screen before starting. A stacked plot is produced by pl2d<(options)>.

Equivalent nonwhitewashed stacked displays and plots may be obtained with dss and pl, respectively.

Label Display

The dssl macro displays a label for each element in a set of stacked spectra. The label is an integer value starting with 1 and ranging to the number of spectra in the display.

If wysiwyg='n', labels can appear at incorrect positions. The positions were empirically determined for a large screen display and are not guaranteed to be correct for all displays.

The following options control the dssl display (more than one option can be entered as long as the options do not conflict with each other):

- 'center', 'left', 'right', 'top', 'bottom', 'above', and 'below' are keywords setting the position of the displayed index relative to each spectrum.
- 'value' is a keyword that produces a display of the values of each array element, instead of an integer index.
- 'list=xxx' produces a display of the values contained in the arrayed parameter xxx.
- 'format=yyy' uses the format yyy to control the display of each label. See the write command for information about formats.

Projection of 2D Data

The proj(exp_number<, 'sum', start, width>) command projects data onto the axis parallel to the screen x-axis, which can be either the f_1 or f_2 axis, depending upon the parameter trace

- In a skyline projection, the data is searched and the maximum intensity at any given frequency becomes the intensity in the projection.
- In a summing projection, the data at each frequency are summed and the result becomes the projection. start and width (in Hz) define the traces to be projected.
- If omitted, the whole data set is projected.

The argument exp_number is the number of the experiment in which the resulting spectrum is stored. If the 'sum' keyword is given, the projection is calculated as the sum of the data; otherwise, as the maximum (skyline projection).

2D Referencing

The macros rl and rll set the direct dimension and first indirectly detected dimension reference lines, respectively. By default, each reference line is set at the cursor position (cr for direct dimension and crl for the first indirect dimension) after taking into account any frequency scaling with the scalesw or scaleswl parameters, respectively.

To set the reference lines to other than the cursor position, include the frequency argument. A number for frequency with no suffix (e.g., rll(0)) indicates the frequency in hertz. A number with the suffix p, d, or k (e.g., rll(7.2p)) indicates the frequency in ppm, decoupler ppm, or kilo, respectively. These suffixes are defined by the unit command. The default definition is multiplication by (respectively) reffrq, dfrq, and 1000. Thus, if you are doing a 2D experiment in which the indirect axis is determined by the decoupler channel (e.g., HMQC or HETCOR experiment), you might enter, for example, rll(10d), which is equivalent to rll(10*dfrq). The refsource1 parameter is used to signify whether the frequency along the indirect dimension is related to the observe transmitter (refsource1='sfrq') or if it is related to the decouple frequency (refsource1='dfrq').

The command crl sets the reference parameters rfl and rfp to zero and sets refpos1='n'. In 2D spectra, this clears the referencing along f_2 . In 2D spectra the command crll can be used to set rfll and rfpl to zero and refpos1='n', therefore clearing referencing along f_1 .

The macro setsw1(nucleus, downfieldppm, upfieldppm) sets the spectral width parameter sw1 for a given spectral window. setsw1 also does referencing.

The macros centersw and centersw1 set the cursor to the center of the spectrum in the directly detected dimension and the first indirect dimension, respectively.

The macros setref1 and setref2 set frequency referencing using the setref macro. Given a nucleus (e.g., tn or 'C13'), setref1(nucleus) calculates the value of rfl1, rfp1, refpos1, and reffrq1; and setref2(nucleus) calculates rfl2, rfp2, refpos2, and reffrq2.

Rotating Homonuclear 2D-J Spectra

The rotate<(angle)> command rotates homonuclear 2D-J data 45° (rotation in frequency-space) to line up multiplets. Use the angle argument to specify other angles.

Symmetrizing Data

The foldt, foldj, and foldcc commands symmetrize data as follows:

- foldt symmetrizes or triangularizes COSY, NOESY, or similar 2D spectra, by "folding" about the diagonal; it requires fn=fn2 and sw=sw1.
- foldj symmetrizes heteronuclear 2D-J or rotated homonuclear 2D-J experiments by "folding" along fl=0 (J=0) axis.
- foldcc symmetrizes 2D INADEQUATE data along the appropriate axis and also applies an automatic dc correction.

Setting Negative Intensities to Zero

The command zeroneg is used for the projection of proton 2D-J spectra at 45° to strip a high resolution proton spectrum down to a list of chemical shifts. zeroneg sets all negative intensities to zero.

Automatic Analysis

The acosy and acosyold commands automatically analyze a COSY data set with fn=fn1 and sw=sw1 (acosyold uses an older algorithm from previous VNMR software versions). Symmetrization of the data with the command foldt is recommended, but not required.

- Select a proper threshold and perform a 2D line listing with the command ll2d('peak', 'volume').
- 2. Display the 2D data with the command dcon, leaving enough room at the left side of the display for the connectivity table.
- 3. Enter a cosy or a cosyold to analyze the data and display the connectivities on the screen.

The command pacosy performs the same analysis as acosy and plots the connectivities on a plotter.

3.12 Interactive 2D Color Map Display

The dconi<(options)> command interactively displays 2D traces and projections, which can also be plotted. The following keywords can be used as optional arguments:

- 'again' means identify current screen mode and redraw the screen in that mode.
- 'avcolor' means use absolute-value color set and display positive peaks only.
- 'dpcon' means display true contour plot.
- 'ds2d' means display "whitewashed" stacked plot.
- 'gray' means use grayscale color set.
- 'linear' means display in linear instead of logarithmic increments.
- 'phcolor' means use phased color set, and display positive and negative peaks.
- 'restart' means activate dconi without redrawing the 2D data set (make sure that 2D data is already displayed).

When dconi in invoked with the 'dpcon' keyword, it draws the 2D data on the screen similar to Figure 23.

At the top of the display (directly below the Permanent menu) is the Main menu for dconi. Below the menu is a graphics window with a 2D display box and, if the 'dcon' option is selected, a color/grayscale adjustment bar to the right of the box. The default display mode is the color map drawn by dcon, but optional keywords 'dpcon' and 'ds2d' can be used to override the default.

Readouts for parameters such as the following appear in the lower part of the screen:

- cr shows the current cursor position.
- crl shows the current cursor position along the first indirectly detected dimension.
- delta shows the cursor difference.
- delta1 shows the cursor difference along the first indirectly detected dimension.
- vs2d shows the vertical scale of the display.
- vsproj shows the vertical scale of the trace or projection.

Table 27 lists commands and parameters associated with the dconi program.




Figure 23. Interactive 2D Contour Display (dconi Program)

<pre>Commands boxes<('graphics' 'plotter')> dconi<(options)> vsadj</pre>	Draw boxes selected by the mark command Interactive 2D contour display Automatic vertical scale adjustment
<pre>Parameters</pre>	Current cursor position
cr {number}	Cursor position along 1st indirectly detected dimension
cr1 {number}	Control display selection for the dconi program
dconi*	Difference of two frequency cursors
delta {pos. number, in Hz}	Cursor difference in 1st indirectly detected dimension
delta1 {pos. number, in Hz}	Gray level window adjustment
grayctr {0 to 64}	Gray level slope adjustment
grays1 {-10 to +10}	Start of plot in directly detected dimension
sp {number, in Hz}	Start of plot in 1st indirectly detected dimension
vs2d {1e-6 to 1e9}	Vertical scale for 2D displays
vsproj	Vertical scale for projections and traces
wp {number, in Hz}	Width of plot
wp1 {number, in Hz}	Width of plot in 1st indirectly detected dimension
* dconi {",'display_program',' display_pro	gram,option1,option2'}

Table 27. Interactive 2D Color Map Display Commands and Parameters

The string parameter dconi also controls the display selection. For example, if the parameter dconi is set equal to the string 'dpcon, pos, 12, 1.2', the dconi command displays twelve positive contours with dpcon, using a spacing of 1.2. For details on using the dconi parameter, refer to the *VNMR Command and Parameter Reference*.

Interactive 2D Display Menus

Upon starting the dconi program, the 2D Display Main Menu is active with the following buttons:

By clicking on the Proj button, the following 2D Display Projection Menu replaces the 2D Display Main Menu:

Hproj(max)	Hproj(sum)	Vproj(max)	Vproj(sum) Plot	Cancel
------------	------------	------------	-----------------	--------

Each of the buttons on these menus is described in the sections below.

Interactive 2D Color Map Display Main Menu

The buttons the Interactive 2D Color Map Display Main Menu function as follows:

Box	The first button is labeled Box or Cursor, depending on the display mode you are in. If labeled Box, you are in the cursor mode, and this button changes the display to he box mode with two pairs of cursors.
Cursor	If labeled Cursor , you are in the box mode, and this button changes the display to the cursor mode with one pair of cursors.
Trace	Selects the trace display mode.
Proj	Displays the Interactive 2D Display Projection Menu, see below.
Expand	The fourth button is labeled Expand or Full depending on the mode you are in. If labeled Expand, you are in the box mode and this button expands the area between the cursors.
Full	If labeled Full, you are in the cursor mode and this button displays the full area.
Redraw	Repeats the last 2D or image display with current parameters.
Plot	Plots the current trace.
Peak	Displays the Interactive 2D Peak Picking Main Menu—112d program (see page 121).
Return	Returns to the previous menu.

Interactive 2D Display Projection Menu

The buttons the 2D Display Projection Menu function as follows:

Hproj (max)	Displays a horizontal projection of the maximum intensity at each frequency.
Hproj (sum)	Displays a horizontal projection of the summed intensity at each frequency.
Vproj (max)	Displays a vertical projection of the maximum intensity at each frequency.
Vproj (sum)	Displays a vertical projection of the summed intensity at each frequency.
Plot	Plots the current projection.
Cancel	Returns to the Interactive 2D Color Map Display Main Menu, see above.

Controlling the Display with the Mouse

The left and right mouse buttons are used to move cursors, the center button to adjust the vertical scale of traces, projections and contour maps, as well as to adjust the threshold in the color bar. The cursors can be used to select regions for expansions of the display. The cursors can also be used to select positions to "mark" using either the mark command or the ll2d('mark') command. Both commands display and record spectral frequencies, maxima, intensities, and volumes. ll2d('mark') is recommended, however, because it also allows for interactive display and editing of mark locations.

Left Mouse Button

The left mouse button adjusts the position of the 2D cursor. The corresponding frequencies are displayed at the bottom of the graphics window. Both the horizontal and vertical cursors move if the left mouse button is pressed within the 2D display box.

Above and below the box, only the vertical cursor can be moved; at the left and the right of the box, only the horizontal cursor. In addition, holding the mouse button down and then moving the mouse moves the cursor with the mouse. This "dragging" mode is not available on the GraphOn terminal for speed reasons.

Center Mouse Button

The function of the center mouse button depends on the location of the cursor:

- If the cursor is within the 2D display box, in gray scale images, pressing the center button sets the point to medium gray. Otherwise, for color map and contour displays, if there is no intensity at that point, the center button changes vertical scale to show intensity at that point. If there is intensity at the point, the center button changes the scale to show no intensity, then changes the parameter vs and redraws.
- If the cursor is near an active trace and active horizontal or vertical projection, pressing the center button changes the vertical scale of trace or projection, so that spectrum goes through the current mouse position.
- If the cursor is near the color/grayscale bar and in the color mode, pressing the center button sets the threshold to remove low intensity peaks. If in the grayscale mode, pressing the center button sets the grayscale intensity (the right button adjusts contrast).

Right Mouse Button

A second cursor pair is displayed with the right mouse button. The second pair can be moved in exactly the same way as the first pair, and is used to select a box within the 2D display. The right mouse button also switches the display into the box mode, the same as clicking on the Box button in the menu.

Changing the Display

The user interactively modifies the display through selecting buttons on the menus, pressing the buttons on the mouse, and moving the mouse. If desired, commands and macros can also be typed in at any time.

Displaying Traces

The Trace button can be used to display a trace at the current position of the first horizontal cursor. The left mouse button still moves the cursor, and different traces will be selected and displayed accordingly.

The vertical scale of the trace can be adjusted with the center mouse button, by placing the mouse arrow above the spectrum at the requested height and then pressing this button. Do not press the center button at this time within the 2D display box.

The trace mode can be left by displaying a box with the right mouse button, or by selecting any other display mode.

Displaying Projections

With the Proj button, you can display projections from the 2D Display Projection Menu. Projections can be made in horizontal and vertical direction, and are available as the sum or maximum of the data. Select one of the four available modes or cancel the operation with the Cancel button. Again, the center mouse button adjusts vertical scales.

Expanding the Display

Once a box is selected, you can click on the Expand button to obtain an expanded display. Alternatively, if you are in the one cursor mode, you can click on the Full button to display the full 2D display. In this case, the two cursor pairs mark the last expanded region and the Full/Expand button toggles between the full and expanded mode.

An alternative way to select an expansion is to type in new values for the parameters sp, wp, sp1, and wp1 (e.g., sp=100 wp=50 sp1=100 wp1=50), then use the Redraw button to redisplay.

Setting the Vertical Scale

The vertical scale can be adjusted in a number of ways. If a peak is expected at a certain position in the spectrum but is not visible, the mouse arrow can be moved to that position and then the *center* mouse button pressed once. This selects a new vertical scale, so that the intensity at that point is by a factor of 2 above the threshold, and the display is redrawn. Be careful in this mode not to queue up several redraw operations.

Adjusting the Threshold

If noise is visible at a certain position in a spectrum, but should be suppressed below the threshold, move the mouse arrow to that position and press the center button. A vertical scale is calculated so that this intensity falls by a factor of 2 below the threshold, and again the spectrum is redrawn. If the peak is visible but is not a factor of 2 above the threshold, clicking on the center button increases vs2d.

On Sun color screens only, the threshold of the 2D display can be adjusted in real time. For color displays, the threshold is adjusted by placing the mouse button on the color bar at the right edge of the display, selecting one of the colors and pressing the center mouse button. All colors below that level are set to black. For grayscale images on the Sun color display, the center of the grayscale is adjusted in the same way. At the same time, the grayscale is expanded by a factor of 4. On monochrome terminals, this mode is not available.

In order to perform threshold adjustment on grayscale images, two new parameters must be created: grayctr, which controls the center of the grayscale, and graysl, which

controls the slope of the grayscale. Enter the parim macro to create grayctr and grays1, or create the parameters by hand as follows:

- To create grayctr, enter create('grayctr','real') setgroup('grayctr','display') setlimit('grayctr',64,0,1).
- To create graysl, enter create('graysl','real') setgroup('graysl','display') setlimit('graysl',10,-10,0.01).

If these parameters do not exist, the interactive display still lets you adjust the grayscale threshold and contrast, but these adjustments are not retained.

Treating 2D Traces as 1D Spectra

After a trace has been selected in the interactive 2D display program, entering the command ds allows the trace to be displayed as if it were a simple 1D spectrum. All standard 1D data manipulations, including line listing, integration, etc., are then accessible for that trace. The command ds(tracenumber) also can be used to display an f_1 or f_2 trace, depending on the value of trace.

3.13 Interactive 2D Peak Picking

The ll2d program is used to automatically or interactively pick peaks in 2D spectra or 2D planes of 3D spectra. The peaks can be displayed on top of the spectrum in the dconi display or can be plotted using the pll2d command. Table 28 lists commands and parameters related to the ll2d program.

Commands	
addpar('ll2d')	Add ll2d parameters to the current experiment
dconi<(options)>	Interactive 2D contour display
112d*	Automatic and interactive 2D peak picking
ll2dbackup<(file)>	Copy current ll2d peak file to another file
parll2d	Create parameters for 2D peak picking
pll2d<(options)>	Plot results of 2D peak picking
* ll2d<(options)><:\$nu	m> ,
ll2d('info'<,#>):\$pe	<pre>ak_number,\$f1,\$f2,\$amplitude,\$volume,\$label,</pre>
\$comment,\$FWHH1,\$F	WHH2,\$f1_min,\$f1_max,\$f2_min,\$f2_max
Parameters	
ins {number}	2D volume value
<pre>ins2ref {number}</pre>	Fourier number scaled volume of a peak
ll2dmode*	Control display of peaks picked by 112d program
th2d {0.0 to 1.0]	Threshold for integrating peaks in 2D spectra
xdiag {number, in Hz}	Threshold for excluding diagonal peaks when peak picking
* 112dmode {4 characters from	'y' and 'n'}

Table 28. Interactive 2D Peak Picking Commands and Parameters

The results of all peak picking operations are stored in a binary file in the ll2d subdirectory of the current experiment directory:

• For 2D spectra, the results are stored in the file peaks.bin.

• For 2D planes of 3D spectra, the results are stored in peaks_f#f#_#.bin, where f#f# denotes the orientation of the plane being picked (e.g., flf3 or f2f3) and the last # denotes the number of the plane.

Binary peak files can be converted to text files for printing or for export to other programs.

For each peak in a peak file, the following information is stored:

- Peak number
- · Interpolated peak frequency in both dimensions
- Interpolated peak amplitude
- Full width at half-height (FWHH) in both dimensions
- Bounds of the peak in both dimension
- Volume of the peak
- 15-character peak label
- 80-character comment

The parameter ins2 adjusts the 2D volume value. Volume is independent of is and vs2d. It is scaled by Fourier numbers for the two dimensions.

The parameter ins2ref is set to the Fourier number scaled volume of the selected peak. The reported volume is *volume**ins2/ins2ref/fn/fn1. If ins2ref is "not used," the sum of all volumes is ins2. The "not used" mode is equivalent to a "normalized" volume mode. If ins2ref is zero or not defined, the reported volumes will be *volume**ns2/fn/fn1.

A typical use of ins2ref would be to position a cursor within a peak region and set ins2ref equal to the scaled volume returned by the command ll2('info'). The reported volume of that peak would then be the value of ins2. This operation is analogous to the 1D integral scheme.

- vs2d shows the vertical scale of the display.
- vsproj shows the vertical scale of the projection or trace.

The options listed below are available for 112d:

- 'peak' is a keyword to find all peaks above the current threshold in the area of the spectrum displayed in dconi (if in cursor mode) or the area defined by the cursors (in box mode). This option gives each peak a number and determines peak frequencies and amplitude.
- 'volume' is a keyword that, for all peaks picked using the 'peak' option, finds the bounds, volume, and FWHH) of the peak in both dimensions.
- 'adjust' is a keyword to adjust the bounds of all peaks in the displayed area so that none overlap, and then to recalculate volumes.
- 'reset' is a keyword to delete all peaks in the spectrum.
- 'read' is a keyword to prompt for a binary peak file name and read in that file.
- 'read', file reads in a binary peak file named file.
- 'readtext' is a keyword to prompt for a text peak file name and read in that file.
- 'readtext', file reads in a text peak file named file.
- 'writetext' is a keyword to prompt for a file name for a peak file and write out a text file with that filename.
- 'writetext', file writes a peak file to a text file with the name given by file.
- 'draw' is a keyword to draw peaks in the peak file to the graphics window.

- 'mark' is a keyword to insert a peak at the current cursor location (in the dconi cursor mode) or to use the area defined by the cursors as peak bounds and calculate the volume in this area (in the dconi box mode). This option assigns these bounds to each peak within this area that does not have its bounds already defined. If a peak without bounds defined does not exist in this area, it finds the highest point in this area, marks it as a peak, and assigns it the bounds defined by the cursors.
- 'unmark' is a keyword to delete the peak nearest the cursor (in dconi cursor mode) or to delete all peak bounds which are completely within the area defined by the cursors (in box mode). The peaks are not deleted in box mode.
- 'unmark', # deletes peak number #.
- 'clear' is a keyword to delete all peaks within the displayed area (in dconi cursor mode) or to delete all peaks within the area defined by the cursors (in box mode).
- 'label' is a keyword to prompt for a 15-character label. The label is assigned to the nearest peak (dconi cursor mode) or to all peaks within the area defined by the cursors (dconi box mode).
- 'label', string executes the 'label' option using the string argument instead of prompting for a label.
- 'label', string, # assigns string to be the label of peak number #.
- 'comment' is a keyword to prompt for a 80-character comment. The comment will be assigned to the nearest peak (cursor mode) or to all peaks within the area defined by the cursors (box mode).
- 'comment', string executes the 'comment' option using the string argument instead of prompting for a comment.
- 'comment', string, # assigns string to be the comment of peak number #.
- 'info' is a keyword to print information to the text window about the peak nearest the cursor.
- 'info', 'total' prints the total number of peaks in the spectrum, or if a return value is requested, returns the total number of peaks in the spectrum.
- 'info', # prints information to the text window about peak number #. If return values are requested, then printing is suppressed and the values are returned in this order: peak_number, f1, f2, amplitude, volume, label, comment, FWHH1, FWHH2, f1_min, f1_max, f2_min, f2_max.
- 'combine' is a keyword to combine all peaks within the area defined by the cursors into a single peak (in dconi box mode only). The individual peaks to be combined are permanently deleted. You may wish to back up the peak file before using this option, since it is not possible to undo combining peaks.
- 'combine', #1, #2,... performs the 'combine' option on the list of peaks with numbers #1, #2,.... If a return value is requested, the value returned is the peak number of the new combination peak.
- 'pos' or 'neg' can be used in addition to 'peak', 'volume', or 'clear' to operate only on either positive or negative peaks.

Display of peaks in dconi is controlled by the global parameter ll2dmode. This parameter has four characters, each of which can take the value 'y' or 'n'. The first character controls display of a "+" to mark the peak maximum, the second controls display of the peak number, the third controls display of the peak bounds, and the fourth controls display of the peak label.

The parameter th2d controls the threshold for integrating peaks, and the parameter xdiag excludes diagonal peaks within xdiag Hz of the diagonal from peak picking. If

these parameters do not exist, a default value is used for each. To use different values than the defaults, create both parameters by entering addpar('ll2d') and then setting the values as usual (the macro parll2d functions the same as addpar('ll2d')).

Interactive 2D Peak Picking Menus

Most of the above options are accessible through a series of user-programmable menus of the dconi program (described in "Interactive 2D Color Map Display," page 113). From the 2D Display Main Menu of dconi, the Peak button brings up the 2D Peak Picking Main Menu of the 112d program with the following buttons:

Auto Edit File	Display Return
----------------	----------------

These buttons provide access to the following menus (the labels on some buttons change depending on what mode you are in):

• 2D Peak Picking Automatic Menu for automatically picking peaks, selected by the Auto button in the 2D Peak Picking Main Menu (Figure 24 shows this menu with a typical 112d screen):

Box Peak Volume Full Both Adju	ıst Reset Return
--------------------------------	------------------

Read Read Text Write Text Backup File Return

• 2D Peak Picking Edit Menu for interactively editing peaks, selected by the Edit button in the 2D Peak Picking Main Menu:

```
Box Mark Unmark Full Clear Combine Label Comment Info Set Int Return
```

• 2D Peak Picking File Menu for manipulating peak files, selected by the File button in the 2D Peak Picking Main Menu:

Box) Peak) Volume)	11) Both] Adjust] Reset] Return)	
F C		
	$\begin{array}{cccccccccccccccccccccccccccccccccccc$	
cr(ppm) de <mark>9.17 </mark> 6.	8.5 8.0 7.5 7.0 6.5 6.0 5.5 5.0 F1 (pom) a(ppm) cr1(ppm) delta1(ppm) specvs vs 8.54 3.89 44.4 376.5	

Figure 24. Interactive 2D Peak Picking (ll2d Program)

• 2D Peak Picking Display Menu for controlling the display of peaks, selected by the Display button in the 2D Peak Picking Main Menu.

Hd Pk Hd Num Hd Box Hd Lbl Sh All Hd All Return

Each of the buttons on these menus is described in the sections below.

2D Peak Picking Main Menu

This menu selects another 2D peak picking menu. The buttons function as follows:

2D Peak Picking Automatic Menu

Auto	Displays the 2D Peak Picking Automatic Menu (see below)
Edit	Displays the 2D Peak Picking Editing Menu (see below)
File	Displays the 2D Peak Picking File Menu (see below)
Display	Displays the 2D Peak Picking Display Menu (see below)
Return	Displays the 2D Display Main Menu (see page 115).

This menu provides automatic peak picking. The buttons functions as follows:

Box	The first button is labeled Box or Cursor, depending on the dconi display mode you are in. If labeled Box, you are in the dconi cursor mode, and this button changes the display to the dconi box mode with two pairs of cursors.
Cursor	If labeled Cursor, you are in the dconi box mode, and this button changes the display to the dconi cursor mode with one pair of cursors.
Peak	Automatically finds peaks in the 2D spectrum. If one cursor is visible (dconi cursor mode), all peaks above the current threshold in the currently displayed region of the spectrum are found and marked. A peak is defined as a data point that is higher than the eight points around it. Once such a point is found, the actual peak location is determined by interpolation in both dimensions.
Volume	Automatically finds the bounds of a peak and the integral of all points within these bounds. The bounds are found by descending down the sides of a peak until the point is reached where the amplitude of a data point is less than th2d times the current threshold. Thus, using a smaller value for th2d will cause 112d to find and integrate a larger area for the bounds of the peaks. The peak volume is calculated by summation of all data points within these bounds. If the bounds of a peak already exist, the volume is recalculated.
Expand	The fourth button is labeled Expand or Full depending on the mode you are in. If labeled Expand, you are in the box mode and this button expands the area between the cursors.
Full	If labeled Full, you are in the cursor mode and this button displays the full area.
Both	Pick peaks and calculate volumes. The Both button does both the peak and volume operations at once.

Adjust	Adjust peak bounds so that none overlap. The Adjust button adjusts
	all peak bounds in the displayed region of the spectrum so that none
	overlap and recalculates peak volumes with the new peak bounds.
Reset	Deletes all peaks that have been found in the current spectrum.
Return	Display the 2D Peak Picking Main Menu (see above).

2D Peak Picking Edit Menu

This menu provides interactive peak editing. The buttons functions as follows:

Box	The first button is labeled Box or Cursor, depending on the display mode you are in. If labeled Box, you are in the cursor mode, and this button changes the display to he box mode with two pairs of cursors.
Cursor	If labeled Cursor , you are in the box mode, and this button changes the display to the cursor mode with one pair of cursors.
Mark	In dconi cursor mode, this button inserts a peak at the current cursor location. In dconi box mode, the cursors are taken as peak bounds and the area inside the cursors is integrated. These peak bounds are then assigned to all peaks within the cursors that do not already have their bounds defined. If a peak without bounds does not exist inside the area defined by the cursors, the highest point within that area is found, marked as a peak, and assigned the bounds defined by the cursors.
Unmark	In dconi cursor mode, this button deletes the peak nearest the cursor. In dconi box mode, this button deletes peak bounds from peaks whose bounds are entirely within the area defined by the cursors.
Expand	The fourth button is labeled Expand or Full depending on the mode you are in. If labeled Expand, you are in the box mode and this button expands the area between the cursors.
Full	If labeled Full, you are in the cursor mode and this button displays the full area.
Clear	In dconi cursor mode, this button deletes all peaks in the area of the spectrum displayed in dconi. In dconi box mode, the Clear button deletes all peaks that are within the area defined by the cursors.
Combine	This button works only in dconi box mode. It combines all peaks within the area defined by the cursors into a single peak. This combination peak is located at the average frequencies of all of the original peaks and has bounds that encompass all of the original bounds of the peaks. The volume of the combination peak is calculated by summation of all data points within its bounds. You may wish to back up the peak file using the Backup File button in the 2D Peak Picking File Menu (see below) prior to using this button, because the original peaks are permanently deleted when the combination peak is created.
Label	Prompts for a 15-character label to be assigned to the peak nearest the cursor (dconi cursor mode) or to all peaks within the area defined by the cursors (dconi box mode). Based on the value of the parameter ll2dmode, this label can be displayed next to the peak in dconi.

Comment	Prompts for an 80-character string to be assigned to the peak nearest the cursor (cursor mode) or to all peaks within the area defined by the cursors (box mode).
Info	Prints the peak file information about the peak nearest the cursor to the text window.
Set Int	Set the value of the peak volume.
Return	Display the 2D Peak Picking Main Menu (see above).

2D Peak Picking File Menu

Read	Prompts for the filename of a binary peak file and reads that file into
	VNMR. When a file is read in, the current peak file (peaks.bin
	for 2D spectra) is overwritten by a copy of the peak file that was read
	in.

- Read Text Prompts for the file name of a text peak file and reads that file into VNMR. When a file is read in, the current peak file (peaks.bin for 2D spectra) is overwritten by a new binary copy of the peak file that was read in.
- Write Text Prompts for a filename to write a text version of the current ll2d peak file.
- Backup Prompts for a file name to copy the current binary peak file. It is a File good idea to do this occasionally when doing a significant amount of interactive peak editing, so that intermediate versions of the peak file can be recovered in the event of an error (such as inadvertently selecting the Clear or Reset button or making a mistake using the Combine button).
- Return Display the 2D Peak Picking Main Menu (see above).

2D Peak Picking Display Menu

Sh Pk	The first button is labeled either Sh Pk (for "show peak") or Hd Pk (for "hide peak"), which is used to select whether or not dconi automatically displays a "+" to mark each peak. If labeled Sh Pk, the "+" is hidden and this button shows a "+" at the location of each peak.
Hd Pk	If labeled Hd Pk, the "+" is shown and this button hides the "+" at the location of each peak.
Sh Num	The second button is labeled either Sh Num (for "show number") or Hd Num (for "hide number"), which is used to select whether or not dconi automatically displays the peak number next to each peak. If labeled Sh Num, the peak numbers are now hidden and this button shows a peak number next to each peak.
Hd Num	If labeled Hd Num, the peak numbers are now shown and this button hides the peak numbers.
Sh Box	The third button is labeled either Sh Box (for "show box") or Hd Box (for "hide box"), which is used to select whether or not dconi automatically displays the peak bounds of each peak. If labeled Sh Box, the box is now hidden and this button shows a box with the area integrated to get the volume of the peak.

- Hd Box If labeled Hd Box, the box is now shown and this button hides this box.
- Sh Lbl The fourth button is labeled either Sh Lbl (for "show label") or Hd Lbl (for "hide label"), which is used to select whether or not dconi automatically displays the peak label next to each peak. If labeled Sh Lbl, the peak labels are now hidden and this button shows a peak label next to each peak.
- Hd Lbl If labeled Hd Lbl, the peak labels are now shown and this button hides peak labels.
- Sh All The fifth button is labeled either Sh All (for "show all") or Hd All (for "hide all"), which is used to make dconi automatically display a "+", the peak number, the peak bounds, and the peak label for each peak. If labeled Sh All, the "+", number, box, and label are now hidden and this button shows "+", number, box, and label for all peaks.
- Hd AllIf labeled Hd All, the "+", number, box, and label are now shown and
this button hides the "+", number, box, and label for all peaks. The
Hd All button is used to make dconi display no peak information.
- Return Display the 2D Peak Picking Main Menu (see above).

Automatic 2D Peak Picking

This section and the next few sections describe techniques for using the 112d program.

Once a 2D spectrum has been Fourier transformed and the threshold has been adjusted to a suitable level in dconi, select the Peak button from the 2D Display Main Menu in dconi to display the 2D Peak Picking Main Menu. Then select the Auto button in that menu to display the 2D Peak Picking Automatic Menu.

If there are not many unwanted peaks visible above the threshold, then with one cursor showing (dconi cursor mode) selecting the Peak button will pick all peaks in the displayed region of the spectrum and write the frequencies and amplitude of each peak to the current peak file. If there are many artifacts or unwanted peaks in the spectrum, it is often better to either expand the display to show only regions of interest or select dconi box mode (two cursors showing) to enclose regions of interest with the cursors before selecting the Peak button. In dconi box mode, only peaks in the region enclosed by the cursors are picked.

As a general rule, it is more efficient to be selective when picking peaks and only pick peaks that you want, rather than picking all peaks and artifacts in a spectrum and later deleting the unwanted peaks.

Interactive Peak Picking or Editing

If unwanted peaks have been picked while using the 2D Peak Picking Automatic Menu, return to the 2D Peak Picking Main Menu and select the Edit button to display the 2D Peak Picking Edit Menu.

Individual peaks can be deleted by placing the cursor near the peak to be deleted and selecting the Unmark button. Groups of peaks can be deleted by entering the dconi box mode and enclosing the peaks to be deleted with the cursors prior to selecting the Clear button. The Clear button can also be used to delete all peaks in the currently displayed region of the spectrum if only a single cursor is visible (dconi cursor mode). A new peak can be interactively marked by placing the cursor in dconi cursor mode on a peak and

selecting the Mark button. In this case the peak location is set to be exactly at the cursor location. Subsequent automatic peak picking on the same peak will not necessarily coincide with peaks marked in this way.

Automatic Integration

Once peaks of interest have been picked, the peak bounds, full width at half-height (FWHH), and volume can be automatically determined using the Volume button in the 2D Peak Picking Automatic Menu. The volume is calculated by summing all data points within the bounds of the peak. The size of a peak's bounds that 112d finds is determined by the current threshold and the value of the parameter th2d. If th2d is set to its maximum value of 1.0, the peak bounds will be selected so that they contain just the portion of the peak which is visible in dconi above the current threshold.

To make 112d integrate a larger area, the value of th2d must be decreased. A value of th2d of 0.5, for example, will cause 112d to find peak bounds that contain all portions of a peak greater than 0.5 times the current threshold in amplitude. Selecting too small of a value for th2d can cause extensive overlap of bounds of neighboring peaks. The Adjust button will attempt to adjust the bounds of all peaks within the displayed region of the spectrum, so that none overlap, and then recalculate peak volumes. In cases of extensive overlap of several peaks, however, this function may not always adjust the peak bounds in an optimal way.

Interactive Integration and Editing

The bounds of peaks can be interactively added or deleted using buttons in the 2D Peak Picking Edit Menu. To delete peak bounds, enclose the bounds of the peak(s) to be deleted entirely with the cursors in dconi box mode and select the Unmark button. This operation deletes the bounds, but not the peaks. To insert peak bounds on a peak without bounds, place the cursors in dconi box mode so they enclose the peak and the area to be integrated for that peak and select the Mark button. This mode can also be used to mark and integrate peaks. If the Mark button is selected when the two cursors in dconi box mode enclose a region where no peaks have yet been picked, it will find the highest point within the area defined by the cursors, interpolate, mark this as a peak, integrate the area inside the cursors, and assign these bounds and volume to the peak.

Peak information from the peak file can be written to the text window using the Info button in the 2D Peak Picking Edit Menu. This button writes out the peak table entry for the peak nearest the cursor in dconi. Adjust the value of the peak by using the Set Int button in the following manner:

- 1. Position the cursor over the peak.
- 2. Press the **Set Int** button.

The program displays the prompt:

Current integral is xx. New value?

3. Type in the value you want to assign to that peak.

Labeling and Commenting Peaks

The Label and Comment buttons in the 2D Peak Picking Edit Secondary Menu can be used to store additional information with a peak. The Label button prompts for a 15-character label to be assigned to the peak nearest the cursor (dconi cursor mode) or all peaks within the area defined by the cursors (dconi box mode). Assigning peak labels or comments to

a group of peaks in dconi box mode has no special significance in the ll2d program the label or comment is simply stored in the peak file record of each peak in the group. The label can be displayed next to the peak in dconi based on the value of the parameter ll2dmode. An 80-character comment can be assigned to a peak or group of peaks using the Comment button. This comment is stored in the peak file with the peak entry and can contain any information desired.

Displaying Peaks in dconi

Peak information— a "+" to mark the peak top, the peak number, a box to mark the peak bounds, and the peak label—will be automatically displayed in dconi based on the value of the parameter 112dmode, as described above. The 2D Peak Picking Display Menu provides buttons to set ll2dmode for easy selection of which peak attributes are displayed in dconi or plotted with pll2d.

Peak File Manipulations

As described above, the default binary peak file for an experiment is stored in the 112d subdirectory of the current experiment directory. The 2D Peak Picking File Menu provides buttons to read, write and backup peak files. The Backup button prompts for a file name to which the current default 112d peak file should be copied. Unless a full UNIX (starting with a "/") is specified, the file is copied to the specified filename in the current working directory. If a full UNIX path is specified, the peak file is copied to the specified file in the specified directory. If no file name is specified at the prompt, the peak file is copied to the default peak filename with .bck appended (peaks.bin.bck for 2D peak files) in the current working directory.

The other buttons in this menu are used to read and write binary and text peak files. Each prompts for the name of a file, which is searched for by the following rules. If a full UNIX path is given, the specified file is read or written, otherwise a "read" searches for the file first in the current working directory and then in the 112d subdirectory of the current experiment directory, while a "write" writes the file to the current working directory.

3.14 3D NMR

VNMR includes full support for 3D NMR, including acquisition, processing, and display. Table 29 and Table 30 list commands and parameters connected with 3D NMR.

Many of the 3D-related macros and parameters—for example, centersw2, cr2, cr12, delta2, dmg2, lp2, lsfid2, phfid2, rfl2, rfp2, rp2, sp2, wp2—are normally used in the same manner as their 1D and 2D counterparts and are not described further in this section.

In a non-arrayed 3D experiment, there are two implicitly arrayed parameters: d2 and d3. d2 is associated with ni and sw1, d3 with ni2 and sw2. The order of these two arrayed parameters is such that d2 is cycled the most rapidly.

In an arrayed 3D experiment, such as a single 3D with "superhypercomplex" data acquisition (States-Haberkorn method applied along both t_1 and t_2), there are, in reality, at least three arrayed elements. By convention, such an arrayed 3D experiment is implemented using four arrayed elements: d3 (t_1 evolution time), phase2, d2 (t_2 evolution time), and phase.

Commands		
addpar('3d')	Add selected 3D parameters to the current experiment	
av2	Select AV mode on 2nd indirectly detected dimension.	
centersw2	Move cursor to center of spectrum in 2nd indirect dimension	
crl2	Clear reference line in 2nd indirectly detected dimension	
dplane*	Display a 3D plane	
dproj<(plane_type)>	Display a 3D plane projection	
dsplanes*	Display a series of 3D planes	
ft3d*	Perform a 3D FT on 3D FID data set (VNMR, UNIX)	
getplane*	Extract planes from a 3D spectral data set	
killft3d(exp_number)	Terminate any ft3d process started in an experiment	
nextpl	Display the next 3D plane	
par3d	Create 3D acquisition, processing, display parameters	
ph2	Select phased mode on 2nd indirectly detected dimension	
plplanes*	Plot a series of 3D planes	
prevpl	Display the previous 3D plane	
pwr2	Select power mode on 2nd indirectly detected dimension	
resetf3	Reset parameters after a partial 3D Fourier transform	
rl2<(frequency)>	Set reference line in 2nd indirectly detected dimension	
set3dproc*	Set 3D processing	
setsw2*	Setspectral width in 2nd indirectly detected dimension	
wft*	Weight and Fourier transform 1D data	
wftld*	Weight and Fourier transform f ₂ for 2D data	
wft2d*	Weight and Fourier transform 2D data	
wftt3	Process f3 dimension during 3D acquisition	
<pre>wti<(element_number)></pre>	Interactive weighting	
<pre>* dplane(<plane_type,>plane_number)</plane_type,></pre>		
dsplanes(start_plane		
ft3d<(<data_dir><,number_files><,'nocoef'><,plane_type>)> (VNMR)</data_dir>		
	f -r <options> (UNIX)</options>	
getplane<(data_directory><,place_directory><,plane_type>)>		
<pre>plplanes(start,stop<,'pos' 'neg'><,number_levels><,spacin g>) </pre>		
<pre>set3dproc<(<'nocef'><,direectory>)> setsw2(nucleus,downfieldppm,upfieldppm):offset</pre>		
<pre>wft<(<options,><'nf'><,start><,finish><,step>)>,</options,></pre>		
wft('inverse',exp_number,expansion_factor)		
wftld(element_number), wftld<(<options,><coefficients>)></coefficients></options,>		
wft2d<(<options,><coefficeints)></coefficeints)></options,>		
(Continued on next page)		

 Table 29. 3D NMR Commands and Parameters (Part 1 of 2)

Assuming that array='phase, phase2' (see below), the order of arrays is such that the phase2 array is cycled the most rapidly, followed by the phase, d2, and d3 arrays.

3D Acquisition

3D data acquisition is accomplished with pulse sequences using the parameter d3, which is incremented according to the parameters ni2 and sw2. This is analogous to d2, which is incremented according to ni and sw1 for 2D NMR (of course, d2, ni, and sw1 are active in 3D as well). In addition, the parameter phase2 is used to control the "mode" of acquisition (hypercomplex, TPPI, or absolute value) in the third frequency domain, just like phase in the second domain. All of these 3D parameters are created with the macro

Parameter	
cr2 {number}	Cursor position along 2nd indirectly detected dimension
d3 {number, in sec}	Incremented delay for 2nd indirectly detected dimension
delta2 {pos. number, in Hz}	Cursor difference in 2nd indirectly detected dimension
dmg2 {'ph2','av2','pwr2'}	Display mode along 2nd indirectly detected dimension
flcoef	Coefficient to construct F1 interferogram
f2coef	Coefficient to construct F2 interferogram
fiddc3d {3-char string}	3D time-domain dc correction
index2 {0, 1 to fn/2}	Projection or 3D plane index selected
lp2 {-3600 to +3600, in deg.}	First-order phase in 2nd indirectly detected dimension
lsfid2 {'n',number}	Number of complex points to left-shift ni2 interferogram
lsfrq2 {number,in Hz}	Frequency shift of the fn2 spectrum, in Hz
ni {number}	Increments in 1st indirectly detected dimension
ni2 {number}	Increments in 2nd indirectly detected dimension
path3d {path}	Path to currently displayed 2D planes from a 3D data set
phfid2 {number}	Phase selection for 3D acquisition
phfid2*	Zero-order phasing constant for ni2 interferogram
plane*	Currently displayed 3D plane type
ptspec3d {3-char string}	Region-selective 3D processing
rfl2 {number,in Hz}	Reference peak position in 2nd indirectly detected dimension
rfp2 {number,in Hz}	Reference peak frequency in 2nd indirectly detected dimension
rfp2 {-360.0 to +360.0, in deg.}	Zero-order phase in 2nd indirectly detected dimension
sp2 {number, in Hz}	Start of plot in 2nd indirectly detected dimension
ptspec3d {3-char string}	Flag for 3D spectral dc correction
scalesw2 {'n', $n > 0$ }	Scale spectral width in 2nd indirectly detected dimension
sw2 {number, in Hz}	Spectral width in 2nd indirectly detected dimension
trace {'f1','f2','f3'}	Mode for <i>n</i> -dimensional data display
wp2 {number, in Hz}	Width of plot in 2nd indirectly detected dimension
* plane {'f1f3','f3f1','f2f3','f3f2','f1f2	
phfid2 {'n',-360.0 to +360.0, in deg	}

Table 30. 3	D NMR	Commands and Parame	ters (Part 2 of 2)
-------------	-------	---------------------	--------------------

addpar('3d') along with other 3D parameters, including fiddc3d for 3D timedomain dc correction, ptspec3d for region-selective 3D processing, and path3d for the path to the currently displayed 2D planes extracted from a 3D data set. (The macro par3d is functionally equivalent to addpar('3d').)

By convention, 3D sequences are described with the first evolution time being known as t_1 , the second evolution time as t_2 , and the time during which data are acquired as t_3 . After transformation, these same dimensions are called the f_1 , f_2 , and f_3 dimensions.

3D Processing

Data processing includes the ft3d command for full 3D processing, governed by the usual parameters to control transform sizes, weighting, phasing, etc., with a "2" at the end of the parameter name signifying the third dimension. Unlike other VNMR commands, ft3d occurs in the background by default; that is, it is run as a separate task by UNIX, leaving VNMR free to continue with other tasks (including 1D and 2D processing of the same data set!). To increase the speed of 3D transforms further, the wftt3 macro allows the software to process one dimension (the acquisition or t₃ dimension) as the data are being acquired. Also, the ft3d software can be configured to run on several computers simultaneously, for

even greater speeds. The killft3d macro terminates any ft3d program that has been started in an experiment.

3D Display

Once the data are processed, the data can be displayed as two-dimensional planes of the 3D data set in any of the three orthogonal directions. Skew planes are not supported, nor are "full 3-dimensional" displays. One command, getplane, extracts the 2D planes from the 3D data set in one or more of the three orientations. After the planes are "extracted" in this manner, they are displayed with the dplane macro. The parameter index2 keeps track of which plane is on display. The macro nextpl displays the next plane from the plane currently on view. Another macro, prevpl, shows the previous plane from the current plane.

The dsplanes (start_plane, stop_plane) macro produces a graphical 2D color or contour map for a subset of 3D planes specified by the arguments. The dconi program is used to display the planes. The plplanes macro is available to plot a series of 3D planes.

The new concept of time-domain frequency shifting can be employed to good use in 3D NMR, where spectra in the indirectly detected directions are often "folded" by accident or by choice. The parameters <code>lsfrq.lsfrql</code>, and <code>lsfrq2</code> cause the frequency of the spectrum to be shifted as part of the Fourier transformation process.

3D Pulse Sequences

No standard and fully documented pulse sequences are provided for 3D NMR in the released software, although a number of sequences will be found in the user library. If you are writing your own sequences, you simply need to write a sequence that includes a d2 and d3 delay (these delays may also be d2/2 or d3/2). If your sequence is to operate in the hypercomplex (or the hyper-hypercomplex) mode, you should use the parameters phase and phase2 to select between the two orthogonal components of the hypercomplex experiment in the relevant domain. To ensure that your experiment is processed correctly using the default processing coefficients, you should write your pulse sequence so that the phase=2 (and phase2=2) experiments leave the receiver unchanged (compared to phase=1) and either increment the phase of the pulse (or pulse sandwich) just prior to the relevant evolution, or decrement the phase of the pulse following evolution by 90 degrees (or for multiple-quantum experiments, by 90/n).

Experiment Setup

Setup is necessary in 3D experiments to position transmitters and decoupler, adjust pulse widths, etc. Just as the setup of 2D experiments can often be assisted by performing "first increment" experiments (i.e., a 1D experiment that represents the first increment of the 2D), so 3D experiments can be assisted not only by 1D setup experiments, but also by "first plane" 2D experiments (not available on *GEMINI 2000* systems). To perform a 2D experiment in the sw1 dimension, set ni2=1 and phase2=1, with ni greater than 1 and phase=1, 2 (or phase=3 for TPPI experiments). This combination of parameters will perform a "normal" 2D experiment, incrementing d2, and the data can be processed with the wft2da command (or its variants).

The "third dimension" 2D experiment is performed by setting ni=1 and phase=1, with ni2 greater than 1 and phase2=1, 2 (or phase2=3, as desired). These parameters will produce a 2D experiment in which d3 is incremented, resulting in a spectral width sw2.

The wft2d command must be given the special argument ni2 to process this data correctly, for example, wft2d('ni2',1,0,0,0,0,0,-1,0). You cannot use the wft2da('ni2') because the wft2da macro does not support this argument.

Notice that when you process a "first plane" 2D experiment, the axes are always labeled f_1 and f_2 because this is considered to be a 2D experiment, and hence the axis labeling corresponds to conventions used in 2D NMR.

When you are finished setting up the 3D experiment, reset ni, ni2, phase, and phase2 to their desired values. Check the value of the parameter array and make sure that array='phase, phase2' and not 'phase2, phase', which will acquire data in the incorrect order. To ensure the correct order, always enter phase before phase2, or simply enter array='phase, phase2'.

Data Processing

Just like processing 2D NMR, the proper processing of 3D NMR requires coefficients to select various components of the data to be combined to form the final data set. There are actually up to 40 coefficients required that are explained in more detail in the *VNMR Command and Parameter Reference*. In normal operation, the coefficients will be transparent to you, just as the 2D coefficients are. The set3dproc command can create a 3D coefficient file for processing 3D FID data under certain conditions.

The ft3d command determines from the values of phase and phase2 what the expected coefficients are, based on whether a hypercomplex ("States-Haberkorn") or TPPI experiment has been performed in a particular dimension. This assumes that the pulse sequence has been written to perform "standard" phase cycling as described above. If your data are reflected along a particular dimension, it is possible (or probable) that different coefficients are required for data processing. In this case, the ft3d('nocoef') form is used to allow you to specify you own coefficients (which are found in a text file named coef in the 3D experiment directory, unlike in ft2d, where they are given as arguments to the command). By default, ft3d calls the make3dcoef macro to create a coefficient file using the f1coef and f2coef string parameter values.

The format for the 3D coefficient file is an extension of that used for 2D coefficients. The coefficient file contains four rows of eight coefficients used to construct the t_2 hypercomplex interferograms and a final row of eight coefficients used to construct the t_1 interferogram. The actual values of the coefficients depends on the order in which the States-Haberkorn components of the 3D FID data set were collected. This order depends in turn on the values of the parameters phase, phase2, and array.

If TPPI phase cycling is used to collect data along one or both of the indirectly detected dimensions, instead of four data sets per (ni,ni2) increment, there are only two or one data sets, respectively, per (ni,ni2) increment. If there are only two data sets per (ni,ni2) increment, the coef file contains four rows of four coefficients that are used to construct the t_2 hypercomplex interferograms, and a final row of eight coefficients that are used to construct the t_1 interferogram. If there is one data set per (ni,ni2) increment, the coef file contains four rows of two coefficients that are used to construct the t_2 hypercomplex interferograms and a final row of eight coefficients that t_2 hypercomplex interferograms and a final row of eight coefficients that are used to construct the t_1 interferograms and a final row of eight coefficients that are used to construct the t_1 interferograms.

Phasing a 3D data set is best accomplished using 2D transforms. In general, the recommended method in writing 3D pulse sequences is to attempt to minimize frequency-dependent phase shifts in f_1 and f_2 . Even so, there are generally small phase shifts that must be dealt with. The following steps are suggested:

1. Set **pmode='full'** to allow full phasing in both dimensions after a 2D transform.

- 2. Adjust **rp** and **lp** on a 1D spectrum (the first increment of the 3D), just as you would for 2D (e.g., by typing wft(1)).
- You now have an f₁f₃ 2D spectrum (with incorrectly labeled axes). Set trace='f1' to adjust the f₁ phase, then set trace='f2' to trim the f₃ phasing. You can now adjust rp1 and lp1 (as well as rp and lp).
- 5. Enter wft2d('ni2',1,1,0,0,0,0,0,0,0,0,0,0,-1,0,0,0,0,0) to adjust f₂ phasing (note that this argument has nine consecutive zeros in the middle and five zeros at the end).
- 6. You now have an f_2f_3 2D spectrum. Set **trace='f1'** to adjust the f_2 phasing (**rp2** and **lp2**), then set **trace='f2'** to trim the f_3 phasing if necessary.

One additional point on phasing. Some pulse sequences are written to result in a 180° phase shift across the spectrum. Remember that in VNMR, the "origin" for phasing is defined as the right edge of the spectrum; however, in "real" terms, the actual origin of phasing (i.e., the zero-frequency point) is at the center of the spectrum. Thus, if you expect a certain lp1 or lp2 value, such as -180° , you should simultaneously use a value of rp1 or rp2 equal to -1p1/2 or -1p2/2 (e.g., 90°).

If you want to adjust the weighting functions for the 3D transform by using the wti command and examine interferograms, you can do so along either the t_1 or t_2 axes. Use the same commands given above to adjust the phasing (the commands with the long series of zeros), but use wftld instead of wft2d.

For the final transformation, the specdc3d parameter controls the dimensions in which a spectral drift correction is performed on the data. A three-letter value of 'ynn' gives drift correction along f_3 (the first letter) but not along f_1 (the second letter) or f_2 (the third letter); this value is probably a good starting point for your efforts.

The pmode parameter is ignored by the 3D transformation; no phasing is possible after the 3D transform.

The 3D transformation process needs to be followed by the process of extracting the 2D planes from the full 3D data set. This can be done separately, with the getplane command, but most often is combined with the ft3d command. In general, and especially for heteronuclear experiments, the f_1f_3 and f_2f_3 planes are the most interesting. The f_1f_2 plane is not only generally less useful, but also is considerably slower to extract from the data. The recommended command to use for 3D transformation, therefore, is ft3d('ff1f3','f2f3'), which performs the 3D transform and extracts the two interesting planes in one step.

Solvent suppression works on t₃ FIDs of 3D spectra just like in the 1D and 2D cases.

Following the transform, set plane='flf3' or 'f2f3' and then use the dproj macro to display the projection of the data on that plane, or dplane(n) to display the nth plane. The resetf3 macro will reset parameters after a partial 3D Fourier transform.

3.15 4D NMR Acquisition

The addpar('4d') macro creates the parameters ni3, sw3, d4, and phase3 that can be used to acquire a 4D data set (the macro par4d functions the same as addpar('4d')).

The parameter ni3 is the number of t_2 increments, sw3 is the spectral width along the third indirectly detected dimension, d4 is the incremented delay, and phase3 is the phase selection for 4D acquisition. Processing and display in 4D is currently not available in VNMR.

Table 31 summarizes 4D acquisition commands and parameters.

Commands addpar('4d') par4d	Add 4D parameters to the current experiment Create 4D acquisition parameters
Parameters d4 {number, in sec}	Incremented delay for 3rd indirectly detected dimension
ni3 {number}	Number of t_2 increments in 4D acquisition
phase3 {number}	Phase selection for 4D acquisition
sw3 {number,in Hz}	Spectral width along the 3rd indirectly detected dimension

 Table 31. 4D NMR Acquisition Commands and Parameters

Chapter 3. Multidimensional NMR

Chapter 4. Multidimensional and Advanced Experiments

Sections in this chapter:

- 4.1 "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136
- 4.2 "COSY—Correlated Spectroscopy," page 139
- 4.3 "COSYPS—Phase-Sensitive COSY," page 139
- 4.4 "DQCOSY—Double-Quantum Filtered COSY," page 141
- 4.5 "HET2DJ—Heteronuclear 2D-J," page 143
- 4.6 "HETCOR—Heteronuclear Chemical Shift Correlation," page 145
- 4.7 "HETCORPS—Absolute-Value and Phase-Sensitive HETCOR," page 147
- 4.8 "HOM2DJ—Homonuclear J-resolved 2D," page 148
- 4.9 "INADEQUATE—Double-Quantum Transfer Experiment," page 150
- 4.10 "MQCOSY—Multiple-Quantum Filtered COSY," page 151
- 4.11 "NOESY—Nuclear Overhauser Effect Spectroscopy," page 153
- 4.12 "ROESY—Rotating Frame Overhauser Effect Spectroscopy," page 156
- 4.13 "TNCOSYPS—COSYPS with Water Suppression," page 158
- 4.14 "TNDQCOSY—DQCOSY with Water Suppression," page 158
- 4.15 "TNMQCOSY—MQCOSY with Water Suppression," page 159
- 4.16 "TNNOESY—NOESY with Water Suppression," page 159
- 4.17 "TNROESY—ROESY with Water Suppression," page 160
- 4.18 "TNTOCSY—TOCSY with Water Suppression," page 161
- 4.19 "TOCSY—Total Correlation Spectroscopy," page 162
- 4.20 "TROESY—Transverse ROESY," page 164:
- 4.21 "HCCHTOCSY Pulse Sequence," page 164
- 4.22 "HMQCTOCSY Pulse Sequence," page 166
- 4.23 "HMQC-TOCSY 3D Pulse Sequence," page 166
- 4.24 "HSQC-TOCSY 3D Pulse Sequence," page 167

In these experiments, each pulse sequence has a macro, usually with the same name as the pulse sequence, that sets up the parameters for the experiment and then displays information on the experiment. The macro retrieves parameters such as pw90, tpwr, dmf, etc. from a central location like /vnmr/probe or \$vnmruser/probe. Other specific parameters, such as mix for NOESY, come from /vnmr/parlib.

It is important that these parameters be correct. The first time the macro for an experiment is entered (e.g., by typing noesy), the system retrieves the default parameters and values.

To change any of these values (e.g., the default mix time in noesy is 0.2 seconds and you want the default to be 0.5 seconds), make the appropriate change in the displayed parameters, and then save the modified parameters either in your user's parlib or in the system /vnmr/parlib. Notice that because files in the directory /vnmr/parlib are available to all users, only the system administrator vnmrl has permission to save the files in this directory.

To view complete listings of each pulse sequence, print out or look at the contents of the files in the directory /vnmr/psglib on your system disk. You can also enter e dps to view a graphical representation.

The 2D pulse sequences HMQC (Heteronuclear Multiple-Quantum Coherence) and HMQCR (HMQC in Reverse Configuration) are described in Chapter 5, "Indirect Detection Experiments," of this manual.

4.1 Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY

The relayh macro sets up parameters for absolute-value COSY, single RELAY-COSY, or double RELAY-COSY pulse sequences. Figure 25 is a diagram of the absolute-value COSY sequence, and Figure 26 is a diagram of RELAY-COSY, single and double.




Figure 25. Absolute-Value COSY Pulse Sequence




Figure 26. RELAY- COSY Pulse Sequences

Applicability

These sequences are available on all systems.

Parameters

relay is the number of relays to be performed:

if relay=0, a normal absolute value COSY is performed.

```
if relay=1, a RELAY-COSY is performed.
```

if relay=2, a double-RELAY-COSY is performed.

pw is a 90° pulse.

pl is a 90° pulse if relay is not equal to 0; if relay is equal to 0, diagonal peaks can be deemphasized by using pulses greater than 90° when doing P-type peak selection or by using pulses less than 90° when doing N-type peak selection.

d1 is a delay set to 1 to 3 times T_1 .

d2 is the evolution time (t_1 period).

at is the acquisition time (t_2 period).

sw is the spectral width (sw1=sw).

v1, v2, and v5 are pulse phases. oph denotes the phase of the receiver.

ni is the number of t_1 increments (set so that sw1/ni equals 6 by default). If speed is essential, ni can be decreased.

tau is the propagation time for long-range COSY (relay=0) and for relayed COSY (relay greater than 0).

Technique

The COSY experiment, as implemented on UNITY *INOVA*, *MERCURY-VX*, *MERCURY*, UNITY *plus*, *GEMINI 2000*, and UNITY systems, can be run several different ways. The experiment macro called is relayh. If the parameter relay is set to 0 (the default), a normal COSY is performed.

The standard setup uses two 90° pulses and is the traditional way to run COSY experiments. The second pulse can be changed to a 45° pulse to decrease the size of the diagonal. This technique is called "COSY-45." Besides decreasing the diagonal, COSY-45 emphasizes active coupling partners relative to passive couplings so if sufficient digital resolution is present, the crosspeaks may show an interesting multiplet structure. To emphasize long-range couplings, increase the setting of the parameter tau to 0.2 seconds.

To set up and acquire:

- 1. Set up by acquiring a 1D proton experiment in expn, where n is 1, 2, 3, etc.
- 2. Narrow the ¹H spectral window to leave approximately 1 ppm on either side of the peaks of interest by using two cursors and entering **movesw**.
- 3. Reacquire the ¹H with the new spectral window.
- 4. Phase the 1 H spectra.
- 5. Move the FID to another experiment (e.g., to move the FID from expl to exp2, enter mf(1,2)).

- 6. Enter **relayh** to modify the parameters for the RELAYH experiment. The relayh macro calculates appropriate weighting functions.
- 7. Turn the spinner off.
- 8. Enter **au** to acquire data.

To Fourier transform:

- 1. Enter wft2d. This performs weighted Fourier transforms in both dimensions and displays the data as a contour map.
- 2. Adjust the threshold and vertical scale.

To interactively adjust weighting during processing:

- 1. Enter **wft(1)** to Fourier transform the first increment of data.
- 2. Enter wti to start interactive weighting.
- 3. Absolute-value data is usually processed using sinebell weighting. This is the default weighting function calculated by the macro setting up the experiment. Adjust the sinebell so that data decays to zero before the end of the window. Adjusting the weighting function on the first increment of data sets the weighting function in the 2D time dimension t₂.
- 4. Enter **wftld** to Fourier transform the t₂ dimension. A contour map of f₂,t₁ is displayed, showing the individual interferograms.
- 5. Click on the **trace** button and choose a trace through one of the horizontal interferograms.
- 6. Enter **wti** to bring up interactive weighting of the interferogram. Adjust the weighting function the same as before.
- 7. Enter **wft2d** to complete the Fourier transformation.
- 8. Enter **foldt** to symmetrize the data.
- 9. You can adjust the vertical scale and threshold to provide a better display in which the cross peaks and diagonal are more easily seen. Adjust the vertical scale by changing the value of the parameter **vs** (e.g., vs=12).

If the vertical scale is very high, you can decrease it by typing vs=vs/10 (or a similar divisor) to drop the vertical scale by that amount. Adjust the threshold by using the mouse to adjust the sliding color scale beside the 2D contour plot.

- 10. Enter **dcon** to display a non-interactive color intensity map, or enter dconi to permit interaction.
- 11. Enter dpcon(12,1.3) to see what your data will look like if plotted with 12 contours spaced 1.3 levels apart. If you wish, try other arguments for dpcon to see how other values of contours and levels will look when plotted.
- 12. Enter plcosy(8,1.5) to plot the COSY or RELAYH data with eight contours spaced 1.5 levels apart. Try other values if you want. The macro plcosy actually takes three arguments: the number of contours, the spacing of levels, and the experiment number that contains the 1D data. If there is no third argument, plcosy assumes that the 1D data is in experiment 1. Therefore, plcosy(8,1.5) plots the contour, retrieves the 1D data, and plots it above the contour.
- 13. Enter **dconi** to make an expansion of the data and redisplay the contour map if it is not present. Expand around the upfield region by using the left mouse button to

determine the lower left corner and the right mouse button to determine the upper right corner, then click the expand button.

- 14. Click on the **expand** button.
- 15. Enter **plcosy** to plot this expanded region. Note that the appropriate 1D region is plotted on the top and to the side of the 2D plot.

Potential Problems

The COSY experiment is very forgiving and usually works (to a greater or lesser extent) with a slightly incorrect 90° pulse width. The most common reason for failure of the experiment is that the value for p1 is incorrect. Make sure that the 90° pulse is correct before beginning the experiment. Remember that pw, the second pulse, may be set to 45° to deemphasize the diagonal and provide a filter for passive couplings.

A second problem that can arise are artifacts caused by pulsing too rapidly. If you symmetrize the spectrum with foldt (normal in COSY), these should be minimized.

4.2 COSY—Correlated Spectroscopy

The cosy macro converts a parameter set to a COSY experiment.

Applicability

COSY is available on all systems.

Parameters

For information on COSY parameters, refer to page 136.

Technique

To set up, acquire, and process, use the technique given in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.

4.3 COSYPS—Phase-Sensitive COSY

The cosyps macro sets up a phase-sensitive COSY pulse sequence. Figure 27 is a diagram of the COSYPS sequence.

Applicability

COSYPS is available on all systems.

Parameters

pw is a 90° pulse.

p1 is 90° if phase is not equal to 0; if phase equals 0, pulses greater than 90° deemphasize diagonal peaks for P-type peak selection.




Figure 27. COSYPS Pulse Sequence

d1 is set to equal 1 to 3 times the value of T_1 .

d2 is the evolution time (t_1 period).

at is the acquisition time (t_2 period).

sw is spectral width (usually sw1=sw, except for phase=3).

ni is the number of t_1 increments (by default swl/ni=12); it may be decreased if speed is essential.

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr.presat does not depend on dm but does depend on dmm and is activated as a part of d1 if presat is greater than 0.presat is available on *MERCURY-VX* and *MERCURY* systems with *GLIDE* (uppercase pulse sequences only).presat is not active on *GEMINI 2000* systems.

phase=1, 2 (suggested value) for phase-sensitive data. phase=0 for 2D data in an *av* display (P-type peaks). phase=1, 2 for 2D hypercomplex data (States-Haberkorn method). phase=3 for 2D TPPI data. For phase=3, remember that sw1=2*sw.

nt is multiple of 8 (minimum, phase=0) or multiple of 4 (minimum, phase=1, 2 or 3).

Notice that for all TPPI experiments, the resulting spectrum appears doubled. Just display and phase one-half of the data.

Phase Cycling

v1 is the phase for first pulse. v2 is the phase for second pulse. oph is the phase for the receiver.

v1 = x -x y -y y -y x -xv2 = x x y y y y x xoph = x -x y -y y -y x -x

These phases are for phase=1. For phase=2, add 90° to v1. For phase=3, add $90*(ix - 1)^\circ$ to v1, where ix is the increment counter. For phase=0, the subcycle of P-type peak selection is mixed into v2.

Technique

To set up and acquire data, use the technique given in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.

To process the data, enter wft(2) to display the first increment with phase=2 in t_1 , phase it to be pure absorptive, then add 90° to rp. Set rp1=90 and lp1=0 initially. If the f_1 phasing is not satisfactory, adjust rp1 and lp1 so that the diagonal peaks are pure absorptive and then add 90° to rp1.

References

Bodenhausen, G.; Freeman, R.; Niedermeyer, R.; Turner, D. L. J. Magn. Reson. 1977, 26, 133–164.

Bachmann, P.; Aue, W. P.; Muller, L; Ernst, R. R. J. Magn. Reson. 1977, 28, 29-39.

4.4 DQCOSY—Double-Quantum Filtered COSY

The dqcosy macro sets up parameters for the DQCOSY pulse sequence. Figure 28 is a diagram of the sequence.




Figure 28. DQCOSY Pulse Sequence

Applicability

DQCOSY is available on all systems.

Parameters

pw is a 90° pulse on the observed nucleus.

d1 is a relaxation delay (1 to 3 times the value of t_1).

d2 is the evolution time (t_1 period).

at is the acquisition time (t_2 period).

sw is the spectral width (usually sw1=sw, except for phase=3).

ni is the number of t_1 increments (set so that sw1/ni=12).

phase=0 for 2D data in an av display (P-type peaks); phase=1, 2 for 2D hypercomplex data (States-Haberkorn method); phase=3 for 2D TPPI data.

sspul='y' activates a homospoil–90–homospoil sequence that precedes d1; this is used to achieve a less oscillatory steady-state for 2D experiments where the recycle time is shorter than t₁. sspul is not active on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems.

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr.presat does not depend on dm, but does depend on dmm, and is activated as a part of d1 if presat is greater than 0. presat is available on *MERCURY-VX* and *MERCURY* systems with *GLIDE* (uppercase pulse sequences only). presat is not active on *GEMINI 2000* systems.

nt is a multiple of 8 (phase=0) (minimum) or multiple of 4 (phase=1, 2, or 3) (minimum).

Phase Cycling

v1 is the phase for first pw pulse. v2 is the phase for second pw pulse. v3 is the phase for third pw pulse.

oph is the phase for receiver.

These phases are for phase=1. For phase=2, add 90° to v1. For phase=3, add 90* $(ix - 1)^{\circ}$ to v1, where ix is the increment counter. For phase=0, the subcycle of p-type peak selection is added in after the basic four-step cycle of v3.

Technique

To set up and acquire data, use the technique given in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.

No phasing should be necessary for processing. For full processing, use menus or enter wft2da. For processing, proceed as follows:

- 1. Enter wft(1) to Fourier transform the first increment of data. There should be no signals in the first increment of a DQCOSY.
- 2. Enter wti to start interactive weighting. Data from DQCOSY is phase-sensitive and usually processed using Gaussian weighting, which is the default weighting function. Adjust the Gaussian weighting so that the data decays to zero before the end of the window. Adjusting the weighting function on the first increment of the data sets the weighting function in the 2D time dimension t₂.
- 3. Enter **wftlda** to Fourier transform the t_2 dimension and display a contour map of f_2, t_1 showing individual interferograms. Click on the trace button and choose a trace through one of the interferograms. Enter **wti** to bring up interactive weighting of the interferogram. Adjust the weighting function as before.

- 4. Enter **wft2da** to complete the Fourier transformation.
- 5. Display and plot using the procedures described in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.

Potential Problems

The DQCOSY experiment is much more sensitive to artifacts than the COSY experiment. DQCOSY experiments are not symmetrized so t_1 noise is a much greater problem. To minimize t_1 noise, make sure that the 90° pulse is correct before beginning the experiment. Another problem that can arise is the presence of artifacts due to pulsing too rapidly. Make sure you set dl to at least 1 to 3 times the T_1 value of the protons in the sample. Run the sample nonspinning and use homospoil pulses (homospoil is not available on *GEMINI 2000* systems).

References

Piatini U.; Sorenson, O. W.; Ernst, R. R. J. Am. Chem. Soc. **1982**, 104, 6800–6801. Rance, M. et al Biochem. Biophys. Res. Comm. **1983**, 117, 479–485.

4.5 HET2DJ—Heteronuclear 2D-J

The het2dj macro sets up parameters for the HET2DJ experiment. Absolute-value (av) mode is required. The experiment can be performed either in a gated or a non-gated mode. Figure 29 is a diagram of the sequence.




Figure 29. HET2DJ Pulse Sequence

Applicability

HET2DJ is available on all systems.

Parameters

pw is the observe 90° pulse (X nucleus).

pp is the proton 90° pulse width supplied from decoupler.

p1 is the observe 180° pulse (if p1=0, 2*pw is used).

tpwr is the power level for observe pulses (systems with linear amplifiers).

pplvl is the power level for proton pulses (systems with linear amplifiers).

sw1 should cover the maximum multiplet pattern.

dm='ynyy' for decoupler gating during the evolution period. dm='ynny' for no decoupler gating during the evolution period (proton flip experiment).

dmm='wcw' (WALTZ) or dmm='fcf' (no WALTZ).

nt is a multiple of 2 (minimum) or a multiple of 16 (maximum and recommended). Note that 1 is possible.

If decoupler gating is used during the evolution period, the actual J value is twice the measured J ($J_{act} = 2*J_{meas}$).

Technique

- 1. A good sample to try the first time is 30% menthol in CDCl₃. This sample is concentrated enough that good data can be obtained in a short period of time.
- 2. Enter jexp2 to join experiment 2 (this can also be done from the Workspace menu).
- 3. Set up a carbon experiment by clicking the C13,CDCl3 button in the setup menu.
- 4. Acquire a ¹³C spectrum.
- 5. Enter mp(2,3) jexp3 to move the parameters from experiment 2 to experiment 3 and join experiment 3.
- 6. Enter **het2dj** to set up the HET2DJ experiment.
- 7. Enter **go** to acquire the data.
- 8. Enter **wft2d** to Fourier transform.
- 9. After transformation, the data is tilted along a 45° angle. To rotate the data, enter **rotate(45.0)**.
- 10. Enter **foldj** to symmetrize the data.
- 11. Display the data the same as described in "COSY—Correlated Spectroscopy," page 139.
- 12. Enter **pcon(8,1.5)** page to plot the data with 8 contours spaced 1.5 levels apart. If you wish, try other values for the arguments to see various plots.

The pcon command actually takes three arguments: positive or negative, the number of contours, and the spacing of levels. Positive or negative refers to the ability to plot exclusively positive or negative contours. The default is to plot both positive and negative contours. Since the HET2DJ data has no negative contours, this option can be ignored.

Potential Problems

The most common reason for failure of the HET2DJ experiment is that the values for pp and pw are incorrect. Make sure that the 90° pulse on the decoupler, pp and the 90° pulse on ^{13}C , pw are correct before beginning the experiment.

Another possible problem is the presence of artifacts due to pulsing too rapidly. This experiment is carbon detected. Make sure you set dl to a value that permits sufficient relaxation, the T_1 value of the protonated carbons.

References

Bodenhausen, G.; Freeman, R.; Turner, D. J. Magn. Reson. 1977, 27, 511.

Freeman, R.; Keeler, J. J. Magn. Reson. 1981, 43, 484-487.

4.6 HETCOR—Heteronuclear Chemical Shift Correlation

The hetcor<(exp_number)> macro sets up parameters for HETCOR pulse sequence, where the optional argument exp_number is the number of the experiment, from 1 through 9, in which a proton spectrum of the sample already exists, (e.g., hetcor(2) specifies experiment 2).

hetcor includes a presaturation option, composite 180° pulses, and simultaneous pulses on transmitter and decoupler rf channels. Figure 30 is a diagram of HETCOR set for decoupled multiplets (hmult='n').




Figure 30. HETCOR Pulse Sequence

Applicability

HETCOR is available on all systems.

Parameters

nt should be a multiple of 4.

dl is 1 to 2 times the proton T_1 recommended; it must not equal 0.

j1xh is an average one-bond X-H coupling constant.

jnxh is an average 2 or 3-bond coupling constant; set it only if long-range correlations are wanted (hmult='y' only in this case).

hmult='y' for H-H multiplets; hmult='n' for decoupled multiplets.

pp is a 90° pulse on protons.

pw is a 90° pulse on observe nucleus.

pplvl is the power level for the proton pulse on the decoupler channel (systems with a linear amplifier on the decoupler rf channel).

tpwr is the power level for the heteronuclear pulse on the observe channel (systems with a linear amplifier on the observe rf channel).

sw1 should equal the spectral width sw used to obtain proton spectrum.

dof should equal the value of tof used to obtain proton spectrum.

presat='y' gives presaturation pulse train; preset='n' for no presaturation presat is available on *MERCURY-VX* and *MERCURY* systems with *GLIDE* (uppercase pulse sequences only). presat is not active on *GEMINI 2000* systems.

rfll, rfpl, wpl, spl should be set to the corresponding 1D parameters from the 1 H spectrum.

sw1 and dof should be set before typing ga.

Technique

The experiment can be run two different ways: (1) the standard setup optimizes for onebond proton-carbon couplings, or (2) setting jnxh to the average two- or three-bond coupling constant emphasizes three-bond (and two- or four-bond) couplings.

- 1. A good sample to try the first time is 30% menthol in CDCl₃. This sample is concentrated enough that good data can be obtained in a short period of time.
- 2. If not in experiment 1 already, join experiment 1 by entering **jexp1** or by using menu buttons.
- Set up a proton in CDCl₃ experiment by clicking the H1,CDCl3 button in the Setup menu. If you are using the 30% menthol sample, set pw=1, gain=0, nt=1; otherwise, acquire the spectrum with the standard parameters.
- 4. Narrow the spectral window to include only the peaks of interest by placing a cursors approximately 1 ppm upfield and 1 ppm downfield of the region containing the peaks of interest and entering movesw.
- 5. Enter **ga** to acquire the spectrum with the new window and phase the data. This spectrum serves as a reference spectrum for the proton section of the HETCOR.
- 6. Enter jexp2 to join experiment 2 and set up a carbon experiment by clicking the C13,CDCl3 button in the Setup menu. Acquire a spectrum. If your are using the menthol or a similar sample, the sample is concentrated enough that one scan should be sufficient to see a carbon spectrum with good signal to noise.
- Place a cursor upfield and a cursor downfield of the carbon signals. Narrow the spectral window by entering movesw. Set nt=16 and reacquire the carbon spectrum. Phase the spectrum. This spectrum will be used as the reference for the ¹³C dimension of the HETCOR.
- 8. Enter mp(2,3) jexp3 to move the ¹³C parameters from experiment 2 to experiment 3 and join experiment 3.
- 9. Enter **hetcor** to set up the HETCOR experiment and show a time estimate for performing the experiment.
- 10. If you are using a concentrated sample such as the menthol sample, you can modify the default parameters to produce acceptable data in a much shorter time: set

nt=16, ni=64, np=1024, fn=1024, fn2=256, and d1=1. Then enter time to check the experiment time. These parameter values should reduce the time to 20 minutes yet provide data with sufficient resolution to assign the spectrum.

- 11. Enter **go** to acquire the data.
- Enter wft2d to Fourier transform the data. Adjust vertical scale and threshold according to the procedures outlined in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.
- 13. Enter **plhxcor(8,1.5)** to plot HETCOR data with 8 contours spaced 1.5 levels apart.

If you wish, try other values to create different plots. plhxcor can have up to four arguments: number of contours, spacing of levels, experiment number for the proton data, and experiment number for the carbon data. If the third or fourth argument is omitted, plhxcor assumes that the proton 1D data is in experiment 1 and that the carbon 1D data is in experiment 2. Therefore, plhxcor plots the contour, retrieves the 1D data, and plots them above the contour. plhxcor also plots expansions with the appropriate 1D plots on the side.

Potential Problems

The most common reason for failure is that the values for pp and pw are incorrect. Make sure that the 90° pulse on the decoupler, pp and the 90° pulse on 13 C, pw are correct before beginning the experiment.

Since both DEPT and HETCOR rely on polarization transfer between proton and carbon, a simple method to troubleshoot a HETCOR is to run DEPT on the sample. If the DEPT does not work, HETCOR probably will not work either. Once the DEPT has been properly calibrated and is working, the same values of pw90, pp (and pplvl if the parameter exists) should be used for the HETCOR.

Another possible problem is the presence of artifacts from pulsing too rapidly. Because of the nature of the experiment, the critical T_1 s are the proton T_1 s (usually 1 to 10 seconds), not the relaxation times of the carbons (which may be much longer). Make sure you set d1 to at least the time of T_1 of the protons in the sample. A common type of artifact is an extra peak exactly between the correlations of two non-magnetically equivalent protons attached to the same carbon.

References

Bax, A.; Morris, G. A. J. Magn. Reson. 1981, 42, 501.
Bax, A. J. Magn. Reson. 1983, 53, 51.
Rutar, V. J. Magn. Reson. 1984, 58, 306.
Wilde, J.; Bolton, P. J. Magn. Reson. 1984, 59, 343–346.

4.7 HETCORPS—Absolute-Value and Phase-Sensitive HETCOR

The hetcorps macro sets up parameters for a HETCORPS (Heteronuclear Chemical-Shift Correlation, Absolute Value and Phase Sensitive) pulse sequence.

Applicability

HETCORPS is not available on MERCURY-VX, MERCURY, or GEMINI 2000 systems.

Parameters

pw is a 90° pulse on the observe nucleus.

tpwr is transmitter power level; only for systems with a linear amplifier on the transmitter channel.

pp is a proton 90° pulse on the decoupler channel.

pplvl is decoupler power level; only for systems with a linear amplifier on the decoupler channel; otherwise, the decoupler is turned to full-power for pulses on systems that have bilevel decoupling capability.

dhp is decoupler power level during acquisition.

dpwr is decoupler power level during acquisition for systems with linear amplifiers.

hmult='n' removes non-geminal proton-proton couplings in F1; hmult='y' preserves all proton-proton couplings in F1.

chonly='y' gives CH only spectrum.

oddeven='y' gives CH and CH₃ positive and CH₂ negative; oddeven='n' gives all positive. oddeven is irrelevant in av (phase=0) spectrum or if chonly='y'.

j1xh is a one-bond heteronuclear coupling constant.

phase=1, 2 gives hypercomplex; phase=0 gives absolute value.

nt=1 is the minimum for hypercomplex; nt=2 is the minimum for absolute value. nt set to a multiple of 2 is recommended for hypercomplex; nt=4 is recommended for absolute value.

Recommendations

HETCORPS cannot and should not be used for long-range correlation. For long-range correlation, use lrhetcor.c.

If ¹³C parameters (after phase correcting a ¹³C spectrum) are moved to set up HETCORPS, there should be no need to phase correct either F1 or F2.

Use gaussian windows.

4.8 HOM2DJ—Homonuclear J-resolved 2D

The hom2dj macro sets up parameters for a HOM2DJ pulse sequence. Absolute-value mode is required. Figure 31 is a diagram of the sequence.

Applicability

HOM2DJ is available on all systems.




Figure 31. HOM2DJ Pulse Sequence

Parameters

pw is a 90° observe pulse.

pl is a 180° observe pulse.

nt is a multiple of 2 (minimum) to a multiple of 16 (maximum and recommended).

d1 is set for 1 to 3 times T_1 (maximum); default is 1.0 seconds.

sw1 should cover maximum multiplet pattern; the default is 50 Hz.

Technique

HOM2DJ is one of the earliest and simplest 2D experiments. Rarely used today, it has been mostly replaced by phase-sensitive COSY experiments such as DQCOSY.

- 1. A good sample to try at first is a 5% heptanone sample in CDCl₃. This sample is concentrated enough that good data can be obtained in a short time. If not there already, join experiment 1 by entering jexp1. Set up a proton in CDCl₃ by using the H1,CDCl3 button in the Setup menu. Acquire a spectrum.
- 2. Narrow the spectral window to include only the peaks of interest by placing a cursors about 1 ppm upfield and 1 ppm downfield of the region containing the menthol peaks and entering movesw. If a peak is present from residual CHCl₃, do not include this peak because the peak will fold in but the filters will remove most of the peak. This greatly decreases the amount of spectrum that is just noise and provides better digital resolution in the spectrum and decreased experiment time.
- 3. Enter **ga** to acquire the spectrum with the new window and phase the data.
- 4. Enter **mp(1,2) jexp2** to move the parameters from experiment 1 to experiment 2 and join experiment 2.
- 5. Enter **hom2dj** to set up the HOM2DJ experiment and display an estimate of the time required to perform the experiment.
- 6. You can modify the default parameters to produce acceptable data in a much shorter period of time. If you are using the heptanone sample, set nt=8, ni=64, np=512, fn=128, fn2=512, and dl=1. This provides an experiment time of 12.5 minutes yet gives sufficient resolution to assign the spectrum.
- 7. Enter **go** to acquire the data.
- 8. Enter wft2d to Fourier transform the data.
- 9. If you wish to interactively process the data, the data is in absolute-value mode and interactive processing is the same as for COSY (see page 136).

- 10. After transformation, the data is tilted along an angle of 45°. Enter **rotate(45.0)** to rotate the data.
- 11. Enter **foldj** to symmetrize the data.
- 12. Display the data contours in the manner described in "Absolute-Value COSY, RELAY-COSY, Double RELAY-COSY," page 136.
- 13. Plot the data using the procedure described in "HET2DJ—Heteronuclear 2D-J," page 143.

Potential Problems

The HOM2DJ experiment is very forgiving and usually works (to a greater or lesser extent) with slightly incorrect pw and pl values. The most common reason for failure of the experiment is that the values for pw are incorrect. Make sure that the 90° pulse is correct before beginning the experiment.

A major problem, inherent to the experiment, is that HOM2DJ results in phase twisted lineshapes that can make interpretation difficult in crowded regions and causes artifacts (which show up as extra peaks) due to strong coupling.

References

W. Aue, J. Karhan, and R. Ernst, J. Chem. Phys. 64:4226 (1976).

K. Nagayama, O, Backmann, K. Wuthrich, and R. Ernst, J. Magn. Reson. 31:133 (1978).

4.9 INADEQUATE—Double-Quantum Transfer Experiment

The inadqt macro sets up parameters for an INADEQUATE (Incredible Natural Abundance Double-Quantum Transfer Experiment) pulse sequence. Figure 32 shows a diagram of the sequence.




Figure 32. INADEQUATE Pulse Sequence

Applicability

INADEQUATE is available on all systems except GEMINI 2000.

Parameters

pw is a 90° pulse on observed nucleus (e.g., carbons) at power equal to tpwr.

tau is set to $1/(4*J_{CC})$.

jcc is the average scalar coupling constant between the two heteronuclei (usually onebond constants).

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

dmf sets the modulation frequency for WALTZ-16 broadband decoupling of protons at power equal to dhp or dpwr throughout the pulse sequence.

nt is a multiple of 8 (minimum: phase=0), a multiple of 128 (maximum: phase=0), a multiple of 4 (minimum: phase=1, 2 or phase=3), or a multiple of 64 (maximum: phase=1, 2 or phase=3).

dm is set to 'y'.

dmm is set to 'w'.

hs is set to 'yn'.

hst is set to 0.01.

phase=0 for 2D absolute-value data, phase=1, 2 for 2D hypercomplex data (States-Haberkorn method), or phase=3 for 2D TPPI data. phase=1, 2 is the suggested value. For phase=3, remember that sw1 must be set to *twice* the desired value.

Technique

Set up the sequence from an existing ¹³C spectra.

If data was acquired with phase=0, process with wft2d, but if data was acquired with phase=1,2 or phase=3, process with wft2da. If phase-sensitive data without f_1 quadrature are desired, set phase=1 and process with wft2da.

This experiment must be performed non-spinning and with VT regulation. To obtain reasonable results within 24 hours, a one-transient ¹³C spectrum should present a signal-to-noise ratio of at least 25:1 (1.5 sec. recycle time, 24 hours total acquisition time, 128 total increments). For 1D spectra, set phase=1 for maximum sensitivity.

4.10 MQCOSY—Multiple-Quantum Filtered COSY

The mqcosy<(level)> macro sets up parameters for a MQCOSY pulse sequence, where the optional argument level is the desired quantum level of filtration (e.g., mqcosy(3)). Figure 33 is a diagram of the sequence.

Applicability

MQCOSY is available on all systems except MERCURY-VX, MERCURY, and GEMINI 2000.

Parameters

pw is the 90° pulse on the observed nucleus.

d1 is the relaxation delay (1 to 3 times the value of t_1).




Figure 33. MQCOSY Pulse Sequence

sw is the spectral width (usually sw1=sw, except for phase=3).

qlvl is the quantum level of filtering desired in the experiment.

phase=1, 2 for 2D hypercomplex data (States-Haberkorn method) or phase=3: 2D TPPI data. phase=1, 2 is suggested. Note: For phase=3, remember that sw1=2*sw (for homonuclear experiments) or that sw1 must be set to *twice* the desired value (for heteronuclear experiments).

sspul='y' activates a homospoil-90-homospoil sequence that precedes d1; this is used to achieve a less oscillatory steady-state for 2D experiments where the recycle time is shorter than t_1 .

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr.presat does not depend on dm but does depend on dmm and is activated as a part of dl if presat is greater than 0.

nt is a multiple of qlvl (minimum) or a multiple of 2*qlvl (recommended).

Technique

Set up the MQCOSY in the same manner that you set up DQCOSY. Set the parameter glvl to the desired multiple-quantum level.

To process, enter wft2da or use the menu buttons to perform full processing.

To interactively process, use the technique given in "DQCOSY—Double-Quantum Filtered COSY," page 141.

References

Piantini, U.; Sorenson, O.; Ernst, R. J. Am. Chem. Soc. **1982**, 104, 6800–6801. Rance, M. et al. Biochem. Biophys. Res. Comm. **1983**, 117, 479–485.

4.11 NOESY—Nuclear Overhauser Effect Spectroscopy

The noesy macro sets up parameters for a NOESY Laboratory Frame Overhauser or 2D Exchange pulse sequence. It can be performed in either a phase-sensitive or absolute-value mode. Either TPPI or the hypercomplex method can be used to achieve f_1 quadrature in a phase-sensitive presentation. No attempt is made to suppress J-cross peaks. Figure 34 is a diagram of the sequence.




Figure 34. NOESY Pulse Sequence

Applicability

NOESY is available on all systems.

Parameters

pw is a 90° pulse on the observed nucleus (power = tpwr).

d1 is the relaxation delay (1 to 3 times the value of T_1).

d2 is the evolution time (t_1 period).

at is the acquisition time (t_2 period).

sw is the spectral width (sw1=sw usually, except for phase=3).

ni is the number of t_1 increments (set up by default so that sw1/ni=12).

mix is the mixing time for magnetization exchange.

phase=1, 2 is the suggested value; use phase=0 for P-type peak selection for av display; phase=1, 2 for 2D hypercomplex data (States-Haberkorn method); phase=3 for 2D TPPI data.

sspul='y' activates a homospoil-90-homospoil sequence that precedes d1; this achieves a less oscillatory steady-state for 2D experiments where recycle time is shorter than T_1 , sspul has no effect on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems.

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr. presat does not depend on dm but does depend on dmm and is activated as a part of d1 if presat is greater than 0. presat is available on *MERCURY-VX* and *MERCURY*

systems with *GLIDE* (uppercase pulse sequences only). presat is not active on *GEMINI* 2000 systems.

dm= 'nnnn' is no decoupler presaturation (unless presat>0); dm= 'nynn': decoupler presaturation during evolution period; dm= 'nnyn': decoupler presaturation during mixing period; dm= 'nnny': homonuclear decoupling during t2; dm= 'nyyn' is recommended when solvent saturation is needed.

nt is a multiple of 8 (minimum) for phase=1, 2 or phase=3; nt is a multiple of 16 (minimum) for phase=0.

Phase Cycling

v1 is the phase for the second pw pulse. v2 is the phase for the first pw pulse. v3 is the phase for the third pw pulse. oph is the phase for receiver.

v1	=	x	х	х	х	х	x	х	х	У	У	У	У	У	У	У	У
		-x	-y														
v2	=	х	-x	х	-x	х	-x	х	-x	У	-у	У	-y	У	-у	У	-у
v3	=	х	х	У	У	-x	-x	-y	-y	У	У	-x	-x	-у	-у	х	x
oph	=	х	-x	У	-y	-x	х	-y	У	У	-y	-x	х	-y	У	х	-x
		-x	х	-y	У	х	-x	У	-y	-y	У	х	-x	У	-y	-x	x

These phases are for phase=1. For phase=2, add 90° to v1. For phase=3, add $90^{\circ}(ix-1)^{\circ}$ to v1, where ix is the increment counter. For phase=0, the P-type selection subcycle is added in after the first eight steps in the phase cycle.

For phase=3, remember that sw1 equals 2*sw (for homonuclear experiments) or that sw1 must be set to *twice* the desired value (for heteronuclear experiments).

If mix is arrayed, then phase=3 is necessary in order to be able to process the data with the wft2dac macro.

Technique

Because NOESY is a phase-sensitive 2D experiment, NOESY spectra may need to be phased before good results are obtained.

- A good sample to try the first time would be a 5% sucrose sample in D₂O. This sample is concentrated enough that good data can be obtained in a relatively short time. If not there already, join experiment 1 by entering jexp1 or by using menu buttons. Set up a proton in D₂O by clicking on the following buttons from the Main Menu: Setup > Nucleus, Solvent > H1 > D2O. Acquire a spectrum.
- 2. Narrow the spectral window to include only the peaks of interest by placing a cursors about 1 ppm upfield and 1 ppm downfield of the region containing the peaks and enter movesw. This greatly decreases the amount of spectrum that is just noise and provides better digital resolution in the spectrum and decreased experiment time. Generally you want to obtain better digital resolution in phase-sensitive experiments than in absolute-value experiments.
- 3. Enter **go** to acquire the spectrum with the new window and phase the data. This spectrum will serve as a reference spectrum for the NOESY.
- 4. Enter **mp(1,4) jexp4** to move the parameters from experiment 1 to experiment 4 and join experiment 4.

- 5. Enter **noesy** to set up the NOESY experiment and display an estimate of the time required to perform the experiment.
- 6. If you are running the sucrose sample, set nt=16, ni=128, np=1024, fn=1024, fn2=1024, and d1=2. The experiment takes about three hours.
- 7. Enter **go** to acquire the data.
- 8. Enter wft(1) to Fourier transform the first increment of data. The spectrum appears inverted, but this is normal. When the spectra are Fourier transformed, the diagonal is below the plane of the 2D, and NOE crosspeaks due to positive NOEs appear positive, and crosspeaks due to negative NOEs are negative. Phase the spectrum, leaving it properly phased but inverted. This sets the phase for the f₂ dimension of the 2D.
- 9. Enter wti to start interactive weighting. NOESY data is phase-sensitive and usually processed using Gaussian weighting, the default weighting function calculated by the setup macro. Adjust the Gaussian weighting so that data decays to zero before the end of the window. Adjusting the weighting function on the first increment of the data sets the weighting function in the 2D time dimension t₂.
- 10. Enter **wftlda** to Fourier transform the t_2 dimension. A contour map of f_{2,t_1} is displayed, showing individual interferograms. Click on trace and choose a trace through one of the interferograms. Enter **wti** to bring up interactive weighting of the interferogram. Adjust the weighting function as before.
- 11. Enter wft2da to complete the Fourier transformation.
- 12. The f_1 dimension may now need phasing. To phase f_1 , click on trace and select a trace at the top (upfield) section of the 2D. Enter **ds** to display the trace, and phase normally using the parameter **rp**.
- 13. Enter dconi to redisplay the contour map with a new rp. Click on trace again, and select a trace at the bottom (downfield) section of the 2D. Enter ds and click on phase. Move the cursor upfield and click. *Do not adjust the phase at this point*. Clicking at this point sets rp and retains the rp value obtained previously. Move downfield and click. Adjust phase normally. This adjusts lp. Enter dconi to redisplay the properly phased 2D.
- 14. Enter **plcosy** to plot the data. The plcosy macro is general and plots all homonuclear correlated data.

Potential Problems

Unlike the COSY experiment, obtaining good NOESY spectra requires proper values of 90° pulse width and a consideration of delay times. Make sure that the 90° pulse is correct before beginning the experiment. If 90° pulse is incorrect, many "COSY type" (i.e. antiphase cross-peaks) appear, complicating the analysis. In small molecules, some "COSY-type" cross-peaks may be unavoidable even when everything is carefully calibrated. Fortunately these can be easily distinguished because of their antiphase nature, i.e., the cross-peaks have both positive and negative components but true NOE peaks are pure absorptive.

Another possible problem is the presence of artifacts due to pulsing too rapidly. Make sure you set dl to at least 1 to 3 times the T_1 of the protons in the sample. The NOESY experiment must be interpreted more carefully than the COSY experiment because cross-pulses arise from COSY interaction as well as dipolar interaction.

Reference

States, D. J.; Haberkorn, R. A.; Ruben, D. J. J. Magn. Reson. 1982, 48, 286–292.

4.12 ROESY—Rotating Frame Overhauser Effect Spectroscopy

The roesy<(ratio)> macro sets up parameters for a ROESY pulse sequence, where the optional argument ratio is the desired value of the parameter ratio used in the sequence (ratio is not used in the ROESY sequence provided with *MERCURY-VX* and *MERCURY*). Either time-shared or a continuous-wave spin lock can be used. Compensation for off-resonance effects is performed. Figure 35 shows the sequence.




Figure 35. ROESY Pulse Sequence

Applicability

ROESY is available on all systems except *GEMINI 2000* systems; however, the ROESY pulse sequence on *MERCURY-VX* and *MERCURY* systems are slightly different from other systems. *MERCURY-VX* and *MERCURY* use a simple CW spinlock, and there are only two parameters that determine the spinlock: slpwr is the power for the CW spinlock and mix is the time. sspul and ratio do not exist in the *MERCURY-VX* and *MERCURY* parameter sets. presat is available on *MERCURY-VX* and *MERCURY* systems with *GLIDE* (uppercase pulse sequences only).

Parameters

ratio is used to adjust the effective spin-lock field strength provided by the transmitter during the spin-lock period. As ratio is increased, the effective spin-lock field strength is decreased. If ratio is set too short, TOCSY peaks begin to appear regardless of the value for pw. If ratio is set too long, the protons are not spin locked in the rotating frame and therefore do not exchange magnetization via cross relaxation in the rotating frame. The parameter ratio should be increased until the protons of interest begin to show sufficient amplitude and phase stability.

mix (the mixing time) is the length of the spin-lock, which can be a continuous wave (ratio=0) or pulsed (ratio is not 0). If pulsed, a series of pulses of length pw and delays of length ratio*pw are applied to form the spin-lock period. The effective spin-lock field strength (B(spinlock)) can be calculated by the relation B (spin-lock) = (x/(360*pw)/(ratio + 1)), where x is the flip angle of the pw pulse.

Note that to achieve resonance offset compensation, this sequence employs two hard 90° pulses on both sides of the spin-lock. If resonance offset compensation is not desired, set rocomp='y'.

p1 is a 90° pulse on the observed nucleus (power is p11v1).

pw is an x° (default is 30°) pulse on the observed nucleus (power is tpwr); note that pw and ratio are used to define the B₁ field for the spin-lock. A typical value of ratio is 10 for a 30° pw pulse.

- d1 is the relaxation delay (1 to 3 times the value of T_1).
- d2 is the evolution time (t_1 period).
- at is the acquisition time (t_2 period).

sw is the spectral width (sw1=sw usually, except for phase=3).

ni is the number of t_1 increments (set up so that sw1/ni equals 12).

phase=1, 2 (suggested value), phase=1, 2 for 2D hypercomplex data (States-Haberkorn method), or phase=3 for 2D TPPI data.

sspul='y' activates a homospoil-90-homospoil sequence that precedes d1; this is used to achieve a less oscillatory steady-state for 2D experiments where the recycle time is shorter than t_1 .

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr. presat does not depend on dm but does depend on dmm and is activated as a part of d1 if presat is greater than 0.

nt is the multiple of 4 (phase=1, 2 or phase=3) (minimum).

Phase Cycling

v1 is the phase for pw pulse. v2 is the phase for the spin lock. oph is the phase for receiver.

v1 = y -y -x x v2 = x x y y -x -x -y -y oph = y -y -x x

These phases are for phase=1. For phase=2, add 90° to v1. For phase=3, add $90^{\circ}(ix-1)^{\circ}$ to v1, where ix is the increment counter.

For phase=3, remember that sw1 equals 2*sw (for homonuclear experiments) or that sw1 must be set to *twice* the desired value (for heteronuclear experiments).

Technique

To set up and acquire, use the technique given in "NOESY—Nuclear Overhauser Effect Spectroscopy," page 153.

To process for normal phase=1, 2, use the technique (except spectra are *not* inverted) given in "NOESY—Nuclear Overhauser Effect Spectroscopy," page 153. If mix is arrayed, phase=3 is necessary to be able to process the data with the wft2dac macro.

Reference

Kessler, et al. J. Am. Chem. Soc. 1987, January.

4.13 TNCOSYPS—COSYPS with Water Suppression

The tncosyps macro sets up parameters for homonuclear correlation (phase-sensitive version) that uses transmitter solvent saturation.

Applicability

TNCOSYP requires a system with a linear amplifier on the observe channel and a T/R switch. It is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 139 for a discussion of COSYPS.

Parameters

satmode determines when the saturation happens. satmode should be set analogously to dm, (e.g., satmode='yyn' or satmode='ynn').

satpwr is the power level during saturation period(s).

satdly is the length of presaturation period (saturation may also occur in d2 as determined by satmode).

sspul='y' selects for trim(x)-trim(y) sequence at start of pulse sequence
(recommended)

4.14 TNDQCOSY—DQCOSY with Water Suppression

The tndqcosy macro sets up parameters for a DQCOSY experiment with transmitter solvent saturation only. The experiment assumes an on-resonance solvent (tof is at solvent position).

Applicability

DQCOSY requires a system with a linear amplifier on the observe channel and a T/R switch. It is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 141 for a discussion of DQCOSY.

Parameters

satmode determines when the saturation happens. satmode should be set analogously
to dm (e.g., satmode='yyn' or satmode='ynn').

satpwr is the power level during saturation period(s).

satdly is the length of presaturation period (saturation may also occur in d2 as determined by satmode).

sspul='y' selects for trim(x)-trim(y) sequence at start of pulse sequence
(recommended)

4.15 TNMQCOSY—MQCOSY with Water Suppression

The tnmqcosy macro sets up parameters for a multiple-quantum filtered COSY experiment with transmitter saturation.

Applicability

TNMQCOSY uses the hardware digital phaseshifter for the transmitter with direct synthesis rf. It uses the software small-angle phaseshifter for the transmitter with old-style rf. It is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 151 for a discussion of MQCOSY.

Parameters

pw is a 90° excitation pulse (at power level tpwr).

phase=1, 2 for hypercomplex phase-sensitive experiment or phase=3 for TPPI phasesensitive experiment. If phase=3, remember that sw1 must be set to twice the desired value.

satmode='ynn' saturates during relaxation delay.

satdly is the saturation time.

satpwr is the saturation power.

sspul='y' selects for trim(x)-trim(y) sequence at start of pulse sequence; sspul='n' selects a normal MQCOSY experiment.

Note: Earlier versions included the following sequence at the beginning of the pulse sequence: *homospoil - 90° pulse - homospoil*. This was used to eliminate both the DQ-like artifacts in the 2D spectrum and the oscillatory nature of the steady-state. This inclusion was selected if sspul='y'.

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

qlvl selects the quantum order for filtering (2, 3, etc.).

nt minimum is a multiple of 2*qlvl; nt maximum is a multiple of 8*qlvl.

References

Piantini, U.; Sorenson, O. W.; Ernst, R. R. J. Am. Chem. Soc. 1982, 104, 6800-6801.

Rance, M.; et al. BBRC, 1983, 117, 479-485.

4.16 TNNOESY—NOESY with Water Suppression

The tnnoesy macro sets up parameters for a 2D cross-relaxation experiment with transmitter saturation. It assumes an on-resonance solvent (tof is at solvent position)

Applicability

TNNOESY requires a system with a linear amplifier on the observe channel and a T/R switch. It is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 153 for a discussion of NOESY.

Parameters

satmode determines when the saturation happens. satmode should be set analogously
to dm (e.g., satmode='yyyn' or satmode='ynyn').

satpwr is the power level during saturation period(s).

satdly is the length of presaturation period (saturation may also occur in d2 and mix as determined by satmode).

sspul='y' selects for trim(x)-trim(y) sequence at start of pulse sequence.

4.17 TNROESY—ROESY with Water Suppression

The throesy macro sets up parameters for a rotating frame NOE experiment. Observe transmitter should be set at solvent position. Saturation, spin lock, and pulses all use the observe transmitter.

Applicability

TNROESY is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 156 for a discussion of ROESY.

Parameters

p1 is the 90° pulse on protons (power level at p11v1).

pw is a small (30°) pulse on protons (active only if ratio is greater than 0). If pw=0, pw is set to p1/3.

pllvl is the power level for the pl pulse.

tpwr is the power level for the spin lock pulse(s).

A cw lock is used if ratio is zero.

 $\verb|phase=1, 2 gives f_1 quadrature by the hypercomplex method (uses f_1 axial peak displacement). \verb|phase=3 gives f_1 quadrature by the TPPI method.$

mix is the mixing time.

sspul='y' selects for trim(x)-trim(y) sequence at start of pulse sequence.

```
rocomp='n' sets no resonance offset compensation
rocomp='y' sets resonance offset compensation (recommended).
```

satmode is the saturation mode. Use analogously to dm, for example,

satmode='nnn', satmode='ynn', or satmode='yyn' (recommended).

satdly is the length of saturation during relaxation delay.

satpwr is the power level for solvent saturation.

nt minimum is a multiple of 2; nt maximum is multiple of 8 (recommended)

d2corr is an empirical correction, in μ s, of d2 (dependent on effective field of spin lock, i.e., tpwr and ratio). It can be determined from the lpl and swl values from a properly phased spectrum by d2corr=(lpl*le6)/(360*swl). Note that the d2corr seems to be dependent on swl; however, it is independent of swl since changes in swl result in corresponding changes in lpl so that their ratio is constant.

Technique

The following procedure finds d2corr so that lp1 is 0, giving better baselines in f₁:

- 1. Run a TNROESY experiment with **d2corr** set either at 0 or at a value found previously. (nt and ni can be smaller, and the spectrum can be transformed early to do step 2.)
- 2. Phase the resulting spectrum in f₁. Determine lp1 and calculate d2corr from the relationship d2corr=(lp1*1e6)/(360*sw1).
- 3. Add this value to the value of d2corr used in step 1.
- 4. Rerun the experiment. **1p1** should be close to zero.
- 5. Note this value for any future experiment with the same value of tpwr and ratio.

Use of any method to make lpl=0 will result in a dc offset of F_1 slices. This should be removed by dc2d('fl') after the 2D transform. Enough noise should be left on the edges (in F_1) to permit this dc correction.

4.18 TNTOCSY—TOCSY with Water Suppression

The thtocsy macro sets up parameters for total correlation spectroscopy (also known as HOHAHA) using transmitter presaturation. It features "clean" TOCSY with optional windowing and MLEV16 + 60° spin lock. The transmitter must be positioned at the solvent frequency.

Applicability

TNTOCSY requires a T/R switch, and linear amplifiers and computer-controlled attenuators on the observe channel. It is not available on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* systems. See page 162 for a discussion of TOCSY.

Parameters

pw is a 90° pulse during mlev periods, at power level tpwr.

pl is a 90° excitation pulse, at power pllvl.

window is a clean-TOCSY window, in µs.

satdly is the length of presaturation.

satmode='yn' for presat control during relaxation delay only, satmode='yy' for presat control during both the "relaxation delay" and d2.

phase=1, 2 for hypercomplex phase-sensitive f_1 quadrature, or phase=3 for TPPI phase-sensitive f_1 quadrature.

sspul='y' gives the trim(x)-trim(y) sequence at beginning of d1 delay; sspul='n' gives the normal d1 delay.

trim is a spinlock trim pulse time (0.002 recommended).

mix is the mixing time (can be arrayed).

nt minimum is a multiple of 2, nt maximum is a multiple of 8 (recommended).

References

Bax, A.; Davis, D. G. J. Magn. Reson. 1985, 65, 355.

Levitt, M.; Freeman, R.; Frenkiel, T. J. Magn. Reson. 1982, 47, 328.

4.19 TOCSY—Total Correlation Spectroscopy

The tocsy macro sets up parameters for the TOCSY pulse sequence. Figure 36 is a diagram of this sequence. TOCSY is also known as the Homonuclear Hartmann-Hahn experiment (HOHAHA).




Figure 36. TOCSY Pulse Sequence

Applicability

TOCSY requires systems having linear amplifiers with computer-controlled attenuators on the observe transmitter channel. It is not available on the *GEMINI 2000*.

Parameters

p1 is the 90° pulse on the observed nucleus (at power equals p11v1)

pw is the 90° pulse on the observed nucleus (at power equals tpwr); note that tpwr is used to define the B_1 field for both the trim (trim) pulses and the MLEV-16 spin lock. The B_1 field at tpwr should be on the order of 15 ppm (in Hz).

trim is the trim pulse time (CW transmitter irradiation).

sl is MLEV-16 spin lock. mix is the duration of the spin lock, recommended mixing time is 30 to 80 ms. mix can be arrayed.

d1 is a relaxation delay (1 to 3 times the value of t_1).

phase=1, 2 (suggested value) for 2D hypercomplex data (States-Haberkorn method); phase=3 for 2D TPPI data. For phase=3, remember that sw1=2*sw (for homonuclear experiments) or that sw1 must be set to *twice* the desired value (for heteronuclear experiments).

sspul='y' activates a homospoil-90-homospoil sequence that precedes d1; This is used to achieve a less oscillatory steady-state for 2D experiments where the recycle time is shorter than t_1 .

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

presat is decoupler presaturation period using a decoupler power specified by dhp or dpwr.presat does not depend on dm but does depend on dmm and is activated as a part of d1 if presat is greater than 0. presat is available on *MERCURY-VX* and *MERCURY* systems with *GLIDE* (uppercase pulse sequences only).

nt is a multiple of 4 (minimum, phase=1, 2 or phase=3).

Phase Cycling

v1 is the phase for p1 pulse v13 is the phase for trim pulses v2 is the master phase for the MLEV-17 spin lock oph is a phase for receiver v1 = y -y -x x v13 = x x y y v2 = y y -x -x -y -y x x oph = y -y -x x y -y -x x

These phases are for phase=1. For phase=2, add 90° to v1. For phase=3, add $90^{*}(ix-1)^{\circ}$ to v1, where ix is the increment counter.

Technique

To set up and acquire, use the technique given in "NOESY—Nuclear Overhauser Effect Spectroscopy," page 153. Set mix to 0.060.

To process simple TOCSY, use the technique (except spectra in the first increment are *not* inverted) given in "NOESY—Nuclear Overhauser Effect Spectroscopy," page 153.

For a series of TOCSY experiments acquired in an interleaved fashion, each with a different mixing time, process the data with the wft2dac(mult1,mult2,...) command. Up to eight interleaved experiments can be processed, which means that mix can be given up to eight different values. Note that if mix is arrayed, phase=3 is necessary in order to be able to process the data with the wft2dac macro.

References

Levitt, M.; Freeman, R.; Frenkiel, T. J. Magn. Reson. 1982, 47, 328.

Bax, A.; Davis, D. J. Magn. Reson. 1985, 65, 355.

Related Macros

The following macro are used with TOCSY (with multiple mixing times):

- ftldac and ft2dac combine 2D FID matrices within the 2D Fourier transform framework.
- wftldac and wft2dac also combine 2D FID matrices within the 2D Fourier transform framework but with weighting added.

Refer to description of these macros in the *VNMR Command and Parameter Reference* for further information, including the arguments used with each macro.

4.20 TROESY—Transverse ROESY

The troesy macro sets up parameters for the TROESY pulse sequence, a transverse cross-relation experiment in a rotating frame.

Applicability

TROESY is not available on MERCURY-VX, MERCURY, or GEMINI 2000 systems.

Parameters

TROESY uses typical ROESY parameters. Note that ratio is not a parameter.

Set slpwr to the spinlock desired region.

slpw is the 90° pulse width at slpwr.

Reference

Shaka, et. al. J. Am. Chem. Soc. 1992, 114, 3157.

4.21 HCCHTOCSY Pulse Sequence

The hochtoosy macro sets up parameters for the HCCHTOCSY pulse sequence, used for sidechain assignments in fully ¹³C-enriched biopolymer.

Applicability

Not supplied with MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

coshape is the decoupler pattern for CO decoupling. fhfdwt1 is a flag to indicate half dwell start in t_1 . fhfdwt2 is a flag to indicate half dwell start in t_2 . tof is set in center of proton spectrum or on H₂O frequency if in H₂O. dof is set in center of the ¹³C aliphatic region. satfrq is the presaturation frequency. pwca is the 90° pulse width for the ¹³C nucleus. pwco is the 90° pulse width for the ¹³C carbonyl decoupling.

satpwr is the low-level ¹H transmitter power for presat.

pwcalvl is the power level for ¹³C pulses.

pwcolvl is the power level for C=O decoupling pulse.

jch is coupling for C-C (set to 40 Hz).

ncyc is the number of cycles through the dipsi loop.

trim is the trim pulse length, in seconds.

dipsipwr is the power level for ¹³C spin lock.

pl is the 90° pulse width for dipsi.

cycletime is an informational parameter (do not enter).

Technique

Optional ¹³C decoupling of carbonyl carbons during t2 uses a 180° pulse done through a shaped SLP (phase-ramped) pulse. The ¹³C transmitter is normally in aliphatic region (do not try for both aliphatic and aromatic hcchtocsy because the bandwidth is too large for good spinlock). The pulse is of length pwco at power pwcolvl and of name coshape.

- 1. Determine the length, nature and power for this 180 in an on-resonance calibration experiment.
- 2. Prepare this pulse by preparing a "template" pulse that has 5 times as many steps as width, in μ s.
- 3. Convolute this to prepare the SLP pulse. For example, if pwcolvl is 38, a sinc pulse does a 180° pulse in 250 μs when on-resonance; therefore, for the most accuracy you would need a sinc.RF file of 250*5 steps from shapelib.

For rectangular pulses, you can run **makehard**(number_steps) to do this. For complex shapes, you can generate it out of **pulsetool** or by a separate program.

If the sinc pulse has the name sincl250.RF, and the distance from the ${}^{13}C$ transmitter is, say, 15000 Hz, enter

convolute('sinc1250','co180_1250us_+15000',250,15000)

to prepare the new shape. This gives you a shape that is used by decshaped_pulse to do a 180 on the carbonyl region when dof is positioned in the aliphatic region. Use pulsetool to verify that your shape is correct.

It is not necessary to do carbonyl decoupling. The only visible effect is CC splitting in f_2 for alpha carbons. Setting pwco=0 is not desirable.

4. Set **dm='nny'** to do ¹³C decoupling during acquisition. One dipsi-3 cycle is 217.33*p1.

For example, if p1 is 36 µs, a single dipsi cycle is 7.8 ms.

5. Set **phase=1,2** and **phase2=1,2** for hypercomplex in t₁ and t₂.

Typical acquisition times are 28 ms for t_1 , 10 ms for t_2 , and 47 ms for t_3 , with 128 complex points for t_1 , 32 complex for t_2 , and 512 real for t_3 .

4.22 HMQCTOCSY Pulse Sequence

The hmqctocsy macro sets up parameters for HMQCTOCSY pulse sequence, with an option to null or invert the direct responses.

Applicability

Not supplied with MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

p1 is the 90° pulse during the clean-TOCSY period.

pllvl is the power level for pulse pl.

window is the off-time for clean-TOCSY. If p1=25, window=50 (in µs) is typical.

pw is the 1 H 90° pulse during the pulse sequence outside of the mlev period.

tpwr is the power level for pw.

mix and trim are the isotropic mixing period, 0.020 and 0.002 are typical values.

mult controls the direct responses:

mult=0 gives a normal HMQCTOCSY (direct and relayed in-phase).
mult=1 nulls the direct responses. dm='nnnn' is best for this.

mult=2 inverts the direct responses.

Technique

Remember to enter **dps** to check things before entering **go**.

4.23 HMQC-TOCSY 3D Pulse Sequence

The hmqctocsy macro sets up parameters for HMQC-TOCSY 3D pulse sequence with presaturation option, written in hypercomplex phase-sensitive mode only. Figure 37 is a diagram of the sequence.

Figure 37. HMQC-TOCSY 3D Pulse Sequence

Applicability

Not supplied with MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

d2 is first evolution time.

d3 is second evolution time.

mix is the TOCSY mixing time.

pwxlvl is the power level for X pulses.

pwx is a 90° X pulse.

j1xh is a X-H coupling constant.

dpwr is the power level for X decoupling.

tpwr is the power level for H pulses.

pw is a 90° H hard pulse.

slpwr is the power level for spinlock.

slpw is a 90° H pulse for mlev17.

trim is a trim pulse preceding mlev17.

phase=1, 2 gives hypercomplex (t1) acquisition; ni is number of t1 increments. phase2=1, 2 gives hypercomplex (t2) acquisition; ni2 is number of t2 increments

satflg='y' is presaturation during satdly.

satfrq is the presaturation frequency.

satdly is the saturation time during the relaxation period.

satpwr is the saturation power for all periods of presaturation with the transmitter.

wdwfctr is multiplication "window" factor of slpw.

nullflg is TANGO nullflg flag for protons not attached to X.

hmqcflg='n' turns off HMQC part of the sequence.

Technique

For F1 x F3, use wft2d(1,0,0,0,0,0,0,-1,0) (i.e., wft2da).

For F2 x F3, use wft2d(1,0,0,0,0,0,-1,0) if hmqcflg='n' (i.e., wft2da) or wft2d(1,0,0,0,0,0,1,0) if hmqcflg='y'.

For 3D processing: create 3D coefficients using **make3dcoef** macro: **flcoef='1 0 0 0 0 -1 0'** and **f2coef='1 0 0 0 0 1 0'**.

4.24 HSQC-TOCSY 3D Pulse Sequence

The hsqctoxySE macro sets up parameters for HSQC-TOCSY 3D pulse sequence with many features: HSQC and TOCSY are each "sensitivity enhanced" as per Rance et. al. The sequence features a dipsi or mlev-17 spinlock option (dipsiflg flag). F1 dimension can be band selected with shaped pulse. (Use of shaped pulse introduces a small phase error along F1. rp1 and lp1 values can be obtained by examining first np x ni plane) The sequence has a gradient to kill unwanted signals during H-X INEPT. Water presaturation can be on- or off- resonance.

Applicability

Not supplied with MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

phase=1,2,3,4 and phase2=1,2,3,4.

dm= 'nnyny' (on during t2 and t3).

flcoef and f2coef are processing parameters used by the make3dcoef macro.

Chapter 5. Indirect Detection Experiments

Sections in this chapter:

- 5.1 "Requirements for Indirect Detection Experiments," this page
- 5.2 "The Basic HMQC Experiment," page 175
- 5.3 "Phase-Sensitive Aspects of the Sequence," page 179
- 5.4 "Cancellation Efficiency," page 179
- 5.5 "Pros and Cons of Decoupling," page 180
- 5.6 "Specifications Testing," page 181
- 5.7 "Using the HMQC and HMQCR Sequences," page 182
- 5.8 "Recabling Single-Broadband Systems," page 183
- 5.9 "Recabling Dual-Broadband Systems," page 183
- 5.10 "Filters for Indirect Detection," page 184
- 5.11 "Tuning the Probe in the Reverse Mode," page 184
- 5.12 "Controlling Transmitter Power in the Reverse Mode," page 185
- 5.13 "Indirect Detection Calibration," page 185
- 5.14 "Typical Experimental Protocol for HMQC Experiments," page 192
- 5.15 "Differences for ¹⁵N Indirect Detection," page 197
- 5.16 "HSQC Experiment," page 197

This chapter describes indirect detection experiments, also known as *heteronuclear multiple-quantum coherence* (HMQC) experiments. Indirect detection experiments show correlations between heteronuclei while detecting high-sensitivity protons. HMQC differs from the more traditional heteronuclear correlation techniques that detect the low-sensitivity heteronucleus (for example, 13 C or 15 N).

5.1 Requirements for Indirect Detection Experiments

Indirect detection experiments have three basic requirements:

- A probe with a proton channel and an X-nucleus channel.
- An rf (radio-frequency) system capable of generating pulses on ¹H and X and possibly decoupling X while observing ¹H.
- One or more pulse sequences that perform the relevant experiments.

Each of these requirements is described in more detail in the following sections.

Probes

The most commonly used probes for indirect detection experiments are the "indirect detection" probes, such as the Varian Indirect•nmrTM probe and triple-resonance probes, such as the Varian Triple•nmrTM probe. Indirect detection probes have a ¹H coil and an X-nucleus coil with the ¹H coil positioned closer to the sample for the highest possible sensitivity of the observed nucleus.

Normal "broadband" probes similarly have a ¹H coil and an X-nucleus coil and can be used for indirect detection. But broadband probes have significantly lower proton sensitivity (about half that of indirect detection probes) and so are not optimum for indirect detection experiments. Nevertheless, broadband probes usually provide some sensitivity improvement over direct detection heteronuclear correlation experiments. Four-nucleus and "Switchable" probes also have a ¹H coil and an X-nucleus coil, with the X coil closer to the sample, and can satisfy the needs for indirect detection experiments.

Probe types are not discussed further in this chapter. See the manual *Getting Started* and the probe installation manual for tuning instructions. When connecting cables to the probe, ignore words like "observe" and "decouple" and think of ¹H (for observe) and X (for decouple) for connections.

RF System

The rf (radio frequency) part of a spectrometer has at least four relevant aspects to indirect detection that exist in various combinations on the various systems:

- · Generation of rf for the transmitter and receiver
- Amplifiers
- Filters
- Routing of rf signals

RF Generation

The first aspect is the most basic—generation of rf for the transmitter and the receiver (specifically the local oscillator or L.O. frequency). Most high-field systems, and a few specially-configured lower-field systems, are "dual-broadband" systems. Because both of the rf channels in the system are broadband, performing an indirect detection experiment involves setting the transmitter nucleus parameter to proton and the decoupler nucleus to the appropriate X-nucleus; the rf is generated correctly.

The remaining systems, including virtually all 200-, 300-, and 400-MHz systems, are "single-broadband" systems, in which the "observe" channel is broadband but the "decoupler" channel can generate only ¹H frequencies. These systems, therefore, must operate in a "reverse" configuration, in which the observe channel (controlled by the transmitter nucleus parameter) generates the X-nucleus decoupling frequency, while the decoupler channel (controlled by the decoupler nucleus frequency) generates the ¹H frequency. Insofar as generating the transmitter frequencies, this process is straightforward, but to allow the decoupler channel to be used as the reference frequency for detection (that is, to be the observe channel), the L.O. frequency must first be generated from the "decoupler" board and then properly connected to the receiver section of the spectrometer. Single-broadband systems may involve changing the cabling, which is discussed later in this manual.

Amplifiers

Amplifiers represent the second aspect of rf relevant to indirect detection. There are two configurations: systems with linear amplifiers and systems with class C amplifiers. All UNITY *INOVA*, *MERCURY-VX*, *MERCURY*, UNITY *plus*, UNITY, broadband *GEMINI 2000*, VXR-500, and late-model VXR-S systems have linear amplifiers. Such amplifiers have computer-controlled output and can sustain power levels appropriate for X-nucleus decoupling indefinitely (*MERCURY-VX*, *MERCURY*, and *GEMINI 2000* broadband systems cannot do presaturation but can do X-nucleus decoupling).

Systems without linear amplifiers have class C amplifiers with power output that cannot be switched rapidly under computer control. Class C amplifiers preclude experiments (such as presaturation) that require power level switching during the sequence. Furthermore, on systems with class C amplifiers, the X-nucleus amplifier cannot be left on for long periods of time at the power levels appropriate for X-nucleus decoupling. If X-nucleus decoupling during acquisition is not desired, these systems pose no problem. If X-nucleus decoupling is desired, however, a limiting acquisition time of 50 to 100 ms is required unless some hardware is modified.

GARP modulation is available on UNITY *INOVA* and UNITY *plus* systems and should be used for X decoupling.

When WALTZ decoupling is used, the maximum power level for decoupling is the level that provides an rf field strength (in Hz) comparable to half the range of expected X shifts. The normal spread of protonated carbons is 150 ppm, which is 15 kHz on a 400-MHz system, and, consequently, a ¹³C 90° pulse of 25 μ s (corresponding to an rf field strength of 8 kHz) is adequate. The somewhat long proton pulses on broadband and switchable probes does not seem to present a problem because indirect detection experiments demand no more proton pulse power than DEPT or HETCOR.

Filters

Filters are the third aspect of rf that is relevant to indirect detection experiments. Filters may be needed on the transmitter, receiver, decoupler, and lock channels. Filters are part of the probe kit shipped with each indirect detection probe that Varian sells.

RF Signal Routing

The routing of rf signals is the fourth and final relevant aspect of rf. Here the major difference between systems arises from the "switchable relays" in the magnet leg. The relays are designed to allow the user to switch between ¹H observation and X-nucleus observation when using a normal switchable probe. Most systems are equipped with these relays, but some are not, and cable routing differs between these two configurations. An additional difference arises on systems with linear amplifiers (mostly UNITY*INOVA*, UNITY*plus*, and UNITY systems) that also have a relay intended for switching the L.O. between two sources (the "observe" and "decouple" slots). The control signal for this relay is wired on all UNITY*Plus*, and UNITY*plus*, and UNITY systems, while the rf cabling to the relay is installed on most, if not all, systems that have a factory-installed L.O. SELECT switch.

All UNITY *INOVA* and UNITY *plus* spectrometers and all single-broadband UNITY spectrometers produced after October 1990 have a configuration that allows full computer-controlled switching between direct and indirect detection. These systems have a relay in the magnet leg that accomplishes the relevant rf signal routing and computer-controlled L.O. selection, instead of the L.O. SELECT switch described above. These systems are

easily identified because most of the connectors in the magnet leg, including those attached to the relay, are N-type connectors.

MERCURY-VX, *MERCURY*, and *GEMINI 2000* broadband systems have relays to enable full computer-controlled switching between direct and indirect detection. The *GEMINI 2000* 1 H/ 13 C system does not do indirect detection.

There is a "catch" with this configuration—the filters used for indirect detection tend to degrade specifications approximately 10% in terms of longer pulse widths and lower signal/noise. The user thus faces a classic trade-off of performance (manually insert filters only when needed but achieve better specs) versus convenience (leave filters in place continuously and achieve worse specs). The convenience factor, of course, is nonexistent if the instrument does anything other than ¹³C and ¹H, because one cannot leave the ¹³C bandpass filter in place on the X line while doing ³¹P, ¹⁵N, or anything else. All standard specifications are given with the indirect detection filters *not* in place.

Pulse Sequences

There are three different data systems that users of UNITY INOVA, MERCURY-VX, MERCURY, UNITY plus, UNITY, GEMINI 2000, Gemini, VXR, and XL spectrometers can have on their spectrometer systems—the Sun-based data system used with the UNITY INOVA, MERCURY-VX, MERCURY, UNITY plus, UNITY, GEMINI 2000, and VXR-S (the Sun-based VXR, originally named the VXR-5000 data system), the data system used with VXR-4000 and Gemini spectrometers, and the V77-200-based data system used with XL spectrometers. This chapter covers only the Sun-based data system.

Combining the Sun-based data system with the two different rf generation types discussed previously (single- and dual-broadband) gives two fundamentally different system configurations and hence two sets of pulse sequences:

- For Sun-based dual-broadband systems, the standard system pulse sequence library contains the macros hmqc and hmqcr for setting up parameters for the HMQC and HMQCR pulse sequences, respectively. An option in hmqc selects another pulse sequence of interest, HMBC (Heteronuclear Multiple-Bond Coherence, the long-range sequence). Another version of HMQC, which uses decoupler 2 for the X nucleus, is also available. A compiled version of each of these sequences is in the library and ready to use; documentation on using these sequences is contained in the pulse sequence listings produced by running the macros.
- The HMQC pulse sequence can be used on single-broadband UNITYINOVA, MERCURY-VX, MERCURY, GEMINI 2000, and UNITYplus. Beginning with VNMR version 5.1 (5.2F for GEMINI 2000), channels 1 and 2 are automatically swapped, if tn='H1' and dn='X' with channel 2 used for observe. On UNITY and VXR-S single-broadband systems, only HMQCR can be used.
- On the UNITY INOVA, MERCURY-VX, MERCURY, GEMINI 2000, and UNITY plus, the sequence S2PUL pulses the decoupler channel in the "reverse" configuration if tn='H1' and dn='X', for testing and calibration purposes. On the UNITY and VXR-S, the pulse sequence S2PULR and its corresponding macro s2pulr are provided for this purpose.

HMQC Pulse Sequence

The hmqc<(isotope)> macro sets up parameters for a HMQC (heteronuclear multiplequantum coherence) pulse sequence. The optional isotope argument is the isotope number of the heteronucleus of interest, for example, hmqc(1) for ¹H (the default is ¹³C). Figure 38 is a diagram of this sequence. The first 2 pwx pulse on the X heteronucleus is a composite 180 consisting of 90(v9) - 180(v1) - 90(v9).




Figure 38. HMQC Pulse Sequence with null<>0 and mbond='n'

Parameters

pw is a 90° pulse on the observed nucleus (protons) at power equal to tpwr.

pwx is a 90° pulse on the heteronucleus at power equal to pwxlvl.

dpwr is the decoupler power level for broadband X-decouphng.

dmf sets the modulation frequency $(4*\gamma B_1)$ at decoupler power (dpwr).

dmm is decoupler modulation mode. For UNITY INOVA, MERCURY-VX, MERCURY, and UNITY plus, dmm='ccg' is recommended; otherwise, set dmm to 'ccw' or 'ccp'.

dm= 'nny' activates heteronuclear broadband decoupling (recommended) during acquisition. Note that dm can be set to either 'nnn' or 'nny', and that the duty cycle for the decoupler should be less than 20%.

j is the average scalar coupling constant between the protons and the heteronucleus (usually one-bond constants). j is 140 for ¹³C or 90 for ¹⁵N. The time Δ , shown in Figure 38, is calculated as 1/2j.

null is a WEFT-like delay used to improve the suppression of the protons connected to ${}^{12}C$ (and not to ${}^{13}C$) that have been inverted by the preceding BIRD pulse. Try a null value of 0.3 for ${}^{13}C$, 1.0 for ${}^{15}N$, and 0 for macromolecules. To optimize, set ss=-8 and array null with nt=1 and phase=1. This selects the value of null that best minimizes the sample's signals (typically 0.2 to 0.7 seconds). If null is set to 0, the BIRD element is omitted from the pulse sequence.

at is the acquisition time (t_2 period).

ni is the number of t_1 increments (set up with default values for either ${}^{13}C$ or ${}^{15}N$).

ss is the number of complete executions of the pulse sequence not accompanied by data collection prior to the acquisition of the real data: if ss is positive, ss steady-state pulses are applied on the first increment only; if ss is negative, -ss transients are applied at the start of each increment.

nt is a multiple of 4 (minimum) or multiple of 8 (recommended).

phase=1, 2 (2D hypercomplex data with hypercomplex-TPPI method) or phase=3 2D TPPI data). phase=1, 2 is suggested. For phase=3, remember that hmqc sets sw1 to *twice* the desired value for heteronuclear experiments.

satflg='yn' gives presaturation during satdly, and satflg='yy' gives
presaturation during satdly and null (not on MERCURY-VX, MERCURY, and
GEMINI 2000).

satfrq=x is the presaturation frequency (using the transmitter), satdly is the length of saturation time during the relaxation period (immediately after dl), satpwr is the power level for presaturation using the transmitter (not on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*).

hs='yn' gives a homospoil pulse at beginning of d1 (length=hst). hs='yy' gives a homospoil pulse at beginning of both d1 and null. For *GEMINI 2000*, homospoil is not supported and hs must be 'nn'.

taumb is a fixed delay associated with the multiple-bond HMQC experiment (taumb=0.055 is recommended).

mbond = 'n' is a normal HMQC experiment. mbond = 'y' is a multiple-bond HMQC experiment (HMBC).

To run HMBC (mbond='y'): (1) set null=0, otherwise, only protons that are both long-range and short-range (one-bond) coupled to a given heteronucleus (13 C, for example) will not be suppressed, (2) set dm='nnn', (3) set taumb, and (4) run the single-bond (HMQC) and multiple-bond (HMBC) experiments with phase=1, 2 or phase=3.

Phase Cycling

The phase cycling is the following:

v1, v2, v3, v4, v5, v9 are phases for pulses. oph is the phase for receiver.

These phases are for phase=1. For phase=2, add 90° to v3. For phase=3, add $90^{\circ}(ix - 1)^{\circ}$ to v3, where ix is the increment counter.

Technique

The usual setup is to place a ¹H bandpass filter between the observe port on the probe and the ${}^{1}\text{H}/{}^{19}\text{F}$ preamplifier, and to place a 250-MHz lowpass LC filter and either a ${}^{13}\text{C}$ bandpass or a ${}^{15}\text{N}$ bandpass filter in the decoupler line just before the probe connection.

The experiment should be performed non-spinning and with VT regulation.

HMQCR Pulse Sequence

The hmqcr macro sets up a HMQCR (heteronuclear multiple-quantum coherence in reverse configuration) pulse sequence. This sequence is normally used in systems with a ¹H only decoupler, *without automatic channel reversing*. For best results, perform the HMQCR experiment non-spinning. Figure 39 is a diagram of the sequence.

Parameters

j is the X-H spin coupling constant.




Figure 39. HMQCR Pulse Sequence

null is nulling time for protons not attached to X (if 0, nulling section of sequence is omitted).

dn and dof set ¹H (observe) frequency.

pw is a 90° pulse for protons.

tn and tof set X nucleus frequency.

pwx is a 90° pulse for X nucleus at power level pwxlvl.

dm='nny' gives software WALTZ-4 decoupling of X during acquisition; dm='nnn' gives coupled spectrum.

pwxw is a 90° pulse for X nucleus decoupling at power level dpwr.

nt=8 for best results; multiple of 2 is minimum.

ss is the number of steady-state transients. If *ss* less than 0, then -ss transients are performed before *each* increment (recommended for calibration/setup experiments, not for 2D).

phase=1, 2 for hypercomplex; phase=3 for TPPI.

5.2 The Basic HMQC Experiment

The essence of the HMQC experiment is the cancellation or elimination of the signals from protons attached to ¹²C, leaving only signals from protons attached to ¹³C, contributing to a ¹³C–¹H chemical shift correlation spectrum. There are three basic, independent mechanisms to generate this discrimination: the spin-echo difference experiment, BIRD (Bilinear Rotation Decoupling) nulling, and presaturation.

Spin-Echo Difference Experiment

The heart of the HMQC sequence can be reduced to a heteronuclear spin-echo difference experiment that looks like Figure 40.




Figure 40. Heteronuclear Spin-Echo Difference Experiment

In Figure 41, *a*, *b*, and *c* represent the protons attached to carbons, where *a* are the protons attached to up-¹³C, *b* are protons attached to 12 C, and *c* are protons attached to down-¹³C. Assume that we are at the resonance frequency of the protons attached to the 12 Cs. In the rotating form, the following steps (shown in Figure 41) occur:




Figure 41. HMQC Pulse Sequence, Showing Movement of Attached Protons

- 1. The first proton 90° pulse places all protons along the y axis.
- 2. After a time $\Delta = 1/(2J)$, the *b* protons are still along the *y* axis, but the *a* protons are along the -x axis and the *c* protons are along the +x axis.

- 3. Next, the 180°_{X} proton pulse places the *b* protons along the -y axis but does not affect the *a* and *c* protons.
- 4. The next pulse has the following effect:
 - a. The $90^{\circ}_X 90^{\circ}_{-X}$ carbon pulse is effectively a null pulse. All rotational directions are maintained.
 - b. The $90^{\circ}_{X}90^{\circ}_{X}$ (= 180°_{X}) carbon pulse reverses the ¹³C, which makes the *a* protons attach to the down-¹³C and the *c* protons attach to the up-¹³C, essentially reversing their rotational direction.
- 5. After another period $\Delta = 1/(2J)$, the following occurs:
 - a. The *a*, *b*, and *c* protons are refocused along the -y axis.
 - b. The *b* protons are still along the -y axis, and the *a* and *c* protons are refocused along the +y axis.

Subtracting the signal resulting from step 5b and 5a, by changing the receiver phase oph, results in cancellation of the *b* protons, while the signal for the *a* and *c* protons doubles.

To create a 2D experiment with information about heteronuclear chemical shifts, we introduce an evolution time t_1 that occurs between the two X-nucleus 90° pulses, as shown in Figure 42.




Figure 42. Evolution Time Added Between X-Nucleus Pulses

In this 2D experiment, which is now a full HMQC experiment, protons attached to ${}^{12}C$ show no different behavior and are still cancelled after two scans. For the ${}^{13}C$ nuclei, however, whether they experience a 180° pulse, a 0° pulse, or something in between, depends on the time between the two 90° pulses and their rate of precession during that time (i.e., their chemical shift). Therefore, this experiment produces a modulation of the intensity of the ${}^{13}C$ -bound protons, and the Fourier transform of that modulation yields the chemical shift of the ${}^{13}C$ bound to that proton.

In this way we detect ¹³C chemical shifts with the intensity of protons, and simultaneously we obtain a correlation of the ¹³C and ¹H chemical shifts. Appropriate variations of the experiment produce long-range coupling information.

BIRD Nulling

The second (optional) type of cancellation that can occur during an HMQC sequence is the so-called BIRD (Bilinear Rotation Decoupling) pulse nulling effect (Summers, Marzilli, and Bax, *JACS*, **1986**, *108*, 4285). A particular sequence of the BIRD pulse, three pulses on the ¹H channel and one on the X channel, inverts the z-magnetization of protons bound

to ¹²C and leaves the z-magnetization of protons bound to ¹³C unaffected. The full sequence is illustrated in Figure 43, where $\Delta = 1/2J_{XH}$.




Figure 43. HMQC with BIRD Pulse Nulling Effects

After the BIRD pulse, a variable waiting period (τ in Figure 43) is inserted, allowing the ¹²C-bound protons to relax back to equilibrium. If τ is adjusted so that the ¹²C-bound protons are approximately at a null, then when the remainder of the pulse sequence (the normal HMQC sequence) is executed, cancellation of the ¹²C-bound protons is enhanced (since those protons had very little magnetization at the start of the HMQC sequence). Obviously, not all protons will have the same relaxation time, so the choice of τ must be a compromise; generally, unless only one proton is involved, the additional suppression from the BIRD nulling will be a factor of two to five.

For systems that exhibit a negative NOE, such as macromolecules, cross-relaxation between the inverted protons on 12 C and the noninverted protons on 13 C will decrease the intensity of the desired proton signal. The extent of this decrease can vary between 0% and 100%. For this reason, omission of the BIRD part of the sequence is advised for macromolecules.

BIRD pulse nulling is also not possible when long-range indirect detection experiments (Heteronuclear Multiple-Bond Coherence, or HMBC) are performed. In this case, protons that have long-range couplings to ¹³C are directly bonded to ¹²C (99% of them, anyway) so that BIRD pulse nulling would lose all intensity in the protons of interest.

Transmitter Presaturation for High-Dynamic Range Signals

When high-dynamic range situations, such as observing signals in H_2O , are involved, HMQC phase cycling and/or BIRD pulse nulling may be insufficient to produce cancellation of the large proton signals. For this reason a third mechanism, presaturation, may be necessary. Since one channel of the instrument is set to an X-nucleus like ¹³C or ¹⁵N, this presaturation must be accomplished with the other channel; that is, the same channel that will be applying observe pulses to the protons. During one or two different periods of the sequence (during the initial delay and during the τ delay), a change in power level and possibly frequency may be appropriate in order to perform the presaturation.

In the pulse sequences discussed in this chapter, presaturation is *only* supported on UNITY*INOVA* and UNITY*plus* systems running VNMR 5.1 or later, and on dual-broadband UNITY and VXR-S systems (presaturation is not supported on fixed-frequency decoupler UNITY, VXR-S, *MERCURY-VX*, *MERCURY*, or *GEMINI 2000* systems). With appropriate pulse sequence modification, single-broadband systems with linear amplifiers could perform these experiments.

5.3 Phase-Sensitive Aspects of the Sequence

The parameter phase, as in other phase-sensitive 2D experiments, controls the f_1 phase detection. For 1D setup experiments or a 2D experiment without quadrature detection in f_1 , use phase=1. For a normal 2D experiment using the States-Haberkorn-Ruben (hypercomplex) method, use phase=1, 2. To acquire data with TPPI, use phase=3, and make sure sw1 is twice the expected range.

The original HMQC experiments, including the experiments found in earlier software distributions from Varian (VNMR versions 2.2, 2.3, 3.1, and 3.2), were quite subject to axial artifacts (at $f_1=0$) when performed using the hypercomplex method. TPPI is preferable because the axial artifacts are shifted to the edge of the spectral region of interest (instead of falling in the center of the region of interest).

A newer technique was invented (Marion, D.; Ikura, M.; Tschudin,R.; Bax, A. J. Magn. Reson. **1989**, *85*, 393) involving only a change of phase cycling that shifts the axial artifacts in a hypercomplex experiment to the edge of the spectrum, giving the hypercomplex version the benefit of TPPI with none of the disadvantages. We have termed this technique, which is applicable to all 2D experiments, the name FAD, for "F1 Axial Displacement." It is also referred to as "States-TPPI." The hmqc and hmqcr macros in VNMR version 4.1 and later include FAD. Once implemented, use of the technique is totally transparent—just perform a standard hypercomplex experiment with phase=1, 2.

5.4 Cancellation Efficiency

Because indirect detection experiments involve cancellation of non-¹³C-bound protons that are two orders of magnitude more intense (assuming unlabeled compounds), cancellation efficiency is critical. Cancellation efficiency, in turn, depends on the fundamental stability of the system rf and the reproducibility of anything else that can affect the signal. While stability is fixed by the instrumentation, you can control a number of operating conditions that can influence the quality of any cancellation experiment (NOE difference experiments are another good example). Some of these conditions are discussed here, roughly in order of importance:

- · Run experiments non-spinning. This is a must.
- Use the highest lock power at which the lock is stable (be sure to shim with a nonsaturating level, however) and keep the lock gain as low as possible, sufficient only to be sure that you don't lose lock during the experiment.
- Use a ²H band-pass filter in the lock line. Interference between X-nucleus decoupling or even X-nucleus pulses can affect the lock and cause field instabilities, limiting the ability to perform cancellation experiments (if you don't have such a filter and want to prove this to yourself, try a short-term experiment in the unlocked mode).
- Use VT regulation, even at room temperature. Large changes in temperature of the environment can affect the VT gas stream. The frequency of peaks in the spectrum and of the lock resonance (which affects all peaks) is temperature-sensitive to some degree. Shimming may also change if the probe temperature varies, which can affect the lineshape.
- Be sure the system is in thermal equilibrium. If you are performing experiments with X-nucleus decoupling, you are applying large amounts of power to the system, which is almost certain to change the temperature of the probe, the sample, or both, even when you are performing VT regulation. The best way to ensure thermal equilibrium is to set up a "dummy" experiment with *identical* conditions (in terms of duty cycle) to your actual experiment, but which runs for perhaps several minutes (easily

accomplished by setting *ni* to a small number). Now, if you queue your real experiment to follow the dummy one, the sample and probe are properly equilibrated.

- Be sure the system is in an NMR steady-state by using steady-state pulses.
- Use a large value of nt. Cancellation improves with larger nt, so the relevant cancellation is that which occurs at nt comparable to what you will be using in an indirect detection experiment (16 to 1024). Do not expect perfect results with nt=2.
- Minimize floor vibration. Where this cannot be fixed by spectrometer placement, an antivibration system should be installed.
- Use a moderate flow of body air through the probe to eliminate "rattling" from turbulent flow.
- Use lengthened pulses (attenuated rf) if you have a rise time or phase glitch problem.

Before beginning an HMQC experiment, you should assess the quality of your reproducibility by performing some simple difference experiments. The standard S2PUL pulse sequence is a good one to use for this purpose. The first pulse of S2PUL, controlled by p1, is held at a constant phase, while the receiver varies in phase. Thus, after four scans with p1 set to the 90° value, pw=0, no signal should be seen. This can be compared to four scans with pw set to the 90° value, p1=0, which produces a full signal. Taking the ratio of these two spectra gives a concrete measurement of your cancellation efficiency, while repeating the null spectrum a number of times gives a measure of the reproducibility of the cancellation. Use this test to assess the value of the various steps and modifications described above, or of other differences (for example, the relative cancellation efficiency of experiments with and without X-nucleus decoupling).

Because rf stability is an issue and rf generation on single-broadband systems differs between observe and decoupler channels, on systems prior to the UNITY*plus* you should repeat the same experiment using S2PULR to pulse the decoupler channel (on *MERCURY-VX, MERCURY,* and *GEMINI 2000* running VNMR 5.2F or later, and on UNITY*INOVA* and UNITY*plus* running VNMR 5.1 or later, use S2PUL with tn='H1' and dn=X instead of both tn and dn set to 'H1'). While the overall cancellation ratio depends on all factors (vibration, lock channel, etc.), the ratio of the cancellation achieved gives you a measure of the relative performance of your decoupler channel in comparison to the observe channel.

5.5 Pros and Cons of Decoupling

It may seem that decoupling of X during acquisition is always better. The advantages are clear—the spectrum is less crowded, with only half as many peaks, and each peak has twice the sensitivity. This is, however, a simplistic view.

The disadvantage of X-nucleus decoupling stems from the need to use large (up to 8 kHz) decoupling fields. This high power can cause significant heating, particularly in lossy samples. As a consequence of sample heating, experiments with X-nucleus decoupling are generally limited to relatively short acquisition times, which in turn may produce less resolution in f_2 as well as less sensitivity for molecules with long T_2 . Furthermore, the heating that does occur frequently produces worse cancellation efficiency. And finally, to prevent the buildup of heat in the sample, the duty cycle of the experiment may need to be limited to 10 to 20%, again possibly reducing sensitivity. For all these reasons, experiments performed without X-nucleus decoupling are perfectly reasonable, and may well be preferable.

If X-nucleus decoupling is desired, it is important to avoid sample heating. This form of sample heating can be non-uniform within the sample and can cause microconvection,

producing poor cancellation. Keep the acquisition time short and the overall duty cycle less than 20%. To lower the decoupler power to tolerable levels, it may be necessary to add a fixed 6-dB attenuator to the X-nucleus channel on systems in which that power is not under computer control; this can have the unavoidable consequence of lengthening the pulse widths of the X-nucleus pulses.

Dual-broadband and *GEMINI 2000* broadband systems have no problem performing modulated decoupling, because the decoupling is being performed by the normal spectrometer decoupling channel. On single-channel broadband UNITY and VXR-S systems, however, the decoupling is being performed by the normal observe channel, and the standard modulation (WALTZ, for example) is not present.

In the sequences described here, broadband decoupling is achieved by using the acquisition computer to provide WALTZ-4 modulation of the X-nucleus channel through explicit software control. This requires hardware looping capability and explicit acquisition, which UNITY and VXR^{-S} systems support. It also imposes some limitation on spectral widths and pulse widths, since the WALTZ-4 sequence (whose length is 6*pw90) must fit in between successive data point samplings (which occur at time intervals of 1/sw).

GARP modulation is available on UNITY INOVA, MERCURY-VX, MERCURY, and UNITY plus systems and should be used for X decoupling.

5.6 Specifications Testing

Before performing any calibrations or specifications testing, connect the cables to the system as if it were a standard broadband probe, and tune both ¹H and X-nucleus channels of the probe just as you would a broadband probe, using tuning procedures described in the indirect detection probe section of the probe installation manual. Now using the standard samples and standard tests (see the acceptance tests procedures manual for the system), measure the pulse width, lineshape, linewidth, spinning sidebands, and ¹H signal/noise.

Note: All specifications on Varian indirect detection probes are guaranteed in the "normal" (i.e., direct observation) configuration. Changes in cabling, filters, and signal routing for indirect detection experiments may affect performance.

Before proceeding for the first time to indirect detection experiments on single-channel broadband systems, one basic test is important: verify that the local oscillator (L.O.) signal from the decoupler board is both enabled and of sufficient amplitude and quality to produce an acceptable NMR signal. The steps below show how to do so:

- 1. Measure the S/N of the system using any sample (the standard S/N test sample is useful but certainly not required).
- 2. Use s2pul with tn='H1' and dn=X (X can be 'C13') to switch to the "reverse" mode, changing the cabling if necessary (described later in this chapter) for UNITY INOVA and UNITY plus running VNMR 5.1 or later, and on MERCURY-VX, MERCURY, and GEMINI 2000 running VNMR 5.2F or later. Use S2PULR on single broadband UNITY and VXR-S.
- 3. Recalibrate the 90° pulse width, and check the setting of the gain parameter to ensure that you are filling the ADC with the signal (using ddff(1) to digitally examine the FID).
- 4. Measure the S/N again in the reverse configuration. The measured value should be *approximately* the same as in the normal mode. Values in the range of $\pm 15\%$ from the "normal" value seem typical. If the S/N drops more than 15%, the L.O. signal on the decoupler board probably needs to be adjusted.

5.7 Using the HMQC and HMQCR Sequences

Figure 44 shows again a diagram of the basic HMQC sequence. The parameters that correspond to this diagram depend on which sequence is used—HMQC or HMQCR. Table 32 lists the associated parameters.




Figure 44. Basic HMQC Pulse Sequence

Parameter	HMQC	HMQCR
¹ H 90° pulse	рw	pw
¹ H 180° pulse	derived from pw	derived from pw
¹ H amplifier power (if appropriate)	tpwr	dpwr
¹ H frequency	tn, tof	dn, dof
¹ H spectral width	SW	SW
X 90° pulse	pwx	pwx
X 90° pulse for WALTZ decoupling	1/(4¥dmf)	pwxw
X 180° pulse	derived from pwx	derived from pwz
X amplifier power for pulses (if appropriate)	pwxlvl	pwxlvl
X amp power for decoupling (if appropriate)	dpwr	tpwr
X frequency	dn, dof	tn, tof
X spectral width	swl	swl
Δ delay	1/(2j) [if j=0, D=0]	1/(2j)
τ delay for BIRD nulling (if null=0, entire BIRD sequence is skipped)	null	null
Coupled experiment	dm='nnn'	dm='nnn'
X decoupling during acquisition	dm='nny'	dm='nny'
Setup experiments	phase=1	phase=1
Hypercomplex experiment	phase=1,2	phase=1,2
TPPI	phase=3	phase=3
Minimum nt possible	2	2
Presaturation and/or multiple-bond correlation	see text	not applicable
Axis parameter for proper ppm on both axes	bq	dp

Table 32. Parameters for HMQC and HMQCR Pulse Sequences

5.8 Recabling Single-Broadband Systems

On systems with computer-controlled switching between normal and reverse modes (described above), the spectrometer switches automatically to reverse mode when you use either HMQCR or S2PULR sequences. If you have such a system, ignore any further recabling instructions in this chapter.

MERCURY-VX, *MERCURY*, and *GEMINI 2000* have automatic switching. If you have a *MERCURY-VX*, *MERCURY*, and *GEMINI 2000* broadband system, skip to "Filters for Indirect Detection," page 184.

To manually switch single-broadband systems to a reverse configuration, perform the following steps:

- 1. Toggle the **L.O. SELECT** switch that enables the local oscillator (L.O.) gate on the decoupler channel of the instrument. This switch is installed as follows:
 - As part of the Indirect Detection Modification Kit 00-990841-00 (see Magnetic Moments Vol. IV, No. 2, p. 16)
 - By service personnel using Varian Service Bulletin MP880002
 - Customer-installed using Magnetic Moments Vol. III, No. 4, p. 7
 - Factory-installed
- 2. Unplug the L.O. cable from the transmitter board (J2X5 on UNITY*INOVA* or UNITY*plus*, J3402 on other systems) and plug it instead into the L.O. connector of the decoupler board (J2X5 on UNITY*INOVA* or UNITY*plus*, J3302 otherwise).
- 3. Reconnecting the cables carrying the transmitter and decoupler signals depends on the system configuration. You may need to try more than one variation to find the correct method for your system; this can be done by either using the HMQCR pulse sequence itself or by using an oscilloscope:
 - Spectrometers shipped after October, 1990, with "N-type" connectors, should require no recabling.
 - The following method should work on earlier UNITY spectrometers: Disconnect the decoupler cable from the Bird wattmeter inside the magnet leg, and connect it directly to the X-decouple input to the probe (via the appropriate filter).
 - The following method has been found applicable on some systems: Reverse the observe and decoupler inputs to the probe; that is, take the cables (and associated filters) normally connected to the 1H/OBS/DEC BNC jack and the NORMAL jack on the inside of the magnet leg and reverse them.

5.9 Recabling Dual-Broadband Systems

Because dual-broadband systems are operated in the "normal" mode, there is little recabling necessary (with the exception of filter placement; see the next section). Both lock and ¹H observe signals follow their normal pathway.

The X channel, however, does have to be recabled. Under conventional (direct observation of the X-nucleus) circumstances, the X channel of the probe will be connected to the NORMAL jack on the inside of the magnet leg; for indirect detection, the X channel of the probe must be connected to the (now mislabeled) 1H/OBS/DEC jack. (The above information applies to 500-MHz and 600-MHz system only; dual-broadband systems at other frequencies require different cabling).

5.10 Filters for Indirect Detection

A number of filters are necessary for optimum performance of indirect detection experiments, regardless of whether a single- or dual-broadband system is used.

In the lock channel line, install a ²H band-pass filter. When the filter is added, expect the lock phase to change. This filter can be left in the system at all times; it will, however, cause a small (about 3 dB) loss in lock sensitivity.

Route the cable attached to the probe connector marked {x} decouple through an X-nucleus band-pass filter. Route the cable attached to the ¹H OBSERVE spigot through a ¹H high-pass filter; route it through a ²H band reject filter as well (not supplied) if 2 kHz lock birds are seen (5 kHz on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*).

Table 33 lists part numbers for the filters supplied with UNITY INOVA, MERCURY-VX, MERCURY, UNITY plus, GEMINI 2000, and UNITY indirect detection probes. If you are performing indirect detection experiments on other probes (Varian broadband or switchable probes, or probes purchased from other vendors), you may need to purchase these filters separately.

 Table 33.
 System Bandpass Filters for Indirect Detection Probes.

Filter	300-MHz System	400-MHz System	500-MHz System	600-MHz System	750-MHz System
15 _N	BE30.4-7.6-9BB	BE40-10-9BB	BE53-15-8BB	BE61-10-8BB	BE77-15-4BB
$2_{\rm H}$	BE46-4.5-6BB	BE61-10-8BB	BE77-3.8-8BB	BE92-9-6BB	BE115-11-6BB
13 _C	BE75-15-8BB	BE109-22-8BB	BE135-35-8BB	BE151-40-8BB	BE188-20-7BB
31 _P	BE135-35-8BB	BE151-40-8BB	BE175-60-8BB	BE240-100-8BB	BE301-46-8BB

5.11 Tuning the Probe in the Reverse Mode

The probe can easily be tuned once the recabling described above is performed. On *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*, use a separate experiment with the standard parameters to tune the probe (see the *NMR Probes Installation* manual for further information).

Tune the ¹H Channel

- 1. Make sure the spectrometer is transmitting the proton frequency by entering tn='H1' su.
- 2. Connect the cable that comes from the ¹H channel of the probe to the TUNE connector on the preamplifier or magnet leg.
- 3. Tune as described in the probe installation manual. After tuning, put the cables back the way they were.

Tune the X Channel

- 1. Enter the nuclei for the X channel (e.g., to use ${}^{13}C$, enter tn='13C' su).
- 2. Connect the cable that comes from the X channel of the probe to the TUNE connector on the preamplifier or magnet leg.

- 3. Reconnect the cables containing the transmitter and decoupler signals as described in step 3 on page 183.
- 4. Tune as described in probe installation manual. After tuning, put the cables back the way they were.

5.12 Controlling Transmitter Power in the Reverse Mode

On a Sun-based system, you will generally find that the parameter dpwr is limited to a value of 49 to prevent excessive decoupling power. When dpwr is used to control the ¹H transmitter power, as it is in the S2PULR and HMQCR sequences, this will be insufficient power. To remove this limitation, you can use the config program to reset the Upper Limit of the coarse attenuator for the decoupler.

5.13 Indirect Detection Calibration

On both single- and dual-broadband systems, calibrations are not the same after cabled for indirect detection as in "normal" mode, and so it is necessary to calibrate the system after recabling. The first system calibration should be done with the standard sample (Part No. 00-968120-96), which contains 1% ¹³CH₃I (as well as a number of other species). This sample also enables running a quick HMQC spectrum to verify overall operation of the system.

- *Note:* Throughout the following instructions, refer to Table 32 to see which parameters control the features in your configuration.
 - The first time you perform this experiment, cable the system for normal ¹³C observe, and obtain a normal ¹³C spectrum in the usual way. The ¹³C signal from ¹³CH₃I is extremely far upfield (22.3 ppm), so depending on your standard parameters you will probably need to either increase your spectral width or move your transmitter offset upfield (to small or more negative numbers) to prevent foldover of the signal. A typical result is illustrated in Figure 45 in which the peak at 22 ppm is the ¹³CH₃I signal. Note that the decoupler should be on. Increase scans for low field systems (for example, 16 for 300 MHz).

Once you have identified the signal 22 ppm (it should be easily visible with a single pulse—this is a 50% enriched sample!), position the cursor near the line and type nl movetof. The software will change the transmitter offset tof so that it is on-resonance for this ¹³C resonance. Make a record of this value of tof for future use. Now, if needed, re-cable the system to the indirect detection configuration.

2. Obtain a normal ¹H spectrum. On the UNITY INOVA, MERCURY-VX, MERCURY, GEMINI 2000, UNITYplus, and dual-broadband UNITY and VXR-S systems, use S2PUL with tpwr set to the desired value (it will also be active in HMQC), tn set to 'H1', and dn to the appropriate X nucleus. On single-broadband UNITY and VXR-S systems, use S2PULR to pulse the decoupler, positioning the spectrum with dn and dof the same as in the indirect detection experiment (and also set tn to the proper X-nucleus at this time). S2PULR pulses the decoupler at a power level controlled by the parameter dpwr.

If you observe a proton spectrum, you have confirmed that you have properly recabled the ¹H channel. Using the movesw command, narrow the spectral width to the region containing the full proton spectrum, as seen below (note that if you are operating in the reverse mode with S2PULR, the commands movesw and




Figure 45. Normal ¹³C Spectrum of ¹³CH₃I

movetof are not appropriate, since they change the value of tof rather than dof. Either recalculate the proper values "by hand," or use a series of commands like rl=tof tof=dof movesw dof=tof tof=rl).

The large three-line pattern (approximately 1:1:1) in the spectrum, illustrated in Figure 46, (centered at 2.2 ppm) represents unenriched CH₃I (the central line) and




Figure 46. Normal ¹H Spectrum of ¹³CH₃I

the enriched ¹³CH₃I (the outer two lines). The peaks centered at 1.5 ppm in the spectrum are a ¹³C enriched impurity peaks, and the peaks near 3.7 ppm are trimethylphosphite, which is also present in the sample (for ³¹P indirect detection experiments). Now calibrate the proton pulse width pw in the usual way.

3. Enter hmgc or hmgcr (as appropriate) to set up the HMQC experiment.

Before you proceed any further, make sure the spinner is off by setting **spin=0** or by using the acqi window. Listen for the spinner air to turn off and check in the acquisition status window that the spinner is off and actually shows 0 Hz.

4. Verify that the proton amplifier power (if appropriate) and pw are what you established in step 2. If you have a linear amplifier on your X channel, set its pulse power level with pwxlvl. Set the position of the ¹³C channel to the on-resonance position measured in step 1. Enter phase=1 ni=1 nt=1 dm='nnn' null=0 ss=0 j=151 spin=0 ai pwx=0, 15 wexp='wft dssh' au.

Two spectra are collected. In the first, no X-nucleus pulses are used (pwx=0), so a normal spectrum should appear (inverted); in the second, X-nucleus pulses are used, and if the value of pwx is correct, the ¹³C satellites of the ¹³CH₃I should appear right side up, as seen in Figure 47. If the two spectra are identical, you have either




Figure 47. HMQC Without and With X-Nucleus Pulses

entered a value for pwx that is exactly 360° (unlikely but you can easily find out), you have not cabled the X channel correctly, or dof might be set wrong.

Note two occurrences in this experiment:

- Because the pulse sequence involves spin echoes whose timing may not be perfectly adjusted, it may not be possible to perfectly phase the first spectrum; the spectrum on the left in Figure 47 is typical. If the spectrum phases poorly, recheck the proton 90° pulse width.
- Because the pulse sequence is very much optimized for a particular coupling constant, the ¹³C satellites of the upfield impurity, which have a coupling of 128 Hz, exhibit a strange phase in the spectrum obtained with a 90° *pwx* pulse on the right. This is a consequence of the sequence and you should not attempt to phase this spectrum "properly;" since you will be unable to do so.
- 5. Assuming you see a difference between the two spectra in step 4, perform an array of values for pwx and select the one for which the ¹³C satellite signals are maximized (make sure to set d1 sufficiently long for this experiment, say, 10 seconds). In the spectra shown in Figure 48, the series on the left represents a broad array of pwx values from 0 to 52 µs in steps of 4 µs. Such a series can be used to get the "big picture" when you have no idea of the proper value of pwx.




Figure 48. Calibration of pwx, Coarse and Fine

Note that unlike a standard pulse calibration, the "null points" in this calibration are *not* the points of interest. The signal starts out inverted at 0°, becomes maximally upright at 90°, inverts at 180°, etc. In the series on the right, *pwx* has been varied from 12 to 16 μ s in steps of 0.5 μ s to more properly select the 90° value. Note also that this calibration will be *severely* impacted if the X-nucleus transmitter is not placed on-resonance during the experiment, even though the nucleus you are observing is ¹H. Note that the view is set to ai (absolute intensity) so that peak heights can be compared.

You can also calibrate pwx with the **pwxcal** macro and sequence (not supplied with the *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*). Refer to section 21.5 "Calibrating 13C (or X) Decoupler Pulse Width with PWXCAL" in the manual *System* Administration for the procedure.

6. For systems without linear amplifiers, skip to step 7.

For systems with linear amplifiers, the following is optional: Increase the power of the amplifier (using **pwxlvl**) by 1 or 2 units and repeat step 5. pwx should decrease as you increase the X-nucleus pulse power. When pwx no longer decreases, the probe is arcing; decrease the power by 3 units to be safe and do not go above this value. Proceed to step 7.

7. Check cancellation. Set nt=4 dl=10 ss=1 and array pwx with the 0° and 90° pulse widths. The two spectra that result should show no signal in the first case and ¹³C satellites only in the second case. If the cancellation is worse when the X-nucleus pulses are present, as in Figure 49, this is an almost certain indication of a lock interference problem. The spectrum in Figure 49 was performed without a ²H band-pass filter on the lock channel. Compare it to the spectrum in Figure 50, a repeat of the same experiment with a ²H band-pass filter.

When you repeat this experiment on a "real" (i.e., natural abundance) sample, the difference between these two spectra will give you confidence that what you are observing are indeed ¹³C satellites and not residual uncanceled signals.

8. Choose either step a or b—on systems other than single-broadband UNITY and VXR-S, perform step a and skip step b below. On single-broadband UNITY and VXR-S systems, skip step a and perform step b.





Figure 49. HMQC 4-Pulse Cancellation without ²H Bandpass Filter





After this step, the carbon parameters need to be calibrated.

a. Systems other than single-broadband UNITY and VXR-S – Modify the parameters to perform two single-frequency, off-resonance decoupled experiments in order to calibrate γB_2 . Set **axis='h' dpwr=43 dm='nny' dmm='c' at=0.256 lb=2**. Enter an array of **dof** values that are ± 4000 Hz from the on-resonance position. h2cal returns γH_2 =xxxhz, pw90 at that power, and coelesc frequency. Set dmf=1/90 or $4*\gamma H_2$.

Obtain two spectra and measure the residual coupling between the two lines. To do this, enter ds(1), place two cursors on the two outer lines, and enter r1=delta. This will save the first splitting in the variable r1. Now enter ds(2), place two cursors on the two outer lines, and enter r2=delta. Now

enter h2cal(r1,r2,151) (151 is the full coupling constant), and the computer will display the value of γB_2 (to obtain a printed copy, enter printon h2cal(r1,r2,151) printoff).

Increase or decrease **dpwr** by 1 or 2 units and repeat step 2. The highest possible value of dpwr that you will be able to use for WALTZ or GARP modulated decoupling may be limited by the maximum decoupler modulation frequency, dmf. For ¹³C indirect detection, this will be a reasonable power level if you wish to decouple the full ¹³C spectrum.

If you wish to decouple only part of the ${}^{13}C$ spectrum, you should decrease γB_2 accordingly; the lower the amount of decoupler power used, the less sample heating and convection that will occur, and the better will be your cancellation in indirect detection experiments. The greater the decoupling power used, the more critical it is to establish a temperature steady-state in the experiment and to keep the duty cycle low (see discussion above in the section "Cancellation Efficiency").

b. Single-broadband UNITY and VXR-S systems – Now that you are convinced you are seeing ¹³C satellites, check your decoupling. Set pwx to the 90° value, pwxw to the same value and tpwr=pwxlvl dm='nnn', 'nny' at=0.256 lb=2 au. These two spectra should now both contain ¹³C satellite signals, in the first case coupled, and in the second case ¹³C-decoupled, as seen in Figure 51.




Figure 51. HMQC Without and With ¹³C Decoupling

Note that these two spectra may not have precisely the same phase, due to the nature of the pulse sequence; this is of no consequence. Further note that the peaks in the decoupled spectrum, even after optimizing the decoupling, will generally only be approximately the same size or slightly larger than the coupled spectrum, although you might have expected a peak twice as large. ¹³C decoupling requires the ability to spread the decoupling power over a very wide range, and even with modern decoupling techniques such as WALTZ, this is typically accompanied by a line broadening in the decoupled spectrum compared to the coupled spectrum.

Once you have verified decoupling is functional by comparing these two spectra, set **dm='nny'** and array **pwxw** to find the value that gives optimum decoupling. This value should be close to pwx but may be a few microseconds longer.

We started with the X-nucleus power the same for pulses (pwxlvl) and decoupler (tpwr). It may be desirable to use less power for decoupling (see "Cancellation Efficiency," page 179). If you wish to do so, you should decrease **tpwr**. For every 3 units (3 dB) tpwr is decreased, pwxw by should be increased by about a factor of $\sqrt{2}$ (1.414), or by a factor of 2 for a 6 dB decrease in tpwr. A constraint is placed on the length of pwxw by the fact that 6*pwxw must be less than 1/sw.

- 9. You have now determined proper values for a number of parameters. A good idea is to save the parameters in one of the following ways:
 - Enter the parameters into the appropriate probe calibration file.
 - Join another experiment and recall the standard parameters by entering rtp('/vnmr/parlib/hmqc13') or rtp('/vnmr/parlib/hmqcr'). Enter the correct values for the various parameters and save the parameter set by using svp to save the parameters in a file of the same name. Now these parameters will be reestablished when you use macros to set up future experiments.
- 10. Perform a simple 2D experiment on the enriched sample. Starting with the experiment in which you have been performing all the calibrations, set ni=128 nt=4 phase=1,2. Move the ¹³C transmitter position 17 ppm downfield (to high frequency) from the on-resonance position for the ¹³CH₃I by entering dof=dof+17d (for hmqc) or tof=tof+17p (for hmqcr). Set the ¹³C spectral width sw1 to 50 ppm by entering sw1=50d (for hmqc) or sw1=50p (for hmqcr). To acquire a coupled 2D spectrum, set dm='nnn'; to obtain a decoupled spectrum, set dm='nny'.

Now acquire the data by entering **go**. Following acquisition, set **fn=2*np fn1=512 gf=.5*at** and **gf1=.5*ni/sw1** and process with **wft2da**; the result should look something like Figure 52.




Figure 52. Coupled HMQC Spectrum of ¹³CH₃I

11. When processing the data, one of the trickiest aspects of this experiment is phasing in the f_1 direction (phasing in f_2 is accomplished simply by performing a 1D transform on the first increment with wft(1) and phasing the spectrum, paying attention only to the ¹³C satellite peaks). In F1, the combination of the usually large spectral width and the pulse in the center of the evolution time produces large negative values for lp1.

Reasonably good starting points for the F1 phase can be calculated according to the following formulas:

```
 lp1= -sw1 \times 360^{\circ} \times (2*rof1+2*pw+4*pwx/p)  (for hmqcr)
  lp1= -sw1 \times 360^{\circ} \times (4*pwx/p)  (for hmqc)
  rp1= -lp1/2
```

5.14 Typical Experimental Protocol for HMQC Experiments

To run HMQC experiments on a "normal" sample, you use a procedure similar to that described above for the standard (enriched) sample. The following instructions mostly emphasize the small differences in operation. Where the instructions are cryptic, refer back to the instructions in the previous section for detail.

A good "normal" sample to use for your first natural abundance sample is the Varian ¹H App Test sample (Part No. 00-968120-93), which is 1% 3-heptanone in CDCl₃.

Throughout the following instructions, refer to Table 32 to understand which parameters control the features in your configuration.

- Insert the sample and, after shimming, leave the spinner off. If you are going to run the experiment at a controlled temperature, enter temp=x su to start the temperature regulation process.
- 2. Set up to obtain a normal carbon spectrum and narrow the spectral width to the appropriate region. In some cases, the ¹³C spectrum will be too weak to observe in a reasonable amount of time. To set the parameters controlling the ¹³C frequency and spectral width if this is the case, you can take two approaches. First, if you have done similar experiments in the past on similar samples, just use the same parameters. Alternatively, if you use the standard command setup('Cl3','CDCl3') (or whatever solvent is appropriate), the spectrum obtained should be properly referenced. Now even if you can't see the peaks in the spectrum, you can apply the appropriate knowledge of the expected chemical shift range to place two cursors where you think the edges of that range will be, and use the movesw command to narrow the spectral width.
- 3. Obtain a proton spectrum using the same cabling you will use for indirect detection (and s2pulr if appropriate) and narrow the spectral width (remember, when using s2pulr you move dof and not tof). Check the calibration of the pulse width by entering pw=4*pw ga. Look only at the signals near the center of the spectrum and see if they produce a null signal. If they are negative, enter pw=pw+0.8 ga; if they are positive, enter pw=pw-0.8 ga; repeat until a good null is found, then enter pw=pw/4.
- 4. Switch to the HMQC experiment, and set the relevant parameters based on the results of steps 2 and 3.
- 5. Enter **phase=1 ni=1 dm='nnn' null=0 ai wexp='wft dssh'**. Set **j** to an appropriate value (normally 140 for C–H), and set **nt** to 4 or more transients, depending on the concentration of the sample (signal to noise needs to be sufficient

to allow you to see the ¹³C satellites). Now set **pwx** to an array of 0 and 90° and enter **au** to acquire two spectra. Proceed only if the two spectra are sufficiently different to give you confidence that the second spectrum is showing you satellite peaks only and not just residual uncanceled intensity of the protons attached to ¹²C.

If you are convinced that you are correctly connected but not happy with the quality of the spectra achieved at this step, skip ahead to step 8 and optimize the null parameter, then return here to check and optimize pwx. In either case, this is a good time to go over the checklist in the section "Cancellation Efficiency," page 179, making sure you have done everything possible to optimize cancellation.

The spectrum in Figure 53 shows the result of this experiment on a sample of 1% 3-heptanone at 300 MHz, using nt=64 null=2.0 and d1=2.




Figure 53. Verifying Cancellation with pwx=0,90

- 6. Now run an array of **pwx** around your expected 90° value, picking the one that gives you the largest satellite signals.
- If you will be decoupling during acquisition, set pwx to its 90° value and enter an array of dm='nnn', 'nny'. These two experiments should show coupled and decoupled spectra, respectively.
- 8. Now, if appropriate, optimize the parameter null. Set nt=1 ss=4 and enter an array of null values with at least one very short value (e.g., 0.001) and one very long value (e.g., 2.0). Because this experiment depends on the relaxation times of the spins involved, you'll also want to set at and dl to the same values you'll be using in the 2D experiment. Now run the array and select the value of null for which either most of the peaks, or the biggest peaks, or the peaks you are most interested in (the criterion is up to you), are approximately zero; remember, no one value of null will be correct for all peaks. Figure 54 shows this experiment run on a sample of 28 mg of gramicidin, with null arrayed over the range of values: 0.001, 0.05,0.1,0.2,0.3,0.4,0.5, and 2.0; examination of the spectra shows clearly how different values of null might be chosen.
- 9. If presaturation is desired (hmqc only, and not on *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*), you can set it up in the following manner. We will need to observe




Figure 54. Optimizing the BIRD Nulling Time

the ¹²C-bound protons, so set **nt=1 dm='nnn'**, and set **null=0** to omit the nulling period (for now at least).

- a. Enter **ga** and a proton spectrum will be observed.
- b. Move the FID to a different experiment, join the different experiment, and retransform the data (e.g., mf(1,2) jexp2 wft).
- c. Set the cursor on a peak that is to be removed by presaturation, and enter nl movetof. Note the value of tof selected and then copy this value back to your original experiment into the parameter satfrq (e.g., jexpl satfrq=x).
- d. Now set **satflg='yn'** and **satdly** equal to a time significant compared with T_1 of the peaks (e.g., satdly=1).
- e. Array satpwr to find the minimum value for which the peak will be removed (e.g., satpwr=10, 7, 4, 1 au). When this is determined, if you wish to use presaturation, set satflg='yy', reset null, and set satpwr to the value determined. If you do not wish to use presaturation, set satflg='nn'.
- 10. Set up the 2D experiment. Set **ni** between 128 and 256, **phase=1,2**, and **nt** to an appropriate number (comparable to what you were using in step 5).
- 11. Phasing in f_2 is accomplished by performing a 1D transform on the first increment with **wft(1)** and phasing the spectrum, paying attention only to the ¹³C satellite peaks. In f_1 , the combination of the usually large spectral width and the pulse in the center of the evolution time produces large negative values for 1p1.

Reasonably good starting points for the f_1 phase can be calculated according to the following formulas:

 $lpl= -swl \times 360^{\circ} \times ((2 \times rofl) + (2 \times pw) + (4 \times pwx/\pi) \text{ (for hmqcr)}$ $lpl= -swl \times 360^{\circ} \times ((4 \times pwx)/\pi) \text{ (for hmqc)}$ rp = -lpl/2

Table 34 lists typical parameters for a quick 2D experiment on the 1% 3-heptanone sample on UNITY and VXR-S systems. Expect to see artifacts in these spectra. The residual

uncanceled signals from protons attached to ${}^{12}C$ show up as stripes parallel to the f₁ axis at the frequency of each ${}^{1}H$ peak. This artifacts will be larger for peaks with long T_{1} , such as solvent peaks (e.g., residual protons on a deuterated solvent) or methyl groups. In Figure 55 they are seen at 2.4 ppm, 1.0 ppm, and 0.9 ppm.

Parameter	300 MHz	400 MHz
Proton spectral width	1000	1000
Carbon spectral width	4000	4000
Number of transients (nt)	16	16
tof	-5732.9	-7060.7
dof	-598.6	-1761.1
Number of increments (ni)	32	32
Mode	Hypercomplex (phase=1,2)	Hypercomplex (phase=1,2)
d1	2.0	2.0
null	2.0	2.0
at	0.256	0.256
gf	0.128	0.128
gfl	0.008	0.008
Solvent	CDCl ₃	CDCl ₃
j	130	130
lpl	-125.2	-125.2
rpl	62.6	62.6

Table 34. Parameter Values for HMQCR on Natural Abundance Sample




Figure 55. Coupled HMQC Spectrum of 3-Heptanone

Axial peak artifacts, which are common, will show up either at $f_1=0$ (the center of the spectrum in f_1) or, as in Figure 55, at the edges of the spectrum in f_1 (if FAD is used), again

Chapter 5. Indirect Detection Experiments

at f_2 frequencies corresponding to each ¹H peak and possibly through the entire spectrum. Another common artifact seen in Figure 55 is the "0,0" artifact in the exact center of the spectrum. Some peaks in the ¹H spectrum, of course, will not appear in the HMQC spectrum, because they represent protons that are not bound to ¹³C (e.g., protons from water or NH groups). This is not the case with 3-heptanone, however.

If you are unfamiliar with HMQC spectra, you may be surprised to see multiplet structures. You should realize that during the detection period we are detecting a normal (i.e., with ¹H-¹H couplings) proton spectrum, albeit a spectrum of only those protons attached to ¹³C. Thus, in Figure 56, we see that the proton attached to the carbon at 37.2 ppm is a quartet (it's adjacent to a CH₃ group), while the proton attached to the carbon at 43.4 ppm is a triplet (it's adjacent to a CH₂ group). In the ¹H spectrum itself, these two groups of protons are heavily overlapped (see Figure 57).




Figure 56. Expansion of Coupled 3-Heptanone HMQC Showing Multiplets




Figure 57. Decoupled HMQC Spectrum of 3-Heptanone

5.15 Differences for ¹⁵N Indirect Detection

Calibrations and operations for ¹⁵N proceed largely along the lines outlined above for ¹³C. In the standard sample, 2% ¹⁵N-benzamide (Part No. 00-968120-97), the ¹⁵N satellite lines are partially obscured by other resonances in the conventional 1D spectrum, and so the ¹⁵N pulse width calibration must be done using multi-transient HMQC experiments as described in "Typical Experimental Protocol for HMQC Experiments," page 192. Be sure to use a J appropriate for NH couplings (90 Hz) in this case.

A step that can often be done in ¹⁵N work of peptides is to make sure that the γB_2 is sufficient to decouple the relatively narrow range of ¹⁵N chemical shifts expected in such samples but no more. This minimizes heating effects and improves cancellation. Typical acquisition times (at) are 0.075 to 0.1 seconds.

5.16 HSQC Experiment

The hsqc macro sets up parameters for the HSQC pulse sequence, a heteronuclear Overbodenhousen experiment using REVINEPT.

Applicability

Not supplied with MERCURY-VX, MERCURY, and GEMINI 2000 systems.

Parameters

sspul='y' selects for *trim(x)-trim(y*) sequence at the start of the pulse sequence; sspul='n' selects a normal experiment.

satmode='yn' gives presaturation during relaxation period (satdly) with the transmitter; satmode='nn' gives no presaturation during relaxation period (satdly); satmode='ny' gives presaturation during only the null period.

satfrq sets the presaturation frequency.

satdly sets the saturation time during the relaxation period.

satpwr sets the saturation power for all periods of presaturation with xmtr.

hs='yn' sets a homospoil pulse (hst) during the d1 relaxation delay.

null is the delay associated with the BIRD nulling.

tpwr is the power level for ¹H transmitter pulses.

pw is a 90° transmitter pulse length for protons (the observed nucleus).

pwxlvl is the power level for X decoupler pulses.

pwx is a 90° decoupler pulse length for X.

jxh is a one-bond heteronuclear coupling constant to X (in Hz).

phase=1, 2 for hypercomplex experiment with F1 quadrature (complex F1-FT).

Chapter 5. Indirect Detection Experiments

Chapter 6. Data Analysis

Sections in this chapter:

- 6.1 "Spin Simulation," this page
- 6.2 "Deconvolution," page 205
- 6.3 "Reference Deconvolution," page 208
- 6.4 "Addition and Subtraction of Data," page 211
- 6.5 "Regression Analysis," page 216
- 6.6 "Chemical Shift Analysis," page 221

6.1 Spin Simulation

VNMR software includes an iterative spin simulation program based on the FORTRAN program LAME, also known as LAOCOON with magnetic equivalence added. LAME calculates the theoretical spectrum for spin-1/2 nuclei, given the chemical shifts and the coupling constants.

Up to eight closely coupled, non-equivalent spins (ABCDEFGH) can be handled. Equivalent spins can be treated by magnetic equivalence factoring to extend the simulation to systems such as A3B2CD3. The X-approximation can be used to handle different types of nuclei. Nuclei are treated as different types if there is at least one spare letter in the alphabet between their groups (e.g., ABD and ABX are both systems using the Xapproximation.) Frequencies, intensities, energy levels and transitions can be listed, and simulated spectra can be displayed and plotted.

Parameters can be adjusted by iteration to approach a given experimental spectrum. For iterative runs, one or several parameters can be kept constant. In addition, one or several parameters can be set equal to each other and held equal during the course of the iteration.

To acquaint yourself with the spin simulation software, we suggest you work through the next section, "Spin Simulation Step-by-Step." Note how the regression menus simplify the procedure. A number of specialized commands and parameters are also available. Table 35 lists these commands and parameters.

For more information about the spin simulation algorithms, refer to the following:

- Bothner-by, A.A. and Castellano, S., J. Chem. Phys., 41, 3863 (1964).
- Emsley, Feeney, and Sutcliffe, eds. 1966. *Progress in Nuclear Magnetic Resonance Spectroscopy*, Vol.1, Chap. 3. Oxford: Pergamon Press.
- Stanley, R.M.; Marquardt, D.W.; and Ferguson, R.C., J. Chem. Phys., 41, 2087 (1964).

Commands	
assign*	Assign transitions to experimental lines
cla	Clear all line assignments
dga	Display group of spin simulation parameters
dla<('long')>	Display spin simulation parameters arrays
dlalong	Long display of spin simulation parameter arrays
dll*	Display listed line frequencies and intensities
dsp<(file)>	Display calculated spectrum
initialize_iterate	Set iterate string to contain relevant parameters
<pre>spinll<('mark')></pre>	Set up a slfreq array
spins<(options)>	Perform spin simulation calculation
<pre>spsm(spin_system)</pre>	Enter spin system
undospins	Restore spin system as before last iterative run
<pre>* dll<('pos'<,noise_mult>)></pre>	<:lines>
assign<('mark')>, assign(transistion_number,line_number)
Parameters	
cla {array of real values}	Calculated transition number
clamp {array of real values}	Calculated transition amplitude
clfreq {real values}	Calculated transition frequency
clindex {array of real values}	Index of experimental frequency of a transition
<pre>iterate {string of parameters}</pre>	Parameters to be iterated
niter {1 to 9999}	Number of iterations
<pre>slfreq {real values}</pre>	Measured line frequencies
slw {0.01 to 1e6}	Spin simulation linewidth
smaxf {-1e10 to 1e10}	Maximum frequency of any transition
sminf {-1e10 to 1e10}	Minimum frequency of any transition
sth {0 to 1.00}	Minimum intensity threshold
svs {0 to 1e10}	Spin simulation vertical scale

Table 35.	Spin Simulation	Commands and Parameters
-----------	-----------------	-------------------------

Spin Simulation Step-by-Step

The simplest way to get acquainted with the spin simulation software is to work through a step-by-step example. The following example is complete with comments to help you understand what you are doing at each step.

1. Click on Main Menu > File > Set Directory > Parent.

The text window displays a list of directories (entries with a backslash as the last character in the name) and files (if any). The status window (near the top of the screen) displays the pathname of the current directory.

- 2. Click on **Parent** as many times as necessary until the status window displays the message: Directory now "/".
- 3. Click the mouse on the directory entry vnmr / in the text window until it turns to inverse video. Then click on the **Change** button.

The text window displays the list of subdirectories and files (if any) in vnmr.

4. Click the mouse on the directory entry fidlib/ in the text window until it turns to inverse video. Then click on the **Change** button.

The text window displays the list of subdirectories and files (if any) in fidlib.

5. Click the mouse on the directory fidld.fid/ until it turns to inverse video.

6. Click on Return > Load > Process > Transform.

The graphics windows displays a spectrum and the opening menu from the interactive spectrum display program (ds) appears.

7. Click on **Next** > **Dscale**.

A scale is displayed under the spectrum.

8. Click the left mouse button near 3.5 ppm and the right mouse button near 3.2. Then click on **Expand**.

The six-line pattern shown will be simulated as an AX2Y system.

9. Click on Th.

A horizontal line for the threshold appears.

- 10. Use the left mouse button to move the threshold line below the tops of the peaks.
- 11. Enter **dll**.

The text window displays a line listing, which will be used later.

With a spectrum now set up, you are ready to use the spin simulation menus to display a simulated spectrum.

1. Click on Main Menu > Analyze > Simulation.

You are now in the Spin Simulation Main Menu.

2. Click on **Spin System** > **other** > **other** > **AX2Y**.

This picks the spin system and initializes its parameters.

3. Click on Set Params.

The spectrum reappears.

4. Click the left button in the center of the six-line pattern and enter **A=cr**.

This sets the chemical shift of spin A to the position of the cursor.

5. Click the left button on the center of the left-most line, the right button in the center of the second left-most line, and enter JAY=delta.

This sets the JAY coupling constant to match the difference frequency.

- 6. Click the right button on the center of the third line, and enter **JAX=delta**. Then click on **Return**.
- 7. Click on Show Params.

This confirms your entry of the spin system parameters.

8. Click on Simulate.

After a brief moment, the simulated spectrum will appear.

If you want to continue with iterative spin-simulation, take the following steps:

1. Enter iterate?.

The status window displays iterate='A, JAX, JAY', which confirms that the iterate parameter was automatically set.

2. Click on assign > auto assign.

The assign macro is executed, which assigns the lines from the dll listing to the lines from the previous simulation.

3. Click on **iterate**.

This performs an iterative optimization and displays the resulting spectrum.

4. Click on list.

The listing contains the values of the A, JAX, and JAY parameters that give the best iterated fit to the experimental spectrum.

Spin Simulation Menus

Almost the entire spin simulation analysis is available through six menus:

Spin Simulation Main Menu

Spin System Show Params Set Params Simulate Original Next

• Spin Simulation First Menu

AB ABC A2B ABCD A2BC A3B other return

• Spin Simulation Second Definition Menu

ABCDE A2BCD A3BC A2B2C A3B2 A3B2C other return

• Spin Simulation Third Definition Menu

AX AXY AX2 AXYZ AX2Y AX3 other return

• Spin Simulation Secondary Menu

list params assign display iterate observe return main

• Spin Simulation Line Assignment Menu



The Spin Simulation Main menu is the opening menu. Each of the other menus is entered through buttons on this menu or submenus of this menu. To open the Spin Simulation Main menu, click on Main Menu > Analyze > Simulation.

Entering a Spin System

The command spsm(spin_system) enables creating a spin system as an alphanumeric string of upper-case letters, and creates and initializes appropriate parameters to describe the various chemical shifts and coupling constants. Chemical shifts are stored in parameters A through Z, and the coupling constants are stored in the parameters starting with JAB and ending with JYZ. Different nucleus types are handled by using letters starting with A for the first type, X for the second, and M for the third.

Spin Simulation Parameters

Spin simulation commonly uses the following global parameters:

• cla and clfreq together make up a table consisting of line numbers assigned by the spin simulation program and the corresponding frequency of a measured line when the intensity of the line is above a threshold value set by the parameter th.

- clamp stores the transition amplitude of calculated transitions when they are above a threshold set by parameter sth.
- clindex is an index of experimental frequency of a transition.
- slfreq is a list of measured line frequencies.
- sminf and smaxf are the minimum and maximum frequency limits for calculation of the final simulated spectrum. These should be set before the calculation is performed. If the Set Params button is used, sminf is initialized to sp, and smaxf is initialized to sp+wp.
- sth is the minimum intensity threshold above which transitions are listed and included in the simulated spectrum. A typical value is 0.05.
- svs is the maximum intensity of calculated transitions.

The command dga displays the file of simulation parameters.

Performing a Spin Simulation

The command spins performs a spin simulation, using the current spin system parameters. The following variations are available:

- spins('calculate', 'energy') puts an energy level table in the output file.
- spins('calculate','transitions') puts a second table of transitions ordered by transition number in the output file.
- spins('iterate') runs in an iterative mode to match experimental and calculated lines.
- spins('iterate','iteration') lists parameters after each iteration in the output file.

The output file is spins.list in the current experiment. This file always includes the calculated transitions ordered by frequency.

Using the current table of transitions and intensities, dsp<(file)> calculates the simulated spectrum (using the current value of parameter slw for the linewidth) and displays the spectrum. dsp can only be used after the spins program has been run. When dsp has a filename as an argument, the spectral information is taken from that file. After the display appears, it may be modified and plotted like any other 1D spectrum.

Iterative Mode

The following commands are used to set up files for the "iterative" mode of spin simulation in which the calculated spectrum approximates an experimental spectrum.

- spins('iterate') performs the simulation in the iterative mode.
- initialize_iterate selects a default value for the parameter iterate that will cause all parameters to be iterated. iterate is a string parameter that contains a list of parameters (separated by commas) to be iterated during iterative spin simulations. Typical value is 'A,B,JAB'. If the Set Params button is used, iterate is initialized to a string containing parameters appropriate to the current spin system. The parameter niter is the maximum number of iterations that will be used by an iterative simulation. If the Set Params button is used, niter is initialized to 20.
- cla clears the file of line assignments used for iteration (matching simulated spectra to actual data). dla displays this file.

- dlalong stores the line assignments in the file spini.la of the current experiment. This command is most useful in more complex problems where the text window is too small for the dla display. dlalong displays the file in the text window.
- spinll copies the list of frequencies from the last line listing by nll or dll (contained in the parameter llfrq) into the simulation line frequency parameter slfreq. spinll also clears the previous line assignments and runs dla.
- spinll('mark') places the line positions in the file markld.out into the parameter slfreq. This is useful if you want to manually assign the lines.
- assign without an argument assigns the nearest calculated transitions to the lines from a dll or nll listing after spinll has placed them in slfreq. For positive lines only, use dll('pos'). All lines may not be assigned and transitions must be greater than sth, which should be 0.05 or greater to prevent assignment of extremely small lines. assign can also be run by clicking on the auto assign button.
- assign('mark') is the same as assign except that the file markld.out is used instead of the dll listing. Use the cursor and the mark button to place the lines to be assigned in the markld.out file. This file is cleared by mark('reset'). Use nl to move the cursor to the center of a selected line.
- assign(t#,l#) assigns a single calculated transition number (t#) to a line from a dll listing (the index is l#). assign(t#,0) removes the calculated transition assignment.
- undospins restores a spin system as it was before the last iterative run. Chemical shifts, coupling constants, and transition assignments are returned to those existing immediately before an iterative spin simulation.

The experimental line listing for assign should be made from the part of the spectrum whose limits are sminf and smaxf. If a frequency in slfreq has already been assigned a line number, a new entry will be made with the same frequency. If the line assignments produced by assign needs to be modified, refer to the description of the spins command in the *VNMR Command and Parameter Reference* for details.

The Spin Simulation Line Assignment menu is especially useful for iterative spin simulation. However, individual assignments are made using the assign command as previously described.

Spin Simulation Files

The spins.list file is an output table made by the spin simulation program. This file can be displayed by clicking on the list button. The following files can exist in the current experiment but, except for the file spini.la, are not normally of interest when spin simulation is run from VNMR:

- spini.la is the current transition assignments for an iterative spin simulation (produced by the dlalong command).
- spini.savela is the transition assignments for iterative spin simulation in a format readable by the macro undospins.
- spini.outpar are the values of the chemical shifts and coupling constants after an iterative spin simulation.
- spini.indata is the line assignment input for the UNIX-level program spins (deleted by spins after iterations are completed).
- spini.inpar is a list of parameters whose values are to be determined by spins('iterate').

- spins.inpar is a list of initial settings of a number of spin simulation parameters.
- spins.outdata is a file of frequencies, amplitudes, and transition numbers from a spin simulation. It is used in calculating the displayed spectrum.
- spins.stat contains constants related to iteration (deleted by the spins program).

6.2 Deconvolution

The VNMR software allows the deconvolution of observed spectra into individual Lorentzian and/or Gaussian lines. Up to 2048 data points from an expansion of an experimental spectrum can be deconvoluted at one time, and up to 25 lines can be fit to this section of the observed spectrum. For each line, the line shape can be defined to be Lorentzian, Gaussian, or a combination of both.

The following parameters are available for each line:

- Frequency (in Hz) of line
- Intensity of line
- Linewidth (in Hz) at half-height of line
- Gaussian fraction of line: 0.0 (completely Lorentzian) to 1.0 (completely Gaussian)

All parameters can be fit at the same time, or selected parameters can be removed from the fit. In addition, a linear baseline correction is always added to the fit to avoid large errors produced by base line offsets.

Deconvolution accomplishes much of its work by means of text files, which are written into the user's current experiment directory:

- fitspec.inpar contains the starting parameters (frequency, intensity, linewidth, and Gaussian fraction) for a subsequent fitting operation.
- fitspec.indata contains the point-by-point intensity of the spectrum in the region of interest (the region that is displayed when the fitting is begun).
- fitspec.outpar contains the final parameters (frequency, intensity, linewidth, and Gaussian fraction) after a fit has been done.
- markld.out contains the result of a mark operation during a spectral display. By using the Use Mark button, this file may be used as an alternative to the last line list in setting up initial guesses for a fitting operation.

For best results, use the following settings:

- Use fn at least 2*np if not larger for adequate digitalization of the line shapes.
- For complex problems, use the macro usemark (see the description below) to set the best possible initial guesses.

Deconvolution Step-by-Step

- 1. Retrieve a sample 1D FID by entering rt('/vnmr/fidlib/fidld').
- 2. Enter fn=65000.

For proper digitalization of the line shape, you should almost always use a larger Fourier transform size than "normal" when the spectrum is to be deconvoluted.

3. Transform the FID by entering wft.

- 4. Using standard spectral manipulation operations, expand the two-line pattern near 8.0 ppm until it fills the center third of the display, with baseline on both sides. Enter ai (or click on buttons to select the ai mode) to select the absolute intensity mode, which is always required for deconvolutions, simulations, etc. Set a threshold that lists exactly two lines, then enter dll.
- 5. Click on Main Menu > Analyze > Deconvolution.

The Deconvolution Menu is displayed.

6. Click on Use Line List.

This produces a line list and a file containing the starting point for the deconvolution. This button automatically measures the linewidth of the tallest line on the screen and uses that as the starting linewidth for the calculation.

7. Click on Fit.

The analysis is performed. This particular example is a 6 parameter fit (2 frequencies, 2 intensities, and 2 linewidths). When the analysis is done, the calculated spectrum is displayed in the graphics window, and the numerical results appear in the text window (click on Flip as necessary to alternate between the two). The numerical output should be similar to this:

```
Number of data points:
                            404
Final chi square:
                            318046.844
Total number of iterations: 19
Successful iterations:
                          19
Digital resolution
                            0.365 Hz/point
ITERATION HAS CONVERGED
Parameters:
line frequency intensity integral linewidth gaussian fraction
     3176.780 125.293
                          423.948
                                    2.154
                                                0.000*
1
     3168.700 135.451
2
                          421.147
                                    1.979
                                                0.000*
```

8. Click on Plot.

The original spectrum, the calculated spectrum, and each of the component lines is plotted automatically, along with the numerical results of the calculation. At the end of this operation, the original spectrum replaces the calculated one.

9. Click on Show Fit to return to viewing the calculated spectrum. Click on Add/Sub to view the original spectrum simultaneously with the calculated one. Select sub from the menu to view the difference between the two.

Performing Deconvolution

Table 36 lists the commands and parameters associated with deconvolution.

The parameter slw is the starting default linewidth for deconvolution calculations. This linewidth is set automatically when deconvolution is operated using the menu mode and is bypassed if the usemark macro has been used in conjunction with two cursor input. Typical value for slw is 1.

The command fitspec performs spectrum deconvolution by fitting experimental data to Lorentzian and/or Gaussian line shapes. fitspec uses as a starting point data in a file fitspec.inpar, which must be prepared prior to performing the calculation. This file contains the frequency, intensity, linewidth, and (optionally) the Gaussian fraction of the line shape. Any number followed by an asterisk (*) is held fixed during the calculation; all other parameters are varied to obtain the best fit.

fitspec creates a file fitspec.indata, which is a text representation of the spectral data (that part of the spectrum between sp and sp+wp). After the calculation is finished,

Commands		
dsp<(file)>	Display calculated spectrum	
fitspec<(<'usell'><,><`setfreq'>)>	Perform spectrum deconvolution (VNMR)	
fitspec	Perform spectrum deconvolution (UNIX)	
mark*	Determine intensity of spectrum at a point	
plfit	Plot deconvolution analysis	
<pre>setgauss(fraction), setgauss(fraction*)</pre>	Set a Gaussian fraction for line shape	
showfit	Display numerical results of deconvolution	
usemark	Use mark as deconvolution starting point	
<pre>* mark<(fl_position)><:intensity>, mark<(left_edge,region_width)><:intensity,integral> mark<(fl_position,f2_position)><:intensity> mark<(fl_start,f1_end,f2_start,f2_end)><:intensity,integral,c1,c2> mark<('trace',<options>)><:intensity,integral,c1,c2>, mark('reset')</options></pre>		
Parameter slw {0.01 to 1e6}	Spin simulation linewidth	

Table 36. Deconvolution	Commands and Parameters
---------------------------------	-------------------------

the results of the fit are contained in a file fitspec.outpar, with a format identical to fitspec.inpar. All lines are set to have a linewidth of slw, and a fixed Gaussian fraction of 0. (Refer to the *VNMR Command and Parameter Reference* for information about keyword options available with fitspec.)

The setgauss macro modifies the output of the last deconvolution (fitspec.outpar) and makes it the input for a subsequent analysis (fitspec.inpar), after first modifying the Gaussian fraction:

- To allow this fraction to vary, use the format setgauss(fraction) where fraction is the Gaussian fraction of the line shape, a number naturally limited from 0 to 1, for example, setgauss(0.4).
- To fix the fraction, use the format setgauss (fraction*), suffix the fraction value (defined the same as above) with an asterisk and enclose the value in single quotes, for example, setgauss('1.0*').

In some cases it will not be possible to produce a line list that is a suitable starting point for a deconvolution. In this case, or in any case, the results of using the Mark button during a previous spectral display (ds program) may be used to provide a starting point. If the mark has been made with a single cursor, the information in the file markld.out contains only a frequency and intensity, and the starting linewidth is taken from the parameter slw. If the mark is made with two cursors, placed symmetrically about the center of each line at the half-height point, markld.out will contain two frequencies, an intensity and an integral. In this case, the starting frequency is taken as the average of the two cursor positions; the starting linewidth is taken as their difference.

Display and Plotting

After a deconvolution, the results are written into a file fitspec.outpar in an abbreviated format. The macro showfit converts these data to an output format more suitable for examination and printing.

The command dsp('fitspec.outpar') displays the theoretical spectrum described by the parameters produced by a deconvolution calculation. The macro plfit produces a

complete output plot of a deconvolution analysis, plotting the observed spectrum, the full calculated spectrum, each individual component, as well as the numerical results of the analysis. It can be invoked with the Plot button in the Deconvolution menu.

Deconvolution Menu

All actions necessary for deconvolution are accessible from the Deconvolution Menu.

This menu is opened by clicking on Main Menu > Analyze > Deconvolution. The Deconvolution Menu can also be opened by the command menu('fitspec').

6.3 Reference Deconvolution

The fiddle program performs reference deconvolution, using a reference signal with known characteristics to correct instrumental errors in experimental 1D or 2D spectra. The main command to start the program can take multiple string and numeric arguments: fiddle(option<,file><,option<,file><,start><,finish><,increment>)

option can be the following keywords:

'alternate'	Alternate reference phase $+ / -$ (for phase sensitive gradient 2D data).
'autophase'	Automatically adjust phase.
'displaycf'	Stop at the display of the correction function.
'fittedbaseline'	Use cubic spline baseline correction defined by the choice of integral regions.
'invert'	Invert the corrected difference spectrum/spectra.
'noaph'	Do not automatically adjust zero order phase of the reference region.
'nodc'	Do not use dc correction of the reference region.
'nohilbert'	Do not use Hilbert transform algorithm; use the extrapolated dispersion mode reference signal unless 'noextrap' is also used as an option.
'normalise'	Keep the corrected spectrum integrals equal to that of the first spectrum.
'readcf'	Read the correction function from file; the argument file must immediately follow 'readcf'.
'satellites'	Use the satellites defined in file in the ideal reference region; file should be in /vnmr/satellites.
'stop1'	Stop at the display of the experimental reference FID.
_	stop at the display of the experimental reference rise.
'stop2'	Stop at the display of the correction function.
'stop2' 'stop3'	
-	Stop at the display of the correction function.
'stop3'	Stop at the display of the corrected FID.
'stop3' 'stop4'	Stop at the display of the correction function. Stop at the display of the corrected FID. Stop at the display of the first corrected FID.

Table 37 list other commands for 1D and 2D variations of the fiddle program.

Commands		
fiddle*	Perform reference deconvolution	
fiddled*	Perform reference deconvolution subtracting alternate FIDs	
fiddleu*	Perform reference deconvolution subtracting successive FIDs from first	
fiddle2d*	Perform 2D reference deconvolution	
fiddle2D*	Perform 2D reference deconvolution	
fiddle2dd*	Perform 2D reference deconvolution subtracting alternate FIDs	
fiddle2Dd*	Perform 2D reference deconvolution subtracting alternate FIDs	
<pre>* (option<,file><,option<,file>><,start><,finish><,increment>)</pre>		

 Table 37. Reference Deconvolution Commands

Reference Deconvolution of 1D Spectra

Only spectra that contain a well-resolved reference signal dominated by a single component (i.e. not a simple multiplet) are suitable for reference deconvolution.

- Fourier transform the raw FID with ft, preferably having zero filled (i.e. set fn >= 2*np). (If there are sinc wiggles, use wft with gf = at*0.6.)
- 2. Set the reference line to the chosen signal using the **rl** command, and then use two cursors on either side of the line to define a region of spectrum that includes all of the reference signal plus a little clear baseline but no other signals. This reference region will be used to define the instrumental line shape.
- 3. Decide what line shape you would like to convert the instrumental line shape to, and set the weighting parameters accordingly. Thus, if you want a 1-Hz wide Lorentzian, set **1b** to 1 and all other weighting parameters to 'n'.

Bear in mind the signal-to-noise ratio penalty for resolution enhancement: if the experimental line is 2 Hz wide and you set lb=0, you get an infinitely sharp line with infinitely poor signal-to-noise. For most purposes, a sensible strategy is to set lb to *minus* the expected *natural* linewidth, and choose gf to give reasonable S/N; this strategy should convert the instrumental line shape to Gaussian. Where the signals of interest are broader than those of the reference, resolution enhancement can easily be obtained by making lb more negative.

- 4. Enter the **fiddle** command to carry out the reference deconvolution and display the corrected spectrum. The integral should remain unchanged, so any resolution enhancement will result in an increase in the amplitude of both signal and noise.
- 5. To save the corrected data, use the option 'writefid' when doing the reference deconvolution. For example, to store the file correctedfid.fid in the current working directory, enter fiddle('writefid', 'correctedfid').

The options 'writecf'<, file> and 'readcf'<, file> respectively write and read the correction function. Therefore, when you perform reference deconvolution on one FID using fiddle with the 'writecf' option and then use fiddle with 'readcf' to process another FID, the first correction function corrects the second FID. Reference deconvolution can be useful for heteronuclear lineshape correction (provided that the spectral widths for the two nuclei are in the ratio of the respective magnetogyric ratios) or for correcting spectra in which a reference signal has been suppressed (e.g., you could correct an INADEQUATE spectrum for lineshape errors by using a correction function derived from the normal carbon spectrum). To correct a series of spectra in an arrayed or 2D experiment, use numeric arguments, as with ft: fiddle(1) corrects spectrum 1, fiddle(2,3) spectra 2 and 3, etc.

Many reference signals have satellites. Like the familiar one-bond ¹³C satellites, for example, TMS has singlet satellite signals from coupling to ²⁹Si and quartet satellites (normally unresolved) from three-bond coupling to ¹³C. For most purposes, ¹³C satellites are small enough to be ignored, but where high accuracy is required or there are stronger (e.g. ²⁹Si) satellites, satellite signals can be included in the specified form of the ideal reference signal by invoking the 'satellites' option.

The /vnmr/satellites directory contains the file TMS with details of the TMS satellite signals. The command fiddle('satellites', 'TMS') allows for the satellite signals when deconvoluting using TMS as a reference.

The format for satellite files is that each line in the file consists of three real numbers in the following order:

- Separation of the satellite line from the parent signal, in Hz (0.5 JXH in the absence of homonuclear coupling).
- Intensity relative to the parent signal (natural abundance divided by the number of satellite lines [usually 2]).
- Isotope shift to high field, in ppm.

For example, the line

3.3 0.023 0.001

would correspond to a pair of satellite signals from a spin-1/2 isotope with an abundance of 4.6%, a coupling to the observed nucleus of 6.6 Hz, and an isotope shift to high field of 0.001 ppm.

Multiple lines in the satellite file account for multiple satellite signals on the parent peak.

To perform corrected-difference spectroscopy, use the command fiddled to produce the corrected difference between successive spectra, which divides arraydim in half. The difference spectrum is written into the second element of the pair. Because the main aim of reference deconvolution here is to optimize the purity of the difference spectrum, the target line shape would normally be chosen to give the best possible S/N; this method corresponds to choosing a target line shape approximately twice the width of the raw experimental signals of interest. The command fiddleu produces corrected differences between successive FIDs and the first FID.

Reference Deconvolution of 2D Spectra

The commands fiddle2d, fiddle2D, fiddle2dd, and fiddle2Dd function in just the same way as the parent fiddle program. Because the principal objective in 2D reference deconvolution is usually the reduction of t1-noise, ideal line shape parameters are normally chosen for optimum S/N ratio rather than resolution enhancement.

To perform 2D reference deconvolution:

- 1. Choose fn (preferably with $fn \ge 2*np$) and fn1.
- 2. Enter **ft** to transform the raw data (as mentioned earlier, if there is significant signal left at the end of at, it might be necessary to use wft with gf set).
- 3. Display the first increment with **ds(1)**, adjust the phase of the reference signal, and use **r1** to select the reference signal.

In earlier versions of fiddle, it was necessary to create a parameter, phinc, to anticipate the changes in the reference signal phase with increasing evolution time.

The current algorithm automatically adjusts the phase (unless you select the 'noaph' option). Deconvolution will set the reference signal phase as a function of t1 to place the reference signal at frequency rfp1 in f1. Therefore, remember to set rfl1 and rfp1 before using fiddle2D or the f1 frequencies might unexpectedly change.

4. Define the reference region with the two cursors, and then enter the command fiddle2D('writefid', <file>) (or fiddle2Dd if a 2D difference spectrum is required, as with corrected HMBC). The 'writefid' option is essential, because fiddle2D alone does not store the corrected time-domain data. If phase-sensitive gradient-enhanced 2D data is to be processed, alternate FIDs will have opposite phase modulations (i.e., the experimental array will alternate N-type and P-type pathways); in such a case, use the 'alternate' option.

After deconvolution is complete, the corrected 2D FID data can be read into an experiment and processed as normal (although if fiddle2Dd has been used, arraydim no longer matches the arrays set and it might be necessary to set the arguments to wft2d explicitly rather than using wft2da).

References

Further information on reference deconvolution can be found in the following literature:

Taquin, J. Rev. Physique App. 1979, 14, 669.

Morris, G. A. J. Magn. Reson. 1988, 80, 547.

Morris, G. A.; Cowburn, D. MRC 1989, 27, 1085.

Morris, G.A., Barjat, H., Horne T.J., Prog. NMR Spectrosc. 1997, 31, 197.

Gibbs, A; Morris, G. A. J. Magn. Reson. 1991, 91, 77.

Gibbs, A.; Morris, G. A.; Swanson, A.; Cowburn, D. J. Magn. Reson. 1983, 101, 351-356.

Rutledge, D. N. Ed. *Signal Treatment and Signal Analysis in NMR*, Chapter 16. Elsevier Science, 1996.

6.4 Addition and Subtraction of Data

The process of addition and subtraction of data is implemented using experiment 5 (exp5) as the add/subtract buffer file (overwriting any data that may have existed in exp5). Tools available for working with the add/subtract experiment include:

- Add/Subtract menu.
- Noninteractive commands.
- Interactive add/subtract program.

Table 38 lists the commands and parameters used with the add/subtract experiment.

Add/Subtract Menu

The Add/Subtract menu in the main menu system provides choices that activate command for clearing the add/subtract experiment (clradd command), activating an interactive

Commands		
add*	Add current FID to add/subtract experiment	
addi	Start interactive add/subtract mode	
clradd	Clear add/subtract experiment	
jexp1,, jexp9999	Join existing experiment	
select*	Select a spectrum or 2D plane without displaying it	
setvalue*	Set value of a parameter in a tree	
spadd*	Add current spectrum to add/subtract experiment	
spmin	Take minimum of two spectra in add/subtract experiment	
spsub*	Subtract current spectrum from add/subtract experiment	
sub*	Subtract current FID from add/subtract experiment	
<pre>* add<(multiplier<,'new'>)>, add('new'), add('trace',index)</pre>		
<pre>select<('next' 'prev' selection)><:index>,</pre>		
select<(<'f1f3' 'f2f3' 'f1f2'><,'proj'><,'next' 'prev' plane>)><:i>		
<pre>setvalue(parameter,value<,index><,tree>)</pre>		
<pre>spadd<(multiplier<,shift>)>, spadd('new'), spadd('trace',index)</pre>		
<pre>spsub<(multiplier<,shift>)>, spsub('new'), spsub('trace',index)</pre>		
<pre>sub<(multiplier<,'new'>)>,</pre>	<pre>sub('new'), sub('trace',index)</pre>	
Parameters		
arraydim {number}	Dimension of experiment	
lsfid {'n',-fn/2 to ni or fn/2}	Number of complex points to left-shift np FID	
phfid {'n',-360 to 360, in deg}	Zero-order phasing constant for np FID	

Table 38.	Add/Subtract	Experiment	Commands and	Parameters
-----------	--------------	------------	--------------	------------

spectral add/subtract mode (addi), performing noninteractive spectral addition (spadd) and subtraction (spsub), and taking the minimum of two spectra (spmin):



Typical entry into this menu is by clicking on Main Menu > Analyze > Add/Sub. The menu can also be opened by entering the command menu('addsub'). The commands activated by the buttons are described in more detail in the next section.

Noninteractive Add/Subtract

The buffer used as for the add/subtract experiment (exp5) is first cleared using clradd, then different FIDs or spectra are added to or subtracted from the accumulating total by the commands add, sub, spadd, and spsub.

Adding and Subtracting FIDs

The add and sub commands add and subtract the last displayed or selected FID to and from the contents of the add/subtract experiment, respectively. An optional argument allows the FID to first be multiplied by a multiplier (the default is 1.0). The parameters lsfid and phfid can be used to shift or phase rotate the selected FID before it is combined with the data in the add/subtract experiment.

A multi-FID add/subtract experiment with FIDs 1 and 2 can be created with the add or sub command using the 'new' keyword as follows:

 Create the add/subtract experiment with a single FID by entering the following commands from some experiment: clradd select(1) add Make the add/subtract experiment contain an array of two FIDs corresponding to the original FIDs 1 and 2 by entering: select(2) add('new')

If you had entered select (2) add in step 2, a single FID that is the sum of the original FIDs 1 and 2 would have been made instead of an array.

The arraydim parameter may need to be updated after constructing a multi-FID add/ subtract experiment. To do this, join exp5 by entering jexp5 and then enter setvalue('arraydim',numFIDs,'processed')

where numFIDs is the number of FIDs in that experiment. (e.g., if 12 FIDs were put into exp5, enter setvalue('arraydim', 12, 'processed')).

Individual FIDs in a multi-FID add/subtract experiment can be added to and subtracted from. The add and sub commands without a 'trace' keyword adds or subtracts from the first FID in the add/subtract experiment. Adding the 'trace' keyword followed by a required index number selects another FID to be the target of the add/subtract.

For example, select(4) add('trace', 6) takes the fourth FID from the current experiment and adds it to the sixth FID in the add/subtract experiment.

When using the 'trace' keyword, that FID must already exist in the add/subtract experiment by using an appropriate number of add('new') or sub('new') commands.

Adding and Subtracting Spectra

Noninteractive spectral addition and subtraction uses the spadd and spsub commands. The last displayed or selected spectrum is added to (spadd) or subtracted from (spsub) the current contents of the add/subtract experiment.

Each spectrum can be optionally multiplied and shifted using the multiplier and shift arguments, respectively. For example, entering spadd(0.75,10) multiplies the spectrum by 0.75 and shifts the spectrum by 10 to the left. A positive value of shift shifts the spectrum being added or subtracted to higher frequency, or to the left. A negative value of shift shifts the spectrum being added or subtracted to lower frequency, or to the right. To shift a spectrum without multiplying it, use a multiplier of 1.0.

A multi-element add/subtract experiment can be created with the spadd or spsub command. Using the keyword 'new' creates a new spectrum in the add/subtract experiment, for example:

- 1. Create the add/subtract experiment with a single spectrum in it by entering **clradd select(1) spadd** from some experiment.
- 2. Enter **select(2) spadd('new')** to make the add/subtract experiment contain an array of two spectra corresponding to the original spectra 1 and 2, respectively.

If you instead entered select(2) spadd, a single spectrum that is the sum of the original spectra 1 and 2 would have been made instead of an array.

Individual spectra in a multi-element add/subtract experiment can subsequently be added to and subtracted from. The spadd and spsub command without a 'trace' keyword adds to or subtracts from the first spectrum in the add/subtract experiment. Adding the 'trace' keyword followed by a required index number selects another spectrum to be the target of the add/subtract. For example, select(4) spadd('trace', 6) takes the fourth spectrum from the current experiment and adds it to the sixth spectrum in the add/ subtract experiment. When using the 'trace' argument, the indexed spectrum must already exist in the add/ subtract experiment by using an appropriate number of spadd('new') or spsub('new') commands. The results can be examined by joining experiment 5 (jexp5) and using the normal spectral display (e.g., ds) and plotting commands.

Interactive Add/Subtract

Interactive add/subtract provides a convenient way to manipulate two different spectra. A "spectrum" may be a 1D spectrum, a trace from a 2D spectrum, or even a spin simulated spectrum. Both horizontal displacement and vertical scale of the two spectra to be added or subtracted are under interactive control. If the spectra can be phased, they can be phased independently. The result can be manipulated using any of the standard software (e.g., the command p1), including further interactive add/subtract with another data set.

Interactive add/subtract is a multi-step process using the command addi. The same as noninteractive add/subtract, interactive add/subtract uses exp5 as an add/subtract buffer file, so no important data should be present in that experiment when you begin the process.

To Start Interactive Add/Subtract

Interactive add/subtract starts exactly the same as noninteractive add/subtract:

- 1. Enter **clradd** to clear the buffer.
- 2. Enter **spadd** to add in the current spectrum as the starting point.
- 3. Enter **addi** to start the interactive add/subtract program. The addi program can also be opened by clicking on the Interactive Mode button in the Add/Subtract Menu or by clicking on the Add/Sub button in the Deconvolution Menu.

Upon opening the addi program, a second spectrum is selected and the interactive process is started.

Spectrum 1, the spectrum selected by the spadd command, appears in the center of the display. Spectrum 2, the spectrum which was active when addi was typed, appears on the bottom. The sum or difference of these spectra appears on top of the screen; when addi is first entered, this spectrum is the sum (1 + 2) by default. The Interactive Add/Subtract menu (described in the next section) also appears.

Displayed at the bottom of the screen is the name of the currently active spectrum (the interactive one). Also displayed there is the current result mode (add, sub, or min).

Interactive Add/Subtract Menu

The Interactive Add/Subtract Menu has the following buttons (the labels on some buttons change depending on what mode you are in):



Each button functions as follows:

The first button is labeled Box or Cursor. When Box appears, you
are in the cursor mode, and clicking on this button changes you to
the box mode with two cursors.

Cursor When Cursor appears, you are in the box mode, and clicking on this button changes you to the cursor mode with one cursor.

Select Expand	Selects whether the current, add/sub, or result spectrum is active. The third button is labeled Expand or Full, depending on whether you are in the box or cursor mode. When Expand appears, you are in the box mode, and clicking on this button expands the area between the cursors.
Full	When Full appears, you are in the cursor mode, and clicking on this button displays the full area.
sp wp	Adjusts the start and width of the active spectrum.
sub	The fifth button is labeled sub, min, or add. When sub appears, clicking on the button makes the result spectrum to be the difference between the current and the add/sub spectra.
min	When min appears, clicking on the button makes the result spectrum to be a minimum intensity of either the current or the add/sub spectra.
add	When add appears, clicking the button makes the result spectrum to be a sum of the current and the add/sub spectra.
save	Saves the result spectrum in the add/sub experiment and returns to the last menu.
return	Returns to the last menu without saving the result.
This menu is n	not user-programmable

This menu is not user-programmable.

To Manipulate Spectra

Manipulation of the spectra is performed using the mouse:

- Left button positions the cursor or pair of cursors. If the sp wp button is selected, the left mouse button adjusts the start of the display.
- Center button changes the vertical scale of the spectrum so that it goes through the current mouse position. If the mouse is positioned at the left edge of the spectrum, the horizontal position of the spectrum is adjusted.
- Right button positions the second cursor relative to the first cursor. If the sp wp button is selected, the right mouse button adjusts the width of the display.

The important points to understand are few. The Select button is used to toggle between different modes of control. When the label at the bottom of the screen reads "active: current," all of the parameters (except wp) control spectrum 2, and spectrum 2 can be phased, scaled, or shifted relative to spectrum 1.

After clicking on the Select button, the label at the bottom of the screen reads "active: addsub," and all of the parameters except wp control spectrum 1.

Clicking on Select again toggles the label to read "active: result," so now parameter changes affect only the sum or difference spectrum. Note that wp always controls all spectra, since differential expansions of the two spectra are not supported. Note also that the colors of the labels change to match the colors of the different spectra.

The sum/difference spectrum displayed on the screen while addi is active is strictly a temporary display. Once all manipulations have been performed, and assuming the sum/ difference is something you wish to perform further operations with (such as plotting), it must be saved into the add/subtract buffer (exp5) by clicking on the save button. At this point spectrum 1, which was in the add/subtract buffer, will be overwritten by the sum or difference spectrum, and addi will cease operation.

In most cases, you enter jexp5 ds to display the difference spectrum on the screen, ready for further manipulation (expansion, line listing, etc.) and plotting. If you wish to continue with the add/subtract process by adding in a third spectrum, display that spectrum in the usual way and enter addi again.

6.5 Regression Analysis

The process of establishing correlations between two or more variables is called *regression analysis* or *correlation analysis*. The established regression or correlation can then be used to predict one variable in terms of the others. Often, paired data indicate that a regression may have a certain functional form, but we do not want to make assumptions about any underlying probability distributions of the data.

This type of problem is often handled by the least squares curve-fitting method. Specific examples of this were used for the analysis of T_1 and T_2 NMR data and for the analysis of kinetics data. Also available within VNMR are tools for fitting arbitrary data to selected functional forms.

The regression process in VNMR takes a set of data pairs from the file regression.inp and attempts to fit a curve to the set. The implemented curves are first, second, and third order polynomials and an exponential in the form:

y = a1 * exp(-x/tau) + a3

There are further possibilities as the original data may be displayed against a choice of linear, squared, or logarithmic scales.

Regression Commands and Menus

Table 39 lists the commands associated with regression analysis.

Commands			
analyze*	Generalized curve fitting in regression mode		
autoscale	Resume autoscaling after limits set by scalelimits		
expfit*	Make least-squares fit to exp. or poly. curve (UNIX)		
expl<(<options,>line1,>)></options,>	Display exponential or polynomial curves		
<pre>pexpl<(<options,>line1,>)></options,></pre>	Plot exponential or polynomial curves		
poly0	Display mean of data in regression.inp file		
rinput	Input data for regression analysis		
scalelimits*	Set limits for scales in regression		
<pre>* analyze('expfit',xarray<,option,option,>)</pre>			
expfit options <analyze>analyze.list</analyze>			
<pre>scalelimits(x_start,x_end,y_start,y_end)</pre>			

Table 39. Regression Commands

To remove the need to determine the correct arguments for the commands involved, especially expl and analyze, regression can be performed almost completely through the following menus.

Regression 1 menu:

x-linear x-squa	re x-log	dp-linear	dp-quad	dp-exp	Next	Rt	
-----------------	----------	-----------	---------	--------	------	----	--

Regression 2 menu:



These menus are accessible by selecting Analyze in the Main Menu, and then selecting Regrs from the Analyze Menu.

The actions of the buttons 1 through 6 on each menu are described below. The next button in the Regression 1 selects the Regression 2 menu, the next button in the Regression 2 selects the Regression 3 menu, selects the next regression menu. The Return button in the Regression 3 menu returns to the Regression 1 menu.

Regression Analysis Step-by-Step

- 1. Write and save the text file regression.inp that contains the data pairs you wish to analyze. The next section describes the format of this file. Create it by one of the following methods:
 - Use the macro rinput. This program displays a series of prompts requesting the axis label titles and the data pairs. When you finish, rinput creates the file regression.inp in the correct format. rinput does not allow you to edit the data; correct errors using a text editor after you complete the program.
 - Use a text editor such as vi or textedit.
 - Create a MAGICAL II macro for this purpose.
- 2. Enter the command **expl('regression')**. Alternatively, select any one of the buttons x-linear, y-linear, x-square, y-square, x-log, and y-log. Each of the buttons also scales for displays and plots either the x or y axis as labeled on the button. In the case of multiple data sets in the input file, data sets may be selected with the command expl('regression1', line#, line#...).

The expl command uses the values in the regression.inp file to display a graph of the data points. It also creates the files analyze.inp (needed by analyze to run the analysis) and expl.out (display information for expl).

You might also want to use the poly0 macro to calculate and display (as horizontal lines) the mean of the data in the file regression.inp.

3. Run the analyze('expfit', 'regression', option, 'list') command, where option is 'poly1', 'poly2', 'poly3', or 'exp', and then enter expl to see the results as a graph. Alternatively, click on one of the buttons dp-linear, dp-quad, dp-cubic, or dp-exp. These buttons include displaying the results using expl.

analyze('expfit', 'regression', option, 'list') calls the UNIX
program expfit, which creates the files analyze.out (used by expl to display
the results) and analyze.list (a table of results). The type of fitting is
determined by analyze option or button you provide:

Fitting	analyze option	Button
linear	'poly1'	dp-linear
quadratic	'poly2'	dp-quad

Chapter 6. Data Analysis

Fitting	analyze option	Button
cubic	'poly3'	dp-cubic
exponential curve	'exp'	dp-exp

The menu system allows immediate display of the regression results and selection of another type of regression if results are not satisfactory. Figure 58 shows quadratic fittings for the data given in the example of the regression.inp file in the next section.




Figure 58. Display of Regression Fittings (expl Program)

To plot the analysis, select the **plot** button in the Regression 3 menu or enter **pexpl** page. To plot selected lines, enter **pexpl** (index#, index#, ...).

Because expl and pexpl set the scale automatically to show all points (if possible), you may want to enter the **scalelimits** macro to set limits for the scales. Entering scalelimits with no argument starts an interactive process in which you are prompted for the four scale limits. You can also enter the limits for the x-axis and y-axis as arguments: scalelimits(x_start, x_finish, y_start, y_finish). The limits are retained as long as an expl display is retained. To return to automatic scaling by expl, enter **autoscale**.

5. To show the results of the analysis in the text window, select the **dp output** button from the Regression 3 menu or enter **cat(curexp+'/analyze.list')**.

Contents of "analyze.out" File

The data input file is analyze.out, except for regression when the input file is regression.inp. The file expl.out saves certain display and plot parameters.

Values can be 2048 points maximum from a data set, 2048 points maximum from all sets displayed/plotted, 8 data sets maximum displayed, and 128 data sets maximum are read.

The following is an example of analyze.out. Numbers **0**, **2**, etc. identify lines in the example and are not part of the actual file:

- exp 7 regression
- ❷ D1 C0 C1 C2

€	1 5 lin	near linear
4	Exponential	Data Analysis
	time	
	amp	
6	NEXT 5	
6	1 -248.962	22.8025 226.157
0	2	4
	3	9
	4	16
	5	25
	6	36
	6	36

The following is a description of the numbered parts of this file:

1 The keyword exp is followed by a number for a curve type from Table 40.

Туре	Function	Functional Form
0	T_{1}/T_{2}	(a0-a2)*exp(-t/a1) + a2
1	Increasing kinetics	a0*exp(-t/a1) + a2
2	Decreasing kinetics	-a0*exp(-t/a1) + a2 + a0
3	Diffusion	a0*exp(-D1*uu) + a2*exp(-a1*D1*uu) where $uu=C0+C1*t+C2*t^2$
4	None	No theoretical curve (use 'link')
5	Linear	a0 + a1*t
6	Quadratic	$a0 + a1*t + a2*t^2$
7	Exponential	a0 * exp(-t/a1) + a2
8	Contact time	(a3 - (a3 - a0) exp(-t/a1)) exp(-t/a2) + a0
9	Cubic	$a0 + a1^{*}t + a2^{*}t^2 + a3^{*}t^3$

Table 40. Curve Types.

The keyword regression, if present, indicates regression output

Ploating point constants *D1*, *C0*, *C1*, and *C2*, if present, are used only with the diffusion function (curve type 3).

• An integer for the number of data sets (curves), followed by an integer for the number of data point pairs in the set. For regression, the words are scale types for the *x* and *y* axes: linear, square, and log.

Q Title line. Use No Title when a title is not desired. Two additional text lines for the x and y axis titles are present in output from regression.

• The keyword NEXT identifies the start of a data set, and the integers that follow give the number of data point pairs in the data set.

6 The first integer specifies the number (usually 1) of the data point symbol used for the data set. The next three integers are the coefficients a0, a1, and a2 (see Table 40) and must all be present, even for functions that do not use all three (e.g., first-order polynomial). If a particular number is not appropriate, put any number there. a3 must be also be present for cubic and contact time functions (curve types 7 and 8).

7 Data point pairs in the set.

Next is an example of regression. inp file for generalized curve fitting:

① time ② amp ③ 0 0 ④ NEXT ⑤ 2.000000 4.000000 3.00000 9.000000 4.000000 16.000000 5.000000 25.000000 6.000000 36.000000

Description of this example:

• Line with text for *x*-axis label displayed by expl('regression').

Line with text for y-axis label (line must not be too long, usually less than 20 characters).
 The first non-blank character must not be a digit.

Subscription Line containing an integer for the number of data sets followed by another integer for the number of pairs per data set. Both values are 0 if the number of pairs is variable.

• A line beginning with the keyword NEXT is inserted at the start of each data set when the number of pairs per peak is variable.

5 The data pairs, listed one pair to a line.

The final example uses the 'file' argument to the expl command:

Description of this example:

1 Keyword exp followed by curve type number.

2 Number of data sets, followed by number of data point pairs.

3 Title.

Data point symbol number, followed by three coefficients.

5 Data point pairs in the set.

Contents of "regression.inp" File

The data input text file regression.inp contains a listing of axis labels and data pairs. The data file can contain up to 128 data sets. Data sets are selected by expl indexes (up to 6, depending upon length of data sets) with a default to the beginning data sets. The analysis is limited to 1024 data points, with the first part of larger data sets selected.

The following is an example of a regression.inp file that shows the format used. Numbers **0**, **2**, etc. identify lines in the example and are not part of the actual file:

0 time

0 temp

6 0 0

```
a
    NEXT
6 0.000000
             5.000000
  1.000000
           9.00000
  2.000000
             12.000000
  3.000000 20.000000
  4.000000
           33.000000
  5.000000
              42.000000
O NEXT
6 0.000000
           20.000000
  1.000000 22.000000
           27.000000
  2.000000
             29.000000
  3.000000
  4.000000
              37.000000
  5.000000
              45.000000
```

Description of this file:

• X-axis label for display by expl ('regression'). The label is optional. If used, the first non-blank character in the label must not be a digit.

Y-axis label for display by expl ('regression'). The label is optional. If used, it must be less than about 20 characters and the first non-blank character in the label must not be a digit.

③ If the number of pairs per data set is fixed, this line contains an integer for the number of data sets, followed by another integer for the number of data pairs per data set. If the number of pairs per data set is variable (as in this example), both integers are set equal to 0.

(9 If the number of pairs per data set is variable, a line with the word NEXT is inserted at the start of each data set.

• Data pairs, one to a line, are listed for each data set, in this order: first pair of first set

second pair of first set third pair of first set

. . .

first pair of second set second pair of second set

6.6 Chemical Shift Analysis

Chemical shifts can be analyzed using the commands pcss and do_pcss listed in Table 41. When you use pcss and do_pcss, the list of chemical shifts that is created is saved in the file pcss.outpar.

Commands	
do_pcss	Calculate proton chemical shifts spectrum
pcss	Calculate and show proton chemical shifts spectrum

Chapter 6. Data Analysis

Sections in this chapter:

- 7.1 "Pulse Shape Analysis," this page
- 7.2 "Pandora's Box," page 228

7.1 Pulse Shape Analysis

The pulsetool program is designed to display and examine shaped rf pulses. The standard pulse template file format is the same as for shaped pulses in /vnmr/ shapelib. Data points are listed as phase amplitude time-count, where phase is in degrees, amplitude is a value between 0 and 1023, and time-count is an integer which describes the relative time duration of the step. The program is started by entering the command pulsetool in a UNIX window. Table 42 summarizes the command and parameters associated with pulse shape analysis.

Table 42. Pulse Shape Analysis Commands and Parameters

Command pulsetool <-shape filepath>	RF pulse shape analysis (UNIX)
Parameters phi theta	Amount of rotation about the Z-axis Declination relative to XY-plane

The amplitude and phase are displayed in the small windows at the top of the display, along with the effective frequency of the pulse, the quadrature components of the pulse, and its Fourier transform. You can select the contents of any of the smaller windows for display in the large graphics window in the center of the screen.

Between the small windows at the top of the display and the large central, graphics window is the control panel, home to a number of buttons that perform various operations or activate the routines described below.

Below the main graphics window is a panel that contains miscellaneous information about the current pulse and display status. The directory file name, pulse name, vertical scale, and vertical reference fields display current information that can be altered by the user.

The Steps, Fourier size, Power factor, and Integral fields are advisory only, and may not be entered or changed by the user. Power factor is calculated when a pulse is loaded, and is the mean square amplitude of the pulse. A square pulse has a power factor of 1. The integral of

the pulse is an attempt to calculate the tip-angle per unit time and B_1 field strength. This number is strictly valid only for pulses that are modulated in amplitude only and can be used to determine the B_1 field required to obtain, for example, a 90° tip with a sinc pulse. To do this, divide the desired tip angle (in revolutions) by the product of the integral value and the pulse length (in milliseconds). The result is the required B_1 field strength, in kHz.

The directory system may be viewed, and pulse files selected for loading through the use of the Files button.

Simulation of the actual response to a pulse, based on Bloch equation calculations, is available by selecting the Simulation button.

A number of standard pulses can be created, with attributes tailored through the Create utility. The data currently displayed in the main graphics window can be saved with the Save button.

Directory and File Operations

In selecting files, both the working directory and the pulse template file name can be specified by direct entry into the Directory and Pulse name fields found in the panel at the bottom of the display (use the Delete key to erase characters, if necessary, and type in the desired name, followed by Return to indicate completion). When the Pulse name field is selected, pressing Return causes the named file to be loaded and displayed.

Alternatively, the Files button causes a popup window to be displayed, listing the contents of the current directory. A trailing slash "/" following a member of the list indicates a subdirectory, and an asterisk "*" an executable file. The Load, Chdir, and Edit buttons operate on an item selected from this listing with the left mouse button:

The Load button causes the selected file to be read, and displayed in the graphics windows. If the file does not correspond to the proper format for pulse template files, an error message is displayed. Comment lines beginning with the pound character "#" are ignored. Descriptive information about the pulse is displayed in the bottom panel—the name of the file, the number of steps in the pulse, the Fourier size required to do the FFT of the pulse, and a "power factor" calculated for the pulse. The power factor is based on the mean square amplitude of the pulse.

The Chdir button changes to and then lists the selected directory.

The Parent button changes to and then lists the parent of the current directory.

The Save button located in the main control panel can be used to save data currently displayed in the main graphics window to a file. When this button is selected, a second button labeled Done appears, along with a type-in field that holds the name of the file that will be created. First, enter an appropriate name, then select the Save button once again to write the file. Once you have entered the Save mode, you can repeat this as many times as you like— display a different attribute in the main window, enter a new file name, and select Save. To exit from this mode, select Done.

The Print button located in the main control panel can create a file that can be used to print the main graphics window on a PostScript printer. When this button is selected, a second button labeled Done appears, along with a type-in field that holds the name of the file that will be created. Selecting the Print button with a appropriate file name in the type-in field writes the file. The file can subsequently be sent to a PostScript printer with the UNIX lp command.

Attribute Selection

The six small graphics windows at the top of the tool initially display the different attributes of the current pulse:

- Amplitude
- Phase
- Effective off-resonance frequency
- · Real and imaginary quadrature components
- Fourier transform

Any of these six windows can be displayed in the large graphics window by clicking in the appropriate small window with either the left or middle mouse buttons:

- The left mouse button causes the large window to be cleared before drawing and sets the clear mode to on.
- The middle mouse button turns off the clear mode and displays the selected attribute, overlaying any current display in the large graphics window.

Repeated selection of the small Fourier transform window will result in the large window cycling through the magnitude of the Fourier transform, the real component, and the imaginary component.

Scale and Reference

The vertical scale can be adjusted either by clicking the middle button inside the boundary of the large graphics window or by manually entering a value in the Vertical scale field of the bottom panel, ending by pressing Return. Using the middle mouse button causes the scale to be adjusted interactively so that the active curve passes through the mouse arrow. Note that no rescaling occurs if the y-value specified with the middle button does not have the same sign as the actual attribute value at that point on the x-axis. A negative value can, however, be entered as a vertical scale if so desired.

The vertical reference controls the vertical position of the active curve on the large graphics window, representing the offset from zero measured in y-axis units. A positive value moves the curve up, and a negative value moves it down. Like the vertical scale, the vertical reference can be adjusted in one of two ways—a value may be entered manually into the Vertical reference field in the bottom panel, or the middle mouse button can be used interactively anywhere in the large graphics window, while simultaneously holding down the Control key. In the second case, the vertical reference is set so that the curve passes through the mouse arrow.

The vertical scale and reference are reset whenever an attribute is selected from any of the small graphics windows. If things get out of hand, use this by reselecting the current small window with the left mouse button.

Cursors

Interactive left, right, and horizontal cursors are available, and display a readout of position at the bottom of the large window when active. The left cursor is activated by clicking the left mouse button inside the large window. When the left cursor is present, the right cursor can be activated by clicking the right mouse button anywhere to the right of the left cursor. At this point, the right mouse button controls the position of the right cursor independently, while the left mouse button moves both cursors in tandem.

When both cursors are active, the control panel button normally marked Full will read Expand, and can be used to display an expanded view of the region selected between the two cursors. (Note that the clear mode will always be set to on after an Expand or Full operation.) The left and right cursors are turned off by clicking the appropriate mouse button in the large window while simultaneously pressing the Control key.

The horizontal cursor is activated with the Thresh button located in the control panel. When this cursor is active, it is controlled interactively with the middle mouse button. The interactive scale and reference functions normally controlled with the middle mouse are not available when the horizontal cursor is present. Select the Scale button in the control panel to turn off the horizontal cursor and reactivate the scale and reference functions (vertical scale and reference can be adjusted even with the horizontal cursor active by direct entry in the appropriate fields in the bottom panel).

Simulation Overview

The simulation routine simulates the effects of an rf pulse by use of the classical model of nuclear spin evolution described by the Bloch equations. T_1 and T_2 relaxation effects are ignored, in which case the evolution of a magnetization vector in the presence of an applied rf magnetic field can be evaluated by multiplication with a 3 by 3 rotation matrix. The simulation consists of repeated multiplication by such a matrix, whose elements are determined at each step by the values of amplitude and phase found in the pulse template file, and by user input values of initial magnetization, B_1 field strength, pulse length, and resonance offset. The simulation is performed over one of three possible independent variables— resonance offset, B_1 field strength, or time, and is determined by the Sweep cycle in the small button panel.

Simulation Parameters

Select the Simulation button in the control panel to activate the Bloch Simulation subwindow. This window consists of a panel containing all required parameters (the pulse length is taken from the value in the bottom panel of the main window) and a small button panel at the bottom of the window. To change the value of any parameter, select it with the left mouse button, then delete the appropriate characters and enter the desired value from the keyboard. Parameters are updated each time the Go button is selected or when the Steps button is selected with Index equal to zero.

The first three parameters in the left hand column describe the starting values for the magnetization components Mx, My, and Mz, whose vector sum must be less than or equal to 1.

The next three fields change to reflect the state of the Sweep cycle, which can be toggled between B_1 , Freq, and Time. When Freq is selected, the first of these fields will read B_{1max} , the value of B_1 at the maximum pulse amplitude. The second and third fields determine the lower and upper off-resonance frequency boundaries. When the Sweep cycle is set to B_1 , these three fields are reversed so that the first determines a constant off-resonance frequency and the remaining two determine the lower and upper boundaries of the maximum B_1 amplitude. Selecting Time will yield a display of the magnetization as a function of progression through the pulse, at the frequency and B1 field strength specified by the parameter values displayed. In this last case, the number of steps in the simulation is taken from the number of points in the pulse template and may not be altered externally. To get finer resolution in the simulation, use a pulse template with a greater number of steps.

The Initialize cycle determines if the magnetization is reinitialized to the values of Mx, My, and Mz, or if the simulation uses the values at each point that were the result of the previous simulation. In this way, the effect of a series of pulses can be evaluated by loading the pulse and performing the simulation with Initialize set to Yes, then loading the next pulse, setting Initialize to No, and selecting Go. Any number of pulses can be concatenated in this fashion. This feature works only for Frequency and B_1 sweep, but not Time.

When Time sweep is selected, the results can also be displayed in the form of a projected three-dimensional coordinate system, showing the path of the magnetization over the course of the pulse. This display is obtained by selecting the 3D button after first selecting the Go button. When the 3D display is active, the left mouse button controls the viewing angle from within the canvas region delineated by the blue corner markers. This viewing angle is described by the two parameters phi (the amount of rotation about the Z-axis) and theta (the declination relative to the XY-plane). A "family" of trajectories can be displayed by first selecting any of the small canvases with the middle mouse button, then selecting the 3D button. Changing either the B_1 field strength or the resonance offset followed by the Go button will result in display of the result without clearing the display. To reactivate the automatic clearing feature, select any of the small canvases with the left mouse button. To see the 3D display drawn in real-time, enter a nonzero integer value in the Time field. The appropriate value depends on the number of steps in the pulse and the type of computer you have. Try a value like 100 for a SPARCstation.

The last parameter in this column determines the number of points at which the simulation will be performed along the y-axis. A larger number will give more detail in the result, but will require proportionally more time to complete.

The Index parameter is a counter that updates the status of the simulation, and cannot be set externally. The value displayed is the number of steps in the pulse template that have been completed.

The Step Inc parameter is used by the Step button, described in the section, "Performing a Simulation," to control the number of intermediate steps to be calculated.

Performing a Simulation

When you have adjusted the parameters to your liking, you will probably want to select the Go button. This does simulation calculations and then displays the results in the first five small graphics windows, replacing (but not destroying!) the pulse information that was displayed there. The Fourier transform information remains unaffected, so that comparisons can be made between this and the exact simulation results.

All of the display functions described elsewhere are active as well, with the simulation data. Additionally, the original pulse data is still present in the background and can be swapped into view with the Display cycle found in the main control panel.

The Step button offers the ability to view the course of the magnetization at intermediate stages through the pulse. When this function is selected, the next Steps Inc steps of the pulse are simulated, starting at the current value of Index. The intermediate result is then displayed in the normal fashion.

During a Go simulation, a small panel containing a Cancel button will pop into view. Use this to stop the simulation if necessary (there may be some delay between selecting the button and the end of the process; it won't do any good to click on Cancel more than once).

Creating a Pulse

The pulse creation routine currently offers the following pulse types:

Square	Hermite 90	Tan swept inversion
Sinc	Hermite 180	Sin/cos 90
Gaussian	Hyperbolic secant inversion	

A file containing the pulse template for any of these pulses can be created from scratch with this utility. Alternatively, pulses can be created for examination only, using the display capabilities of pulsetool. Each pulse is generated with user-definable parameters appropriate for the pulse in question.

When the Create button is selected, a menu of pulse types appears. Hold the right mouse button down on the Create button, select one of the pulses in the resulting menu, and release the mouse button. If you decide you don't like any of the possibilities, move the mouse arrow out of the menu area and release the button. When a pulse type is selected, a small window appears with a brief description of the characteristics of the pulse and a set of changeable attributes whose values you may alter if so desired. The number of steps in the pulse is limited to powers of 2 and can be set by clicking the left mouse button, or by holding the right mouse button down and selecting the desired value from the resulting menu. All other attributes, which vary depending on the pulse type, can be altered from their default values by first selecting the appropriate field with the left mouse button, deleting with the Delete key, and typing in the desired value (pressing Return is not required).

At this point, you may select one of the three buttons at the bottom of the window: Preview, Execute, or Done:

- Preview uses the attributes as they appear on the screen to create a pulse that is loaded internally into pulsetool. All pulsetool features can then be used to examine and evaluate the new pulse. Any previous pulse information is deleted.
- Execute uses the attributes as they appear on the screen to create a pulse, which is written to a standard UNIX file. The name of the file is taken from the file name field in the Create window and written into the current directory, listed in the Directory field in the bottom panel. If a file of the same name already exists, you are asked to confirm your request. If, for any reason, the program is unable to write to the named file, an error message appears. This is generally symptomatic of not having write permission in the current directory.

Currently, there is no convenient way for a user to add new pulse types to those listed above. Suggestions for those pulse types that should be included in the future are welcomed. However, any user-created shaped pulse may be examined using the Files button.

7.2 Pandora's Box

Pandora's Box (Pbox) software creates shape pattern files for experiments involving shaped rf pulses, composite pulses, decoupling and mixing patterns, adiabatic rf sweep waveforms, and pulsed field gradient shapes. The goal of Pbox is to simplify generation and use of different waveforms in NMR experiments to a level where the user does not need to be an expert in selective excitation. Pbox makes the use of complex waveforms as simple as using ordinary rectangular pulses. Indeed, not only does Pbox provide all the necessary parameters (pulse width, power, dmf, dres, etc.) when the shape files are created, but this information can be extracted at any time from Pbox shape files by macros or directly within pulse sequences. More than 160 different shapes are available from the Pbox library.

Getting Started

The simplest mode of operation is from Pbox menus:

- 1. Enter **ds** to display a spectrum.
- 2. Click on **Pbox**.

The Pbox menu system appears:

These buttons have the following actions:

Wavelib	Selects the wavelib directory and provides file system utility for selection of shapes.
90	Prompts for cursors to be placed around the signals to be excited and displays the Pbox menu for the definition of 90° excitation pulse shapes.
180	Prompts for cursors to be placed around the signals to be excited and displays the Pbox menu for definition of 180° inversion/refocusing pulse shapes.
Homo-dec	Prompts for cursors to be placed around the signals to be homodecoupled and displays the Pbox menu for definition of decoupling waveforms.
Het-dec	Displays the Pbox menu for the definition of decoupling waveforms. By default, it assumes ¹ H decoupling over 10 ppm
Mixing	Displays the Pbox menu for the definition of mixing waveforms. By default, it assumes bandwidth of 10 ppm (^{1}H) .
Grad	Displays the Pbox menu for the definition of gradient shapes.
Quit	Exits Pbox and displays the 1D Interactive Display menu.

The following steps are typical for creating an excitation pulse:

- 1. Click on **Pbox** > 90.
- 2. Select an excitation band using cursors.
- 3. Click on e-Burp1 > Close > Name > Close.

For help in understanding Pbox menus, click on the Help button.

Calibrating the RF Field

To obtain the pulse calibration numbers in the Pbox output, provide the rf field calibration data ref_pwr and ref_pw90 in the input. If menus are used for the first time in the current experiment, Pbox prompts you to provide the necessary information. Therefore, before waveform creation, make sure the rf field has been calibrated and you know the length of the 90° pulse at a given power level.

If the spectrometer amplifiers are linear, which is standard on Varian NMR spectrometer systems, it does not matter at what power level the rf field is being calibrated. However, for maximum accuracy, do the calibration close to the field used in the experiment. An estimate of the rf field is obtained by providing approximate calibration data and using cal as an output file name. No waveform is created in this case, and only the calibration results appear in the output.

Creating Waveforms from Macros

Pbox macros provide useful tools for customizing NMR experiments. The simplest way to create a shape is using the pxshape macro. For example, a single band excitation pulse using the E-BURP-1 shape, covering 400 Hz, and shifted off-resonance by -880 Hz from the carrier frequency (middle of the spectrum) can be created and stored in the alpha.RF file as follows:

pxshape('eburp1 400.0 -880','alpha.RF')

The following steps are necessary to create multiply-selective pulses. If the spectrum of interest is on the screen, use the cursors.

- 1. Enter **opx('hadamard.RF')** to open the Pbox.inp file and write the file header.
- 2. Select an excitation band using cursors.
- 3. Enter selex('rsnob').
- 4. Select the second excitation band using cursors.
- 5. Enter selex('rsnob').
- 6. Repeat steps 2 to 5 as many times as needed.
- 7. Enter **cpx** to close the Pbox.inp file.
- 8. Enter **dshape** to display the last created shape.

If an experimental spectrum is not available, the following slightly different set of macros are used:

- 1. Enter **opx('myshape')** to open Pbox and provide a file name.
- 2. Enter setwave('sech 400.0 -880.0') to select first band at -880 Hz.
- 3. Enter setwave('sech 400.0 1240.0') to select second band at 1240 Hz.
- 4. Enter cpx(ref_pw90, ref_pwr) to close Pbox.
- 5. Enter **dshape** to display the shape file.

The pbox_pw and pbox_pwr macros are used to load the parameters of the last created shape file into the current experiment:

pbox_pw:selpw
pbox pwr:selpwr

Alternatively, the calibration data is directly retrieved from the shape file provided as an argument to the pbox_dmf and pbox_dres macros:

```
pbox_dmf('ccdec.DEC'):dmf
pbox_dres('ccdec.DEC'):dres
```

where ccdec.DEC is the name of the decoupling shape file.

The excitation profile of shaped pulses is conveniently verified using the Pbox Bloch simulator:

- 1. Enter **opx** to open Pbox.
- 2. Enter setwave('iburp2 400.0 -880') to select first band at -880 Hz.
- 3. Enter setwave('iburp2 400.0 1240.0') to select second band at 1240 Hz.
- 4. Enter **pbox_rst** to reset par-s and write comments.
- 5. Enter **pboxpar('name', 'test.DEC')** to define the output file name.

- 6. Enter **pboxpar('bsim', 'y')** to activate the Bloch simulator.
- 7. Enter **cpx** to close Pbox.
- 8. Enter **dshape** to display the shape file.
- 9. Enter **dprofile('z')** to display inversion (Mx) profile.

In the vast majority of cases, you don't have to tell Pbox whether you are creating a 90° excitation pulse or 180° inversion pulse, or even whether an .RF, .DEC, or .GRD type of waveform is needed, because this information is stored with the corresponding wave file in the wavelib directory. Pbox can be forced to change the waveform type by you simply providing the required extension to the output shape file name. Wave files are modified, if necessary, by copying it into your wavelib and editing the text file as required. See "Pbox Macro Reference," page 241, for a more complete description of macros.

Creating Waveforms from UNIX

It is sometimes more convenient to create waveforms from the UNIX shell:

> Pbox

I

I

The name of the output shape file is passed as the first argument:

```
> Pbox filename
```

The input data are typically stored in the Pbox.inp file in your vnmrsys/shapelib directory and are modified using standard text editors. Alternatively, most of the necessary data can be provided as arguments to the Pbox command. For example,

```
> Pbox myfile -w 'eburp1 480 -1200' -p 40 -1 104
```

generates an E-BURP-1 excitation pulse covering 480-Hz-wide band and shifted -1200 Hz off-resonance using for calibration 104 μ s long pw90 at 40 dB power level and stored in myfile.RF. Note that the name of the output shape file is always passed as the first argument.

Several other options are accepted by Pbox; for example, -b activates the Bloch simulator, -c calibrates the waveform without creating a shape file, and -o prints out the available options. (see "Pbox UNIX Commands," page 258, for further information).

Pbox File System

All the information about the waveform to be created (e.g. calibration data, output file name, excitation band definition) is stored in the Pbox.inp text (ASCII) file in the user directory vnmrsys/shapelib. This file is generated whenever Pbox menus or macros are used. You can also create it by using one of the standard text editors.

Any shape file can consist of one or several shaped pulses that are combined into a single waveform. Each excitation band is defined by a wave definition string (a string of wave variables enclosed between delimiters {and}). The number of wave definition strings in a single Pbox.inp file is unlimited. In order to simplify the input file format, the wave variables are entered without names in a strongly predefined order:

sh bw(/pw) ofs st ph fla trev d1 d2 d0 wrap php

The following list describes each of the variables.

sh	Shape name as stored in wavelib
bw(/pw)	Bandwidth in Hz, or pulsewidth in sec, or both
ofs	Offset from transmitter offset or carrier in Hz

I

I

st	Spin status (0 for Mz or 1 for Mxy)
ph	Phase or phase cycle
fla	Flip angle
trev	Time reversal flag
d1	Prepulse delay
d2	Postpulse delay
d0	Delay before the pulse
wrap	Wraparound parameter
php	Phase pivot

The order of parameters has been chosen such that the importance of parameters is decreasing and rarely used parameters can be omitted or defaulted by assigning a value of n (not used). The following examples are valid wave definition strings.

{qsneeze}	q-SNEEZE pulse applied on resonance, the pulse length will be internally defaulted to 5 ms.
{G3 800}	G3 pulse covering bandwidth of 800 Hz and applied on-resonance.
{sech 400/0.05 -1200}	50 ms long hyperbolic secant pulse covering 400Hz and shifted off-resonance by -1200 Hz.
{WURST2 2k/5m 12k n t5}	5 ms long WURST-2 decoupling pulse covering 2 kHz and shifted off-resonance by 12 kHz uses t5 phase cycle.
{eburp1 450 0.0 n 180} {eburp1 450 820 1 0.0}	Two E-BURP-1 pulses mixed in a single waveform, both covering 450 Hz wide band. The first pulse is applied on-resonance with a phase of 180°. The second pulse is shifted to 820 Hz of-resonance, has zero phase and is a de-excitation pulse (status 1). By default such a pulse is time reversed.

A set of Pbox parameters can be used to define the waveform to be generated. The syntax of the Pbox.inp file is straightforward, parameter=value, for instance, name= myshape.RF, or simply name=myshape. The following list describes Pbox parameters and their default values (see "Pbox VNMR Parameters," page 236, for more details):

name=Pbox	Shape file name, the extension is optional.
type=r	Shape type, r - RF, d - DEC, g - GRD.
dres=9.0	As in VNMR, degrees. The default value is stored in wavefile.
steps=200	Minimum number of steps (< 64k). The default value is stored in a wave file.
maxincr=30	Max phase incr, degrees (<<180)
attn=i	Attenuation, i (internal), e (external) or d (nearest dB step)
sfrq=0	Spectrometer frequency, MHz
refofs=0	Reference offset, Hz (/ppm)
sucyc=d	Super Cycle, d (default), n (no), name as in wavelib/ supercycles. The default value is stored in wavefile.
reps=2	Amount of reports (0-4)
stepsize=n	Size of a single step (ms)
wrap=0	Wraparound parameter (0-1)
header=y	Shape header, y (yes) n (no) i (imaging)
bsim=n	Bloch simulation, y (yes), n (no), a (add), s (subtract), 200 (time in sec)
T1=n	Relaxation time T1 (sec)

T2=n	Relaxation time T2 (sec)
dcyc=1	Duty cycle (0 - 1).
sw=0	Spectral width (Hz)
padl	Zero padding; number of zero steps at the beginning of a waveform
pad2	Zero padding; number of zero steps at the end of a waveform
ptype=selective	pulse type (for imaging only).

The number and order of input parameters is optional and not important.

1

I

You can redefine the internally defaulted Pbox parameters by entering the default values in the .Pbox_globals file.

Parameters describing software and hardware limitations are also pre-defined internally and can be redefined by the user in the .Pbox_globals file that is stored in user's home directory. The following list describes global parameters and their default values.

shdir=\$HOME/vnmrsys/shapelib/	Default shape directory
wvdir=/vnmr/wavelib	Default wave directory
maxst=65500	Maximum number of steps in waveform
defnp=100	Default number of steps
minpw=0.2	Minimum step length, in µs
minpwg=2.0	Minimum gradient step length, in µs
drmin=1.0	Minimum dres
maxamp=1024.0	Maximum amplitude
maxgr=32767.0	Maximum gradient amplitude
amres=1.0	Amplitude resolution
phres=0.1	Phase resolution, in degrees
tmres=0.05	Time resolution, in µs
dres=9.0	Default dres
maxpwr=63	Maximum power level, in dB
minpwr=-16	Minimum power level, in dB
maxitr=5	Maximum number of iterations
maxdev=2.0	maximum deviation, in percent
cmpr=y	Waveform compression
minstps=64	Minimum steps in Bloch simulation
pw=0.005	Default .RF and .DEC pulse length, in sec
pwg=0.001	Default .GRD pulse length, in sec

The parameters of individual shapes—Gaussian, E-BURP-1, or hyperbolic secant pulse, etc.—are stored in the wavelib directory, which has several subdirectories, such as excitation, inversion, refocusing. Every individual shape is defined by a set of parameters that can be grouped in several categories.

Wave definition parameters are the following:

amf	Amplitude modulation function
fmf	Frequency modulation function
pmf	Phase modulation function
su	Default supercycle
fla	Default flip angle on resonance
pwbw	Pulsewidth to bandwidth product
pwbl	Pulsewidth to B1max product

pwsw	Pulsewidth to sweepwidth product
adb	Adiabaticity on resonance
ofs	Offset of excitation bandwidth
dres	Default tipangle resolution, in degrees
dash	Dash variable
wf	Window function
st	Default status
dutyc	Duty cycle
c1	Constant
c2	Constant
c3	Constant
steps	Default number of steps

Wave truncation parameters are the following:

min	Minimum truncation threshold (0 to 1)
max	Maximum truncation threshold (0 to 1)
left	Truncation from left (0 to 1)
right	Truncation from right (0 to 1)
cmplx	Flag, retain real (1), imag (-1) or complex(0) part of wave
wrap	Wraparound factor (0 to 1)
trev	Time reversal flag (yes = 1, no = 0)
srev	Frequency sweep reversal flag (0 to 1)
stretch	Stretching factor (>= 0)
dcflag	dc correction, y/n

Additional parameters are usually data matrices, such as Fourier coefficients or square wave parameters e.g. length, phase, amplitude, etc. These matrices are listed without parameter names. The size of the data matrix given is defined by:

cols	Number of columns
rows	Number of rows

Pbox incorporates the following amplitude modulation (AM) functions:

sq	Square (constant amplitude)
sqa	Square wave amplitude modulation (used for "composite" pulses)
gs	Gaussian
lz	Lorentzian
sch	sech (hyperbolic secant)
hta	tanh (hyperbolic tangent)
tra	Triangular amplitude (ramp)
SC	Sinc function
csp	Cosine power
wr	WURST (wideband uniform rate smooth truncation)
sed	Seduce-1, mixture of sech and sin
db	Quadrupolar
ata	Amplitude mod for CA atan frequency sweep pulse
exa	Exponential amplitude
tna	Tangential amplitude

- fs Fourier Series
- ft Inverse Fourier Transform

Pbox incorporates the following frequency modulation (FM) functions:

ls	Linear sweep (chirp)
tns	Tangential sweep (tan)
ht	Hyperbolic tangent sweep (tanh)
lzs	Constant adiabaticity Lorentzian sweep
ca	Constant adiabaticity (CA) sweep (frequency modulated frame)
cas	Constant adiabaticity sweep (phase modulated frame)
CS	Cosine / sine pulse frequency sweep
cs2	CA cosine square frequency sweep
CCS	CA cosine frequency sweep
sqw	Squarewave phase modulation
fsw	Frequency switch (step function)
fslg	Frequency switched as per Lee-Goldburg

Pbox Equation Evaluator

Pbox equation evaluator enables you to write your own expressions for amplitude, frequency, or phase modulation functions. The syntax for such equations is based on the "C" programming language. An example of a user-defined shape is provided in the wavelib/usr directory.

Reserved Shape Names

There are several shape names that are reserved for internal use and do not appear in wavelib.

NULL	A null shape is used to generate shapes with zero amplitude. Such shapes are useful for waveform generator tests and can be convenient in conjunction with creating complex waveforms that include delays.
RDC, rdc	Radiation damping compensation shape is a fictitious shape that is used to compensate other shapes for radiation damping. It has three active parameters:
	pw/bw entry is used to define the radiation damping rate (in Hz) or time constant (in seconds)
	ofs entry is used to define the frequency of the radiation damping source, e.g. water signal
	st entry is used to define the position of the radiation damping source (0 for Mz, 1 for Mxy, etc.).
	If the RDC shape is defined, it is active during the whole waveform.
a-N	Internal phase cycle that alternates the phase of a waveform between $-N$ and N, where N denotes the phase value. For example, a-45 is a two-step phase cycle with phase alternating between +45° and -45°. It can be used in conjunction with other phase cycles, for instance in nested phase cycles.

A-N	Internal phase cycle that alternates the phase of a waveform between 0 and N, where N denotes the phase value. For example A-180 is a two-step phase cycle
	with phase alternating between 0° and 180°. It can be used in conjunction with other phase cycles, for instance, in nested phase cycles.
R-N,r-N	Internal repetitive phase cycle that simply repeats the waveform N times. For example, r-4 repeats the same waveform four times. It can be used in conjunction with other phase cycles, for instance, in nested phase cycles.

Pbox VNMR Parameters

The following list describes Pbox VNMR parameters.

name	Name and extension of the output shape file. If the extension is not given, the shape type is set according to the type parameter. The default name is internally set as Pbox. This can be changed in the .Pbox_globals file.
type	 Shape type. Allowed values are r (.RF type), d (.DEC) or g (.GRD). If the shape type is not defined and the shape file is given without an extension, the shape file type is determined from the wave file according to the following criteria: type is set to r if pwbw > 0.0. type is set to d if dres > 0.0. type is set to g otherwise.
dres	Corresponds to dres parameter in VMNR. Active only with .DEC files.
steps	 Defines the required number of steps in the waveform. The default number of steps is stored with each individual shape in the corresponding wave file. This number can be overridden by Pbox if it is smaller than the internally calculated minimum number of steps, which is necessary to maintain the functionality of the waveform. This number is defined according to the following criteria: By the minimum number of steps necessary for adequate representation of the waveform (as in wave file). If the waveform is shifted off-resonance, by the Nyquist condition (see maxincr).
maxincr	Maximum phase increment. By default, set to 30°. This number is active only if the waveform is shifted off-resonance or the shape itself is frequency modulated (e.g., adiabatic sweeps). In order to satisfy the Nyquist condition, maxincr should not exceed 180°, otherwise the waveform gets folded back. In fact, the degradation of performance and interference with sidebands can be observed even with a maxincr of greater than 90°. Therefore, a maxincr of less than 90° is recommended

attn	Fine attenuation mode, which uses the following allowed values:			
	i	(Internal), default. Fine attenuation is implemented by internally		
		rescaling the waveform within the amplitude range set by maxamp		
		(0 to 1023).		
	е	(External) Fine attenuation is implemented by externally rescaling the		
		waveform using linear modulators. The internal amplitude is set to maxamp (1023.0) and the required fine attenuator setting is produced in		
		the output.		
	d	Attenuate to the nearest dB step by changing the pulse width, which will		
		affect the excitation bandwidth typically within 5%, which is tolerable		
		in most applications. The internal amplitude is set to maxamp (1023.0)		
	4.5i	Internally attenuate to a given (4.5 kHz) B1 field strength by adjusting the pulse length.		
	4.5e	Externally attenuate to a given (4.5 kHz) B1 field strength by adjusting the pulse length.		
	45I	Internally attenuate, keeping course power level at a given (45 dB) power level.		
	45E	Externally attenuate (with fine power), keeping course power level at a given (45 dB) power level.		
	45d	Attenuate to a given (45 dB) power level by changing the pulse width. The internal amplitude is set to maxamp (1023.0).		
sfrq	Spectro	ometer frequency in MHz.		
refofs	Referen	nce offset, usually 0.0. Can be specified if the excitation bands are shifted		
	by or re	eferenced to some frequency. Units: Hz, kHz, or ppm (if sfrq is defined).		
sucyc	Super of	cycle. Allowed values are n (no), d (default) or any name of a super cycle		
		in the wavelib/supercycles directory. By default, it is internally set		
		uper cycles can be nested by separating the names with a comma, for		
	exampl super c	e, t5,m4 represents 5 step TPG super cycle nested in four step MLEV-4 ycle.		
reps		s level of reporting. Allowed values are 0-4: 0=silent, 1=single line, mum, 3=medium, 4=maximum. The default is 2.		
stepsize		agth of a single step in a waveform. The default units are μ s. Note that size disables the maxingr parameter.		
bscor		s correction for Bloch-Siegert effect in multiple band excitation, inversion		
		cusing pulses. Allowed values are y (yes) or n (no, default). Active only if		
		nber of bands is two or more. Reduces the rf interference effects (see M.		
		, L.M.K. Vanderseypen and I.L.Chuang, Abstracts of the 41st ENC, p. 268,		
	Asilom	ar 2000).		
wrap		ound parameter. It allows wrapping around the waveform. The allowed are between 0 and 1.0.		
header		ile header. Allowed values are y (yes, default), n (no shape file header) and		
		ging). Information required for imaging systems is stored in the shape file		
	header.			
bsim		simulator. Performs Bloch simulation for the given waveform at the moment		
		eform generation. Allowed values are y (yes), n (no, default), a (add to the		
		is simulation), \mathfrak{s} (subtract from the previous simulation) and any positive limiting the simulation time in seconds. The default maximum length of		
	-	tion is internally set to 60 seconds and can be redefined in		
		box_globals file. Note, that Bloch simulator can also be externally		
		ed, e.g., from menus or using the dprofile macro.		
T1	-	udinal relaxation time, T1 in seconds. Can be required by some waveforms LURP pulses). Optional for the Bloch simulation.		
Т2		ersal relaxation time T2, in seconds. Can be required by some waveforms		
	(e.g. SLURP pulses). Optional for the Bloch simulation.			

I

dcyc	Duty cycle. Usually required for homonuclear decoupling applications. Only values between 0.0 and 1.0 are active. Outside these boundaries dcyc is reset to 1.0 (default).
SW	Spectral width. Required for Bloch simulation. Also recommended for H-H homo-decoupling applications. If given, the step size of waveform is set equal to the dwell time $(1/sw)$. sw helps to make sure that excitation sidebands are kept outside the spectral window.
ref_pw90	Reference 90° pulse width (in μ s) at ref_pwr. Required for calibration of waveforms. If set to 0.0, the maximum B1 field intensity (in kHz) is reported instead of the power setting.
ref_pwr	Reference power level (in dB steps). See ref_pw90.
ptype	Pulse type. Only necessary with imaging header. By default, set to selective.
pad1	Zero padding at the beginning of shapefile. Adds a specified number of steps of zero amplitude to the beginning of a shape file (gate is closed).
pad1	Zero padding at the end of shapefile. Adds a specified number of steps of zero amplitude to the end of a shape file (gate is closed).

Wave String Variables

A reminder is given in Pbox.inp files generated by menus and macros because these parameters appear without names. The wave string variables are listed as they appear in the reminder.

sh	Shape name as in wavelib.
bw/pw	Bandwidth and/or pulsewidth. For most waveforms, only one of the two parameters is required. Pbox distinguishes between bw (in Hz), which is always greater than 1.0, and pw (in sec), which is always less than 1.0. It is up to you which of the two parameters to provide for input, because they are mutually related via the pw*bw product, which is stored with each individual shape in wavelib. Some waveforms (e.g., adiabatic sweep pulses) can require both bw and pw. In such cases, both variables can be provided in a single string using the "/" separator. For example, {WURST2 200.0/0.05} denotes a 50-ms long WURST-2 pulse covering 200 -Hz-wide band. Alternatively, units can be used for clarity, e.g., {WURST2 0.2k/50m}. If the sfrq parameter is defined, bandwidth can also be specified in ppm, e.g., {WURST 20p/5m}.
ofs	Offset of the center of the excitation band in Hz with respect to the carrier frequency (middle of the spectrum). Note that if the sfrq spectrometer frequency, (in MHz) is defined, ofs can also be specified in ppm. In order to specify ofs in terms of absolute frequency, the reference offset refofs (i.e., chemical shift value of carrier frequency) must be defined. For instance, {WURST2 20p/5m 170p} sfrq=225.0 refofs=55p.
st	Spin status. Defines whether the waveform is used for excitation $(st=0)$, refocusing $(st=0.5)$ or de-excitation $(st=1)$, which, in turn, defines whether the wave starts with phase defined by ph $(st=1)$, the ph occurs in the middle of the pulse $(st=0.5)$, or the pulse ends with phase ph (status 0). In addition, the waveforms are time reversed if status is 1, as required for proper de-excitation. Undesired time reversal can be undone using the trev parameter. Furthermore, if several waves of different width are generated, they are bound to the beginning $(st=1)$, middle $(st=0.5)$, or end $(st=0)$ of the waveform. The spin status of the first wave is also used by Bloch simulator as the starting magnetization.

ph	Phase in degrees or phase cycle (super cycle). Usually phase is externally set in the pulse program and this parameter is not required. You can also apply any phasecycle (super cycle) from wavelib/supercycles. The difference between this phase cycle and the sucyc parameter is that phase cycling is carried out before waveform mixing and is therefore independent of other Super cycles, whereas sucyc is applied to the final (mixed) waveform. In this way, several waves of different width can be independently phase cycled and use different super cycles.
fla	Flip angle, in degrees. Usually, fla is defined in the wave file and there are very few applications where intermediate flip angles are required.
trev	Time reversal flag (see st). Allowed values are y (yes) and n (no, default).
d1	Prepulse delay, in seconds. Normally not required. If defined, it disables the internal wave shifting according to the spin status.
d2	Postpulse delay, in seconds. Normally not required. If defined, it disables the internal wave shifting according to the spin status.
d0	Pre-dl delay, in seconds. Essentially repeats dl. It is used only for convenience, e.g., if internal duty cycle is defined in shape parameters in wavelib. If set to 'a', the wave is appended to the previous wave.
wrap	Wraparound parameter. Can take values between 0 and 1.0.
php	Phase pivot. It determines the point in the shaped pulse when phase takes the value set by the ph (phase) parameter. The php parameter takes values between 0 (beginning of the pulse) and 1 (end of the pulse). By default, if $st=0$, then php is set to php = 1. If $st=1$, then php=0.

Creating Waveforms Using Menus

A set of menus is provided for convenience of pulse shaping from the VNMR window:

Pbox	Pbox180	Pbox180a	Pbox180b
Pbox180r	Pbox90	Pbox90a	Pbox90b
PboxDec	PboxDec2	PboxGrad	PboxHoDec
PboxHoDeca	PboxLib	PboxMix	PboxOpt
PboxOpt2	PboxSol	PboxWva	PboxWvb

To enter Pbox menus, do the following steps:

- 1. Enter **ds** in the VNMR input window.
- 2. Click on **Pbox**.

For example, to create a simple excitation pulse, such as an E-BURP-1 shaped pulse, do the following steps:

- 1. Click on **Pbox** > **90**.
- 2. Select an excitation region using cursors
- 3. Click on e-Burp1 > Close > Name > Close.

The Name button is optional and can be omitted. The shape files are stored under a default name (usually Pbox.RF).

To create a multiply-selective pulse, do the following steps:

- 1. Click on **Pbox** > **90**.
- 2. Select an excitation region using cursors.

- 3. Click on e-Snob.
- 4. Select the next region using cursors.
- 5. Click on e-Snob.
- 6. Repeat steps 4 and 5 as many times as required.
- 7. Click on Close > Close.

To create 4-ply selective excitation pulses for phase encoding (e.g., in experiments using the Hadamard transform), use the following sequence:

- 1. Click on **Pbox** > **90**.
- 2. Click on **Options** > **Phase** > **Return**.
- 3. Select an excitation region using cursors
- 4. Click on e-Burp1.
- 5. Repeat steps 1 through 4 as required.
- 6. Click on Close > Name > Close.

To create inversion and refocusing pulses:

- 1. Click **Pbox** > **180**.
- 2. Select an inversion band using cursors.
- 3. Click on **iBurp2** > **Close** > **Bloch**, enter **y**, and click on **Close**.

In the last example, a Bloch simulation is activated and the inversion profile is displayed on top of the spectrum. A relaxation-sensitive simulation on a refocusing pulse can be carried out by using the following menus:

- 1. Click on **Pbox** > **180** > **Refoc**.
- 2. Select a refocusing band using cursors.
- 3. Click on reBurp > Close > T1 > T2 > Bloch, enter y, and click on Close.

If the spectrum is not available, as is frequently the case with homo-decoupling in indirectly detected experiments, the bandwidth and offset can be manually entered, using the Options button. For example, off-resonance SEDUCE-1 decoupling centered on carbonyls can be created as follows:

- 1. Click on Pbox > Homo-dec > Options > Offset > Bandwidth > Return
- 2. Click on Seduce1 > Close > Name > Close.

For homodecoupling in the directly detected dimension, it is necessary to specify the decoupler duty cycle. For adiabatic decoupling, the J-coupling must also be specified:

- 1. Click on **Pbox** > **Homo-dec** > **Adiabatic** (set J = 10 Hz for H-H decoupling).
- 2. Select decoupling band using cursors.
- 3. Click on WURST-2 > Close > Dutycyc > Name > Close.

Heteronuclear decoupling and mixing sequences usually do not need to be shifted offresonance and the menu sequence can be very simple:

- 1. Click on **Pbox**.
- 2. Click on Mixing > bw (Hz) > DIPSI-3 > Name > Close.

A default supercycling is usually assigned to all waveforms in the directories wavelib/ decoupling and wavelib/mixing. However, the supercycle can be changed or assigned to other (nondecoupling) shapes either by using the Phase button in the Options menu or using the **Sucyc** button in the Close menu. It is necessary to click the **Sucyc** button several times in order to create a nested supercycle:

- 1. Click on Pbox > Het-dec > bw (Hz) > WALTZ16 > Name.
- 2. Click on Sucyc, enter t5, click on Sucyc, enter m4, and click on Close.

It is even simpler to create shaped gradients:

• Click on Pbox > Grad > h-Sine.

The Other button in different menus or the Wavelib button provide access to all wavelib items. The Wavelib sequence of menus is slightly different from the usual menus. For example, a multiply selective adiabatic inversion pulse using sech (hyperbolic secant) shape is created as follows:

- 1. Click on Pbox > Wavelib > (select inversion/) Ch Dir > (select sech).
- 2. Click on Set Shape.
- 3. Select inversion region using cursors.
- 4. Click on Set Wave.
- 5. Repeat step 4 as required.
- 6. Click on Close > Name > Close.

Information about various shapes is also accessed through the Wavelib menu:

- 1. Click on **Pbox** > **Wavelib**.
- 2. Select decoupling/ and click on Ch Dir.
- 3. Select MPF7 and click on i?
- 4. Repeat step 3 as required.

Pbox Macro Reference

Although most of needs for generating selective pulses can be satisfied by using Pbox menus, a set of macros is provided for those who prefer macros over menus. The following table lists the macros in the order of decreasing importance. For additional information on Pbox macros, refer to the manual *VNMR Command and Parameter Reference*.

opx	Opens Pbox, writes the Pbox.inp file header, and resets parameters r1-r7 and n1-n3.
selex	Defines the excitation band from the position of cursors in the graphics window and reports them to the user. It also sets r1 to excitation bandwidth and r2 to offset. selex uses the pbox_bw and putwave macros.
срх	Calls the Pbox command, which generates the specified waveform as defined by the Pbox.inp file.cpx also checks if parameters ref_pwr and ref_pw90 exist in the current experiment and puts their values into the Pbox.inp file. If the parameters do not exist, cpx creates them and asks the user for parameter magnitudes.
setwave	Sets up a single excitation band in the Pbox.inp file. An unlimited number of waves can be combined by reapplying setwave.

I

putwave	Sets up a single excitation band in the Pbox.inp file. An unlimited number of waves can be combined by reapplying putwave.
pxshape	Plots real (X) and imaginary (Y) components of a shaped pulse. Any type of waveform (.RF, .DEC or .GRD) can be plotted.
pboxpar	Adds a parameter definition to the Pbox.inp file.
pboxget	Extracts calibration data from the file shapefile.ext generated by Pbox or, if the file name is not provided, from the pbox.cal file containing parameters of the last created Pbox shape file. Note that the parameter is not changed by this macro if it was set to 'n' (not used)!
pbox_pw	Extracts pulse length from the file shapefile.RF generated by Pbox or, if the file name is not provided, from pbox.cal file containing parameters of the last created Pbox shape file.
pbox_pwr	Extracts the power lever from the file shapefile.ext generated by Pbox or, if the file name is not provided, from the pbox.cal file containing parameters of the last created Pbox shape file. Note that the parameter will not be changed by this macro if previously set to 'n' (not used).
pbox_pwrf	Extracts the fine power lever from the file shapefile.ext generated by Pbox or, if the file name is not provided, from the pbox.cal file containing parameters of the last created Pbox shape file. Note that the parameter will not be changed by this macro if it was previously set to 'n' (not used).
pbox_dmf	Extracts the dmf value from the file shapefile.DEC created by Pbox or, if the file name is not provided, from the pbox.cal file containing parameters of the last created Pbox shape file.
pbox_dres	Extracts the dres value from the file m shapefile.DEC created by Pbox or, if the file name is not provided, from the pbox.cal file containing parameters of the last created Pbox shape file.
pbox_name	Extracts name of the last shape file generated by Pbox and stored in the pbox.cal file. Note, that the file name extension is not stored explicitly and is not provided by this macro.
dshape	Displays real (X) and imaginary (Y) components of a shaped pulse. Any type of waveform (.RF, .DEC or .GRD) can be displayed.
pshape	Plots real (X) and imaginary (Y) components of a shaped pulse. Any type of waveform (.RF, .DEC or .GRD) can be plotted.
dshapef	Displays the real (X) and imaginary (Y) components of last generated shaped pulse, stored in pbox.fid file.
dshapei	Interactively displays the real (X) and imaginary (Y) components of last generated shaped pulse, stored in pbox.fid file
dprofile	Displays the X, Y, and Z excitation (inversion) profile for a pulse shape generated by the Pbox software.
pprofile	Plots the X, Y, and Z excitation (inversion) profile for a pulse shape that has been generated with the Pbox software. If a shape name is not provided, the last simulation data stored in shapelib/Pbox.sim are plotted.
pph	Prints out the shape file header (i.e., all lines starting with #).
pbox_bw	Defines the excitation band from the position of cursors in the graphics window and reports them to the user. It also sets r1 to excitation bandwidth and r2 to offset. This macro is used mainly in Pbox menus and macros.
pbox_bws	Defines the excitation band from the position of cursors in the graphics window and reports them to the user. It also sets $r1$ to excitation bandwidth and $r2$ to offset. Note, the left cursor should be placed on the left side of the excitation band and the right cursor on resonance of the solvent signal. This macro is mainly used in Pbox menus and macros.

pbox_rst	Resets r1=0, r2=0, r3=0, r4=0, n2= 'n', n3= '', and adds some
	standard comment lines to the Pbox.inp file. This macro is used in menus
	and other Pbox macros.
pbox_files	This macro is used only in conjunction with Pbox file menus.

Pbox PSG Statements

Pulse programs can be significantly simplified using a set of Pbox PSG statements. Because the calibration parameters are stored in the shape file header in a predefined form, they can be directly retrieved from the Pbox shape file within the pulse program. Storage reduces the number of statements in the pulse program and also the number of created and displayed parameters in the current experiment, especially in the case of complex multidimensional pulse sequences. Because parameters like dres, dmf, and selpw are normally not changed, there is no need to create and display them.

The Pbox psg statements are provided as Pbox_psg.hinclude file, which is normally stored in the /vnmr/psg directory. Compare the two DPFGSE pulse programs in Listing 1 and Listing 2.

Note that power, fine power, and pulse length of shaped pulses are internally adjusted by the pbox_pulse statement making pulse program code compact and easy to read. Fine power is set only if the external attenuation mode (attn=e in Pbox.inp) was used to generate the shape file. Although, you can still have access to all experimental parameters via the getval statements, in the majority of applications, these parameters do not need to be changed. For example, because of the high tolerance of adiabatic pulses to B1-inhomogeneity, possible incorrect calibration of this pulse is not a concern, and the only parameter requiring adjustment in the previous example is power level for homonuclear decoupling.

Pbox PSG statements can also be used to generate the shaped gradient:

pbox_grad(pattern, amp, gof1, gof2)

where gof1 and gof2 are (recovery) delays. The gradient pulse width is retrieved from the .GRD file shape header (defaulted to 1 ms by Pbox).

Shaped gradients for spectrometers with no PFG waveform generators is another feature available with the Pbox_psg.h include file, and the syntax is almost the same (except loops is substituted by number of steps, np):

shapedgradient(pattern, width, amp, channel, np, wait);

The shapedgradient macro in the /vnmr/psg/macros.h file needs to be redefined to implement shaped gradients on spectrometers with no gradient WFG as follows:

```
#define
shapedgradient(pulsefile,pulsewidth,gamp0,which,loops, \
    wait_4_me)shaped_gradient((int)(loops), \
    (double)(gamp0),(double)(pulsewidth),which)
```

Standard Pbox PSG Statements

shape

Declares shape parameters in VNMR pulse program

Description: Declares shape parameters in VNMR pulse program. The shape structure contains a set experimental parameters required to implement a shaped pulse or decoupling in a pulse sequence. The shape structure has the following members

```
Listing 1. Text of DPFGSE Pulse Program, Version A
```

```
/* dpfgse.c - version A */
#include <standard.h>
#include <Pbox_psq.h>
                                /* Pbox psg-functions included */
static int ph1[2] = \{0, 2\};
static int ph2[8] = \{0, 0, 2, 2, 1, 1, 3, 3\}
pulsesequence()
{
   double gzlvl1 = getval("gzlvl1"),
           gzlvl2 = getval("gzlvl2");
/* Declare & retrieve shape parameters */
   shape grdpat = getGshape("grdpat"), /* .GRD shape */
                                          /* .RF shape */
          pwpat = getRshape("pwpat"),
          decseq = getDshape("decseq"); /* .DEC shape */
   char Hdecflg[MAXSTR];
  getstr('Hdecflg', Hdecflg);
          decseq.pwr = getval('decpwr'); /* optional */
  settable(t1, 2, ph1);
  settable(t2, 8, ph2);
/* Pulse sequence starts here */
   status(A);
delay(d1);
   status(B);
obspower(tpwr);
rgpulse(pw, t2, rof1, rof2);
                                      /* hard ninety */
pbox_grad(grdpat, gzlvl1, 0.0, d2);
                                     /* Pbox shaped gradient */
pbox_pulse(pwpat, t1, rof1, rof2); /* Pbox selective 180 pulse */
pbox_grad(grdpat, gzlvl1, 0.0, d2);
                                     /* Pbox shaped gradient */
pbox_grad(grdpat, gzlvl2, 0.0, d2);
pbox_pulse(pwpat, t1, rof1, rof2);
pbox_grad(grdpat, gzlvl2, 0.0, d2);
   status(C);
setreceiver(t2);
   if (Hdecflg[0] == 'y') /* Pbox homo-decoupling */
homodec(decseq);
}
```

- name, pw, pwr, pwrf, dmf, and dres. In .GRD and .RF type waveforms, the dmf and dres parameters are not active and their values internally are set to zero.

Examples: shape grdpat, pwpat, decpat;

getRshape Retrieves shape parameters from .RF shape file header

Syntax: getRshape("parname"); escription: Retrieves shape name, pulse length, power and fi

Description: Retrieves shape name, pulse length, power and fine power from a .RF type shape file header provided that parameter parname exists in the current experiment. Note that shape parameters can be altered if necessary.

Listing 2. Text of DPFGSE Pulse Program, Version B

```
/* dpfgse.c - version B */
#include <standard.h>
static int ph1[2] = {0,2};
static int ph2[8] = \{0, 0, 1, 1, 2, 2, 3, 3\};
pulsesequence()
{
   double
            selpwr = getval("selpwr"),
            selpw = getval("selpw"),
            gt1 = getval("gt1"),
            gzlvl1 = getval("gzlvl1"),
            gt2 = getval("gt2"),
   gzlvl2 = getval("gzlvl2"),
            gof1 = 0.05*d2,
            gof2 = 0.95*d2,
            dutyc = getval("dutyc"),
            decpwr = getval("decpwr"),
            decres = getval("decres"),
            decdmf = getval("decdmf"),
            dpls90 = 1.0/decdmf,
            decdel = dutyc/sw,
            acqdel = (1.0-dutyc)/sw;
            decflg[MAXSTR], decseq[MAXSTR];
   char
   getstr("pwpat", pwpat);
   getstr("decflg", decflg);
   getstr("decseq", decseq);
   getstr("satmode", decseq);
   satpwr = getval("satpwr");
   satdly = getval("satdly");
   if ((dutyc > 0.2) && (decpwr > 49.0))
   {
      printf("decpwr too high! Aborting...\n");
      abort(1);
   }
   if( satpwr > 10.0)
   {
      printf("satpwr too high! Aborting...\n");
      abort(1);
   }
   settable(t1, 8, rec);
   settable(t2, 2, phb);
      status(A);
   if (satmode[0] == 'y')
   {
   obspower(satpwr);
   rgpulse(satdly, zero, rof1, rof2);
   delay(d1-satdly);
```

Listing 2. Text of DPFGSE Pulse Program, Version B (continued)

```
}
   else
   {
      delay(d1);
   }
      status(B);
  obspower(tpwr);
  rgpulse(pw, t1, rof1, rof2); /* hard 90 */
  obspower(selpwr);
  delay(gof1);
  rgradient('z',gzlvl1);
  delay(gt1);
  rgradient('z',0.0);
  delay(gof2);
  shaped_pulse(pwpat, selpw, t2, rof1, rof2);
  delay(gof1);
  rgradient('z',gzlvl1);
  delay(gt1);
  rgradient('z',0.0);
  delay(gof2);
  delay(POWER_DELAY);
  delay(gof1);
  rgradient('z',gzlvl2);
  delay(gt2);
  rgradient('z',0.0);
  delay(gof2);
  shaped_pulse(pwpat, selpw, t2, rof1, rof2);
  delay(gof1);
  rgradient('z',gzlvl2);
  delay(gt2);
  rgradient('z',0.0);
  delay(gof2);
      status(C);
  setreceiver(t1);
  if (decflg[0] == 'y')/* homo-decoupling */
   {
      delay(alfa + 1/(1.3*fb) - PRG_START_DELAY - PRG_STOP_DELAY);
      txphase(zero);
      obsprgon(decseq, dpls90, decres);
      initval(np/2.0, v14);
      starthardloop(v14);
                                   /* explicit acquisition */
         acquire(2.0, acqdel);
         rcvroff(); xmtron();
         delay(decdel); xmtroff(); rcvron();
      endhardloop();
      obsprgoff();
  }
}
```

	Examples:	<pre>(shape) pwpat=getRshape("pwpat"); pwpat.pwr=getval("selpwr");</pre>
	getDshape	Retrieves shape parameters from .DEC shape file header)
	Syntax:	getDshape("parname");
	Description:	Retrieves shape name, pulse length, dres, dmf, power and fine power from a .DEC type shape file header provided that parameter parname exists in the current experiment. Note that shape parameters can be altered if necessary.
	Examples:	<pre>(shape) decpat=getDshape("decpat"); decpat.dmf=getval(dmf2);</pre>
	getGshape	Retrieves shape parameters from .GRD shape file header
	Syntax:	getRshape("parname");
	Description:	Retrieves shape name and pulse length from a .GRD type shape file header provided that parameter parname exists in the current experiment. Note that shape parameters can be altered if necessary.
	Examples:	(shape) grdpat=getGshape("grdpat"); grdpat.pwr=getval("gt1"); /* optional */
	getRsh	Retrieves shape parameters directly from .RF shape file header
I		getRsh("shname");
I	Description:	Retrieves pulse length, power and fine power directly from the shname.RF shape file header.
	Examples:	<pre>(shape) fbshape=getRsh("H2Osinc"); (shape) fbshape=getRsh(pwpat);</pre>
	getDsh	Retrieves shape parameters directly from .DEC shape file header
I		getDsh("shname"); getDsh(dseq);
•	Description:	Retrieves pulse length, dres, dmf, power and fine power directly from shname.DEC shape file header.
	Examples:	<pre>(shape) decpat=getDsh("COdec"); (shape) decseq=getDsh(dseq);</pre>
	getGsh	Retrieves shape parameters from .GRD shape file header
I	Syntax:	getGsh("shname");
	Description:	Retrieves shape pulse length directly from the shapename.GRD shape file header.
	Examples:	(shape) grdpat=getGsh("waterg");
	pbox_pulse	Sets pulse width and fine power, performs a shaped pulse
	Syntax:	<pre>pbox_pulse(shpat,v1,rof1,rof2); shape shpat; codeint v1; /* phase */ double rof1, rof2;</pre>

Description: Sets pulse length, power and fine power (if used) and performs a shaped pulse on the transmitter channel. Note, that power is not automatically reset to the previous magnitude.

Examples: pbox_pulse(H2Oshape,v9,rof1,rof1); pbox_pulse(getRshape("pwpat"), oph, rof1, rof2); pbox pulse(getRsh(pwpat),zero,rof1,rof1);

pbox_decpulse Sets power, performs a shaped pulse on first decoupler

Syntax: pbox_decpulse(decsh,v1,rof1,rof2); shape decsh; codeint v1; /* phase */ double rof1, rof2;

- Description: Sets power, fine power (if used), and pulse duration and performs a shaped pulse on the first decoupler channel. Note, that power is not automatically reset to the previous magnitude.
 - Examples: pbox_decpulse(decsh,v5,rof1,rof1);

pbox_dec2pulse Sets power, performs a shaped pulse on second decoupler

Syntax:	<pre>pbox_dec2pulse(dec2pat,v1,rof1,rof2);</pre>				
	shape dec2pat;				
	codeint v1; /* phase */				
	double rof1, rof2;				
Description:	Sets power, fine power (if used), and pulse length. Performs a shaped pulse on the second decoupler channel. Note, that power is not automatically reset to the previous magnitude.				
Examples:	<pre>pbox_dec2pulse(dec2sh,zero,rof1,rof2);</pre>				

pbox_dec3pulse Sets power, performs a shaped pulse on third decoupler

Syntax:	<pre>pbox_dec3pulse(dec3pat,vi,rof1,rof2);</pre>					
	shape	dec3pa	at;			
	codeint	vi;	/*p	hase*/		
	double	rof1,	rof2;			

Description: Sets power, fine power (if used), and pulse length. Performs a shaped pulse on the third decoupler channel. Note, that power is not automatically reset to the previous magnitude.

Examples: pbox_dec3pulse(dec3sh,v2,rof1,rof2);

pbox_simpulse Sets power, performs simultaneous shaped pulse

Applicability:	Not applicable on <i>MERCURY</i> .
Syntax:	<pre>pbox_simpulse(shpat,decpat,vi,vj,rof1,rof2)</pre>
	shape shpat, decpat;
	codeint vi, vj; /* phases */
	double rof1, rof2;
Description:	Sets power and performs a simultaneous shaped pulse on the transmi

itter and the first decoupler channels. The pulse lengths and power levels are internally set within the Pbox statement. Note, that both transmitter and decoupler power levels are not automatically reset to the previous magnitude.

Examples: pbox_simpulse(shpat, decpat, zero, v4, rof1, rof2);

pbox_sim3pulseSets power, performs a simultaneous shaped pulse

- Description: Sets power and performs a simultaneous shaped pulse on the transmitter, the first and the second decoupler channels. The pulse lengths and power levels are internally set, within the Pbox statements. Note, that both transmitter and decoupler power levels are not automatically reset to the previous magnitudes.
 - Examples: pbox_sim3pulse(pwpat,decpat,dec2pat,zero,v5,three, \
 rof1,rof2);

pbox_sim4pulseSets power, performs a simultaneous shaped pulse

Syntax:	<pre>pbox_sim4pulse(pwpat,decpat,dec2pat,dec3pat, \ vi,vj,vk,vn,rof1,rof2); shape pwpat,decpat,dec2pat,dec3pat codeint vi,vj,vk,vn; /* phases */ double rof1,rof2;</pre>		
Description:	Sets power and performs a simultaneous shaped pulse on four RF channels. The pulse lengths and power levels are internally set, within the Pbox statements. Note, that both transmitter and decoupler power levels are not automatically reset to the previous magnitudes.		
Examples:	<pre>pbox_sim4pulse(shpat,decpat,dec2pat,dec3pat, \zero,two,v2,v9,rof1,rof2);</pre>		
pbox_xmtron	Sets power, dmf, and dres, initiates a programmable rf irradiation		
Applicability:	Not applicable on <i>MERCURY</i> .		
Syntax:	<pre>pbox_xmtron(mixpat); shape mixpat;</pre>		
Description:	Sets power, dmf, and dres and initiates a programmable rf irradiation via the observe transmitter channel. The gating of the observe transmitter (xmtron) is automatically executed within the pbox_xmtron statement.		
Examples:	<pre>pbox_xmtron(mixpat); pbox_xmtron(getDshape("decpat"));</pre>		
pbox_xmtroff	Terminates a programmable rf irradiation		
Applicability:	Not applicable on <i>MERCURY</i> .		
Syntax:	<pre>pbox_xmtroff();</pre>		
Description:	Terminates a programmable rf irradiation via the observe transmitter channel. Note that the power level is not automatically reset to its previous magnitude.		
pbox_decon	Sets power, \mathtt{dmf} and \mathtt{dres} , initiates a programmable rf irradiation		

Applicability: Not applicable on MERCURY.

I

I

	Syntax:	pbox_decon(decsh); shape decsh
	Description:	Sets power, dmf, and dres and initiates a programmable rf irradiation via the first decoupler channel. The gating of the first decoupler [decon()] is internally executed within the pbox_decon statement.
Ι	Examples:	<pre>pbox_decon(decsh); pbox_decon(getDshape("dseq")); pbox_decon(getDsh(dseq));</pre>
	pbox_decoff	Terminates a programmable rf irradiation
	Applicability:	Not applicable on <i>MERCURY</i> .
	Syntax:	<pre>pbox_decoff();</pre>
	Description:	Terminates a programmable rf irradiation via the first decoupler channel. Note that the power level is not automatically reset to the previous value.
	pbox_dec2on	Sets power, dmf and dres, initiates programmable rf irradiation
	Syntax:	pbox_dec2on(dec2sh); shape dec2sh;
	Description:	Sets power, dmf and dres and initiates a programmable rf irradiation via the second decoupler channel. Gating of the second decoupler [dec2on()] is internally executed within the pbox_dec2on statement.
	Examples:	<pre>pbox_dec2on(dec2sh); pbox_dec2on(getDshape("dseq2")); pbox_dec2on(getDsh(dseq2));</pre>
	pbox_dec2off	Terminates a programmable rf irradiation
	Syntax:	<pre>pbox_dec2off();</pre>
	Description:	Terminates a programmable rf irradiation via the second decoupler channel. Note that the power level is not automatically reset to the previous value.
	pbox_dec3on	Sets power, dmf, and dres, initiates a programmable rf irradiation
	Syntax:	<pre>pbox_dec3on(dec3sh); shape dec3sh;</pre>
	Description:	Sets power, dmf, and dres, and initiates a programmable rf irradiation via the third decoupler channel. Gating of the third decoupler [dec3on()] is internally executed within the pbox_dec3on statement.
I	Examples:	<pre>pbox_dec3on(dec3sh); pbox_dec3on(getDshape("dseq3")); pbox_dec3on(getDsh("WALTZ16"));</pre>
	pbox_dec3off	Terminates a programmable rf irradiation
	Syntax:	<pre>pbox_dec3off();</pre>
	Description:	Terminates a programmable rf irradiation via the third decoupler channel. Note that the power level is not automatically reset to the previous value.

pbox_spinlock Sets power and phase, executes a programmable spin lock

Applicability: Not applicable on *MERCURY*.

```
Syntax: pbox_spinlock(mixshape,mixtime,mixphase);
    shape mixshape;
    double mixtime;
    codeint mixphase;
```

Description: Sets power, fine power, dmf, dres, calculates the required number of loops and executes a programmable spin lock via the observe transmitter channel.

```
Examples: pbox_spinlock (mixsh,mix,zero);
    pbox_spinlock(getDshape("mixpat"),mix,zero);
    pbox_spinlock(getDsh(mixpat),mix,v2);
```

pbox_decspinlock Sets power and phase, executes a programmable spin lock

Applicability: Not applicable on MERCURY.

Syntax:	<pre>pbox_decspinlock(mixshape,mixtime,mixphase);</pre>		
	shape	mixshape;	
	double	mixtime;	
	codeint	mixphase;	

- Description: Sets power, fine power, dmf, dres, calculates the required number of loops and executes a programmable spin lock via the first decoupler channel.
- Examples: pbox_decspinlock(mixsh,mix,zero); pbox_decspinlock(getDshape("mixpat"),mix,zero); pbox_decspinlock(getDsh(dipsi2),mix,v2);

pbox_dec2spinlock Sets power and phase, executes a programmable spin lock

Syntax:	<pre>pbox_dec2spinlock(mixshape,mixtime,mixphase);</pre>		
	shape	mixshape;	
	double	mixtime;	
	codeint	mixphase;	
Description:	Sets power, fine power, dmf, dres, calculates the required number of loops and executes a programmable spin lock via the second decoupler channel.		
Examples:	<pre>pbox_dec2spinlock(mixsh,mix,zero); pbox_dec2spinlock(getDshape("mixpat"),mix,zero); pbox_dec2spinlock(getDsh(mixpat),mix,v2);</pre>		

pbox_dec3spinlock Sets power and phase, executes a programmable spin lock

	C		
	codeint	mixphase;	
	double	mixtime;	
	shape	mixshape;	
Syntax:	pbox_dec3spinlock(mixshape,mixtime,mixphase		

- Description: Sets power, fine power, dmf, dres, calculates the required number of loops and executes a programmable spin lock via the third decoupler channel.
- Examples: pbox_dec3spinlock(mixsh,mix,zero);
 pbox_dec3spinlock(getDshape("mixpat"),mix,zero);
 pbox_dec3spinlock(getDsh(mixpat),mix,v2);

pbox_grad	Sets gradient length and performs shaped z gradient pulse
Applicability:	Systems with no WFG on PFG module (see page 243). On <i>MERCURY</i> , installation is not required.
Syntax:	<pre>pbox_grad(grdpat,gzlvl,gof1,gof2); shape grdpat; double gzlvl, gof1, gof2;</pre>
Description:	Sets the gradient length and performs a shaped 'z' gradient pulse. Unlike the standard shapedgradient statement, the gradient is turned off automatically. gof1 and gof2 are pregradient and postgradient (recovery) delays.
pbox_xgrad	Sets the gradient length and performs shaped x gradient pulse
Applicability:	Systems with no WFG on PFG module (see page 243). Not applicable on <i>MERCURY</i> .
Syntax:	<pre>pbox_xgrad(grdpat,gw,gxlvl,gof1,gof2); shape grdpat double gxlvl, gof1, gof2 gw;</pre>
Description:	Performs a shaped 'x' gradient pulse of length gw. Unlike the standard shapedgradient statement, the gradient is automatically turned off.gof1 and gof2 are pregradient and postgradient (recovery) delays.
pbox_ygrad	Sets the gradient length and performs shaped y gradient pulse
Applicability:	Systems with no WFG on PFG module (see page 243). Not applicable on <i>MERCURY</i> .
Syntax:	<pre>pbox_ygrad(grdpat,gw,gylvl,gof1,gof2); shape grdpat; double gylvl,gof1,gof2,gw;</pre>
Description:	Performs a shaped 'y' gradient pulse of length gw. Unlike the standard shapedgradient statement, the gradient is automatically turned off. gof1 and gof2 are pregradient and postgradient (recovery) delays.
pbox_zgrad	Sets the grad length and performs a shaped z grad pulse
Applicability:	Systems with no WFG on PFG module (see page 243). On <i>MERCURY</i> , installation is not required.
Syntax:	<pre>pbox_zgrad(grdpat,gw,gzlvl,gof1,gof2); shape grdpat double gzlvl,gof1,gof2</pre>
Description:	Performs a shaped 'z' gradient pulse of length gw. Unlike the standard shapedgradient statement, the gradient is automatically turned off.gof1 and gof2 are pregradient and postgradient (recovery) delays.
Miscella	aneous PSG Statements

homodec Set power, dutyc, dmf, and dres, then executes irradiation

Applicability: Not applicable on MERCURY.

Syntax:	homodec(decseq) shape decseq
Description:	Sets power, dmf, and dres, and executes a programmable rf irradiation via the observe transmitter during the acquisition in a time-shared mode.
Examples:	<pre>homodec(Hdecsh); homodec(getDshape("Hdecpat")); homodec(getDsh("wurst20"));</pre>
pfg_pulse	Perform a rectangular z gradient pulse
Syntax:	pfg_pulse(gzlvl,gw,gof1,gof2); double gzlvl, gw, gof1, gof2
Description:	gw is gradient pulse width, gof1 and gof2 are pregradient and postgradient (recovery) delays.
Examples:	<pre>pfg_pulse(gz2lvl,gt2,rof1,grec);</pre>
presat	Performs a long rf pulse (presaturtion)
Syntax:	presat();
Description:	Performs a long rf pulse (presaturtion) via the observe transmitter. Uses the standard presat parameters satpwr, satdly, and dl. presat() is executed only if satdly is greater than 0.0 (satflg and satfrq are not active).
Examples:	<pre>presat();</pre>
pre_sat	Performs a long rf pulse (presaturtion)
Syntax:	pre_sat();
Description:	Performs a long rf pulse (presaturtion) via the observe transmitter. It uses the standard presat parameters satmode and satfrq in addition to satpwr, satdly, and dl.
Examples:	<pre>pre_sat();</pre>
setlimit	Sets a (safety) limit for the given parameter
Syntax:	<pre>setlimit(name, parameter, limit); char name[MAXSTR]; double parameter, limit;</pre>
Description:	Sets a (safety) limit for the given parameter. The execution of the pulse program
	is terminated if the parameter magnitude exceeds the given (safety) limit.

Pulse Shaping "On-Fly"

It is often convenient to create pulse shapes from within a pulse program (on-fly), i.e. during acquisition (go), which enables you to keep the experiment setup very simple and independent of spectrometer frequency. In "C," pulse shaping is implemented using the system() statement. For example:

system("Pbox sh.RF -w \"esnob 20p 170p\" -attn e -refofs 55p -sfrq 150.05 -ref_pwr 51 -ref_pw90 60.2");

I

A more convenient and flexible approach is to construct the command line first and then execute it using the system() statement. This method enables you to retrieve the necessary parameters from the experiment. The following example shows how to create and use a wide-band adiabatic inversion pulse with the maximum efficiency:

```
...
char cmd[MAXSTR];
double pwxlvl, pwx;
...
if((getval("arraydim") < 1.5) || (ix==1)) /* execute only once */
{
    sprintf(cmd, "Pbox ad180.RF -w \"cawurst-10 %.1f/%.6f\" -s
    1.0 -0",1.0/pwx, 20.0*pwx);
    system(cmd);
}
...
    decpower(pwxlvl);
    decshaped_pulse("ad180", 20.0*pwx, zero, rof1, rof1);
...</pre>
```

Pbox_psg.h include Pulse Sequence Statements

For even more convenience and readability, the following pulse sequence statements are provided in the Pbox_psg.h include file.

opx	Initiates Pbox within pulse program					
Syntax:	opx(shname); char shname[MAXSTR];					
Description:	Initiates Pbox from within a pulse program. Similar to the opx macro, except the Pbox.inp file is not created.					
Examples:	<pre>opx("water");</pre>					
setwave	Writes wave string into Pbox buffer					
Syntax:	setwave(wvstr); char wvstr[MAXSTR];					
Description:	Writes a wave definition string into a Pbox buffer. Used for pulse shaping from within a pulse program (on fly). Similar to the setwave macro, except the Pbox.inp file is not created.					
Examples:	<pre>setwave("eburp1 200.0 -1.2k"); setwave("esnob 20p 170p");</pre>					
putwave	Writes wave string into Pbox buffer					
Syntax:	<pre>putwave(sh,bw,ofs,st,pha,fla); char sh[MAXSTR]; double bw,ofs,st,pha,fla;</pre>					
Description:	Writes a wave definition string into a Pbox buffer. Used for pulse shaping from within a pulse program (on fly). Similar to the putwave macro, except the Pbox.inp file is not created. Unlike the macro, a full set of arguments is required. Zero can be used to request a default value.					
Examples:	putwave("eburp1",200.0,-1200.0,1.0,90.0,0.0);					

	pboxpar	Writes a Pbox parameter into Pbox buffer
	Syntax:	<pre>pboxpar(parname, parval); char parname[MAXSTR]; double parval;</pre>
	Description:	Writes a Pbox parameter into a Pbox buffer. Used for pulse shaping from within a pulse program (on fly). Note that parval is a real number.
	Examples:	<pre>pboxpar("stepsize",2.0); pboxpar("sfrq",dfrq);</pre>
I	pboxSpar	Writes a Pbox string parameter into Pbox buffer
	Syntax:	<pre>pboxpar(parname,parval); char parname[MAXSTR],parval[MAXSTR];</pre>
	Description:	Writes a character type Pbox parameter into a Pbox buffer. Used for pulse shaping from within a pulse program (on fly). Note that parval is a string variable.
	Examples:	<pre>pboxpar("attn","e"); pboxpar("sucyc","t5,m4"); pboxpar("refofs","55p");</pre>
	pboxUpar	Writes a Pbox parameter with units into Pbox buffer
	Syntax:	<pre>pboxpar(parname,parval,units); char parname[MAXSTR],units[MAXSTR]; double parval;</pre>
	Description:	Writes a numeric Pbox parameter with units into a Pbox buffer. Used for pulse shaping from within a pulse program (on fly). Note that units is a string variable.
	Examples:	<pre>pboUxpar("attn",45,"d"); pboxUpar("refofs",55,"p");</pre>
	срх	Executes Pbox within pulse program
	Syntax:	<pre>cpx(ref_pw90,ref_pwr); double ref_pwr,ref_pw90;</pre>
	Description:	Executes Pbox from within a pulse program. Similar to the cpx macro.
	Examples:	<pre>cpx(ref_pw90,ref_pwr); cpx(compH*pw90,tpwr);</pre>
	pbox_get	Retrieves shape parameters within pulse program
	Syntax:	<pre>pbox_get();</pre>
	Description:	Retrieves shape parameters from the pbox.cal file and sets pbox_pw, pbox_pwr, pbox_pwrf, pbox_dres and pbox_dmf from within a pulse program. These parameters can then be assigned to appropriate variables.
	Examples:	<pre>pbox_get(); selpw=pbox_pw; selpwr=pbox_pwr;</pre>

Chapter 7. Pulse Analysis

isarry	Checks whether given parameter is arrayed
Syntax:	isarry(parname); char parname[MAXSTR];
Description:	Returns 1 if parameter is arrayed and zero otherwise.
Examples:	<pre>isarry("selpwr");</pre>
FIRST_FID	Flag indicating the first FID in an arrayed experiment
Description:	Indicates the first experiment in an arrayed or multi-dimensional experiment. Useful if a particular pulse sequence statement needs be executed only once, at the beginning of an arrayed experiment.
Examples:	if(FIRST_FID) opx("Pbox");
pbox_shape	Generates a shape file and returns shape parameters
Syntax:	<pre>shape pbox_shape(shn, wvn, pw_bw, ofs, rf_pw90 \</pre>
	rf_pwr); char shn[MAXSTR], wvn[MAXSTR];
	double pw_bw, ofs, rf_pw90, rf_pwr;
Description:	Generates a shape file from a given set of input parameters.
Parameters:	shn is the shape file name.
	wvn is the wavelib waveform name.
	pw_bw is the pulse length in seconds or, alternatively, excitation bandwidth in Hz.
	ofs is the excitation offset from the carrier frequency, in Hz.
	rf_pw90 is the reference 90° pulse width, in ms.
	If necessary, other parameters, such as number of steps or attenuation level, must be predefined using the pboxpar() and related statements. This function returns a shape structure containing all the necessary experimental parameters. It is used for pulse shaping from within a pulse program (on fly).
Examples:	offC7=pbox_shape("offC7","sinc90",118*dfrq,0.0,pwC, \
	<pre>pwClvl); pwC7=offC7.pw; rf7=offC7.pwrf</pre>
pboxAshape	Generates a shape file and returns shape parameters
Syntax:	<pre>shape pboxAshape(shn, wvn, bw, pws, ofs, rf_pw90, \ rf_pwr); char shn[MAXSTR], wvn[MAXSTR];</pre>
	double bw, pws, ofs, rf_pw90, rf_pwr;
Description:	Generates a (adiabatic) shape file from a given set of input parameters.
Parameters:	shn is the shape file name.
	wvn is the wavelib waveform name.
	bw is the excitation bandwidth in Hz.
	pws is the pulse length in seconds.
	ofs is the excitation offset from the carrier frequency, in Hz.
	rf_pw90 is the reference 90° pulse width, in ms.

Note that both pulse duration and excitation bandwidth must be provided, as generally required for adiabatic waveforms. If necessary, other parameters, such as number of steps or attenuation level, must be predefined using pboxpar() and related statements. This function returns a shape structure containing all the necessary experimental parameters. It is used for pulse shaping from within a pulse program (on fly).

```
stC80 = pbox_shape("stC80","sech",80*dfrq,0.001,0.0, \
pwC,pwClvl);
pwC80=stC80.pw; rfst=stC80.pwrf;
```

shonfly.c Sequence

The following two examples show a simple shonfly.c sequence. You must check the arraydim and ix parameters to avoid excessive pulse shaping in arrayed experiments.

Version A

```
/* shonfly.c - version A, allows to array ofs */
#include <standard.h>
#include <Pbox_psg.h>
static double selpw=0.0, selpwr=0.0;
static char sharr[MAXSTR];
pulsesequence()
{
   double
            n=0.0,
                                     /* frequency offset */
             ofs = getval("ofs"),
             ref_pwr = getval("ref_pwr"),
             ref_pw90 = getval("ref_pw90");
   char
             repflg[MAXSTR];
   getstr("pwpat", pwpat);
                                 /* Retrieve the variables */
   getstr("repflg", repflg);
/* create a shaped pulse using Pbox macros only when needed */
 if((getval('arraydim') < 1.5) || (ix==1) || isarry("ofs"))
{
     sprintf(sharr, "%s_%ld", pwpat, ix);
     opx(sharr);
                                         /* open Pbox */
     putwave("gsneeze", 200.0, ofs, n, n, n);
     pboxpar("stepsize", 10.0);
                                        /* stepsize in us */
     if (repflg[A] == 'y') pbox_par("reps", "2");
                         pbox_par("reps", "0"); /* silent mode */
     else
                                                  /* close Pbox */
     cpx(ref_pw90, ref_pwr);
     pbox_get();
                                    /* retrieve shape data */
     selpwr = pbox_pwr;
                                        /* use Pbox power */
                                         /* use Pbox pw
     selpw = pbox_pw;
                                                           */
}
/* THE PULSE PROGRAM STARTS HERE */
status(A);
  obspower(selpwr);
  delay(d1);
status(B);
```

```
shaped_pulse(sharr, selpw, oph, rof1, rof2);
status(C);
}
```

Version B

```
/* shonfly.c - version B, allows to array ofs */
#include <standard.h>
#include <Pbox_psg.h>
static char sharr[MAXSTR];
pulsesequence()
{
double reps=0.0,
                                        /* silent mode */
        ofs = getval("ofs"),
                                        /* frequency offset */
        ref_pwr = getval("ref_pwr"),
       ref_pw90 = getval("ref_pw90");
char
       repflg[MAXSTR], cmd[MAXSTR];
                                 /* Retrieve the variables */
getstr("repflg", repflg);
if (repflg[A] == 'y') reps = 2.0;
/* create a shaped pulse using Pbox macros only when needed */
if((getval("arraydim") < 1.5) || (ix==1) || isarry("ofs"))
{
   sprintf(sharr, "shtst_%ld", ix);
   sprintf(cmd, "Pbox %s -w \"qsneeze 200.0 %.2f\" -s 10",
                 sharr, ofs);
   sprintf(cmd, "%s -p %.0f -l %.2f %.0f", cmd, ref_pwr,
                 ref_pw90*1.0e6, -reps);
   system(cmd);
}
/* THE PULSE PROGRAM STARTS HERE */
status(A);
 delay(d1);
status(B);
 pbox_pulse(getRsh(sharr), oph, rof1, rof2);
status(C);
}
```

Although version B is less readable, it is clearly more compact and efficient.

Pbox UNIX Commands

The Pbox program is always executed when a shaped pulse is created. Any of the Pbox parameters can be used as an argument followed by the parameter value. The arguments

and shortcuts listed in Table 43 are available. Note that the output filename is optional and is always the first argument.

Command	Parameter	Action
Pbox*	-b time	Activate Bloch simulator, opt=a (add), s (subtract), or time in sec.
	-C	Calibrate only, do not create a shape file.
	-f file	Set name of the output file.
	-h wave	Print wave file header.
	-i wave	Print wave file parameters.
	-l ref_pw90	Length (in µs) of reference pw90 pulse.
	-0	List options.
	-p ref_pwr	Reference power level (dB).
	-r file	Reshape Pbox pulse.
	-s stepsize	Define the length (in μ s) of a single step in the waveform.
	-t wave	Print shape title from wave file.
	-u userdir	Set user home directory.
	-w wavestr	Set wave definition string.
	-v	Run in verbose mode. Also print Pbox version.
	-x	Prints all Pbox parameters.
	-value	Sets reps to value.
Pxsim		Used in Pbox menus and macros for simulation of excitation profiles of shaped pulses.
Pxfid		Used by dshape and dshapei to format shape file into a FID-format text file.
Pxspy		Converts alien shapes (.RF, .DEC and .GRD) into Pbox compatible file format. Essentially converts a time-domain shape file into (frequency-domain) Fourier coefficients, which can be used to create a wave file in the wavelib directory.
Pbox sel Pbox -w	shape -wc 'eburp .RF -w 'eburpl 4 'eburpl 200 -120 .RF -w 'esnob 20	01 450 -1280.0' -1 220 -800' 'eburp1 420 1200' 00' -attn e -pl 45 54.2 -b p 170p' -sfrq 150.02 -refofs 55p -refpwr 45 ∖

 Table 43. Pbox Commands and Parameters (continued)

Chapter 7. Pulse Analysis

Chapter 8. Variable Temperature Operation

Sections in this chapter:

- 8.1 "Startup," this page
- 8.2 "Operating Procedures," page 262
- 8.3 "Temperature-Related Command," page 264
- 8.4 "Operating Recommendations," page 265
- 8.5 "VT Controller Safety Circuits," page 266
- 8.6 "VT Interlock Parameters," page 267

This chapter describes startup and operation of the optional variable temperature (VT) unit. A VT unit is available for all Varian NMR spectrometers to vary the temperature of the sample. Table 44 lists commands and parameters associated with VT operation.

Table 44.	Variable	Temperature	Unit O	peration	Commands	and Parameters
-----------	----------	-------------	--------	----------	----------	----------------

```
Command
acqmeter<('host')>
                                                     Open the Acqmeter window
dqs
                                                     Display group shims, automation parameters
ga<('nocheck')>
                                                     Submit experiment to acquisition and FT result
go<(<'acqi><,'nocheck'><,'nosafe'>)>
                                                     Submit experiment to acquisition
sethw('vt','reset')
                                                     Reset VT controller (not on GEMINI 2000)
su
                                                     Submit a setup experiment to acquisition
tempcal(solvent)<:temperature>
                                                     Temperature calculation
Parameter
                                                     Lock and spin interlock
in {'n','w','y'}
                                                     Preacquisition delay
pad*
temp \{-150 \text{ to } +2000, \text{ in deg. } C\}
                                                     Sample temperature
tin {'n', 'w', 'y'}
                                                     Temperature interlock
vtc {0 to 50, in deg. C}
                                                     Variable temperature cutoff point
vttype {2,0}
                                                     Variable temperature controller type
vtwait {number, in sec}
                                                     Variable temperature wait time
* pad {0 to 4095 (GEMINI 2000), 0 to 8190 (others), in seconds}
```

8.1 Startup

On systems equipped with the VT unit, the parameter temp, set by the user, changes the internal probe temperature. A thermocouple senses the temperature, which the VT controller continuously displays on the front panel. The controller compares the user-requested value with the current probe temperature and changes the heater current accordingly. The VT controller then reports the temperature of the gas flow and status to the spectrometer through a serial port at the rear of the console. The vtc parameter (for

"variable temperature cutoff"), also set by the user, determines the temperature below which the gas is cooled. The values of the temp and vtc parameters are shown in the "SPECIAL" section of the parameter display called by the dgs macro.

Starting up the VT unit takes the following steps:

- 1. Refer to the installation manual to make sure that the hardware is installed and connected properly, and that the VT controller is calibrated correctly.
- 2. If the system power has been off or the VT unit has been disconnected from the probe, reset the VT controller by pressing the POWER switch to turn the unit off, then pressing POWER again to turn it on. The VT controller also can be reset from the Temperature Control window (see page 264) or by entering the VNMR command sethw('vt','reset').

CAUTION: For VT and probe operation, use either dry nitrogen gas or air. A mixture of nitrogen gas and air can cause spikes in the baseline adjacent the large peaks in the spectra. For temperatures above 100°C, the use of air as the VT gas is not recommended. Such use will destructively oxidize the heater element and the thermocouple.

3. Use dry nitrogen gas if the requested temperature is over 100°C or below the dew point or 0° C, whichever is higher. Otherwise, air may be used as the VT gas. If the requested temperature is below -40°C, dry nitrogen gas is recommended for cooling the bearing, spinner, and decoupler. This prevents moisture condensation in the probe and spinner housing.

The source of heating or cooling gas is not automatically selected. To use nitrogen, you must attach a nitrogen gas source to the VT system. The same is true when using air. The VT system only selects the routing of the gas flow.

- 4. Use the flow control meter on the magnet leg to adjust the flow to about 10 LPM (as shown on the flow gauge).
- 5. A sample that can be handled at ambient temperature can now be placed in the probe, NMR lock obtained, and field homogeneity adjusted. Samples that cannot be handled at ambient temperature should wait until the system reaches the requested temperature.

8.2 Operating Procedures

The following procedures are recommended:

1. (Highland VT units) The power switch is located on the back panel. The heater on/ off switch is located on the front panel.

(Oxford VT units) Although the VT unit should be left on at all times, it's a good idea to check that the VT controller power is on (front panel button pressed in with the light on).

If the VT controller is off and you cannot turn it on, run the config program and check that the VT Controller label is set to Present.

- 2. Enter the acquisition parameters for the experiment as usual.
- 3. Enter **acqmeter** to open the Acqmeter display. Click in the display and select VT in the popup menu. To change the look of the display, select Properties from the popup menu.

- 4. Enter the desired temperature value(s) using the temp parameter.
- 5. If the temp parameter is an array, set the pad parameter for preacquisition delay that allows sufficient time for the sample to equilibrate after a temperature change. The system will then wait pad seconds in between each temperature before starting data acquisition. Delays of several minutes are optimum because the sample will take longer to equilibrate than it takes the VT controller to stabilize the heating/ cooling gas at the set point.
- 6. Set the vtc parameter to a temperature near the ambient VT gas temperature (normally vtc is correct and need not be changed). Based upon the value of vtc compared to the value of temp, the system will route the VT gas flow for either heating or cooling.
- 7. Start temperature control by entering a **su**, **go**, or **ga** command. These commands act as follows:
 - If su is entered, the temperature control and acquisition hardware controls are set and the sample temperature changed to the desired temperature. The experiment will not be started when the desired temperature is reached. After waiting pad seconds, go or ga must be issued before data acquisition will begin.
 - If a go or ga is entered, the same actions as su occur, except that after reaching the desired temperature, the system waits pad seconds then begins the pulse sequence and data acquisition. The pad time delay will occur every time the temperature is changed under program control.

CAUTION: Do not use aromatic, ketone (including acetone), and chlorinated solvents in the coolant bucket. Such coolant media attack the standard polystyrene bucket. Another type of container must be used (not supplied by Varian).

After entering su, go, or ga, the selection of the VT gas routing occurs, and the VT controller begins to control the gas temperature in the probe at the requested value of temp. The temperature readout will begin to change and the VT indicator light will begin flashing. At this time, if the requested temperature is below ambient, add coolant liquid to the coolant bucket.

CAUTION: Operating the system with the coolant bucket filled with liquid nitrogen and with the value of temp greater than the value of vtc results in the condensation of liquid nitrogen inside the exchanger coil tube. If the exchanger coil is then warmed above -210°C or if nitrogen gas is passed through the coil (when temp is less than vtc), very cold liquid nitrogen is forced through the transfer line and into the probe. This will cause a sudden pressure surge in the transfer lines and probe as the liquid nitrogen boils, and it can blow the flexible connector apart. If the liquid nitrogen reaches the glass components of the probe and sample tube, the glass will probably break. Instrument damage can be avoided by following these precautions:

Do not immerse the exchanger coil in liquid nitrogen when no nitrogen gas is flowing through the coil.

Do not stop the VT nitrogen gas flow while the exchanger is immersed in liquid nitrogen. Arrayed VT experiments that have a temperature range from above *vtc* to below *vtc* should be set up starting at the lowest temperature and ending at the highest temperature. When the experiment passes the *vtc* crossover, remove the liquid nitrogen coolant.

To avoid water in the exchanger when the low temperature experiment is complete, warm up the exchanger by removing it from the liquid nitrogen and maintain a flow of dry nitrogen until room temperature is reached.

When the temperature reaches the value requested by temp (it may initially overshoot), the VT indicator light stays on steadily. A sample that could not be handled at ambient temperature can now be transferred into the probe. The VT readout is the temperature of the cooling/heating gas and may be different from the true sample temperature. The exact temperature of the sample is correctly determined by a calibration curve that must be constructed for each probe, and must include flow rate and equilibration time. Refer to the VT installation manual for the NMR calibration method.

CAUTION: Before running sealed samples at elevated temperatures, check the samples in an oven at a temperature higher than the highest temperature during the experiment. When heated, volatile materials in a sealed tube can build up high pressures. If the tube ruptures while in the probe, the glass components and insert coil will probably be destroyed. For the same reason, do not insert into the probe any sealed sample of volatile material that must be kept cold to avoid excessive pressure buildup. Be sure the probe has equilibrated to a safe, cold temperature before inserting the sample.

8.3 Temperature-Related Command

For systems equipped with the optional variable temperature controller, the temp command opens the Temperature Control window, shown in Figure 59.

The Temperature Control window can be used for the following purposes

- Turn off temperature control.
- Set temperature control on at a specified temperature in degrees C.
- Enable temperature control from within an experiment using the temp parameter and the su, go, ga, or au commands and macros.
- Alternatively, turn off experiment control of the temperature and allow only the Temperature Control window (and the command sethw) to set the temperature.
- Reset the temperature controller when the temperature cable is reconnected to a probe.

If the temperature is controlled only through the Temperature Control window, two actions (to be taken after a temperature error) can be selected:

- Display a warning but continue acquisition.
- Stop acquisition and display a warning.

If experiment control of temperature is selected, the two previous selections appear faded because they are inoperative, and the selection of the action to be taken after a temperature error is provided by the parameter tin.




Figure 59. Temperature Control Window

8.4 Operating Recommendations

The following recommendations should help achieve better VT performance.

- The spectrometer system was designed and tested with a VT gas flow rate of about 10 LPM. Sizable deviation from this rate can result in significant inaccuracy in temperature calibration and reduce the attainable temperature limits.
- Initial cool-down of the exchanger and transfer tubing after the coolant is added increases the initial time required to reach regulation (about 5 to 10 minutes for -40°C with liquid nitrogen). Because this may be longer than the vtwait parameter, an su command is the best way to start up.
- Below -40°, using dry nitrogen gas for the spinner and bearing air supply avoids moisture and frost buildup on the spinner housing and turbine. Should this happen, the sample spinning can become erratic or stop altogether.
- Every sample has some vertical temperature gradient. Minimize the gradient by *not* filling the sample tube more than about 25 to 32 mm (1 to 1.25 in.), by inserting a vortex plug or glass wool plug in the tube just above the sample solution, and by entering the liquid column to the probe coil center lines. The plug reduces refluxing of the solvent in the upper portion of the tube. Any mass movement, such as refluxing or convection, can seriously degrade resolution and lock stability.
- Above 100°C, use dry nitrogen gas to reduce heater and thermocouple oxidation.
- High-power decoupling adds heat to the sample. The increase in temperature depends on the dielectric of the solution and the power level. Under these conditions, the temperature accuracy under VT control is significantly affected. If necessary, reduce decoupler power and use a more efficient decoupling mode, such as WALTZ-16 or GARP.
- Overnight or long-term unattended VT operation at low temperatures is hampered by the fact that the usual coolant, liquid nitrogen, provides only about 1 to 2 hours of operation on a single fill of the coolant bucket. Some other coolant that lasts longer can

be used if the operating temperature does not require the low temperature of liquid nitrogen. A common alternative is a mixture of dry ice and acetone. Another option is a fluid such as isopropyl alcohol or ethylene glycol cooled indirectly by a refrigerating device. Do not use aromatic, ketone, and chlorinated solvents (including acetone) in the coolant bucket. Such coolant media attack the standard polystyrene bucket.

- The ability of the VT unit to achieve temperature stability is directly affected by the stability of the room temperature. The VT unit compensates for about 80% of external changes (leaving 20% uncompensated for). Thus if the temperature of the room changes by 1°, the sample temperature will change by about 0.2°, which will not be reflected in a change in the numerically displayed temperature. For best results, the room temperature should be made as stable as possible. Any cycling of the temperature due to air conditioning or heating should be accomplished with the shortest possible cycle time and the minimum possible temperature variation.
- High-stability and independence from room temperature can be achieved if the VT controller is equipped with an optional cold junction (CJ) compensator. With the high-stability feature, the VT controller is no longer compensated for room temperature changes, but instead receives its reference from the cold junction devices mounted in the magnet leg. As the CJ compensator reduces the room temperature influences on the system, the influences of the VT gas supply become more apparent. For optimum performance of the CJ compensator, the flow and temperature of the VT gas supply must be as stable as possible.
- A possible setup to help stabilize the VT gas supply is to run the VT gas through a heatexchanger coil in a water bath at a regulated temperature. For best results, use an ice bath to cool down the VT gas to between 5°C and 10°C, and keep the flow as stable as possible for experiments below 40°C. Generally, for best performance of the VT controller and heater in the probe, the VT gas supply temperature should be a minimum of 10°C below the set temperature.
- For exact determination of sample temperature, a temperature calibration curve must be made for each probe used. All data, such as gas flow, must be noted. Samples of ethylene glycol are used for high-temperature calibration, and samples of methanol are used for low-temperature calibration.
 - a. After bringing the sample to the desired temperature and allowing sufficient time for equilibration, obtain a spectrum.
 - b. Display two cursors and align them on the two resonances in the spectrum.
 - c. If the sample is ethylene glycol, enter tempcal('e'); if the sample is methanol, enter tempcal('m').
 - d. The temperature is calculated and displayed based on the difference frequency between the cursors.

8.5 VT Controller Safety Circuits

The VT controller includes safety circuits to avoid damage to the heating element and probe. The following error conditions produce an error code in acqstatus[2]:

- Open circuit in the thermocouple circuit.
- Open circuit, short circuit, or over-temperature at safety sensor.
- Short circuit or software/microprocessor failure at the output transistor.

Over-temperature at the safety sensor initially turns off the heater. If this method fails to correct the condition within 5 seconds, either the gas flow has been interrupted or an output

transistor failure has occurred, whereupon a protective relay operates, isolating the heater from the control electronics. Failure of any of the sensors also results in this relay operating.

Once the protective relay has operated, the output will remain off. A power-down and power-up cycle of the VT controller is required to release the relay.

The over-temperature circuit can be inadvertently tripped if the VT is started at a below ambient temperature and the temperature is increased greater than 70°C. If the circuit is tripped, reset it by turning the VT off and on, then change to the desired temperature in 50°C steps.

Excessive heat requirements that cause the current to remain near the maximum can also trip the second circuit. Therefore, when using liquid nitrogen for cooling and when operating from 0°C through +25°C, reduce the gas flow rate to between 8 and 9 LPM. Reset will also occur if the VT cable is removed from the probe while the VT is on.

Refer to the VT installation manual for system failure analysis.

8.6 VT Interlock Parameters

The tin (temperature interlock) and vtwait (VT wait time) parameters check VT operation and stop the experiment if temperature regulation is lost. The tin parameter functions much the same as the in (lock and spin interlock) parameter:

- If tin='y', the VT regulation light is monitored during the course of the experiment, and if it starts to flash (regulation lost), the current data acquisition is stopped. The acquisition does not resume automatically if regulation is regained.
- If tin='w', the VT regulation light is monitored during the course of the experiment, and if it starts to flash (regulation lost), a warning is generated but acquisition is not stopped.
- If tin= 'n', the temperature interlock feature is turned off.

For both tin='y' and tin='w', the lost regulation causes werr processing to occur, thus providing a user-selectable mechanism to respond to VT failure.

The interlock operation does not apply to the cases when VT regulation is temporarily lost as a result of a programmed temperature change in an experiment where temp is an array. Also if tin='y', a maximum limit is imposed on the time that the system waits for regulation to be established. This limit is determined by the vtwait parameter and is independent of the pad parameter. If regulation is not established after the vtwait time (normally set to 180 seconds), the system displays the message VT FAILURE and does not proceed with the experiment. If the regulation problem is later corrected, the experiment can be resumed with the go and ga commands.

The VT gas flow has no sensor or interlock. If gas flow stops, the heater is protected by an internal temperature limit sensor that turns off the heater current before the element overheats. Because a loss of gas flow will result in a loss of regulation, any experiment in progress is stopped if tin='y'. Only the sample is left unprotected if VT gas stops.

Although tin and vtwait are not part of a parameter display group, each can be checked and changed in value the same as other parameters (e.g., by entering tin? to check the value or by entering tin='y' to set the value).

CAUTION: Do not run unattended a sealed sample of highly volatile materials that must be kept cold to avoid excessive pressure buildup. The undetected loss of VT gas or exchanger coolant could result in the rupture of the sample tube and damage to the probe components. Chapter 8. Variable Temperature Operation

Chapter 9. Carousel, SMS, and NMS Automation

Sections in this chapter:

- 9.1 "Carousel Autosampler," this page
- 9.2 "SMS Autosampler," page 280
- 9.3 "NMS Autosampler," page 285
- 9.4 "General Automation Tasks For All Sample Changers," page 289
- 9.5 "Changing Sample Changers or Serial Ports," page 295
- 9.6 "Using Gradient Autoshimming with Automation," page 295
- 9.7 "Automation Run Description," page 295
- 9.8 "Customizing the Sample Entry Window," page 300
- 9.9 "Automated Data Acquisition," page 301
- 9.10 "Automated Data Processing," page 308
- 9.11 "File Structures in an Automation Run," page 308

See Chapter 10, "VAST Accessory Operation," for coverage of the VAST sample changer. This chapter is limited to the SMS, Carousel, NMS, and the ASM-100 sample changers.

The SMS autosampler, Carousel autosampler, NMS, and ASM-100 automatic sample changers mechanically manage the process of removing and inserting samples from a magnet as part of spectrometer operation. With a minor extension of "normal" operating procedures, however, the sample changer can also be integrated into the process of experiment queueing.

The VNMR software can handle as many as 9999 data sets simultaneously in separate experiment files labeled exp#, where # is the experiment number. Without a sample changer, these experiment files can be used to queue up a series of NMR experiments (up to nine) on the same sample. By adding a sample changer, these experiments can instead be used for experiments on up to 100 different samples; the system identifies the sample by location number (the loc parameter) as one of the parameters in the experiment.

The sample changer is most useful in fully automatic, unattended operation, called *automation mode*, or an *automation run*.

9.1 Carousel Autosampler

This section is organized as follows:

- "Configuring VNMR for the Carousel," page 270
- "Checking Out the Carousel," page 271
- "Mounting and Removing the Carousel," page 273
- "Adjusting the Eject Air," page 274

- "Loading and Unloading Samples," page 275
- "Running NMR on One Sample at a Time," page 276
- "Running Automated NMR on Up to Nine Samples," page 277
- "Inserting Samples Manually with the Carousel Attached," page 278
- "Carousel Error Codes and Recovery," page 279

A typical use scenario for the Carousel autosampler is as follows:

- Remove the carousel from the driver (Figure 60 is an illustration of the carousel and driver).
- Load samples and turbines.
- Place loaded carousel onto the driver.
- Use the enter program to set up the automation run.
- Enter autogo to begin the automation run.

If you have not used the Carousel autosampler before or if you want to reacquaint yourself with it, follow the steps in the next three sections.

After you are familiar with the carousel and have attached it to the driver, go to "Adjusting the Eject Air," page 274.




Figure 60. Carousel Autosampler Carousel and Driver

Configuring VNMR for the Carousel

Before you can use the carousel, you must select the carousel in the VNMR configuration window and edit the enter.conf file.

- 1. Log in as **vnmr1**.
- 2. Enter config in the VNMR input window.
- 3. In the VNMR Configuration window:
 - Set Sample Changer to Carousel.
 - Set **Sample Changer Serial Port** to **Port A** or **Port B**, depending on which port is used to connect the sample changer.
 - Click Exit and Save.
- 4. Using a text editor, such as vi, open the file /vnmr/asm/enter.conf.
- 5. Find the line set loc(max) at about line 22 in the file. Change the number, if necessary, so that the line is the same as the following:

```
set loc(max) 9
```


6. Save the changes.

VNMR is now configured to use the carousel.

Checking Out the Carousel

Inspecting the carousel before mounting it on the driver helps you understand how it works.

1. With the sensor arm locked (See Figure 61 and Figure 62), look into the carousel tubes. The sensor arm is locked when the arm is positioned between position 9 and position 1 and when the arrow on the knurled knob is pointing to the left or right. The pin of the knurled knob should be in the hole located between positions 9 and 1 on the upper white disk of the carousel (see Figure 62).






Figure 62. Lock Pin on Sensor Arm Disengaged and Engaged

In *locked mode*, the retaining disk appears in all nine tubes, preventing samples from falling through the bottom. In "run mode," the retaining disk leaves a gap in the tube over the upper barrel, allowing the sample to drop into the magnet.

- If the sensor arm is locked, disengage the lock pin by pulling out the knurled knob and rotating it 90°. The arrow on the knurled knob now points up or down. See Figure 62. Notice how the sensor arm is free to rotate.
- 3. Move the sensor arm to the right, from between positions 9 and 1 to directly over position 9. When the carousel is mounted to the driver, the arm only moves between the run and locked positions, shown in Figure 61.
- 4. Look through the tubes.

Notice that the white retaining disk inside the carousel has a cutout that aligns with the sensor arm. This cutout allows one sample to drop into the magnet, while the other eight samples rest on the retaining disk.

5. Look at the bottom of the carousel. You see a connector and a groove that fit into the connector and alignment bar on the top of the driver, as shown in Figure 63. Notice the orientation of the connector and the alignment groove and compare this to the driver connector and alignment bar.




Figure 63. Bottom View of Carousel, Showing Connector and Groove

- 6. Move the sensor arm to between position 9 and position 1, and then engage the lock pin by rotating it $90^{\circ\circ}$ and allowing it to fall into the hole in the upper white disk.
- 7. Place the empty carousel on the driver according to the instructions in the next section, "Mounting and Removing the Carousel."

After the empty carousel is mounted, it is ready to be used in walk-up mode, as described in "Inserting Samples Manually with the Carousel Attached," page 278.

Mounting and Removing the Carousel

To Mount the Carousel

This procedure describes how to mount the carousel on the driver.

- 1. Make sure the sensor arm on the carousel is locked in between positions 9 and 1. The groove in the bottom of the carousel (see Figure 63) should be 90° (perpendicular) to position 9.
- 2. Make sure the driver has position 9 over the upper barrel. If not, use the manual index button to rotate the driver.
- 3. Align the numbers on the carousel to the numbers on the driver.
- 4. Push down on the carousel until it seats. This engages the electrical connector and drive pin (see Figure 63).
- 5. At the console, turn on the eject air.
- 6. Disengage the lock pin by pulling out the knurled knob and rotating it 90°. The arrow on the knurled knob will point up or down.
- 7. Move the sensor arm to the right, until it is over position 9. This is the run position. You will notice some resistance as you move the arm.

The carousel is now ready to use for NMR experiments.

To Remove the Carousel

This procedure describes how to remove the carousel from the driver. Position 9 must be over the upper barrel before the carousel is removed.

- **CAUTION:** Samples could fall through and break if the carousel is not in the locked mode. Before removing the carousel, rotate the carousel so that position 1 is over the upper barrel. The sensor arm must be in the locked position, as shown in Figure 61, with the lock pin engaged and the arrow on the knurled knob pointing to the left or right.
 - 1. After an automated run is finished, turn on the eject air at the console.
 - 2. Use the manual index button to rotate the carousel, until position 9 aligns with the upper barrel.
 - 3. Move the sensor arm to the left, from run to lock mode (see Figure 61), until the arm stops between position 9 and position 1.

In the lock position, all the samples rest on the retaining disk.

4. Engage the lock pin by pulling the knurled knob and turning it 90°. The lock pin fits into the hole in the white upper disk of the carousel.

The locking pin keeps the sensor arm from moving out of the lock position.

5. Remove the carousel by pulling up on the carousel while holding down the driver, rocking gently if necessary. The upper barrel may lift slightly. The sensor arm provides a convenient, central gripping point.

CAUTION: Without the carousel mounted, the eject air pressure is too high and samples may launch out of the magnet when ejected. The procedure in "Inserting Samples Manually with the Carousel Attached," page 278, is recommended for nonautomated work. Otherwise, lower the eject air pressure if you intend to use the magnet without the carousel. Refer to the next section, "Adjusting the Eject Air."

Adjusting the Eject Air

This procedure describes how to adjust the eject air pressure so that the sample and turbine float where the optical sensor on the sensor arm can *see* the sample. Figure 64 shows the optical sensor and the correct sample floating height. The empty carousel should be mounted on the driver as described in the procedure for mounting the carousel on the previous page.

Variation in supply air pressure can cause failures. Adjust the eject air to the lowest level expected from the supply air during a 24-hour period that maintains an acceptable turbine float height for the sensor.




Figure 64. Optical Sensor and Proper Sample Floating Height

 Make sure the air hoses for body and VT air are attached to the probe and flow meters are set for about 11 cfm.

Some probes require more eject air than others. If the pressure required to eject the sample and float the turbine adequately exceeds 65 psi, check for the guide hole at the bottom of the probe.

If your probe has a guide hole in the bottom for a capacitor or inductor stick and if your experiment does not require the stick, plug the guide hole with the provided rubber stopper.

- 2. In VNMR, turn on the eject air by typing **e** and then pressing Return.
- 3. With the knurled knob disengaged (see Figure 62), move the sensor arm to the left so that it clicks between positions 9 and 1 (see Figure 61).
- 4. Place a sample in position 9 of the carousel.
- 5. Move the sensor arm back in front of position 9.
- 6. Look through the round hole in the position 9 tube of the carousel.

The sample should be floating high enough for the optical sensor to see the turbine, as shown in Figure 64. In other words, the turbine should cover at least half of the inside of the round hole.

• If the turbine is too low, turn up the eject air at the source pressure regulator. If eject air is still insufficient, air may be leaking through a tuning capacitor or inductor guide hole in the bottom of the probe. Use the provided rubber stopper to plug the hole.

- If the turbine is too high, turn down the eject air at the source pressure regulator.
- 7. After the sample is floating at the correct height, insert the sample using VNMR by typing **i** and pressing Return.

The sample should drop into the probe. If the sample drops too rapidly, adjust the SLO-DROP air in the magnet leg.

8. Eject the sample again and recheck the height. Adjust the eject air again if necessary.

The Carousel autosampler can now be used for automation. You may have to lower the bearing air to compensate for the higher eject air. Proceed to the next section, which describes how to load samples into the carousel.

Loading and Unloading Samples

The Carousel autosampler provides two methods for loading and unloading samples. The only difference between the two procedures is whether the carousel is left on the driver or removed.

Do not mix turbine types within the carousel, and be sure that eject air is on before the sensor arm is moved from the locked to the run position.

WARNING: Removing the carousel while reaching over a high-field magnet (such as 400-, 500-, 600-, or 750-MHz) could cause injury if you to lose your balance and fall. Take care when removing the carousel from a high-field magnet.

To Load and Unload Samples with the Carousel Installed

This procedure describes loading and unloading samples from the carousel while it is still attached to the driver. If you want to remove the carousel before loading and unloading samples, use the next procedure.

- 1. Enter **e** to turn on the eject air.
- 2. Move the sensor arm to the left, until it stops between position 9 and position 1.
- 3. Remove the sample from position 9. Then, remove the rest of the samples.
- 4. With the eject air still on, insert the samples into the carousel in any order. Be sure to remember which samples are in which positions.
- 5. Move the sensor arm back to the right, until it aligns with position 9.

The sample in position 9 will float. Now is a good time to check sample float height. Adjust eject air if necessary, as described in "Adjusting the Eject Air," page 274.

The Carousel autosampler can now be used for NMR.

To Remove the Carousel and Load Samples

This procedure describes removing the carousel then loading or unloading samples from the carousel.

- 1. Enter **e** in VNMR to turn on the *eject* air.
- 2. If position 9 is not already aligned over the upper barrel, use the manual index button to index the carousel to position 9, or enter **loc=9** change in VNMR.

- 3. Move the sensor arm to the locked position and engage the lock pin as follows. Carousel must be in the lock position when removed; otherwise, samples could fall through the carousel and remounting the carousel would be extremely difficult.
 - a. Move sensor arm to the left, until it stops between position 9 and position 1.
 - b. Engage the lock pin by rotating the knurled knob 90°, until the arrow points left or right.
- 4. Remove the carousel by pulling up on the carousel while holding down the driver, rocking gently if necessary. The upper barrel may lift slightly. The sensor arm provides a convenient, central gripping point.
- 5. Place the carousel on a table.
- 6. Remove the samples from the carousel.
- 7. Load the new samples into the carousel.
- 8. Mount the carousel on the driver. Match the sample position numbers on the carousel to the numbers on the driver. Press down on the carousel until it seats.
- 9. Enter e in VNMR to turn on the eject air.
- 10. Disengage the lock pin and move the sensor arm to the right, until it aligns with position 9.

The sample in position 9 should float at the correct height. If not, adjust the eject air as described in "Adjusting the Eject Air," page 274.

The Carousel autosampler can now be used for NMR.

Running NMR on One Sample at a Time

You can run NMR on any individual sample in the carousel using the loc parameter and the change command.

- 1. Make sure the carousel is loaded with nine samples and is installed on the carousel driver, as described in previous sections.
- 2. Enter **e** in VNMR to turn on the eject air.
- 3. Disengage the lock pin and then move the sensor arm to the right until it stops in front of position 9.
- 4. Check that the eject air pressure is sufficient to float the turbine/sample to the correct height. If not, adjust the eject air as given in "Adjusting the Eject Air," page 274.
- 5. Using the manual index button, rotate the carousel through all positions. Check the turbine/sample float height for each position. Stop when position 9 is over the upper barrel.
- 6. Enter **loc=1** change to change to position 1.

The carousel does the following:

- The message expl: Experiment started appears on the Sun computer.
- The carousel rotates a half position (in this case between positions 9 and 1).
- The carousel rotates to position 1 and inserts the sample.
- The message expl:Setup complete appears on the Sun computer.
- 7. Run the desired NMR experiment.
- 8. Enter **loc=2** change to change to position 2. Run an experiment if desired.

Running Automated NMR on Up to Nine Samples

To run a series of NMR experiments on some or all of the samples in the carousel, load up to nine samples into the carousel and use the enter program to set up the automation run.

Begin the automation run with the sample in position 9 inserted into the magnet. When the automation run begins, the Carousel ejects the sample from location 9, rotates to location 1, and inserts the sample from location 1.

After an automation run finishes, the last sample location run is left in the magnet.

- 1. Prepare your samples and have a list of samples and experiments ready.
- 2. Load up to 9 samples as described in "Loading and Unloading Samples," page 275.
- 3. If not already done, insert the sample from location 9.

When the automation run begins, the Carousel ejects the sample from location 9, rotates to location 1, and inserts the sample from location 1.

- 4. Enter **cd** to change to your home directory.
- 5. Type enter in the VNMR input window.

You are prompted for a directory name that stores the information set by the enter program.

You can also use the command enter('abc') creates a directory named abc. In the abc directory, there is a file named abc, which contains experiment information. Also in the directory is a directory named abc.macdir, which contains *GLIDE*-related information for an automation run.

The Sample Entry Form window (see Figure 65) appears.

			Sai	mple Entry Forr	n		· 🗆				
	Sample Number										
123456	1 2 3 4 5 6 7 8 9										
	User identification										
vnmr1											
			So	lvent Select	ion						
CDC13	D20	Benz	zene	DMSO							
Acetone	cetone Cyclohexane Toluene Methanol										
			E×pe	riment Selec	ction						
H1			C13			F19					
P31			HC			H-gCOSY					
H-TOCSY			H-gHSQC			H-gCOSY-gHSQC					
H-gCOSY-HSQ	C-gHMBC		H-gCOSY	-HSQC-gHMBC-	HSQCTOXY	Y H-COSY-C-DEPT-HETCOR					
H-COSY-C-AP	Т		H-TOCSY	-NOESY		H-TOCSY-ROESY					
H-TOCSY-HMQ	C										
Text											
Customi:	ze Parameters	3	Ac	dd Entry	E×i	t and Save	Quit				
Number of sa	amples submit	ted:	0								

Figure 65. Sample Entry Form Window for the Carousel Autosampler

- 6. Fill in the Sample Entry Form window.
 - **Sample Number** select one or more sample locations on which to run experiments.

- User Identification select the user.
- Solvent Selection select the solvent for the chosen locations(s).
- **Experiment Selection** select one or more experiments for the chosen locations(s).
- Text enter information about the experiment, if desired.
- To customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles, click the **Customize Parameters** button.

Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

• When finished with this location (or locations), click Add Entry.

When you are done adding entries, click Exit and Save.

- 7. To start the automation run, enter **autogo**.
 - When the system asks Location of automation queue, enter the directory name you used in step 5.
 - When the system asks Location of automation data, accept the default or enter a new directory name.

You can also enter autogo('MySamples') or autogo('MySamples','AutoRun_621').

The data is stored in the automation directory, with the usual automation-style file names (as specified by autoname if desired). The acquisition takes place in background within VNMR, allowing the user complete freedom to process other dataset in any desired experiment location.

- 8. If any error messages appear, refer to "Carousel Error Codes and Recovery," page 279 for an explanation.
- 9. To monitor the automation run, enter **status** to bring up the Status window. Refer to "Monitoring an Automation Run," page 291 for more details.
- 10. Return later to retrieve your results.

For each experiment, a plot is created and a FID file is saved in the automation directory (e.g., AutoRun_621).

The FID files and the fid extension are the data files collected during the acquisition run. Normally, as data is collected, the FIDs are placed into the automation directory with the name xxyy.fid, where xx is the sample location number and yy is an "experiment number" on that particular sample (1 through n, where n experiments have been run on the same sample).

Inserting Samples Manually with the Carousel Attached

Since the eject air is higher than normal, the carousel should remain attached to the driver, even when empty. Therefore, the Carousel autosampler provides a mode that allows you to insert and eject samples when the carousel is not being used for automation. In other words, if you want to run one sample at a time, as in normal operation, you should do so with the carousel in place.

This procedure describes how to insert and eject samples, one at a time, through the carousel. The carousel must be on the magnet with no sample in position 9.

1. Enter **e** in VNMR to turn on eject air.

- 2. Make sure no sample is in position 9, and make sure position 9 is aligned over the upper barrel.
- 3. If the sensor arm is over position 9, move it to the left until it stops between position 9 and position 1, as shown in Figure 66. You will notice some resistance moving the arm. Do not engage the locking pin.




Figure 66. Manually Loading and Unloading Samples Through Position 1

- 4. Place your sample and turbine in position 9.
- 5. With eject air still on, rotate the sensor arm back to position 1.

The sample can be inserted and ejected, as in normal, nonautomated operation.

6. To remove the floating sample from the magnet, rotate the sensor arm to the left until it stops between position 9 and position 1. Then lift the sample straight up and out.

Carousel Error Codes and Recovery

This section defines the error codes and suggests ways to recover from the errors. The following messages are common to all Varian autosamplers. The accompanying descriptions and remedies, however, refer to the Carousel autosampler.

Sample changer arm unable to move sideways during retrieve.

Index air is too low. Check and adjust the index air regulator to increase the air pressure to 38 psi.

Sample tray size is not consistent.

System is configured with wrong tray size. Open the CONFIG window and make sure Sample Tray Size is set to 9.

Invalid sample number during retrieve.

Location value too high for Carousel and traymax configuration incorrect. Only positions 1 through 9 are valid for the Carousel. Open the CONFIG window and make sure Sample Changer is set to Carousel (or Sample Tray Size is set to 9).

Sample out of range during automatic retrieve.

Sample tube is floating too low. Slightly increase the system eject air one of the following ways:

- Increase the cooling air flow.
- Make sure the body air and VT air hoses are connected to the probe and the flow meters are set to 11 cfm.
- Increase supply air pressure.
- Plug the capacitor/inductor guide hole in the bottom of the probe with the provided rubber stopper.

Robot arm failed to find home position during retrieve.

Carousel is not in an indexed position. Adjust the air adjust knob on the driver so that the carousel rotates properly. The micro switch should insert fully into the recess on the underside of the white disk of the driver.

- If the carousel rotated too far, turn the air adjust knob 1/16 of a turn at a time in the clockwise direction. Test using the manual index button (press and hold the manual index button and wait for motion).
- If the carousel does not rotate far enough, turn the air adjust knob 1/16 of a turn in the counter-clockwise direction. Test using the manual index button (press and hold the manual index button and wait for motion).

Air supply to sample changer failed during retrieve.

No tube or sample is detected. Verify that a sample is in the magnet or carousel.

- If no sample is present, do an eject and place a sample in the carousel. Then, continue to the next sample.
- If a sample is present, check the eject air and adjust if necessary, as described in "Adjusting the Eject Air," page 274.
- If the sample sticks in the lower part of the barrel, be sure the carousel is fully seated on the driver.

9.2 SMS Autosampler

This section is organized as follows:

- "Configuring VNMR for the SMS Autosampler," this page
- "Preparing Sample Tubes," page 282
- "Running NMR on One Sample at a Time," page 282
- "Running Automated NMR," page 282
- "SMS Error Codes and Recovery," page 284

The SMS (Sample Management System) autosampler accessory (illustrated in Figure 67) is an electromechanical, robotic device designed to automatically insert and retrieve samples from the bore of the magnet. The system is designed for completely automatic, unattended operation and can handle from 1 to 100 samples.

When installed, the robot is attached to a table next to the magnet, where the robot is calibrated and "taught" the magnet and sample rack positions.

Two sample tray configurations are available:

- A 50-sample tray accommodates 5-mm and 10-mm sample tubes.
- Two 50-sample trays accommodate one hundred 5-mm and 10-mm sample tubes.

Each sample position is individually encoded so that the computer can address a position either randomly or sequentially. In addition, sample locations are staggered to improve the




Figure 67. SMS Autosampler

identification of each sample position. Sample zero is a default position accessible by the autosampler for various error recovery scenarios, including power failures.

The autosampler and the acquisition system communicate through a standard RS-232C serial link at 9600 baud. Time between sample changes is about 35 seconds for a complete cycle—retrieve the sample from the magnet, return it to the sample tray, locate and pick up the next sample, and deliver it to the magnet.

Configuring VNMR for the SMS Autosampler

Before you can use the SMS, VNMR must be configured as follows:

- 1. Log in as vnmr1, and enter **config**.
- 2. In the VNMR configuration window, make the following changes:
 - Sample Changer to SMS 50 Sample or SMS 100 Sample.
 - Sample Changer Serial Port to Port A or Port B, depending on which one of the Sun serial ports is used for the sample changer.
 - Click Exit and Save.
- To get the 100 sample positions (for dual 50-sample trays) to appear in the Sample Entry Form window, use a text editor, such as vi, to edit the file /vnmr/asm/enter.conf.
- Find the line set loc(max) at about line 22 in the file. Change the number, if necessary, so that the line is the same as the following: set loc(max) 100
- 5. Save the changes.

VNMR is now configured to use the SMS autosampler.

Preparing Sample Tubes

Every automation run starts with sample preparation. Tube lengths are constrained by the geometry of the sample changer to be either 7 in. or 8 in. tubes.

Do not attempt to use longer (e.g., 9 in.) tubes because the top of such a tube may come into contact with the bottom of a tube being held by the gripper arm.

Positioning the sample in the spinner and adjusting the height of the liquid in the tube are both done the same way as in manual operation, but position and sample height are much more critical in automated operation. The more precise your tube position and the more you are careful to adjust each sample to exactly the same height, the fewer problems you will have with automatic locking and shimming, and the less time you will have to take for these operations to occur.

If at all possible, fill tubes to a height of 50 mm (about 2 inches).

Running NMR on One Sample at a Time

You can run NMR on any individual sample in the sample tray using the loc parameter and the change command.

- 1. Prepare the samples you want to run and place them in the sample tray. It might be helpful to write down the location of each sample.
- 2. Enter **loc=#** change, where # is the location of the sample you want to run.

For example, to insert the sample at location 3, enter **loc=3 change**.

The SMS does the following:

- The message expl: Experiment started appears on the Sun computer.
- The SMS arm retrieves the sample from the specified location.
- The eject air turns on, the SMS inserts the sample, and the sample inserts into the probe.
- The message expl:Setup complete appears on the Sun computer.
- 3. Run the desired NMR experiment.
- 4. Enter **loc=# change**, where # is the location of the next sample you want to run.

The previous sample is removed from the magnet and returned to its location. The next sample is retrieved from the tray and inserted into the magnet.

Run an experiment if desired.

Running Automated NMR

To run a series of NMR experiments on some or all of the samples in the sample tray(s), load up to 50 samples for a single 50-sample tray, or up to 100 samples for dual 50-sample trays, and use the enter program to set up the automation run.

- 1. Prepare the samples you want to run and place them in the sample tray. It might be helpful to write down the location of each sample.
- 2. Enter **cd** to change to your home directory.
- 3. Type **enter** in the VNMR input window.

You are prompted for a directory name to store the information set by enter.

You can also invoke the command enter (*directory*) to create a directory. For example, enter('abc') creates a directory named abc. In abc, a file named abc contains experiment information. Also in the abc directory, the directory abc.macdir contains *GLIDE*-related information for an automation run.

The Sample Entry Form window (see Figure 68) appears.

	Sample Entry Form																			
	Sample Number																			
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	
21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	
41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	
61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	
81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100	
vnm	User identification																			
								So	olve	nt S	elec	tion	1							
CDC	13		D20			Ber	izene	2	DMS	50										
Ace	tone		Cyc	lohe	exane	Tol	uene	2	Me	than	51									
								Е×р	erim	nent	Sele	ctio	on							
H1							C13	3						F1	9					
P31							HC							H-	gCOS	SY				
H-T	OCSY						H-g	gHSQC	2					<u>H-</u>	gCOS	SY-gł	HSQC			
H-g	COSY	-HSQ	C-g⊦	IMBC			H-8	gCOS1	r–HS	QC-g	нмвс	-HSQ	сто×	Y H-	COSY	'-C-I	DEPT	-HET	COR	
H-C	OSY-	C-AP	Т				H	FOCSY	r-NO	ESY				<u>H</u> -	TOCS	SY-RO	DESY			
H-T	OCSY	-HMQ	С																	
Text	Text Customize Parameters Add Entry Exit and Save Quit																			
Numb	per d	of sa	amplo	es s	ubmi	tted	: 0													

Figure 68. Sample Entry Form Window (100 samples) for the SMS Autosampler

- 4. Fill in the Sample Entry Form window.
 - **Sample Number** select one or more sample locations on which to run experiments.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen location(s).
 - **Experiment Selection** select one or more experiments for the chosen location(s).
 - Text enter information about the experiment, if desired.
 - To customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles, click the **Customize Parameters** button.

Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

• When finished with this location (or locations), click Add Entry.

When you are done adding entries, click Exit and Save.

5. To start the automation run, enter **autogo**.

- When the system asks Location of automation queue, enter the directory name you used in step 3.
- When the system asks Location of automation data, accept the default or enter a new directory name.

```
You can also enter autogo('MySamples') or autogo('MySamples','AutoRun_621').
```

The data is stored in the automation directory, with the usual automation-style file names (as specified by autoname if desired). The acquisition takes place in background within VNMR, allowing the user complete freedom to process other dataset in any desired experiment location.

- 6. If any error messages appear, refer to "SMS Error Codes and Recovery" for an explanation.
- 7. To monitor the automation run, enter **status** to bring up the Status window. Refer to "Monitoring an Automation Run," page 291, for more details.
- 8. Return later to retrieve your results.

For each experiment, a plot is created and a FID file is saved in the automation directory (e.g., AutoRun621).

The FID files the fid extension and are the data files collected during the acquisition run. Normally, as data is collected, the FIDs are placed into the automation directory with the name xxyy. fid, where xx is the sample location number and yy is an "experiment number" on that particular sample (1 through n, where n experiments have been run on the same sample).

SMS Error Codes and Recovery

This section contains a list of the error messages returned by the spectrometer when there is a problem. Under each error message is an explanation and correction.

General Operation Errors

The following errors indicate a problem with the general operation of the autosampler.

• ERROR 96

SMS autosampler is plugged into the RS-232C but not turned on.

• ERROR 97

SMS autosampler is not plugged into the RS-232C connector.

• ERROR 98

If you receive ERROR 98 while executing loc=x or LOC=X change command or while attempting an automation run, reinitialize the System V Controller by turning its power off and then on. Wait for the beep to indicate that it is ready. Try the operation from the console.

• ERROR 99

SMS autosampler is not responding to issued commands.

Sample Removal Errors

The following errors indicate a problem with removing samples from the sample tray.

• ERROR 01 NO SAMPLE

The gripper has been commanded to close with no sample because the fingers close completely. This may be a normal indication, for example, when executing R0 (retrieve sample 0) to check the magnet for a sample. If a fault occurs, the sample switch may need adjustment.

• ERROR 05 INVALID SAMPLE NUMBER

The sample number exceeds the positions available on the tray. This error message appears if you attempt to load sample 0.

• ERROR 07 GRIPPER ABORT

The gripper is compressed and unable to operate, usually because it has contacted another surface.

• ERROR 13 ILLEGAL COMMAND (shutdown?)

An illegal character has been received. This is a normal response to any command except E (energize) or H (help) when in the shutdown mode.

Sample Loading Errors

If a sample removal error occurs during sample loading, the same error message will appear on the monitor. However, the error number will be different; the error number will equal the corresponding sample removal error number plus 20. For example, if a gripper abort error occurs during sample loading, error number 27 (7 + 20) will appear on the monitor.

9.3 NMS Autosampler

This section is organized as follows:

- "Configuring VNMR for the NMS Autosampler," page 286
- "Running NMR on One Sample at a Time," page 286
- "Running Automated NMR," page 287

The Nano Multisampler (NMS) automates the steps required to position a 4-mm rotor into the bore of a magnet for NMR analysis. The rotor is spun at a magic angle of 54.7° in a stator within the body of a Nano probe. The probe is subsequently lifted into the bottom bore of an NMR magnet.


The NMS is comprised of the following main units:

- *Probe elevator* Attaches to the shim coil flange at the bottom of the magnet, in place of a probe. This unit lifts the probe into the magnet bore and lowers the probe out of the magnet bore.
- *Probe flange* Fits onto the probe shield and then snaps into the elevator.
- *NMS Controller* Sits on the floor near the magnet. The NMS controller provides serial (RS-232) communications between the NMS and the NMR host computer. The NMS controller also provides signal and pneumatic connections to the probe elevator and the carousel.

• *Carousel* – Sits on the floor aligned with the tilted probe. The carousel holds a 48-sample rack and replaces the sample in the probe with a new one from the rack.

The probe elevator and NMS controller can be used together to provide a semi-automated means of inserting a Nano probe into the NMR magnet. The addition of the carousel provides an automated means for loading each of 48 rotors into a probe.

The NMS window, shown in Figure 69, provides general controls for installing and operating the NMS sample changer.



Before Using NMS

Figure 69. NMS Window

Before the NMS autosampler is used, the probe elevator must be attached to the shim coil and the carousel must be properly positioned. *No sample can be in the probe when you begin using the NMS autosampler.*

Configuring VNMR for the NMS Autosampler

Before you can use the NMS, make sure VNMR is configured as follows:

- 1. Log in as vnmr1, enter config.
- 2. In the VNMR configuration window, make the following selections:
 - Sample Changer to NMS.
 - Sample Changer Serial Port to Port A or Port B, depending on which one of the Sun serial ports is used for the sample changer.
 - Click Exit and Save.
- 3. To get the proper number of sample positions to appear in the Sample Entry Form window, use a text editor, such as vi, to edit the file /vnmr/asm/enter.conf.

Find the line set loc(max) at about line 22 in the file. Change the number, if necessary, so that the line is the same as the following: **set loc(max) 48**

4. Save the changes.

VNMR is now configured to use the NMS.

Running NMR on One Sample at a Time

You can run NMR on any individual sample in the sample tray using the loc parameter and the change command.

- 1. Prepare the samples you want to run and place them in the sample tray. It might be helpful to write down the location of each sample.
- Enter loc=# change, where # is the location number of the sample you want to run. For example, to insert the sample at location 3, you would enter loc=3 change.

The NMS does the following:

- The message expl: Experiment started appears on the Sun computer.
- The probe elevator lowers the probe out of the magnet and tilts the probe.
- The suction cup arm on the carousel unit removes the sample from the probe stator and returns it to the tray, the sample tray rotates to the desired location, and the suction cup arm retrieves the desired sample and places it into the probe stator.
- The message expl:Setup complete appears on the Sun computer.
- 3. Run the desired NMR experiment.
- 4. Enter **loc=# change**, where # is the location of the next sample you want to run.

The previous sample is removed from the magnet and returned to its location. The next sample is retrieved from the tray and inserted into the magnet.

Run an experiment if desired.

Running Automated NMR

To run a series of NMR experiments on some or all of the samples in the sample tray, load up to 48 samples, and use the enter program to set up the automation run.

- 1. Prepare the samples you want to run and place them in the sample tray. It might be helpful to write down the location of each sample.
- 2. Enter **cd** to change to your home directory.
- 3. Type **enter** in the VNMR input window.


You are prompted for a directory name to store the information set by enter.

You can also invoke the command enter (*directory*) to create a directory. For example, enter('abc') creates a directory named abc. In abc, a file named abc contains experiment information. Also in the abc directory, the directory abc.macdir contains *GLIDE*-related information for an automation run.

The Sample Entry Form window (see Figure 70) appears.

- 4. Fill in the Sample Entry Form window.
 - **Sample Number** select one or more sample locations on which to run experiments.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen locations(s).
 - Experiment Selection select one or more experiments for the chosen locations(s).
 - Text enter information about the experiment, if desired.
 - To customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles, click the **Customize Parameters** button.

— Sample Entry Form										
Sample Number 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 48										
vnmr1	User identification									
			So	lvent Select	ion					
CDC13 D20 Benzene DMS0										
Acetone	Cyclohexane	Tolu	uene Methanol							
			E×pe	eriment Selec	ction					
H1		1	C13			F19				
P31			нс			H-gCOSY				
H-TOCSY			H-gHSQC			H-gCOSY-gHSQC				
H-gCOSY-HSQ	IC-gHMBC		H-gCOSY	-HSQC-gHMBC-	HSQCTOXY	H-COSY-C-DEPT-HETCOR				
H-COSY-C-AF	Ϋ́Τ		H-TOCSY-NOESY			H-TOCSY-ROESY				
H-TOCSY-HMQ	iC									
Text Customize Parameters Add Entry Exit and Save Quit Number of samples submitted: 0										



Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

• When finished with this location (or locations), click Add Entry.

When you are done adding entries, click Exit and Save.

- 5. To start the automation run, enter **autogo**.
 - When the system asks Location of automation queue, enter the directory name you used in step 3.
 - When the system asks Location of automation data, accept the default or enter a new directory name.

You can also enter autogo('MySamples') or autogo('MySamples', 'AutoRun_621').

The data is stored in the automation directory, with the usual automation-style file names (as specified by autoname if desired). The acquisition takes place in background within VNMR, allowing the user complete freedom to process other dataset in any desired experiment location.

- 6. If any error messages appear, refer to "Carousel Error Codes and Recovery," page 279, for an explanation.
- 7. To monitor the automation run, enter **status** to bring up the Status window. Refer to "Monitoring an Automation Run," page 291, for more details.
- 8. Return later to retrieve your results.

For each experiment, a plot is created and a FID file is saved in the automation directory (e.g., AutoRun_621).

The FID files and the fid extension are the data files collected during the acquisition run. Normally, as data is collected, the FIDs are placed into the automation directory with the name xxyy.fid, where xx is the sample location number and yy is an "experiment number" on that particular sample (1 through n, where n experiments have been run on the same sample).

9.4 General Automation Tasks For All Sample Changers

The procedures in this section (listed below) and general automation tasks that can apply to all sample changers:

- "Preparing and Initiating an Automation Run," this page
- "Setting Up an Automation Run for Multiple Users," page 290
- "Monitoring an Automation Run," page 291
- "Using Sample Changers in Continuous Walkup Mode," page 293
- "Adding Samples to an Automation Run in Progress," page 294

Preparing and Initiating an Automation Run

- 1. Prepare your samples and have a list of samples and experiments ready.
- 2. Insert the samples in the carousel or sample tray (an automation run should not be in progress).
- 3. Enter **cd** to change to your home directory.
- 4. Type **enter** in the VNMR input window.

You are prompted for a directory name to store the information set by enter.

You can also invoke the command enter(*directory*) to create a directory. For example, enter('abc') creates a directory named abc. In abc, a file named abc contains experiment information. Also in the abc directory, the directory abc.macdir contains *GLIDE*-related information for an automation run.

The Sample Entry Form window appears (an example is shown in Figure 70).

- 5. Fill in the Sample Entry Form window:
 - **Sample Number** select one or more sample locations on which to run experiments.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen peak(s).
 - **Experiment Selection** select Autoscout, Autowettocsy, or both for the chosen peaks.
 - Text enter information about the experiment, if desired.
 - **Customize Parameters** Click this button to customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles.

Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

- When finished with this peak (or peaks), click Add Entry.
- 6. Repeat step 5 for each sample.

- 7. Click on the Exit and Save button when you are finished.
- 8. Enter autogo('mydata', 'myauto') to start the automation run.

The command autogo('mydata', 'myauto') creates a new automation directory (called myauto) in the users home directory and submits the experiment information contained in the mydata directory to the new myauto directory.

- 9. To monitor the automation run, enter **status** to bring up the Status window. Refer to "Monitoring an Automation Run," page 291, for more details.
- 10. Return later to retrieve your results.

For each experiment, a plot is created and a FID file is saved in the automation directory (e.g., myauto).

The FID files and the fid extension are the data files collected during the acquisition run. Normally, as data is collected, the FIDs are placed into the automation directory with the name xxyy.fid, where xx is the sample location number and yy is an "experiment number" on that particular sample (1 through n, where n experiments have been run on the same sample).

Setting Up an Automation Run for Multiple Users

Automation runs can be set up for multiple users.

The sample entry process might typically occur during the course of the day (assuming a common case of open access operation during the day, and sample changer run at night). One or more operators can use the enter program to create a file of information about the samples and experiments they wish to run.

All the information created by individual users must reside in a single file before the automation run is started. To get this final, single file, multiple users could do one of the following:

- Enter information into a single file consecutively.
- Enter information into different files simultaneously. The system operator must eventually merge the various files into a single file using a text editor.

Below are two examples of how this could be implemented.

Example of Using a Common File

- The system operator vnmrl creates an empty file in a location where everyone can access it. A simple way to create a file is to log in as vnmrl and enter the command touch /vnmr/samples. Also enter chmod a+w /vnmr/samples to give the file write access to all. Notify each user about the name of the file and give them a time to use it.
- 2. To add information to the agreed upon text file, each user runs the enter program by entering enter ('/vnmr/samples'), clicking on the Add Entry and then the Exit and Save buttons.
- 3. The sample changer operator supplies the name of the common file as the first argument to autogo, as shown in step 8 of the procedure "Preparing and Initiating an Automation Run," page 289.

Example of Using UNIX Mail to Gather Files

- 1. Each user enters the command **enter mydata** from the UNIX level, on any computer on the network that has VNMR software on it.
- 2. Next, the user enters sample information as noted in step 5 in "Preparing and Initiating an Automation Run," page 289, and then mails a message to the sample changer operator, like this: mail bill Subject: sample changer my sample changer info is in file /home/jim/mydata

Press Control-D to end the message.

3. The sample changer operator enters the following command (on one line): cat /export/home/jim/samples \ export/home/george/samples \ export/home/expts > /export/home/vnmr1/asm

This command combines files from all users and must be typed on a single line. Note that the order in which experiments are entered into this file will determine the order in which they are run.

- 4. The sample changer operator, using a UNIX text editor such as vi, edits the file just created, changing the sample location numbers to correspond to actual location numbers as samples are entered into the tray (as an alternative, the operator might preassign particular locations to certain users, and require them to have entered the correct numbers themselves in step 1).
- 5. The sample changer operator supplies the name of the common file as the first argument to autogo, as shown in step 8 of the procedure "Preparing and Initiating an Automation Run," page 289.

Monitoring an Automation Run

During an automation run, you can check to see if an experiment is queued, active, or completed.

• To monitor an automation run, enter **status** in the VNMR input window. The Sample Status window appears, as shown in Figure 71.

At the top of the Sample Status window, a table summarizes how many samples are queued, how many are active, etc. The color coding in the table identifies the status of each sample in the list below. For example, if the row in the table for queued samples is blue, each entry in the list below displayed in blue indicates a queued sample.

- To see detailed information of a particular entry, **click on the entry to highlight it**. The detailed information appears in the Log information section of the Status window.
- To retrieve the selected entry into the currently-joined experiment, click on **Retrieve Data**. This button is active only if data is present in the selected entry.




Figure 71. Sample Status Window (status Program)

• To sort the entries or find a particular entry in the Status window, click on **Find Entry**. The locate window appears (see Figure 72).

You can sort the list chronologically, by user, by location, or by status. Additionally, if you click on the label of one of the sort types, a pull-down menu is displayed so that you can further refine the sort. For example, if




Figure 72. Locate Window (status Program)

you click on the **by User** label, a list of user names appears. Select the user you want to put at the top of the list, and automatically the entries in the list change position according to your selection.


Using Sample Changers in Continuous Walkup Mode

You can use a sample changer throughout the day in a continuous walkup mode. This mode enables you to add samples to the queue throughout the day.

- 1. Prepare your sample and place it in an available location in the carousel or in a sample tray.
- 2. Enter walkup in the VNMR input window.

The Sample Entry Form window (see Figure 73) appears.

🗖 Sample Entry Form 🔽 🗖								
	Sample Number							
1 2 3 4	1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20							
21 22 23 24	21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40							
41 42 43 44	45 46 47 48 4	19 5	0				_	
			 Lise	r identif	icatir			
vnmr1								
			So	lvent Sel	ectior	1		
CDC13	DC13 D20 Benzene DMSO							
Acetone	Cyclohexane	Tolu	Jene	Methanol				
Experiment Selection								
H1			C13				F19	
P31			HC		H-gCOSY			
H-TOCSY			H-gHSQC		H-gCOSY-gHSQC			
H-gCOSY-HSQ	C-gHMBC		H-gCOSY-HSQC-gHMBC-HSQCTOXY		H-COSY-C-DEPT-HETCOR			
H-COSY-C-APT			H-TOCSY-NOESY		H-TOCSY-ROESY			
H-TOCSY-HMQ	H-TOCSY-HMQC							
Text								
Customize Parameters Add Sample Priority Sample Quit								
Number of samples submitted: O								



- 3. Fill in the Sample Entry Form window.
 - **Sample Number** select the location(s) to which you added samples. The sample locations that are grayed out are already in the day's queue.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen location(s).
 - **Experiment Selection** select one or more experiments for the chosen location(s).
 - Text enter information about the experiment, if desired.
 - To customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles, click the **Customize Parameters** button.

Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

• When finished with this location (or locations), click **Add Sample** or **Priority Sample**.

Add Sample – appends your requested samples and experiment to the current enterQ file.

Priority Sample – puts your requested samples and experiment at the head of the enterQ file. You can configure this button to request a password.

Quit – exits the program without adding samples to enterQ.

The walkup macro creates a new daily automation directory named auto_dd.mm.yy, where dd is the current day of the month, mm is the month, and yy is the year (e.g., auto_01.04.99). The automation directory for each day is saved in the directory specified by the global parameter globalauto. If the parameter and/or directory do not exist, the walkup macro creates them.

An alternative is to call the walkup macro from a menu button. A prototype asm menu with the buttons Walkup, Status, and Return has been put in menulib. The Status button is useful because it runs the status program to monitor the status of an automation run and it can retrieve data into VNMR for viewing. The asm menu can be made accessible from the Main menu by removing the comment symbols on the lines near the bottom of the menulib/main file.

The walkup macro also permits a fully continuous automation queue in which only the sample being run and any queued locations are inaccessible (grayed out). To enable this mode of operation, change the following line in the file /vnmr/asm/auto.conf:

```
From -
set exList {expl/sampleinfo locQ}
To -
set exList {expl/sampleinfo}
```

Adding Samples to an Automation Run in Progress

You can add samples or experiments to an automation run already in progress.

- 1. Prepare your sample and place it in an available location in the carousel or in a sample tray.
- 2. Type enter in the VNMR input window.

A Sample Entry Form window appears, similar to Figure 73.

- 3. Fill in the Sample Entry Form window.
 - **Sample Number** select the location(s) to which you added samples. The sample locations that are grayed out are already in the day's queue.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen location(s).
 - **Experiment Selection** select one or more experiments for the chosen location(s).
 - Text enter information about the experiment, if desired.
 - **Customize Parameters** Click on this button to customize any of the experimental parameters, including spectral width, number of scans, delays, and pulse angles. Click on the item you want to customize. A new window opens and provides appropriate parameters to change.

• When finished with this location (or locations), click Add Sample or Priority Sample.

Add Sample – appends your requested samples and experiment to the current enterQ file.

Priority Sample – puts your requested samples and experiment at the head of the enterQ file. You can configure this button to request a password.

• Quit – exits the program without adding samples to enterQ.

9.5 Changing Sample Changers or Serial Ports

To change to a different sample changer or to change to a different serial port, you must stop and restart acqproc by entering su acqproc twice:

```
>su acqproc
```

>su acqproc

9.6 Using Gradient Autoshimming with Automation

Gradient autoshimming requires the Automated Deuterium Gradient Shimming module to be installed.

Before you can run gradient shimming in automation, you must set up gradient shimming and map the shims with deuterium as described in the section on gradient autoshimming in the section "Deuterium Gradient Shimming," page 383.

After the shims are mapped, set wshim='g' in the parameter sets to be used with the automation run.

You might also wish to customize the macros setlk and gmapz for gradient shimming with different solvents. See "Example of Customizing a Macro," page 304, for more information.

9.7 Automation Run Description

Table 45 lists the commands and parameters related to automation run operation.

In an automation run, the user starts by using the enter program to provide information about the nature of up to 100 samples and the experiments to be performed on them.

After the information is entered and the automation run started, the user can run the status program to find out which experiment have been completed, which are in progress, and which are still in the queue. As part of an automation run, a single directory is created to handle the data collection. You can choose where the acquired data sets are to be distributed.

The enter program keeps track of which locations in the sample tray have been assigned, which are not assigned, and allows entry of information about each sample in a simple format. The enter program can be considered the "tray manager." The enter program requires the user to input information about the sample and the experiments to be run on a particular sample. If a method is already defined for performing a typical experiment, the name of the macro that describes that method (for example, hl or hc) can be entered or selected, and that experiment or experiments will be run.

Commands	
	A1 / 1.1/1 1/1
aa	Abort acquisition with error
auto<(automation_dir>	Prepare for an automation run
auto_au	Controlling macro for automation
<pre>autogo<(file<,automation_dir>)></pre>	Start automation run
<pre>autogo<(file<,automation_dir>)></pre>	Start automation run
autosa	Suspend current automation run
enter<(file<, <config_file>)></config_file>	Enter sample information for automation run (VNMR)
enter <file> <config_file></config_file></file>	Enter sample information for automation run (UNIX)
gilson	Open the Gilson Liquid Handler window
halt	Abort acquisition with no error
<pre>status<(directory<config_dir>)></config_dir></pre>	Display status of sample changer (VNMR)
status directory <config_dir></config_dir>	Display status of sample changer (UNIX)
walkup	Walkup automation
Parameter	
auto {'y','n'}	Automation mode active
autoname	Prefix for automation data file
gilpar {array of 4 integers}	VAST sample changer values
loc {0, 1 to traymax}	Location of sample in tray
saveglobal {array of parameter names}	Save selected parameters from global tree

Table 45. Automation Mode Commands and Parameter
--

Sometimes, however, a particular sample requires other conditions for data acquisition; this is not a problem because enter also allows a complete line of keyboard entry to describe the experiment to be run; for example, setup('H1','dmso') nt=4 pw=10 would be a legitimate entry. You can customize the selection of what information is requested and in what manner the information is requested. See the manual *VNMR User Programming* for details.

During the automation run, three files in the automation directory keep track of the status of different experiments:

- enterQ file holds the entries for the experiments waiting to be submitted.
- psgQ holds the experiment in the process of being submitted.
- doneQ file holds both active and completed experiments.

The status window automatically displays the contents of these files and enables the user to closely examine individual files. As the automation run proceeds, the status program updates its display to reflect the current status of the run.

An automation run uses its own experiment files, exp1 through exp4, which reside in the automation directory; thus the "automation" experiments are separate and distinct from any experiments in the files of individual users. After data is accumulated, you can choose where the collected data is stored. Data can be saved in the following locations:

- · In the automation directory
- · In subdirectories of the automation directory
- At any arbitrary directory

Basic Automation Run

An automation run involves three basic steps:

1. Enter information about the samples to be run during the sample changer run.

The enter program provides a mechanism to define essential information about the sample and the experiment to perform. Depending on the specific details of the enter program configuration, sample and experiment information might be selected by pressing a button or typing in information.

The sample entry process might typically occur during the course of the day (assuming a common case of open access operation during the day, and sample changer run at night). One or more operators can use the enter program to create a file of information about the samples and experiments they wish to run.

If multiple users are involved, the situation becomes slightly complicated because all the information must in the end reside in a single file. Either different users could enter information into a single file consecutively or they could enter information into different files simultaneously, with the system operator required to then merge these various files into a single file before the automation run is started. Since the file that is created is a simple text file, the latter task is easily accomplished with standard UNIX tools.

2. Create a directory for your "automation directory" to hold the various "queues," the automation experiments, and perhaps the data collected during the run.

To do this, use the auto program. This can be done before or after step 1, or even automatically when the automation run is started.

3. Enter autogo to start the run.

The autogo command starts the automation run, submitting the first experiment and holding subsequent experiments in a queue, ready to be submitted one at a time as the previous experiment terminates.

Note that the experiment queueing facility is not used, and that only one experiment is typically submitted at a time. This fact makes it easier to change the order of experiments in the queue once the automation run is in progress, to add experiments to the queue, etc.

Sometimes at the start of an automation run, there is a sample already in the magnet. Often, the system does not know the tray location to which that sample should be returned. The result is that the sample is returned to position zero in the sample tray. The sethw command provides a convenient method of identifying the location the sample should be returned to (e.g., sethw('loc', 22) causes the sample to be returned to location 22 in the sample tray).

Automation Behind the Scenes

An automation run consists of two main items: a file containing information on the samples to be run and a directory containing the data along with supporting automation files. Samples are entered using the enter program, which is specifically designed for automatic sample entry. You should enter the proper information. Once the information is entered, it is saved as a simple text file named by you. If you examine this file, you will see that the last item of each entry, Status, is set to queued.

An automation run does not begin until the command autogo is issued. You can enter autogo at the command line in VNMR to initiate automation (see the VNMR Command and Parameter Reference for further details about autogo).

Several things happen when autogo is issued. First, a directory is created to hold the automation data. The directory path is contained in the VNMR global variable autodir. Next, several subdirectories and files are created within this main data directory. Automation creates a new set of experiments (expl to exp4) to temporarily hold the data

during acquisition. The experiments reside in \$autodir. Other critical files created are enterQ, doneQ, psgQ, locQ, and sampleinfo.

To examine what the system is doing during an automation run, open a console window. This window echoes each event as it happens. Experiments to be performed are copied by autogo from the file created by the enter program to the enterQ file. As each sample is submitted for acquisition, the information about the sample is removed from the enterQ file and placed in the sampleinfo file. After submission, the information is written into the doneQ file.

The psgQ file holds information generated by a pulse sequence, e.g., S2PUL.psgQ is used by the automation system to start an acquisition.

The locQ file contains loc values that have been used during the automation run. This file is created by the auto_au macro. During walkup operation, the enter program uses the locQ file to determine which sample locations are available for another user.

The sampleinfo file contains information on the sample currently being run. This file is accessed by the acquisition macros to retrieve needed information. The sampleinfo file is placed in the appropriate experiment directory. The doneQ file contains information about *both* active and complete experiments.

The enterQ file is the critical file for sample management. Once an automation run is started, the only way to add more samples, delete samples, and so forth, is to directly modify the enterQ file. enterQ can be updated manually according to the procedures given later in this chapter. When the enterQ file contains no entries, the system recognizes that automation is complete.

While an Automation Run is in Progress

Automation mode is a separate mode of operation from the normal multiexperiment mode and must be kept distinct from it. An automation run can only be initiated when all other acquisitions are finished. Likewise, all data acquisition of an automation run must be finished (or at least paused) before the system can be returned to the multiexperiment mode of data acquisition. You should realize that these statements apply only to acquisition. Data processing is possible using any experiment. Data can be recalled, processed, and plotted in the normal way. Because the automation run acts in effect as a separate user, you will have no problems trying to plot at the same time as the automation run.

Also, by using the status command (see Figure 71), data that is part of the automation run and already finished (including 2D runs that are still in progress but for which some of the FIDs have already been acquired) can be recalled using the mouse and then processed or reprocessed. Of course, any processing that has been specified to take place as part of the automation run will occur at the appropriate time separately and concurrently with any processing you may be doing.

Interacting with the acquisition itself during an automation run is more limited than in multiexperiment mode. You cannot start an acquisition, but can abort the current acquisition using the Abort Acquisition button or by issuing the aa command. Depending on what you have used for the werr parameter, it is likely that an abort of an acquisition in progress terminates the acquisition with no processing whatsoever. Instead of entering aa, you can enter halt, in which case the normal end of experiment (wexp) processing occurs.

As the automation run progresses, experiments are removed from the enterQ file, placed the sampleinfo file, and then placed into the doneQ file as they are active or completed. All processing messages that occur during the automation run appear in the console window, and thus the progress of the automation run can be easily monitored by leaving the console window open. The status command can be used to monitor this progression.

When an Automation Run is Finished

At the finish of an automation run, the system automatically returns to the "normal" mode and normal data acquisition can be initiated with the go or au commands.

If you reach the end of the automation run and wish to continue the run with more samples, use the enter program to enter new samples into a new sample file (or perhaps you already did this while the automation run was in progress). Then enter autora.

The autoname parameter controls the file names to be used and can use the value of VNMR parameters as part of the file name. For a complete description of the autoname parameter, see the *Command and Parameter Reference* manual.

Parameters for Automation

Basic parameters for automation are recalled from the directory stdpar. This directory can be owned by the user who started autogo or it can reside in the /vnmr directory, in which case the directory is not alterable by anyone other than vnmr1.

stdpar should contain an entry for each nucleus. For example, the stdpar/H1.par file is retrieved by the macro h1 and used to run the proton NMR spectrum. To a significant extent, acquisition, processing, and plotting can be controlled by parameters preset in stdpar. Some parameters in stdpar that directly affect functioning of automation (not including general parameters such as pw) are listed in Table 46 with descriptions of values.

Certain global parameters are also important for the correct submission of additional acquisitions. The saveglobal global parameter saves an array of these global parameter names. Whenever a go, ga, or au is entered, the parameters listed in the array are copied to the current experiment parameter tree, and an underscore is appended to the parameter names (e.g., loc becomes loc_). Some of the global parameters saved are loc, lockpower, lockphase, lockgain, z0, lkof, gilpar, pkpick, and parstyle. Additional global parameter names can be added to the saveglobal array. The saveglobal parameter is also saved as saveglobal_, even though it is not an explicit member of the array.

Whenever any conditional processing occurs, the saveglobal_list of parameters is copied back into the global tree. By doing this, the processing macros have access to the values used to initiate the acquisition. Whenever a go, ga, or au is executed, any parameters listed in the experiment's saveglobal_parameter are deleted before the new parameters from the global saveglobal parameter are copied.

Variable Temperature Control During Automation

Variable temperature operation during an automation run proceeds in a straightforward way. The sequence of events is the following:

- 1. The current sample is removed from the magnet by the sample changer.
- 2. The new temperature is set and the temperature starts to change, at a rate of no more than 12 degrees per minute.
- 3. When the new temperature is reached, the system checks if the tin parameter is set to 'y'. If so, the system waits for the temperature controller to achieve regulation.

Parameter	Value	Description	
alock	's'	Autolock only on sample insert.	
	'Y'	Lock using parameters based on the solvent. Lock phase and lock gain are not adjusted (quick lock mode).	
	'u'	Do not lock (does not apply to MERCURY or GEMINI 2000).	
	'a'	Autolock searches for lock resonance and adjusts parameters lockpower and lockgain.	
wshim	's'	Shim using shim method specified by the parameter method.	
	'g'	Automatic shimming using gradient shimming is done only at beginning of the first experiment, following the sample change. method is ignored. This value is available only in automation; it is not used with go, ga, or au.	
method	'z1z2'	Shim gradients Z1 and Z2 with a criterion of medium to medium, total time 60 seconds or less (default method).	
gain	'n'	System adjusts gain before acquisition.	
spin	'n'	System does not regulate spinning before acquisition.	
	20	System tries to adjust spin speed to 20 Hz before starting acquisition. If speed does not regulate and in='y', ejects the sample and proceeds to the next sample. If in='w', only a warning message is added to the log file and acquisition continues with the unregulated speed.	
intmod	'off'	No integrals are plotted (typical for c13 parameter set).	
	'partial'	Integrals are plotted for each peak separately (typical for h1 parameter set). The region command determines integrals.	
pltmod	'off'	Suppress plotting. If an error occurs when processing (e.g., aph failure), the processing macro sets pltmod to off and no plotting occurs.	
	'fixed'	Plot spectrum using stdpar parameters. This allows setting all plotting-related parameters when stdpar is saved.	

Table 46.	Basic Parameters	s for Automation
Indic 10	Duble I ulumeters	, ioi riatomation

The VNMR parameter vtwait determines how long the system waits for the VT to regulate. The default is 300 seconds (5 minutes). If you are making a large temperature change, or have the VT controller slow rate set to very slow, you might need to change this default. If the VT does not regulate before this time and tin='y', the system will move to the next sample.

- 4. The sample changer inserts the sample into the magnet.
- 5. If spinning is requested, the sample spin regulation system starts.
- 6. If a pad delay is requested, that delay occurs. This delay allows the sample temperature to come to equilibrium.
- 7. The sequence then proceeds in the normal fashion, doing locking, shimming, and receiver gain adjustment, as requested.

9.8 Customizing the Sample Entry Window

The choices displayed by the enter program can be customized. The *users*, *solvents*, and *experiments* are specified in text files in the /vnmr/asm directory. Each text file is specified as pairs of lines of entries. The first line of each pair is the label that will be

displayed in the Sample Entry Form window by the enter program. The second line is what will be written into the text file generated by the enter program when that particular labeled item is selected.

Listing 3 shows the default contents of the experiments text file. The solvents file and user file are similarly constructed.

Listing 3. Contents of Default /vnmr/asm/experiments File

```
Н1
AuHexp
C13
AuCexp
F19
AuF
P31
AuP
HC
AuHexp(solvent, `CARBON`)
H-qCOSY
AuHexp(solvent,`gCOSY`)
H-TOCSY
AuHexp(solvent, `TOCSY`)
H-qHSOC
AuHexp(solvent,`gHSQC`)
H-gCOSY-gHSQC
AuHexp(solvent,`gCOSY`,`gHSQC`)
H-gCOSY-HSQC-gHMBC
AuHexp(solvent,`qCOSY`,`HSQC`,`qHMBC`)
H-qCOSY-HSOC-qHMBC-HSOCTOXY
AuHexp(solvent,`gCOSY`,`HSQC`,`gHMBC`,`HSQCTOXY`)
H-COSY-C-DEPT-HETCOR
AuHexp(solvent, `COSY`, `CARBON`, `DEPT`, `HETCOR`)
H-COSY-C-APT
AuHexp(solvent, COSY`, CARBON`, APT`)
H-TOCSY-NOESY
AuHexp(solvent, `TOCSY`, `NOESY`)
H-TOCSY-ROESY
AuHexp(solvent, `TOCSY`, `ROESY`)
H-TOCSY-HMOC
AuHexp(solvent, `TOCSY`, `HMQC`)
```

9.9 Automated Data Acquisition

As previously stated, when entering information about the sample, experiments are most easily begun by selecting the name of a macro from the enter program. Any user-created macro, with associated arguments (if any), can also be selected by the enter program, so that any particular experiment needed can be specifically programmed. A complete description of automation macros appears in the *VNMR Command and Parameter Reference*.

In the enter program, these macros are typically used as a simple name and without any of the optional arguments that may be used if the macro is used "by itself." The solvent,

which is required by each of these macros when run in the multi-experiment mode, is automatically added by the enter program.

Note that multiple entries in the enter file may refer to the same sample. Thus it is not necessary to use a "combined" experiment like hc to run a proton and carbon experiment on a particular sample. Instead, you can fill out one entry with hl and then a second for the same sample with cl3.

All of the standard single-experiment automation macros set the wexp parameter to 'procplot'. All of the standard combined-experiment automation macros set the wexp parameter to 'autolist'. None of the standard automation macros include a call to au (this is a change from earlier operation of these macros because they used to include the call to au). Therefore, customization in the enter program is possible by entering the standard macro followed by the changes (e.g., hl nt=4 selects the standard hl experiment but runs it with four transients).

During an automation run, the experiment information from the enterQ file is placed in expl and in the autodir directory. The auto_au macro is then called. This macro reads the sampleinfo file and sets the solvent and loc parameters, fills in the text, and executes the command defined by MACRO. After that, auto_au examines the value of the wexp parameter. If wexp is set to 'procplot', then auto_au calls au. If wexp is set to 'autolist', then auto_au inserts 'auto' as the first argument to autolist and calls au('wait'). If wexp is set to anything else, auto_au does not call au. As a result, any existing automation macros with a built-in call to au will still work as long as wexp is not set to 'procplot' or 'autolist'.

When au is executed, the global parameters listed in the saveglobal parameter are copied into the experiment parameter list. Each parameter name is appended with an underscore (_). Whenever any conditional processing occurs (e.g., wexp), the previously saved parameters are returned to the global tree. The conditional processing macros therefore have access to the correct values of global parameters, such as loc and the lock parameters.

Optimizing Acquisition Macros

Acquisition in automation is usually initiated by macros such as h1 and c13. By understanding the order of events in acquisition, automation can be optimized.

The macro first retrieves the appropriate stdpar file. Important parameters such as pw90 and tpwr are retrieved from the probes file. The .def files determine default sw, nt, etc. In addition, the parameters that directly affect automation (alock, wshim, etc.) should be set appropriately. These parameters are listed in Table 46 with explanations of reasonable values. Once the stdpar file has been retrieved, the macro calls setlk. After calling setlk, acquisition begins.

Listing 4 is the text of the AuHexp macro.

Data is processed by the macro procplot, which is called by autolist. The procplot macro processes and plots both 1D and 2D data. A full explanation of procplot and associated macros is provided in the *VNMR Command and Parameter Reference*.

Customizing Macro Operation

Extensive opportunity exists for customizing the operation of existing macros—this customization can always occur on either the systemwide level (in the /vnmr directory) or at the user level (in the user's personal vnmrsys directory). The starting parameter sets

Listing 4. Contents of AuHexp Macro

```
if (\$ = 0) then \$olv=\$1 else \$olv=solvent endif
setup('H1',$solv)
Autoclrexp
wexp=''
Autosetgpar
if (auto <> 'y') then
  explist='PROTON','glidewexp'
else
  explist='PROTON' setlk(solvent)
endif
if (\$ \# > 1) then
  $arg=''
  \$x = 2
  repeat
    format(\$x,0,0): \$arg
    if typeof('$'+$arg) then
      AutoAddEXP(${$x}, 'nodialog')
    endif
    $x=$x+1
  until $x > $#
endif
macro=$0 pw(45) seqfil='s2pul'
nt=8 clear(2) d1=1 wbs='' wnt=''
setsw(14,-2)
Automacrodir('make')
Autosetwexp
dg
```

can be modified to contain parameters of your choosing. And finally, the macros themselves (h1, c13, setlk, etc.) can be modified. An obvious recommendation is to make such changes at the single-user level first, and then determine their acceptability before making the changes at the system-wide level.

Other possibilities for customizing include macros used during an automation run:

- hlp, cl3p, fl9p, and p31p for processing.
- hregions to select integral regions in a proton spectrum.
- react for error processing.
- getsn to get a signal-to-noise estimate and testsn (which uses the testct macro and the sn parameter) for signal-to-noise checking.
- cleanexp to remove old experiment files and directories.
- cptmp to copy experiment data into an experiment subfile.
- setlk for setting up automatic locking choices (see example in next section).

Table 47 lists commands and parameters related to customizing macro operation.

A number of commands (au, halt, resume, etc.) are important in constructing macros for sample changer operation and are more or less unique to that mode in that they are rarely if ever used at other times. Note that jexp cannot be used in an automation macro. A

Commands	
au*	Submit experiment to acquisition and process data
auto_au	Controlling macro for automation
c13p	Process 1D carbon spectra
<pre>cleanexp<(file1<,file2,>)></pre>	Remove old files and directories from an experiment
cptmp<(file)>	Copy experiment data into experiment subfile
f19	Process 1D fluorine spectra
fixup	Adjust parameter values selected by setup macros
getsn:current_sn,predicted_sn	Get signal-to-noise estimate of a spectrum
hlp	Process 1D proton spectra
halt	Abort acquisition with no error
hregions	Select integral regions in proton spectrum
p31p	Process 1D phosphorus spectra
react<('wait')>	Recover from error conditions during werr processing
resume	Resume paused acquisition queue
setlk(solvent)	Set up lock parameters
testct	Check ct for resuming signal-to-noise testing
testsn	Test signal-to-noise of a spectrum
<pre>* au<(<'nocheck' ><,><'next'><</pre>	,><'wait'>)>
Parameters	
auto {'y','n'}	Automation mode active
autodir {string}	Automation directory absolute path
autoname {string}	Prefix for automation data file
lastlk {string}	Last lock solvent used
sn {number}	Signal-to-noise ratio

Table 47. Commands and Parameters for Customizing Macro Operation

number of parameters are also important in constructing macros for sample changer operation. These parameters are rarely if ever used at other times.

Example of Customizing a Macro

The macro setlk can be customized for each system to enhance autolocking and autoshimming. Because many of these enhancements are system specific, various lines in setlk are commented out and the system macro is nonfunctional until modified by the system administrator. As a result, it is advisable to explicitly set the values for the z0, lockpower, and lockgain parameters as discussed in item 3 below. Listing 5 provides an example of a modified setlk macro.

For each sample, setlk permits the system to perform the following adjustments:

 Retrieve solvent-based shim sets. Starting from a shim set for each solvent has proven useful on many automation systems. To retrieve a set of shims for each solvent, remove the quotes around the following line in setlk: rts(\$1):\$e

Alternatively, separate shims for each solvent may not be necessary but specific shims for only one or two solvents may be desirable. The following statements show how to implement this scheme:

```
format($1,'lower'):$solv
if ($solv='d2o')
  then rts('d2o')
  else rts('cdcl3')
endif
```

In this example, if the solvent is d20, then specific shims for d20 are called; otherwise, cdcl3 shims are used as a starting point for all organic solvents.

2. Choose a specific shim method based upon the T_1 of the individual solvents. Quotes should be removed from the following lines in setlk:

```
if ($solv='acetone')or($solv='cd3od')or($solv='cd2cl2')
then method='longt1'
else if ($solv='c6d6')or($solv='cdcl3')
then method='medt1'
else if ($solv='dmso')or($solv='d2o')
then method='shrtt1'
else method='zall'
endif
endif
endif
```

3. Explicitly set z0, lockpower, and gain. This is a useful technique to improve the autolock reliability. This is especially apparent for the first sample of an automation run. If the lock power has been low, this first sample may not be run because of an autolock failure.

To use this feature effectively, lock and shim on a sample in each one of the solvents that you will normally use. Record the values of z0, lockpower, and gain that produce a stable lock.

Remove the quotes from the following lines in setlk and input the appropriate values of z0, lockpower, and gain. Notice that the example below shows some of the entries with the z0 values inserted. The number represents the value of z0 as shown by the lock display in acqi. Replace the numbers with the appropriate values for your system.

```
if $solv='cdcl3'
 then Z0=128 lockpower=37 lockgain=41
  else if $solv='d2o'
   then Z0=35lockpower=37 lockgain=48
   else if $solv='acetone
      then Z0=12 lockpower=21 lockgain=39
      else if $solv='dmso'
        then lockpower=33 lockgain=36
        else if $solv='c6d6'
          then lockpower=20 lockgain=39
          else if $solv='cd3od'
            then lockpower=25 lockgain=35
            else lockpower=30 lockgain=40
          endif
        endif
      endif
    endif
endif
```

In addition, the macro gmapz can be customized to run gradient shimming on different solvents. Uncomment the appropriate lines in the following section of the macro:

```
"uncomment this section for automation (wshim='g')"
if (auto='y' and tn='lk') then
   $solv='' format(solvent,'lower'):$solv
"add other solvents with long T1"
   if ($solv='acetone') then
      d1=d1*2
```

```
endif
"add other solvents with weak signal"
    if ($solv='cdcl3') then
        nt=nt*4
    endif
endif
```

Listing 5 shows a more extensively customized setlk macro that handles multiple probes and provides solvent-based parameter sets, all based upon the system probe files.

Listing 5. Modified setlk Macro

```
exists('probe','parameter','global'):$ep
if ($ep > 0.5) and (probe <> '') then
 $ep = $ep
else
 $ep = 0
endif
if ($# < 1) then
 $solv = solvent
else
 $solv = $1
endif
"if solvent-based shim sets are desired for new solvents, store"
"appropriate shim sets (if you use more than one probe you should "
"have a shims library for each) then remove quotes from the next line "
format($solv,'lower'):$solv
if $solv = 'benzene' then $solv = 'c6d6' endif
if ($ep) then "probe exists & has a value "
 $shimfile = userdir+'/shims/'+probe+'/'+$solv
 exists($shimfile,'file'):$eshims
 if ($eshims > 0.5) then
   rts($shimfile)
 else
   $shimfile = '/vnmr/shims/'+probe+'/'+$solv
   exists($shimfile,'file'):$eshims
   if (\$eshims > 0.5) then
     rts($shimfile)
   else
     write('error','Sorry, no shims for this probe and solvent')
   endif
 endif
else
 rts($solv) "This happens if no probe-specific shim directories"
endif
"if shim method is to be determined by solvent, then remove the"
" quotes in next lines (you must create the methods given)"
if
($solv='acetone')or($solv='cd3od')or($solv='cd2cl2')or($solv='cd3cn')
 then method='longt1'
else if ($solv='c6d6')or($solv='cdcl3')or($solv='thf')
 then method='medt1'
else if ($solv='dmso')or($solv='d2o')
 then method='shrtt1'
else method='zall' endif endif endif
```

Listing 5.	Modified	setlk Macro	(continued)
------------	----------	-------------	-------------

```
"if you prefer to explicitly set lockpower/lockgain for each solvent,"
"the lines below must be updated to correspond to"
"the non-saturating values for the spectrometer in use. this"
"can speed up the automatic locking process."
"z0 can also be set (for weak locks), along with lockpower/lockgain"
if ($ep > 0.5) then
 if probe = 'asw305' then
   if $solv='cdcl3' then z0=415 lockpower=25 lockgain=32 else
   if $solv='d2o' then z0=536 lockpower=15 lockgain=29 else
   if $solv='acetone' then z0=630 lockpower=12 lockgain=24 else
   if $solv='dmso' then z0=629 lockpower=20 lockgain=26 else
   if $solv='c6d6'
                     then z0=420 lockpower=20 lockgain=24 else
   if $solv='cd3od' then z0=575 lockpower=24 lockgain=26 else
   lockpower=20 lockgain=25 endif endif endif endif endif
  else if probe = 'pfgid' then
   if $solv='cdcl3' then z0=413 lockpower=14 lockgain=29 else
   if $solv='d2o'
                     then z0=534 lockpower=12 lockgain=24 else
   if $solv='acetone' then z0=628 lockpower=6 lockgain=19 else
   if $solv='dmso' then z0=630 lockpower=11 lockgain=24 else
   if $solv='c6d6'
                     then z0=416 lockpower=10 lockgain=20 else
   if $solv='cd3od' then z0=573 lockpower=24 lockgain=24 else
   lockpower=24 lockgain=40 endif endif endif endif endif
  else if probe = 'sw231' then
    "values in this section are for the AutoSW probe "
   if $solv='cdcl3' then z0=452 lockpower=26 lockgain=29 else
   if $solv='d2o'
                      then z0=592 lockpower=15 lockgain=29 else
   if $solv='acetone' then z0=702 lockpower=12 lockgain=20 else
   if $solv='dmso' then z0=690 lockpower=20 lockgain=26 else
   if $solv='c6d6' then z0=452 lockpower=20 lockgain=24 else
   if $solv='cd3od' then z0=630 lockpower=24 lockgain=20 else
   lockpower=20 lockgain=25 endif endif endif endif endif
  else if probe = '10mm' then
   if $solv='cdcl3' then z0=414 lockpower=26 lockgain=46 else
   if $solv='d2o' then z0=533 lockpower=22 lockgain=42 else
   if $solv='acetone' then z0=696 lockpower=16 lockgain=30 else
   if $solv='dmso' then z0=603 lockpower=20 lockgain=42 else
   if $solv='c6d6'
                     then z0=415 lockpower=20 lockgain=40 else
   if $solv='cd3od' then z0=570 lockpower=24 lockgain=40 else
   lockpower=24 lockgain=40 endif endif endif endif endif
  endif endif endif "close probe-specific section"
else "probe variable does not exist or is not set"
  "values currently in this section are for AutoSW probe "
  if $solv='cdcl3' then z0=415 lockpower=26 lockgain=29 else
 if $solv='d2o' then z0=539 lockpower=15 lockgain=29 else
 if $solv='acetone' then z0=629 lockpower=10 lockgain=20 else
 if $solv='dmso' then z0=631 lockpower=20 lockgain=26 else
 if $solv='c6d6'
                    then z0=415 lockpower=20 lockgain=24 else
 if $solv='cd3od' then z0=570 lockpower=24 lockgain=20 else
 lockpower=20 lockgain=25 endif endif endif endif endif
endif "end of z0/power/gain section"
  " alock = 'y' "
```

9.10 Automated Data Processing

In an automation run, any processing done as part of the experiment is governed by four parameters: werr, wbs, wnt, and wexp.

werr contains the name of the macro that runs if an error is detected during the experiment (receiver overflow, lost lock, etc.).

- wbs contains the name of the macro that runs after each block of the acquisition, assuming that no other processing is occurring. wbs processing is not done, for example, if data processing is still going on for a previous sample. In automation mode, the basic use of wbs processing is for experiments in which it is desired to test the signal-to-noise ratio periodically in order the stop the experiment at a point when the desired ratio has been reached (discussed under the descriptions of cl3 and testsn in the VNMR Command and Parameter Reference).
- wnt applies to processing that must be done after each FID of a multi-FID experiment; therefore, wnt processing is rarely used in automation mode.
- wexp normally specifies a processing macro that performs data processing appropriate to the particular experiment that has been performed. This processing may include setting up and starting further experiments, based on the results of the first experiment. The most important processing to be done on each sample is selected by the wexp parameter.

If all four of these parameters are set to the null string (wexp=''), no data processing whatsoever is done, and the experiment is simply acquired and stored. Normally, however, you will want automatic data processing.

Whenever wbs, wnt, wexp, or werr processing occurs, the acquisition condition that initiated the processing is available from the parameter acqstatus. Many of the codes are related to sample changer operation. Refer to the manual *Getting Started* for more information about using acquisition codes. The codes are listed in the description of acqstatus in the *VNMR Command and Parameter Reference*.

9.11 File Structures in an Automation Run

This section describes the files and directories created in an automation run.

The command enter('abc') creates a directory named abc (referred to as the *automation setup directory*). The abc directory contains a file named abc, which lists experiment information, as follows:

```
varian> cat /export/home/vnmr1/abc
 SAMPLE#:
          2
    USER:
          vnmr1
   MACRO: h1
 SOLVENT: acetone
   TEXT:
          ethyl-vanillin
 USERDIR:
          /export/home/vnmr1/vnmrsvs
   DATA:
  STATUS:
          Queued
_____
 SAMPLE#:
          1
   USER: vnmr1
   MACRO:
          c13
 SOLVENT:
          cdcl3
    TEXT:
          menthol
```

USERDIR: /export/home/vnmr1/vnmrsys DATA: STATUS: Queued

Also in the automation directory (in this case, abc) is a directory named abc.macdir, which contains *GLIDE*-related information for an automation run, including the following files:

```
varian> ls /export/home/vnmrl/abc/abc.macdir
curloc loc1_AuHexp loc2_AuCexp
```

For the automation run, the autogo command uses the information from the *automation setup directory* created by the enter program (abc in our example) and places information about current and completed experiments in a new directory, called the *automation directory*. For example, if autogo('abc','abcl') is entered, information is copied from the abc directory into the abcl directory created by autogo. varian> ls /export/home/vnmr1/abcl

autoinfo enterQ enterQ.macdir doneQ exp1 exp2 exp3 exp4 psgQ sampleinfo

The experiment files (exp1, exp2, etc.) are used for data acquisition and data processing. The enterQ, psgQ, doneQ, and sampleinfo files contain information about the experiments to be run. The autoinfo file is currently unused. The information from the file (e.g., abc) in the automation setup directory is copied to enterQ. As each sample is submitted for acquisition, the information about that sample is copied to sampleinfo. After the experiment is finished, the information about completed sample is moved from sampleinfo to doneQ. Information generated by pulse sequences is stored in psgQ.

A sampleinfo file typically looks like the following, containing information about the sample that is just about to be run, for possible use by the macros that are being executed: varian> cat /export/home/vnmr1/abc1/exp1/sampleinfo

varian> cac	/ CAPOI C/ HOMC/ VHMILI/ ADCI/ CAPI/ BA
SAMPLE#:	1
USER:	vnmr1
MACRO:	c13
SOLVENT:	cdcl3
TEXT:	menthol
USERDIR:	/export/home/vnmr1/vnmrsys
DATA:	
STATUS:	Queued

The sampleinfo file is placed in the appropriate experiment directory and is saved with the data. After the automation run has run for a while, it might look like this:

varian> cd	/home/vnmr	1; ls abc	1		
0101.fid	autoinfo	enterQ	exp2	exp4	
0201.fid	doneQ	expl	exp3	3 psgQ	lab0401.fid

The new files, all containing the fid extension, are the data files collected during the acquisition run. Normally, as data is collected, the FIDs are placed into the automation directory with the name xxyy.fid, where xx is the sample location number and yy is an "experiment number" on that particular sample (1 through n, where n experiments have been run on the same sample). If the user who started autogo sets the parameter autoname, then autoname is interpreted to determine the directory and file name where the data will be stored. The autoname parameter controls the version number attached to the file name and uses the value of VNMR parameters as part of the file name.

The data files have the following internal structure:

varian> cd abc1; ls 0101.fid
fid log procpar sampleinfo text

Of these files, a key file that relates to automation is the log file, which maintains a log of everything that happens during the run, for you to examine at a later time.

```
varian> cat 0101.fid/log
======> Wed Dec 25 17:05:18 1996 <==== New Exper: '0101'
17:05:17
        Experiment Submitted.
17:05:18
          Experiment Started ID = 1
  FID
         Time
               Error
  ___
          ____
                  _____
17:09:06
         Experiment Completed ID = 1
The doneQ file might look like this at this time:
varian> cat doneQ
  SAMPLE#:
           2
    USER: vnmr1
    MACRO: h1
  SOLVENT: acetone
    TEXT: ethyl-vanillin
  USERDIR: /home/vnmr1/vnmrsys
     DATA:
          /home/vnmr1/slptest/0201
   STATUS: Complete
SAMPLE#: 1
    USER: vnmr1
    MACRO: c13
  SOLVENT: cdcl3
    TEXT: menthol
  USERDIR: /home/vnmr1/vnmrsys
     DATA: /home/vnmr1/slptest/0101
   STATUS: Error
_____
```

Another key file is /vnmr/acqqueue/acquisitioninfo, which contains information about the state of the acquisition computer after the last sample. This file is a text file that might look like this:

```
varian> cat /vnmr/acqqueue/acquisitioninfo
nt 16 ct 16 scale 0 np = 29952
gain 36 lockphase 256 lockpower 22 lockgain 53
```

Finally, the file /vnmr/acqqueue/lastlk contains information about the last lock solvent:

varian> cat /vnmr/acqqueue/lastlk
last lock solvent was cdcl3

Chapter 10. VAST Accessory Operation

Sections in this chapter:

- 10.1 "Using the VAST Accessory," page 311
- 10.2 "Solvent Suppression in VAST," page 317
- 10.3 "Processing, Displaying, and Plotting VAST Data Sets," page 323
- 10.4 "Using CombiPlate to Analyze Data," page 331
- 10.5 "Vast Process, Display, and Plot Macros," page 337
- 10.6 "Preparing the Hardware and Configuring VNMR," page 341
- 10.7 "Calibrating Volumes and Flow Rates," page 345
- 10.8 "Acquiring Data on Standard Test Samples," page 354
- 10.9 "Evaluating Carryover," page 355
- 10.10 "VAST Interface Description," page 355
- 10.11 "Customizing the enter Window for VAST," page 362
- 10.12 "Files that Control VAST Operation," page 363
- 10.13 "Writing VAST Protocols," page 363

The VAST (versatile automatic sample transport) accessory uses the VNMR automation capability of UNITY *INOVA* and *MERCURY-VX* NMR spectrometers to take a sample from a vial or well on the Liquid Handler (see Figure 74) and move it into the Microflow probe. When the sample is in place, NMR data can be acquired. When the NMR experiment is finished, the sample is removed and either discarded into a waste container or returned to the source container. After removing the sample, the liquid path is then washed to minimize cross-contamination or carryover between samples.

The VAST accessory comprises the following:

- Gilson 215 Liquid Handler, Rheodyne (Gilson 819) injector and Valco valves
- VAST software interface to VNMR (sold separately, requires 6.1B or later)
- PEEK tubing, PEEK connectors, and signal cables
- Varian Microflow probe (sold separately)

10.1 Using the VAST Accessory

The VAST accessory is highly mobile and can be easily moved between spectrometers or moved away from the spectrometer when not in use. Therefore, a certain amount of setup might be required each time the VAST accessory is to be used.

After setting up, most of your interaction with the VAST autosampler will occur with normal VNMR automation controls, such as the enter program and with the loc parameter in the VNMR input window.



Figure 74. Gilson 215 Liquid Handler

To Prepare VAST for Use

- 1. Make sure the VAST accessory hardware is set up and that VNMR is configured, as described in "Preparing the Hardware and Configuring VNMR," page 341.
- 2. Calibrate the Z position of the arm (if necessary). Z position calibration is described in "Calibrations," page 361.
- 3. Make sure an appropriate waste container is properly installed on the rinsing station.
- 4. Prepare an appropriate solvent and make sure the inlet line to the syringe assembly is inserted into your solvent reservoir. Be sure to have enough solvent to complete the automation run.
- **CAUTION:** Make sure the sample/solvent you are using is compatible with, and does not react with, the sample container. For example, chloroform solvents will rapidly dissolve polystyrene microtiter plates. Microtiter plates can be made of polystyrene, polypropylene, polycarbonate, glass, and many other materials.
 - 5. In the VNMR input window, enter **gilson** to open the Liquid Handler window shown in Figure 75.

-	Gilson Liqui	id Handler	•				
SAMPLE Def. Rad	ck Def. Main Contro	Calibrations					
Select parameters Default 1 🛁							
Probe Volume:	400.0	Probe Slow Vol:	100.0				
Sample Volume:	350.0	Probe Slow Rate:	0.5				
Push Volume:	50.0	Probe Fast Rate:	2.0				
Keep Sample:	yes no	Sample Extra Vol:	25.0				
Number of Rinses:	1	Rinse Extra Vol:	25.0				
Rinse Delta Vol:	10.0	Additional Paramet	ers				
Save Now							

Figure 75. Gilson Liquid Handler Window

Field	Default Value*
SAMPLE Def. pane (Figure 103)	
Probe Volume	400
Sample Volume	350
Push Volume	50
Keep Sample	no
Number of Rinses	0
Rinse Delta Vol	10.0
Sample Well Rate	4.0
Sample Height	20
Sample Depth	NOSEEK
Needle Rinse Volume	1000 (for a 1000 µL syringe)
Needle Rinse Rate	16
Probe Slow Vol	100
Probe Slow Rate	0.3
Probe Fast Rate	0.3
Sample Extra Vol	25
Rinse Extra Vol	25
Mix Volume	0
Mix Time	0
Mix Flow Rate	0
Mix Height	0
Rack Def. pane (Figure 104)	
Settings for racks 1 to 5 must set, severe damage to the need	reflect the racks that exist on the system. If incorrectly ille (probe) can result.
Main Control pane (Figure 105)	
Air Valve:	Off
Plunger:	Up
Syringe Volume:	Volume of the syringe on the Liquids Handler
Arm Z Scale (mm):	Set at installation

Check the default settings in the Liquid Handler window. The default values are listed in Table 48.

Table 48. Typical Default Values in Liquid Handler Windows

* These default values are designed for the most viscous solvents and will work for D_2O .

- 6. Click the **Rack Def.** tab and configure the window for the racks to be used.
- 7. Click on the Main Control tab and then click on Return Home.

This may take a few minutes. To speed up the process, press the STOP button on the front of the Liquid Handler, lift the needle, and slide the arm toward the left.

No values at this time.

8. Prime the pump by clicking on the **Prime pump** button in the **Main Control** pane.

Calibrations pane (Figure 106)

The needle should automatically move over to the rinse station before the priming starts. Make sure to prime until you no longer see bubbles in the inlet line to the syringe module.

The priming routine runs for about 3 minutes and uses the solvent attached to the inlet line.

9. Clear the NMR probe of any leftover sample or solvent by setting **Air valve** to **ON**.

Let this run for about two minutes or until you hear air hissing out of the inject port. Be prepared to wipe up any sample or solvent that seeps out of the inject port.

- 10. Calibrate **Probe Volume**, as described in "Calibrating Volumes and Flow Rates," page 345.
- Acquire data on standard test samples, as described in "Acquiring Data on Standard Test Samples," page 354.

This verifies that the VAST hardware and software are operational.

- 12. If desired, calibrate Sample Volume and the flow rate parameters as described in "Calibrating Volumes and Flow Rates," page 345.
- 13. Prepare your samples and transfer them to the sample wells in a sample rack. Be sure to note where each sample is located in the sample tray.

The VAST accessory is now ready for use with the VNMR automation software. Refer to the following procedures in this section:

- "To Set Up NMR Experiments for VAST," page 314
- "To Change Samples with VNMR," page 314
- "To Shut Down a VAST System," page 317

To Set Up NMR Experiments for VAST

If solvent suppression is not needed, only routine parameters such as tpwr, pw90, d1, and sw need to be optimized. The usual commands (go, ga, au) can be used to start acquisition.

If solvent suppression is needed, use WET and Scout, as described in section 10.2 "Solvent Suppression in VAST," page 317. Processing, displaying, and plotting of data obtained with VAST is described in section 10.3 "Processing, Displaying, and Plotting VAST Data Sets," page 323.

To Change Samples with VNMR

You can change samples manually or automatically.

Manually

Sample location is defined with a number (loc=3) or with a string (vloc='G7'), using the loc or vloc parameters.

To change samples and not run acquisition:

• In the current zone and rack—for example, location 3—enter:

loc=3 change

You can also set vloc to a value like 'A3' or 'G7' if you know the microtiter-plate address.

• In a different zone or rack—for example rack 1, zone 1, location 3—enter: loc=3 vzone=1 vrack=1 change

To change samples and run acquisition (ga, go, au, or su):

- In the current zone and rack—for example, location 3—enter: loc=3 ga
- In a different zone or rack—for example rack 1, zone 1, location 3—enter: loc=3 vzone=1 vrack=1 ga or

```
vloc='G7' ga
```

Automatically with the enter Program

To change samples automatically with the enter program, use the following steps.

1. Type **enter** in the VNMR input window.

The enter program prompts you for a file name that stores the information set by the enter program. You can also use the command enter('file').

The Sample Entry Form window appears, similar to Figure 76.

— Sample Entry Form 🗾 👘	
Rack Number 1 2 3 4 5	
Zone Number	
2 one Number	
Sample Number	
H1 G1 F1 E1 D1 C1 B1 A1	
H2 G2 F2 E2 D2 C2 B2 A2	
H3 G3 F3 E3 D3 C3 B3 A3	
H4 G4 F4 E4 D4 C4 B4 A4	
H5 G5 F5 E5 D5 C5 B5 A5	
H6 G6 F6 E6 D6 C6 B6 A6	
H7 G7 F7 E7 D7 C7 B7 A7	
H8 G8 F8 E8 D8 C8 B8 A8	
H9 G9 F9 E9 D9 C9 B9 A9	
H10 G10 F10 E10 D10 C10 B10 A10	
H11 G11 F11 E11 D11 C11 B11 A11	
H12 G12 F12 E12 D12 C12 B12 A12	
User identification	
vnmr1	
Solvent Selection	
None D20 CH3CN CH3OH	
DMSO CDC13 CH2C12	
Experiment Selection	
Autoscout Autowettocsy AutoCOSY	
AutoPFG-COSY AutoDQCOSY AutoPFG-MQCOSYps	
AutoNOESY AutoROESY AutoPFG-HMQC	
AutoPFG-HMQCps AutoPFG-HSQC	
Selected Experiment autoscout	
Text	
Add Entry Exit and Save Quit	If you did not configure the racks you will see this message
Racks have not been configured. Use gilson	
Number of samples submitted: 0	Click on Quit and configure
	the racks.

Figure 76. VAST Sample Entry Form Window (enter program)

- 2. Select the desired samples in the window.
- 3. Click Add Entry.
- 4. Click Exit and Save.
- 5. Enter autogo

You are prompted for the file set up by the enter program and for the name of a directory in which the automation data is to be stored.

Enter the file name from step 1 and the name of a directory. You can also include the file and directory in the command by entering

```
autogo('file','directory').
```

The automation run begins.

The autogo command by itself uses a different set of parameters than the **Start Run** button on the LC-NMR pane, which is described below.

Automatically from the LCNMR/STARS Pane

A preferred alternative to using the enter program is to use the LCNMR/STARS pane (shown in Figure 77), which has the buttons **Sample Entry** and **Start Run**. This pane is automatically converted to a VAST Experiment Setup pane when you initialize VAST using the menu buttons or by typing VAST1D on the command line. Unlike other sample changer software, these buttons use the current parameter set as the starting conditions for the subsequent Autoscout acquisition, which is described in "To Set Up NMR Experiments for VAST," page 314.

Multiple Frequency Solv				VAST Experiment Setup	Text entry	Acq & Obs
WET Center(ppm) Solv 1 1.98	0.0	Width 80.0 Auto	o peak find	Number of transients Steady state transients		Decouplers
■ Solv 2 4.14	1082.2	80.0 Aut	o peak find	Reference Frequency	1.95 ppm	Sequence
□ Solv 3 2.60	309.0		et Solv 3	WET 13C Dec Freq	-11640	Flags & Cond.
□ Solv 4 7.54	2783.64 Hz	80.0 Se	et Solv 4	Start of Plot -0.00 Width of Plot 8.99	Autoshim n Equil. Delay 5	Process
				Plot Mode Fixed	Integrals Off	Process2
📕 C13WET		T	rial WET			Display
📕 Composite observe pul	Composite observe pulse					
Solvent subtraction	One	All		Sample	Entry	display2
Acqiuisition		Processing		Directory VAST_data_1a		LCNMR/STARS
Acquisition time	1.82048	📕 Line Broadeni	ing 0.5 🔺 💌			Text
Repetition rate (sec)	1.93	👅 Complex Fouri	ier No 32k 🔺 🛨	Start	Run	
Spectral Width	9000.9					Spare
Transmitter offset	-1258.3	Process		Stat	us	Setup EXP

Figure 77. LC-NMR Pane for VAST

The LC-NMR pane makes it easier to use automation.

- 1. In VNMR click Main Menu > Setup > Flow-NMR > Initialize VAST.
- 2. Click the LCNMR/STARS tab.
- 3. Click the Sample Entry button.

A Sample Entry window pops up.

- 4. Fill in the Sample Entry Form window.
 - Peak Number select one or more peaks on which to run experiments.
 - User Identification select the user.
 - Solvent Selection select the solvent for the chosen peak(s).
 - Experiment Selection select Autoscout, Autowettocsy, or both for the chosen peaks.
 - Text enter information about the experiment, if desired.
 - When finished with this peak (or peaks), click Add Entry.

When you are done adding entries, click Exit and Save.

The sample specified by the Queue name file is run, and the data is stored in the specified directory.

After the automation run is started, you can click the Status button to bring up a status window.

Using the LC-NMR pane replaces the commands enter, autogo, status. Unlike the autogo command, however, the Start Run button saves the current parameter set (including solvent suppression parameters) and uses this for all subsequent acquisitions called for by the autoscout macro.

To Shut Down a VAST System

The following procedure to shut down the VAST system consists of: setting the loc parameter to the appropriate value and then flushing the sample out of the probe with the air valve.

- 1. Set the loc parameter to 0 or 'n' as follows:
 - To extract a sample from the NMR probe after a manual run (needle is in the injector port):

loc=0 change

- To disable the Gilson Liquid Handler so samples are not changed inadvertently: loc='n'
- 2. In the VNMR input window, enter gilson to open the Liquid Handler window.
- 3. Click on the Main Control tab in the Liquid Handler window.
- 4. Clear the NMR probe of any leftover sample or solvent by setting Air valve to ON.

Let this run for about two minutes or until you hear air hissing out of the injector port. Be prepared to wipe up any sample or solvent that seeps out of the injector port.

10.2 Solvent Suppression in VAST

Solvent suppression is invariably the first step needed for any VAST experiment. After the initial solvent suppression setup is done, however, you only need to click the Trial WET button in the VAST pane to perform further optimizations.

This section contains the following:

- "Setting Up Solvent Suppression," page 317
- "Troubleshooting Solvent Suppression," page 322
- "Evaluating Solvent Mixture Equilibration," page 322
- "Solvent Suppression: Background Information," page 322

Setting Up Solvent Suppression

- Retrieve a parameter set appropriate for VAST by using either the menu system (Setup -> Flow-NMR -> Initialize VAST) or through the rtp command.
- 2. Regulate the probe temperature (preferably to a value close to ambient). The value of temp is displayed in the Acq & Obs pane.
- 3. Turn the lock off. Set lockgain and lockpower to zero.

Turning off lock eliminates *lock pull* during each spectrum of a solvent gradient experiment and *lock jumps*. *Lock pull* is caused by the constantly changing HOD frequency during the signal averaging. *Lock jumps* occur if the concentration of the deuterated species drops below the detection threshold (see the Noise reject parameter).

The sudden loss of lock, with its commensurate sudden change in field, is annoying during Trial WET setup operations. Even more critically, it can disrupt solvent suppression during a VAST run. Although turning the lock off is not mandatory, data is measurably better unless you are running a long experiment.

- 4. Tune the probe and adjust shimming if necessary.
- 5. Check if the global parameters ref_pw90 and ref_pwr exist:

```
ref_pw90?
```

```
ref_pwr?
```

or click on the Sequence tab and look at the VAST sequence pane, see Figure 78.

	Multiple	Frequency	Solvent Suppr	ession						Acq & Obs
	Active	Search	👅 WET	👅 Filter	Center(ppm)	Offset	WET Width		WET (no peak find)	
	👅 Solv 1	F	I	H	1.95	0.0	80.0	Set		Decouplers
	■ Solv 2	F	۳ 2	F	2.77	411.171	80.0	Set	Trial WET	Sequence
	📕 Solv 3	F	X 3		2.57	309.021	80.0	Set	ITIAI WEI	Flags & Cond.
	👅 Solv 4		□ 4	F	7.52	2783.64	80.0	Set		
	📕 Solv 5		F 5	μ.	3.95	1000.0	80.0	Set	process	Process
	🗆 Solv 6		□ 6		1.95	0.0	0.0	Set		Process2
	🗆 Solv 7		□ 7		1.95	0.0	0.0			Display
	WET RF Puls	ses (Calc.)	WET Gradien	t Pulses		Hz	Hz	VAST Inform	mation	
	WET pulse	21557.7	Strength	24000	WET Decoupl	ing		Location:	0	display2
	Power	-2	Time	0.002	👅 C13WET			(Location:) AO	LCNMR/STARS
rof pw00	WET shape	pbx	Recovery	0.002	Dec Freq	-11640.		Zone :	1	Text
ref_pw90	Reference I	Power 46	Null delay	0.002	Power	35		Rack:	1	
ref_pwr —	Reference p	pw90 16.5			Mod Freq	750		Protocol:		Spare
-	📕 Composit	e Shape comp	_3g		ON/OFF stat	us nnn		Noise reje	ct 2	Setup EXP

Figure 78. VAST Sequence Pane

6. If ref_pw90 and ref_pwr do not exist, create them:

```
create('ref_pw90','pulse')
create('ref_pwr','real')
```

- 7. If they do exist, make sure they are set correctly for the probe being used (e.g., ref_pwr=44 and ref_pw90=25). The ref_pwr used must be within the linear portion of the amplifier's operational range. These values are used to calculate the power levels for the shaped pulses. Generally, pw90's determined at power levels that are 12 to 20 db below the level used for the observe pw90 are within the linear operating range of the amplifier.
- 8. Click on the LCNMR/STARS tab to open the VAST pane see Figure 79.
- 9. Click on WET and select: Solv 1, Solv 2, Solv 3, or Solv 4 to set the number of NMR signals (up to 4 signals using this pane see step 13 for suppression of up to 7 signals) you intend to suppress.

The first NMR signal to be suppressed must be suppressed using Solv 1. If the spectrum contains two or more NMR signals that are to be suppressed, select Solv 2, Solv 3, or Solv 4 until the appropriate number of signals have been chosen. Solv 2 does not have to be selected to use Solv 3, or Solv 4. A more detailed discussion on how to determine the which signal suppression option to select when 2 or more signals are to be suppressed is given in step 13.

Multiple Frequency Solv	ent Suppres	sion	VAST Experiment Setup	Text entry	Acq & Obs
Center(ppm)	Offset		Number of transients	32	
■ Solv 1 1.95	0.0	Auto peak find	Steady state transients	2	Decouplers
👅 Solv 2 2.77	411.2	Auto peak find	Reference Frequency	1.95 ppm	Sequence
_ Solv 3 2.57	309.0	Set Solv 3	WET 13C Dec Freq	-11640	Flags & Cond.
_ Solv 4 7.52	2783.64	Set Solv 4	Start of Plot -7.05	Autoshim n	
	Hz		Width of Plot 18.01	Equil. Delay 5	Process
👅 WET			Plot Mode Fixed	Integrals Off	Process2
C13WET		Trial WET			Display
👅 Composite observe pu	lse		Queue name queue_1a		1. 1. 0
Solvent subtraction			Sample H	intry	display2
Acquiisition		Processing	Directory VAST_data_1a		LCNMR/STARS
Acquisition time	1.82047	📕 Line Broadening 🛛 0.5 🔺 💌			Text
Repetition rate (sec)	1.93	👅 Complex Fourier No 32k 🔺 👻	Start 1	Run	-
Spectral Width	9000.9				Spare
Transmitter offset	-150.99(Process	Statu	IS	Setup EXP

Figure 79. VAST Pane

The value of the Center and Offset for Solv 1 is determined automatically by searching over a wide band of frequencies.

- If Solv 2 is not selected, the tallest peak in the spectrum is suppressed.
- If Solv 2 is selected, the high field signal is assigned to Solv 1 (this may or may not be the tallest signal in the spectrum) and then the region of the spectrum that is to the low field side of Solv 1 is searched for the tallest signal in that region.

Center and Offset are determined automatically and displayed. These values are non changeable on this panel. These values can, however, be entered in the corresponding Sequence pane.

For Solv 3, or Solv 4 you must supply either a frequency for Center (in ppm) or a value for the offset from the carrier for Offset (in Hz). These values can be entered manually, or they can be extracted from a spectrum by placing a cursor at the appropriate location in the displayed spectrum and pushing the appropriate Set Solv button. When either Solv 3, or Solv 4 or both are selected, a narrow range of frequencies, about the set frequency is searched. The frequency range of the search is discussed in step 13.

Three additional NMR signals can be suppressed for a total of 7. The Sequence pane provides control of the suppression of these additional signals - see step 13.

10. Select the C13WET box to suppress the ¹³C satellites of the solvent resonances.

This assumes the values for Dec Freq, Power, and Mod Freq (located in the Sequence pane) have been properly set. The ON/OFF status is typically nnn, although nny may be useful if the size of the residual 13C satellites are a problem.

- 11. Click on the Acq & Obs tab to open the Acq & Obs pane see Figure 80.
- 12. The **obs pulse** and **power** should be set to value that produce pw90 that is as short as practical. In general, the shorter the better, although 3.0 or 4.4 microsecond pw90s are particularly useful. see step 14.
- 13. VAST Multiple Frequency Suppression Options.
 - a. Click on the Sequence tab to open the Sequence pane for suppression of more than four signals.Select up to 7 signals using the LCNMR/STARS pane. One (must be Solv 1) or two (must be Solv 1 and Solv 2) frequencies can be searched for automatically. Five additional frequencies can be searched for over a narrow range based on the value of the center frequency. In the Sequence pane, up to 7 signals can selected. Solv 1 must be active and search selected if two or more signals will be suppressed.

Acquisition			Observe C	hannel		Samp	le		Acq & Obs
Spect width	8999.89	Hz	Nucleus	HI		Date	February 10,	2000	Decouplers
Acq. time	1.82047	sec	Spect Freq.	500.618	MHz	File	lc1d		Doorderoro
Acquired complex pts	16384		offset	-1270.7	Hz	Solvent	hone		Sequence
Recycle delay	0.001	sec	obs pulse	3.0	us 🛁	temperature	20	<u> </u>	Flags & Cond.
Transients	8		power	63					
Steady state	2								Process
			Calibration	:		Channel	: Nucl.	Status	Process2
			pi/2 pulse	: 3	us	Observe	: H1		Display
			power	63		Decoupler	: 013	- nnn	
						Decoupler2	:	- nnn	display2
						Decoupler3	:	- mm	LCNMR/STARS
						Gradient Ty			Text
						Grautent Ty	de: um		Spare
									Setup EXP

Figure 80. VAST Acq & Obs Pane

Active	Search	WET	Filter	Center(ppm)	Offset	WET Width		WET (no peak find)	Acq & Obs
		I 1		1.95	0.0	80.0	Set	the (no point sind)	Decouplers
	F		2 🔳	2.77	411.171	80.0	Set	1	Sequence
Solv 3	F	X 3		2.57	309.021	80.0	Set	Trial WET	Flags & Con
Solv 4	F		4 🔳	7.52	2783.64	80.0	Set		
Solv 5		F 5	.	3.95	1000.0	80.0	Set	process	Process
⊥ Solv 6			6 🗆	1.95	0.0	0.0	Set		Process2
Solv 7		□ 7		1.95	0.0	0.0			Display
WET RF Puls	es (Calc.)	WET Gradie	ent Pulses		Hz	Hz	VAST Inform	ation	display2
WET pulse	21557.7	Strength	24000	WET Decoupli	ing		Location:	0	ursprayz
Power	-2	Time	0.002	📕 C13WET			(Location:)	A0	LCNMR/STAF
WET shape	pbx	Recovery	0.002	Dec Freq	-11640.		Zone :	1	Text
Reference Po	ower 46	Null delay	7 0.002	Power	35		Rack:	1	
Reference p	90 16.5	1		Mod Freq	750		Protocol:		Spare
Composite	Shape comp_	3g		ON/OFF statu	ıs nnn		Noise rejec	t 2	Setup EXE

Composite shape Noise reject/

Figure 81. VAST Sequence Pane and Frequency Suppression Options

Both Solv 1 and Solv 2 (see step 9 for description of Solv 1 and Solv 2 frequency search) use algorithms for finding the signal frequency that are different from the algorithm used by Solv 3 through Solv 7 to locate the signal frequency. If the search option is active for any one of four signals, Solv 3 through Solv 7, the search is restricted to an approximate range of + or - 0.2 ppm about the value displayed in the Sequence pane for the Center frequency. The Center frequency can be either entered manually or set by placing the cursor over the signal and pressing any one of the SET buttons. The value of Offset is calculated automatically.

- b. Specify **acquisition and processing options** for each signal. Each **Solv** selection, see Figure 81, has four choices:
- Search for a signal.
- Suppress a signal using WET.
- Apply post acquisition signal suppression filtration.
- Specify the Wet Width of the signal suppression.

Each choice may be applied individually or in any combination. The choices are only active if the Active button for a given Solv is RED, see Figure 81. Currently the software limits post acquisition filtration to a maximum of 4 signals (more then four then filter buttons to may be set to active in the tcl

pane). If five or more filter buttons are active, only the first four active Solv with filtration selected will be active.

Wet Width is typically set to 80 Hz (for each solvent frequency; this corresponds to about a 20 msec pulse). You can adjust these entries to larger or smaller values to generate a larger or smaller bandwidths of suppression for each signal (as desired). This value of Width determines the duration of the WET pulse, which influences the power of the WET pulse (both parameters being displayed in the Sequence pane).

Pressing the **Wet** (**no peak find**) button will execute scan using the parameters as they appear in the Sequence window. WET suppression and or post acquisition Filtering will only be applied to those signals that are Active and have either or both of these options selected.

Pressing the **Trial WET** button will execute a search for each active signal with search enabled and apply all other options that are enabled for each active solvent.

- 14. Select Composite observe pulse in the LCNMR/STARS pane (Figure 79) if you can achieve a pw90 of 3.0 or 4.4 μs. This will produce better looking data. Comp shape in the Sequence pane (Figure 81) must also be properly set (for pw=3.0 use comp_3pos and for pw=4.4 use comp_44g).
- 15. In the **Sequence** pane, the parameters under WET Gradient Pulses are typically set as follows:

Label	Value	Parameter Equivalent
Strength	24000	gzlvw
Time	0.002	gtw
Recovery	0.002	gswet
Null delay	0.001	dz

Because the solvent resonances are often left upside down in relation to the peaks of interest, the Null delay can be optimized to a value between 0 and 0.005 to catch the signals as they pass through zero (in analogy to the WEFT experiment).

16. Set Noise reject to 3 in the Sequence pane, see Figure 81.

The Noise reject parameter affects the determination of the frequencies for solvents 2 through 7, but not solvent 1. Set Noise reject to a higher or lower number as appropriate:

- To maximize the reliability of peak-finding for large solvent resonances, set Noise reject to a larger number (5 to 10).
- To detect smaller resonances, set Noise reject to a smaller number (1 to 3).

Noise reject does not accept values of less than 1.

- If Noise reject is set too small, the Scout Scan[™] incorrectly detects an incorrect frequency to suppress (e.g., a noise spike or a tall sample resonance).
- If Noise reject is set too large, the Scout Scan[™] does not detect the peak you intend to suppress (in which case Scout Scan[™] sets the frequency to the last value used, which might not be accurate for the current sample), although this is a useful way to suppress a fixed frequency or a very tiny signal.
- 17. Click on the **Trial WET** button in the **VAST** pane, Figure 79, to perform further optimizations. This starts the Scout Scan process. Two spectra are produced:

- The first spectrum is a one-transient no-solvent-suppression 1D spectrum (full spectral width). The NMR software spends a few seconds analyzing this data, creating the appropriate shapes, and setting up parameters.
- The second spectrum, a signal-averaged data set, is acquired using the parameter optimization that was done after the initial acquisition completed. The data is automatically Fourier transformed and displayed as the second spectrum.

Troubleshooting Solvent Suppression

If the suppression does not look as good as it should, the cause might be one of the following:

- PFG gradients are not properly connected, turned on, or configured (check the values of gradtype and pfgon).
- Probe is not well rinsed out.
- Shims are not optimized. Note that the shims can change significantly between acetonitrile and methanol solvent mixtures.
- Parameters ref_pwr and ref_pw90 are incorrect for the probe currently in use.

Evaluating Solvent Mixture Equilibration

The NMR chemical shifts of the solvent signals change as a function of mobile phase composition. You should always allow the pumped LC solvent mixture to equilibrate before starting an NMR run. One way to evaluate this is by doing the following:

- 1. In the VAST pane, press the Trial WET button.
- 2. After Trial WET finishes, check the value of transmitter offset as well as the **Offset** values for any selected **Solv 2**, **Solv 3**, and **Solv 4** entries.
- 3. Repeat 1 and 2 until these values remain constant. The system is equilibrated when these numbers do not change after you run Trial Wet.

The offsets for Solv 1 and Solv 2 are usually determined by the Scout Scan. They are not directly changeable by the user in the VAST pane, but manual entry is possible in the Sequence pane.

Solvent Suppression: Background Information

For VAST, the VAST pane and the Sequence pane are of the most interest:

• Virtually all routine VAST operations are available from the VAST pane.

WET Experiments

Using a modified form of the WET (Water Eliminated through Transverse gradients) experiments, the peaks of (dis)interest are excited with a selective 90° pulse, using a WET waveform that has been convoluted with one or more frequencies (SLP). This pulse is then followed with a gradient to dephase the undesired signals, and the process is repeated multiple times. The WET shape has been found to be a suitable choice for VAST applications.

WET Shapes

The WET shape used (set to pbx in the Sequence pane) is calculated on-the-fly by the Pbox software package, which is now resident within VNMR. This shape is recalculated many times during an VAST experiment. For this reason, the single shape name pbx is used by the pulse sequence although the shape is recalculated for each use.

While a number of basic shapes can be used to apply selective pulses to solvent lines, there are always trade-offs. Longer pulses are more selective, but require more time (which is of the essence, particularly in flowing systems). WET (based on SEDUCE-1) is a shape taken from the literature that has been found to be very effective in pulsed field gradient solvent suppression experiments.

One of the two shaped pulses comp_3g.RF or comp_44g.RF is used when the Composite observe pulse box in the VAST pane is selected. The user can define which of these is used in the Comp shape field in the Sequence pane.

Related Parameters

Solvent suppression is achieved not only by the pulse sequence, but also (in part) by software processing. Software processing is most frequently accomplished by using wftlc after the parameters ssfilter and ssntaps (and perhaps ssorder) have been set to reasonable values (typically, ssfilter and solvent suppression:LC-NMR ssfilter is set to about 80, ssntaps is set to about 251, and ssorder is not used).

Another parameter of interest is sslsfrq. This parameter shifts the location of the filter notch (in Hz from tof). Often, the quality of solvent suppression can be improved somewhat by setting sslsfrq to a small positive number (e.g., sslsfrq=3). This parameter can also be arrayed to produce multiple filter notches within the same spectrum. Selection of which signal to which the notch filter is applied is made in the **Sequence** pane. Any four of the available seven signals may be selected. NOTE: the filter will be applied only to the first four signals if more then four are selected.

10.3 Processing, Displaying, and Plotting VAST Data Sets

Spectra obtained during a VAST automation run can be processed and displayed in different ways:

- Individually, using standard VNMR 1D processing.
- As a pseudo 2D data set.
- As an array of 1D spectra.
- Using Combiplate

Processing the spectra as a pseudo 2D data set or as an array of 1D spectra facilitates the inspection of all the spectra at the same time.

Combiplate provides a graphical representation of the samples based predefined conditions.

Creating a Pseudo 2D Data Set

The individual spectra from a VAST automation run can be merged into a pseudo 2D data set. Two related display macros are used to create the pseudo 2D data set:

- vastglue
- vastglue2

The criteria for selecting which macro to use is determined by the value of the parameter autoname.

- If autoname=' ' (*the default and preferred value*) at the beginning of an automation run, construct the pseudo 2D data set with vastglue.
- If autoname='<user defined>' at the beginning of the automation run, construct the pseudo 2D data set with vastglue2. The special case, autoname='<user defined>', is not as user friendly as the autoname=''. The macro vastglue2 must be hard coded to accommodate the user defined name.

vastglue

The macro vastglue accepts either of two pairs of arguments or no argument. The two pairs of arguments for vastglue are: rack, zone and glue order, plate.

You must provide two arguments or no argument. The default, rack=1 zone=1, is specified as vastglue(1,1). Where the argument rack, zone, or the default is used, the glue order is determined from the doneQ file. This is the *default and preferred* gluing option. Part of a doneQ file is shown Figure 82.

The order in which each spectrum was acquired is specified by the entry on the line DATA. In this example, the micro titer plate was defined with each micro titer well having a number (1 through n). The first two digits are the well number and the second two digits specify the experiment. In an automation run, the first spectrum is the *autoscout* spectrum and the second spectrum is the optimized *wet solvent suppressed* spectrum. The macro "glues" together the even numbered spectra to create the pseudo 2D data set that can be saved using svf('filename').

The doneQ file, Figure 82, is written by the automation macro into the same directory as the spectra. Before starting a VAST automation run two directories were specified, one in the Queue name field and one in the Scout directory field, Figure 83. The Scout directory field specifies where the spectra doneQ, and other files are to be stored. If you are moving data sets you should move this entire directory, not just the spectra. In the other directory are a Queue file and directory used by automation software during the automation run.

SAMPLE#: 1
RACK: 1
ZONE: 1
USER: vnmr1
MACRO: autowet
SOLVENT: D2O
TEXT: A Micro titer plate mix ACDMT
USERDIR: /export/home/vnmr1/vnmrsys
DATA: 0101
STATUS: Complete
·
-
SAMPLE#: 1
RACK: 1
ZONE: 1
USER: vnmr1
MACRO: autowet
SOLVENT: D2O
TEXT: A Micro titer plate mix ACDMT
USERDIR: /export/home/vnmr1/vnmrsys
DATA: 0102
STATUS: Complete
_
SAMPLE#: 2
RACK: 1
ZONE: 1
USER: vnmr1
MACRO: autowet
SOLVENT: D2O
TEXT: A Micro titer plate mix ACDMT
USERDIR: /export/home/vnmr1/vnmrsys
DATA: 0201
STATUS: Complete
SAMPLE#: 2
RACK: 1
ZONE: 1
USER: vnmr1
MACRO: autowet
SOLVENT: D2O
TEXT: A Micro titer plate mix ACDMT
USERDIR: /export/home/vnmr1/vnmrsys
DATA: 0202
STATUS: Complete
SAMPLE#: 3

Figure 82. Partial doneQ File

Sample E Scout directory five_sam			LCNMR/STARS
Queue name queue_1a	r		display2
Reference Frequency	2.44		Display
Steady state transients	2	F	Process2
Number of transients	8		Process

Figure 83. Queue Name and Scout Directory Fields

The glue order, plate option allows the glue order to be defined by the user when the data set is glued rather than using the doneQ order. This method will produce a pseudo 2D data that can be displayed with as a contour map, color map, or stacked plot. The association of the well number with the spectrum is lost and renders the display macro, dsvast, less useful (although still functional). The plat_glue program (see "Defining a Custom Display Order with plate_glue" on page 327) allows the user to redefine the order in which the spectra in a pseudo 2D data set are displayed. Since this is a post-gluing operation, the association between the well and the spectrum can always be determined.

vastglue2

The vastglue2 macro artificially reconstructs a 2D data set from a series of 1D data sets having similar filenames. It is crucial to ensure that the format of the file names of each of the 1D data sets is identical. The vastglue2 macro reads each 1D file in succession, and adds it to the previous data, but in a 2D format. It assumes that the wells will be named using a predefined format, for example:

autoname='filename_R%RACK:%_Z%ZONE:%_S%SAMPLE#:%_'

This definition must be hard coded in the macro by the user. A different version of the macro must be written for each variation of the autoname format.

Creating a pseudo 2D DATA Set Using vastglue

This procedure creates a pseudo 2D data set for rack 1 zone 1. See "vastglue" on page 324 for other rack and zone combinations.

- 1. Start VNMR.
- 2. Join any experiment other than exp5
- 3. Change directories to the directory containing the VAST experiment FIDs.
- 4. Type vastglue.

If the data was acquired other then in rack1 zone 1, recall the first spectrum and then type vastglue(rack, zone).

The gluing process takes place in exp5. Once the process is finished the glued data set is returned to the original experiment.

- 5. The pseudo 2D data set can be saved using the standard command, svf('filename').
- 6. If the data set was acquired with fn>65k, set fn to a value \leq 65 k before processing the data set.

There are several glue order templates provided with VNMR, the default is to use the done Q if no glue order is specified. Glue templates are found in ~/vnmr/user_templates/glue.

Processing, Displaying, and Plotting Glued VAST Data

Process the glued data similarly to a spectrum from a 2D experiment.

 Set the display and display orientation, enter: f full trace='f2'

```
2. Set the 'f1' axis label:
    setwell
    setvalue('sw1',1,'processed')
    wf1=100 dconi
```

This example is for a 96 well data set. The setwell macro changes the label of the f1 axis to "well number" and setvalue('swl', 1, 'processed') labels each trace with the well number corresponding to that sample's spectrum. The number of traces displayed is set using wf1=100 dcon.

3. Enter pmode='full'

If pmode='full' is set before wftlda, the 2D spectrum can be phased as needed and redisplayed using dconi without transforming all the spectra.

- 4. Enter **wftlda** to process the glued data.
- 5. Display the data set as a color map, a contour plot, or as a stacked plot.

The standard 2D display options, vs2d, trace, projection, etc. all apply. For more flexible data presentation, use the processing macros dsvast, dsvast2d, and vastget and the plotting macros plvast, plvast2d and plvastget.

- To generate a dss stacked display, enter vastget.
- To display on spectrum, enter vastget('well#').
- To display selected spectra, enter:

vastget('well#','well#')

where 'well#', 'well#' is an arbitrary list of one or more sample wells. The names of the wells are defined as A1 through H12 and are the names attached to each spectra by dsvast or plvast. The spectra are displayed or plotted as a dss stacked display.

- To display spectra 1 through arraydim, use the dsvast macro.
- 6. Plot the data in the same way you would plot a true 2D spectrum.

For more flexible data presentation, use the plotting macros plvast, plvast2d and plvastget.

- Plot an individual spectrum with pl.
- Plot a series of spectra with fixed intervals or all the spectra using:

```
plww(<start,finish,step> or <'all'>)
```

• To plot spectra 1 through arraydim, use the plvast macro.

VAST display macros, dsvast and vastget, and their corresponding plotting macros, plvast and plvastget, present or extract information from the wftlda or wft processed pseudo 2D data. The macro dsvast (or plvast) default to display (or plot) spectra 1 through arraydim (the number of spectra in the pseudo 2D data set) and labels each spectrum with the name of the sample well. A preset naming convention is used: rows A through H and columns 1 through 8. The spectra are labeled A1 through H8. These macros accept one or two arguments to determine the format of the display (dsvast) or plot (plvast):

dsvast<(display order),(number of columns displayed)>

plvast<(display order),(number of columns displayed)>

- If no arguments are given the default will display will 8 columns of spectra.
- If the arraydim of glueorderarray is ≤ 8 , one row of spectra is displayed.
- If dsvast or plvast is provided with one argument, the argument is defined as the number of columns to be displayed. If two arguments are provided, the first argument, display order, must be the name of the file containing the glue order and the second argument is the number of columns displayed.

Defining a Custom Display Order with plate_glue

To define a display order that is different from the order defined in the doneQ file, use the plate_glue program.

1. Start plate_glue from a terminal window by entering:

plate_glue &

The plate glue window appears, similar to Figure 84.

2. Select the wells you want to display.

A well turns gray when selected. As you select more wells, a red line joins the grayed out wells.Wells do not have to be adjacent to each other.

- To select a well, click on it.
- To deselect a well, right-click on a selected (grayed) well.
- To select a row, partial row, column, partial column, or diagonal grouping of wells, click and hold the left mouse button on the first well, drag to the last well in the selection, and release the mouse button.




Figure 84. plate_glue Window and e4x4 glue order.

3. After the display order is defined, click the save button and save your custom glue order file in

the ~/vnmrsys/templates/glue directory with a file name of your choice. In this example the glue order file name is e4x4.

The display and plotting macros, dsvast and plvast, look in ~/vnmrsys/templates/glue for the glue order file.

Examples of Plots of a VAST Data Set

The plots in this section represent the data from a 96-well titer plate in zone 1, rack 1. Figure 85 shows an array of 96 1D spectra plotted using plvast with no arguments.




Figure 85. Array of 1D Spectra from a VAST Data Set

Figure 86 shows a plot of 16 spectra from the lower right corner. The plate_glue program was used to create a glue file called e4x4 and the spectra were plotted using plvast('e4x4', 4). This produced a 4x4 matrix using the glue file e4x4 and specifying the number of columns as 4.




Figure 86. Plot of a Subset from a 96-Well Sample Plate




Figure 87 shows a pseudo 2D contour plot of the VAST data set.

Figure 87. Contour Plot of a VAST Data Set

The pseudo 2D in Figure 87 was created using vastglue with no arguments. The fl axis of the plot is modified to display the well numbers by typing the following on the VNMR command line:

```
setwell
```

```
setvalue('sw1',1,'processed')
wf1=100 dconi
```

Pseudo 2D data sets processed with wftlda can be displayed and plotted using the same tools that are used for a typical 2D data set. The plot in Figure 87 was created using pcon('pos',8).

Figure 88 shows the same data plotted as a stacked spectra using the pl2d macro.

Figure 89 shows the spectrum of the sample in well F6, extracted with vastget('F6') and plotted with plvastget('F6') pap pscale page.

Summary of VAST Display and Plot Options

To display or plot VAST microtiter plate data, you have the following options:

- Stack of 96 spectra in the output bin of the plotter. This is the normal output from an automation run.
- 8 x 12 matrix display of the full spectrum of each compound (similar to Figure 85). This display gives you the feeling of the entire plate at once, including perhaps approximate reaction yields, but little specific structural information. Use the vastglue, dsvast, and plvast macros.
- 8 x 12 matrix display of one specific spectral expansion of each compound, showing an interesting region (e.g., 6.5 to 9 ppm or 0.5 to 3 ppm). This display gives you more specific information about the chemical functionality and diversity available in the plate. Use the vastglue, dsvast, and plvast macros.




Figure 88. Stacked Plot of a VAST Data Set




Figure 89. Spectrum of the Sample in Well F6

- 8 x 12 matrix display of one specific spectral expansion, where the expansion only includes a reference (control) resonance (e.g., TMS). This display allows you to calibrate concentration or to verify the quality of each injection. Use the vastglue, dsvast, and plvast macros.
- 8 x 12 matrix display of the integral of one specific spectral expansion of each compound, showing an interesting region (e.g., reference resonance, 6.5 to 9 ppm, or 0.5 to 3 ppm). This display can give a better picture of quantitation, especially if the

peak amplitudes are affected by differing NMR resolution for each well. Use the vastglue, dsvast, and plvast macros.

- To show a contour plot of the entire plate (similar to Figure 87), enter trace='f2' wftlda. You can look at each spectrum using the *trace* display to scan through each *increment* as desired. You can expand around either a given chemical shift range, or a certain number of wells (a row or column, depending upon the orientation of the plate). Use the vastglue, dsvast, and plvast macros.
- To view region intensities by colors and color densities, use the CombiPlate window, which is described in the next section.

10.4 Using CombiPlate to Analyze Data

CombiPlate allows you to view integral region intensities by colors and color densities. The CombiPlate window provides a map of the microtiter plate, which allows you to click on individual sample wells to bring up a spectrum in the VNMR display window.

Preparing VNMR Data For Analysis Using CombiPlate

The information displayed by *CombiPlate* is obtained from a pseudo 2D data set that was previously created and saved.

1. Create the pseudo 2D data set.

See "Creating a Pseudo 2D Data Set" on page 323 for instructions on creating a pseudo 2D data set. The data set should be saved with all the spectra properly phased and referenced.

- 2. Change to the directory containing the VAST pseudo 2D data set to be used to construct the *CombiPlate* display.
- 3. The load the pseudo 2D spectra into the current experiment and process the spectra (*wft dsvast* is the most convenient processing and display method). Use any experiment except experiment 5. Experiment 5 is used by various macros during the processing steps. The pseudo 2D set will supply the data to be displayed by CombiPlate.
- 4. Start CombiPlate from the VNMR command line by typing: combi_preproc.
- 5. The macro will process the data in exp5 and display exp5 when it starts the processing. When the processing is completed, a spectrum representing the sum of all spectra is displayed in exp5. Use this spectrum is used to set the integral regions. Adjust the phase if necessary. Only the regions for the integration need to be set at this time so perfect phasing is not necessary. Adjust the amplitude of the integral as needed. Tilt and Level do not need to be critically adjusted.
- 6. Define up to three integral regions using the sum of spectra spectrum created in step 4.
- 7. Join the original experiment and move the parameters from exp5 to the original experiment.
- 8. Type mp(5, 'original exp')
- 9. Type *dsvast* to display all the spectra then, using *vastget*('*well*') (each well is defined [A->H][1->8] e.g. F8), select one having the widest range of resonances.

- 10. Carefully **phase** this spectrum, **check** the **integral Lvl/Tlt**. If necessary correct the integral Lvl/Tlt to give proper integration. DO NOT change the integral regions.
- 11. Type nlivast(number of wells analyzed). The number of wells analyzed must be an even number. The file ASCII text file, integ.out, will be created in the current experiment's directory (~/vnmrsys/exp(#)/integ.out). CombiPlate uses the data in integ.out to create the display. If you intend to continue the analysis at a later time and you will be running CombiPlate from a terminal window, copy this file to the directory containing the VAST data. Renaming the file, although not required, is suggested. Each time nlivast is run it creates another integ.out file. Any previous copies of this file in the current active directory are overwritten.


Data Analysis Using CombiPlate And VNMR

Complete the section, "Preparing VNMR Data For Analysis Using CombiPlate," page 331 before continuing with the instructions here.

1. On the VNMR command line type:

combishow(1,2,3).

This assigns region (1) as red, region (2) as green, and region (3) as blue and starts the CombiPlate display using the data in ($\sim/vnmrsys/exp(#)/integ.out$). The initial colors displayed are from a synthetic data included with CombiPlate and may not be correct, see Figure 90 and "Checking And Fixing The Color Map," page 336. Correcting the colors is *not* necessary at this time.



CombiPlate displays synthetic data on startup.The color map, as shown here, may not be correct. This is normal.

The integral regions defined ' prior to running dlivast

Color intensity control for each integral region


Figure 90. CombiPlate Started by combishow From vnmr


- 2. Click the OK button on the Welcome! window to close when you finish reading its contents.
- 3. **Click** the **Display** button to load the data into the CombiPlate window. Depending upon the values of the integrals you may or may not have any colors displayed.

- If no colors are displayed set all the *log(100% color intensity)* slider bars to (log(100%)= 2.000) and **Click** the **Display** button.
- 4. Set the color intensity for each region.

The integral regions defined in step 6 of "Preparing VNMR Data For Analysis Using CombiPlate," page 331 are displayed in the scrolling text box at the at the bottom of the CombiPlate window and in each color field. Color intensity for each region is controlled independently. Each region has two preset values, 0% and 100% and a variable. If the variable option is selected the intensity of the color can be related to the integral area.

- If the *log(100% color intensity)* equals the largest integral value the color intensity is related to the integral area.
- If *log(100% color intensity)* is set to a value equal to or less then the smallest positive non zero integral value for that region the intensity is not related to integral area. The color is binary, on or off. In this case the color is turned on (a non zero positive integral) or it is off (integral is either zero of negative).
- a. Set all the log(100% color intensity) slider bars to (log(100%)=-2.000). This will set up a binary condition. The display in Figure 91 is set up to turn on the color of a given region only if there is a positive integral for that region.





- b. Click on the Variable radio button.
- c. Click the Display button.
- d. To use the intensity of the color as a measure of the integral area of a region, move the slider bar for that regions to the right. Click on the Display button. If all the color for that region diappears, you moved the bar to far move it back to the left. Select a value that provides a display with both the smallest integral value and largest integral value represented. Repeat for each region as desired. The display does not automatically update as the slider bar is moved. You must Click on the Display button to view all changes.

- 5. Individual spectra can be examined by selecting and clicking on the cell of interest. In Figure 91 cell c6 has been selected. The cell is white and the label "Show values" is displayed. The spectrum and inset spectral regions are displayed in the vnmr window. This feature is active only when vnmr is active and CombiPlate is started from combi_show on the vnmr command line.
- 6. Save the CombiPlate display by clicking the **Save** button.

In the CombiPlate save file window, Figure 92, you specify the directory and file name and which region or regions the file is to contain. A data set that was analyzed with all three regions active can be saved in three

SAVE DATA
Current directory: /export/home/vnmr1
New Directory:
Enter Filename: combidata
📕 Save red 📕 Save green 📕 Save blue
OK Cancel

Figure 92. CombiPlate Save File Window

independent files, on for each region by selecting one region. You must select one or more colors. If you do not select any colors a message "No colors selected" is displayed and no further actions takes place.

Analyzing Data Using CombiPlate Without VNMR.

When CombiPlate is started from a terminal window, clicking on a cell will display the integral information for that well in the scroll box at the lower right corner of the CombiPlate window. No spectra will be displayed, even if VNMR is running.

• Complete the instructions in "Preparing VNMR Data For Analysis Using CombiPlate," page 331 before continuing.

CombiPlate loads a synthetic data set that will give it an evenly graded color scheme running from green in the top left to yellow (top right), light blue (bottom left) and white (bottom right). If you do not see a smooth, even gradation of color, there are other programs running that are using colors requested by CombiPlate. This makes the colors unavailable to CombiPlate. See "Checking And Fixing The Color Map," page 336.

Analyzing VNMR Generated DATA

1. From a shell tool or terminal window UNIX, enter:

combiplate&



- Click on the Load button. An Input screen window appears with two input file options, From VNMR and Reload. Click on From VNMR to open the File Browser.
- 3. **Click** on the VNMR file containing the output from the macro *nlivast* such as *integ.out*. **Click** the **open** button to load the file. The *Update* button updates the browser.



The next *Input screen* opens, see Figure
 93. Before the selected input file can be loaded, information on the number of row and the starting and ending columns must be provided.

The default values assume that the microtiter plate is in a portrait orientation and the cell in the upper left is H1 with columns labeled H->A and rows labeled 1->12 (top to bottom). If the number of rows and columns in the data set are different (fewer then 8 columns or 12 rows) enter




Figure 93. CombiPlate Column and Row Input Screen

the correct values in the appropriate fields. Remember to press the return key after each entry or the change will not take affect.

- 5. Click the OK button to load the file and the row and column data into CombiPlate.
- 6. The File is loaded. The integral regions defined in step 6 of "Preparing VNMR Data For Analysis Using CombiPlate," page 331 are displayed in the scrolling text box at the at the bottom of the CombiPlate window in the scroll box, see Figure 90.
- 7. **Enter** the number of the integral region in the field entry box that you want associated with this color, see Figure 94. CombiPlate *does not* automatically assign




Figure 94. CombiPlate Field Box

a region to a color. Each integral region or field must be associated with a color. You must press the ENTER key after a value is placed in each field.

- 8. **Click** the **Display** button to load the data into the CombiPlate window. Depending upon the values of the integrals you may or may not have any colors displayed.
 - If no colors are displayed set all the *log(100% color intensity)* slider bars to (log(100%)= 2.000) and **Click** the **Display** button.
- 9. Set the color intensity for each region.

Color intensity for each region is controlled independently. Each region has two preset values, 0% and 100% and a variable. If the variable option is selected the intensity of the color can be related to the integral area.

- If the *log(100% color intensity)* equals the largest integral value the color intensity is related to the integral area.
- If *log(100% color intensity)* is set to a value equal to or less then the smallest positive non zero integral value for that region the intensity is not related to integral area. The color is binary, on or off. In this case the color is turned on (a non zero positive integral) or it is off (integral is either zero of negative).
- a. Set all the log(100% color intensity) slider bars to (log(100%)=-2.000). This will set up a binary condition. The display in Figure 91 is set up to turn on the color of a given region only if there is a positive integral for that region.
- b. Click on the Variable radio button.
- c. Click the Display button.

- d. To use the intensity of the color as a measure of the integral area of a region, move the slider bar for that regions to the right. **Click** on the **Display** button. If all the color for that region diappears, you moved the bar to far move it back to the left. Select a value that provides a display with both the smallest integral value and largest integral value represented. Repeat for each region as desired. The display does not automatically update as the slider bar is moved. You must **Click** on the **Display** button to view all changes.
- 10. Click on any cell to view the integration data.

Analyzing CombiPlate DATA

- 1. From a shell tool or terminal window UNIX, enter: combiplate
- Click on the Load button. An Input screen window appears with two input file options, From VNMR and Reload. Click on From VNMR to open the File Browser.
- 3. **Click** on a CombiPlate display file. **Click** the **open** button to load the file. The *Update* button updates the browser.
- 4. The next *Input screen* opens, see Figure 93. Click the OK button to load the file.
- 5. Continue by following the instructions in "Analyzing VNMR Generated DATA," page 334 from step 6 to the end.

Checking And Fixing The Color Map

When you start CombiPlate (see Figure 95) the plate that is displayed is filled with synthetic data. This will give it an evenly graded color scheme running from green in the top left to yellow (top right), light blue (bottom left) and white (bottom right).

							C	CombiPlate
								√ ∿ 0% Field selected for Red: 0
h1	g1	f1	e1	d1	c1	b1	a1	◆ 100% Norichia log(100%) has a value of
h2	g2	f2	e2	d2	c2	b2	a2	◆ Variable 109(100%) has a value of 1.200
h3	g3	fЗ	e3	dЗ	c3	b3	a3	-2.000 0.000 2.000 4.000 6.000
h4	q4	f4	e4	d4	c4	b4	a4	♦ 0% Field selected for Green: 0
	q5	f5	e5	d5	c5	b5	a5	◆ Variable log(100%) has a value of 1.200
	ļ							
h6	g6	f6	e6	d6	C6	b6	a6	-2.000 0.000 2.000 4.000 6.000
h7	g7	f7	e7	d7	с7	b7	a7	O% Field selected for Blue: 0 O
h8	g8	f8	e8	d8	c8	b8	a8	◆ Variable log(100%) has a value of 1.200
h9	g9	f9	e9	d9	c9	b9	a9	-2.000 0.000 2.000 4.000 6.000
h10	g10	f10	e10	d10	c10	b10	a10	Current cell
h11	q11	f11	e11	d11	c11	b11	a11	Weight – Red Green Blue Blue Raw Data
h12	q12	f12	e12	d12	c12	b12	a12	
n12	y12	112	elZ	ulz	CIZ	012	alz	
Loa	ad	Sa	we		Displa	xy	Clear	r Help Exit

Figure 95. CombiPlate Window for VAST Data Analysis

If the color gradation is *not smooth*, you have run out of color map entries. The color coding of the results may be displayed incorrectly in the *CombiPlate* window. A warning will be issued whenever *CombiPlate* thinks that it has run out of colors. This warning is in red in the text window.

- 1. Fixing the color map.
- 2. Close any applications that may be color intensive such as:
 - Acrobat Reader
 - Netscape
 - Frame Maker
 - VNMR

Applications running on the SUN computer all use colors form the color map and these colors are not available to CombiPlate. As a result, displays dependent upon gradations in color may be incorrect if a color is not available.

- 3. Exit *CombiPlate*
- 4. Restart CombiPlate.

The desired color scheme should now be displayed.

5. If VNMR is running and you run out of colors, you may wish to save results inspection later when VNMR is not active.

10.5 Vast Process, Display, and Plot Macros

combiplate	View a color map for visual analysis of VAST microtiter plate (U)
Syntax:	(From UNIX) combiplate
Description:	Opens the CombiPlate window, which provides a map of microtiter plate, allowing data to be viewed from individual sample wells. The window enables viewing integral region intensities by colors and color densities.
combishow	Display regions as red, green, and blue in CombiPlate window (M)
Syntax:	combishow(r,g,b)
Description:	Displays integral regions shown on the spectrum as red (r) , green (g) , and blue (b) in the CombiPlate window. CombiPlate reads the regions automatically. 1, 2, or 3 integral regions can be designated. At least one integral region must be specified. Combishow displays spectra associated with individual wells.
dsvast	Display VAST data in a stacked 1D-NMR matrix format (M)
Applicability:	Systems with the VAST accessory.
Syntax:	dsvast<(display order,number of columns displayed)>
Description:	dsvast will arrange and display the traces from a reconstructed 2D data set (see (see vastglue) as an array of 1D spectra in a matrix of 1D spectra. If no arguments are provided, the number of rows and columns will be determined by the periodicity of the display order based on the doneQ. For example, if a block of 96 spectra (typical for a microtiter-plate) have been acquired using VAST

automation, the spectra will be displayed in a matrix 8 rows and 12 columns with the well label using the format $(A \rightarrow H)(1 \rightarrow 12)$.

The spectra can be plotted using the macro plvast.

Arguments: display order is optional and its default value is the glue order as listed in glueorderarray. A display order can be defined using the plate_glue program.

number of columns displayed. The default value of is deduced by examining the periodicity of the requested display order. The number of columns displayed can entered as the second argument or as the first argument if the default display order is used.

dsvast2d Display VAST data in a pseudo-2D format (M)

Applicability: Systems with the VAST accessory.

Syntax: dsvast2d(number)

- Description: If an array of 1D spectra have been acquired (in particular if a block of 96 spectra has been acquired using VAST automation, especially in a microtiterplate format), and if these spectra have been glued into a reconstructed 2D data set (see vastglue), this macro will arrange and display them (on the screen) in a convenient pseudo-2D format (almost like an LC-NMR chromatogram).Well labels are not attached to the spectra and spectra are plotted with 8 spectra per row.
- Arguments: The default is to display all the spectra (from 1 through arraydim) with 8 columns (spectra) and 12 rows. An optional argument dsvast2d(number) allows one to specify that only spectra from *1* through *number* should be plotted. The number of spectra displayed is rounded up to the nearest multiple of 8.

nlivast Produces a text file of integral regions without a sum region (M)

Applicability:	Systems with	VAST accessory.
----------------	--------------	-----------------

Syntax: nlivast(last)

- Description: Using predefined integral regions from the spectra for each well, nlivast writes a text file, integ.out, containing the integrals of the regions. The file is written into the current experiment. Does not add an additional region that is the sum of all the defined regions for each well (see dlivast).
- Arguments: last is the number of the last well. The default is 96.

nlivast2 Produces a text file with normalized integral regions (M)

Applicability: Systems with VAST accessory.

Syntax: nlivast(well#)

- Description: Using predefined integral regions from the spectra for each well, nlivast2 writes a text file, integ.out, containing the integrals of the regions. The file is written into the current experiment. Integrals are normalized to the integral specified by the argument well. The macro nlivast2 does not add an additional region that is the sum of all the defined regions for each well (see dlivast). All of the spectra are integrated.
- Arguments: well is the number of the reference sample well. The default reference is well 96.

nlivast3	Produces a text file with normalized integral regions (M)
Applicability:	Systems with VAST accessory.
Syntax:	nlivast(well#)
Description:	Using predefined integral regions from the spectra for each well, nlivast3 writes a text file, integ.out, containing the integrals of the regions. The file is written into the current experiment. Integrals are referenced to the integral specified by the argument well. The integral of spectrum from the sample specified by well is set to 1000. The macro nlivast3 does not add an additional region that is the sum of all the defined regions for each well (see dlivast). All of the spectra are integrated.
Arguments:	well is the number of the reference sample well. Reference integral set to 1000. The default reference is well 96.
plate_glue	Define a glue order for plotting and display (U)
Applicability:	Systems with VAST accessory
Syntax:	plate_glue
Description:	In a Unix terminal or shell window type plate_glue. The glue order is determined by clicking on the wells to be displayed. Save the glue order file in the user 's vnmrsys/templates/glue directory.
plvast	Plot VAST data in a stacked 1D-NMR matrix format (M)
Applicability:	Systems with the VAST accessory.
Syntax:	plvast<(display order,number of columns plotted)>
Description:	plvast will arrange and plot the traces from a reconstructed 2D data set (see vastglue)as an array of 1D spectra in a convenient format (as a matrix of 1D spectra). If no arguments are provided, the number of rows and columns will be determined by the periodicity of the display order. For example, if a block of 96 spectra, as is typical for a microtiter-plate, have been acquired using VAST automation, the spectra will be plotted in a matrix 8 rows and 12 columns.
	The default is to plot the spectra from 1 through arraydim (the number of spectra in the 2D data set). An optional argument (plvast(##)) allows one to specify that only spectra from 1 through ## should be plotted.
Arguments:	display order is optional and its default value is the glue order as listed in glueorderarray.
	number of columns plotted. The default value of is deduced by examining the periodicity of the requested display order. The number of columns plotted can entered as the second argument or as the first argument if the default display order is used.
plvast2d	Plot VAST data in a stacked pseudo-2D format (M)
Applicability:	Systems with the VAST accessory.
Syntax:	plvast2d<(number)>
Description:	If an array of 1D spectra have been acquired (in particular if a block of 96 spectra has been acquired using VAST automation, especially in a microtiter- plate format) and if these spectra have been glued into a reconstructed 2D data set (see vastglue), plvast2d will arrange and plot them (on the plotter) in a convenient pseudo-2D format (almost like an LC-NMR chromatogram). Well

labels are not attached to the spectra and spectra are plotted with 12 spectra per row.

Arguments: number specifies that only spectra from 1 through number should be plotted. The default is to plot all the spectra (from 1 through arraydim).

vastget Selects and displays VAST spectra (M)

Applicability: Systems with VAST accessory.

Syntax: vastget(<well>, <well>, ...)

Description: Selects and displays the spectra from any arbitrary well or wells using the well label(s) as arguments. The spectra are displayed in a dss stacked plot.

Arguments: well is the well label from which you want to select and display spectra. The wells are labeled [A->H][1-8].

vastglue Assemble related 1D data sets into a 2D (or pseudo-2D) data set (M)

Applicability: Systems with the VAST accessory.

Syntax: vastglue(<rack,<zone>)

vastglue(<glue order>,<plate>)

- Description: Used to artificially reconstruct a 2D data set from a series of 1D data sets having similar filenames. It is crucial to ensure that the format of the file names of each of the 1D data sets is identical. vastglue reads in each 1D file, in succession, and adds it to the previous data, but in a 2D format. It assumes that file names are of the format obtained when using the default setting of autoname (autoname=''). If autoname has been redefined, use a macro like vastglue2. Save the resulting reconstructed 2D data set in the normal manner using svf.
- Arguments: rack is the rack number; the default is 1. If you enter a rack number, you must also enter a zone number.

zone is the zone number; the default is 1. If you want to specify a zone number, you must enter a rack number.

glue order is the specific glue order to be defined based on the order defined in a plate_glue file. If glue order is specified, you can provide a plate number as the second argument and used with the glue order argument.

vastglue2 Assemble related 1D data sets into a 2D (or pseudo-2D) data set (M)

Applicability: Systems with the VAST accessory

Syntax: vastglue2<(number)>

Description: Used to artificially reconstruct a 2D data set from a series of 1D data sets having similar filenames. It is crucial to ensure that the format of the file names of each of the 1D data sets is identical.vastglue2 reads in each 1D file, in succession, and adds it to the previous data, but in a 2D format. It assumes that file names are of the format obtained using a non-default setting of autoname (autoname='filename_R%RACK:%_Z%ZONE:%_S%SAMPLE#:%_'). This definition must be hard coded into the macro by the user. If autoname has not been redefined, use a macro like vastglue. Save the resulting reconstructed 2D data set in the normal manner using svf.

Arguments: number is used to specify that only spectra from 1 through number are to be glued. The default is to glue all the spectra stored in the current directory that have the proper file name format (from 1 through arraydim).

10.6 Preparing the Hardware and Configuring VNMR

This section describes how to prepare the VAST hardware and configure VNMR for VAST operation. Use the procedures in this section if the Liquid Handler is being installed for the first time or if the Liquid Handler is being reconnected after it has been moved.

- Connecting the Transfer Tube
- Connecting the Air Tubing
- Connecting Signal and Power Cables
- Configuring VNMR for VAST

Connecting the Transfer Tube

This procedure describes how to connect the transfer tube between probe IN port and the injector valve (Gilson 819).

1. Determine the length of the transfer tube (between injector valve and probe), and then calculate the volume for this length using:


$$V = \pi r^2 l = \frac{\pi}{4} D^2 l$$


where V is the volume in cubic inches, r is the radius of the tube in inches, l is the length in inches, and D is the inside diameter (I.D.) in inches. To convert the volume from cubic inches to μ L, multiply V by 16387. Some examples are shown below.

Length	0.010-in. I.D.(Blue)	0.020-in. I.D.(Red)
5 ft	77 μL	308 µL
10 ft	154 μL	616 µL
20 ft	310 µL	1240 µL

- 2. Connect one end of the transfer tube (blue PEEK tubing) to the IN port on the probe.
- 3. Use a union to connect the other end of the transfer tube to the stainless steel tube coming out of port 5 of the Rheodyne injector valve.

Figure 96 shows the connection and part numbers for each part.





4. Place the inlet tube of the Liquid Handler (with filter on the end) into the solvent container.

The total volume from the injector port to the probe may be reduced by making a direct connection between the probe transfer tubing (going to the IN port on the NMR probe) and the injector port. This will bypass the Rheodyne injector valve.

Connecting the Air Tubing

This procedure describes how to connect the air tubing between the air regulator, the Valco valve, and the OUT port on the Microflow probe.

- 1. If the Microflow probe is not installed, install it now as described in the manual *Microflow NMR Probes Installation*.
- 2. Connect the Tygon (HDPE) tubing between an air regulator and port 3 of the Valco valve, as shown in Figure 97.
- 3. Use a PEEK connector to connect Teflon tubing between OUT port on the probe and port 2 of the Valco valve.
- 4. Make sure port 1 is open and insert a plug into port 4. Ports 5 and 6 are not used.

Connecting Signal and Power Cables

This procedure describes how to connect the signal and power cables.

- Make sure the GSIOC cable is connected between the GSIOC connector on the back of the injector valve and the GSIOC connector on the back of the Liquid Handler. Refer to Figure 98.
- 2. Make sure the Remote Switching cable is connected. Connect one end of the cable 01-905126-00 to the Remote Switching cable on the back of the Valco valve. Connect the other end to the OUTPUT connector on the back of the Liquid Handler, as follows:
 - a. On the OUTPUT connector, locate relay 2, which is the third set of two pins from the top.
 - b. Connect the red wire to the top pin.
 - c. Connect the white wire to the bottom pin.
- 3. Connect the RS-232 cable between the RS-232 connector on the back of the Liquid Handler and serial port A or B on the back of the Sun computer.
- 4. On the back of the Liquid Handler (below the RS-232 connector), verify that SW1 is set to 2 and SW2 is set to 6, as shown in Figure 98.
- 5. On the back of the injector valve (Gilson 819), verify that UNIT ID is set to 9, as shown in Figure 98.
- 6. Make sure the power is off in each of the following units:
 - Liquid Handler
 - Valco valve
 - · Injector valve

Connect the power cord for each unit into the power strip (with the power switch off), and connect the power strip to a power outlet.

7. Turn on the power to the units:




Figure 97. VAST Air Connections




Figure 98. VAST Signal Cable Connections

- Power strip
- Liquid Handler
- Valco valve
- Injector valve

Configuring VNMR for VAST

This procedure describes how to use the VNMR configuration window to configure VNMR for the VAST accessory.

- 1. Log into VNMR as vnmr1.
- 2. Enter **config** in the VNMR input window.
- 3. Set the following values in the VNMR Configuration window:
 - Set Sample Changer to VAST.
 - Set **Serial Port** to the Sun computer serial port (**A** or **B**) that is connected to the Liquid Handler.

- 4. Click Exit and Save in the VNMR Configuration window.
- 5. Enter **admintool** and verify the settings for the serial port used with the Liquid Handler. Make sure Template is set to Terminal-Hardwired and that Bidirectional *is not* set under Options, which is only shown if Detail is set to More or Expert.

To make changes to the serial port configuration, log in as root and use admintool to make the appropriate changes.

6. In the VNMR menus, select Main->Setup->LCNMR->Initialize VAST.

10.7 Calibrating Volumes and Flow Rates

This section provides procedures for calibrating the following parameters for VAST accessory operation:

- Probe Volume
- · Sample Volume
- Probe Slow Rate or Probe Slow Rate and Probe Slow Vol
- Calibrate needle rinse and Z position

Also in this section are some guidelines for calibrating the XYZ positions for the Liquid Handler arm.

To Calibrate Probe Volume

Probe Volume is the total volume of fluid contained in the transfer tube and the probe. You calibrate Probe Volume to determine the optimum amount of fluid necessary to adequately fill the flow cell of the Microflow probe. Figure 99 illustrates where the different volumes are in a VAST system.




Figure 99. Microflow Probe and Transfer Tube Volumes

Calibrating Probe Volume requires a comparison of several data sets to one of two control data sets. The first control data set is acquired after using a syringe to completely fill the Microflow probe with D_2O . The second control data sets is acquired after filling the probe

with the Liquid Handler (with either D_2O or a prepared sample) by setting the Probe Volume parameter larger than needed. Finally, several data sets are acquired (with either D_2O or a prepared sample), each with a slightly smaller value of Probe Volume, until the optimal probe volume is found.

One of two methods can be used to find the optimal Probe Volume:

- Using D_2O For Varian installation engineers and users who do not need to have Sample Volume calibrated. This method uses D_2O to minimize chemical handling, but D_2O does not allow you to do the Sample Volume calibration.
- Using a prepared sample For users who want to calibrate Sample Volume. This method, which is described in "To Calibrate Sample Volume," page 350, requires preparing a sample, such as 1% CH₃CN or CH₃OH in D₂O.

Use the following procedures to calibrate the Probe Volume parameter:

- "Obtaining a Control Data Set Using a Syringe," next
- "Verifying the Control Data Sets Using the Liquid Handler," page 347
- "Finding the Optimum Value of the Probe Volume Parameter," page 349

Start with a sample that provides an observable 1 H NMR signal. The default is D₂O. The sample should not contain any dissolved or undissolved solids.

Obtaining a Control Data Set Using a Syringe

You obtain the first data set after using the syringe to completely fill the Microflow probe with D_2O .

- Enter gilson to open the Liquid Handler window. Select the Main Control pane and make sure Air valve is set to OFF.
- 2. Disconnect the high-pressure air valve from the Microflow probe OUT connector. Connect a waste line to the probe OUT connector. See Figure 100.
- 3. Use the syringe and manually fill the Microflow probe with D_2O . Inject the D_2O until it begins to drip out of the waste line.

Refer to the manual *NMR Probes Installation* for details on filling and flushing the Microflow probe.

- 4. After the probe is full:
 - Tune the probe
 - Lock and shim
 - Acquire an NMR spectrum

Use a parameter set appropriate for determining pw90 on the Microflow probe. Use a 90° pulse and a sufficiently long d1 delay.

5. Expand around the resulting water resonance (about a 100-Hz window) and adjust the integral reset regions to obtain an integral of the residual HOD resonance.




Figure 100. Injecting Solution into a Microflow Probe

6. Save the file and write down the lock level.

You will eventually compare the resolution, lineshape, peak height, peak area, and lock level of future data sets to this data set.

7. Using the syringe, flush the probe with air to remove the D_2O solution.

Refer to the manual *NMR Probes Installation* for instructions on flushing the Microflow probe.

8. Disconnect the waste line from the probe OUT port and reconnect the air line.

Verifying the Control Data Sets Using the Liquid Handler

With this procedure, you use the Liquid Handler to verify the control data set. You can use either D_2O or a prepared sample.

Using D_2O – Read these points and go to step 1 below:

- Change the rinse solvent for the Liquids Handler to the same D₂O used above. Currently, the Liquid Handler is supplied with a waste container for the needle rinse station, but it is not supplied with any additional solvent reservoir containers.
- With this method, you use the same D₂O solution used for both the control data sets and the Probe Volume calibration (and for the rinse solvent).

Using a Prepared Sample – Read these points and go to step 1 below:

- With this method, you prepare a sample of your choice and fill a container in a Liquid Handler rack (instead of using just the rinse solvent). Although this requires the end user to make up a sample, such a sample is also required for the Sample Volume calibration, described in the procedure "To Calibrate Sample Volume," page 350.
- The test is easiest if the ¹H spectrum of the sample contains an isolated, tall, sharp resonance; a suggested sample is 1% CH₃CN or CH₃OH in D₂O (with D₂O as the rinse solvent), although almost any sample with a sufficiently strong signal (signal-to-noise greater than 100:1) can be used.
- Do not try to run the test by observing the HOD line of D₂O samples stored in open containers because such samples absorb too much H₂O from the atmosphere during the test to allow accurate peak areas to be measured.
- 1. Place the Liquid Handler inlet line (fritted) into the supplied bottle of D_2O .

Cover the opening of the D_2O bottle to prevent the absorption of H_2O .

- 2. If you are using a prepared sample, do the following substeps. If you are using D₂O, skip to step 3:
 - a. Prepare your sample and transfer it to the sample wells in a Liquid Handler sample rack.
 - b. Make sure each sample well contains enough sample to run the calibration. The total volume should be more than Probe Volume. Remember to write down the locations where you place the samples; you will use these locations in later steps.
- 3. Flush the Microflow probe with air:
 - a. Enter **gilson** in the VNMR input window to open the Liquid Handler window.
 - b. Click on the Main Control pane.
 - c. Set Air valve to ON.

- d. Wait 1 to 2 minutes to ensure that the Microflow probe is empty.
- e. Set Air valve to OFF.
- 4. Click on **Prime pump** in the **Main Control** pane of the Liquid Handler window. Wait for the priming routine to finish.
- 5. Switch to the **SAMPLE Def.** pane and set the following values:
 - Probe Slow Rate to 0.3
 - Probe Fast Rate to 0.3
 - Probe Slow Vol to 250
 - Number of Rinses to 0
 - · Keep Sample to No
- 6. Set **Probe Volume** to a value large enough to ensure that the active region of the Microflow probe is completely filled. Use one of the following values as a starting point.
 - More than 440 μ L for a 60- μ L probe
 - More than 600 μ L for a 120- μ L probe

The value needed also depends upon the length and diameter of the transfer tube that connects the Microflow probe IN port to the Liquid Handler inject port. The default transfer tubing is 0.010-in. I.D. and 10-ft. long.

This value must also be constrained to be somewhat smaller than the available syringe volume (Sample Volume must be less than SyringeVolume *minus* the absolute value of Sample Extra Vol *plus* Retrieve Extra Vol).




Figure 101. Microflo w Probe Flow Cell

- 7. Set the **Sample Volume** value as follows:
 - If you are using D₂O, set Sample Volume to 0. This value ensures that all of the fluid will come from the solvent bottle and not from the sample containers on the racks.
 - If you are using a prepared sample stored in a container in a rack, set Sample Volume to 150 µL for a 60-µL probe or to 250 µL for a 120-µL probe.
- 8. Click the Save Now button in the SAMPLE Def. pane.
- 9. Make sure all racks and vials are properly placed on the bed of the Liquid Handler so that the robotic arm will not accidently get caught on anything.
- **WARNING:** Avoid personal injury by using caution around the Liquid Handler work area. The Liquid Handler is capable of operating at high speeds. Leaning into the Liquid Handler work area, placing your hand or arm in its path, or wandering within range of the arm movements could result in serious physical injury.
- **CAUTION:** Avoid damage to the Liquid Handler and other equipment by clearing the work area around the Liquid Handler. The Liquid Handler arm moves rapidly around its work area.
 - 10. Enter the **loc** parameter as follows:

- If you are using D₂O, enter **loc=1** change. Note that the location used is unimportant since all the fluid comes from the rinse solvent container.
- If you are using a prepared sample stored in a container in a rack, enter **loc=# change**, where *#* is the location of the prepared sample.

Verify that the sample changed properly. If so, go to the next step.

- 11. Enter ga.
- 12. Compare the resulting spectrum to the control spectrum.

The spectrum you just obtained should be as good as the first control spectrum in all aspects (lock level, resolution, lineshape, peak height, and peak area/integral).

- If the peak height and peak area differ, but the lineshape and lock level are as good, it could be because of a different amount of H₂O absorbed from the atmosphere—if so, save this spectrum as control #2 for future comparisons.
- If the resolution is worse, or especially if the lock level is lower, either the current Probe Volume is not large enough to fill the NMR probe or something else is wrong that needs to be corrected before continuing.
- If the two spectra compare favorably and you have saved the second spectrum, enter: loc=0 change.
- 14. After the sample change finishes, manually turn on the air valve to ensure that the Microflow probe is completely empty before continuing with the rest of this procedure.

If part of the sample is not being recovered (for example, if Probe Fast Rate or Probe Slow Rate is too large), you can still perform the Probe Volume calibration but you need to perform the air-valve routine between each spectrum to blow the probe empty.

Finding the Optimum Value of the Probe Volume Parameter

The ideal Probe Volume value should be just slightly bigger than the minimum required volume (adding 10 to 20 μ L of extra volume) to ensure reliable filling during routine operation. Write down this optimized Probe Volume value should for future use.

Because the value obtained depends upon the probe and the transfer tubing being used, both of these variables should be documented. Table 49 lists typical Probe Volume values for various probe and transfer tube combinations.

Microflow Probe Flow Cell Volume, Tube Diameter	Probe Volume with 10-ft 0.01-in. ID transfer line	Probe Volume with5.5-ft 0.01-in. ID transfer line
60 μL, 0.005/0.01 in.	410 μL	320 (+20 µL)
60 µL, 0.01/0.02 in.	440 μL	350 μL
120 µL, 0.01/0.02 in.	520 μL	450 μL
240 µL, 0.01/0.02 in.	820 μL	750 μL

Table 49. Approximate Probe Volume Values for the Various Probes and Transfer Tubes

Note: The dead volume of a 10-ft length of 0.010-in. tubing is approximately 155 µL.

To find the optimum probe volume:

1. Acquire a series of spectra (on different sample locations) using different values of the Probe Volume parameter.

Use at least 10 different values for Probe Volume, in 10 μ L steps, decreasing from the Probe Volume value used in the previous procedure (this should be greater that the corresponding value in Table 49). Use the largest expected value of the Probe Volume and work towards smaller values. When you run spectra with different values of the Probe Volume, make sure the resulting value of the Push Volume remains positive.

Remember to click the **Save Now** button in the SAMPLE Def. pane after each change. It is always a good precaution to blow the probe empty once or twice during the series of spectra.

If the Flow Rate parameters have not been calibrated, the probe should be blown empty after each sample (the command loc=0 change can be used for this).

Invoke each change with loc=# change, where # is location of the prepared samples.

Do not shim between spectra. All measurements should be reproducible, ideally with minimal hysteresis, regardless of whether measurements range from maximum-to-minimum or minimum-to-maximum values.

2. Make a chart of the Probe Volume value versus the quality of the spectrum—lock level, resolution, lineshape, peak height, and peak area/integral.

As long as the Probe Volume is large enough to fill the active region of the Microflow probe, the spectral quality will be good. If the spectral quality is not good, the Probe Volume is too small. If the active region is not filled, the spectral quality will be poor and will rapidly deteriorate as the Probe Volume gets smaller. Also, any samples in the future will be hard to shim.

- 3. Choose the smallest value of Probe Volume that still generates good data. Write down this Probe Volume value as well as the probe and transfer tube geometry (length and diameter) used.
- 4. After you find the optimal Probe Volume, use the same sample to calibrate Sample Volume as described in the next section "To Calibrate Sample Volume."

To Calibrate Sample Volume

Sample Volume can be equal to Probe Volume, or it can be a portion of Probe Volume, where the remainder of Probe Volume is filled with Push Volume. In other words, the sum of Sample Volume and Push Volume equals Probe Volume. Portioning Probe Volume into Sample Volume and Push Volume allows you to use less sample because you fill the dead volume of the transfer tube with Push Volume. Using a gas bubble to push the sample is undesirable.

You calibrate Sample Volume to find the optimum (minimum) volume of sample to fill the flow cell of the Microflow probe. An accurate calibration of Sample Volume depends upon an accurate value of Probe Volume. You should verify or recalibrate Probe Volume immediately before calibrating Sample Volume.

As a general rule, the larger the value of Push Volume, and the "stickier" the solvent, the lower the NMR sensitivity—DMSO is "sticky," D₂O is not; other organic solvents are intermediate. As another general rule, the bigger the difference in magnetic susceptibility between the sample and the rinse solvent, the bigger the chance that NMR resolution will suffer for a given Sample Volume. If the NMR resolution suffers when a push solvent is used, then you probably need to either increase Sample Volume or match the magnetic susceptibilities of the two solvents better; the line broadening effect usually cannot be corrected by shimming.

The Sample Volume calibration requires some visible difference in the ¹H NMR spectra between the sample and the rinse solvent. The sample needs to have some additional ¹H NMR resonances, which are not present in the rinse solvent. You can use 1% CH₃OH, CH₃CN in D₂O, or the same prepared sample used for calibrating Probe Volume. This calibration cannot be done with D₂O alone (see "To Calibrate Probe Volume," page 345).

- Note: Do not reshim during this calibration. Any given set of conditions should produce the same spectral quality reproducibility—if not, it might be a good idea to enter loc=0 change and then blow the probe empty with the air valve between each run. If this is necessary, it suggests that the values of Probe Slow Rate or Probe Fast Rate are too large.
 - 1. If not already done, prepare your sample and transfer it to the sample wells in a Liquid Handler sample rack.

Make sure the sample wells contain sufficient amounts of the sample to accomplish the calibration.

- 2. Set the following values in the SAMPLE Def. pane of the Liquid Handler window:
 - Sample Volume to 150 uL for a 60 uL probe or to 250 uL for a 120 uL probe
 - Keep Sample to No
 - Number of Rinses to 0
 - **Probe Volume** to the value calibrated by the procedure "To Calibrate Probe Volume," page 345.
- 3. Click the Save Now button in the SAMPLE Def. pane.
- 4. Obtain a spectrum on the first sample and expand around a representative sharp peak (other than HOD) in the ¹H NMR spectrum.
- 5. Save this data set and note the peak height (and lineshape).
- 6. Obtain additional spectra with decreasing values of Sample Volume (Probe Volume remains fixed) and look at the peaks heights (or S/N).
- 7. Plot the peak heights (or S/N) as a function of the Sample Volume.

The NMR peak heights should decrease as the Sample Volume parameter decreases, even when the Probe Volume parameter is kept constant. (the NMR peak heights should drop to zero when Sample Volume is 0).

The exact nature of the curve you plot depends upon the solvent being used (DMSO generates a curve, while D_2O produces more of an angle or *knee* at the intersection of two straight lines), but this chart allows you to determine how much sample must be used to obtain the percentage of the available signal-to-noise. See Figure 102.

For a 60-uL probe in D_2O with Probe Volume set to 350, the signal intensity when Sample Volume is 150 should be 75 to 95% of the signal intensity obtained when Sample Volume is 350.

To Calibrate Flow Rate Parameters

This section describes how to calibrate the parameters Probe Slow Rate, Probe Slow Vol, and Probe Fast Rate.

- If these parameters are set too large, the sample is not completely removed from the probe and subsequent samples are significantly contaminated.
- If these values are set too small, time is wasted by running the system at a slower than optimal speed.




Figure 102. Finding the Optimum Sample Volume

The values obtained are very dependent upon the following:

- · Sample solvent
- Rinse solvent
- Probe flow cell volume
- Tubing used, length, I.D., etc.
- Pressure of the gas delivered to the Valco valve.

Choose one of the three methods described below. In all three methods, Probe Slow Rate and Probe Fast Rate must both be small enough (giving slow enough flow rates) that the following are true:

- No leaks or drips occur at either the syringe vent port or the injector port when pushing
- No cavitation bubbles appear in the syringe when retrieving solvents (or in the Gilson tubing connecting the needle to the syringe)
- When Keep Sample is set to Yes, all of the sample is recovered

Method 1: Fast Setup, Slow Runs

This method consists of just setting the parameters to be so slow that the solvent will never leak or cavitate. This method is fast to set up (since no calibrations are done), but slow in operation (since no parameters are optimized). This method is not recommended for high-throughput operation, but may be acceptable for testing or open-access operation.

- 1. Set the following values in the SAMPLE Def. pane of the Liquid Handler window:
 - Probe Slow Vol equal to Probe Volume or larger.
 - Probe Slow Rate to 0.3
 - Probe Fast Rate to 0.3

Larger values might be possible for your system, but 0.3 mL/min should work for the most viscous solvent in the most narrow tubing.

If Probe Slow Vol equals Probe Volume, setting Probe Fast Rate has no effect, but you should set it to a small value just for safety.

2. Click the Save Now button.

Method 2: Calibrating Probe Slow Rate for Faster Runs

This method consists of performing a simple calibration to allow the Liquid Handler to run faster. This flow rate value must be determined for the most viscous sample in the system. This might be the rinse solvent, a push solvent, or a sample, but it must be the most viscous sample you expect to encounter!

- 1. Set the following values in the SAMPLE Def. pane of the Liquid Handler window:
 - · Probe Slow Vol equal to Probe Volume or larger.
 - Probe Slow Rate to 0.3
 - Probe Fast Rate to 0.3
 - Number of Rinses to 1
- 2. Click the Save Now button.
- 3. Run samples with the solvents in question repeatedly, increasing the Probe Slow Rate in 0.1 mL/min steps each time, while watching the clear glass syringe barrel on the Liquid Handler.

If bubbles start appearing in the syringe barrel or needle tubing when the sample is being withdrawn from the probe, Probe Slow Rate is too large (the bubbles typically get bigger as the Probe Slow Rate value gets larger). Any cavitation bubbles formed must disappear *before* the needle (probe) is pulled up out of the injector port; otherwise, sample is left behind in the probe (all the bubbles always disappear the moment the needle [probe] is pulled out of the injector port).

4. Use the largest value of Probe Slow Rate that still allows complete sample recovery. Table 50 lists some typical flow rates for common solvents.

Solvent	Typical Flow Rate*
CH ₃ CN	4.4mL/min for 0.01/0.02**
50:50 CH ₃ CN:D ₂ O	1.5 mL/min for 0.01/0.02**
D ₂ O	1.0 mL/min for 0.01/0.02** 0.7 mL/min for 0.005/0.01**
DMSO	0.5 mL/min for 0.01/0.02**

 Table 50.
 Typical Flow Rates for Common Solvents

* Flow rates depend primarily upon the inner diameter of the inlet tubing.

** Probe inlet tubing diameter (in inches/outlet tubing diameter [in inches]).

Method 3: Calibrating Probe Slow Rate and Probe Slow Vol for Fastest Runs

This method consists of performing a more complex, three-parameter calibration. Although this calibration allows the system to run with the fastest possible throughput, it needs to be reoptimized for every different solvent and plumbing geometry. The nature of the calibration precludes a simple step-by-step description.

After the majority of a sample (or rinse) is withdrawn from the Microflow probe, the remainder can be withdrawn faster. The remainder is withdrawn at a speed determined by the Probe Fast Rate parameter. This calibration consists of balancing Probe Slow Vol and Probe Slow Rate with Probe Fast Rate.

A faster Probe Slow Rate or a smaller Probe Slow Vol necessitates a slower Probe Fast Rate.

The criteria for optimization are as follows:

- · No cavitation bubbles in the syringe or needle tubing
- · The full volume of recovered samples is returned to the original sample containers
- Sample insertion and removal occur as fast as reasonably possible
- Probe Slow Vol is greater than Sample Volume

To Calibrate XYZ Positions of the Arm

You have control over the XYZ coordinates of the rinse station, injector port, and sample racks on the Liquid Handler. You can make changes to these units through the Calibrations pane of the Liquid Handler window See "VAST Interface Description" on page 355.

The X and Y positions of all three units will probably never need recalibration after installation. You might possibly need to calibration the Z position for the following:

- *Rinse station* You can adjust the Z position of the needle to place the needle up to 45 mm down inside the rinse station for an inside-and-outside needle rinse. See "To Calibrate XYZ Positions of the Arm" on page 354.
- *Injector port* You might have to optimize the Z position after changing the needle to ensure a leak-free injection.
- Sample racks No calibration of Z position needed.

10.8 Acquiring Data on Standard Test Samples

After the Probe Volume is calibrated, you can acquire data on some standard test samples to verify that all the VAST hardware and software are operational.

1. Fill four different sample locations with a standard sample.

This test does not require a specific sample. You can use D_2O , the prepared sample, sucrose, or methyl p_hydroybenzoate.

- 2. Set the following values in the SAMPLE Def. pane of the Liquid Handler window:
 - Set **Probe Volume** to the value determined in the section "To Calibrate Probe Volume," page 345.
 - Sample Volume to 150 µL for a 60-µL probe or to 250 µL for a 120-µL probe
 - Keep Sample to Yes
 - Number of Rinses to 1
 - Make sure the **Push Volume** is positive
- 3. Click the Save Now button.
- 4. Acquire spectra both manually and with the enter program.
 - a. Acquire spectra manually:

loc=# ga

where # is one of the four locations filled in step 1.

b. Acquire the same spectra using the enter program to obtain ¹H NMR spectra of the four samples.

10.9 Evaluating Carryover

The amount of carryover from one sample to the next can never be zero in a flow probe. The amount of carryover is affected by the probe geometry, the direction(s) of solvent flow, the number of rinses, the effectiveness of sample recovery, and in particular by the solvents used. Multiple rinses reduce carryover, but the additional time consumed by more that one rinse is probably prohibitive.

To measure the amount of carryover:

- 1. Select a solute that generates an NMR singlet (e.g., CH₃OH, CH₃CN, TSP, or DSS) and make a concentrated solution of this solute (1 to 10%) in your solvent of choice (e.g., D₂O).
- 2. Place this sample in location 1, then place pure solvent in locations 2 and 3.
- 3. Set up proper parameters in the VNMR parameter set and in the Liquid Handler window.
- 4. Enter loc=1 change.
- 5. Enter **ga** to obtain a ¹H NMR spectrum. Expand the spectrum around the region of interest, and set the vertical scale, integral resets, and integral scale.
- 6. Save the data and plot the data:

```
svf(filename)
```

```
where filename is a name you choose (e.g., svf('mydata')).
```

```
region vp=12 pl dpir page
```

- 7. Enter loc=2 ga.
- 8. Use the resulting spectrum to quantitate the amount of carryover remaining from sample 1.
- 9. Enter **loc=3** ga to repeat the process for the sample in location 3.

With one rinse, the carryover averages about 1%, but the value you obtain can be highly variable (0 to 8%) depending upon the solvents used and other specifics of the experiment.

10.10 VAST Interface Description

The user interface on the front panel of the Gilson 215 Liquid Handler is a display and a STOP button.

The VAST interface to VNMR includes the LC-NMR pane (shown in Figure 77), which provides controls for setting up VAST experiments and the Gilson Liquid Handler window shown in Figure 75.

Otherwise, the VAST accessory is set up using the Liquid Handler window in VNMR, which provides the following four panes for setting up the VAST accessory:

- SAMPLE Def. (shown in Figure 103) For defining sample volumes and fluid-control parameters.
- Rack Def. (shown in Figure 104) For defining which racks and containers are available to the Liquid Handler.
- Main Control (shown in Figure 105) For setting up and controlling the hardware.
- Calibrations (shown in Figure 106) For calibrating arm positions.

Each of these panes is described in a separate section below. Table 48 lists the default values for the fields in the Liquid Handler window.

SAMPLE Def.

The Sample Def. pane (shown in Figure 103) contains most of the routinely used fluidcontrol parameters. These parameters allow you to

- Adjust sample volume
- Specify whether to keep the sample
- · Perform probe rinses
- Specify flow rates

Clicking the **Save Now** button saves the values in the Liquid Handler panes and applies them to the next injection. Clicking the Additional Parameters button expands the screen and provides more fields, as shown in Figure 103.

-		Gilson	Liquid Handler		•
SAMPLE Def. R	ack Def.	Main Control	Calibrations		
Select parameters	PAK 60u10	in2765 D20 2 🛁			
ProbeVolume:	340.0		ProbeSlowVol:	250	
SampleVolume:	140.0		ProbeSlowRate:	1.5	
Push Volume:	200.0		ProbeFastRate:	2.5	
SampleKeepFlag:	yes n	5	SampleExtraVol:	25.0	
NumRinses:	0		RinseExtraVol:	25.0	
RinseDeltaVol:	10.0		Additional Paramete	rs	
SampleWellRate:	4.0		SampleMixVol:	0	
SampleHeight:	20		SampleMixRetrievVol:	0	
SampleDepth:	NOSEEK		SampleMixRetrievRate:	d	
NeedleRinseVolume:	500				
NeedleRinseRate:	16.0				
Save Now					

Figure 103. Sample Definition Pane

After the values are set and saved in this window, all NMR activities (typically controlled by either the enter program or the loc command) will use the same values for every sample, until the parameters in the Liquid Handler window are changed again.

The SAMPLE Def. window contains the following fields:

Select parameters The file names of the available sets of Gilson parameters. The files are selected by a pull-down menu. The file name of an active file can be changed by double-clicking on the current name.

ProbeVolume	 Total volume, in μL, held in the probe flow cell plus the volume held in the transfer tubing. In other words, this is the volume of liquid that the Liquid Handler will need to fill the probe flow cell, the dead volume in the tubing, etc. These values can vary widely, depending upon the length and inner diameter of the transfer tubing (see Table 49); however, typical values are as follows: For a 60-μL Microflow probe, the default value of Probe Volume is about 440 μL (with 10 ft of 0.010-in. I.D.
	transfer tubing).
	 For a 120-μL Microflow probe, the default value of Probe Volume is about 520 μL.
SampleVolume	The amount of sample, in μ L, needed to fill the probe flow cell. This is the amount that will be removed from the sample container. This value depends only on the Microflow probe that is used and not on the transfer tube. This volume can be calibrated to obtain the best NMR sensitivity per amount of sample injected.
	The default value is to set Sample Volume equal to Probe Volume; however, smaller sample quantities can be used at the expense of sensitivity. If this value is set too small, the NMR sensitivity, and possibly the NMR resolution, decreases. If this value is too large, you will consume more sample than is necessary.
	 For a 60-μL Microflow probe, a typical value of Sample Volume is about 140 μL.
	 For a 120-μL Microflow probe, a typical value of Sample Volume is about 240 μL
PushVolume	The amount of solvent, in μ L, used to "push" the sample into the probe flow cell. Push Volume is automatically calculated Probe Volume minus Sample Volume). A push volume allows you to minimize the amount of sample used to create the Probe Volume.
SampleKeepFlag	Sets whether or not the sample is kept after NMR is finished.
	• Set to Yes to have samples returned to the vials on the sample racks.
	• Set to No to have samples placed in waste. The used sample is flushed into the rinse station container.
NumRinses	The number of NMR probe rinse cycles performed between sample changes. Typical values are either 0 or 1, where 1 is recommended. A needle rinse is performed between each sample change, but a probe rinse is only performed if a value is entered in the Number of Rinses field. More rinses (e.g., Number of Rinses set to 2 or more) reduces carryover, but takes longer.

	An NMR probe rinse cycle consists of the following actions: After the previous sample is aspirated from the probe and dispensed to waste or to the sample vial, the syringe aspirates from the solvent reservoir a volume equal to Probe Volume + Rinse Extra Vol. This volume is injected (dispensed) into the probe, and then immediately withdrawn from the probe and dispensed into the waste container.
RinseDeltaVol	The amount of solvent, in μ L, in addition to the Probe Volume, injected into the probe during a rinse cycle. This volume ensures that the flow cell is thoroughly rinsed. A typical default is 10 μ L; any value over 2 μ L seems sufficient.
ProbeSlowVol	The volume, in μ L, removed from the NMR probe at the Probe Slow Rate. An easy default is to set Probe Slow Vol equal to Probe Volume.
	• If Probe Slow Vol is too small, the sample might not be completely removed from the NMR probe.
	• If Probe Slow Vol is too large, turnaround times might be longer than necessary.
	After the Probe Slow Vol is removed, the rest of the Probe Volume is removed at the faster rate (Probe Fast Rate).
ProbeSlowRate	The flow rate, in mL/min., used to remove the initial volume specified by Probe Slow Vol.
	• If this rate is too fast, bubbles will appear in the syringe when the sample is withdrawn (because of the vacuum that is formed) and not all of the sample is recovered when the sample is removed from the probe.
	• If this rate is too slow, turnaround times might be longer than necessary.
	Values typically range from 0.5 for DMSO to 1.0 for D_2O to 1.5 for CH ₃ CN but depend greatly on the sample and hardware.
ProbeFastRate	The flow rate, in mL/min., used to inject samples into the NMR probe and to remove the remaining volume (after Probe Slow Vol) from the NMR probe. Some factors that might limit this speed are solvent composition (viscosity), NMR probe volume, and tubing diameter.
	• If the Probe Fast Rate is too fast, the pressure relief valve next to the syringe leaks when sample is pushed into the NMR probe, and bubbles can appear in the syringe when the sample is withdrawn (because of the vacuum that is formed). As a result, not all of the sample is recovered when the sample is removed from the probe.
	• If this rate is too slow, turnaround times might be longer than necessary.

SampleExtraVol	The extra amount of sample, in μ L, removed from the probe when the sample is aspirated from the probe. This volume helps ensure that all the sample is recovered and that the probe is clear before the next sample is injected. A typical default is 25 μ L; any value over 5 μ L seems to be sufficient.
RinseExtraVol	The extra amount of solvent, in μ L (in addition to Probe Volume and Rinse Delta Vol), that is removed from the probe during a rinse cycle. This volume ensures that more solvent is removed than injected, thus clearing the probe. A typical default is 25 μ L; any value over 5 μ L seems to be sufficient.
Save Now	This button applies the indicated values and saves the current settings to the Sun computer.
Additional Parameters	This button expands the SAMPLE Def. pane and provides the following parameters:
SampleWellRate	The rate, in mL/min, at which samples are aspirated from, and dispensed to, the sample container. A typical default is 4.
SampleHeight	The height (in units of 0.1 mm) above the rack at which the needle is placed before aspirating the sample. A value of zero should place the needle against the metal rack; a value of 25 (2.5 mm) is commonly used (but only if Sample Depth is set to NOSEEK. Because this may coat the outside of the needle with sample, the needle rinse station should be recalibrated to wash both the inside and outside of the needle.
SampleDepth	The depth of the liquid sample in the sample container. This value is used to allow the needle to follow the sample amount down as sample is withdrawn from the container. This parameter is commonly set to NOSEEK (especially when the Sample Height is set to 25) so the needle does not move down as the sample is aspirated.
NeedleRinseVolume	The volume of rinse solvent used to wash the needle. The default (and maximum usable value) is to set this equal to the syringe volume.
	Values that are too small could lead to increased carryover from one sample to the next.
	Values that are too large increase rinse solvent usage and increase the turnaround time.
NeedleRinseRate	The rate at which the syringe aspirates and dispenses rinse solvent through the needle. The default (and maximum usable rate) is 16 (mL/min).
Mix Volume	Defined but not yet implemented.
Mix Time	Defined but not yet implemented.
Mix Flow Rate	Defined but not yet implemented.
Mix Height	Defined but not yet implemented.

Rack Def. Pane

Use the Rack Def. pane (shown in Figure 104) to set up sample racks on the Gilson Liquid Handler.

The settings for racks 1 to 5 must reflect what racks exist on the system; otherwise, severe damage to the needle (probe) can result.

		Gilson Liquid Handler	•
SAMPLE Def.	Rack Def.	Main Control Calibrations	
	Rack 1 205h — Zones:	Rack 2 Rack 3 Rack 4 202 Rack 5 Zones: Zones: Zones: Zones: Zones: Left & Back ◇ Right & Back Left & Front ◇ Right & Front 	
	() () () Number	2 3 1 2 3 1 4 7 1 6 7 5 6 6 5 4 2 5 8 2 5 8 8 9 7 8 9 3 6 9 3 4 9 locations with: Numerals (1-96) Names (A1-H12)	

Figure 104. Rack Definition Window

Main Control

Use the Main Control pane (shown in Figure 105) to control the Liquid Handler hardware.

-	Gilson Liquid Handler	
SAMPLE Def. Ra	ack Def. Main Control Calibrations	
Í Í		
	Air valve: ON Plunger: UP OFF DOWN	
	Reset arm Prime pump Return Home	
	Syringe Volume (ul): 500 🛁	
	Arm Z scale (mm): 121 💻	

Figure 105. Main Control Window

This pane contains the following fields:

Air valve	Sets the Air valve to ON or OFF. The default is OFF, which allows the valve to be triggered by the VNMR automation software during normal operation. The valve provides gas pressure to the probe outlet port for pressure-assisted sample withdrawal. The gas pressure works in conjunction with gravity and the syringe to speed sample withdrawal and reduce cavitation. For normal operation, set Air valve to OFF. To provide continuous gas pressure to the probe (for example
Plunger	to empty the probe), set Air valve to ON. Moves the syringe plunger up or down. During normal operation, Plunger is set to UP, which is the default. When replacing the syringe, select DOWN; after the new syringe is in place, select UP.
Reset arm	Resets the Liquid Handler controller after an error or fault. If an error or fault occurs while the VAST autosampler is running, click on the Reset arm button to reset the Liquid Handler controller.
Prime pump	Moves the needle to the rinse station and runs the pump priming cycle, which takes about 3 minutes. Use this button to prime the pump after switching to a new solvent or after the system has been moved. When priming the pump, make sure no bubbles are in the inlet tubing to the barrel.
Return Home	Returns the arm to the predefined home position.
Syringe Volume	The syringe volume, in μ L. This number must match the value printed on the syringe module.
Arm Z scale (mm)	Sets the arm Z scale value to accommodate the Liquid Handler probe (needle) module installed on the system. The scale from 0 to 220 mm is printed on the side of the vertical arm assembly. This value is normally set at installation; however, if you replace or adjust the position of the Liquid Handler vertical arm assembly, be sure to update the Arm Z scale value.

Calibrations

Use the Calibrations pane (shown in Figure 106) to calibrate the various positions of the Liquid Handler arm and probe (needle). The X, Y, and Z adjustments work in real time, so you must watch the probe move as you make adjustments.

- 1. In the Cal. Selections menu, select rinse, sample, or inject.
- 2. Use the X and Y adjustment arrows to center the Liquid Handler probe over the position you are calibrating.
- 3. Use the Z adjustment arrows to calibrate the height of the Liquid Handler probe as follows:
 - *rinse* For an *inside* wash, lower the needle until the tip is centered directly above the center of the rinse station. For an *inside and outside* wash, lower the

-		Gilson	Liquid Ha	ndler			F	
SAMPLE Def.	Rack Def.	Main Control	Calibratio	ins				
				Cal. :	Selections	:	1	
		×1:	×5:	×10:	×50:	×100:		
	X Position :			$\square \square$	$\square \square$	$\square \square$		
	Y Position :			$\square \vdash$	$\square \square$	$\square \vdash$		
	Z Position :			$\square \vdash$	$\square \square$	$\square \vdash$		
						Save		

Figure 106. Calibrations Window

needle until the tip is inserted several centimeters into the rinse station (depending on how much of the outside of the needle is contaminated; refer to the Sample Height parameter).

- *rack* Lower the needle below the level of the sample. Keep in mind that the needle lowers as it draws sample out of the container, unless the NOSEEK option is used for Sample Depth.
- *inject* Lower the needle 0.1 mm (x1) at a time until the needle makes snug contact with the injector port. If the needle begins to bow, back off until the needle is straight.
- 4. After a position is calibrated, click on **Save**. Then select another position, until all three positions are calibrated.

The values of X, Y, and Z for the center of the injector manifold, and the values of X and Y for the racks, are written into the file /vnmr/asm/racksetup. The XYZ coordinates of the rinse station are stored only in the Liquid Handler memory.

If bubbles appear in the transfer line (between the needle and the syringe) as the sample is withdrawn, this usually indicates that the needle is not calibrated to go deep enough into the inject port.

10.11 Customizing the enter Window for VAST

The file entervast.conf, located in /vnmr/asm, controls the appearance of the enter window, when Sample Changer is set to VAST in the VNMR configuration window.

The entervast.conf file uses information in the files experiments.vast, solvents.vast, users.vast and protocol.vast to determine which buttons are displayed. If you want to define protocols for injecting and removing samples, you will also need to uncomment the *infield* line in entervast.conf. The experiments.vast file typically points to the autoscout macro (among others).

After you click Exit and Save in the enter program, the text file to be used by autogo is written out. This text file (the enter queue) can be edited, if desired, before

running autogo. For further information on customizing the enter program, refer to the manual *VNMR User Programming*.

An example of an autoname value that is especially useful for VAST is autoname='sucrose_R%RACK:%_Z%ZONE:%_S%SAMPLE#:%_'

10.12 Files that Control VAST Operation

The files that control VAST automation are located in the directory /vnmr/asm.

current	contains information that indicates the current sample in the magnet. Additionally, the sampleInfo file (a link to samp0 to samp9 files) is appended into this file.
default	contains all the possible parameters, set to appropriate default values.
info	directory that contains the sample definition files.
racks	contains the supported Gilson rack definition files (e.g., code_205.grk). This directory also contains the file (rackInfo) used by the gilson program that defines the graphical characteristics of the racks.
racksetup	contains the injector and all racks and alignment values.
samp0 to samp9	are in the info directory and contain the parameters and values as set by the user (using the gilson program) to override any defaults. The gilson program modifies the parameter values in these files. These parameters define the conditions for injecting and removing a sample from the probe.
tcl	contains the standard set of Tcl scripts provided by Varian NMR.
protocols	a directory to hold VAST protocols.

10.13 Writing VAST Protocols

This section describes the elements for programming custom VAST protocols. Some familiarity with the Tcl programming language will be useful. Sample protocols are in the /vnmr/asm/protocols directory.

gWriteDisplay	Display message on Gilson display (8 characters maximum)
Syntax:	WriteDisplaymessage
Description:	Displays the string message on the LED display of Gilson 215.
Arguments:	message is a string of 8 characters maximum.
Examples:	gWriteDisplay "Rinse"
	set msg [format "Samp: %d" \$SampleNumber]
	gWriteDisplay "\$msg"
gMoveZ2Top	Move needle to further most upward position
Syntax:	MoveZ2Top
Description:	Moves the Gilson 215 needle to the fully retracted position.
Examples:	gMoveZ2Top

gMove2RinseStation Move to Rinse Station

Syntax: gMove2RinseStation

Description: Moves the Gilson 215 needle to the Gilson's Rinse Station and lowers the needle. Note the Rinse Station position is maintained within the Gilson's internal EEPROM.

Examples: gMove2RinseStation

gFlush Draw solvent from reservoir, expel through needle

Syntax: gFlush Volume InFlow OutFlow

Description: Draws the given volume from the solvent reservoir at the given flowrate and expels the same volume of solvent through the Gilson needle at given flowrate.

When used with an InFlow greater than zero and OutFlow of zero, gFlush is asynchronous and allows the next command to proceed while the syringe is filling. gFlush with InFlow as zero and OutFlow with a greater than zero value waits until the syringe is filled. However, if you use other commands that test or perform actions with the syringe then follow gFlush with the command gStopTestAll. This command ensures that any active command is complete before proceeding..

Arguments: Volume is the volume in μ l.

Inflow is the flow rate in ml/min in which solvent is drawn from solvent reservoir; if set to zero, then no solvent in drawn.

OutFlow is the flow rate in ml/min in which the solvent is expelled through the Gilson needle. If Zero, then no solvent is expelled.

Examples: gFlush \$NeedleRinseVolume \$NeedleRinseRate \$NeedleRinseRate gFlush \$RinseVol \$MaxFlow \$ProbeFastRate gFlush \$RinseVol \$MaxFlow 0 gStopTestAll SetCurrentVolume [gCurrentSyrVol]

Related: gStopTestAll Wait for all axis motion and syringe operation to complete

gMove2InjectorPortMove to Injector port

Syntax: gMove2InjectorPort

Description: Moves the Gilson 215 needle to the Gilson's Injector Port and lowers the needle. Examples: gMove2InjectorPort

gDelayMsec Delays TCL instruction execution

Syntax: gDelayMsec delaytime

Description: Delays Tcl instruction execution for the given number of milliseconds.

Arguments: delaytime is the time in milliseconds to delay Tcl execution.

Examples: gDelayMsec 500

gDelayMsec 1000 (sets a 1 second delay)

gSetContacts Sets output contacts on Gilson 215 Syntax: gSetContacts contact# state Description: Sets the Gilson's output contacts to connected or disconnected. Arguments: contact# is the contact number (1 to 4) to connected or disconnected. state is 1 for connected or 0 for disconnected. Examples: gSetContacts 2 1 switches the gas valve to put pressure on NMR Probe. gSetContacts 2 0 switches the gas valve off. Draw in liquid through Gilson needle gAspirate Syntax: gAspirate volume flow Zspeed Description: Draws the specified volume in through Gilson needle at the given flowrate, with needle lowering at the speed given. The speed depends on the flowrate and sample well dimensions. Arguments: volume is the volume in μ l to draw. flow is the flowrate in ml/min. Zspeed is the speed in mm/sec to lower Gilson needle while drawing in liquid through the needle. If Zspeed is set to zero, the needle is not lowered. Examples: gAspirate 250 4 0.0 gAspirate \$SampleVolume \$SampleWellRate \$SampleZSpeed gDispense Expel liquids though Gilson needle Syntax: gDispense volume flow Zspeed Description: Expels the specified volume out through Gilson needle at the given flowrate, with needle lowering at the speed given. Speed depends on flowrate and sample well dimensions. Arguments: volume is the volume in μ l to draw. flow is the flowrate ml/min. Zspeed is the speed in mm/sec to lower Gilson needle while expelling liquid through the needle. If Zspeed is zero, the needle is not lowered. Examples: gDispense 250 4 0.0 gDispense \$SampleVolume \$ProbeFastRate 0.0 gMoveZ Move Gilson needle to a Z position Syntax: gMoveZ Z_location Description: Move Gilson needle to the given Z location. Arguments: Z_location is the Z location in tenths of a millimeter. Examples: qMoveZ \$SampleZTop gMove2Sample Move Gilson needle to sample well location Syntax: gMove2Sample Rack Zone Sample Description: Move Gilson needle to the sample well location. Arguments: Rack is the rack location 1 to 7.

Zone is the rack zone, which is rack dependent.

Sample is the sample tray number.

Examples: gMove2Sample \$RackLoc \$SampleZone \$SampleNumber

gInjector2Load Switch Rheodyne injector valve to the load position

- Syntax: gInjector2Load
- Description: Switches the Rheodyne injector valve to the load position.

Examples: gInjector2Load

gInjector2Inject Switch Rheodyne injector valve to the Inject position

Syntax: gInjector2Inject

Description: Switches the Rheodyne injector valve to the inject position.

Examples: gInjector2Inject

gCurrentSyrVolReturn current liquid volume in syringe

Syntax:	SetCurrentVolume [gCurrentSyrVol]
Description:	Returns the current volume in μ l in Gilson syringe pump.
Arguments:	$\texttt{CurrentVolume}$ is the current volume in μl in Gilson syringe pump.
Examples:	SetCurrentVolume [gCurrentSyrVol]

gAspirateAsync Draw in liquid through Gilson needle

Syntax: gAspirateAsync volume flow Zspeed

- Description: Draws the specified volume in through Gilson needle at the given flowrate, with needle lowering at the speed given. This function returns immediately, so other operation can be accomplished while aspirating. The speed depends on flowrate and sample well dimensions. Use gStopTestAll before performing any axis moves.
- Arguments: volume is the volume in μ l to draw.

flow is the flowrate in ml/min.

Zspeed is the speed in mm/sec to lower Gilson needle while drawing in liquid through the needle. If Zspeed is zero the needle is not lowered.

Examples: gAspirateAsync 250 4 0.0

gAspirateAsync \$SampleVolume \$SampleWellRate \$SampleZSpeed (starts aspiration operation)

gSetContacts relay state (while still aspirating switch a relay)

gStopTestAll (now wait for syringe operation to complete before moving)

gilMoveZ2Top (move needle up)

gMove2InjectorPort (move to injector port)

gDispenseAsync Expel liquids though Gilson needle

Syntax: gDispenseAsync volume flow Zspeed

Description:	Expels the specified volume out through Gilson needle at the given flowrate, with needle lowering at the given speed. This function returns immediately, so other operation can be accomplished while dispensing. Speed will depend on flowrate and sample well dimensions
	Use gStopTestAll before performing any axis moves.
Arguments:	volume is the volume in μ l to draw.
	flow is the flowrate in ml/min.
	Zspeed is the speed in mm/sec to lower Gilson needle while expelling out liquid through the needle. If Zspeed is zero the needle is not raised.
Examples:	gDispenseAsync 250 4 0.0
	gDispenseAsync \$SampleVolume \$SampleWellRate \$SampleZSpeed (start dispensing operation)
	gSetContacts relay state (while still aspirating switch a relay)
	gStopTestAll (now wait for syringe operation to complete before moving)
	gilMoveZ2Top (move needle up)
	gMove2RinseStation (move to rinse station)
gStopTestAll	Wait for all axis motion and syringe operation to complete
Syntax:	gStopTestAll
-	Return when all axes and syringe operations have completed.
Examples:	gStopTestAll
	gAspirateAsync \$SampleVolume \$SampleWellRate \$SampleZSpeed (start aspiration operation)
	gSetContacts relay state (while still aspirating switch a relay)
	gStopTestAll (now wait for syringe operation to complete before moving)
	gilMoveZ2Top (move needle up)
	gMove2InjectorPort (move to injector port)
gZSpeed	Returns the Z speed (mm/sec) to follow the liquid while aspirating
Syntax:	set zspeed [gZSpeed rackloc zone sample flowrate]
Description:	Return when all axes and syringe operations have completed.
Arguments:	zspeed is the variable that contains the zspeed.
	rackloc is the location of rack on the Gilson bed.
	zone is the sample zone.
	sample is the sample number.
	flowrate is the flowrate for sample (ml/min).
Examples:	set SampleZSpeed [gZSpeed \$RackLoc \$SampleZone \$SampleNumber \$SampleWellRate]
gMix	Mix a Sample in place
Syntax:	
	mixflow

Chapter 10. VAST Accessory Operation

Description:	Mixes a sample in the given sample position by aspirating the sample into the syringe then dispensing back into the sample well. This is done for the volume, flow, and number of times given.
Arguments:	<pre>rackloc is the location of rack on the Gilson bed. zone is the sample zone. sample is the sample number. height is the height above bottom of tube to place needle during mixing operation. mixvolume is the volume that is aspirated and dispense for mixing. mixtimes is the number of time the aspirating and dispensing are done. mixflow is the flowrate for mixing sample (ml/min).</pre>
Examples:	gMix \$RackLoc \$SampleZone \$SampleNumber \$SampleHeight \$MixVolume \$MixTimes \$MixFlow
gTransfer	Transfer sample from one well to another
Syntax:	gTransfer racksrc zonesrc samplesrc heightsrc flowsrc rackdst zonedst sampledst heightdst flowdst
Description:	Transfer the source sample to the destination sample location.
Arguments:	<pre>racksrc is the sample source location of rack. zonesrc is the sample source zone. samplesrc is the sample source number. heightsrc is the source height above bottom of tube to place needle. flowsrc is the source flowrate for sample (ml/min.) rackdst is the sample destination location of rack. zonedst is the sample destination zone. sampledst is the sample destination number. heightdst is the destination height above bottom of tube to place needle. flowdst is the destination flowrate for sample (ml/min).</pre>
Examples:	gTransfer \$XferLoc \$XferZone \$ZferNumber \$XferHeight \$XferFlowRate \$RackLoc \$SampleZone \$SampleNumber \$SampleHeight \$SampleWellRate
gTubeX	Returns the X axis location of the Gilson needle
Syntax:	CurrentXAxisPosition [gTubeX]
Description:	Returns the current X axis position in 0.1 mm. This is the left-to-right direction when facing the Gilson 215 Liquid Handler. The greater the number, the further the needle is to the right of the liquid handler.
Arguments:	CurrentXAxisPosition is set to the current X position in 0.1 mm.
Examples:	CurrentXPosition [gTubeX] - A value of 2310 is 231.0 mm or 23.10 cm.
gTubeY	Returns the Y axis location of the Gilson needle
Syntax:	CurrentYAxisPosition [gTubeY]
Description:	Returns the current X position in 0.1 mm. This is the front-to-back direction when facing the Gilson 215 Liquid Handler. The greater the number, the further the needle is toward the front of the liquid handler.
Arguments:	CurrentYAxisPosition is set to the current Y position in 0.1 mm.

Examples: CurrentYPosition [gTubeY] - A value of 435 is 43.5 mm or 4.35 cm.

Lower needle until liquid detected or limit reached

gMoveZLQ

Syntax: gMoveZLQ \$TRAYBOTTOM

Description: Seeks the liquid level until it detects or reaches a specified Z limit. The greater the number, the futher down the limit. Specifying a limit beyond the tray bottom can drive the needle through the bottom of the tray! The sample requires electrical conductivity for liquid detection to work. Therefore liquid detection will not work for most nonpolar liquids. In addition liquid detection usually requires large sample volumes, greater than ml.

Arguments: TRAYBOTTOM is the Z axis limit in.01 mm.

Examples: gMoveZLQ 1000 - Seeks liquid until Z 100.0 mm is reached.

gMove2LiqLevel Move needle to the liquids level

Syntax: gMove2LiqLevel rackloc zone sample height depth

Description: Moves needle down into sample based on height and depth.

Arguments: rackloc is the location of rack on the Gilson bed.

zone is the sample zone.

sample is the sample number.

height is the height above the tube bottom to initially place the needle. depth is the depth that the needle is extended into the sample. If depth is equal to NOSEEK then no further action is taken. If depth is not equal to NOSEEK then the needle is moved down from the detected liquid level or the height to the depth. The needle for safety reason will not be placed below the tube bottom or above the tube top as specified the Gilson rack definition file.

Examples: gMove2LiqLevel \$RackLoc \$SampleZone \$SampleNumber \$SampleHeight \$SampleDepth

gMove2LiqLevel 1 1 10 30 NOSEEK

This example places needle 3 mm from the bottom of the sample tube.

gMove2LiqLevel 1 1 10 180 100

This example moves needle down attempting to detect liquid level using the Gilson liquid sensor until the liquid is detected or needle is 18-mm above the sample tube bottom. Then the needle is extended down another 10-mm (depth) into the sample.

gRackZoneSequenceOrder Sample sequencing order by zone

Syntax: gRackZoneSequenceOrder rackloc zone startloc pattern

Description: Defines the sample sequencing order for the rack and zone.

Arguments: rackloc is the location of rack on the Gilson bed. zone is the sample zone. startloc is the position of sample number 1 (NW, NE, SW, SE). pattern is the sequencing pattern (HST, HZZ, VST, VZZ). HST is the horizontal straight. HZZ is the horizontal zig-zag. VST is the vertical straight.

VZZ is the vertical zig-zag.

Examples: gRackZoneSequenceOrder \$RackLoc \$SampleZone NW HST

gRackSequenceOrder Sample sequencing order by rack

Syntax: gRackSequenceOrder rackloc startloc pattern

Description: Defines the sample sequencing order for all zones on the rack.

- Arguments: rackloc is the location of rack on the Gilson bed. startloc is the Position of Sample number 1 (NW, NE, SW, SE). pattern is the Sequencing pattern (HST, HZZ, VST, VZZ). HST is the horizontal straight. HZZ is the horizontal zig-zag. VST is the vertical straight. VZZ is the vertical zig-zag.
- Examples: gRackSequenceOrder \$RackLoc NW HST

gRackLocTypeMap Maps the type of rack to a location on the Gilson 215 bed

Syntax: gRackLocTypeMap rackloc racktype

- Description: Maps the type of rack to a location on the Gilson 215 bed. These calls are created by the Gilson program and place in the /vnmr/asm/info/racks file.
- Arguments: rackloc is the location of rack on the Gilson bed. racktype is the Gilson rack type (205, 205H, 505, etc).
- Examples: gRackLocTypeMap 1 205 This example shows that in position one there is a 205 rack.

ResumeAcq Allows acquisition to start while Tcl script still runs

Syntax: ResumeAcq

Description: Sends a signal to console that sample change is complete, this allows acquisition to continue even while the Gilson is still performing other operations. Any operation must be complete by the time the next sample is to be changed.

Otherwise, unpredictable results will occur.

Examples: ResumeAcq

gPuts

Diagnostic output for Tcl scripts

Syntax: gPuts "'string' \$parameter \n"

- Description: Prints function for diagnostic output. Output is inhibited when Tcl parameter debug is zero (default). Output is enabled when Tcl parameter debug is 1. The default value is set within the /vnmr/asm/info/default file and can be changed here to enable output.
 - Examples: gPuts "My parameter value is \$myparameter n" outputs My parameter value is 42.

Chapter 11. PFG Modules Operation

Sections in this chapter:

- 11.1 "Configuring the Software," this page
- 11.2 "PFG Amplifier Operation," on page 372
- 11.3 "Shimming PFG Systems," on page 374
- 11.4 "Setting Up Software for Imaging Pulse Sequences," on page 374
- 11.5 "Homospoil Gradient Type," on page 375
- 11.6 "Gradient Shimming," on page 376

This chapter covers operation of the Varian Performa pulsed field gradient (PFG) modules. The Performa I, Performa II, Performa III, and Performa XYZ PFG modules add new capabilities to high-resolution liquids experiments on most systems. Only the Performa I is available on the *GEMINI* and *GEMINI 2000*. The *MERCUYRY* series spectrometers support all Z axis single gradient options. Single-axis systems apply a gradient in B₀ at programmed parts of the pulse sequence. This gradient can perform several functions, including solvent suppression and coherence pathway selection. Triple-axis systems can also perform solvent suppression and coherence pathway select as well as imaging.

The gradient subsystem produces an intense gradient of up to 65 gauss/cm for a time and then returns to spectroscopic conditions quickly. The relevant parameters are: gradient strength and stability, duty cycle of the gradient system, and time to recover to spectroscopic conditions.

The PFG module installation manual covers installation, calibration, and test procedures for PFG modules.

11.1 Configuring the Software

When the PFG module is installed, values are set in the CONFIG window for PFG by selecting from the choices provided. Only choices appropriate for the system are provided in the CONFIG window. The CONFIG window is opened by entering the config command. The configuration process is described in the manual *VNMR and Solaris Software Installation*.

From the choices made in the CONFIG window, the parameter gradtype is set as a string of three characters (e.g., gradtype='nnp'):

- The first character is the gradient for the X axis, the second for the Y axis, and the third for the Z axis.
- Each axis is 'n' (None choice in CONFIG window), 'w' (WFG+GCU), 'l' (Performa I), 'p' (Performa II/III), 'q' (Performa II/III + WFG), 't' (Performa XYZ), 'u' (Performa XYZ + WFG), 's' (SIS (12 bit), or 'h' (Homospoil). WFG

Chapter 11. PFG Modules Operation

stands for the waveform generation, and GCU stands for the gradient compensation unit. Homospoil is functional only on the Z axis.

Standard PFG operation with the single-axis probe uses the following parameter settings:

- gradtype='nnl' is for Performa I hardware with only the Z channel available.
- gradtype='nnp' is for Performa II hardware with only the Z channel available.

Performa XYZ operation with the triple-axis probe uses the following parameter settings:

- gradtype='ttt' is for Performa XYZ hardware with all three channels available.
- gradtype='nnt' is for Performa XYZ hardware with only the Z channel available.

After the software is configured, the system is ready to use. Table 51 summarizes the commands and parameters associated with PFG.

Command		
creategtable	Generate system gradient table	
grecovery	Eddy current testing	
<pre>setgcoil<(file)></pre>	Assign sysgcoil configuration parameter	
updtgcoil	Update gradient coil	
Parameters		
boresize {number, in cm}	Magnet bore size	
gcal {number, in G/cm-DAC}	Gradient calibration constant	
gcoil {string}	Current gradient coil	
gmax {number, in gauss/cm}	Maximum gradient strength	
gradtype*	Gradients for x, y, and z axes	
gxmax,gymax,gzmax*	Maximum gradient strength for each axis	
pfgon*	PFG amplifiers on/off control	
sysgcoil {string}	System gradient coil	
<pre>trise {number, in sec}</pre>	Gradient rise time	
* gradtype (3-character string from 'n','w','l','p', 'q', 's','t','u','h')		
gxmax, gymax, gzmax {number, in gauss/cm)		
pfgon (3-character string from 'n	' and 'y')	

Table 51. Pu	ilsed Field	Gradients	Commands	and Parameters
--------------	-------------	-----------	----------	----------------

11.2 PFG Amplifier Operation

The PFG current amplifier is left on during the experiment, eliminating the need for current blanking. Although the PFG amplifiers are quiet, they do produce a small amount of quiescent current, resulting in a shift of the Z1 shimming. The Performa XYZ PFG amplifier produces a small amount of quiescent current into each of the X, Y, and Z gradient coils, resulting in a change to the X1, Y1, and Z1 shimming.

Although leaving the amplifier on to have the shimming stable is a reasonable procedure, it may be disabled by the global string parameter pfgon to check for noise sources or to change probes. pfgon is a three-character string with the first character controlling the X channel, the second controlling the Y channel, and the third controlling the Z channel. The value for each channel is y or n (e.g., pfgon='nny' turns on the PFG amplifier on the Z channel only). A su or go command must be sent to activate pfgon.

It is useful to translate the gradient control in DAC units to G/cm by a constant that represents G/cm-DAC units. The parameter gcal, a user global real-valued parameter, makes this translation.

The amplifier system is well-behaved during power up and other exceptional conditions. For example, if the console is on and the current amplifier is subsequently switched on, the unit warms up but its output is disabled by the internal logic. During the setup, the interface transmits a zero to the unit. The output is enabled whenever a su or go operation is requested with pfgon='nny'. The enable logic in the current driver is edge-triggered so that the output current is not enabled unless an explicit command is sent.

If the amplifier was on and the console power cycled or rebooted, the AP bus control of the enable is off, disabling the current. An su then restores normal functioning. If an experiment requiring the gradient is started without the su command, the experiment fails.

As an extra safety precaution with Performa II systems, pressing the reset button on the interface board zeroes the main DAC. Similarly, a reset aborts zero but does not disable the output current. Using pfgon='nnn', or setting the amplifier standby switch disables amplifier output.

PFG amplifiers have a series of lights to enable you to diagnose operation:

• Performa I amplifier lights:

POWER	Power is on.
ON	Amplifier is active.
ACK	AP response.
DATA	Flickers when data is sent to it, indicating the interface to the amplifier from the system is working.
ERR	Duty cycle exceeded or internal overload.
LOAD	Amplifier sees a short or open circuit. Check probe and probe connections.

• Performa II and Performa III amplifier lights:

POWER	Power is on.
ON	Amplifier is active.
DATA	Flickers when data is sent to it, indicating the interface to the amplifier from the system is working.
HITEMP	fault: internal thermal overload.
ERROR	Duty cycle exceeded.
WARMUP	Warming up (operation not recommended).
LOAD	Amplifier sees a short or open circuit. Check probe and probe connections.

• Performa XYZ amplifier lights:

POWER	Power is on.
ON	Amplifier is enabled.
DATA	Flickers when data is sent to it, indicating the interface to the amplifier from the system is working.
PULSE	Gradient pulse is active.
ERROR	Duty cycle exceeded, open load, or blown coil fuse.

11.3 Shimming PFG Systems

The procedures in this apply to the Performa I, Performa II, and the Performa XYZ systems. Once in operation, leave the amplifier on while using the gradient system, to allow the amplifier to reach a long-term equilibrium.

Performa I and Performa II

- 1. Enter **pfgon='nnn' su** to turn off the amplifier. Verify a drop in the lock level from the small dc zero current from the amplifier.
- 2. Shim the system to the desired level.
- 3. Enter **pfgon='nny' su** to turn on the amplifier. The shimming changes from the small dc offset current.
- 4. Adjust Z1 to restore the homogeneity. The lock level should have identical stability on the meter.

This two-stage approach is not strictly necessary, but it does separate any problems that might arise.

Performa XYZ

- 1. Prepare the amplifier by moving the switch from STANDBY position to ON.
- 2. Enable the amplifier by entering **pfgon='yyy' gradtype='ttt' su** (or **gradtype='uuu'** for Performa XYZ with waveshaping). The yellow RUN lights turns on. Shim the system to the desired level.

11.4 Setting Up Software for Imaging Pulse Sequences

After the gradient configuration has been selected, to make use of the imaging, obliquing, and variable angle pulse sequences and pulse sequence statements, calibrate the gradients, create a gradient table, and set the configuration parameter sysgcoil to the desired gradient coil. Each of these steps is explained in the next sections.

Calibrating the Gradients

Follow the steps in the "PFG Test Procedure" section in the PFG installation manual to calibrate the gradients. Calibrate each axis and save the values for the maximum gradient strength in gauss/cm. These values will be needed when creating a gradient table.

Creating a Gradient Table

A gradient table should be created for each gradient coil/amplifier combination. Run the creategtable macro, then enter the following information:

- Maximum gradient strength of each axis
- Gradient rise time
- Gradient bore size

Setting the System Gradient Coil

Once a gradient table has been created, it can then be configured in the system. Any number of gradient tables can be created, although the config program can only handle 30. Configuring a table in the system means setting the system gradient coil parameter, sysgcoil to the name of the gradient table. To set sysgcoil, use the CONFIG window or the setgcoil macro.

- In the CONFIG window, select the Gradient window to initially configure the gradients. Notice the field for Imaging Gradient Coil. Select the desired table from the list of choices. If the name does not exist, it means a table has not been created. For a list of configuration items, see the manual *VNMR and Solaris Software Installation*.
- To configure a gradient table using setgcoil, enter setgcoil(file), where file is the file name of the table (e.g. setgcoil('tbl')). This updates the values of sysgcoil and gcoil as well as related parameters in the current experiment.

After a system gradient coil has been configured, any experiments that have a gcoil parameter will automatically be updated whenever that experiment is "joined" or new parameters are retrieved. When the gcoil parameter is updated, all the parameters corresponding to the entries in the gradient table file (gxmax, gymax, gzmax, trise, and boresize) are also updated. These parameters are then used by the pulse sequence. gxmax, gymax, and gzmax, in particular, are used by the obliquing PSG elements to convert gauss/cm to DAC units. To create and set a gcoil parameter for an experiment parameter set, type updtgcoil.

11.5 Homospoil Gradient Type

Starting with VNMR 6.1A, it is possible to use homospoil (room temperature Z1 shim coil) as a general gradient type. It does not require the use of a pulsed field gradient module and thus is available on systems without PFG. Homospoil gradients are implemented only on the Z axis.

When homospoil is switched on in a pulse sequence, the shim current is set to maximum for a given period of time. Homospoil control within a pulse sequence is done in the following manner:

- To use homospoil as a quick homogeneity spoil, use hsdelay. This is the traditional homospoil method, and is usually done at the beginning of a relaxation recovery delay (e.g., hsdelay(dl)). The parameter gradtype is ignored. See the *User Programming manual* for details of how to use hsdelay.
- To use homospoil as a general gradient type, first select the homospoil gradient type. Enter **config** and under **Gradients** select **Homospoil** (this sets gradtype='nnh'). The parameter pfgon is ignored, since a separate gradient amplifier is not needed. Homospoil is then triggered by gradient statements such as rgradient('z',gzlvl1). If the value of gzlvl1 is non-zero, homospoil is switched on; if the value of gzlvl1 is zero, homospoil is switched off. Only one sign and strength of gradient current is available during a pulse sequence, and is set by hardware.

Homospoil gradients may be switched on only for a limited period of time, usually 20 ms. This time limit is determined by hardware in UNITY*INOVA* and UNITY*plus* systems (see Table 52 for system configurations). Check your pulse sequences to ensure this time limit is not exceeded.

System	Shim Supply	Homospoil Time Limit
UNITYINOVA	Varian 14	20 ms/200 ms ^a
UNITYINOVA	Varian 18 to 40	20 ms/200 ms ^{a,b}
UNITYINOVA	RRI Ultrashims	20 ms/200 ms ^a
UNITY <i>plus</i>	Varian 13	20 ms/200 ms ^c
UNITY <i>plus</i>	Varian 18 to 40	20 ms/200ms ^{a,b}
UNITY <i>plus</i>	RRI Ultrashims	20 ms/200 ms ^a
UNITY/VXR-S	Oxford 18	Weak homospoil ^d
MERCURY-VX, MERCURY	Varian 14	No time limit ^e
GEMINI 2000	Varian 13	No homospoil

a. Hardware upgrade to 200 ms with the Automated Deuterium Gradient Shimming module is required for compatibility with ²H gradient autoshimming.

b. Hardware adjustment required for both ¹H and ²H gradient autoshimming. Adjust homospoil potentiometer resistor labeled HOMO (blue square) on front of Z0/Z1 board to maximum in either direction for maximum homospoil gradient strength.

c. Hardware upgrade required for both ¹H and ²H gradient autoshimming. Homospoil gradient strength is too weak for autoshimming without hardware upgrade.

d. Hardware upgrade not available. Homospoil gradient strength is too weak for gradient autoshimming.

e. A homospoil time limit of 20 ms is set by software for hsdelay.

The behavior of homospoil gradients is quite different from that of a pulsed field gradient. The gradient strength is much weaker than the traditional PFG, and the recovery time is much longer because of eddy currents. The strength and recovery of the gradient depends on the shim coils and system hardware. Typically, these gradients are suitable only for profile-type experiments and unsuitable for gradient coherence-selection experiments such as GCOSY and GNOESY. For all gradient experiments, pulsed field gradients are preferred if available.

Homospoil gradients are suitable for ¹H and ²H gradient shimming on some systems (see Table 52 for system configurations). The Automated Deuterium Gradient Shimming module (see *Getting Started* and *Accessories Installation* manuals) is required on UNITY*INOVA* and UNITY*plus* systems to upgrade the homospoil hardware for compatibility with deuterium gradient shimming.

11.6 Gradient Shimming

Gradient shimming provides rapid, automatic shim adjustment and applies to all systems with PFG or homospoil and gradient shimming software installed.

Proton gradient autoshimming with PFG is available on all systems configured with a PFG accessory. Deuterium gradient shimming is only available on UNITYINOVA, UNITYplus, *MERCURY-VX*, and *MERCURY* systems and not on *GEMINI 2000*, UNITY, and VXR-S systems because lock sample and hold capability is required to perform deuterium gradient

shimming. A hardware upgrade is also required to perform deuterium gradient shimming. System configuration requirements are summarized in Table 53.

System	Gradient Shim Availability
UNITY INOVA, UNITY plus	¹ H or ² H with PFG or Homospoil ^a
MERCURY-VX, MERCURY	¹ H or ² H with PFG or Homospoil ^b
UNITY, VXR-S, GEMINI 2000	¹ H with PFG only

 Table 53. Gradient Shim Availability

a Automated deuterium gradient shimming module required for deuterium gradient shimming with PFG or homospoil.

b Automated liquids/solids spinner controller and automated deuterium gradient shimming module required for deuterium gradient shimming.

Configuring Gradients and Hardware Control

- 1. Confirm that PFG or homospoil gradients are installed on your system. See the previous sections in this chapter.
- 2. Confirm that the gradients are active by checking that gradtype and pfgon are set appropriately for your system. Use config to change settings if necessary.
- 3. If you have the Ultra•nmr shim system, enable control of the shims from the Acquisition window, as described in the section "Shimming Using the Ultra•nmr Shim System" in the manual *Getting Started*.

Gradient Shimming Method

The full gradient shimming method consists of these steps:

- 1. Map the shims.
- 2. Perform autoshimming.

The shims must be mapped before autoshimming is used. Mapping the shims is necessary when a new probe is installed, but can be repeated at any time.

Table 54 summarizes gradient shimming commands and parameters.

CAUTION: Spinning the sample during gradient shimming can cause motion artifacts.

Mapping the Shims

Mapping the shims is necessary after installing a new probe. 90% H_2O is recommended for first time shimming on ¹H and 1% H_2O /99% D_2O is recommended for ²H.

- 1. Insert a sample and find lock.
- 2. Stop sample spinning (**spin=0**). Disable sample changer control (**loc='n'**).
- 3. Adjust lock power, lock gain, and lock phase. Make coarse shim adjustments on Z1, Z2, X1, and Y1.
- 4. Use **s2pul** to find the 90° pulse for tn='H1' or tn='lk', as appropriate.

Commands	
dg2	Display group of parameters.
<pre>gmapshim<('files' 'quit')></pre>	Run gradient autoshimming, get files and parameters, quit.
gmapsys*	Enter Gradient Shimming System menu, make shimmap.
gmapz<(mapname)>	Get parameters/files for gmapz pulse sequence.
* gmapsys<'shimmap'<,'auto'	'manual' 'overwrite' mapname>
Parameters	
d2	Incremented delay for 1st indirectly detected dimension.
d3	Incremented delay for 2nd indirectly detected dimension; arrayed to two values
gradtype*	Gradients for x, y, and z axes
gzlvl {DAC value}	Pulsed field gradient strength
gzsize {integer,1 to 8}	Number of z-axis shims used by gradient shimming
gzwin {0 to 100}	Percentage of spectral window used by gradient shimming
pl	First pulse width—If > 0 , it is used between the gradient pulses as a 180 refocusing pulse, and the gradients have the same sign.
pfgon{'nny' if on}	PFG amplifiers on/off control
wq	Pulse width; it can be <90° if p1=0.
* gradtype {3-character string fro	m'n','w','l','p','q','s','t','u','h'}

Table 54.	Gradient	Shimming	Commands	and	Parameters
-----------	----------	----------	----------	-----	------------

5. Enter **gmapsys** to display the Gradient Shimming System menu.

Standard parameters are retrieved from gmapz.par the first time gmapsys is entered, or if a shimmap was previously made, parameters are retrieved from the current shimmap. If desired, enter **gmapz** to retrieve standard parameters from gmapz.par.

- To set parameters for a particular gradient and nucleus, enter gmapsys and click on Set Params > Gradient, Nucleus, and then click on the appropriate button (Pfg H1, Pfg H2, Homospoil H1, or Homospoil H2). Next, set pw as follows:
 - For PFG, set **pw** to the 90° pulse or less.
 - For homospoil, set **pw** to 90° pulse and p1 to 180° pulse.
- 7. Test the parameters. Enter **gmapsys** and click on **Set Params** > **Go**, dssh.

You should see two profile spectra. If you don't, check that the gradients are active and check pw, tpwr, and gain.

8. To make a shimmap, again enter **gmapsys** and click on **Shim Maps** > **Automake Shimmap**. Enter a mapname (any string valid for a file name) at the prompt.

Starting Gradient Shimming

To start shimming as a system administrator, enter **gmapsys** and then click on **Autoshim on Z**. This button starts gradient shimming using current parameters, and displays the curve fit and shim adjustments for each iteration.

Quitting the Gradient Shimming System Menu

Enter **gmapsys** and click on **Quit** to exit from the gmapsys menu system. This also retrieves the previous parameter set and data, including any data processing done on the previous data set.

General User Gradient Shimming

For the general user, gradient shimming can be run from outside gmapsys from any experiment. Any one of the following methods is recommended for routine use:

- Click on **Main Menu** > **Setup** > **Shim** > **Gradient Autoshim on Z**. Parameters are retrieved from the current mapname, which is displayed at the start of shimming, and the spinner is automatically turned off. The curve fit and shim adjustments are not displayed. The previous parameter set and data are retrieved when shimming is finished. This button only functions after a shimmap is made.
- Enter gmapshim. This performs the same action as clicking on Gradient Autoshim on Z.
- Within automation parameter sets, use **wshim='g'** (UNITY*INOVA, MERCURY-VX*, UNITY*plus*, and *MERCURY* only).

To stop gradient shimming before it is completed, use one of the following methods:

- Click on Main Menu > Setup > Shim > Quit Gradient Shim. Quitting aborts the experiment and retrieves the previous parameter set and data.
- Enter gmapshim('quit'). This performs the same action as Quit Gradient Shim.
- Abort the acquisition with aa and click on Cancel Cmd. Then enter gmapshim('quit') to retrieve previous data set and parameters.

How Gradient Shimming Works

The basis of gradient shimming is differential phase accumulation from shim gradients during an arrayed delay. The phase is spatially encoded by a pulsed field gradient. Figure 107 shows the gradient shimming pulse sequence.




Figure 107. Gradient Shimming Pulse Sequence

The gradient shimming pulse sequence in Figure 107 is shown with p1=0, in which case pw can be set to a small flip angle. If p1>0, the pulse field gradients are both set to the same sign, and p1 should be set to $180\square^\circ$ and pw to 90° , so that rf inhomogeneities are refocused. p1=0 is usually sufficient for most cases.

Phase accumulation from all gradients present is as follows:

 $\phi = z Gz(-at/2 + t) + dG(d3 + at/2 + 3*d2 + t)$

where t is the time during acquisition at, Gz is the z-axis pulsed field gradient strength, and dG is the sum of the shim gradient fields, shown as being on during relevant times in the pulse sequence.

The effect of the shim gradients dG can be isolated by arraying d3 and taking the difference in the phases:

 $\Delta \phi = \phi 2 - \phi 1 = dG^*(d3[2] - d3[1])$

For example, at a particular point, $\Delta \phi$ can be $2\pi * 100$ Hz * 10 ms, or 2π radians. Thus, a pair of profiles with different d3 values can be used to calculate the B₀ field along z.

The effect of any one shim gradient can be isolated by arraying the shim value, represented by dG, and taking the difference in the phase differences:

 $\begin{aligned} \Delta(\Delta \phi) &= \Delta \phi 2 - \Delta \phi 1 = dG2^*(d3[2] - d3[1]) - dG1^*(d3[2] - d3[1]) \\ &= (dG2 - dG1)^*(d3[2] - d3[1]) \end{aligned}$

Therefore, two pairs of profiles can be used to map out the effect of a shim. By arraying all the shim values, a set of phase difference maps or shim field maps can be constructed for a given shim set. Shimming can then be performed by constructing a background field map for the starting shim values ($\Delta \phi$) and fitting the result to the shim field maps. The calculations are quite fast, so the entire shimming process is usually limited by the data acquisition time, typically taking only a few minutes.

In practice, the phase is calculated from $\phi = \arctan(x, y)$ from the real and imaginary values at each point in the spectrum, and $\Delta \phi$ is calculated from the difference in the phases of a pair of spectra with d3 arrayed. Figure 108 shows an example of mapping the z1 shim.




Figure 108. Mapping the z1 Shim

References

Van Zijl, P. C. M., et al. J. Magn. Reson. 1994, 111 (Series A), 203–207.
Sukumar, S., et al. J. Magn. Reson. 1997, 125 (Series A), 159–162.
Barjat, H., et al., J. Magn. Reson. 1997, 125 (Series A), 197–201.

How Making a Shimmap Works

Automake Shimmap first runs an experiment that calibrates gzwin and tof to set the spectral window. Next, it runs an experiment with the shims arrayed to map the shims, and processes the experiment when done. Coarse shims are used if present. The parameters and data for the shimmap are stored in the file userdir + '/gshimlib/shimmaps/' + mapname + '.fid'. These parameters are retrieved the next time gradient shimming is run if the gradient shimming system menu is exited.

Displaying the Shimmap

After the shims are mapped, display the shimmap by entering **gmapsys** and clicking on **Display > Display Shimmap**. The shimmap display is a multicolored plot of the shimmap, with Z1 as #1 and Z2 as #2, and so on (see Figure 109).




Figure 109. Shimmap Plot

The shimmap is specific to the probe used, and can also be dependent on sample volume for small volumes. The shimmap shows the actual field dependence of the shims, except for a dc offset added for display purposes. Good signal-to-noise in the shimmap is needed for the shimming to work well. Poor signal-to-noise might result in incorrectly set shims.

Calibrating gzwin (optional)

The parameter gzwin is the percentage of the spectral window used in calculating the field maps. gzwin should be adjusted only when making a new shimmap. If this parameter is not calibrated correctly, you may see excess noise data at the edge of the shimmaps, which corresponds to the region in the profile spectrum where the signal goes to zero. It is normal to have a few noise data points at the edge of the shimmap, but if it is more than a few data points (greater than 25% of the window), gzwin may be miscalibrated. This can occur if there is low signal-to-noise or if gzwin has not previously been calibrated for the current parameter set. If the gain is too high, "wings" will appear on the sides of the spectra and

Chapter 11. PFG Modules Operation

may result in miscalibrated gzwin. This can also occur if there are multiple chemical shifts in the presence of a weak gradient.

Automatic Calibration of gzwin

Clicking the Automake Shimmap button causes the following actions:

- 1. Enter **gmapsys** and click on **Set Params** > **Find gzwin**, which calibrates gzwin and sets tof to center the window used for calculation.
- 2. Click on **Return** > **Shim Maps** > **Make Shimmap**, which makes the shimmap with the current values of gzwin and tof.

You may click through these steps separately to see if gzwin is calibrated correctly. The box cursors at the end of step 1 should be at either edge of the profile.

Manual Calibration of gzwin

Manual calibration of gzwin may be used to avoid noise spikes in the spectrum, or other artifacts. To manually calibrate gzwin, do the following:

- 1. Click on **Set Params** > **Go, dssh**. Wait until the experiment is done.
- 2. Enter **ds**, and set the box cursors near the edges of the profile.
- 3. Enter gmapsys, and click on Set Params > Calculate gzwin.
- 4. Click on Return > Shim Maps > Make Shimmap.

The parameter gzwin should be adjusted only when making a new shimmap. The calibrated value of gzwin is saved when the new shimmap is saved at the end of the mapping experiment. The same value of gzwin must be used in shimming as in making a shimmap, and should not be adjusted when shimming.

Shimmap Files and Parameters

The parameters and shimmap files saved under a mapname are retrieved when that mapname is retrieved. When reinserting a probe, reload the shimmap for that probe. If you are unsure if the shimmap is correct, make a new shimmap, which typically only takes a few minutes. The last parameters and files used are automatically retrieved the first time gmapsys is entered. If gmapsys is entered again, the parameters are not retrieved. Gradient shimming uses the current parameters after the pulse sequence is loaded (seqfil='gmapz').

Standard parameters can be loaded before making a shimmap by entering gmapz or by using the **Gradient,Nucleus** menu button. Parameters and files can also be explicitly loaded and distributed, as described in the following subsections:

Loading a Shimmap

To change shimmaps as a system administrator, do the following:

- 1. Enter gmapsys.
- 2. Click on Shim Maps > Shimmap Files > Cd to Userdir.
- 3. Select a file.
- 4. Click on **Load Shimmap** (loads the shimmap files gshim.list and gshim.bas from gshimlib/shimmaps/mapname.fid into gshimlib/data, but does

not load the parameters) or click on **Load Shimmap & Params** (loads shimmap files and parameters).

The general user can also change shimmaps by entering gmapshim('files') from any experiment to display the Gradient Autoshimming Files menu, and then clicking on either the Load Shimmap or Load Shimmap & Params buttons.

Distributing a Shimmap

The system administrator can copy a shimmap file from vnmrsys/gshimlib/ shimmaps into the directory /vnmr/gshimlib/shimmaps so that the file is accessible to all users. To copy files, do the following steps:

- 1. Log in as **vnmr1**.
- 2. Enter gmapsys.
- 3. Click on Shim Maps > Shimmap Files > cd to Userdir.
- 4. Select a file.
- 5. Click on Copy to Systemdir.
- 6. Become the new user to be given access to gradient shimming.
- Enter gmapsys.
- 8. Click on Shim Maps > Shimmap Files > cd to Systemdir.
- 9. Select a file.
- 10. Click on Load Shimmap.

How Automated Shimming Works

The shims must be mapped before gradient autoshimming is used. See "Mapping the Shims," on page 377 for details.

When gradient shimming is run from the **gmapsys** menu, the curve fit plot is displayed for each iteration. The plot shows the raw data as #1 and the curve fit as #2 (see Figure 110).

Shim adjustments for each iteration are also displayed (see Figure 111) and have converged when the rms error number is less than 1.0. Gradient shimming continues until convergence or until a maximum of 10 iterations are reached.

If a shim goes out of range, the shim is set to maximum and shimming continues with the remaining shims. If convergence is then reached, shimming is tried once more with all Z shims and continues unless a shim goes out of range again.

If the parameter gmap_zlz4 is set to 'y', then if gzsize is greater than 4, shimming is done first on Z1–Z4 and then proceeds with all shims specified by gzsize. Gradient shimming takes longer and goes through more iterations, but this may avoid the problem on some systems where a high-order shim (i.e. Z5, Z6) goes out of range because it contains impurities from lower-order shims. This parameter may be set at any time while shimming from gmapsys. In order to use this parameter in user autoshimming, set it before making a shimmap, or in the corresponding parameter set in gshimlib/shimmaps.

Deuterium Gradient Shimming

Deuterium gradient shimming is feasible for most deuterated solvents for which lock solvent has a single, strong deuterium resonance with sufficient signal.




Figure 110. Curve Fit Plot

-	e 5mm_Tria z 4 gzsiz	—		rms err	1.892
Shim	Offset	Old	New	Diff	Error
z1	800	-9405	-9269	-136	48
z2	800	-3118	-3104	-14	13
z3	3200	-4356	-4321	-35	37
z4	-3200	4049	4885	-836	104
z5	-3200	13443	14537	-1094	322
zб	3200	-15619	-12568	-3051	467
z7	3200	0	0	0	0
z8	3200	0	0	0	0

Figure 111. Display of Shim Adjustments for Each Iteration

The automated deuterium gradient shimming module is required to run deuterium gradient shimming. If present, this module automatically holds the lock at its current value and switches the transmitter cable to pulse the lock coil when an experiment is run with tn='lk'. The module is strongly recommended for all users who wish to run deuterium gradient shimming in automation. See Table 53 for system configurations supported.

The system administrator must make a shimmap on deuterium before deuterium gradient shimming can be used. Follow the procedure "Mapping the Shims," on page 377, using the deuterium signal for all steps, (use tn='lk' for step 4, and select **Pfg H2** at step 6). The transmitter power (tpwr) should be kept low to avoid probe arcing, with a 90° pulse greater than about 200 µs.

The recommended parameters for different solvents are shown in Table 55.

The deuterium parameters are saved for future use when the shimmap is saved, and are used the next time gradient autoshimming is run.

		Gain		
Solvent	nt	dl (sec)	Inova	Mercury
deuterochloroform	8-32	3	36	18
dmso-d6	4-16	3	28	10
D ₂ O	1-4	3	24	6
deuterobenzene	1-4	3	24	6
deuteroacetone	1-4	6-12	24	6

 Table 55.
 Deuterium Parameters

Note: Actual parameters might vary, depending on solvent concentration, probe, and system hardware.

Homospoil Gradient Shimming for ¹H or ²H

It is also possible to use the Z1 room temperature shim as a homospoil gradient, instead of using a pulsed field gradient or PFG. Use of this option is recommended only if a PFG amplifier or probe is not available. For details on how to configure a homospoil gradient, see the section "Homospoil Gradient Type" in the chapter on PFG modules operation in the manual *User Guide: Liquids NMR*. The system administrator must make a shimmap using homospoil before homospoil gradient shimming can be used. Follow the procedure in "Mapping the Shims," on page 377.

At step 4 select tn as appropriate, set tof on resonance, and then find the 90° pulse. If deuterium is used, tpwr should be kept low, with a 90° pulse greater than about 200 µs.

At step 6, select either **Homospoil H1** for proton parameters or **Homospoil H2** for deuterium parameters, as appropriate. Homospoil gradients must be configured at this step (use config or set gradtype='nnh').

At step 7, further testing of the gradient shimming parameters for homospoil can be done as follows:

- 1. Enter gmapsys, and click on Set Params > Go,dssh.
- Using the 90° pulse from step 4, calibrate the 90° and 180° pulses to obtain an echo. Enter df to display the FID. You should see an echo forming in the middle of the FID.
- 3. If needed, adjust **sw** so that the gradient covers at least 10% of the spectral window. Increase **np** to 512 to improve Hz/point resolution. However, **np** should be adjusted so that at is not longer than the homospoil time limit (20 ms on standard UNITY *INOVA* and UNITY *plus* and 200 ms with the automated deuterium gradient shimming module). The acquisition time (at) should also be shorter than T_2 . Set **d2=at/4**.

Once all the parameters are set, click on **Automake Shimmap** (step 8). The parameters are saved when the shimmap is done and are used the next time gradient autoshimming is run.

To use homospoil deuterium gradient shimming with different solvents, set the parameter gmap_findtof='y'. This should be done before making the shimmap, or may be set in the corresponding parameter set in gshimlib/shimmaps. Gradient shimming will then perform a calibration to find tof before autoshimming starts. This takes an extra iteration in

the shimming procedure, but is essential for homospoil deuterium gradient shimming to function correctly in automation.

Full Deuterium Gradient Shimming Procedure for Lineshape

The automated deuterium gradient shimming module must be installed to use this procedure.

- 1. Insert the appropriate lineshape sample (chloroform in acetone-d6) and find lock. Turn off spinning and disable sample changer control. Adjust lock power, lock gain, and lock phase as necessary. Do quick shimming on z1, z2, x1, y1 (use z1c, z2c, if present).
- 2. Find the 90° pulse on 2 H as follows:
 - a. Enter s2pul tn='lk' tpwr=42 pw=200.
 - b. Enter **ga** and wait for acquisition to finish. You should see only a single line. Place the cursor on resonance and enter **nl movetof**.
 - c. Enter **array('pw', 20, 100, 100)**. Enter **ga** and wait for acquisition to complete.
 - d. Set pw90 to the first maximum.
- 3. Set up gradient shim parameters. Enter **gmapsys** and click on **Set Params** > **Gradient**, **Nucleus**.
 - a. If you have PFG, click on **Pfg H2**, set **pw** to one half the 90° pulse found in step 2 and set **d1=6**.
 - b. If you have homospoil, click on **Homospoil H2**, set **pw** to the 90° pulse, **p1** to the 180° pulse, and **d1=6**.
- 4. Enter **gmapsys**, click on **Set Params** > **Go**, **dssh**, and wait for acquisition to complete. You should see two profile spectra.
- 5. Map the shims. Enter **gmapsys** and click on **Shim Maps** > **Automake Shimmap**. At the mapname prompt, enter a file name. Wait for acquisition to complete and the message to be displayed: shimmap done!
- Perform shimming on z's. Enter gmapsys and click on Autoshim on Z. Wait for the acquisitions to complete, and the message to be displayed: Gradient Autoshimming on Z done! N iterations.
- 7. Enter **acqi** and adjust lock phase. Shim only on low-order nonspins (x1, y1, xz, yz, etc.). Do not shim on z's (z1, z2, etc.).
- 8. Perform shimming on z's. Repeat step 6.
- 9. Measure proton lineshape. Turn on spinner if appropriate, and enter **acqi** to make fine shim adjustment. Shim on all shims as necessary.

Setting Up Automation

You must have the deuterium gradient shimming module installed to perform deuterium gradient shimming in automation.

- 1. Find your 90° pulse on deuterium:
 - Enter s2pul tn='lk' tpwr=42 pw=200 gain=20.
 - Enter ga.

- Array pw to find the 90° pulse.
- 2. Set the parameters:
 - Enter gmapsys.
 - Click on **Set Params**.
 - Click on Gradient, Nucleus.
 - If you have pulsed field gradient, click on **Pfg H2** to select the parameter for deuterium with pfg. If you have homospoil gradients, select **Homospoil H2** instead.
 - Set tpwr and pw for 90° pulse and set gain as appropriate for your solvent.
 - Click on **Return**, then **Go**, **dssh** to test the parameters.

You should see two top-hat profile spectra. Adjust gain and nt to see good signal-to-noise with no ADC overflow.

The recommended parameters for different solvents are listed in Table 55.

- 3. Make a shimmap for a particular solvent using parameters for good signal-to-noise. If the solvent you most often use has a weak signal, make the map on a solvent with more signal. To make the map, click on **Shim Maps**, then **Automake Shimmap**. When prompted, enter a map name.
- 4. Test autoshimming by entering gmapsys, then click on Autoshim on Z.
- 5. When you are satisfied that autoshiming works well for your particular solvent, open a shell window and edit the gmapz macro. Go to the bottom of the macro and uncomment the section appropriate for your solvent (the solvent on which the shimmap was made). Use parameters as in step 2. Repeat step 2 for all solvents of interest.
- 6. Edit the setlk macro and uncomment the section setting lockpower, lockgain, and z0 for different solvents (you must first determine the appropriate values for different solvents).
- 7. Change the appropriate parameter sets in /vnmr/stdpar to use wshim='g' spin='y' spin=20 (or another appropriate spin value). Use alock='y' if you use the setlk macro.
- *Note:* Use config to set the gradient type to Homospoil or PFG for use in automation. Setting gradtype from the command line does not set the gradient type for automation.

Suggestions for Improving Results

Calibrate the 90° pulse and adjust tpwr, pw, and gain to optimize signal-to-noise. Reduce gain if ADC overflow occurs, which may appear as wings on the profile.

For solvents with long $T_{1,}$ set **d1** to 3 to 5 times T_{1} , or use a small flip angle for pw. Stimulated echoes may otherwise result, which may appear as excess noise or a beat pattern in the spectrum, or as secondary echoes in the FID (use df to observe this).

The phase encode delay d3 is arrayed to two values, the first of which is zero. The second value can be increased for better signal-to-noise in the phase maps, up to about the point where the amplitude of the second profile is half that of the first (about $2/3 T_2$ without radiation damping; radiation damping can be severe in water ¹H). However, longer d3 values increase the phase excursion, and can make it difficult to shim large shim corrections (especially Z1). Typical ¹H values are 5 to 30 ms, and typical ²H values are 30 to 200 ms.

Chapter 11. PFG Modules Operation

If the shims are far off when making a shimmap, the second value of d3 might be too small. If this problem occurs, decrease the second value of d3 to temporarily one-half to onequarter its value.

When reinstalling a probe, make sure it is in the same vertical position in the magnet barrel as when the shimmap was made. If you are unsure, make a new shimmap, which typically takes only a few minutes.

Alternate between z-axis gradient shimming and shimming the low-order x- and y-axis shims by other methods (e.g., on lock level). The z-axis shims account for the majority of sample volume changes (changes in height), and the x- and y-shims are relatively insensitive to change in height. Evaluate shimming for a particular application, since the ideal lineshape may vary with the application.

The high-order shims can sometimes be set off-scale during shimming. This may occur if the sample is short, or if the sample is improperly seated in the probe, or if the high-order shims are weak or other effects. In such cases, the off-scale shim is set to maximum, and shimming continues with lower-order shims. Superior results can be obtained in some cases by first setting **gzsize=4** and clicking on Autoshim on Z to shim on z1-z4, and then shimming the low-order transverse shims, and then increasing gzsize and shimming again. This may also be done using $gmap_z1z4='y'$. On a short sample it also can be useful to remap the shims.

Some shim systems may need additional time when running the shim mapping experiment to allow the shims to settle. The added time is especially noticeable on some systems for Z4. To account for added time, lengthen the d1 delay or add dummy scans in between each array element (e.g., ss=-2). Decreasing the amount a shim is offset also allows the shim to settle more quickly. Enter gmapsys('vi') to edit the values in the Offset column, and then enter gmapsys('shimmap', 'manual') to map the shims with user-defined offsets. A new mapname may also be set using gmapsys('vi').

Coarse shims are used on systems on which they are available. To use fine shims on these systems, enter gmapsys('vi') to edit the entries in the shim column (e.g., change zlc to zl), and then enter gmapsys('shimmap', 'manual') to map the shims.

For samples in H_20 , the water protons provide sufficient signal for shimming. For samples other than water, deuterium gradient shimming is strongly recommended if there is sufficient deuterium signal. Proton gradient shimming can be made to work in samples other than water if there is sufficient proton signal and the signal is well-resolved (does not overlap with other strong resonances). Gradient shimming can also be done on a water sample of equal height of the sample of interest, and then the sample of interest can be inserted.

For further information, refer to the entries for gmapshim, gmapsys, and gmapz in the *VNMR Command and Parameter Reference*.

Gradient Shimming Menus

The top-level gradient shimming menu is the Gradient Shimming System menu, which is opened by entering gmapsys.

Set Params Shim Maps Autoshim on Z Set Shims Display Quit

From this menu and submenus, there are buttons to open the following menus:

Gradient Shimming Setup menu

Go, dssh Gradient, Nucleus Find gzwin Calculate gzwin Return

• Gradient Nucleus Parameter Setup menu

Pfg H1 Pfg H2 Homospoil H1 Homospoil H2 Return

• Gradient Shimming Map menu

Automake Shimmap Make Shimmap Shimmap Files Current Mapname Save As Return

- Gradient Shimming Files menu
 Cd to Userdir Copy to Userdir Load Shimmap Load Params Rename Return
- Gradient Shim Setting menu

Old Shims New Shims Min Shims Return

• Gradient Shimming Display menu

dgs List Shims Display Shimmap Display Fit Show Record Plot Return

• Gradient Shimming Plot menu

Print dgs Print Shims Plot Shimmap Plot Fit Plot b0 Return

Two other gradient shimming menus are available for the general user:

• Automated Shimming menu. To open, click on Main Menu | Setup | Shims.

Lock Autoshim Lock Autoshim z1z2 Gradient Autoshim on Z Return

• Gradient Autoshimming Files menu. To open, enter gmapshim('files').

Cd to Userdir Copy to Userdir Load Shimmap Load Params Rename Return

A description of each menu (in the same order as above), the button or command for opening the menu, and the action of each button in the menu is given in the following sections. On all menus, the left mouse button is used to click on the desired choice in the menu. The other mouse buttons are not active.

Gradient Shimming System Menu

The Gradient Shimming System menu is the top-level menu for gradient shimming setup and shim mapping. To open this menu, enter gmapsys.

Button	Description
Set Params	Enter Gradient Shimming Setup menu.
Shim Maps	Enter Gradient Shimming Map menu.
Autoshim on Z	Start gradient shimming using current parameters. Fit is displayed with each iteration.
Set Shims	Enter Gradient Shim Setting menu.
Display	Enter Gradient Shimming Display menu.
Quit	Return to previous experiment.

Gradient Shimming Setup Menu

The Gradient Shimming Setup menu is used for calibrating gradient shimming parameters. To open this menu, click on Set Up in the Gradient Shimming System menu.

Button	Description
Go, dssh	Run and display profile spectra.
Gradient,Nucleus	Enter Gradient, Nucleus parameter setup menu.
Find gzwin	Run a calibration experiment to set parameters tof and gzwin, to optimize the spectral window for a given gzlvl. Pulse length and gain should be correctly set before using this button.
Calculate gzwin	Determine tof and gzwin from the profile spectrum using cursor positions.
Find tof	Run a calibration experiment to set tof. Pulse length and gain should be correctly set before using this button
Add Params	Add parameters for gradient shimming.
Return	Return to the Gradient Shimming System menu.

Gradient Nucleus Parameter Setup Menu

The Gradient Nucleus Parameter Setup menu is used for setting parameters according to type of gradient and nucleus. Recabling may be required to observe 1 H or 2 H nucleus. Check the value of gradtype when done. To open this menu, click on **Gradient,Nucleus** in the Gradient Shimming Setup menu.

Button	Description
Pfg H1	Set parameters for pulsed field gradient on ¹ H.
Pfg H2	Set parameters for pulsed field gradient on ² H.
Homospoil H1	Set parameters for homospoil gradient on ¹ H.
Homospoil H2	Set parameters for homospoil gradient on ² H.
Return	Return to the Gradient Shimming Setup menu.

Gradient Shimming Map Menu

The Gradient Shimming Map menu is used for making and retrieving shimmaps. To open this menu, click on **Shim Maps** in the Gradient Shimming System menu.

Button	Description
Automake Shimmap	Find gzwin and then run shimmap experiment.
Make Shimmap	Run shimmap experiment.
Shimmap Files	Enter Gradient Shimming Files menu.
Current Mapname	Show current mapname used for gradient shimming.
Save As	Save current parameter set and shimmap files. Enter mapname when prompted.
Return	Return to the Gradient Shimming System menu.

Gradient Shimming Files Menu

The Gradient Shimming Files menu is used for loading and copying shimmaps. To open this menu, click on **Shimmap Files** in the Gradient Shimming Map menu.

Button	Description
Cd to Userdir	Change directory to userdir shimmap.
or	
Cd to Systemdir	Change directory to systemdir shimmap.
Copy to Userdir	Copy selected shimmap into user shimmap directory.
or	
Copy to Systemdir	Copy selected shimmap into system shimmap directory.
Load Shimmap	Load files for selected shimmap.
Load Shimmap & Params	Load files and parameters for selected shimmap.
Rename	Rename selected shimmap under a new mapname.
Return	Return to the Gradient Shimming Map menu.

Gradient Shim Setting Menu

The Gradient Shim Setting menu sets shim values. To open this menu, click on **Set Shims** in the Gradient Shimming System menu.

Button	Description
Old Shims	Set shims from before first iteration of gradient shimming.
New Shims	Set shims from last iteration of gradient shimming.
Min Shims	Set shims from minimum rms err in last set of gradient shimming iterations.
Return	Return to the Gradient Shimming System menu.

Gradient Shimming Display Menu

The Gradient Shimming Display menu is used for displaying items associated with gradient shimming. To open this menu, click on **Display** in the Gradient Shimming System menu.

Button	Description
dgs	Display shim parameters.
List Shims	Display current mapname and last shim changes.
Display Shimmap	Display current shimmap used by gradient shimming.
Display Fit	Display shim fit from last iteration of gradient shimming.
Show Record	Display record of shim adjustments from last gradient shimming run.
Plot	Enter Gradient Shimming Plot menu.
Return	Return to the Gradient Shimming System menu.

Gradient Shimming Plot Menu

The Gradient Shimming Plot menu is used for printing and plotting shim parameters and files. To open this menu, click on **Plot** in the Gradient Shimming Display menu.

Button	Description
Print dgs	Print shim parameters.
Print Shims	Print current mapname and last shim changes.
Plot Shimmap	Plot current shimmap.
Plot Fit	Plot last iteration shim fit.
Plot b0	Plot b0 shim plot.
Return	Return to the Gradient Shimming Display menu.

Automated Shimming Menu

The Automated Shimming menu is used for running automated gradient shimming after a shimmap is made. To open this menu, click on **Main Menu** | **Setup** | **Shim**.

Button	Description
Lock Autoshim	Start lock autoshim using current method.
Lock Autoshim z1z2	Start lock autoshim on z1z2.
Gradient Autoshim on Z	Start gradient autoshim on z-axis shims
Return	Return to Acquisition Parameter Setup menu.

Gradient Autoshimming Files Menu

The Gradient Autoshimming Files menu is used for loading shimmaps. To open this menu, enter gmapshim('files').

Button	Description
Cd to Userdir	Change directory to userdir shimmaps.
or	
Cd to Systemdir	Change directory to systemdir shimmaps.
Current Mapname	Show current mapname used for gradient shimming.
Load Shimmap	Load files and parameters for selected shimmap.
Load Shimmap & Params	Load files and parameters for selected shimmap.
Rename	Rename selected shimmap under a new mapname.
Return	Return to Main menu.

Chapter 12. PFG Modules Experiments

Sections in this chapter:

- 12.1 "GCOSY—PFG Absolute-Value COSY," this page
- 12.2 "GHMQC—PFG HMQC," page 394
- 12.3 "GHMQCPS—PFG HMQC, Phase Sensitive," page 394
- 12.4 "GHSQC—PFG HSQC, Absolute Value or Phase Sensitive," page 395
- 12.5 "GMQCOSY—PFG Absolute-Value MQF COSY," page 396
- 12.6 "GNOESY—PFG NOESY," page 397
- 12.7 "GTNNOESY—PFG TNNOESY," page 398
- 12.8 "GTNROESY—PFG Absolute-Value ROESY," page 398
- 12.9 "PFG Selective Excitation," page 400

This chapter describes pulse sequences used with Varian PFG (Pulsed Field Gradient) modules. Each pulse sequence has a macro on the system that converts a 1D S2PUL parameter set into a parameter set ready to run the PFG experiment named.

The exact numbers used for gradient levels depends on the gradient amplifier and probe. For this reason, values listed are expressed in terms of gauss/cm.

12.1 GCOSY—PFG Absolute-Value COSY

The gcosy macro sets up parameters for the pulse sequence GCOSY (Pulsed Field Gradient Absolute-Value Correlated Spectroscopy). GCOSY is available on all systems with the PFG module.

Parameters

Table 56 lists the parameters used in GCOSY.

qlvl	1 (quantum selection level)
gzlvl1	Gradient amplitude (-32768 to +32767; use 3 gauss/cm)
grise	Gradient rise and fall time (in seconds, 0.00001)
gstab	Optional delay for stability (in seconds)
gt1	Gradient duration (in seconds, 0.002)
phase=1	Selects echo N-type coherence selection (default)
phase=2	Selects anti-echo P-type coherence selection

Table 56. GCOSY Parameters

Processing

```
Process N-type data with wft2d(1,0,0,1). Process P-type data with wft2d(1,0,0,-1). The optional 't2dc' argument to wft2d can be useful.
```

12.2 GHMQC—PFG HMQC

The ghmqc macro sets up parameters for the pulse sequence GHMQC (Pulsed Field Gradient Heteronuclear Multiple-Quantum Correlation). GHMQC is available on all systems with the PFG module.

Parameters

Table 57 lists the parameters used in GHMQC.

grise	Gradient rise and fall time (in seconds; 0.00001).
gstab	Optional delay for stability (in seconds).
gt1	Gradient duration (in seconds, 0.001).
gt2	Gradient duration (in seconds, 0.001).
gt3	Gradient duration (in seconds, 0.001).
gzlvl1*	Gradient amplitude (-32768 to +32767).
gzlvl2*	Gradient amplitude (-32768 to +32767).
gzlvl3*	Gradient amplitude (-32768 to +32767).
j	1 JXH (in Hz, 140 typical for 1 H- 13 C).
pwx	Decoupler pulsed pw90.
pwxlvl	Decoupler pulse power level.
* For more inform	nation, see the text in the section that contains this table.

Table 57. GHMQC Parameters.

gzlvl1, gzlvl2, and gzlvl3 and their times (gt1, gt2, and gt3) can eventually be fixed in their relationship (i.e., 2:2:-1, 0:4:-3, etc.).

For ¹³C, try gzlvl1=10 (gauss/cm), gt1=0.002, gzlvl2=10 (gauss/cm) gt2=0.00, gzlvl3=5 (gauss/cm), gt3=0.002, array gzlvl3 for maximum signal.

Processing

If gzlvl3 is the same sign as gzlvl1 and gzlvl2 (N-type data), process with wft2d('t2dc').

If gzlvl3 is opposite in sign to gzlvl1 and gzlvl2 (P-type data), process with wft2d('t2dc', 'ptype'). The sequence sets three gradients, all separately.

12.3 GHMQCPS—PFG HMQC, Phase Sensitive

The ghmqcps macro sets up parameters for GHMQCPS (Pulsed Field Gradient HMQC, Phase-Sensitive Version). Note that GHMQCPS does not do gradient HMBC. Features of this pulse sequence include:

- Two gt1's precede and follow d2, instead of follow and precede X 90's.
- TANGO pulse followed by a homospoil gradient precede the actual sequence.

- sspul option included.
- FAD on both first X 90 and first X 180.
- Two-step phase cycling on X 90's.

GHMQCPS is available on all systems, except *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*, with the PFG module.

Processing

Use wft2d(1,0,1,0,0,1,0,-1) with rp and lp the same as in s2pul spectrum. rp1 and lp1 will be close to zero. (wft2dnp can also be used with the "usual" 45° shift in f1 and f2).

Recommendations

For small molecule in H_2O and ${}^{13}C$:

gzlvl1=10000 gt1=0.001 gzlvl3=5025 gt3=0.001 gstab=0.0005 hsgpwr=10000 hsgt=0.005 nt=1 per inc. (FADed axials, though substantially reduced compared to original ghmqcps). nt=2 per inc. (no significant axials or FADed axials even with high concentration samples, recommended). nt=4 per inc. (more improvement in axial suppression).

12.4 GHSQC—PFG HSQC, Absolute Value or Phase Sensitive

The ghsqc<(nucleus)> macro sets up parameters for the pulse sequence GHSQC (Pulsed Field Gradient Heteronuclear Single-Quantum Correlation), absolute value or phase sensitive. The optional argument nucleus is ¹³C or ¹⁵N (e.g., ghsqc(15N)). GHSQC is available on all systems, except *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*, with the PFG module.

Parameters

Table 58 lists the parameters used in GHSQC.

Gradients need to maintain an approximately 4:1 (13 C) or 10:1 (15 N) ratio. The exact value can be determined from a ratio of the gyromagnetic ratios of the nuclei. The ratio can be controlled with either amplitude or time.

For ¹³C, try gzlvl1=20 (gauss/cm), gt1=0.002, gzlvl2=20 (gauss/cm), gt2=0.0005.

For ¹⁵N, try gzlvl1=30 (gauss/cm), gt1=0.0025, gzlvl2=15 (gauss/cm), gt2=0.0005.

le d'auto	$C_{\text{restant time for } \mathbf{Y}}$ such that $(t_{\text{rest}}, 2, \dots, 1)$
bigt	Constant time for X evolution (try 3 milliseconds).
dmf	Decoupler modulation frequency for 1st decoupler; controlled by dpwr.
dpwr	Power level for decoupling.
f1180	y/n flag for 1/2 dwell starting t_1 evolution delay.
gt1	First gradient duration.
gt2	Second gradient duration.
gzlvl1*	First gradient power level.
gzlvl2*	Second gradient duration.
j	One-bond heteronuclear coupling constant (140 for ${}^{13}C$, 90 for ${}^{15}N$).
nt	Number of transients; works with nt=1, but 2 or higher improves data.
phase=1	Use to generate an absolute-value data set.
phase=1,2	Use to select N,P-type selection (to be sorted later).
pwx	Pulse width for hard decoupler pulses.
pwxlvl	Decoupler power level for hard decoupler pulses.
satmode	Saturation mode; 'y' for transmitter presaturation.
satdly	Presaturation delay used if satmode='y'.
satfrq	Frequency desired for presaturation.
satpwr	Presaturation power.
* For more information, see the text in the section that contains this table.	

Table 58. GHSQC Parameters

Processing

Use wft2d('t2dc', 1, 0, 1, 0, 0, 1, 0, -1) for phase=1, 2 (phase up the first increment). The t2dc argument is optional. If the peaks appear to be reversed along the F_1 axis, multiply the last half of the coefficients by -1 to produce wft2d(1,0,1,0,0, -1,0,1). This sequence depends on the sign of the gradients used and on the order of phase (phase=1, 2 versus phase=2, 1).

Use wft2d(1,0,0,1) for phase=1.

12.5 GMQCOSY—PFG Absolute-Value MQF COSY

The gmqcosy macro sets up parameters for the pulse sequence GMQCOSY (Pulsed Field Gradient Multiple-Quantum Filtered Correlated Spectroscopy). GHMQOSY is available on all systems with the PFG module.

Parameters

Table 59 lists the parameters used in GMQCOSY.

Gradient levels for organic samples: try gzlvl1=10 (gauss/cm), gt1=0.003. Gradient levels for H2O samples: try gzlvl1=10 (gauss/cm), gt1=0.008

Processing

Process N-type data with wft2d(1,0,0,1). Process P-type data with wft2d(1,0,0,-1). The optional 't2dc' argument to wft2d can be useful.

qlvl	2 (quantum selection level for DQF COSY).	
grise	Gradient rise and fall time (in seconds, 0.00001).	
gstab	An Optional delay for stability (in seconds).	
gt	Gradient duration (in seconds, 0.001).	
gzlvl1*	Gradient amplitude (use 10 to 15 gauss/cm because it is multiplied by qlvl+1).	
phase=1	Selects echo N-type coherence selection (this is the default)	
phase=2	Selects anti-echo P-type coherence selection.	
* For more information, see the text in the section that contains this table.		

Table 59. GMQCOSY Parameters

12.6 GNOESY—PFG NOESY

The gnoesy macro sets up parameters for the pulse sequence GNOESY (Pulsed Field Gradient Nuclear Overhauser Effect Spectroscopy), absolute value or phase sensitive. GNOESY is available on all systems with the PFG module.

Parameters

Table 60 lists the parameters used in GNOESY.

Table 60. GNOESY Parameters

grise	Gradient rise and fall time (in seconds; 0.00001).
gstab	Optional delay for stability (in seconds).
gtl	Gradient duration (in seconds, 0.003).
gt2	Gradient duration (in seconds, 0.012).
gt3	Gradient duration (in seconds, 0.003).
gzlvl1*	Gradient amplitude (-32768 to +32767; use 30 gauss/cm).
gzlvl2*	Gradient amplitude (-32768 to +32767; use 30 gauss/cm).
gzlvl3*	Gradient amplitude (-32768 to +32767; use 30 gauss/cm).
phase=1	Selects echo N-type coherence selection (phase=1 is default);
phase=2	Selects antiecho P-type coherence selection;
phase=1,2	Selects phase-sensitive acquisition (N,P).
* For more information	on, see the text in the section that contains this table.

For organic samples, try gzlvl=gzlvl2=gzlvl3=10 (gauss/cm), gt1=gt3=0.003, gt2=0.012.

For H_2O samples, try gzlvl1=gzlvl2=gzlvl3=30 (gauss/cm), gt1=gt3=0.003, gt2=0.012.

Processing

Process N-type data with wft2d(1,0,0,1); process P-type data with wft2d(1,0,0,-1). The optional 't2dc' argument to wft2d can be useful. Process phase-sensitive data (phase=1,2) with wft2d(1,0,0,1,0,1,1,0).

12.7 GTNNOESY—PFG TNNOESY

The gtnnoesy macro sets up parameters for the pulse sequence GTNNOESY (Pulsed Field Gradient TN Nuclear Overhauser Effect Spectroscopy), absolute value or phase sensitive (includes optional presaturation). GTNNOESY is available on all systems, except *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*, with the PFG module.

Parameters

Table 61 lists the parameters used in GTNNOESY.

grise	Gradient rise and fall time (in seconds; 0.00001).
gstab	Optional delay for stability (in seconds).
gt1	Gradient duration (in seconds, 0.003).
gt2	Gradient duration (in seconds, 0.012).
gt3	Gradient duration (in seconds, 0.003).
gzlvl1	Gradient amplitude (-32768 to +32767; use 30 gauss/cm).
gzlvl2	Gradient amplitude (-32768 to +32767; use 30/gauss/cm).
gzlvl3	Gradient amplitude (-32768 to +32767; use 30 gauss/cm).
phase=1	Selects echo N-type coherence selection (phase=1 is the default);
phase=1,2	Selects phase sensitive acquisition (N,P).
phase=2	Selects antiecho P-type coherence selection;
satdly	Length of saturation during relaxation delay.
satmode	Saturation mode. Use analogously to dm (i.e., satmode='nnn', satmode='ynn', or satmode='yyn' (recommended)).
satpwr	Power level for solvent saturation

 Table 61. GTNNOESY Parameters

For organic samples, try gzlvl1=gzlvl2=gzlvl3=10 (gauss/cm), gt1=gt3=0.003, gt2=0.012.

For H_2O samples, try gzlvl1=gzlvl2=gzlvl3=30 (gauss/cm), gt1=gt3=0.003, gt2=0.012.

Processing

Process N-type data with wft2d(1,0,0,1). Process P-type data with wft2d(1,0,0,-1). The optional 't2dc' argument to wft2d can be useful.

Process phase-sensitive data (phase=1, 2) with wt2d(1, 0, 0, 1, 0, 1, 1, 0).

12.8 GTNROESY—PFG Absolute-Value ROESY

The gtnroesy macro sets up parameters for GTNROESY pulse sequence (including optional presaturation). GTNROESY is available on all systems, except *MERCURY-VX*, *MERCURY*, and *GEMINI 2000*, with the PFG module.

Parameters

Table 62 lists the parameters used in GTNROESY.

d2corr*Empirical correction, in μ s, of d2 (dependent on effective field of spir lock, i.e., tpwr and/or ratio).griseGradient rise and fall time (in seconds, 0.00001).gstabOptional delay for stability (in seconds).gt1Gradient duration (in seconds, 0.0015).gt2Gradient duration (in seconds, 0.0015).gzlvl1Gradient amplitude (-32768 to +32767; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to +32767; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at p11v1).p1lv1Power level for the p1 pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=1, 2For f1 quadrature by hypercomplex method (uses f1 axial peak displacement)		
griseGradient rise and fall time (in seconds, 0.00001).gstabOptional delay for stability (in seconds).gt1Gradient duration (in seconds, 0.0015).gt2Gradient duration (in seconds, 0.0015).gzlv11Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).pllvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f ₁ quadrature by hypercomplex method (uses f ₁ axial peak		
gstabOptional delay for stability (in seconds).gt1Gradient duration (in seconds, 0.0015).gt2Gradient duration (in seconds, 0.0015).gzlv11Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).pllvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f ₁ quadrature by hypercomplex method (uses f ₁ axial peak		
gt1Gradient duration (in seconds, 0.0015).gt2Gradient duration (in seconds, 0.0015).gzlv11Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to $+32767$; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).p1v1Power level for the p1 pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f ₁ quadrature by hypercomplex method (uses f ₁ axial peak		
gt2Gradient duration (in seconds, 0.0015).gzlvl1Gradient amplitude (-32768 to +32767; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to +32767; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at p1lv1).p1lv1Power level for the p1 pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f1 quadrature by hypercomplex method (uses f1 axial peak		
gzlvl1Gradient amplitude (-32768 to +32767; use 3 gauss/cm).gzlv12Gradient amplitude (-32768 to +32767; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).pllvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f ₁ quadrature by hypercomplex method (uses f ₁ axial peak		
gzlv12Gradient amplitude (-32768 to +32767; use 3 gauss/cm).mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).p1lv1Power level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f1 quadrature by hypercomplex method (uses f1 axial peak		
mixMixing time.ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).p1lvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f1 quadrature by hypercomplex method (uses f1 axial peak		
ntSet to 1 (min.), multiple of 2 (max.), multiple of 8 (recommended).p190° pulse on protons (power level at pllvl).p1lvlPower level for the p1 pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1, 2For f1 quadrature by hypercomplex method (uses f1 axial peak		
p190° pulse on protons (power level at pllvl).p1lvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1,2For f1 quadrature by hypercomplex method (uses f1 axial peak		
pllvlPower level for the pl pulse.phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1,2For f1 quadrature by hypercomplex method (uses f1 axial peak		
phase=1Selects P-type coherence selection (phase=1 is the default);phase=2Selects N-type coherence selection.phase=1,2For f1 quadrature by hypercomplex method (uses f1 axial peak		
phase=2Selects N-type coherence selection.phase=1,2For f_1 quadrature by hypercomplex method (uses f_1 axial peak		
phase=1, 2 For f_1 quadrature by hypercomplex method (uses f_1 axial peak		
phase=3 For f_1 quadrature by the TPPI method.		
pw Small (30°) pulse on protons (active only if ratio > 0). If pw=0, pw set to $p1/3$		
ratio tau/pw (uses CW transmitter lock if ratio is zero).		
rocomp='n' Sets no resonance offset compensation.		
rocomp='y' Sets resonance offset compensation (recommended).		
satdly Length of saturation during relaxation delay.		
satmode Saturation mode. Use analogously to dm (i.e., satmode='nnn', satmode='ynn', or satmode='yyn' (recommended)).		
satpwr Power level for solvent saturation		
<pre>sspul='y' Selects for trim(x)-trim(y) sequence at start of pulse sequence.</pre>		
tpwr Power level for the spin lock pulse(s).		
* For more information, see the text in the section that contains this table.		

 Table 62. GTNROESY Parameters

d2corr can be determined from the lpl and swl values of a properly phased spectrum by d2corr=(lpl*le6)/(360*swl). Note that d2corr seems to be dependent on swl. They are actually independent because changes in swl result in corresponding changes in lpl, which makes their ratio constant.

Do the following steps to find d2corr so that lpl=0, which gives better baselines in f_1 :

- 1. Run a gtnroesy experiment with d2corr set either at 0 or at a value found previously (nt and ni can be smaller, and the spectrum can be transformed early to do step 2).
- 2. Phase the resulting spectrum in f_1 . Determine lp1 and calculate d2corr from the relationship d2corr=(lp1*le6)/(360*sw1).
- 3. Add this value to the value of d2corr used in step 1.
- 4. Rerun the experiment—lp1 should be close to zero.
- 5. Note this value for any future experiment with the same value of tpwr and ratio.

Use of any method to make lpl=0 will result in a dc offset of f_1 slices. This should be removed by dc2d('fl') after the 2D transform. Enough noise should be left on the edges (in f_1) to permit this dc correction.

Processing

Process N-type data with wft2d(1,0,0,1); process P-type data with wft2d(1,0,0,-1). The optional 't2dc' argument to wft2d can be useful.

12.9 PFG Selective Excitation

The selexcit macro sets up parameters for PFG selective excitation with a presat option. Selective excitation is based on double PFG spin-echo, also known as excitation sculpting. Available on all systems, except *GEMINI 2000*, with the PFG module.

Parameters

Table 63 lists parameters used in PFG selective excitation.

gt1	Echo gradient during DPFGSE; typical value is 0.00025.
gt2	Echo gradient during DPFGSE; typical value is 0.0005.
gzlvll	Echo gradient during DPFGSE; typical value is 8000.
gzlvl2	Echo gradient during DPFGSE; typical value is 5000.
hsgpwr	Homospoil gradient during DPFGSE.
pwshape	Selective proton, shape; values are typically created by Pbox.
selpw	Selective proton, pulse width; values are typically created by Pbox.
selpwr	Selective proton, 180 power; values are typically created by Pbox.

Table 63. Selective Excitation Parameters

Presaturation option is included (using satflg, satdly, satfrq, and satpwr).

Reference

Shaka, et.al J. Am. Chem. Soc. 1995, 117, 4199.

Symbols

* notation in deconvolution, 206
+ mark by peak, 124
.Pbox_defaults file, 233
.Pbox_globals file, 233
[] (brackets) notation, 32

Numerics

10-mm sample tubes., 281 1D experiments, 71 1D pulse sequences, 71 1D water suppression, 85 2D Display Main Menu, 115, 121 2D Display Projection Menu, 115 2D DOSY experiment, processing, 53 2D DOSY sequences, 49 2D DOSY spectroscopy, 49 2D Exchange experiment, 153 2D FID matrices, 164 2D NMR, 87 2D Peak Picking menus Automatic, 121, 122, 125, 126 Display, 122, 124, 127 Edit, 121, 123, 125, 126 Edit Secondary, 126 File, 121, 124, 127 Main, 121, 122 3D DOSY experiments, 56 processing, 58 3D NMR, 127-132 4D data acquisition, 132 4-nucleus probes, 170 4-ply selective inversion pulses creation, 240 50-sample tray, 281 5-mm sample tubes, 281

Α

aa command, 298 Abort Acquisition button, 298 aborting acquisition in automation mode, 298 absolute-value color set, 113 Absolute-Value COSY pulse sequence, 136 absolute-value mode, 89, 107 absorption mode, 90 acetone in coolant bucket, 263 acosy command, 113 acosyold command, 113 Acqmeter display, 262 acqueue directory, 310 acqstatus parameter, 266, 308 acquiring data on standard test samples VAST, 354 acquisition block processing, 308 codes, 308 macros, 302 mode, 128 order of events, 302 status window, 91

acquisition computer state, 310 acquisitioninfo file, 310 add button, 215 add command, 212 add FID to add/subtract experiment, 212 Add/Sub button, 214 add/subtract buffer file, 211 add/subtract experiment, 211 Add/Subtract Menu, 211 addi command, 214 addition and subtraction of data, 211 additive weighting constant, 95 addpar macro, 88, 103, 104, 121, 132 adept command, 79 adiabatic decoupling, 240 adiabatic rf sweep waveform creation, 228 Adjust button, 123, 123, 126 air conditioning cycling, 266 Air valve, 361 alfa parameter, 97 amplifier lights on PFG amplifiers, 373 amplifier operation for PFG, 373 amplifiers used with indirect detection, 171 amplitude functions, writing expressions for, 235 amplitude modulated data, 102 analyze command, 44, 217 analyze.inp file, 43 analyze.list file, 217 analyze.out file, 217, 218 angled brackets (<>) notation, 27 apt macro, 72 APT pulse sequence, 72 aptaph macro, 73 Arm Z scale (mm), 361 aromatic solvents in coolant bucket, 263 array for a parameter, 30 array macro, 32 array of 1D spectra, VAST, 323 array of data, 30 array parameter, 32, 131 arraydim parameter, 34, 213 arrayed 1D experiment, 87 arrayed 2D and 3D, 90 arrayed 3D data set, 105 arrayed temperature operation, 267 arrayed VT experiments, 264 arraying order and precedence, 32 asize macro, 107 ASM-100 sample changer, 269 assign macro, 201, 204 assign nearest calculated transitions, 204 Attached Proton Test, 72 attn parameter, 237 au macro, 302, 303 auto assign button, 204 Auto button, 122 auto command, 297 auto_au macro, 298, 302 autodept macro, 79 autodir variable, 297 autogo command, 297, 297 VAST automation, 316

autoinfo file, 309 autolocking enhancements, 304 Automake Shimmap button, 381 automated shimming, 383 Automated Shimming menu, 389, 392 automatic 2D peak picking, 118, 125 analysis of COSY data set, 113 automatic gain, 33 automatic sample changer accessory, 269 automatic teller machine (ATM) cards caution, 25 automation files that control VAST operation, 363 VAST, 316 automation acquisition codes, 308 automation directory, 296, 308 automation file, 297 automation mode, 269, 298 automation run, 269 actions when finished, 299 adding samples during progress, 294 entering samples, 290, 297 experiment files, 296 interacting with acquisition, 298 multiple users, 290, 297 parameters, 299 starting, 297 starting a run, 296 autoname parameter, 299, 309 autoscale macro, 36, 218 autoscaling, 36 Autoshim on Z button, 378 autoshimming, enhancements, 304 av1 command, 107 awc parameter, 95 awc1 parameter, 95 awc2 parameter, 95 axial artifacts, 179 axis labeling, 131 axis parameter, 109

В

background operation, 129 Backup button, 127 Backup File button, 124 backwards linear prediction, 106 baseline correction, 53, 98, 107 baseline spikes, 262 bc command, 98, 107 bearing air supply, 265 bilinear rotation decoupling, 177 binary peak file, 119, 119, 124, 127 binom macro, 73 BINOM pulse sequence, 73 binomial water suppression, 73 biopolymer, 164 bipolar pulse pair stimulated echo experiment, DOSY. 50 bipolar pulse pair stimulated echo INEPT experiment, DOSY, 52 BIRD pulse nulling effect, 177 Bloch equation calculations, 224

Bloch Simulation subwindow, 226 Bloch simulator, 230, 237 boresize parameter, 375 Both button, 122 Box button, 115, 122, 123, 214 brackets ([]) notation, 32 broadband probes, 170 Bruker FID data, 103 bs parameter, 33 bscor parameter, 237 bsim parameter, 237 bw parameter, 238

С

c13p macro, 303 C13WET VAST, 319 calculate simulated spectrum, 203 theoretical spectrum, 199 calculated transitions file, 203 calculated transitions intensity, 203 calculating field maps, 381 calfa macro, 97 calibrate flow rate parameters, VAST, 351 probe volume, control data set, VAST, 346 probe volume, control verification using liquid handler, 347 probe volume, finding optimum value, 349 probe volume, VAST, 345 pulse width, 35 sample temperature curve, 264 sample volume, VAST, 350 volumes and flow rates, VAST, 345-354 XYZ coordinates of the rinse station, 354 calibrations for carbon indirect detection, 185 calibrations for nitrogen indirect detection, 197 Calibrations pane, gilson window, 355, 361 Cancel button, 115, 117 cancellation efficiency, 179 Carousel autosampler, 269 adjusting eject air, 274 adjusting sample float height, 274 error codes, 279 loading and unloading samples, 275 locked mode, 272 mounting the carousel, 273 removing the carousel, 273 run mode, 272 using the carousel, 270 Carr-Purcell Meiboom-Gill T2, 74 carryover evaluating, 355 cautions defined, 23 cdc command, 98 center command, 107 centersw macro, 112 centersw1 macro, 112 centersw2 macro, 127 cexp macro, 29 cfpmult macro, 98

change bar, 28 change samples automatically from the LCNMR/STARS pane, VAST, 316 manually, VAST, 314 using the enter program, VAST, 315 VAST. 314 chemical shift analysis commands, 221 chemical shifts, 202 chlorinated solvents in coolant bucket, 263 CJ compensator, 266 cla macro, 203 cla parameter, 202 clamp parameter, 203 class C amplifier systems, 171 cleanexp macro, 303 clear add/subtract buffer, 212 file of line assignments for iteration, 203 Clear button, 123 clfreq parameter, 202 clindex parameter, 203 clradd command, 212 coef file, 131 cold junction (CJ) compensator, 266 color intensity and CombiPlate, 333, 335 color intensity map, 110 color/grayscale adjustment bar, 113 combine all peaks within an area, 123 Combine button, 123 combine peaks, 120 CombiPlate analyzing data, 336 and integral region colors, 333 data analysis and VNMR, 332 data for analysis without VNMR, 334 fixing the color map, 336 preparing VNMR data for analysis with, 331 using for data analysis, 331-337 combiplate, command in UNIX window, 337 combishow, 332, 337 commands chemical shift analysis, 221 DOSY, 54 comment assigned to peak, 120 Comment button, 124, 126 common coefficients for wft2d processing, 102 common solvents, typical flow rates for, 353 completed transients, 91 complex FT processing, 103 complex interferograms, 99 composite pulse creation, 228 configure VNMR for the VAST accessory, 344 contour plot, 110 control data set for probe volume calibration, 346 control temperature, 265 conventions used in manual, 27 convertbru command, 103 coolant bucket, 263, 265 copying current binary peak file, 124 correct spectral drift, 99 corrected-difference spectroscopy, 210 Correlated Spectroscopy, 139 correlation analysis, 216

COSY acquisition procedure, 89 data automatic analysis, 113 data symmetrization, 112 pulse sequence, 139 COSY experiment, 137 phasing, 95 cosy macro, 89, 139 cosyps macro, 139 COSYPS pulse sequence, 139 coupling constants, 202 cpmgt2 macro, 74 CPMGT2 pulse sequence, 74 cptmp macro, 303 cpx macro, 230, 241 cr parameter, 112 cr1 parameter, 112 cr2 parameter, 127 create 2D parameters, 88 2D peak picking parameters, 121 3D coefficient file, 131 3D parameters, 104, 128 4D parameters, 132 experiment, 29 shaped rf pulses, 228 waveforms, 228 creating a shape file, 256, 256 creating a pseudo 2D data set from a VAST data set, 323 using vastglue, step by step example, 325 credit cards caution, 25 crl command, 112 crl1 command, 112 crl2 macro, 127 crof2 macro, 97 ct parameter, 33, 91 current experiment, 30 current FID number, 91 current transition assignments, 204 cursor move to center of spectrum, 112 position, 112 Cursor button, 115, 122, 123, 214 custom display order, VAST data, 327 customizing macro operation, 302 cycled difference NOE experiment, 75 cyclenoe macro, 75

D

d0 parameter, 239 d1 parameter, 239 d2 parameter, 88, 90, 91, 127, 128, 130, 239 d2pul macro, 76 D2PUL pulse sequence, 76 d3 parameter, 104, 127, 128, 130, 130 d4 parameter, 132 da command, 31 da display, 88 dac_p2 parameter, 40 dac_x parameter, 42 daslp parameter, 104

data processing creating a pseudo 2D display form a VAST data set, 32 FDM, 66 glued VAST data sets, 326-327 VAST data sets, 323–327 data processing type, 91 data systems for indirect detection, 172 data tables, 90 Dbppste experiment, 50 Dbppsteinept experiment, 52 dc command, 97, 98 dc correction, 98 dc2d command, 99, 107 dcg parameter, 98 dcon command, 110 dconi command, 113, 130 dconi parameter, 114 dconn command, 110 dcyc parameter, 238 Deconvolution menu, 208 deconvolution of spectra into individual lines, 205 decoupler as transmitter, 76 channel, 170 decoupling patterns creation, 228 decoupling power control, 185 defining integral regions, 41 delete all peak bounds, 120 all peaks, 123 experiment, 29 peak bounds, 123 peak nearest cursor, 120 peak nearest the cursor, 123 spectra from file, 36 delexp macro, 29 dels command, 36 delta2 parameter, 127 DEPT data processing, analyzing, plotting, 79 dept macro, 77 DEPT pulse sequence, 77 dept.out file, 79 deptgl macro, 79 DEPTGL pulse sequence, 79 deptproc macro, 80 deshapei macro, 242 detected spectral axis, 109 detection time, 88 df2d command, 111 dg, 88 dg2 macro, 105 dga command, 203 Dgcstecosy experiment (AV mode), 57 Dgcstehmqc experiment (AV mode), 57 DgcsteSL experiment, 51 dglp macro, 106 dgs macro, 262 diagonal peaks exclusion, 120 diagonally arrayed parameters, 32 diffusion coefficient, 43 measurements, 38 ordered spectroscopy, 46

projection, showing, 56 digitization of lineshapes, 205 dispersion tails, 93 display 3D plane, 132 3D processing parameters, 105 arrayed parameter values, 31 connectivities, 113 experiment library, 29 exponential curves, 36 file of line assignments for iteration, 203 full 2D display, 117 full screen, 107 grid lines, 110 image, 110 line assignments file, 204 next plane, 130 non-whitewashed stacked spectra, 111 peak bounds, 124 peak label, 125 peak marked with +, 124 peak number, 124 peaks on top of spectrum, 118 polynomial curves, 36 previous plane, 130 projection of data on plane, 132 second cursor, 116 shaped rf pulses, 223 simulated spectrum, 203 simulation parameter file, 203 spectra with whitewashing, 111 spectrum from array, 33 stacked 3D spectra, 111 subset of 3D planes, 130 theoretical spectrum, 207 trace at cursor position, 117 traces, 118 Display button, 122 display limits, 107 display mode parameter, 107 Display Shimmap button, 381 Distortionless Enhancement by Polarization Transfer, 77 dla macro, 203, 204 dlalong command, 204, 204 dll command, 204 dmg parameter, 107 dmg1 parameter, 107 dmg2 parameter, 127 doff parameter, 83 doneQ file, 296, 298, 298, 309, 310 Doneshot pulse sequence, Tcl-Tk acquistion panel, DOSY experiments, 46-?? dosy macro, 54 DOSY tools, 47 dot1 command, 36 double PFG spin-echo, 400 Double RELAY-COSY pulse sequence, 136 Double-Quantum Filtered COSY, 141 Double-Quantum Transfer Experiment, 150 downfield peak shift, 104 dpcon command, 110 dpconn command, 110

DPFGSE pulse program, 243 DPFGSE pulse program listing, 244, 245 dplane macro, 130, 132 dprofile macro, 231, 242 dproj macro, 132 dpwr parameter, 185 dqcosy macro, 141 DOCOSY pulse sequence, 141 dres parameter, 236 drift correction, 132 drift correction calculation, 99 dry nitrogen gas, 262, 265 ds command, 33, 95, 106, 118, 239 ds2d command, 111 ds2dn command, 111 dshape macro, 259 dshapef macro, 230, 242 dshapei macro, 259 dsp command, 203, 207 dsplanes macro, 130 dss command, 34, 111 dssa command, 34 dssan command, 34 dssh command, 34 dsshn command, 34 dssl macro, 34, 111 dssn command, 34 dsvast, 326, 337 dsvast2d, 326, 338 dsww command, 111 dual-broadband system, 170 recabling for indirect detection, 183 dutycycle variable, 40

Ε

E-BURP-1 shaped pulse, 239 echo experiment, 41 Edit button, 122, 125 eject air adjustment, 274 elements gradient stimulated echo, 47 energy level table, 203 enter command, 290, 295, 297, 297, 300, 301, 302 enter programming selections, 300 enter window customizing, VAST, 362 enterQ file, 296, 298, 298, 298, 302, 309 error codes, 279 error messages, 284 error processing, 303, 308 ethylene glycol for VT calibration, 266 evaluating carryover, 355 evolution time, 88 exchanger coil, 263 excitation bandwidth parameter, 256 excitation offset parameter, 256, 256 excitation sculpting, 400 exclude diagonal peaks, 120 exp1 command, 45 Expand button, 115, 117, 122, 215 expanded display, 117 experiment

library, 29 method, 295 parameter list, 302 queueing, 269 status on sample changer, 291 experiment 5 (exp5), 211 experiment, indicating first, 256 experimental frequency index of a transition, 203 experimental line listing limits, 204 experiments 2D DOSY, setting up, 49 3-D DOSY, 56 3D-DOSY, setting up, 56 Dbppste, 50 Dbppsteinept, 52 Dgcstecosy (AV mode), 57 Dgcstehmqc (AV mode), 57 DgcsteSL, 51 DOSY, 46 oneshot DOSY, 52 expfit command, 217 expl command, 36, 217 expl command, 45 expl.out file, 218 expladd macro, 45 explib command, 29 Exponential Analysis Menu, 37 exponential curves analysis, 36 exponential signal change, 38 expressions, writing, 235 extract 2D planes from 3D data set, 130, 132

F

F1 Axial Displacement, 179 f₁ axis, 88 F1 phase detection, 179 F1, F2, and F3 dimensions, 129 f19p macro, 303 f1coef parameter, 131 $f_2 axis, 88$ f2coef parameter, 131 FAD technique, 179 fbc macro, 53 FDM (Filter Diagonalization Method), 66-68 fid file extension, 278, 284, 289, 290, 309 fiddc3d parameter, 129 fiddle command, 208, 209 Fiddle program, 208 fiddle program, 53 fiddle2D command, 210 fiddle2d command, 210 fiddle2Dd command, 210 fiddle2dd command, 210 fiddled command, 210 fiddleu command, 210 FIDs 2D color intensity map, 111 move between experiments, 30 processing after each FID, 308 varying one or more parameters, 30 File button, 122 Filter Diagonalization Method (FDM), 66-68 filter part numbers, 184

filters for indirect detection, 171 filtration for indirect detection, 184 final simulated spectrum frequency limits, 203 Find, 292 Find Entry button, 292, 292 first domain phasing, 96 first evolution time, 129 first experiment, indicating, 256 first point distortion, 96 FIRST_FID flag statement, 256 fitspec command, 206, 206 fitspec.indata file, 205, 206 fitspec.inpar file, 205, 206 fitspec.outpar file, 205, 207, 207 fixed Gaussian fraction, 207 fla parameter, 239 flag indicating first FID, 256 flammable gases warning, 24 flow control meter, 262 flow rate parameter calibration, VAST, 351 flow rate setup faster runs, limited by the solution with the viscosity, 353 optimized for maximum through put of samples with constant viscosity, 353 slow runs, 352 fn parameter, 104 fn1 parameter, 104 fn2 parameter, 104 foldcc command, 112 folding about the diagonal, 112 folding along J=0 axis, 112 foldj command, 112 foldt command, 103, 112 forwards linear prediction, 106 Fourier transformation programs, 101 fp command, 36, 36 fpmult parameter, 96, 97 fpmult1 parameter, 96 fpmult2 parameter, 96 frequencies, correcting, 53 frequency functions, writing expressions for, 235 frequency referencing, 112 frequency shift, 104 frequency-dependent phase shifts, 131 ft command, 33, 95, 96, 98 ft1d command, 101 ft1da macro, 102 ft1dac and ft2dac commands, 164 ft2d command, 96, 98, 101, 102, 104 ft2da macro, 101, 102 ft3d command, 129, 131, 132 full 2D display, 117 Full button, 115, 117, 122, 215 full command, 107 full screen with room for traces, 107 full width at half-height (FWHH), 119 fullt command, 107

G

g_array parameter, 41

g_max parameter, 41 g_min parameter, 41 g_steps parameter, 41 g2pul pulse sequence, 46 g2pulramp pulse sequence, 46 ga command, 263 gain parameter, 33, 305 GARP modulation, 171, 181 Gaussian apodization constant, 95 Gaussian function shift, 95 Gaussian functions, 93 Gaussian lines, 205 Gaussian lineshapes, 206 gaussian macro, 94 gcal parameter, 373 gcosy macro, 393 GCOSY pulse sequence, 393 GCU (gradient compensation unit), 372 generating a shape file, 256 generating a shape file, 256 generating weighting values, 94 getDshape statement, 247 getGshape statement, 247 getplane command, 130, 132 getRshape statement, 244 getsn macro, 303 gf parameter, 93, 94, 95 gf1 parameter, 93, 94, 95 gf2 parameter, 95 gfs parameter, 95 gfs1 parameter, 95 gfs2 parameter, 95 ghmqc macro, 394 GHMQC pulse sequence, 394 ghmqcps macro, 394 GHMQCPS pulse sequence, 394 ghsqc macro, 395 GHSQC pulse sequence, 395 Gilson 819 valve Unit ID, 342 gilson command, 312, 312, 317 glass wool plug, 265 global parameters, 299 globalauto parameter, 294 glue order, 324 gmapshim macro, 379, 383 gmapsys macro, 378, 378, 382, 383, 383, 388 gmapz macro, 305 gmqcosy macro, 396 GMQCOSY pulse sequence, 396 gnoesy macro, 397 GNOESY pulse sequence, 397 go command, 263 grad_cw_coef parameter, 42 grad_p_coef parameter, 42 grad_p1 parameter, 40 grad_p2 parameter, 40 grad_x parameter, 42 Gradient Autoshim on Z button, 379 Gradient Autoshimming Files menu, 383, 389, 392 gradient calibration constants, 42 gradient compensated stimulated echo w/ spin lock experiment, DOSY, 51

Gradient Nucleus Parameter Setup menu, 389 gradient power levels, 42 gradient pulses, typical values for WET VAST, 321 gradient shimming, 379 pulse sequence, 380 Gradient Shimming menus Display, 389, 391 Files, 389, 391 Map, 389, 390 Plot, 389, 392 Setting, 389, 391 Setup, 388, 390, 390 System, 378, 388, 389 gradient stimulated echo element, 47 gradient tables, 374 gradtype parameter, 371 grayctr parameter, 117 grayscale center control, 117 grayscale color set, 113 grayscale images threshold adjustment, 117 grayscale slope control, 118 graysl parameter, 117 grid lines, 110 grid macro, 110 gshimlib directory, 383 gtnnoesy macro, 398 GTNNOESY pulse sequence, 398 gtnroesy macro, 398 GTNROESY pulse sequence, 398 gxmax parameter, 375 gymax parameter, 375 gzmax parameter, 375 gzsize parameter, 388 gzwin parameter, 381

Η

h1 macro, 303 h1p macro, 303 Hadamard transform, 240 halt command, 298, 303 hcchtocsy macro, 164 HCCHTOCSY pulse sequence, 164 Hd All button, 125 Hd Box button, 125 Hd Lbl button, 125 Hd Pk button, 124 HdNum button, 124 header parameter, 237 heat-exchanger coil in a water bath, 266 heating cycling, 266 helium contact with body, 24 helium gas flowmeters caution, 26 het2dj macro, 143 HET2DJ pulse sequence, 143 hetcor macro, 145 HETCOR pulse sequence, 145 hetcorps macro, 147 HETCORPS pulse sequence, 147 Heteronuclear 2D-J, 143 heteronuclear 2D-J data folding, 112 heteronuclear 2D-J experiment, 90, 101, 102, 103 Heteronuclear Chemical Shift Correlation, 145

heteronuclear chemical shift correlation, 92 Heteronuclear Chemical-Shift Correlation, Absolute Value and Phase Sensitive, 147 Heteronuclear Multiple-Bond Coherence, 172, 178 Heteronuclear Multiple-Quantum Coherence, 169, 172 Heteronuclear Multiple-Ouantum Coherence in Reverse Configuration, 174 heteronuclear Overbodenhousen experiment using **REVINEPT**, 197 heteronuclear shift correlation experiment, 98 heteronuclear spin-echo difference experiment, 175 high-power amplifiers cautions, 26 high-power decoupling and VT operation, 265 high-sensitivity protons, 169 high-temperature VT calibration, 266 HMBC pulse sequence, 172, 174, 178 HMOC experiments, 169, 394 parameters, 182 phasing, 107 pulse sequence, 172, 175, 182, 192 hmqc macro, 172, 179, 186 HMOCR parameters, 182 pulse sequence, 172, 174, 183 hmqcr macro, 172, 174, 179, 186 HMQC-TOCSY 3D sequence, 166 hmqctocsy macro, 166, 166 HMQCTOCSY pulse sequence, 166 HOHAHA pulse sequence, 162 HOM2DJ experiment, 149 hom2dj macro, 148 HOM2DJ pulse sequence, 148 homodec statement, 252 homonuclear 2D-J experiment, 92 homonuclear correlation phasing, 107 Homonuclear Hartmann-Hahn experiment, 162 Homonuclear J-resolved 2D, 148 homospoil gradient shimming, 385 homospoil gradient type, 375 horizontal axis, 109 Hproj (max) button, 115 Hproj (sum) button, 115 hregions macro, 303 hsqc macro, 197 HSQC-TOCSY 3D sequence, 167 hsqctoxySE macro, 167 hypercomplex method, 90, 90 hyper-hypercomplex mode, 130

I

```
ice bath to cool VT gas, 266
il parameter, 33
image macro, 110
imageprint macro, 110
implicitly arrayed parameter, 90
Improving Results, suggestions, 387
INADEQUATE data folding, 112
INADEQUATE pulse sequence, 150
inadqt macro, 150
Incredible Natural Abundance Double-Quantum
Transfer Experiment, 150
```

index of experimental frequency of a transition, 203 index2 parameter, 130 indicating first experiment, 256 indirect detection experiments, 169 indirect detection probe, 170 Indirect•nmrTM probe, 170 indirectly detected axis, 109 inept macro, 80 INEPT pulse sequence, 77, 80 Info button, 124, 126 info_n file, 42 initialize VAST, 316 initialize_iterate macro, 203 injector valve, 342 ins parameter, 119 Insensitive Nuclei Enhanced by Polarization Transfer, 80 insert peak at cursor position, 120, 123 inssref parameter, 119 instrumental errors, 208 integral resets, 98 intensity of calculated transitions, 203 intensity threshold for simulated spectrum, 203 interactive 2D peak picking, 118 add or delete peak bounds, 126 display 2D traces and projections, 113 parameter arraying, 32 spectra add/subtract, 214 weighting parameters, 94 Interactive 2D Color Map Display Main Menu, 115 Interactive 2D Display Projection Menu, 115 Interactive Mode button, 214 interferogram left-shift, 104 interferogram phase rotation, 104 interferograms, 87 interleaved acquisitions, 33 intmod parameter, 34, 98 inversion pulse creation, 240 inversion-recovery T_1 experiment, 35 is parameter, 119 iterate parameter, 203 iterative mode of spin simulation, 203 iterative spin simulation, 199

J

jexp command, 303 jumpret macro, 81

Κ

ketone solvents in coolant bucket, 263 killft3d macro, 130 kind macro, 38 kinds macro, 38 kinetics data analysis, 216 kinetics studies, 37 kini macro, 38 kinis macro, 38

L

L.O. frequency, 170 L.O. SELECT switch, 171 Label button, 123, 126 label for a peak, 120 LAME program, 199 LAOCOON program, 199 last lock solvent, 310 lastlk file, 310 lb parameter, 93, 94, 95 lb1 parameter, 93, 94, 95 lb2 parameter, 95 LCNMR/STARS tab, VAST, 316 least squares curve-fitting method, 216 left command, 107 left shift interferogram, 104 lfs (low-frequency suppression), 104 line broadening factor, 95 line number table, 202 line positions, 204 line shapes, correcting, 53 linear amplifier systems, 171 linear prediction, 105 processing control, 103 linewidth for deconvolution calculations, 206 linewidth of simulated spectrum, 203 Liquid Handler, 355 air connections, diagram, 343 connecting signal and power cables, 342 signal cable connections diagram, 344 Liquid Handler arm XYZ coordinates, adjusting, 361 Liquid Handler window, 312, 317 typical default values, 313 liquid nitrogen, 265, 267 list button, 204 list of line frequencies, 204 literature references, 92 ll2d command, 116, 118 ll2d directory, 118, 127 112d file, 124 ll2dmode parameter, 120, 123, 127 llfrq parameter, 204 loc parameter, 269 loc, as used with VAST, 314 local oscillator frequency, 170 local oscillator signal, 181 lock interlock, 267 lock pin, 272, 272 locked mode, 272 lockpower parameter, 305 locQ file, 298 log file, 310 long-range HMQC pulse sequence, 172 Lorentzian lines, 205 Lorentzian lineshapes, 206 low-temperature VT calibration, 266 lp parameter, 96, 97 lp1 parameter, 96, 106, 107 lp2 parameter, 96, 127 Isfid parameter, 212 lsfid1 parameter, 104 lsfid2 parameter, 127 lsfrq parameter, 104, 130

lsfrq1 parameter, 104, 130 lsfrq2 parameter, 104, 130

М

macromolecules, 178 macros acquisition, 302 customizing macro operation, 302 magnet quench warning, 24 magnetic equivalence factoring, 199 magnetic media caution, 25 magnitude-mode 2D experiment, 102 magnitude-mode transform, 102 Main Control, gilson window, 355, 360 make3dcoef macro, 131 making a shape file, 256 making a shape file, 256 manual sample change VAST, 314 mapping the shims, 377 Marion and Wuthrich, 90 Mark button, 123, 126, 207 mark button, 204 mark command, 116, 205 mark1d.out file, 204, 205, 207 match experimental and calculated lines, 203 maximum intensity, 109 maxincr parameter, 236 measured line frequencies, 203 metal objects warning, 23 methanol for VT calibration, 266 Microflow probe connecting air supply from Valco valve, 342 microflow probe and transfer tube volumes, 345 microflow probe cell, 348 min button, 215 mix parameter, 156 mixing patterns creation, 228 modulated decoupling, 181 move last FID and parameters, 30 parameters between experiments, 30 mp command, 30mqcosy macro, 151 MQCOSY pulse sequence, 151 multiexperiment mode, 298 multi-FID add/subtract experiment, 212 multiple arrays of parameters, 31 multiple frequency suppression options VAST, 319 Multiple-Quantum Filtered COSY, 151 multiply first point, 96 multiply selective adiabatic inversion pulse creation, multiply selective pulse creation, 239

Ν

name parameter, 236 names, reserved shape, 235 NeedleRinseRate, 359 NeedleRinseVolume, 359 negative intensities, setting, 112 nested supercycle, 241 nextpl macro, 130 ni dimension, 93 ni parameter, 87, 88, 94, 127, 128, 130 ni2 parameter, 94, 104, 127, 128, 130 ni3 parameter, 132 niter parameter, 203 nitrogen contact with body, 24 nitrogen gas caution, 262 nitrogen gas flowmeters caution, 26 nitrogen gas for VT use, 262 nlivast, 338 nlivast2, 338 nlivast3, 339 nll command, 204 nm2d command, 110 NOE difference experiment, 82, 82 NOEDIF experiment, 82 noedif macro, 82 NOESY data folding, 112 data symmetrization, 112 phasing, 107 pulse sequence, 153 noesy macro, 153 noise reject parameter VAST, 321 noise suppression, 117 non-whitewashed stacked spectra, 111 notational conventions, 27 npoint parameter, 36 nt parameter, 33 nt_array parameter, 41 nt first parameter, 41 N-type data, 400 N-type experiment, 102 N-type peak selection, 101 N-type peaks, 102 Nuclear Overhauser Effect Spectroscopy, 153 Nuclear Overhauser Enhancement difference experiment, 82 null string, 308 number of increments, 87 number of iterations, 203 NumRinses, 357

0

observe channel, 170 observe pulse and power VAST and solvent suppression, 319 off-resonance SEDUCE-1 decoupling, 240 ofs parameters, 238 oneshot DOSY experiment, 52 operating the amplifier for PFG, 373 optical sensor, 274 opx macro, 230, 241 order of 2D experiments, 91 origin for phasing, 132 output transistor failure on VT controller, 266 over-temperature circuit, 267 oxidation of heater and thermocouple, 265

Ρ

p31p macro, 303 pa command, 108 pa1 command, 107, 108 pacemaker warning, 23 pacosy command, 113 pad parameter, 37, 263 pad parameter, 300 padept command, 80 paired data, 216 Pandora's Box, See Pbox par2d macro, 88 par3d macro, 104 par4d macro, 132 parameters array format, 30 array order and precedence, 32 Dbppste, 51 Dbppsteinept, 53 Dgcstecosy, 57 Dgcstehmqc, 58 DgcsteSL, 51 FDM, 67 interactive arraying, 32 iterated, 203 jointly arrayed, 32 move between experiments, 30 move saved display, 30 multiple arrays, 31 oneshot DOSY, 52 returning shape, 256, 256 parim macro, 118 parlib directory, 71 parll2d macro, 121 parlp macro, 103, 106 path3d parameter, 129 Pbox, 228 macros reference, 241 menus, 229 miscellaneous PSG statements, 252 parameter list, 232, 236 PSG statements, 243, 243 pulse calibration numbers, 229 wave definition strings, 232 Pbox command (UNIX), 231, 258 Pbox equation evaluator, 235 Pbox parameter, writing a, 255 Pbox.inp file, 230, 231, 231 pbox_bw macro, 242 pbox_bws macro, 242 pbox_dec2off statement, 250 pbox_dec2on statement, 250 pbox_dec2pulse statement, 248 pbox_dec3off statement, 250 pbox_dec3on statement, 250 pbox dec3pulse statement, 248 pbox_decoff statement, 250 pbox_decon statement, 249 pbox_decpulse statement, 248 pbox_dmf macro, 230, 242 pbox_dres macro, 242 pbox_drex macro, 230 pbox_files macro, 243 pbox_grad function, 243

pbox_grad statement, 252 pbox_name, 242 Pbox psg.h include file, 243 pbox_pulse function, 243 pbox_pulse statement, 247 pbox_pw macro, 230, 242 pbox_pwr macro, 230, 242 pbox pwrf macro, 242 pbox_rst macro, 230, 243 pbox_sim3pulse statement, 249 pbox_sim4pulse statement, 249 pbox_simpulse statement, 248 pbox_xgrad statement, 252 pbox_xmtroff statement, 249 pbox_xmtron statement, 249 pbox_ygrad statement, 252 pbox_zgrad statement, 252 pboxget macro, 242 pboxpar macro, 230, 242 pcon command, 110 pcss.outpar storage file, 221 Peak button, 115, 121, 122, 125 peak comment, 124 peak file information, 124 peak files, 119 peak height, 36 peak intensities, 36 peak label, 123, 126 peak volume, 122 peak(s) selecting solvent for suppression, VAST, 316 peak2d command, 109 peaks.bin file, 118, 124, 124 peaks_f#f#_#.bin file, 119 peek tubing, I.D. and volume calculations, 341 perform spin simulation, 203 Performa I PFG module, 371 amplifier lights, 373 parameter settings, 372 shimming, 374 Performa II PFG module, 38, 371 amplifier lights, 37 parameter settings, 372 reset button, 373 shimming, 374 Performa III PFG module, 371 amplifier lights, 373 Performa XYZ PFG module, 371 amplifier lights, 373 amplifier-produced quiescent current, 372 parameter settings, 372 shimming, 374 pexpl command, 36, 45, 218 pexpladd macro, 45 PFG Absolute Value COSY, 393 PFG Absolute Value MQF COSY experiment, 396 PFG Absolute-Value ROESY experiment, 398 PFG HMQC experiment, 394 PFG HMQC, Phase Sensitive, 394 PFG HSQC experiment, 395 PFG NOESY experiment, 397 PFG Selective Excitement, 400 PFG systems experiments, 393

gradient subsystem, 371 user interface, 372 PFG TNNOESY experiment, 398 pfg_pulse statement, 253 pfgon parameter, 372, 372, 374, 374 pge macro, 41 pge parameter, 46 pge parameter set, 40 PGE pulse sequence, 39, 40 pge_calib macro, 42 pge_data macro, 42 pge_output macro, 45 pge_process macro, 42 pge_process macro, 42 pge_results macro, 43, 45 pge_setup macro, 41, 42 pgeramp pulse sequence, 46 ph command, 95 ph parameter, 239 ph1 command, 95, 107 phase angle mode, 107 Phase button, 106 phase correction, 106 phase cycling, 103 phase modulation functions, writing expressions for, phase parameter, 90, 130, 130, 131, 179 phase parameters, 96 phase rotate interferogram, 104 phase shift the receiver, 103 phase2 parameter, 127, 130, 131 phase3 parameter, 132 phased color set, 113 phased data, 89 phases, correcting, 53 phase-sensitive 2D NMR, 89 phase-sensitive COSY pulse sequence, 139 phase-sensitive HETCOR, 147 phase-sensitive mode, 107 phase-twist, 102 phasing corrections, 107 phfid parameter, 212 phfid1 parameter, 104 phfid2 parameter, 127 phi parameter, 227 php parameter, 239 pi3ssbsq macro, 94 pi4ssbsq macro, 94 pl command, 34, 111 pl2d command, 111 plane parameter, 132 plate_glue, 327, 339 pldept macro, 80 plfit macro, 207 plgrid macro, 110 pll2d command, 118, 127 plot 2D traces and projections, 113 connectivities, 113 deconvolution analysis, 208 exponential curves, 36 grid lines, 110 image, 110 limits, 107

non-whitewashed stacked spectra, 111 polynomial curves, 36 regression analysis, 218 results of 2D peak picking, 118 series of 3D planes, 130 spectra, 34 spectra with whitewashing, 111 Plot button, 115, 115, 208 plplanes macro, 130 Plunger, 361 plvast, 326, 339 plvast2d, 326, 339 plvastget, 326, 326 pmode parameter, 95, 107, 131, 132 point-by-point intensity of the spectrum, 205 polarization transfer experiment, 79 poly0 macro, 217 polynomial curves analysis, 36 polynomial fitting, 98 power mode, 107 power switch (VT unit), 262 pph macro, 242 pprofile macro, 242 pre_sat statement, 253 preacquisition delay, 37, 263 predicting variables, 216 preparing VAST for use, 312 presat macro, 85 presat statement, 253 presaturation mechanism, 178 prevpl macro, 130 Prime pump, 361 probe directory, 71, 135 fast rate, 351 guide holes, 274 impedance mismatch, 46 pulsed Z-gradient, 38 slow rate, 351 slow vol, 351 temperature, 261 tuning in reverse mode, 184 probe volume, 314 calibration, VAST, 345 control data set, obtaining using a syringe, 346 finding optimum value, 349 initial estimates, 348 VAST, 345 verifying using liquid handler, 347 ProbeFastRate, 358 ProbeSlowRate, 358 ProbeSlowVol, 358 ProbeVolume, 357 proc parameter, 103 proc1 parameter, 91, 103, 106 proc2 parameter, 105, 106 processing 2D DOSY experiments, 53 3D-DOSY experiments, 58 processing stored 2D data, 88 processing type, 103 processing, displaying, and plotting from glued data sets, 326-327 macros for VAST data, 337

VAST data sets, 323–325 procplot macro, 302 profile-type experiments, 376 Proj button, 115, 115, 117 proj command, 111 project data onto axis, 111 Projection Menu, 117 prosthetic parts warning, 23 protective relay on VT controller, 267 proton 2D-J spectra, projecting, 112 pseudo 2D data set creating from a VAST data set, 323 pseudo macros, 94 pseudo-2D, 340 pseudo-echo weighting, 93 psglib directory, 71, 136 psgQ file, 296, 298, 309 pshape macro, 242 ptspec3d parameter, 129 P-type data, 400 P-type experiment, 102 ptype parameter, 238 P-type peak selection, 101 pulse creation routine, 228 pulse length parameter, 256, 256 pulse sequence variations, 46 pulse sequences 3D DOSY, 56 Dgcstecosy, 56 Dgcstehmqc, 56 DOSY, 46, 49 pulse shaping, 255, 256, 257 pulse template file, 223 pulse width calibration, 35 Pulsed Field Gradient experiments, 393 pulsed field gradient shape creation, 228 pulsed gradient experiments, 38 pulsed Z-gradient probe, 38 pulsetool command, 223 Pulsetool program, 223 pure 2D absorptive lineshapes, 95 pure absorption spectra, 102 Push Volume, 350 PushVolume, 357 putwave macro, 242 pw parameter, 238 pwr1 command, 107 pwxcal macro, 188 Pxfid command (UNIX), 259 pxshape macro, 230, 242 Pxsim command (UNIX), 259 Pxspy command (UNIX), 259

Q

queueing experiments, 269

R

r macro, 106 Rack Def., gilson window, 355, 360 Rack Definition, 313 radio-frequency emission regulations, 26 ramp length, 46 rate of VT gas flow, 265 react macro, 303 Read button, 124 Read Text button, 124 reading binary peak file, 124 text peak file, 124 real FT processing, 103 real-time 2D, 92 Redraw button, 115, 117 ref_pw90, 318 VAST, 318 ref_pw90 parameter, 238 ref_pwr, 318 VAST, 318 ref_pwr and ref_pw90 calibration data, 229 ref pwr parameter, 238 reference 90 degree pulse width, 256, 256 reference deconvolution, 53, 208 reference literature, DOSY-related, 64 referencing, clearing in 2D spectra, 112 refocusing pulses creation, 240 refofs parameter, 237 refpos1 parameter, 112 refrigerating device, 266 region-selective 3D processing, 129 regionx_results file, 45 regression, 216 regression analysis, 216 regression.inp file, 216, 217, 218, 220, 220 relaxation times measurement, 36 relaxation-sensitive simulation, 240 RELAY-COSY pulse sequence, 136 relayh macro, 136 relief valves warning, 25 removable quench tubes warning, 25 removing systematic errors, 53 reps parameter, 237 reserved shape names, 235 Reset arm, 361 Reset button, 123 reset parameters, 132 resetf3 macro, 132 resolution enhancement, 93, 95 restore spin system to before last iterative run, 204 resume command, 303 retaining disk, 272 Return Home, 361 Return key, 28 returning a shape structure, 256, 257 shape parameters, 256, 256 reverse configuration, 170, 172 review papers, DOSY, 65 **REVINEPT**, 197 rf generation for transmitter and receiver, 170 rf signal routing for indirect detection, 171 rf system of spectrometer, 170 rfl parameter, 112 rfl2 parameter, 127 rfp parameter, 112 rfp2 parameter, 127 Rheodyne Injector Valve bypassing, 342

connecting transfer tubing, 341 right command, 107 rinput macro, 217 RinseDeltaVol, 358 RinseExtraVol, 359 roesy macro, 156 ROESY phasing, 107 ROESY pulse sequence, 156 room temperature stability, 266 rotate command, 112 rotate homonuclear 2D-J data, 112 rotated homonuclear 2D-J data folding, 112 Rotating Frame Overhauser Experiment, 156 rp parameter, 96 rp1 parameter, 96, 106, 106 rp2 parameter, 96, 127 RS-232C serial link, 281 rt command, 34 rubber stopper, 274 run mode, 272

S

s macro, 30 s2pul macro, 85, 181 S2PUL pulse sequence, 85, 172, 180 s2pulr macro, 86, 172, 192 S2PULR pulse sequence, 86, 172, 183 safety circuits on VT controller, 266 safety precautions, 23, 25 safety sensor for VT controller, 266 sample adjusting liquid height, 282 change cycle time, 281 location, 269 positioning in spinner, 282 positions, 281 protocols, 363 removal errors, 285 spinning, 265 temperature, 264 sample changers, 269 sample preparation, 282 Sample Def, gilson window, 356 SAMPLE Def. field definitions, 356 gilson window, 355 sample entry form window VAST, 315 sample tray, 281 sample location, 297 sample locations, 295 Sample Volume, 348 sample volume, calibrate, 350 sample zero, 281 SampleDepth, 359 SampleExtraVol, 359 SampleHeight, 359 sampleinfo file, 298, 298, 309 SampleKeepFlag, 357 SampleVolume, 357 SampleWellRate, 359 satellite signals, 210

save button, 215 saved display parameters, 30 saveglobal parameter, 299, 302 sb parameter, 94, 95 sb1 parameter, 94, 95 sb2 parameter, 95 sbs parameter, 94, 95 sbs1 parameter, 94, 95 sbs2 parameter, 95 sc parameter, 34, 107 sc2 parameter, 34, 107 scalelimits macro, 36, 218 Scout Scan VAST, 321 sdp macro, 56 sealed samples at elevated temperatures, 264 sealed samples caution, 267 second cursor pair, 116 second evolution time, 129 SEDUCE-1 decoupling, 240 Select button, 215 selecting composite observe pulse VAST, 321 selecting solvents for suppression VAST, 318 selective excitation, 228, 400 selective Fourier transformation, 102 selex macro, 230, 241 selexcit macro, 400 sensitivity enhancement, 95 Set Params button, 203, 203 set3dproc command, 131 setgauss macro, 207 setgcoil macro, 375 sethw command, 262, 297 setlimit statement, 253 setlk macro, 302, 303, 303, 304 setLP1 command, 105 setref macro, 112 setref1 macro, 112 setref2 macro, 112 setsw1 macro, 112 setting up 2D-DOSY experiments, 49 3D DOSY experiments, 56 setvalue command, 213 setwave macro, 230, 241 sfrq parameter, 237 Sh All button, 125 Sh Box button, 124 Sh Lbl button, 125 sh parameter, 238 Sh Pk button, 124 shape file name parameter, 256, 256 shape file, generating, 256 shape file, generating a, 256 shape information, 241 shape macro, 242 shape names, reserved, 235 shape parameters, returning, 256 shape statement, 243 shape structure, returning a, 256, 257 shaped gradients creation, 241 shaped pulses, 223

shaped rf pulses creation, 228 shapedgradient macro, 243, 243 shapelib directory, 231 shaping pulses, 256, 257 shift frequency of spectrum, 130 shim field maps, 380 shimmaps displaying, 381 distributing to users, 383 loading, 382 shimming PFG systems, 374 ShNum button, 124 showfit macro, 207 showing diffusion projection, 56 sidechain assignments in fully ¹³C-enriched biopolymer, 164 sign of frequencies, 103 sign reversals along f_1 , 103 signal-to-noise estimate, 303 simple excitation pulse, 239 simulated spectrum frequency limits, 203 sine macro, 94 sine window function, 94 sinebell macro, 94 shift, 95 time period, 95 weighting, 93 sinebell-squared weighting, 93 sinesq macro, 94 sine-squared window function, 94 single-broadband systems, 170 recabling for indirect detection, 183 single-user sample changer operation, 289 skyline projection, 111 slfreq parameter, 203, 204 SLP, 165 slpwr parameter, 156 slw parameter, 203, 206, 207, 207 smaxf parameter, 203, 204 sminf parameter, 203, 204 SMS sample changer, 269 sn parameter, 303 solids high-power amplifiers caution, 26 solvent mixture equilibration VAST, 322 solvent subtraction 2D transform, 104 solvent suppression retrieving parameters for, VAST, 317 VAST, 317 VAST acquisition and processing options, 320 VAST post acquisition signal suppression, 320 wet width, VAST, 321 solvent suppression background information VAST, 322 solvent suppression problems VAST, 322 sp parameter, 109, 117 sp wp button, 215 sp1 parameter, 109, 117 sp2 parameter, 127 spadd command, 213 specdc3d parameter, 132 spectral

display, 33 drift correction, 132 plot, 34 spectral width, setting the, 112 spectroscopy, 2-D DOSY, 49 spectrum deconvolution, 206 move cursor to center, 112 spin lock, 156 spin simulation, 199 spin simulation algorithms, 199 Spin Simulation Line Assignment Menu, 204 spin-echo difference experiment, 175 spini.indata file, 204 spini.inpar file, 204 spini.la file, 204, 204 spini.outpar file, 204 spini.savela file, 204 spinll macro, 204, 204 spinner air supply, 265 spins command, 203, 203, 204 spins.inpar file, 205 spins.list file, 203, 204 spins.outdata file, 205 spins.stat file, 205 spline fitting, 98 spsm command, 202 spsub command, 213 sqcosine macro, 94 sqsinebell macro, 94 ssfilter and solvent suppression VAST, 323 ssfilter parameter, 104 ssorder parameter, 104 st parameter, 238 Standard Two-Pulse in Reverse Configuration, 86 Standard Two-Pulse Sequence, 85 Standard Two-pulse with Decoupler as Transmitter, 76 start of chart, 107 starting point for a deconvolution, 207 States, Haberkorn, and Ruben, 90, 102 status command, 295, 298 stdpar file, 302 stdpar parameters, 299 steps parameter, 236 stepsize parameter, 237 sth parameter, 203, 203, 204 stim parameter, 40 stimulated echo experiment, 40 store line assignments, 204 su command, 263, 265 sub button, 215 sub command, 212 subexperiments, 30 subtract FID from add/subtract experiment, 212 sucyc parameter, 237 summing projection, 111 Sun-based VXR-5000 data system, 172 supercycling, 241 superhypercomplex data acquisition, 127 supply air pressure variations, 274 suppression of 13C satellites, VAST, 319

svf command, 34 svs parameter, 203 sw parameter, 238 sw1 parameter, 87, 127, 128 sw2 parameter, 104, 127, 128, 130 sw3 parameter, 132 switchable probes, 170 switchable relays, 171 symbol, 32 symmetrize data, 112 Syringe Volume, 361 sysgcoil parameter, 375 systematic errors, removing, 53

T

 T_1 and T_2 NMR data analysis, 216 t1 command, 36 t_1 domain, 93 T_1 experiment, 36 T_1 exponential time, 35 T1 parameter, 237 t1s command, 36 t2 command, 36 T_2 experiment, 36 T_2 measurement, 74 T2 parameter, 237 t₂ time, 88 $t\bar{2}s$ command, 36 TANGO pulse, 394 Tcl-Tk acquisition panel, Doneshot pulse sequence, Tcl-Tk process panel, 2D_DOSY pulse sequences, 55 Tcl-Tk process2 panel, 3D_DOSY pulse sequences, 60 temp command, 264 temp parameter, 30, 261, 263, 263 tempcal command, 266 temperature calibration curve, 266 temperature control, 263 temperature control window, 265 temperature interlock parameter, 267 temperature limit sensor, 267 temperature readout, 263 temperature regulation, 267 testct macro, 303 testsn macro, 303 text peak file, 119, 124 th parameter, 110, 110, 202 th2d parameter, 120, 122, 126 theoretical spectrum for spin 1/2 nuclei, 199 thermocouple, 261, 266 theta parameter, 227 third indirectly detected dimension, 133 threshold for integrating peaks, 120 time periods in the 2D experiment, 88 Time Proportional Phase Incrementation, 90 time-domain frequency shifting, 130 tin parameter, 267, 299 TMS satellite signals, 210 tncosyps macro, 158 tndqcosy macro, 158 tnmqcosy macro, 159

tnnoesy macro, 159 throesy macro, 160 tntocsy macro, 161 tocsy macro, 162 TOCSY phasing, 107 TOCSY pulse sequence, 162 tools. DOSY. 47 Total Correlation Spectroscopy, 162 TPPI experiment, 90 TPPI method, 90, 91 TPPI phase cycling, 131 Trace button, 115, 117 trace parameter, 94, 107, 109, 110, 111, 132 tramp parameter, 46 transfer tube calculating the volume, VAST, 341 transform non-arrayed 2D data, 104 transition amplitude of calculated transitions, 203 Transverse ROESY, 164 trev parameter, 239 Trial WET VAST, 317, 321, 322 triangularize 2D spectra, 112 Triple•nmr[™] probe, 170 triple-resonance probe, 170 trise parameter, 375 troesy macro, 164 TROESY pulse sequence, 164 true contour plot, 113 true contour plot display, 110 truncation wiggles, 93 two-component analysis, 45 type parameter, 236 typical flow rates for common solvents, 353

U

undospins macro, 204, 204 UNIT ID of Gilson 819 valve, 342 Unmark button, 123, 126 upfield peak shift, 104 Use Mark button, 205 usemark macro, 205, 206

V

Valco valve connecting air tubing, 342 variable temperature, See VT VAST air connections, 343 autosampler, 269, 355 data analysis, 337 files that control operation, 363 Hardware installation and configuration, 341-362 interface to VNMR, 355 microtiter plate, 337 multiple frequency suppression options, 319 Operation and data processing, 311-341 sample entry form window, 315 sequence pane, 318 set up NMR experiments, 314 shut down, 317

signal cable connections, 344 solvent mixture equilibration, 322 solvent suppression problems, 322 solvent suppression setup, 317 solvent suppression, retrieving parameters for, 317 transfer tube volume. 341 WET shapes, 323 writing protocols, 363 VAST data sets, examples of plots 96 1D spectra plotted using plvast, 328 using a custom display order, 328 vastget, 326, 340 vastglue, 324, 340 vastglue2, 325, 340 versatile automatic sample transport, 311 vertical scale adjustment, 117 vertical temperature gradient, 265 vloc as used with VAST, 314 Volume button, 122, 126 vortex plug, 265 Vproj (max) button, 115 Vproj (sum), 115 vs parameter, 110 vs2d parameter, 110, 119 VT accessory, 261 VT Acqmeter display, 262 VT controller, 264 VT Controller label, 262 VT controller power, 262 VT controller safety circuits, 266 VT cutoff, 262 VT experiment warning, 24 VT FAILURE message, 267 VT gas flow rate, 265 VT indicator light, 263, 264 VT nitrogen gas flow, 263 VT operation, 261 VT operation with a sample changer, 299 VT regulation light, 267, 267 vtc parameter, 261, 263 vtwait parameter, 265, 267 VXR-4000/Gemini data system, 172

W

WALTZ decoupling, 171 WALTZ-4 modulation, 181 warnings defined, 23 water in VT exchanger, 264 water suppression, 73, 85 wave definition string, 231 wave string variables, 238 waveform generation, 228 wavelib directory, 231 wavelib waveform name parameter, 256, 256 wbs parameter, 308 wc parameter, 34, 107 wc2 parameter, 34, 107 weighting values generation, 93 werr parameter, 298, 308 WET VAST, 318

WET experiments solvent suppression, VAST, 322 WET shapes VAST, 323 Wet Width VAST, 320 wexp parameter, 302, 308 WFG (waveform generator), 371 wft command, 33, 95, 98 wft1d command, 101, 132 wft1da macro, 102 wft1dac and wft2dac commands, 164 wft2d command, 92, 98, 101, 102, 131, 132 wft2da macro, 102, 102, 104 wft2dac macro, 10 wftt3 macro, 129 whitewashed stacked plot, 113 whitewashing, 111 width of chart, 107 wnt parameter, 308 Workspace button, 30 Workspace Menu, 30 wp parameter, 109, 117, 215 wp1 parameter, 109, 117 wp2 parameter, 127 wrap parameter, 237, 239 Write Text button, 124 writing a numeric Pbox parameter, 255 writing expressions, 235 writing VAST protocols, 363 wshim parameter, 379 wti command, 94, 132 wtia command, 95

X

X-approximation, 199 X-axis label, 221 xdiag parameter, 120 X-nucleus decoupling, 180 X-nucleus decoupling frequency, 170 XYZ coordinates of the liquid handler arm, 361 XYZ coordinates of the rinse station, 354

Y

Y-axis label, 221

Ζ

z0 parameter, 305 Z1 room temperature shim as a gradient, 385 Z1 shimming shift, 372 zeroneg command, 112 zero-order phase rotation, 104 zfs (zero-frequency suppression), 104