

Table of Contents

UNITY INOVA

Acceptance Tests Specifications

UNITY INOVA NMR Spectrometer Systems

Pub. No. 87-195329-00, Rev. B0795

Chapter 1 Introduction	7
1.1 Introduction	7
Chapter 2 Consoles	9
2.1 Specifications	10
Amplitude and Stability Specifications for ¹ H Observers	11
Pulse-to-Pulse Phase Stability Specifications	12
Lock with ¹ H and ¹³ C and ¹⁵ N	13
Quadrature Impedance	14
Homonuclear Decoupling Specifications	15
WALTZ ¹ H Decoupling Specifications without Programable Decoupling	16
WALTZ ¹ H Decoupling Using Programable Decoupling Specifications	17
Shaped Pulses Specifications	18
Variable Temperature Operation Specifications	19
Magnet Drift Specifications	20
2.2 Uncontracted Custom Console Specifications	21
Temperature Accuracy for VT System	22
Stability Calibration for High-Stability VT	23
Homospin Demonstration	24
Chapter 3 Liquids Probes Specifications	25
3.1 90° Pulse Width Specifications and Performance Characteristics	26
About the Specifications	27
About Test Samples	28
3.2 Resolution Specifications	31
About the Specifications	32
About Test Samples	33
3.3 Spinning Sidechain Specifications	34
About the Specifications	35
About Test Samples	36
3.4 Sensitivity Specifications	39
About the Specifications	40
About Test Samples	41

varian 

nuclear magnetic resonance instruments

3.1 90° Pulse Width Specifications

The 90° pulse width acceptance test specifications for Varian liquids probes are listed in Table 6. Probes are categorized by frequency, and specifications are listed under the appropriate nucleus.

About the 90° Pulse Width Specifications

All 90° pulse width values are in microseconds and are less than or equal to the values listed in the table. Acceptance test specifications are achieved using the Varian UNITYINOVA Acceptance Test Procedures.

Blank spaces (shown as “—”) in the specifications column indicate that Varian does not specify a value or that it is not applicable to specify a value. Asterisks (*) in the specifications column indicate that the 90° pulse width was determined using the indirect method. For custom specifications that have been purchased, please refer to the Custom Specifications Form for the acceptance test specifications. Refer to the *NMR Probes Installation Manual* for the power-handling capability of your probe.

All specifications for Bioproton•nmr probes are achieved using the Ultra•nmr shims.

About Test Samples

Table 5 lists the samples used to achieve the 90° pulse width specifications. Note that the 4-mm sample tube is used only for Nano•nmr probes. The 4-mm sample tube has a 4-mm outside diameter (O.D.) and a 40- μ L sample volume. Note also that doped D₂O can be used for the ¹H pulse width test as an alternative to the sample specified in Table 5.

Table 5. Samples for 90° pulse width tests

Nucleus	Sample Tube (mm)	Test Sample	Sample Part Number
¹ H	5	0.1% ethylbenzene in CDCl ₃	00-968120-70
¹ H	10	0.1% ethylbenzene in CDCl ₃	00-968123-70
¹ H	4	0.1% ethylbenzene in CDCl ₃	00-993143-99
¹⁹ F	5	0.05% trifluorotoluene in benzene-d ₆	00-968120-82
³¹ P	5	0.0485 M triphenylphosphate in CDCl ₃	00-968120-87
³¹ P	10	0.0485 M triphenylphosphate in CDCl ₃	00-968123-87
¹³ C	5	40% p-dioxane in benzene-d ₆ (ASTM)	00-968120-69
¹³ C (indirect)	5	1% iodomethane- ¹³ C in CDCl ₃	00-968120-96
¹³ C	10	40% p-dioxane in benzene-d ₆ (ASTM)	00-968123-69
¹⁵ N	5	90% formamide in DMSO-d ₆	00-968120-83
¹⁵ N (indirect)	5	2% benzamide- ¹⁵ N in DMSO-d ₆	00-968120-97
¹⁵ N	10	90% formamide in DMSO-d ₆	00-968123-83
²⁹ Si	5	25% hexamethyldisiloxane in benzene-d ₆	00-968120-98

Table 6. 90° pulse width specifications

Probe	90° Pulse Width Specifications (μs)						
	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
200 MHz							
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	15	25	12	10	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	15	25	—	18	43	—	—
5-mm ¹ H/ ¹⁹ F proton	10	10	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	20	20	—	15	—	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	15	—	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	20	—	—	—
10-mm ¹⁰³ Rh- ⁹ Be broadband	Available by request.						
300 MHz							
5-mm ¹ H{ ¹⁵ N- ³¹ P} indirect detection	12	—	—	—	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	15	20	12	10	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	15	25	—	15	50	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ²⁹ Si Auto•nmr	15	20	—	10	—	18	18
4-mm (40 μL) ¹ H Nano•nmr	Available by request.						
5-mm ¹ H/ ¹⁹ F proton	10	10	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	20	20	—	15	—	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	15	—	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	20	—	—	—
10-mm ¹⁰³ Rh- ¹⁵ N broadband	Available by request.						
400 MHz							
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	10	—	—	12*	25*	—	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	10	—	—	12*	25*	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	10	—	—	15*	45*	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	10	—	—	15*	45*	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	20	30	20	15	—	—	—
4-mm (40 μL) ¹ H Nano•nmr	Available by request.						
5-mm ¹ H/ ¹⁹ F proton	10	10	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	25	25	—	15	—	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	15	—	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	20	—	—	—
10-mm ¹⁰³ Rh- ¹⁵ N broadband	Available by request.						
500 MHz							
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	10	—	—	12*	25*	—	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	10	—	—	12*	25*	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	10	—	—	15*	45*	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	10	—	—	15*	45*	—	—

3.1 90° Pulse Width Specifications **Table 6.** 90° pulse width specifications (Continued)

Probe	90° Pulse Width Specifications (μs)						
5-mm $^1\text{H}\{^{13}\text{C}/^{31}\text{P}\}$ Triple•nmr	10	—	40*	15*	—	—	—
5-mm $^{13}\text{C}\{^1\text{H}/^{15}\text{N}\}$ Triple•nmr	30	—	—	9	25	—	—
5-mm Proton•nmr	10	10	—	—	—	—	—
4-mm (40 μL) ^1H Nano•nmr	4.5	—	—	—	—	—	—
10-mm Bioproton•nmr	18	—	—	—	—	—	—
5-mm $^1\text{H}/^{19}\text{F}/^{15}\text{N}-^{31}\text{P}$ Switchable•nmr	12	12	18	12	—	—	—
5-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	—	—	15	—	—	—
5-mm $^{15}\text{N}-^{13}\text{C}$ broadband	—	—	—	15	—	—	—
10-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	—	—	25	—	—	—
10-mm $^{15}\text{N}-^{13}\text{C}$ broadband	—	—	—	20	—	—	—
10-mm $^{103}\text{Rh}-^{15}\text{N}$ broadband	—	—	—	—	30	—	—
600 MHz	^1H	^{19}F	^{31}P	^{13}C	^{15}N	^{29}Si	^{39}K
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ Indirect•nmr	10	—	—	14*	30*	—	—
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ PFG Indirect•nmr	10	—	—	14*	30*	—	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ Triple•nmr	10	—	—	15*	45*	—	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ PFG Triple•nmr	10	—	—	15*	45*	—	—
5-mm $^{13}\text{C}\{^1\text{H}/^{15}\text{N}\}$ Triple•nmr	30	—	—	9	25	—	—
5-mm Proton•nmr	10	10	—	—	—	—	—
10-mm Bioproton•nmr	20	—	—	—	—	—	—
5-mm $^{15}\text{N}-^{31}\text{P}$ Broadband•nmr	—	—	15	9	15	—	—
750 MHz	^1H	^{19}F	^{31}P	^{13}C	^{15}N	^{29}Si	^{39}K
Available by request.	Available by request.						

* 90° pulse width determined using the indirect method.

Table 9. Resolution and lineshape specifications (14-shim and 18-shim)

Probe	CHCl ₃ Sample Part No.	Resolution and Lineshape (Hz) for Standard Shim Systems (50.0%/0.55%/0.11%)		
		¹ H Spin (CHCl ₃)	¹ H Non-spin (CHCl ₃)	¹³ C Spin (ASTM)
200-MHz				
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	00-968120-76	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	00-968120-76	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹ H/ ¹⁹ F proton	00-968120-76	0.4/6.0/12.0	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	00-968120-76	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
300-MHz				
5-mm ¹ H{ ¹⁵ N- ³¹ P} indirect detection	00-968120-99	0.4/6.0/12.0	—	—
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	00-968120-99	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	00-968120-99	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ²⁹ Si Auto•nmr	00-968120-99	0.4/6.0/12.0	—	0.2/3.0/5.0
4-mm (40 µL) ¹ H Nano•nmr		Available by request.		
5-mm ¹ H/ ¹⁹ F proton	00-968120-89	0.4/6.0/12.0	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	00-968120-99	0.4/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
400-MHz				
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	00-968120-89	0.45/6.0/12.0	—	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	00-968120-89	0.45/6.0/12.0	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	00-968120-89	0.45/6.0/12.0	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	00-968120-89	0.45/6.0/12.0	—	—
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	00-968120-89	0.45/6.0/12.0	—	0.2/3.0/5.0
4-mm (40 µL) ¹ H Nano•nmr		Available by request.		
5-mm ¹ H/ ¹⁹ F proton	00-968120-89	0.45/6.0/12.0	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	00-968120-99	0.45/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0

Table 10. Resolution and lineshape specifications (23-shim)

Probe	CHCl ₃ Sample Part No.	Resolution and Lineshape (Hz) for Enhanced Shim Systems (50.0%/0.55%/0.11%)		
		¹ H Spin (CHCl ₃)	¹ H Non-spin (CHCl ₃)	¹³ C Spin (ASTM)
400-MHz				
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/8/16	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/8/16	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/8/16	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/8/16	—
5-mm 4-nuc- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	00-968120-89	0.45/6.0/12.0	—	0.2/3.0/5.0
4-mm (40 μL) ¹ H Nano•nmr		Available by request.		
5-mm ¹ H/ ¹⁹ F proton	00-968120-89	0.45/6.0/12.0	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	00-968120-99	0.45/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ³¹ P broadband	—	—	—	0.2/3.0/5.0
500-MHz				
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹³ C/ ³¹ P} Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹³ C{ ¹ H/ ¹⁵ N} Triple•nmr	—	—	—	0.2/3.0/5.0
5-mm Proton•nmr	00-968120-89	0.45/6.0/12.0	—	—
4-mm (40 μL) ¹ H Nano•nmr	00-993143-99	0.6/8.0/13.0	—	—
10-mm Bioproton•nmr	00-968123-89	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P Switchable•nmr	00-968120-89	0.45/6.0/12.0	—	0.2/3.0/5.0
5-mm ¹⁵ N- ³¹ P broadband	00-968120-89	—	—	0.2/3.0/5.0
5-mm ¹⁵ N- ¹³ C broadband	00-968120-89	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ³¹ P broadband	00-968120-89	—	—	0.2/3.0/5.0
10-mm ¹⁵ N- ¹³ C broadband	00-968120-89	—	—	0.2/3.0/5.0
600-MHz				
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	00-968120-89	0.45/6.0/12.0	0.65/10/20	—
5-mm ¹³ C{ ¹ H/ ¹⁵ N} Triple•nmr	—	—	—	0.2/3.0/5.0
5-mm Proton•nmr	00-968120-89	0.45/6.0/12.0	—	—
10-mm Bioproton•nmr	00-968123-89	—	—	—
5-mm ¹⁵ N- ³¹ P Broadband•nmr	—	—	—	0.2/3.0/5.0
750-MHz				
Available by request.		Available by request.		

Table 14. Spinning sidebands specifications

Probe	Spinning Sidebands	
200 MHz	^1H	^{13}C
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{31}\text{P}$ Auto•nmr	1%	1%
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{15}\text{N}$ Auto•nmr	1%	1%
5-mm $^1\text{H}/^{19}\text{F}$ proton	1%	—
5-mm $^1\text{H}/^{19}\text{F}/^{15}\text{N}-^{31}\text{P}$ switchable	1%	1%
5-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
10-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
300 MHz	^1H	^{13}C
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ indirect detection	1%	1%
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{31}\text{P}$ Auto•nmr	1%	1%
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{15}\text{N}$ Auto•nmr	1%	1%
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{29}\text{Si}$ Auto•nmr	1%	1%
4-mm (40 μL) ^1H Nano•nmr	1%	—
5-mm $^1\text{H}/^{19}\text{F}$ proton	1%	—
5-mm $^1\text{H}/^{19}\text{F}/^{15}\text{N}-^{31}\text{P}$ switchable	1%	1%
5-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
10-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
400 MHz	^1H	^{13}C
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ Indirect•nmr	1%	—
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ PFG Indirect•nmr	1%	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ Triple•nmr	1%	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ PFG Triple•nmr	1%	—
5-mm 4-nucleus- $^1\text{H}/^{19}\text{F}/^{13}\text{C}/^{31}\text{P}$ Auto•nmr	1%	1%
4-mm (40 μL) ^1H Nano•nmr	1%	—
5-mm $^1\text{H}/^{19}\text{F}$ proton	1%	—
5-mm $^1\text{H}/^{19}\text{F}/^{15}\text{N}-^{31}\text{P}$ switchable	1%	1%
5-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
10-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	1%
500 MHz	^1H	^{13}C
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ Indirect•nmr	1%	—
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ PFG Indirect•nmr	1%	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ Triple•nmr	1%	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ PFG Triple•nmr	1%	—
5-mm $^1\text{H}\{^{13}\text{C}/^{31}\text{P}\}$ Triple•nmr	1%	—
5-mm $^{13}\text{C}\{^1\text{H}/^{15}\text{N}\}$ Triple•nmr	—	1%
5-mm Proton•nmr	1%	—

Table 16. Sensitivity (S/N) specifications

Probe	Sensitivity (S/N) Specifications						
	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
200 MHz	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	50:1	50:1	30:1	40:1	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	50:1	50:1	—	30:1	6:1	—	—
5-mm ¹ H/ ¹⁹ F proton	85:1	85:1	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	35:1	35:1	35:1	40:1	6:1	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	45:1	50:1	7:1	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	135:1	150:1	20:1	—	—
10-mm ¹⁰³ Rh- ⁹ Be broadband	Available by request.						
300 MHz	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
5-mm ¹ H{ ¹⁵ N- ³¹ P} indirect detection	135:1	—	—	—	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	100:1	90:1	50:1	80:1	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ¹⁵ N Auto•nmr	90:1	90:1	—	70:1	8:1	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ²⁹ Si Auto•nmr	100:1	90:1	—	80:1	—	80:1	80:1
4-mm (40 µL) ¹ H Nano•nmr	Available by request.						
5-mm ¹ H/ ¹⁹ F proton	175:1	175:1	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	80:1	80:1	70:1	80:1	12:1	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	90:1	100:1	15:1	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	280:1	300:1	40:1	—	—
10-mm ¹⁰³ Rh- ¹⁵ N broadband	Available by request.						
400 MHz	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	350:1	—	—	—	—	—	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	350:1	—	—	—	—	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	350:1	—	—	—	—	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	350:1	—	—	—	—	—	—
5-mm 4-nucleus- ¹ H/ ¹⁹ F/ ¹³ C/ ³¹ P Auto•nmr	135:1	135:1	60:1	100:1	—	—	—
4-mm (40 µL) ¹ H Nano•nmr	Available by request.						
5-mm ¹ H/ ¹⁹ F proton	250:1	250:1	—	—	—	—	—
5-mm ¹ H/ ¹⁹ F/ ¹⁵ N- ³¹ P switchable	100:1	100:1	100:1	120:1	15:1	—	—
5-mm ¹⁵ N- ³¹ P broadband	—	—	110:1	150:1	20:1	—	—
10-mm ¹⁵ N- ³¹ P broadband	—	—	350:1	450:1	60:1	—	—
10-mm ¹⁰³ Rh- ¹⁵ N broadband	Available by request.						
500 MHz	¹ H	¹⁹ F	³¹ P	¹³ C	¹⁵ N	²⁹ Si	³⁹ K
5-mm ¹ H{ ¹⁵ N- ³¹ P} Indirect•nmr	450:1	T	—	—	—	—	—
5-mm ¹ H{ ¹⁵ N- ³¹ P} PFG Indirect•nmr	450:1	T	—	—	—	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} Triple•nmr	450:1	T	—	—	—	—	—
5-mm ¹ H{ ¹³ C/ ¹⁵ N} PFG Triple•nmr	450:1	T	—	—	—	—	—

Table 16. Sensitivity (S/N) specifications (Continued)

Probe	Sensitivity (S/N) Specifications						
5-mm $^1\text{H}\{^{13}\text{C}/^{31}\text{P}\}$ Triple•nmr	450:1	T	—	—	—	—	—
5-mm $^{13}\text{C}\{^1\text{H}/^{15}\text{N}\}$ Triple•nmr	—	—	—	250:1	—	—	—
5-mm Proton•nmr	450:1	450:1	—	—	—	—	—
4-mm (40 μL) ^1H Nano•nmr	90:1	T	—	—	—	—	—
10-mm Bioproton•nmr	850:1	T	—	—	—	—	—
5-mm $^1\text{H}/^{19}\text{F}/^{15}\text{N}-^{31}\text{P}$ Switchable•nmr	275:1	275:1	90:1	180:1	—	—	—
5-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	—	180:1	190:1	25:1	—	—
5-mm $^{15}\text{N}-^{13}\text{C}$ broadband	—	—	—	250:1	35:1	—	—
10-mm $^{15}\text{N}-^{31}\text{P}$ broadband	—	—	380:1	380:1	45:1	—	—
10-mm $^{15}\text{N}-^{13}\text{C}$ broadband	—	—	—	650:1	85:1	—	—
10-mm $^{103}\text{Rh}-^{15}\text{N}$ broadband	—	—	—	—	80:1	—	70:1
600 MHz							
	^1H	^{19}F	^{31}P	^{13}C	^{15}N	^{29}Si	^{39}K
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ Indirect•nmr	575:1	T	—	—	—	—	—
5-mm $^1\text{H}\{^{15}\text{N}-^{31}\text{P}\}$ PFG Indirect•nmr	625:1	T	—	—	—	—	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ Triple•nmr	575:1	T	—	—	—	—	—
5-mm $^1\text{H}\{^{13}\text{C}/^{15}\text{N}\}$ PFG Triple•nmr	625:1	T	—	—	—	—	—
5-mm $^{13}\text{C}\{^1\text{H}/^{15}\text{N}\}$ Triple•nmr	—	—	—	340:1	—	—	—
5-mm Proton•nmr	600:1	600:1	—	—	—	—	—
10-mm Bioproton•nmr	1050:1	—	—	—	—	—	—
5-mm $^{15}\text{N}-^{31}\text{P}$ Broadband•nmr	—	—	160:1	365:1	35:1	—	—
750 MHz							
	^1H	^{19}F	^{31}P	^{13}C	^{15}N	^{29}Si	^{39}K
Available by request.	Available by request.						

